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- [1] S.K. Singh, S. Srivastav, A.K. Yadav, S. Srikrishna, G. Perry, Overview of Alzheimer's disease and some therapeutic approaches targeting abeta by using several synthetic and herbal compounds, *Oxid. Med. Cell. Longevity* 2016 (2016) 7361613.
- [2] J. Mattila, H. Soininen, J. Koikkalainen, D. Rueckert, R. Wolz, G. Waldemar, J. Lotjonen, Optimizing the diagnosis of early Alzheimer's disease in mild cognitive impairment subjects, *J. Alzheimer's Disease* 32 (2012) 969–979.
- [3] C.P. Ferri, M. Prince, C. Brayne, H. Brodaty, L. Fratiglioni, M. Ganguli, K. Hall, K. Hasegawa, H. Hendrie, Y. Huang, A. Jorm, C. Mathers, P.R. Menezes, E. Rimmer, M. Sczufca, I. Alzheimer's Disease, Global prevalence of dementia: a Delphi consensus study, *Lancet* 366 (2005) 2112–2117.
- [4] K. Blennow, H. Hampel, CSF markers for incipient Alzheimer's disease, *Lancet Neurol.* 2 (2003) 605–613.
- [5] B. Van Broeck, M. Timmers, S. Ramael, J. Bogert, L.M. Shaw, M. Mercken, J. Slemmon, L. Van Nueten, S. Engelborghs, J.R. Streffer, Impact of frequent cerebrospinal fluid sampling on Abeta levels: systematic approach to elucidate influencing factors, *Alzheimer's Res. Therapy* 8 (2016) 21.
- [6] Lin P, Chen W, Yuan F, et al. An UHPLC–MS/MS method for simultaneous quantification of human amyloid beta peptides A β 1-38, A β 1-40 and A β 1-42 in cerebrospinal fluid using micro-elution solid phase extraction[J]. *Journal of Chromatography B*, 2017, 1070: 82-91.

For Research Use Only. Not for use in Diagnostics Procedures.

AB Sciex is operating as SCIEX.

© 2019. AB Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

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