



# Quantitation Method for Nitrofurans Metabolites in Milk using SCIEX Triple Quad™ 3500 System

M.Chandrasekar<sup>1</sup>; Aman Sharma<sup>1</sup>; Santosh Kapil<sup>1</sup>; Anoop Kumar<sup>1</sup>; Manoj Pillai<sup>1</sup> and Jinaru Stahl-Zeng<sup>2</sup>  
<sup>1</sup>SCIEX, 121, Udyog Vihar, Phase IV, Gurgaon - 122015, Haryana, India, <sup>2</sup>Darmstadt, Germany

## Overview

A liquid chromatography tandem mass spectrometry (LC-MS/MS) method for quantification of Nitrofurans metabolites in milk was developed. The method presented adequate linearity with correlation coefficients above  $r \geq 0.99$  for both analytes in the dynamic range of 0.50–20.0  $\mu\text{g}/\text{kg}$ , with average accuracies for matrix based recovery were in the range 85%–120%. The results qualified the method for the quantification and confirmation of the analytes in milk at concentrations lower to the established Minimum Required Performance Limit (1.0 $\mu\text{g}/\text{kg}$ ).

## Introduction


Nitrofurans are synthetic chemotherapeutic agents which have a broad spectrum of bacteriostatic activity. Nitrofurans mainly inhibit the enzymes involved in the carbohydrate metabolism. These bacteriostatic compounds are prohibited in livestock production by EU which is listed in Annex IV of EC Council Regulation 2377/90. No MRLs have been established for Nitrofurans hence it is necessary to have sensitive confirmatory analytical methods for the detection of nitrofurans residues in food commodities. Further the detection of Nitrofurans has been shown to be difficult as they are quickly metabolized. Therefore the analysis of the protein bound, solvent extractable metabolites of Nitrofurans have been reported as the ideal choice of analysis. Analytically, residues are checked only for marker metabolites of the 4 nitrofurans chemicals, in particular: 3-amino-2-oxazolidinone (AOZ) for furazolidone, 3-amino-5-methylmorpholino-2-oxazolidinone (AMOZ) for furaltadone, 1-aminohydantoin (AHD) for nitrofurantoin and Semicarbazide (SEM) for Nitrofurazone (Figure 2).



Figure 1. SCIEX Triple Quad™ 3500

In general the study of nitrofurans metabolites in food samples requires incubation period for derivatization with nitrobenzaldehyde for 16hr at 37°C in dark. The quantitative and confirmatory determination of nitrofurans metabolites was performed by liquid chromatography/electrospray ionization tandem mass spectrometry (LC/ESI-MS/MS) in positive ion mode, according to European Decision 2002/657/EC. The MRPL for nitrofurans metabolites (individual) is 1.0 $\mu\text{g}/\text{kg}$  as per RMP/EU/2016-17.

The present application note describes a method which is sensitive and selective enough to meet the global guidelines analyze the nitrofurans metabolites in milk using SCIEX Triple Quad™ 3500 LC-MS/MS System.



**Figure 2. Structures of Nitrofuran, Nitrofuran metabolites and Nitrophenyl derivatives.**

## Materials and Methods

### Chemicals

Nitrofuran metabolites Standards were purchased from Clearsynth and 2-Nitrobenzaldehyde was purchased from Sigma Aldrich  $\geq 99\%$  Purity. All other chemicals used were of LC-MS grade.

### Sample Preparation

Milk sample (3 ml) was mixed with 1ml of HCl (0.1M) and 50mM of 2-Nitrobenzaldehyde (0.3ml), vortexed and incubated on ultrasonic bath for 16hr added 0.6ml of 1M K<sub>2</sub>HPO<sub>4</sub> solution and 10 ml of ethyl acetate, vortexed, followed by centrifugation at 4000 rpm, the supernatant was evaporated to dryness reconstituted with 1ml of Methanol: water (5:95) and 10 $\mu$ l is used for LC-MS/MS analysis.

## Experimental Conditions

### LC Conditions

LC separation was performed on a Shimadzu instrument using Zorbax Eclipse Plus C18(150 X 4.6)mm 5.0 $\mu$  and a fast gradient of 1mM Ammonium acetate(Mobile Phase A) and Methanol(Mobile Phase B) at a flow rate of 0.4ml/min (Table 1).

Time (min)	Mobile phase A%	Mobile phase B%
0.01	95	5
0.50	45	55
3.50	45	55
4.00	95	5
12.00	Controller	Stop

**Table 1: Mobile Phase Gradient**

### MS/MS Conditions

The SCIEX Triple Quad™ 3500 was operated in Multiple Reaction Monitoring (MRM) mode. The TurboV™ source was used with an Electrospray Ionization (ESI) probe in positive ionization mode at 5500 ion spray voltage. Two selective MRM transitions were monitored and ion ratio was calculated automatically by software for compound identification (Table 2). Analyst® 1.6.2 Software was used for method development and data acquisition. LC-MS/MS data was processed using the MultiQuant™ Software version 3.0.2

Compound	Precursor ion	Product ion Quantifier	Product ion Qualifier
AOZ	236.0	104.0	78.0
AMOZ	335.0	291.1	128.2
SEM	209.0	166.0	192.0
AHD	249.1	134.0	104.0

**Table 2: MRM transitions**

## Results and Discussions

### Sensitivity, Reproducibility, Linearity and Accuracy

The developed method showed signal-to-noise ratio > 23 for all the analytes with sample extracted at a level of 1.0  $\mu$ g/kg (Spiked) which meets the regulatory criterion (Figure 3)




Figure 3: Chromatogram (1.0 µg/kg) showing signal to noise

Matrix based Calibration curve was plotted, found linear in the range of 0.50 µg/kg (ppb) to 20.0µg/kg (ppb) and correlation regression co-efficient  $r > 0.98$  for both quantifier and qualifier ions by applying weighing factor of 1/X<sup>2</sup> (Table 4).

Repeatability at three levels (1/2 MRPL, MRPL, 1.5MRPL) were evaluated for 6 injections and %relative standard deviation (%CV) was observed to be less than 10 (Table 3). Accuracies observed were in the range from 85% to 120%.

Analyte	Repeatability			Recovery (n=6)		
	½ MRPL (0.5ppb)	MRPL (1.0ppb)	1.5MRPL (1.5ppb)	½ MRPL (0.5ppb)	MRPL (1.0ppb)	1.5MRPL (1.5ppb)
AOZ	3.41	5.04	1.72	98.83	105.00	107.89
AMOZ	5.98	2.89	3.49	94.37	92.82	101.67
SEM	8.70	8.84	3.30	92.37	102.75	108.33
AHD	7.12	5.86	8.00	92.57	86.93	96.56

Table 3: Repeatability (%CV) and recovery statistics and in Milk

Analyte	Calibration Range (ppb)	Linearity (r)	CC $\alpha$	CC $\beta$
AOZ	0.5 -20	0.9994	0.54	0.57
AMOZ	0.5 -20	0.9994	0.57	0.61
SEM	0.5 -20	0.9964	0.59	0.65
AHD	0.5 -20	0.9992	0.58	0.63

Table 4: Summary of CC $\alpha$ , CC $\beta$  and linearity in milk Sample

Decision limit (CC $\alpha$ ) and detection capability (CC $\beta$ ) were calculated for AOZ, AMOZ, SEM and AHD derivatives of Nitrofurans in milk samples. The calculation was based on using linear regression model analyzing spiked milk samples at below MRPL level (Van Loco et al, 2007).

The calculated value of CC $\alpha$  and CC $\beta$  are given in Table 4. The decision limit (CC $\alpha$ ) and detection capability (CC $\beta$ ) of all the compounds were well below the MRPL.




Figure 4: Matrix based calibration curve AOZ, AMOZ, SEM and AHD

## Conclusion

The method and data acquired here gives sensitive and accurate solution for the quantitation and confirmation of Nitrofuran metabolites in Milk samples by LC-MS/MS. The SCIEX™ 3500 system provides good sensitivity and selectivity for this analysis, allowing maximum output for the analysis of a bigger batch of samples in a short time period. Automatic ion ratio calculation in MultiQuant™ software can be used for confirmation of compound. The method showed acceptable accuracies (85%-120%), linearity with  $r > 0.99$  for both quantifier and qualifier, repeatability (%CV) observed was less than 10. The method allows high throughput, selective, rapid and sensitive LC-MS/MS identification and quantitation of banned Nitrofuran metabolites meeting EU MRPL of 1.0 µg/kg.

7. J. Van Loco, A. Janosi, S. Impens, S. Fraselle, V. Cornet, J.M. Degroodt, Calculation of the decision limit (CC $\alpha$ ) and the detection capability (CC $\beta$ ) for banned substances: The imperfect marriage between the quantitative and the qualitative criteria. *Analytica Chimica Acta* Volume 586, Issues 1–2, 14, Pages 8–12 (2007)

## References

1. FAO-Dairy Production and products: milk production facts, Food and Agriculture Organization of United Nations, Rome.-2013 <http://fao.org/agriculture/dairy-gateway/meat-production/en>.
2. Method validation and quality control procedures for pesticide residue analysis in food and feed: SANCO/12495/2011
3. FSSAI manual of methods of analysis of food. [http://fssai.gov.in/Portals/0/Pdf/15Manuals/ANTIBIOTIC S%20AND%20RESIDUES.pdf](http://fssai.gov.in/Portals/0/Pdf/15Manuals/ANTIBIOTIC%20AND%20RESIDUES.pdf)
4. Schreiber: 'MultiQuant™ Software Version 3.0 - Improving Data Quality and Processing Throughput with Better Peak Integration, Quantitative and Qualitative Compound Review for the Analysis of Food, Drinking Water, and Environmental Samples' Application Note AB SCIEX (2013) #8160213-01
5. Neha Bhasin, Prasanth Joseph, Praveen Sharma, Manoj Pillai, Jens Dahlmann, André Schreiber and Christopher Borton; Targeted Multi-Residue LC-MS/MS Method for Sulfonamide and Nitroimidazole Antibiotics in Honey, Poster, AOAC conference (2015)
6. Tomasz Śniegocki, Andrzej Posyniak, And Jan Żmudzki, Determination of Nitrofuran Metabolite Residues In Eggs By Liquid Chromatography-Mass Spectrometry. *Bull Vet Inst Pulawy* 52, 421-425 (2008)

AB Sciex is doing business as SCIEX.

© 2017 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-02-5032-A



**Headquarters**

500 Old Connecticut Path | Framingham, MA 01701 USA  
Phone 508-383-7700  
[sciex.com](http://sciex.com)

**International Sales**

For our office locations please call the division  
headquarters or refer to our website at  
[sciex.com/offices](http://sciex.com/offices)