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1. Acoustic dispensing-mass spectrometry: the next high throughput bioanalytical platform for early drug discovery
DOI: 10.4155/bio-2017-4980
2. Operational Modes and Speed Considerations of an Acoustic Droplet Dispenser for Mass Spectrometry
DOI: 10.1021/acs.analchem.0c02999
3. Ultrahigh-Throughput ESI-MS Sampling Pushed to Six Samples per Second by Acoustic Ejection Mass Spectrometry
DOI: 10.1021/acs.analchem.0c01632
4. Ultrahigh-Throughput and Chromatography-Free Bioanalysis of Polar Analytes with Acoustic Ejection Mass Spectrometry
DOI: 10.1021/acs.analchem.0c03006
5. Ultra-High-Throughput Acoustic Droplet Ejection-Open Port Interface-Mass Spectrometry for Parallel Medicinal Chemistry
DOI: 10.1021/acsmmedchemlett.0c00066
6. Acoustic Ejection Mass Spectrometry for High-Throughput Analysis
DOI: 10.1101/2020.01.28.923938
7. Acoustic Droplet Ejection and Open Port Interface for Rapid Analysis of Metabolic Stability Assays
DOI: 10.1016/j.xphs.2020.07.025
8. Emerging Chromatography-Free High-Throughput Mass Spectrometry Technologies for Generating Hits and Leads
DOI: 10.1021/acsmmedchemlett.0c00314
9. Acoustic Ejection/Full-Scan Mass Spectrometry Analysis for High Throughput Compound Quality Control
DOI: 10.1177/2472630320967625
10. Direct Analysis from Phase-Separated Liquid Samples using ADE-OPI-MS: Applicability to High-Throughput Screening for Inhibitors of Diacylglycerol Acyltransferase 2
DOI: 10.1021/acs.analchem.0c04312
11. Fluid Dynamics of the Open Port Interface for High-Speed Nanoliter Volume Sampling Mass Spectrometry
DOI: 10.1021/acs.analchem.1c01312
12. Acoustic Ejection Mass Spectrometry: A Fully Automatable Technology for High-Throughput Screening in Drug Discovery
DOI: 10.1177/24725552211028135
13. Ultra high-throughput mass spec for life sciences research
DOI: <https://www.nature.com/articles/d42473-020-00307-5>
14. High-Throughput Analysis from Complex Matrices: Acoustic Ejection Mass Spectrometry from Phase-Separated Fluid Samples
DOI: 10.3390/metabo11110789

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Accelerated pharmacokinetic profiling and metabolite monitoring Using Echo[®] MS System

Acoustic Ejection Mass Spectrometry

Neil Devenport¹, Neil Walsh¹, Simon Wood²

¹SCIEX UK, ²Cyprotex Discovery UK

The generation of clearance profiles and calculation of intrinsic clearance (CL_{int}) values is an important stage in pre-clinical drug discovery.¹ The pharmacokinetic (PK) assessment of potential drug candidates is typically performed by a time-course incubation with human liver microsomes (HLM) before LC-MS analysis. The *in vitro* CL_{int} can be scaled to *in vivo* and used to predict human clearance which can be then used to inform dosing studies, as well as provide insight into enzymes responsible for metabolism of drug candidates.

As the process is performed during the early stage of drug discovery, the number of potential candidates can still be in the thousands, or even tens of thousands, and CL_{int} values will need to be calculated for each target. This will require multiple time-points for each target leading to a large analytical demand. This creates a significant bottleneck using current LC-MS solutions.

The Echo MS System using Acoustic Ejection Mass Spectrometry (AEMS) provides a sample introduction method capable of analyzing 1 sample per second and up to 3 samples per second when multiplexed.² With the Echo MS System, very small, reproducible droplets can be ejected from high-density well plates, and this ejection process is compatible with the wide range of solutions and complex matrices used in drug discovery. The ejected droplets are then captured by the Open Port Interface (OPI) device which uses a laminar flow of carrier

solution that dilutes and transports the sample to the electrospray source, where it is ionized by conventional electrospray (Figure 2).² This process forms a contactless sample introduction method with minimal carry-over and minimal matrix effects.

Here, the utility of this high-throughput solution was explored for use in candidate screening. The goal was to highlight the ability of the Echo MS System to alleviate the analytical bottleneck of traditional LC-MS solutions through speed and quantitative performance, and show how that correlates with current Cyprotex methodology.

Key features of the Echo[®] MS System

- Rapid sample analysis using Acoustic Droplet Ejection and Open Port Interface:
 - One sample per second for a single analyte, up to three samples per second in multiplex mode
 - Highly reproducible sample injection
 - Extremely low potential for carryover
- Broad compound coverage using electrospray ionization featuring the OptiFlow[®] Turbo V Ion Source for robust and efficient ionization
- Industry proven high sensitivity and quantitative robustness using the SCIEX Triple Quad[™] 6500+ LC-MS/MS System
- Ease of operation using SCIEX OS Software
- Purpose built for incorporation into HTS workflows
 - Sample plate tray can be accessed by most robotic arms
- Open software API for incorporation of platform into existing automated HTS environments


Figure 1. Microsomal incubation of buspirone. Example data using the Echo MS System were obtained from a microsomal incubation of buspirone with data acquired from 5 timepoints, a control and a blank reference, all completed in under 90 seconds. Sampling was done in a serpentine manner from the 384-well plate.

Methods

Sample preparation: Microsomal incubations were performed in human liver microsomes (0.1M phosphate buffer, 0.5 mg/mL) initiated with NADPH (1mM) for 39 traditional small molecule APIs (1μM) and quenched in acetonitrile at 0, 5, 15, 30 and 45 time-points. These same samples were split for LC-MS analysis and Echo[®] MS System analysis.

Acoustic ejection: The Echo MS System was operated using 70% acetonitrile with 2mM ammonium formate and 0.1% formic acid at a flow rate of 450 μL/min. A total ejection volume of 25 nL was employed for each acoustic event at a frequency of 1.5 Hz.

Mass spectrometry: The SCIEX Triple Quad 6500+ System was operated in MRM mode monitoring a single transition with a total analytical cycle time of 100 msec. The ESI source was operated at 500 °C with a capillary voltage of 5500 V and the GS1/GS2 gas flows at 90 and 50 psi, respectively.

Data processing: SCIEX OS Software 1.6 was employed for the optimization of the target MRM transitions, automatic integration of the MS peaks and the automatic creation of the results file and associated .txt file for LIMS integration.

Results

The Echo MS System is specifically calibrated to eject droplets of 2.5 nL per acoustic event. The effective ejection volume can be increased by firing multiple 2.5 nL droplets at a high frequency (20 Hz), which the MS system detects as a single peak. The process allows the user to optimize analytical sensitivity versus the matrix effects, because of the “larger”


Figure 2. System schematic. Schematic of the Echo MS System showing the acoustic ejection of sample droplets and their capture and transfer to the OptiFlow Source by the OPI.


Figure 3. An example AEMS droplet ladder for lidocaine ranging from 1-10 droplets per MS measurement, equivalent to ejection volumes of 2.5-25 nL. The calibration line shown in insert with linear R^2 value >0.998 .

ejection volumes. Figure 3 shows this process using lidocaine as the target analyte.

The “droplet ladder” shown in Figure 3 shows the typical response obtained from the Echo MS System producing peak widths of approximately 400 msec at 50% peak height. These sharp peaks enable the 1 sample per second acquisition rates.

To assess the quantitative performance of the Echo MS System, solvent standard calibration lines for each of the 39 analytes were created across a target concentration range of 0.01-100nM. The AEMS analysis of the solvent standard calibration line was performed using an ejection volume of 25 nL. This will allow maximum analyte sensitivity from the non-complex solvent


Figure 4. Concentration curve for Midazolam. (Top) Midazolam ejection peaks across the concentration range of 100-0.05nM ran from high to low with $n=6$ technical replicate ejections from each well position. (Bottom) Linear regression for midazolam ejections show good quantitative performance ($R^2 >0.99$).

Table 1. Quantitative performance of the Echo MS System for six representative marketed drugs (0.01-100nM).

Compound name	MRM transition (Q1-Q3)	Calibration range (nM)	Regression value (R2)	Accuracy >LLOQ/@LLOQ	%CV >LLOQ/@LLOQ
Acetaminophen	152-110	1-100	0.994	108% / 97%	6% / 7%
Codeine	300-165	0.5-100	0.996	105% / 94%	3% / 18%
Diphenhydramine	256-167	0.1-100	0.991	95% / 98%	2% / 17%
Midazolam	326-291	0.05-100	0.996	93% / 100%	3% / 11%
Zolpidem	308-235	0.01-100	0.998	98% / 105%	2% / 20%

standard solution. Figure 4 shows example data for midazolam across the calibration range.

Following the quantitative testing, the clearance profiles for each of the target analytes were determined. In all cases, a single MRM transition was monitored with no internal standard used. The timepoints provided were 0, 5, 15, 30, 45 min. The peak

area values obtained from the analysis of T0 were set to 100% parent response, and the following time-points were provided as % parent remaining. Figure 5 provides three clearance regression profiles showing a range of clearance profiles.

Using the gradient of the clearance profiles, the elimination rate constant (k) and compound half-life ($t_{1/2}$) can be determined. The CL_{int} values for each of the target APIs can be calculated using the following equation

$$t_{1/2} = \frac{0.693}{k}$$

$$CL_{int} (\mu L / min / mg) = \frac{0.693}{t_{1/2}} \times \frac{\mu L \text{ of incubation}}{mg \text{ microsomes}}$$

Where; k = calculated elimination rate constant

A strong correlation ($r^2= 0.985$) of calculated CL_{int} is observed between the Echo MS System data and LC-MS data (Figure 6).


Figure 5. Clearance profiles. Clearance profiles for prazosin (top), erythromycin (center) and amodiaquine (bottom) are shown. The rapid acquisition approach using Echo MS System is shown as the blue trace vs. the same samples run with standard LC-MS analysis shown in the green trace.


Figure 6: Correlation analysis of Echo MS System data with LC-MRM data. Blue dotted line represents correlation, solid orange is the line of unity, dashed orange line is +/- 2-fold of unity. CL_{int} values below the assay level of quantification for the assay (<10.3 $\mu L/min/mg$ equivalent to $t_{1/2}$ of 3*incubation time) have been excluded.


Figure 7. Incubation time course. Echo MS System data monitoring hydroxy-imipramine over the incubation time course with n=6 replicate measurements per time point.

The data spans a range of clearance values and is scattered around the line of unity.

In addition to the clearance information, the ability to monitor suspected or known metabolites in a rapid manner might also be of benefit during late-discovery phase. The ability of the Echo MS System to provide this data was tested by repeating the measurement of the imipramine incubation with the additional MRM transition of the known hydroxy metabolite (297.1→86.1). The MS data produced from this scan is provided in Figure 7.

This metabolism data can then be related back to the clearance of the parent imipramine molecule (Figure 8) providing the analyst with DMPK information in rapid time frames suitable for the high through-put environment.

Conclusions

The data presented here confirms the viability of the use of the Echo MS System for early discovery phase studies providing DMPK data in seconds as opposed to the hours required when using traditional LC-MRM based scanning methodologies. The data generated from the Echo MS System was quantitative, with excellent %CV and accuracy values. The technology also provides clearance information with good agreement between AEMS and traditional LC-MS, as well as high-throughput metabolite monitoring, allowing decisions on candidate progression to be made quickly and accurately.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries. Echo and Echo MS are trademarks or registered trademarks of Labcyte, Inc. in the United States and other countries, and are being used under license. The images shown may be for illustration purposes only and may not be an exact representation of the product and/or the technology. Plates are available from Beckman Coulter Life Sciences. Beckman Coulter® is being used under license. © 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-12894-A. AB SCIEX™ is being used under license.


Figure 8. Clearance of imipramine. Clearance of imipramine parent (blue) overlaid with detection of hydroxy-imipramine metabolite (green) across the incubation time scale.

References

1. Chao P, Uss AS, Cheng KC. (2010) Use of Intrinsic Clearance for Prediction of Human Hepatic Clearance, *Expert Opinion on Drug Metabolism and Toxicology*, **6(2)**, 189-98.
2. Rapid MS/MS analysis with Acoustic Ejection Mass Spectrometry (AEMS). SCIEX technical note, RUO-MKT-02-11385-A.


Ultra-fast, accurate and simultaneous quantification of ritonavir and lopinavir in human plasma

Sensitive and robust results using the Echo® MS System

Rahul Baghla, Rolf Kern
SCIEX, USA

Protease inhibitors (PIs) are a class of anti-viral drugs that prevent viral replication by selectively binding to viral proteases and inhibiting their function. The development of PI-based therapies has been of enormous benefit to people infected with HIV. Used in combination with other drugs, PIs have dramatically reduced the number of people who become ill from HIV-related opportunistic infections or who die from AIDS.

Unfortunately, the effectiveness of protease inhibitors can fade over time. Mutations during viral replication can result in viruses that produce new, different proteases that are not targeted by current PI therapies¹. The best way to avoid this drug resistance is to reduce or stop HIV replication. With less HIV replication, there is less of a chance of a new strain that is resistant to anti-HIV drugs. To keep HIV levels as low as possible, PIs are typically taken in combination with at least two other anti-HIV drugs. Such combination therapies are referred to as highly active antiretroviral therapy (HAART)².

Lopinavir and ritonavir are two protease inhibitors that are often used as part of a fixed-dose combination, and serve as the model compounds in this study.

Acoustic Ejection Mass Spectrometry (AEMS), as implemented in the Echo® MS System with a SCIEX Triple Quad™ 6500+ LC-MS/MS system, offers clear benefits for quantification of lopinavir and ritonavir in human plasma. Requiring minimal sample preparation and no chromatographic separation, it provides high sample throughput without sacrificing robustness or reproducibility. The Echo® MS System combines Acoustic Droplet Ejection technology with an Open Port Interface (OPI) to deliver nanoliter sample volumes with minimal carryover. The small sample size reduces matrix suppression and facilitates easy sample preparation.


Figure 1. High reproducibility in quantification. 6 replicates at each concentration level of lopinavir and ritonavir in extracted human plasma, along with 6 replicates at quality control levels, showed excellent %CVs, below 15% at all levels.

Key features of the Echo® MS System for high-throughput bioanalysis quantification

- Ultra-fast 3 sec/sample analyses provide high throughput for quantifying large number of samples
- Liquid-liquid extraction (LLE) sample preparation ensures excellent quantitative accuracy
- Exceptional sensitivity of the Echo® MS System yields high quality data at all concentration levels

Experimental details

Sample preparation: Lopinavir and ritonavir (Sigma Aldrich) were spiked into human plasma (BioIVT) samples in the range of 0.5 ng/mL to 250 ng/mL each. Samples were processed using a liquid-liquid extraction method as follows. 0.5 mL of di-isopropyl ether was added to 0.1 mL aliquots of spiked plasma. Samples were vortexed for 1 minute followed by centrifugation at 12,000 rpm for 5 minutes. 0.4 mL of supernatant liquid was collected and dried under a nitrogen stream. Samples were reconstituted in 100 μ L of 25% v/v methanol in water and transferred to a 384-well plate for analysis by AEMS. Plates are available from Beckman Coulter Life Sciences.

Acoustic ejection method: Methanol with 0.1% v/v formic acid was used as carrier solvent at a flow rate of 425 μ L/min. 50 nL sample volumes were ejected into mass spectrometer for analysis.

Mass spectrometry: The Echo[®] MS System included a SCIEX Triple Quad 6500+ LC-MS/MS system and was controlled by SCIEX OS Software 1.6.10. The optimized MS parameters are listed in Table 1.

Table 1. Optimized MS parameters.

Name	Q1 / Q3 (m/z)	DP (V)	CE (V)	CXP (V)
Ritonavir	721.3/296.1	30	28	10
Lopinavir	629.3/447.2	30	22	10
Source parameters	Value	Source parameters	Value	
Curtain gas (psi)	25	CAD gas (psi)	10	
Ion source gas 1 (psi)	90	Ion spray voltage (V)	5000	
Ion source gas 2 (psi)	70	Source temperature (°C)	350	


Figure 2. Calibration curve for quantification of ritonavir in human plasma using the Echo[®] MS System. Good linearity and reproducibility were observed for the samples ejected in six replicates at each concentration level.


Figure 3. Calibration curve for quantification of lopinavir in human plasma using the Echo[®] MS System. Good linearity and reproducibility were observed for the samples ejected in six replicates at each concentration level.


Figure 4. XICs for transition for ritonavir and lopinavir using the Echo[®] MS System. Blank, 0.5 ng/mL and 1.5 ng/mL (A,B,C for ritonavir and D,E,F for lopinavir, respectively) are shown.

Results

Samples for the calibration curves, along with quality control samples, all analyzed in six replicates, demonstrated the high reproducibility of AEMS when combined with liquid-liquid extraction. Excellent %CVs were achieved across all concentration levels with no interference in blank human plasma samples (Figure 1). Even at the extremely short analysis time of 3 seconds per sample, the method yielded LLOQs of 0.5 ng/mL for both lopinavir and ritonavir. As summarized in Tables 2 and 3, the assay accuracy was 85.31–112.34% for ritonavir and 85.51–113.52% for lopinavir. Accuracies and %CVs were well within standard acceptance criteria for all tested samples. The

Table 2. Quantification summary for ritonavir analysis by the Echo® MS System.

Actual Concentration (ng/mL)	Calculated Concentration (ng/mL)	Accuracy (%)	CV (%)	N
0.50	0.49	97.37	8.15	6
1.00	1.06	106.14	10.63	6
1.50	1.48	98.51	5.86	6
6.00	5.56	92.66	7.51	6
12.00	13.24	110.31	6.64	6
24.00	26.72	111.34	6.49	6
120.00	134.81	112.34	6.86	6
240.00	215.70	89.88	6.58	6
400.00	341.23	85.31	4.46	6
500.00	480.76	96.15	6.49	6
Quality Control				
0.50 (LLOQ QC)	0.46	91.68	14.55	6
1.50 (LQC)	1.45	96.36	6.91	6
240.00 (MQC)	248.37	103.49	7.12	6
400.00 (HQC)	429.18	107.30	10.52	6

Table 3. Quantification summary for lopinavir analysis by the Echo® MS System.

Actual Concentration (ng/mL)	Calculated Concentration (ng/mL)	Accuracy (%)	CV (%)	N
0.50	0.50	100.02	12.02	6
1.00	1.01	101.05	9.88	6
1.50	1.46	97.35	8.33	6
6.00	5.68	94.62	5.79	6
12.00	13.62	113.52	3.26	6
24.00	26.49	110.37	3.89	6
120.00	133.63	111.36	9.46	6
240.00	225.31	93.88	8.23	6
400.00	342.02	85.51	4.23	6
500.00	466.56	93.31	5.24	6
Quality Control				
0.50 (LLOQ QC)	0.49	98.51	14.61	6
1.50 (LQC)	1.41	94.31	12.57	6
240.00 (MQC)	242.70	101.13	9.02	6
400.00 (HQC)	416.81	104.20	12.68	6

calibration curve covered approximately 3 orders of magnitude (0.5–250 ng/mL) for both analytes and displayed linearity with regression coefficients (r^2) of 0.9934 for ritonavir and 0.9946 for lopinavir using a weighting of $1/x^2$ (Figures 2 and 3).

Samples were analyzed without internal standards for this assay. Use of labeled internal standards might further improve these results.

Conclusions

- The Echo® MS System produced very sensitive, accurate and reproducible results for the simultaneous quantitative analysis of lopinavir and ritonavir in human plasma
- The very short analysis time (3 sec/sample) enabled rapid generation of quantitative data for high numbers of samples
- The assay showed great reproducibility without using labeled internal standards, but use of labeled internal standards is recommended to further improve these results

References

1. Koehn J, Ho R.JY. (2021) Novel Liquid Chromatography-Tandem Mass Spectrometry Method for Simultaneous Detection of Anti-HIV Drugs Lopinavir, Ritonavir, and Tenofovir in Plasma. *[Antimicrobial Agents and Chemotherapy](#)*, **58** (5), 2675-2680.
2. Mishra TD, Kurani H, Singhal P, Shrivastav PS, (2012) Simultaneous Quantitation of HIV-Protease Inhibitors Ritonavir, Lopinavir and Indinavir in Human Plasma by UPLC-ESI-MS-MS, *[Journal of Chromatographic Science](#)*, **50** (7), 625-635.
3. Jaiswal S, Sharma A, Shukla M, Lal J, (2017) Simultaneous LC-MS-MS Determination of Lopinavir and Rifabutin in Human Plasma, *[Journal of Chromatographic Science](#)*, **55** (6), 617-624.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Echo and Echo MS are trademarks or registered trademarks of Labcyte, Inc. in the United States and other countries, and are being used under license. Other trademarks and/or registered trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries.

© 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-13058-A. AB SCIEX™ is being used under license.

Headquarters

500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales

For our office locations please call the division headquarters or refer to our website at sciex.com/offices

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1. Rapid MS/MS analysis with Acoustic Ejection Mass Spectrometry (AEMS) - Using the SCIEX Echo™ MS System to break bottlenecks in quantitative mass spectrometry throughput. SCIEX technical note RUO-MKT-02-11385-A.

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Rapid and High Throughput Screening of Multiple Components in Fermentation Broth by SCIEX Echo™ MS System

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Key Words : SCIEX Echo™ MS System, Fermentation Broth, High Throughput Screening, Synthetic Biology

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1. Rapid quantitative analysis of fermentation broth samples to assess efficiency of engineered yeast strain turnover. SCIEX technical note RUO-MKT-02-11695-A.

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Ultra-fast quantification of alprazolam in human plasma extracts

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Mackenzie Pearson, Rahul Baghla, Rolf Kern
SCIEX, USA

Benzodiazepines are often used to treat anxiety, insomnia and central nervous system disorders. They belong to the class of depressant drugs, acting on GABA receptors to dampen neural pathway excitation and consequently provide a calming effect.¹ They are among the most common prescribed psychiatric medications. However, benzodiazepines are often accompanied by other drugs of abuse and detected alongside opioids in many overdose cases and fatalities.²

Alprazolam is the most commonly prescribed benzodiazepine for anxiety and panic disorders.³ As such, rapid, high-throughput detection of alprazolam in plasma is critical to forensic and clinical toxicology. LC-MS analysis strategies can be time consuming to develop and implement as they typically require lengthy and costly solid-phase extraction (SPE) sample preparation. Acoustic Ejection Mass Spectrometry (AEMS) offers a clear alternative for analysis of alprazolam in plasma. Requiring minimal sample preparation and no chromatographic separation, it provides increased sample throughput without sacrificing robustness or reproducibility.

The Echo® MS System with a SCIEX Triple Quad™ 6500+ Mass Spectrometer provides sensitive, high-throughput quantification of alprazolam. The system combines Acoustic Droplet Ejection technology with an Open Port Interface (OPI) to deliver nanoliter sample volumes with minimal carryover. The small sampling size reduces matrix suppression and facilitates easy sample preparation.⁴


Figure 1. Highly reproducible quantification. 6 replicates at each concentration level of alprazolam in extracted human plasma showed excellent % accuracy with %CVs below 15% at all levels.

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- Ultra-fast 2 sec/sample analyses provide high throughput for quantifying large numbers of samples
- Liquid-liquid extraction (LLE) sample preparation ensures excellent quantitative accuracy
- Exceptional sensitivity of the Echo® MS System yields high quality data at all concentration levels

Experimental details

Sample preparation: Alprazolam was spiked into human plasma in the range of 0.5 ng/mL to 250 ng/mL along with 50 ng/mL of internal standard (alprazolam D5). Samples were processed using a liquid-liquid extraction method as follows. 1 mL of di-isopropyl ether was added to 0.3 mL aliquots of spiked plasma. Samples were vortexed for 1 minute followed by centrifugation at 12,000 rpm for 5 minutes. 0.8 mL of supernatant liquid was collected and dried under a nitrogen stream. Samples were reconstituted in 100 µL of 25% v/v methanol in water and transferred to a 384-well plate for analysis by AEMS. Plates are available from Beckman Coulter Life Sciences.

Acoustic ejection method: Methanol was used as carrier solvent in this assay at 425 µL/min. 50 nL sample volumes were ejected into the mass spectrometer for analysis.

Mass Spectrometry: The Echo[®] MS System includes a SCIEX Triple Quad 6500+ System and was controlled by SCIEX OS Software 1.6.10. The optimized MS parameters are listed in Table 1.

Table 1. Optimized MS parameters.

Name	Q1 / Q3 (m/z)	DP (V)	CE (V)	CXP (V)
Alprazolam	309.2 / 205.1	55	61	9
Alprazolam D5	314.2 / 210.0	55	55	9
Source parameters	Value	Source parameters	Value	
Curtain gas (psi)	25	CAD gas (psi)	9	
Ion source gas 1 (psi)	90	Ion spray voltage (V)	5000	
Ion source gas 2 (psi)	70	Source temperature (°C)	350	

Data processing: Data processing was performed using SCIEX OS Software. A calibration curve was generated with six replicates at each concentration level to evaluate ejection reproducibility and accurately determine the lower limit of quantification (LLOQ). The calibration curve is shown in Figure 2 and a quantification summary is listed in Table 2.


Figure 2. Calibration curve for quantification of alprazolam in human plasma using the Echo[®] MS System. Good linearity and reproducibility were observed for the samples ejected in six replicates at each concentration level.

Results

Samples analyzed in six replicates demonstrated the solid reproducibility of AEMS when combined with liquid-liquid extraction. Excellent %CVs were achieved across all concentration levels with no interference in blank human plasma samples (Figure 3). Even at the extremely short analysis time of 2 seconds per sample, the method yielded an LLOQ of 0.5 ng/mL for alprazolam in human plasma. As summarized in Table 2, the assay accuracies of 92.67–105.55% and %CVs were well within standard acceptance criteria for all tested samples. The calibration curve covered 3 orders of magnitude (0.5–250 ng/mL) and displayed linearity with a regression coefficient (r) of 0.9971 using weighting of 1/x² (Figure 2).

Table 2. Quantification summary for alprazolam analysis by the Echo[®] MS System.

Actual Concentration (ng/mL)	Calculated Concentration (ng/mL)	Accuracy (%)	CV (%)	N
0.50	0.50	100.63	8.04	6
1.50	1.50	100.02	10.78	6
2.50	2.38	95.02	3.87	6
5.00	5.04	100.87	8.36	6
10.00	10.56	105.55	3.33	6
25.00	25.70	102.78	4.10	6
50.00	46.33	92.67	4.31	6
100.00	101.41	101.41	5.67	6
200.00	201.35	100.67	2.24	6
250.00	250.94	100.38	4.83	6


Figure 3. XICs for transition for alprazolam. Double blank and blank from extracted plasma (A and B), LLOQ at 0.5 ng/mL (C) and 1.5 ng/mL (D).

Conclusions

- The Echo[®] MS System produced very sensitive, accurate and reproducible results for the quantitative analysis of alprazolam in human plasma.
- The extremely short analysis times (2 sec/sample) enabled rapid generation of quantitative data for high numbers of samples.

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Echo and Echo MS are trademarks or registered trademarks of Labcyte, Inc. in the United States and other countries, and are being used under license. Other trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries.

© 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-13059-A. AB SCIEX[™] is being used under license.

References

1. Votaw VR, Geyer R, Rieselbach MM, McHugh RK. (2019) The epidemiology of benzodiazepine misuse: A systematic review. *Drug Alcohol Depend.* **200**: 95-114.
2. Moore TJ, Mattison DR, (2017) Adult Utilization of psychiatric drugs and differences by sex, age, and race. *JAMA Intern. Med.* **177**, 274–275.
3. Grohol J. (2016) *Top 25 Psychiatric Medication Prescriptions for 2013*.
4. Rapid MS/MS analysis with Acoustic Ejection Mass Spectrometry (AEMS) - Using the SCIEX Echo[®] MS System to break bottlenecks in quantitative mass spectrometry throughput. [SCIEX technical note RUO-MKT-02-11385-A](#).

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AEMS ultra-fast detection of the concentration of 17 drugs in sewage

Using the SCIEX Echo™ MS system to break bottlenecks in quantitative mass spectrometry throughput at the rate of one sample per second

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Li Zhiyuan¹, Chen Junmiao¹, Sun Xiaojie¹, Liu Bingjie¹, Guo Lihai¹

¹ SCIEX China

Key Words : Echo™ MS, AEMS, ultra-fast, quantitative, sewage, drug

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