



Electron activated dissociation (EAD)

A step change in fragmentation technology

EAD is a new paradigm for mass spectrometry. Learn how EAD provides more informed decision-making, accelerates research and development and improves productivity for routine analytical applications.



The Power of Precision

SCIEX WHITE PAPER

Contents

2	Contents
3	EAD fragmentation explained
4	Executive summary
5	Problem statement
6	The Goldilocks principle of MS/MS
7	When CID is not enough
8	EAD fragmentation in action for large molecules
9	Challenges of isomeric amino acid residues
10	Challenges of identifying and localizing glycosylation
11	N-glycosylation
12	O-glycosylation
13	Fast and comprehensive LC-MS/MS glycan characterization
14	Challenges of identifying and localizing phosphorylation
15	Challenges of disulfide-linked peptide characterization
16	EAD fragmentation in action for small molecules
17	Complete characterization of lipids in a single experiment
18	Challenges with metabolite identification
19	Challenges with small polar metabolites
20	Challenges of pesticide identification
21	Summary
22	Learn more about the Zeno revolution
23	References
24	SCIEX Now support network



EAD fragmentation explained



Executive summary

EAD is a groundbreaking approach for tandem mass spectrometry (MS/MS) applications. While the industry-standard collision induced dissociation (CID) has proven to be an invaluable tool for MS/MS experiments, the data produced using CID can leave gaps in our understanding of:


- Molecular structure
- Information on positional isomers
- Characterization of larger molecules, including post-translational modifications [PTMs]
- Structural elucidation of entire compound classes that can be difficult or even impossible to fragment

Electron-based fragmentation has been demonstrated to provide vital information that is key for complete molecule identification and characterization. Its utility extends beyond CID, providing new and crucial information even for traditionally difficult compounds.

Now, with the advent of tunable EAD, a range of different free-electron based fragmentation mechanisms are provided within one device. The availability of both low and high energy electron fragmentation within the same design extends the utility of the approach for both small [singly-charged] and large [multiply-charged] molecules alike.

Further enhancements serve to significantly increase the sensitivity and speed of fragmentation over traditional electron-based fragmentation mechanisms. This makes making the process compatible with fast separations using UHPLC and with the analysis of complex mixtures containing a range of compound masses and abundances.

Continual investments in the development of MS/MS technologies are central to delivering tools and workflows that enable the characterization of an increasingly comprehensive suite of compound classes, molecular structures and sample types. As such, EAD has proven to be a game changer that can transform MS/MS experiments and enable the acquisition of essential fragmentation data. This white paper provides an overview of the challenges encountered with current MS/MS approaches and the significant advantages that can be realized using new tunable EAD fragmentation. Examples show the power of EAD for structural elucidation of small molecules and metabolites, differentiation and quantification of isomeric compounds, identification and localization of protein modifications and complete characterization of lipids.



Problem statement

MS/MS has evolved as one of the most valuable analytical tools available to modern scientists. Introduced over 50 years ago, MS/MS enables the identification, structural elucidation and quantitation of compounds by breaking them apart and using the pieces to infer details about the original parent molecule.¹

Most MS/MS applications today use CID, for inducing fragmentation. CID is the underlying framework supporting most quantitative assays and is responsible for the identification and structural elucidation of countless numbers of compounds. But as with any approach, there are certain limitations to its applicability. With CID, those limitations manifest as insufficient fragmentation of specific molecule classes, sizes and chemistries that inhibit their characterization or selective quantification. As a result, there is clearly a need for improved fragmentation mechanisms to address the shortcomings of CID.


EAD addresses those shortcomings by providing a rich assortment of new and different diagnostic fragment ions, and by enabling the dissociation of intractable compounds previously difficult or impossible to fragment using CID. EAD encompasses a range of electron-based fragmentation mechanisms that differ by the kinetic energy of the irradiating

electron beam and the charge states of the precursor ions dissociated. Thus, EAD provides the flexibility for fragmenting both singly- and multiply-charged ions and can be applied to the characterization of both small molecules and larger biomolecules. Because EAD fragmentation is both fast and highly sensitive, structural elucidation of low-level compounds and variants is possible, even during fast chromatographic separations. Workflows that traditionally use CID are amenable to EAD, such as in-depth analysis of complex biological mixtures, with EAD now delivering new information that can clarify molecular structures.

EAD can be used for the analysis of a wide range of molecules possessing vastly different chemistries and molecular weights. Accordingly, many applications can benefit from EAD, including those within life science research, drug development, clinical studies and environmental, food and forensic analysis. The new information that EAD provides allows scientists to make faster, more informed decisions, accelerating research and development and improving productivity for routine analytical applications.

Key benefit

The ZenoTOF 7600 system with EAD combines high sensitivity with high information content to enable the characterization of important, and previously difficult to characterize metabolites.



The Goldilocks principle of MS/MS

There are a variety of MS/MS scan types for qualitative studies in use today, but the most widely used is undoubtedly the product ion scan. The general approach for a product ion MS/MS scan involves selecting precursor ions based upon their mass, fragmenting those ions and then analyzing the array of resulting product ions. This is accomplished using two stages of mass filtering or scanning with a fragmentation event in between. The masses of the product ions and the mass differences between product ions are used in determining the structure or sequence of the original parent molecule. The nature of specific chemical bonds can also be surmised from the energy required for their fragmentation.

For MS/MS to be effective, the fragmentation must be interpretable and reproducible. Additionally, the energy transferred for fragmentation must be “just right” — that is, enough to break the precursor into diagnostic product ions, but not so much that it obliterates the parent molecule into pieces so small that they lose specificity, destroying more informative higher mass fragment ions [Figure 1]. Even with the right amount of energy, the fragmentation can still be lacking in many ways. With refractory compounds, there can be a failure to produce enough fragmentation, leaving bare areas within the spectrum with no fragment ions for interpretation. Specific diagnostic fragment ions can be missing that are precisely the key for discerning fine details in the structure, such as isomers. Important side chains and modifications can be cleaved off from the parent molecule, leaving no indication of its location within the molecule. Alternatively, all of the fragmentation energy can be directed to the side chain or modification itself, again leaving little indication of its attachment point.

As a result, many different mechanisms have been used for fragmenting ions over the years, including the use of photons, electrons, atoms, molecules and even solid surfaces, with some methods more compatible than others with specific mass analyzers.² But overall, the most common and widespread techniques use atoms, molecules or electrons to impart energy for fragmentation.







Figure 1. Collision induced dissociation [CID] spectra at different collision energies. As the collision energy is increased during MS/MS, the fragments created will change with large fragments typically observed with lower CE values (top) and smaller fragments observed with higher CE values (bottom).



When CID is not enough



CID is the underlying mechanism used for most research and commercial assays in drug development, testing (food, water and environmental), forensic studies, omics research (proteomics, lipidomics and metabolomics), biomarker experiments and clinical assays. But CID does have its limitations with certain molecules and classes of compounds.


Electron-based fragmentation mechanisms can address many of the shortcomings of CID with enhanced qualitative data that can greatly extend the capabilities of MS/MS for structural elucidation. One drawback of most commercially available electron-based fragmentation devices is that they support either small molecule analysis or large molecule analysis, but not both. Electron transfer dissociation (ETD) and electron capture dissociation (ECD) both require multiply-charged precursor ions that capture low-energy electrons to induce fragmentation. Conversely, electron impact excitation of ions from organics (EIEIO), and other higher energy electron fragmentation techniques, fragment singly-charged ions.

Thus, in cases where both CID and ETD are available, a third fragmentation technique is used to fill the void, such as ultraviolet photodissociation (UVPD). For simplicity and general applicability, an ideal electron-based fragmentation device would enable a range of EAD electron energies and remove any requirement for a reagent (as in ETD).

Recently, a tunable EAD device was demonstrated on a quadrupole time-of-flight (QTOF) mass spectrometer.³ The nature of the device allows simultaneous capture of ions and electrons. The ability to tune the device allows the kinetic energy of the electron beam to be adjusted to provide the best fragmentation for the molecule of interest. Large molecule applications using low-energy electron dissociation are possible, as well as small molecule applications better suited to fragmentation using higher electron energies. Thus, both low-energy and high-energy electron voltage fragmentation experiments are possible in addition to CID.



EAD fragmentation in action for large molecules



Proteins as a class come with a rich assortment of varying structures and modifications that often exist as extremely complex and heterogeneous mixtures within biological fluids. As a result, the acquisition of MS/MS data is imperative for their characterization. CID has been used extensively for determining the structure and sequence of large biomolecules - such as peptides, saccharides and oligonucleotides-but full characterization of biomolecules and their modifications can be difficult or impossible using CID. EAD can provide a more complete picture. With EAD, fragmentation of large, multiply-charged precursor ions can be induced by the capture of lower energy electrons. Different fragment ions are produced than those typically observed with CID.⁷ For example, with peptide fragmentation, CID typically produces “b” and “y” ions while EAD produces “c” and “z” ions. These ions enable sequencing of the peptide amino acid backbone through examination of the mass differences between sequential ions within a series.



Challenges of isomeric amino acid residues




Read technical note >


Some amino acids are identical in mass — such as aspartic acid/isoaspartic acid and leucine/isoleucine — making their differentiation impossible using low-energy CID. Using EAD, however, these isomers can be identified from additional fragment ions that are produced. For example, with leucine/isoleucine, a secondary fragment ion called the w-ion, caused by further fragmentation of the backbone z-ion, can be used for identification. With leucine, the z-ion loses 43 Da from the side chain, while for isoleucine, the z-ion loses 29 Da from the side chain [Figure 8].

An example of EAD for leucine/isoleucine differentiation is shown in Figure 8. At the top, EAD clearly indicates the identity of two leucine residues within this peptide sequence through the loss of 43 Da from the z6 and z13 ions. At the bottom, the loss of 29 Da from the z5 ion identifies an isoleucine within this peptide sequence.


Differentiation of leucine and isoleucine using EAD



Leucine



Isoleucine



Key takeaway

EAD enables the differentiation of isomeric amino acids that otherwise cannot be differentiated using conventional low-energy CID.

Figure 8. EAD clearly indicates the identity of two leucine residues within this peptide sequence through the loss of 43 Da from the z6 and z13 ions [top]. At the bottom, the loss of 29 Da from the z5 ion identifies an isoleucine within this peptide sequence.



Challenges of identifying and localizing glycosylation



Protein glycosylation is another important PTM that impacts biological activity. Glycosylation can affect protein folding, protein stability, solubility and cell adhesion, is an important regulatory process. While CID of glycopeptides can provide information about the structure of the glycan modification, fragment ions that show loss of the modification still attached to the backbone residue are typically missing.

However, this specific fragment ion is necessary for localizing the attachment point of the glyco group on the peptide backbone. In contrast, EAD spectra often contain this fragment ion, along with other peptide backbone fragment ions, making localization straightforward.

Additionally, the mass difference between fragment ions on both sides of the attachment point enable calculation of the glycan molecular weight. Thus EAD and CID are quite complementary for glycopeptide analysis (Figure 12 and Figure 13).

Zeno EAD DDA of an N-glycopeptide




Figure 12. EAD provides both peptide sequence and glycan localization (c9++ ion).

Zeno CID DDA of an N-glycopeptide




Figure 13. CID only shows glycan loss or peptide backbone information, not both in the same MS/MS spectrum in a descriptive manner; localization is not possible.



N-glycosylation



N-linked glycosylation is directed by the specific consensus sequence NXS/T along the peptide backbone, which is defined as an asparagine [N], along with any residue except proline [X] and either a serine [S] or threonine [T]. This makes localization of attachment points easier as it limits the number of N residues that comply with the attachment rules.

Figure 14 shows the assignment of an N-glycosylation site on a peptide from the enzymatic digestion of alpha-2-macroglobulin using an error-tolerant database search [Mascot], a workflow typically used in protein identification and characterization experiments. Near-complete c- and z-series fragment ions pinpoint asparagine at position 1424 [N1424] as the site of glycosylation with a Hex[5] HexNAc glycan attached.

Zeno EAD DDA of an N-glycopeptide from plasma




Figure 14. EAD of an N-glycosylated peptide. EAD LC-MS/MS analysis of an N-linked glycopeptide from the tryptic digestion of alpha-2-macroglobulin shows near-complete c- and z-series fragment ions that identify and localize the type and site of modification as Hex[5] HexNAc on N1424 through the loss of N+1914 Da.



O-glycosylation



In contrast to N-linked glycosylation, no unambiguous consensus sequence exists for O-linked glycosylation, which can theoretically occur on any serine [S] or threonine [T] residue. Additionally, O-linked glycosylation can occur as dense clusters with many serines and/or threonines modified in close proximity to one another.

Figure 15 contains the assignment of an O-glycosylation site on a peptide generated from the enzymatic digestion of inter-alpha-trypsin inhibitor using a Mascot error-tolerant search. The nearly-complete c- and z-series fragment ions identify T653 as the site of glycosylation with a Hex[1] HexNAc[1] NeuAc[1] glycan attached [656 Da].

Zeno EAD DDA of an O-glycopeptide from plasma




Figure 15. EAD of an O-glycosylated peptide. EAD LC-MS/MS analysis of an O-linked glycopeptide from the tryptic digestion of inter-alpha-trypsin inhibitor shows near-complete c- and z-series fragment ions that identify and localize the type and site of modification as Hex[1] HexNAc[1] NeuAc[1] on T653 through the loss of 656 Da.



Fast and comprehensive LC-MS/MS glycan characterization




Alpha-2-HS-glycoprotein, or fetuin, is a protein that contains both N- and O-linked glycosylation. Tryptic digestion of this protein produces a complex mixture of glycopeptides, including a large 61-residue peptide with four O-linked glycosylation sites that could be modified in a wide variety of glycan combinations. The difficulty of analyzing O-linked glycopeptides and fragmenting a large 61-residue peptide that is just one of many within a complex mixture makes the characterization of fetuin glycosylation particularly challenging.

In a recent study, EAD was demonstrated for the site-specific analysis of bovine fetuin O-linked glycosylation.¹¹ While previous studies have partially characterized fetuin O-linked glycopeptides using infusion and long ETD accumulation times,¹² the recent study used UHPLC and fast EAD on a ZenoTOF 7600 system to fully characterize 57 O-glycopeptides with different glycosylation patterns and partially characterize another 22 O-glycopeptides.

Figure 11 shows an example for one O-linked glycopeptide from fetuin that contains two glycans. Here, EAD enables both the identification of the glycans and the localization to specific residues along the peptide backbone. The fetuin O-glycopeptide LC-MS/MS study is a nice demonstration of the speed with which newer EAD approaches can now be performed.

Today, fast fragmentation and MS/MS acquisition on newer instrumentation, such as the ZenoTOF 7600 system, allow EAD reaction times on the order of 10-30 ms.^{13,14} Importantly, high resolution is maintained, which is extremely important for structural elucidation, particularly for larger mass compounds. Practically speaking, this means that data dependent LC-MS/MS experiments using EAD are now compatible with fast UHPLC chromatography and can be used for the in-depth analysis of complex biological mixtures.

EAD fragmentation of an o-glycopeptide



With the Zeno trap




Figure 11. O-glycosylation requires comprehensive fragment coverage to confirm the glycan and the site localization. Higher energy EAD [hot EAD] for a complex O-glycopeptide shows near-complete sequence coverage and both glycans localized in the center portion of the peptide. [Right] Activation of Zeno trap for significant sensitivity improvements provides high peptide cleavage coverage even at fast spectrum accumulation times.



Challenges of identifying and localizing phosphorylation



Another area where EAD excels is in the characterization of PTMs on peptides. PTMs can be quite labile and fragmentation using CID can cause cleavage of these labile groups from the backbone peptide, thereby losing their attachment point and making their full characterization difficult. Alternatively, EAD fragmentation can retain these modifications. Thus, while CID can provide basic sequence information and indicate a modification exists somewhere within the peptide, EAD can identify both the type and site of a modified residue, confirming exactly which residue along the backbone is modified.⁶

Figure 10 shows the LC-MS/MS EAD spectrum of a phosphopeptide acquired using a ZenoTOF 7600 system. The location of the phosphorylated residue is easily identified and confirmed with multiple c- and z-ion series fragments and 100% sequence coverage along the peptide backbone. Because the phosphoserine residue is lost as an intact unit, rather than simply losing the phosphate group from the serine (as often happens

with CID), both the identification of serine as the modified amino acid and the location as the 10th residue from the N-terminus are provided.

Previous studies using external trapping devices for electron-based fragmentation have shown they can increase sensitivity without increasing the probability of adverse side reactions, such as neutralization, loss of larger fragment ions and production of charge reduced species.⁹ As shown in Figure 10, the spectrum acquired using the Zeno trap shows increased sensitivity versus the spectrum acquired without the Zeno trap. This added sensitivity is particularly useful when characterizing post-translationally modified peptides, such as phosphopeptides, as their abundances are typically far lower than unmodified peptides.

[Read technical note >](#)

Zeno EAD of a phosphopeptide




Figure 10. EAD spectrum of a phosphopeptide from bovine NDUFA10 subunit of Complex I [P34942] using a 7 eV electron energy. Extensive c- and z-series fragment ions are produced with 100% sequence coverage, enabling the site and type of modification to be deduced. Site differentiating ions are labeled. The spectrum on the top uses the Zeno trap to increase sensitivity while the mirror spectrum on the bottom shows the comparison without the use of the Zeno trap.

Key takeaway

Phosphorylation, which is an important regulatory modification on proteins can be both identified and localized using EAD.



Challenges of disulfide-linked peptide characterization



[Read technical note >](#)

Disulfide bonding through cysteine residues plays an important role in protein structure and function. However, identifying and locating disulfide bonds, including the peptide sequences that surround both cysteine residues, can be a difficult task using traditional CID LC-MS/MS approaches. With CID, often the peptides will fragment along their backbones until reaching the cysteine residues, providing only partial information for the peptide sequences. This also makes determining the exact location of the cysteines and disulfide bond along the peptide backbone more difficult to identify.

Figure 9 compares the LC-MS/MS spectrum obtained using EAD versus CID for a disulfide-bound peptide from the hinge region of a monoclonal antibody (mAb). Fragmentation of this precursor using CID produces both b- and y-type fragment ions along both peptide backbones, up to the cysteine residue.

At this point, fragmentation stops, as CID typically does not cleave the disulfide bond. This leaves the residual peptide backbones still connected and an absence of data for complete sequencing. In contrast, EAD fragmentation shows a preference to cleave the disulfide bond, leading to superior sequence coverage that enables the identification of the linked peptides and locations of cysteine residues.

Key takeaway

EAD fragmentation of disulfide-linked peptides exhibits a preference to cleave the disulfide bond, leading to superior sequence coverage and peptide identification.

EAD and CID comparison of a disulfide peptide






Figure 9. EAD and CID LC-MS/MS of a disulfide-bound peptide. The CID spectrum (top) contains fragment ions up to the cysteine residue, leaving a lack of information for the complete peptide sequences. In contrast, EAD (bottom) dissociates the disulfide-bond linkage, providing greater sequence coverage along both peptide backbones.



EAD fragmentation in action for small molecules

Research and development labs are continuously challenged to characterize small molecules. CID is typically the first choice for MS/MS experiments but can sometimes leave a puzzling absence of information when it is most needed. CID can produce limited or non-selective fragmentation, leading to inconclusive spectra and sub-par quantitative assays. With small molecules, the compounds that make up different molecule classes — such as pesticides, metabolites or lipids — can be extremely diverse in size, polarity and solubility. This problem is often compounded by sensitivity challenges.



Complete characterization of lipids in a single experiment




Lipids are a complex group of compounds with subtypes that share a similar high-level structure. For example, triglycerides consist of a glycerol group bonded to three long hydrocarbon chains with additional functional head groups attached in some cases. Small but meaningful differences between lipid species — such as the location of a single double bond along the hydrocarbon chain—can have important implications for health and disease. In one experiment, EAD provides all of the information for complete lipid characterization that normally requires multiple technologies and experiments (Figure 3). Complete characterization of lipids involves the identification of:

- Head group
- Backbone
- Regioisomerism
- Double bonds


[Read technical note >](#)

Single-experiment lipid characterization




Fragmentation patterns

1. EAD produces diagnostic fragments of:
 - Head group
 - PC (SM), PE, PS, PI, PA, PG



2. Backbone (glycerol/sphingoid)
 - Regioisomerism (sn-1, sn-2, etc.)



3. Chain structure
 - Length and double bonded positions
 - Double-bond stereochemistry




Figure 3. Single-experiment lipid characterization. The complete de novo identification and characterization of the lipid PC 16:0 / 18:1(n-9:cis) in a single experiment.



Challenges with small polar metabolites



The numbers and types of small molecules that make up the metabolome for any given biological system are vast and complex. Although these compounds span a range of different structures, sizes, chemistries and abundances, the metabolome is also composed of a multitude of compounds that differ only slightly in their makeup. As a result, metabolomics studies are particularly challenged in trying to identify and quantify the vast range of very different compounds and differentiate between extremely similar compounds within one sample.

The CID and EAD spectra for 3,5-cAMP are shown in Figure 4. The power of EAD is deeper, more rich fragmentation spectra

for structural elucidation/identification. For cAMP, utilizing EAD resulted in an increase from 6 fragment ions to 50 fragment ions.

In addition to increasing in fragment coverage, the alternative selectivity with EAD can be utilized for quantification (Figure 4). While in this case, the two isomers of cAMP are fully resolved chromatographically, for 3,5-cAMP, a fragment ion with nominal mass of m/z 165 is generated and when it is extracted as an XIC, it shows no interference from the 2,3 isomer (Figure 5). The unique diagnostic fragment ions produced by EAD leave no doubt as to the precursor identity.

Comparison of CID and EAD




Figure 4. CID spectrum (top) of 3,5-cAMP, EAD spectrum (bottom).

Quantitative extraction of EAD fragment ions




Figure 5. EAD extracted ion chromatograms (XICs) of 3,5-cAMP (top). Unique fragment ion for nominal mass 165 for 3,5-cAMP (bottom)



Challenges of pesticide identification



There are over 1,000 pesticides in use around the world today with many different chemistries and toxicities. Human toxicity will depend upon several factors, including the chemical nature of the pesticide, its metabolism within the body, exposure levels and use. Because of the inherent poisonous nature of pesticides, they can also be harmful to the environment if they come in contact with, susceptible plants and animals. Additionally, contact with sunlight, water and other chemicals can cause certain pesticides to degrade, leading to even more toxic breakdown products within the environment.

While approved pesticides have acceptable residual levels that are set by regulatory agencies, there are no acceptable levels for banned pesticides. Testing for residual pesticides is an important task to help ensure the safety of human and animal food products, field crops [such as cotton and grass] and the environment.

Structure and CID spectrum for cyhalofop-butyl




Figure 6. The CID spectrum is shown for cyhalofop-butyl with one major and two minor fragment ions generated.

Because there are so many different types of pesticides, multi-residue testing is a common workflow for pesticide analysis. Library matching is typically used to identify and confirm specific pesticides. Figure 6 and Figure 7 compare the EAD and CID LC-MRM^{HR} spectra for the fungicide cyhalofop-butyl. While CID only produces 2 fragment ions, EAD produces a much richer spectrum with over 60 fragment ions that can be used for vastly improved specificity for library matching and ion ratios.

Qualitative and quantitative EAD




Figure 7. The EAD spectrum [top] is shown with multiple (>60) fragment ions observed. Using EAD with many more choices for targeted MRM^{HR} is shown in the XIC plots [left], compared with only one minor and one major fragment for targeted MRM^{HR} using CID [right].



Summary

EAD is an exciting and innovative new technique for MS/MS. The development of a fast, sensitive and tunable EAD device enables the acquisition of a richer data set versus traditional CID approaches. With EAD, critical information is provided that is often missing from CID MS/MS experiments. This information leads to the characterization of previously intractable compounds, greater clarity of molecular structures, differentiation of isomers, localization of modifications and identification of low-level variants. Entire application areas and compound classes are now accessible due to the benefits that MS/MS can provide through EAD.

The electron energies used in the novel tunable EAD device on the ZenoTOF 7600 system can be adjusted to suit the precursor of interest. Thus, both singly-charged precursors (such as lipids and metabolites) and multiply-charged precursors, (such as peptides) can be analyzed. The ability to rapidly adjust the electron energy makes EAD compatible with fast UHPLC chromatography. Additionally, the sensitivity provided by the Zeno trap of the ZenoTOF 7600 system enables in-depth analysis of complex mixtures containing compounds with a wide variety of abundances, molecular weights and chemistries. These developments will allow more scientists and researchers to take advantage of the benefits of electron-based fragmentation for structural elucidation and quantitative assays.

ZenoTOF 7600 system

Ready to learn more?

VISIT

sciex.com/zenorevolution

to get your inside look at
EAD in action with the
ZenoTOF 7600 system.



Learn more about the Zeno revolution

The cover features a green-tinted photograph of a Zeno trap mass filter. The text is overlaid in white. The title 'Zeno trap' is in a large, bold font. Below it, the subtitle 'Defining new levels of sensitivity without compromise' is in a smaller font. At the bottom, there is a white rounded rectangle containing the text 'SCIEX WHITE PAPER' and a link 'Access here' with a right-pointing arrow.

Zeno trap

Defining new levels of sensitivity
without compromise

SCIEX WHITE PAPER

[Access here](#)

The cover features a green-tinted photograph of a mass spectrometer's detector showing a complex pattern of ions. The text is overlaid in white. The title 'Data independent acquisition (DIA)' is in a large, bold font. Below it, the subtitle 'Innovations in DIA-MS' is in a smaller font. At the bottom, there is a white rounded rectangle containing the text 'SCIEX WHITE PAPER' and a link 'Access here' with a right-pointing arrow.

**Data independent
acquisition (DIA)**

Innovations in DIA-MS

SCIEX WHITE PAPER

[Access here](#)



References

1. Bruins, A.P., Covey, T.R., Henion, J.D. [1987] Ion spray interface for combined liquid chromatography/atmospheric pressure ionization mass spectrometry, *Anal. Chem.*, 59 [22], 2642-2646.
2. Chernushevich, I.V., Merenbloom, S.I., Liu, S., Bloomfield, N. [2017] A W-geometry ortho-TOF MS with high resolution and up to 100% duty cycle for MS/MS, *J. Am. Soc. Mass Spectrom.*, 28, 2143-2150.
3. Covey, T.R., Lee, E.D., Bruins, A.P., Henion, J. D. [1986]. Liquid chromatography/mass spectrometry. *Anal. Chem.*, 58 [14], 1451A-1461A.
4. Whitehouse, C.M., Dreyer, R.N., Yamashita, M., Fenn, J.B. [1985] Electrospray interface for liquid chromatography and mass spectrometers, *Anal. Chem.*, 57, 675-679.
5. Fenn, J.B., Mann, M., Meng, C.K., Wong, S.F., Whitehouse, C.M. [1989] Electrospray ionization for mass spectrometry of large biomolecules, *Science*, 246, 64-71.
6. Douglas, D.J., French, J.B. [1992] Collisional focusing effects in radio frequency quadrupoles, *J. Am. Soc. Mass Spectrom.*, 3:4 398-408.
7. Chernushevich, I.V. [2000] Duty cycle improvement for a quadrupole time-of-flight mass spectrometer and its use for precursor ion scans, *Eur. J. Mass Spectrom.* 6, 471-479.
8. Loboda, A.V., Chernushevich, I.V. [2009] A novel ion trap that enables high duty cycle and wide m/z range on an orthogonal injection TOF mass spectrometer, *J. Am. Soc. Mass Spectrom.*, 20, 1342-1348.
9. Antonoplis, A. et al. Over 40% more proteins identified using Zeno MS/MS. SCIEX technical note, RUO-MKT-02-13286-B.
See technical note here >
10. Hunter, C.H. Large-scale, targeted, peptide quantification of 804 peptides with high reproducibility using Zeno MS/MS. SCIEX technical note, RUO-MKT-02-13346-A.
See technical note here >
11. Barbetti, F., Yang, C., D'Addonna, D., Klaas, C. [2020] Pesticide residues screening and quantification analysis in olive oil using an Orbitrap Exploris 240 HRMS, Thermo Fisher Scientific application note, 65901.
See technical note here >
12. Di Lorenzo, R.A et al, Highly sensitive quantification and selective identification of pesticides in food with Zeno MRMHR. SCIEX technical note RUO-MKT-02-13902-A.
See technical note here >
13. Baba, T., Ryumin, P., Duchoslav, E., Chen, K., Chelur, A., Loyd, B., Chernushevich, I. [2021] Dissociation of biomolecules by an intense low-energy electron beam in a high sensitivity time-of-flight mass spectrometer. *J. Am. Soc. Mass Spectrom.* 32, 8, 1964-1975.
14. Chernushevich, I.V., Merenbloom, S.I., Liu, S., Bloomfield, N. A W-geometry ortho-TOF MS with high resolution and up to 100% duty cycle for MS/MS. *J. Am. Soc. Mass Spectrom.*, 28, 2143-2150 [2017].
15. Zeno trap: Defining new levels of sensitivity. [2021] SCIEX white paper, RUO-MKT-19-13373-B.
16. Wang, Z. et al. High-throughput proteomics of nanogram-scale samples with Zeno SWATH DIA. BioRxiv preprint: <https://www.biorxiv.org/content/10.1101/2022.04.14.488299v1.full>
17. Qualitative flexibility combined with quantitative power: Using the ZenoTOF 7600 system, powered by SCIEX OS software [2021] SCIEX technical note, RUO-MKT-02-13053-B. [edited]



SCIEX Now support network

SCIEX Now

- Manage your instruments
- Submit and manage support cases, track status and view history
- Access online training courses and articles
- Manage software licenses linked to your registered instruments
- View and report critical instrument statistics when connected to the StatusScope remote monitoring service
- Be a part of the SCIEX community by submitting questions and comments
- Receive notifications from SCIEX with content

→ CONTACT SCIEX NOW

SCIEX Learning Hub

Learning Hub success programs provide LC-MS and capillary electrophoresis (CE) training customized to meet your exact needs.

With a selection of training methods and certifications available, you can build a mass spectrometry learning program that is most suited to your lab and users.

Starting with a clear understanding of your desired learning outcomes, we help you improve lab productivity and consistency by designing and delivering a program that is focused on knowledge advancement and retention.

→ FIND OUT MORE

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use, Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures. Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries [see sciex.com/trademarks]. © 2023 DH Tech. Dev. Pte. Ltd. RU0-MKT-19-13373-B

Headquarters
500 Old Connecticut Path
Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please
call the division headquarters
or refer to our website at
sciex.com/offices

