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
提高二级灵敏度对非靶向代谢组学工作流程的影响.....	8
SCIEX ZenoTOF™ 7600 系统 Zeno Trap 功能在代谢组学分析中的应用.....	12
Improved metabolite identification using the ZenoTOF 8600 system to analyze NIST SRM 1950 plasma by DDA analysis.....	16
Untargeted data-dependent acquisition (DDA) metabolomics analysis using the ZenoTOF 7600 system.....	23
SWATH DIA and library driven compound identification revealed upregulation of flavonoids in Ocimum basilicum L. flowers.....	32
Analysis of untargeted metabolomics data from an untargeted Zeno data-dependent acquisition (DDA) workflow using SCIEX OS software.....	38
MEExplorer Ultimate 新一代非靶向代谢组学数据分析软件助力痛风疾病研究.....	48
MS-DIAL software parameters for processing untargeted metabolomics data acquired on the ZenoTOF 7600 system.....	53
全球天然产物社会分子网络系统 (GNPS) 结合 SCIEX 高分辨质谱系统进行天然产物研究.....	59
Zeno MS/MS 技术在大鼠尿液靶向代谢组学研究中的应用.....	63
利用 SCIEX ZenoTOF™ 7600 质谱建立玉米中小分子代谢物的分析方法.....	68
基于高分辨质谱 X500R QTOF 系统油菜籽的代谢组和脂质组学分析.....	75

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Large scale targeted metabolomics assay for quantitative plasma profiling.....	82
Significant sensitivity increases provides 30% more polar metabolites quantified in plasma.....	87
ME/TED: 新一代靶向代谢组学数据分析软件及其应用.....	91
Targeted metabolomics in rat urine using Zeno MS/MS.....	95
液相色谱串联质谱法测定 20 种常见氨基酸的非生化方法.....	100
色氨酸代谢物的 LC-MS/MS 定量分析.....	103
LC-MS/MS 法测定两种甲基化修饰组氨酸.....	106

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Quantitative and qualitative bile acid analysis on the ZenoTOF 8600 system using EAD..... 207

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1. Qualitative flexibility combined with quantitative power. SCIEX technical note, RUO-MKT-02-13053-A.
2. Automated Targeted Screening of Hundreds of Metabolites -Using an Accurate Mass Metabolite Spectral Library with SCIEX OS Software. SCIEX technical note, RUO-MKT-02- 2201-B.
3. What is Zeno-On-Demand and what does it do? SCIEX community post, RUO-MKT-18-13325-A.

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Zeno Trap Application in Metabolomics

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1. Impact of increased MS/MS sensitivity on the untargeted metabolomics workflow. SCIEX technical note, RUO-MKT-02-13434-ZH-A.
2. Complete structural elucidation of lipids in a single experiment using electron activated dissociation (EAD). SCIEX technical note, RUO-MKT-02-13050-ZH-B.
3. Zeno trap white paper. SCIEX technical paper, RUO-MKT-19-13373-B.

Improved metabolite identification using the ZenoTOF 8600 system to analyze NIST SRM 1950 plasma by DDA analysis

Paul RS Baker¹, Suya Liu², Robert Proos¹, Jason Causon², Rebekah Sayers³, Sahana Mollah¹, Elliott Jones¹,

¹ SCIEX, USA, ²SCIEX, Canada, ³SCIEX, UK,

This technical note introduces a novel high-resolution mass spectrometer, the ZenoTOF 8600 system, which was used to identify metabolites in human plasma using an untargeted mass spectrometry approach.

Untargeted metabolomics using a data-dependent acquisition (DDA) scan mode is a common approach to studying the metabolites in biological samples. The experiment is based on a cycle of identifying potential precursor molecules via an MS1 scan and selecting the top candidates for product ion analysis to generate spectra for database matching to identify compounds. To improve the identification of metabolites, two instrument qualities are equally essential: high sensitivity and fast acquisition speed. High sensitivity increases the likelihood of triggering high-quality product ion spectra for low-concentration analytes, and a faster acquisition allows for more DDA cycles during the experiment, improving the chance that the low-level metabolites will be analyzed.

The ZenoTOF 8600 system is a high-resolution QTOF mass spectrometer with a fast analytical speed (~100 MS/MS events/s), and an improved front end to enable a high ion influx paired with a novel optical detector to increase the instrument's sensitivity ~10-fold compared to the ZenoTOF 7600 system, as determined by signal-to-noise [S/N]. Here, human plasma was analyzed using the ZenoTOF 7600 and 8600 instruments to evaluate the improved identification of analytes in untargeted metabolomics analysis on human plasma using reversed-phase (RP) chromatography (Figure 1).

Key features of untargeted metabolomics analysis using the ZenoTOF 8600 system

- Using untargeted analysis, the SCIEX ZenoTOF 8600 system identified 1.4-fold more metabolites than the ZenoTOF 7600 system in human plasma extract (274 vs. 203, analyzing 1 μ L plasma equivalents)
- The ZenoTOF 8600 system demonstrated ~10-fold greater sensitivity at the TOF MS and TOF MS/MS levels than the ZenoTOF 7600 system


Figure 1. A comparative study to evaluate the identification of metabolites using DDA analysis. Human plasma extracts were analyzed on the ZenoTOF 7600 or 8600 system using a top-10 DDA experiment in the positive ion mode. Data were processed using MS-DIAL 5.5 software to identify reference-matched metabolites. Due to its greater sensitivity, data acquired on the ZenoTOF 8600 system enabled the identification of ~1.4-fold more metabolites than the ZenoTOF 7600 system.

Introduction

Metabolomics research has heavily relied on untargeted MS-based analytical methods [i.e., DDA] to detect and potentially quantify small biomolecules within a sample [1,2]. High-resolution mass spectrometry (HRMS) analysis using a DDA scan mode is the primary tool for this type of experiment, which is based on an initial survey scan, such as a TOFMS scan, to select potential precursor molecules for product ion analysis. Individual MS/MS scans are acquired for each chosen precursor, and the cycle repeats itself throughout the whole experiment. Algorithms are in place to minimize repetitive measurements of the same compound, but multiple MS/MS spectra for an individual compound are typically acquired. The time for each cycle depends on the number of dependent MS/MS scans, which must be balanced with the chosen chromatographic separation strategy and the speed and sensitivity of the instrument to maximize coverage.

The primary criterion for DDA precursor ion selection is intensity—ions will be chosen in decreasing order from the most intense peak within the survey scan spectrum, provided it has not been excluded based on the experimental exclusion criteria. A plasma-derived metabolomics sample typically has hundreds of precursor ions in each survey scan, so the DDA duty cycle must be as short as possible, the exclusion criteria algorithm carefully balanced, and the instrument must be sensitive for low-level analytes to be present in the spectra for optimal untargeted analysis. The ZenoTOF 8600 system can execute an MS/MS scan in < 10 ms, enabling a scan rate of ~100/s. This rate for a top-10 DDA experiment results in a duty cycle of 600 ms [100 ms for the survey scan and 500 ms for the 50 dependent MS/MS scans]. That means for a 6-second-wide eluting peak, there is a strong chance that an MS/MS spectrum would be acquired for any given peak, if it is detected at the precursor ion level. This rate also suggests that for a 20 min gradient, there could be > 40,000 MS spectra acquired. This is not the case. For the ZenoTOF 7600 system, analyzing 0.2 μ L plasma equivalents, the number is ~9,000. This means that provided the DDA exclusion criteria were optimally set, the detection capacity of the instrument exceeds the number of detected precursors. It is this aspect of untargeted metabolomics that can benefit from improved instrument sensitivity.

Two instrument qualities are essential for the accurate and comprehensive identification of metabolites using DDA: a high speed of analysis and high sensitivity. Better sensitivity increases the likelihood of triggering product ion analysis for low-concentration analytes, and a faster analysis speed allows for more DDA cycles, which improves the chance that the low-level metabolites will be analyzed. The ZenoTOF 8600 system is a high-resolution QTOF mass spectrometer with fast analytical speed, a reconfigured front end to increase ion flux, and a novel optical detector that combine to increase the instrument's sensitivity ~10-fold compared to the ZenoTOF 7600 system. It is hypothesized that the increased instrument sensitivity will result in an increase in the number of triggered MS/MS events.

MKT-34922-A

Furthermore, the MS/MS spectra acquired from the low-level metabolites will be of better quality due to the better sensitivity. Both improvements are expected to increase the number of reference-matched metabolites.

The NIST SRM 1950 human plasma standard was used in these studies to provide a common reference for comparison. Recently, a comprehensive metabolomic/lipidomic profile of the NIST 1950 standard was performed using high-resolution NMR spectroscopy [700 MHz], direct injection tandem mass spectrometry (DI-MS/MS), liquid chromatography-tandem mass spectrometry (LC-MS/MS), and inductively coupled plasma mass spectrometry (ICP-MS) [3]. Using these diverse methods, the authors identified and quantified ~ 500 metabolites [lipid molecular species were also quantified, bringing the target list to 1058]. The mass spectrometry methods were targeted, which is advantageous in terms of quantitative sensitivity, but they lacked the ability to characterize known- and unknown-unknown compounds. Untargeted analysis is better suited to this type of discovery-based analysis.

A significant challenge associated with untargeted metabolomics is data processing. As is the case for lipidomics data processing, robust interpretation of untargeted metabolomics data is still, effectively, in its early development. For the results presented here, MS-DIAL 5.5 software was used to match DDA-derived MS/MS spectra to its compound database. Because of the diverse programs available to process DDA data and the data themselves having been acquired uniquely despite the sample being the same, direct comparisons of data presented here with literature-based studies should be approached with caution and not necessarily considered an “apples-to-apples” comparison.

Here, human plasma was analyzed using the ZenoTOF 7600 and 8600 instruments to compare the extent and quality of analyte identification in untargeted metabolomics analyses. Due to its improved sensitivity, the ZenoTOF 8600 system identified ~ 1.4-fold more analytes than the ZenoTOF 7600 system under similar analytical conditions and generated superior quality scores, as defined by the MS-DIAL 5.5 software match score.

Materials and methods

Sample preparation: NIST SRM 1950 plasma was extracted using 4 volumes of ice-cold methanol, vortexing for 10 s, centrifuging at 15,000 \times g for 10 min, and decanting the supernatant. The supernatant was dried using a speedvac, and the metabolites were resuspended in water to a final concentration of 1 μ L extract = 0.1 μ L plasma equivalents. The supernatant was directly analyzed by high-performance liquid chromatography electrospray ionization tandem mass spectrometry (HPLC-ESI-MS/MS).

Chromatography: Samples were analyzed using an Exion LC system with a Kinetex F5 column (2.1 × 150 mm, 2.6 μm, Phenomenex). A simple linear gradient from 0 to 95% B was used with standard reversed phase mobile phases at a flow rate of 200 μL/min. Mobile phase A was 0.1% formic acid in water, and mobile phase B was 0.1% formic acid in acetonitrile. The HPLC rinse solvent was water/methanol/iso-propanol/acetonitrile [1:1:1:1, v/v]. The flow rate was 0.200 μL/min, and the gradient conditions are shown in **Table 1**. A 1 to 10 μL injection volume was used [0.1 to 10 μL plasma equivalents], and the column temperature was maintained at 30°C throughout the analysis. The total runtime was 23 min.

Table 1 Chromatographic gradient

Time (min)	Mobile phase A (%)	Mobile Phase B (%)
0	100	0
2	100	0
14	5	95
16	5	95
18.5	100	0
23	100	0

Mass spectrometry: Metabolomics analysis on human plasma extracts was performed on two instruments: a ZenoTOF 7600 system and a ZenoTOF 8600 system, both equipped with an Optiflow Turbo V ion source and an electrospray ionization (ESI) probe. Instrument calibration was maintained using the automated calibrant delivery system (CDS), which calibrated every five samples with an ESI calibration solution specific for the positive ionization mode. DDA experiments were performed using CID-based fragmentation in the positive ion mode.

The systems were configured for CID-based DDA experiments to select the top 10 most abundant ions for fragmentation. Dynamic background subtraction (DBS) with a mass tolerance of 50 mDa was

Table 2. Instrument parameter settings

Parameter	ZenoTOF 7600 system	ZenoTOF 8600 system
Curtain gas (CUR)	40 psi	40 psi
Ion source gas 1 (GS1)	40	40
Ion source gas 2 (GS2)	60	60
CAD gas (CAD)	7	7
Source temperature (TEM)	500	500
Ion spray voltage (IS)	4000 V	4000 V
Declustering potential (DP)	80 V	80 V
Collision energy (CE)	30 V	30 V
Zeno threshold	100,000 cps	100,000 cps
CE (V) - TOFMS	10 V	12 V
CE (V) - TOFMS/MS	30 V	30 V
Declustering potential (DP)	40 V	40 V
TOFMS accumulation time (Acc)	100 ms	100 ms
TOFMS/MS accumulation time (Acc)	50 ms	50 ms
Time bins to sum	6	6

MKT-34922-A

applied to both experiments to minimize noise and maximize the MS/MS quality. Once a precursor ion was selected and fragmented, it was dynamically excluded from candidate selection for 6 s. The TOF MS accumulation time was set to 100 ms, and the accumulation time for the dependent TOFMS/MS analysis was 50 ms. The TOF MS and TOFMS/MS instrument parameter settings are shown in **Table 2**. A detailed description of the ZenoTOF 7600 system instrument parameters and their relevance to metabolomics DDA experiments has been previously published [9]; the parameter descriptions therein are applicable to the ZenoTOF 8600 system.

To test the ZenoTOF 8600 system's sensitivity improvement, a high-resolution multiple reaction monitoring (MRM^{HR}) scan mode was employed. From the untargeted data, 14 analytes were selected to cover a wide range of masses and peak intensities. The same plasma sample was used, and n=5 replicates were run using the MRM^{HR} scan mode. Instrument parameter settings were the same as those used for the DDA experiments (**Table 2**).

Data processing: All DDA data were processed using MS-DIAL 5.5 software [4]. Optimal parameter settings used for MS-DIAL software processing are presented in a recent SCIEX technical note [5], with some exceptions: the DOT product score and the reverse DOT product score thresholds were set to 500, and the TOFMS threshold was set to 100 counts per second (cps). Metabolite coverage determined from data generated by each instrument is given as reference-matched identifications; the total number of spectra acquired far exceeds these numbers. (For a detailed explanation of match scores, see reference 6.) Quantitative data acquired using the MRM^{HR} scan mode were processed using the Analytics module in SCIEX OS software; Qualitative data were visualized using the Explorer module.

Results

Quantitative performance of the ZenoTOF 8600 system

The ZenoTOF 8600 system was designed with an optical detector that can process more ion current than the detector used in the ZenoTOF 7600 system. This ability, along with the upgraded front-end ion optics, imparts greater sensitivity to the instrument. To measure the impact of improved sensitivity on metabolomics analysis, two different approaches were used: a comparison of the peak areas [and their respective S/N ratios] from a small set of metabolites as measured in the TOFMS dimension, and a targeted MRMhr experiment to measure sensitivity improvements in the MS2 dimension.

Figure 2 shows the combined TICs [TOFMS in blue and dependent product ion scans in pink] for the ZenoTOF 7600 and 8600 systems [top and bottom panels, respectively]. A cursory inspection of each instrument's TICs indicates an approximately 10-fold increase in ion intensity for the ZenoTOF 8600 system, reflecting the improved sensitivity of the instrument. Because the two instruments perform at

the same data acquisition speed, the improved sensitivity is likely responsible for the significant increase in the number of features acquired using the ZenoTOF 8600 system compared to the ZenoTOF 7600 system [25413 vs. 11002 for 1 μ L human plasma equivalents; Fig. 1]. The improved sensitivity also enabled more triggered dependent scans on the ZenoTOF 8600 system [6064 vs. 3932 MS/MS spectra acquired using 1 μ L human plasma equivalents; Fig. 1].

The increased peak intensities acquired on the ZenoTOF 8600 system do not necessarily mean the instrument is more sensitive. Sensitivity improvements should be reported in terms of a relative increase in the S/N. Because DDA experiments are not quantitative at the TOF MS/MS level, sensitivity was assessed by comparing the TOF MS XIC peak intensities and their concordant S/N from metabolites detected on each instrument. An evaluation of the peak areas with S/N calculations for 5 example metabolites is presented in Table 3. The S/N values were calculated using SCIEX OS software using the Explorer data analysis module and represent the average value for each analyte from n=5 injections. The noise region from each TOFMS XIC used to calculate the S/N was selected using the following criteria: the "blank"


Figure 2. Untargeted [DDA] analysis of metabolites in human plasma on the ZenoTOF 7600 and 8600 systems using DDA analysis. The top panel shows a combined TIC for TOFMS (blue) and dependent product ion scans (pink) acquired on the ZenoTOF 7600 system. The lower panel shows the same for data acquired on the ZenoTOF 8600 system. The overall intensities of the TICs from the ZenoTOF 8600 system are > 10-fold higher than those acquired on the ZenoTOF 7600 system. The number and intensities of the triggered product ion scans are significantly higher using the ZenoTOF 8600 system—6064 vs. 3932 dependent MS/MS scans on the ZenoTOF 7600 analyzing 1.0 μ L plasma equivalents. The 1.4-fold increase in the number of dependent scans contributed to the greater number of identified metabolites.

Table 3. TOFMS-derived S/N ratios for 5 example metabolites analyzed by DDA analysis

ID	m/z	8600 ZenoTOF system		7600 ZenoTOF system		Fold increase [8600 vs. 7600]
		Peak area avg	S/N avg	Peak area avg	S/N avg	S/N
Arginine	175.1185	1.18E+06	14493	6.611E+05	3316	4.4
Citruline	176.1037	1.37E+05	2700	2.884E+04	222	12.1
Hypoxanthine	137.0462	1.85E+06	2122	2.215E+04	46	46.5
Pantothenate	220.1174	2.32E+04	866	9.556E+03	95	9.2
Hippuric acid	180.0655	6.40E+04	773	1.372E+04	98	7.9

MKT-34922-A

Table 4. MRMhr-derived peak areas [with respective S/N ratios] for 12 example metabolites

Metabolite	ZenoTOF 8600 system		ZenoTOF 7600 system		Average fold-improvement using the ZenoTOF 8600 system	
	Peak area avg (n=5)	S/N avg (n=5)	Peak area avg (n=5)	S/N avg (n=5)	Peak area	S/N
Tyrosine	1.66E+07	1220	2.53E+06	236	6.5	5.2
Naproxen	1.54E+06	2311	9.32E+04	351	16.5	6.6
Piperine	2.81E+06	4514	2.41E+05	1073	11.6	4.2
Octonoyl carnitine	1.33E+06	2765	2.89E+04	111	46.0	25.0
Methionine	3.36E+06	3060	4.60E+05	832	7.3	3.7
Bilirubin	2.05E+07	153222	1.03E+07	6724	2.0	22.8
Thyroxine	2.68E+05	791	1.70E+04	172	15.8	4.6
Sphinganine-1-phosphate	4.95E+04	234	3.59E+03	6	13.8	42.4
Gabapentin	6.34E+04	111	9.44E+03	11	6.7	10.4
Cortisone	3.93E+04	365	3.51E+03	19	11.2	18.8
Indole-3-carbinol	1.54E+05	43	4.05E+03	4	38.0	9.7
Proline	9.73E+07	12424	1.44E+07	2454	6.7	5.1

Annotation: Caffeine
Adduct type: [M+H]⁺
RT [min]: 6.836(ref=4.000)(diff=2.84)
m/z: 195.08835(ref=195.08765)(diff(mDa)=0.705)
Peak height: 276030
Peak area: 1453274
Formula|Ontology: C₈H₁₀N₄O₂ | Xanthines
InChIKey: RYVVLZVUVUJUGH-UHFFFAOYSA-N
Comment: NA

Annotation: Arginine
Adduct type: [M+H]⁺
RT [min]: 1.628(ref=-1.000)(diff=2.63)
m/z: 175.11848(ref=175.11895)(diff(mDa)=0.474)
Peak height: 139620
Peak area: 1193768
Formula|Ontology: C₆H₁₄N₄O₂ | L-alpha-amino acids
InChIKey: ODKSPYDXXOIFQJN-UHFFFAOYSA-N
Comment: NA

Annotation: 3-Formylindole
Adduct type: [M+H]⁺
RT [min]: 6.249(ref=-1.000)(diff=7.25)
m/z: 146.05912(ref=146.06004)(diff(mDa)=0.926)
Peak height: 110846
Peak area: 636873
Formula|Ontology: C₉H₇NO | Indoles
InChIKey: OLNJUISKQQNIM-UHFFFAOYSA-N
Comment: NA




Figure 3. Spectral matching by MS-DIAL 5.5 software. Presented here are three examples of metabolites identified in human plasma by MS-DIAL 5.5 software. More information is available within the software's user interface, but highlighted here are the acquired MS/MS spectra (in blue) compared to the reference spectra (inverted, in red).

chromatographic region was > 1 min in length and within 2 min of the targeted peak elution time, and the region needed to be devoid of any extraneous isobaric peaks. The S/N of the TOFMS XIC peaks was improved by an average of 16-fold [range = 4.4- to 46-fold] on the ZenoTOF 8600 system compared to the ZenoTOF 7600 system.

To compare the sensitivity of the two instruments at the MS/MS level, an MRMhr experiment was performed on 12 representative metabolites that cover a broad mass range and diverse peak intensities and areas. Table 4 shows an overall increase in peak areas and calculated s/n values for all compounds, with an average peak

area increase of 15.5-fold [range of 2 to 46-fold] and an average S/N value increase of 13.2 [range of 3.7 to 42-fold]. These data indicate that the ZenoTOF 8600 system is more sensitive than the ZenoTOF 7600 system for metabolomics studies.

Identification of metabolites using MS-DIAL 5.3 software

The improved sensitivity of the ZenoTOF 8600 system and the increased number of features should translate to better coverage of the metabolome from untargeted metabolomics experiments. In general, better sensitivity allows for detecting less prominent

fragment ions, improving the quality of the product ion spectra, and leading to a better match and more confident metabolite identification. To test this hypothesis, DDA experiments were performed on the ZenoTOF 7600 and 8600 systems to compare the numbers of identified metabolites in human plasma lipid extracts.

The MS-DIAL user interface provides a comprehensive view of the processed DDA data. **Figure 3** shows an abridged example of data that is generated for reference-matched metabolites. Of note is the comparison of experimental data [lower panels, blue spectra] to reference data [red spectra]. The enhanced sensitivity of the latter instrument resulted in a 136% increase [~ 1.4 -fold increase] in coverage of the human plasma metabolome in these experiments [274 vs. 203 reference-matched IDs]. These identified metabolites have match scores ranging from 2.2 to 1.1. Importantly, in all the comparative experiments performed on the ZenoTOF 7600 and 8600 systems, the latter instrument consistently generated more features, and more MS/MS dependent scans were triggered (**Figure 1**). This observation is likely due to the enhanced sensitivity in both the MS1 and MS/MS modes of the ZenoTOF 8600 system. The higher MS1 sensitivity logically leads to more candidate ions for the dependent product ion scan, and the enhanced MS/MS sensitivity will impart better fragment intensities and a better overall quality of MS/MS spectrum.

In summary, the ZenoTOF 8600 system is shown to generate superior untargeted metabolomics results compared to the ZenoTOF 7600 system. Although it is challenging to find true equivalent data sets in the metabolomics literature, a study using NIST 1950 plasma analyzed in the positive ion mode, with analytes resolved using RP chromatography, reported 72 distinct metabolites were found using MS-DIAL software for data processing [7], which is less than found here using either instrument. As mentioned above, it is problematic to do a numbers comparison, especially in metabolomics, where there is little standardization and many different software platforms by which the data are analyzed. In the data presented here, the strongest conclusions should be made from the direct comparison of the two instruments' performances, analyzing the same sample under the same experimental conditions and data processing parameters.

Conclusions

- In comparative studies, the ZenoTOF 8600 system is ~ 10 -fold more sensitive than the ZenoTOF 7600 system at the TOFMS and TOF MS/MS levels of analysis
- The SCIEX ZenoTOF 8600 system identified ~ 1.4 -fold more metabolites than the ZenoTOF 7600 system in human plasma extract [275 vs. 203]
- Using the same analytical method, the improved sensitivity of the ZenoTOF 8600 system generated an increased number of features

and dependent product ion scans compared to the ZenoTOF 7600 system [6064 vs. 3932 for 1 μ L human plasma equivalents]

References

1. Broeckling, C. D., Beger, R. D., Cheng, L. L., Cumeras, R., Cuthbertson, D. J., Dasari, S., Davis, W. C., Dunn, W. B., Evans, A. M., Fernández-Ochoa, A., Gika, H., Goodacre, R., Goodman, K. D., Gouveia, G. J., Hsu, P. C., Kirwan, J. A., Kodra, D., Kuligowski, J., Lan, R. S., Monge, M. E., ... Mosley, J. D. [2023]. Current Practices in LC-MS Untargeted Metabolomics: A Scoping Review on the Use of Pooled Quality Control Samples. *Analytical chemistry*, 95(51), 18645–18654. <https://doi.org/10.1021/acs.analchem.3c02924>
2. Defosse, E., Bourquin, J., von Reuss, S., Rasmann, S., & Glauser, G. [2023]. Eight key rules for successful data-dependent acquisition in mass spectrometry-based metabolomics. *Mass spectrometry reviews*, 42(1), 131–143. <https://doi.org/10.1002/mas.21715>
3. Rupasri Mandal, Jiamin Zheng, Lun Zhang, Eponine Oler, Marcia A. LeVatte, Mark Berjanskii, Matthias Lipfert, Jun Han, Christoph H. Borchers, and David S. Wishart. *Analytical Chemistry* **2025** 97 (1), 667–675. <https://doi.org/10.1021/acs.analchem.4c05018>
4. Takeda, H., Matsuzawa, Y., Takeuchi, M., Takahashi, M., Nishida, K., Harayama, T., Todoroki, Y., Shimizu, K., Sakamoto, N., Oka, T., Maekawa, M., Chung, M. H., Kurizaki, Y., Kiuchi, S., Tokiyoshi, K., Buyantogtokh, B., Kurata, M., Kvasnička, A., Takeda, U., Uchino, H., ... Tsugawa, H. [2024]. MS-DIAL 5 multimodal mass spectrometry data mining unveils lipidome complexities. *Nature communications*, 15(1), 9903. <https://doi.org/10.1038/s41467-024-54137-w>
5. Ozbalci, C., Baker, P. R., & Sayers, R. MS-DIAL software parameters for processing untargeted metabolomics data acquired on the ZenoTOF 7600 system. [MKT-30320-A SCIEX-Technical-Note metabolomics ZenoTOF DDA parameters-FINAL 1-11-2024.pdf](#)
6. Tsugawa, H., Cajka, T., Kind, T., Ma, Y., Higgins, B., Ikeda, K., Kanazawa, M., VanderGheynst, J., Fiehn, O., & Arita, M. [2015]. MS-DIAL: data-independent MS/MS deconvolution for comprehensive metabolome analysis. *Nature methods*, 12(6), 523–526. <https://doi.org/10.1038/nmeth.3393>
7. Barbier Saint Hilaire, P., Rousseau, K., Seyer, A., Dechaumet, S., Damont, A., Junot, C., & Fenaille, F. [2020]. Comparative Evaluation of Data Dependent and Data Independent

MKT-34922-A

Acquisition Workflows Implemented on an Orbitrap Fusion
for Untargeted Metabolomics. *Metabolites*, 10(4), 158.

<https://doi.org/10.3390/metabo10040158>

Untargeted data-dependent acquisition (DDA) metabolomics analysis using the ZenoTOF 7600 system

Source, compound and experimental parameters to perform DDA experiments using collision-induced dissociation

Paul RS Baker and Robert Proos
SCIEX, USA

This technical note details instrument settings and best practices for untargeted metabolomics using the ZenoTOF 7600 system. The ion flow path of this instrument is different from any other QTOF instrument available. Therefore, untargeted metabolomics method development requires unique instrument settings to achieve success. Here, appropriate parameter settings are provided and these settings are explained in the context of how the device functions and is designed. Understanding how the hardware is used to generate untargeted metabolomics data can dramatically improve data quality and increase the number of metabolites detected and identified during data processing.

Ideally, untargeted metabolomics experiments are designed to detect and potentially quantify all biomolecules in a sample. These data represent a snapshot of the metabolic state of the organism studied. A typical results file for DDA metabolomics analysis is shown in Figure 1. In practice, the quality of the data and the extent of metabolite coverage depend on the parameter settings used on the instrument and the experimental method. For targeted analysis experiments, in which internal and primary reference standards are typically available, instrument parameter

settings can be optimized empirically. In contrast, generalized parameter settings are required for global analysis to maximize coverage of the diverse compounds present in the biome.

Key features of untargeted metabolomics analysis using the ZenoTOF 7600 system

- The unique instrument design that enables electron-activated dissociation (EAD) and the use of the Zeno trap requires specific instrument parameter settings for untargeted metabolomics
- Parameter setting suggestions and their detailed explanations are provided to optimize metabolomic coverage and spectral quality
- The untargeted metabolomics method described here leverages the speed and sensitivity of the ZenoTOF 7600 system to maximize throughput while maintaining data quality


Figure 1. TIC from a 15 min untargeted metabolomics experiment using the IDA scan mode. The speed and sensitivity of the ZenoTOF 7600 system enable high-throughput analysis while maintaining high MS/MS spectral quality.

Methods

Sample preparation: NIST SRM 1950 plasma was extracted using 4 volumes of methanol. After centrifugation to separate the precipitated proteins, the supernatant was analyzed directly.

Chromatography: Samples were analyzed using an ExionLC system with a Kinetex F5 column (2.1 × 150 mm, 1.7 μm, Phenomenex). A simple linear gradient from 0 to 95% B was used with standard reversed phase mobile phases (A = 0.1% formic acid in water and B = 0.1% formic acid in acetonitrile) with a flow rate of 200 μL/min. A 1 μL injection volume was used and the column temperature was maintained at 30°C throughout the analysis. The total run time was 15 min.

Mass spectrometry: The extracted sample was analyzed on the ZenoTOF 7600 system equipped with the OptiFlow Turbo V ion source. Data were collected using DDA methods designed to vary from top 20 to top 100 candidate ions, with dynamic background subtraction (DBS), and exclusion criteria of 6s after 3 occurrences, both activated and deactivated. The TOF MS accumulation time was set to 100 ms, and a collision energy (CE) of 30 V was used. A 5 ms accumulation time was used for TOF MS/MS experiments.

Data processing: All data were analyzed using the Analytics module in SCIEX OS software. The MS/MS spectra were compared to the library spectra using version 2.0 of the SCIEX All-in-One and SCIEX Accurate Mass Metabolite spectral libraries and the NIST 17 spectral library.

Instrument calibration

The first step in developing an untargeted metabolomics method is to ensure that the ZenoTOF 7600 system is operating at peak performance. Instrument performance depends on several factors but can be distilled down to 2 requirements: 1) the instrument must be correctly calibrated, and 2) it must be free from contamination. Below are steps that should be taken before any analytical experiments are performed.

1. Open the MS Tuning tab
2. Every 1-2 weeks, run EI Background Reduction to help clean the EAD cell, even if the EAD fragmentation function has not been used. The precursor ions travel through the EAD cell to get to the CID portion of Q2 and, subsequently, to the accelerator, TOF and detector. If the EAD cell is dirty, it may affect the sensitivity and resolution of the system in CID mode.
3. Every day, run MS Quick Check for the polarity that will be used. When using the TwinSpray ion source, there is no need to change any source parameters from the default

settings. However, if using the OptiFlow Turbo V ion source, the source temperature should be changed to 500°C.

Important note: Changing this value in the box does not change the temperature setting in real time. To implement this change, proceed to Step 2-Achieve Stable Signal, change the value from 200°C to 500°C, then click in another of the setting boxes. Press Stop, click on Step 1, then click Next to proceed back to Step 2. The value will then be appropriately set to 500°C. The system might take a couple of minutes to equilibrate to the new temperature.

4. During the MS Quick Check, observe the intensity of the MS/MS product ion when the Zeno trap is deactivated (m/z 520 in the negative ion mode or m/z 494 in the positive ion mode). Ideally, this signal must be above 1e5 counts per second (cps). If the Zeno trap is activated during this step, the signal will be variable because it depends on the MS-level sensitivity as a function of different ion transmission control (ITC) values and therefore will not be consistent over time. Save the tuning reports to a dedicated folder so that the signal of this third ion can be tracked as a measure of overall instrument performance.
5. If the mass accuracy calibration fails (for example, is >2 ppm) but is outside the tolerance range (for example, within 4 ppm), the system is likely suitable to use. The in-run calibration employs a mass tolerance range of 5 ppm.
6. If the sensitivity or resolution fails during calibration, perform the ADC Initialization and run EI Background Reduction, followed by TOF Tuning. These functions all appear in the MS Tuning tab. Below are the steps to take in response to failures during calibration. Note: If there is a need to proceed to steps (e) and (f), these steps can be performed in any order, depending on timing constraints and familiarity with venting and cleaning the QJet. A flow chart is presented in Figure 2 that summarizes these actions.

- ADC Initialization.** ADC initialization re-establishes proper communication between the detector and the computer. Once reset, rerun MS Quick Check. If it passes, then start sample analysis. If it fails, proceed to step (b).
- EI Background Reduction.** The Background Verification Level process will almost always fail the first time. Press next and perform the Background Reduction Scan. The signal should decrease significantly over the first few minutes. After 5-10 minutes, press Next to proceed to the final check. Save the tuning report, including the optional steps. Re-check MS Quick Check. If it passes, then start sample analysis. If it fails, proceed to step (c).
- TOF Tuning.** Select TOF tuning and allow the process to proceed. Re-run MS Quick Check. If it passes, then start sample analysis. If it fails, proceed to step (d).
- Zeno Calibration.** After TOF Tuning, the Zeno trap should be re-calibrated by clicking Zeno Calibration. Once this test is complete, run MS Quick Check again. If the instrument fails this last check, please clean the QJet.
- QJet Cleaning.** If the instrument sensitivity test still fails, vent the system and clean the curtain plate, the orifice plate and the QJet. SCIEX has an official QJet cleaning procedure and detergent, bags and brushes. Do not perform any cleaning or maintenance of the instrument beyond the QJet. Once vented, the QJet is cleaned and the instrument is reassembled, the system will need time (2 - 8 min) to re-achieve the suitable operating vacuum pressure in the TOF chamber. After restarting and allowing the pressure to stabilize, run ADC Initialization and run MS Quick Check.
- Detector Optimization.** If the instrument fails the sensitivity test after re-running MS Quick Check, then run Detector Optimization. This test varies the voltage on the detector to ensure it is held at the correct voltage. As detectors age, the voltage required to maintain sensitivity increases. There is a voltage limit and the detector must be replaced when the limit is reached. In this situation, the analytical quadrupoles or other ion path components might also be contaminated and deep cleaning will need to be performed by a service engineer.

ZenoTOF 7600 system hardware configuration

The functions and mechanisms of several features on the ZenoTOF 7600 system will be explained here to better


Figure 2. Flowchart depicting the correct order of MS tuning procedures to follow if the initial MS Quick Check experiment fails at any stage of the run. At each stage of the process, re-run MS Quick Check to see if the issue is resolved. If it is, proceed with sample analysis. If not, continue to troubleshoot by performing the next tuning procedure in the list and re-testing via MS Quick Check.



Figure 3. Assembly diagram of the ZenoTOF 7600 system. The integrated MS/MS assembly supports both CID and EAD fragmentation. Fragment ions can be trapped by the Zeno trap, which increases the instrument duty cycle to >90%.

understand how specific parameter settings are set to optimally run untargeted metabolomics in a DDA experiment using the information-dependent acquisition (IDA) scan mode. Figure 3 shows the assembly of the primary ion rail components, through

which the ions flow from left to right. As with all LC-ESI MS/MS instruments, ions are generated in the source, the ion beam is created and refined in the interface region, which comprises the

Figure 4. Diagram of the ionization process in the electrospray electrode. Note that both positive and negative ions are produced during the process. In this example, a positive voltage is applied, causing negative ions to migrate to the edges of the electrode and the positive ions to migrate to the center. The spray is not homogeneous, so it is important to consider the position of the electrode and the drying gas (GS2) when tuning the source parameters to maximize signal intensity.

QJet and Q0 region, and ions are selected for product ion analysis in Q1. In the ZenoTOF 7600 system, ions leave Q1 and enter the EAD cell, which can trap and fragment ions or transfer precursor ions to Q2 for CID-based fragmentation. In either case, fragment ions enter the Zeno trap, where they are collected and sent to the TOF accelerator in a mass-dependent manner. EAD-based fragmentation will not be discussed in the context of untargeted metabolomics in this technical note.

The Zeno trap functions by trapping ions generated by EAD- or CID-based fragmentation. The Zeno trap is not used during TOF MS scans because the ion current is already attenuated via ion transmission control (ITC) and the likelihood of detector saturation is high. Once the ions are trapped ($t < 1$ ms), they are sequentially scanned from the trap to the TOF accelerator in a mass-dependent manner. High molecular weight ions are scanned out first and are followed by lower molecular weight ions, such that they all arrive at the TOF accelerator at the same time. This eliminates the mass intensity bias that results in other instruments only taking a portion of the fragment ions. Once the ions are accelerated into the TOF, the Zeno trap fills and repeats the cycle. This process improves the instrument duty cycle and ~90-95% of all fragment ions hit the detector.

RUO-MKT-02-15367-A

A key parameter to consider regarding the Zeno trap is the Zeno threshold, which is a measure of the precursor ion intensity. If the precursor ion intensity is below this value, the Zeno trap is used and fragment ions from that precursor ion are trapped and sent to the TOF accelerator. If the precursor ion intensity is above the threshold value, the Zeno trap is not used. In this case, the fragment ions that hit the detector without the Zeno trap activated are scaled according to the set Zeno calibration values. Therefore, the visualized data are corrected to what the intensity would have been if the Zeno trap had been used.

Parameter settings for DDA experiments

The parameter settings for untargeted metabolomics on the ZenoTOF 7600 system will differ from those used on other SCIEX instruments. The primary reason for this is that the ion path is different from other instruments due to the incorporation of the EAD cell and the Zeno trap. For example, the EAD cell and the Zeno trap occupy a portion of Q2 and therefore make the collision cell shorter. Consequently, the collision energy (CE) parameter setting should be higher than that used on other SCIEX instruments to compensate for the shorter collision cell dimensions.

Figure 5 shows a screenshot of the parameter settings that were optimized for untargeted metabolomics on the ZenoTOF 7600 system in the positive ion mode. For the experiments performed in the negative ion mode, the absolute value for each setting is the same and further discussion points are provided below. In the sections below, many of the critical parameters are discussed, starting with source parameters.

Spray voltage. Typically, the ESI voltage is set at or near its maximum. It is essential to recognize that during ionization, both positive and negative ions are produced (conservation of charge). If a positive voltage is applied (for example, in positive ion mode), positive ions will migrate toward the center of the effluent stream and negative ions will migrate toward the walls of the electrode (Figure 4). When the spray is created, the plume is not homogeneous. Negative ions will be at the outside edges of the plume and positive ions will be at its center. This segregation of the ions makes it essential to always adjust the position of the probe on the ZenoTOF 7600 system. The position of the OptiFlow Turbo V ion source does not need to be adjusted, however, because its geometry is already optimized. Other key source parameters include the nebulizing gas (GS1), which affects droplet size, and the drying gas (GS2), which evaporates

The screenshot displays the configuration interface for a mass spectrometer experiment. It is organized into several sections:

- Experiment:** Includes a dropdown menu set to 'IDA'. Parameters include Polarity (Positive, V) and Spray voltage (5500, V).
- TOF MS:** Parameters include TOF start mass (70, Da), TOF stop mass (1000, Da), Accumulation time (0.1, s), Declustering potential (40, V), DP spread (0, V), Collision energy (10, V), and CE spread (0, V).
- IDA Criteria:** Set to 'Small molecule'. Parameters include Maximum candidate ions (40), Intensity threshold exceeds (100, cps), and checkboxes for Dynamic background subtraction (checked) and Exclude former candidate ions (unchecked). It also includes 'For' and 'After' settings.
- TOF MSMS:** Parameters include Fragmentation mode (CID), Precursor ion (100, Da), Q1 resolution (Unit), TOF start mass (50, Da), TOF stop mass (1000, Da), Accumulation time (0.005, s), Declustering potential (40, V), DP spread (0, V), Collision energy (30, V), CE spread (0, V), and Zero pulsing (checked).

Figure 5. User interface screencap showing the specific recommended settings for global untargeted metabolomics experiments in the positive ion mode. The experimental set up includes parameter settings for the source, TOF MS, TOF MS/MS and important considerations for IDA criteria. These considerations include whether to use DBS and/or exclusion criteria. For the negative ion mode, the absolute value of each setting is the same.

the solvent to generate gas-phase ions. Subtle changes in either or both parameters can affect instrument sensitivity. In samples with heavy matrix contamination, as is common in metabolomics experiments, one can slightly lower the ESI voltage to reduce the ionizing species that come from the matrix. With a high voltage setting, highly ionizable species in the matrix may crowd out the lower ionizable species and suppress the analyte signal.

Declustering potential. The declustering potential (DP) is the voltage potential between the orifice and the front of the QJet. It acts as the motive force to draw ions into the instrument and helps disrupt any solvent clusters surrounding the ions. If the DP setting is too high, it can lead to source fragmentation. For both the TOF MS and the dependent TOF MS/MS portions of the DDA experiment, the DP should be set to 40 V in untargeted metabolomics experiments. Ion intensity does not change significantly with changes to DP but a higher setting (~60 V) might be preferred for larger compounds. In contrast, small molecules such as amino acids, are better analyzed at a lower DP (~20 V). The 40 V setting is a middle ground and gives the best results for a diverse set of analytes.

The parameters below affect the ion flow through and fragment generation by the ZenoTOF 7600 system. Differences between TOF MS and TOF MS/MS experiments are explained.

Collision energy. The collision energy (CE) setting affects the acceleration of ions through Q2. During CID, CE is responsible for imparting enough energy to the ions to fragment when they collide with nitrogen in the collision cell. However, CE also serves as a motive force for ions in non-MS/MS experiments. For the TOF MS survey scan in the DDA method, note that the CE is set to 10 V, even though no fragmentation is desired. As ions move along the ion path, they cool down. This is notable at the entrance to Q2, where nitrogen gas is present and

Table 1. Effects of varying IDA criteria on the features identified, MS/MS events triggered, library hits, and cycle time.

Experiment	DBS	Exclusion list	Peaks (features)	Peaks w/ MS/MS spectra	Formula Finder hits	Library hits with >70% purity score	% MS/MS with library match	Average cycle time (s)
Blank	Yes	No	3781	2319	2499	71	3.06	255
POS IDA Top 30	Yes	No	6811	4435	4200	216	4.87	321
POS IDA Top 40	Yes	No	6065	4184	3771	204	4.88	327
POS IDA Top 50	Yes	No	5576	3978	3524	221	5.56	376
POS IDA Top 60	Yes	No	4879	3524	3011	198	5.62	392
POS IDA Top 70	Yes	No	4404	3294	2724	208	6.31	420
POS IDA Top 80	Yes	No	4025	3060	2528	188	6.14	439
POS IDA Top 90	Yes	No	3825	2918	2429	190	6.51	446
POS IDA Top 100	Yes	No	3641	2838	2314	192	6.77	464
POS IDA Top 30	Yes	Yes	5363	3614	3349	210	5.81	308
POS IDA Top 40	Yes	Yes	4924	3510	3070	219	6.24	336
POS IDA Top 50	Yes	Yes	4644	3398	2902	216	6.36	357
POS IDA Top 60	Yes	Yes	4346	3264	2734	230	7.05	384
POS IDA Top 70	Yes	Yes	4127	3134	2583	209	6.67	391
POS IDA Top 80	Yes	Yes	4028	3114	2566	197	6.33	423
POS IDA Top 90	Yes	Yes	3804	2937	2423	199	6.78	436
POS IDA Top 100	Yes	Yes	3555	2783	2254	197	7.08	451
POS IDA Top 30	No	No	5225	1808	3254	216	11.95	624
POS IDA Top 40	No	No	4718	2139	2958	216	10.10	777
POS IDA Top 50	No	No	4376	2319	2729	227	9.79	928
POS IDA Top 60	No	No	4073	2485	2591	196	7.89	1082
POS IDA Top 70	No	No	3820	2554	2435	203	7.95	1235
POS IDA Top 80	No	No	3641	2616	2326	203	7.76	1380
POS IDA Top 90	No	No	3336	2584	2123	212	8.20	1515
POS IDA Top 100	No	No	3278	2685	2122	172	6.41	1651
POS IDA Top 30	No	Yes	5120	3291	3220	264	8.02	629
POS IDA Top 40	No	Yes	4765	3341	2978	246	7.36	773
POS IDA Top 50	No	Yes	4366	3306	2793	244	7.38	915
POS IDA Top 60	No	Yes	4089	3219	2596	218	6.77	1051
POS IDA Top 70	No	Yes	3816	3095	2454	223	7.21	1181
POS IDA Top 80	No	Yes	3573	3042	2241	217	7.13	1320
POS IDA Top 90	No	Yes	3406	2915	2172	194	6.66	1449
POS IDA Top 100	No	Yes	3189	2776	2063	198	7.13	1587

RUO-MKT-02-15367-A

collisional cooling is significant. By adding a small amount of CE voltage, the ions are pushed through Q2 without significant fragmentation and arrive at the TOF accelerator.

The CE voltage is set to 30 V for CID-based MS/MS analysis. This value gives similar MS/MS spectra to the SCIEX metabolomics library collected using SCIEX TripleTOF systems (4600, 5600 and 6600 series). These instruments were run at 35 V with a collision spread of 15-20 V. Note that on the ZenoTOF 7600 system, the CE spread is set to 0, meaning that all the spectra will be acquired at the designated 30 V. The main reason for this, compared to the instrument settings using the TripleTOF systems, is that it takes approximately 15-20 ms to enable the CE spread function. The ZenoTOF 7600 system significantly improves the instrument duty cycle, sensitivity and detection of low m/z fragments. It was more advantageous to shorten the accumulation time and decrease the cycle time to improve coverage, than it was to use a collision energy spread.

Accumulation time. For global metabolomics experiments, accumulation time is one of the most important parameters to consider. In general, increasing accumulation time will improve the signal-to-noise ratio (S/N), which improves the quality of the MS and MS/MS spectra. Higher-quality spectra generally have increased library hits. However, if the accumulation time is too long, the cycle time of the DDA experiment is extended and fewer MS/MS spectra are acquired. For the MS TOF experiments, it was empirically determined that the best balance is achieved using 0.100 s (100 ms) for untargeted metabolomics experiments. At the MS/MS level, it was determined that 5 ms is the best accumulation time. Although the latter setting might seem too small, the ZenoTOF 7600 system was designed to be fast (133 CID MS/MS per second, 133 Hz). The 8- to 10-fold improvement in the instrument duty cycle from the Zeno trap compensates and enables high-quality MS/MS spectra at high-speed acquisition rates. In general, with this instrument, faster accumulation times are better and there do not appear to be significant gains in spectral quality at accumulation times >10 ms.

Zeno pulsing. The function of the Zeno trap is explained above. By clicking the Zeno pulsing box, the Zeno trap will be activated, provided that the precursor ion intensity is below the set Zeno threshold. For untargeted metabolomics experiments, it is recommended that the Zeno threshold is set to 20,000 cps.

IDA criteria. One of the more challenging aspects of any DDA experiment is selecting the optimal settings that control which precursor ions are selected for subsequent MS/MS analysis. In SCIEX DDA experiments, this group of parameters is controlled by the IDA criteria panel (Figure 4, blue box). To determine the

best IDA criteria parameter settings, a series of experiments was performed to vary the maximum number of candidate precursor ions per cycle ("maximum candidate ions"), whether DBS was used and whether exclusion criteria were enforced (Table 1). In all experiments, the TOF MS and TOF MS/MS instrument parameter settings were set, as shown in Figure 4.

The comprehensive dataset shown in Table 1 is separated into 4 sections, divided by the darker lines within the table. Each section shows 8 experiments that were performed in the positive ion mode, in which the maximum number of candidate ions was varied from the top 30 to 100 candidates. The 4 sections differ by whether DBS and/or exclusion criteria were applied. Each variable and its effects on the data are discussed below. For this type of experiment, no advanced IDA criteria parameters are used.

1. **Maximum candidate ions.** Ideally, DDA experiments are designed to maximize the coverage of unknown molecules while still maintaining high MS/MS spectral quality. Balancing coverage and spectral quality are most affected by the number of candidate ions chosen. Table 1 shows the numbers of peaks (features) detected, peaks detected that generated MS/MS spectra, Library and Formula Finder hits and the average cycle time as a function of candidate ions. Figure 6 shows a comparison of the MS/MS spectral quality for tyrosine acquired with 30, 70 and 100 candidate ions selected. From these data, the best number of candidate ions was determined to be between 30 and 50. As the number of candidate ions increased above 50, the number of features identified decreased.
2. **Dynamic background subtraction.** DBS is a function that improves the detection of precursor ions in a DDA experiment. The algorithm creates an extracted ion chromatogram (XIC) of the candidate ion over the preceding data points. It takes the first derivative of the curve and determines whether the candidate ion is at the apex of the first derivative. If it is at the apex, the MS/MS-dependent scan will be triggered. This algorithm is designed for small molecules and minimizes redundant MS/MS collection. Additionally, it reduces background contamination from triggering dependent MS/MS scans.
3. **Exclusion criteria.** Exclusion criteria are parameters that limit redundant MS/MS spectra from individual precursor ions. At least 1 MS/MS scan should be triggered near or at


Figure 4. Example MS/MS spectra from the DDA data obtained in Table 1. The MS/MS spectral quality is not affected by the use of DBS or exclusion criteria. However, DBS or exclusion criteria are useful to limit false increases in features that result from redundant peak analysis.

the apex of the peak. If a product ion scan is triggered from a candidate ion, that precursor ion will be ignored for the exclusion period, which should be entered as half of the peak width in seconds. Typically, this exclusion is triggered after 1 or 2 triggered MS/MS scans. This feature works well for large- and small-molecule applications.

The effects of DBS and exclusion criteria on global metabolomics results are detailed in Table 1. The use of DBS resulted in more features with MS/MS spectra, however, the number of library matches from the data did not increase correspondingly. As the MS/MS spectral quality was not impacted (Figure 6), this outcome might be attributed to the spectral libraries used not including the lower intensity analytes. The use of exclusion criteria without DBS resulted in slightly more library matches, as additional features included MS/MS data. When no filtering is used for precursor selection, the most intense ions will be continually selected for the MS/MS scans at

the expense of less intense ions, except at the highest number of candidate ions. In this case, the fewest peaks with MS/MS spectra and the fewest library matches are observed. In general, using either DBS or exclusion criteria in the acquisition method will increase both the number of peaks with MS/MS data and the number of library matches.

Conclusions

- Source parameters, including the probe position, drying and nebulizing gases and the electrospray voltage on the ZenoTOF 7600 system must be optimized for untargeted metabolomics experiments. The optimization of these parameters can be affected by the flow rate and the matrix composition of the sample.
- Due to the unique hardware configuration of the instrument, parameter settings, such as collision energy and

accumulation, must be adjusted to accommodate diverse metabolites

- The ZenoTOF 7600 system is ideal for untargeted metabolomics due to its speed, sensitivity, and ability to perform complementary EAD experiments

SWATH DIA and library driven compound identification revealed upregulation of flavonoids in *Ocimum basilicum* L. flowers

Identification of differentially regulated metabolites using SWATH DIA

Marialuca Maldini¹, Kranthi Chebrolu², Paola Montoro³, Luigi d'Aquino⁴ and Christie Hunter²
¹SCIEX, Italy, ²SCIEX, US, ³University of Salerno, Fisciano (SA), ⁴ENEA, Portici (NA)

Indoor plant growing is increasing in urban environments to deal with increasing soil loss and to shorten food chains. Light is among the most important environmental factors to be addressed in indoor cultivations, since it is the main driver for basic plant processes such as growth and development but also for directional movement, circadian rhythms, seed germination, etc. Different plant species are differently affected by light conditions (spectra, photon flux, light / dark ratio) and different phenological stages display different optima about light conditions. A better understanding of the effect of different light wavelengths on plant growth and development in indoor environment is, therefore, needed to better drive plant cultivation.¹ Metabolomics has served as a powerful tool to elucidate gene functions, enzyme activity and provide a detailed



Figure 1. Deconvoluted MS/MS spectrum matching with library spectrum. (Top) The blue colored spectrum is the deconvoluted MS/MS spectrum from the SWATH acquisition data which is overlaid with raw MS-MS spectrum (pink) for parent mass 464.0876 at 5.9 min. (Bottom) Grey colored spectrum represents MS/MS spectra of Hyperoside from NIST library. The fit and purity score of 99 indicate the confidence of the identification for hyperoside from the library match.

blueprint of overall plant health. In this study, basil (*Ocimum basilicum* L.) was grown from seedling to flowering stages in a new concept microcosm (EP 3 236 741 B1) in which the effects of different light spectrums on the growth and metabolome of different plant parts are investigated.^{2,3} SWATH DIA, a data independent acquisition workflow, was used to generate a comprehensive signature of the metabolome of each sample.

In the previously proposed metabolomic workflow, the untargeted identification step was avoided by reduction of total number of features to a smaller feature list that are significantly different across biological samples.² Later, a complete compound annotation of a selected feature list is carried out through targeted identification using library spectral match.

Here, an untargeted data processing pipeline was used where all unknown features present in data that match with library spectrum are identified prior to performing multivariate statistical analysis. High confident identification of features utilizes MS¹, MS/MS and deconvoluted spectrum for library spectral matching (Figure 1).

In the current workflow, metabolites identified with high confidence in SCIEX OS software are plotted in MarkerView software to visualize statistical differences through PCA-DA and volcano plots. A great similarity in identified features in the proposed workflow (SCIEX OS to MarkerView software) to the conventional workflow (MarkerView to SCIEX OS software) illustrates a robust data processing workflow.²

Key features of data acquisition and processing workflow

- Metabolic features from high-resolution, accurate mass MS and MS/MS data acquired using SWATH DIA were identified by leveraging spectral libraries in SCIEX OS software.
 - Natural products HR-MS/MS and NIST libraries were used to identify key metabolites with high confidence
 - Identified metabolite panel were passed through to MarkerView software for multivariate statistical analysis including PCA and T-test
 - Graphical displays allow easy determination of differentiating metabolites

Methods

Sample preparation: The plants were cultivated for approximately 60 days under either white (W) or blue/red (B/R) light. The top leaves (apical, Ap L), middle leaves (median, Me L) and inflorescences (flower) were harvested and extracted with a solution of 1:1, ethanol/water. A Sep-Pak C18 was used to extract chlorophyll.²

Chromatography: An ExionLC system equipped with an HSS T3 column (100×1.0mm, particle size 1.8µm, Waters) was used for metabolite separation at a flow rate of 0.45 mL/min.²

Mass spectrometry: Data were acquired in positive and negative modes using a TripleTOF 6600 system using SWATH DIA. Details of the acquisition strategy are summarized in this technical note.²

Data processing: SWATH DIA data were processed in SCIEX OS software 1.6 using the Natural Products HR-MS/MS, Accurate Mass High Resolution and NIST spectral libraries. Feature identification step was conducted using a non-targeted workflow through Analytics section in SCIEX OS software. Sample data (.wiff files) were uploaded in Analytics, and a new processing method was created using a representative quality control sample. In contrast to the first technical note, here the metabolites were first identified, and then identified metabolites were used in the statistical analysis to find differentiating features. For this workflow type, non-targeted screening was selected, and components selection can be left blank (Figure 2). In the Library Search section, the Smart Confirmation Search is selected for the library search algorithm. Purity was selected for how the results would be sorted. Then, all the libraries used in compound identification were selected.

The mass tolerance for both MS and MS/MS was set to 0.02 Da (Figure 2). Finally, the polarity box was checked which enables the software to recognize the mode in which data is acquired. Under flagging rules, Mass Error (ppm), Library hit score and Formula finder score boxes were selected, and appropriate weighting percent are input in the boxes while the percent difference isotope ratio is optional.

Figure 2. Building a processing method for non-targeted metabolomics workflow for SWATH DIA data. The left panel depicts the selection of workflow type and selection of a representative sample for building a processing method. In the middle panel, the components section can be left empty or only first row should be filled with a known endogenous compound. The right panel shows the parameter settings, where smart confirmation search is selected. The results are sorted by purity of the spectrum. All libraries that are relevant for the sample types are selected, and precursor and fragment mass tolerances are selected as 0.02 Da. Finally, by checking the polarity box, software recognizes the polarity mode in which SWATH DIA data is acquired.

RUO-MKT-02-14667-A

In Formula Finder section, since these are plant samples, only naturally occurring compounds option was chosen because the data acquired is on plant samples. In Max element section, S (sulfur) was de-selected because sulfur compounds are uncommon in plants. Under non-targeted peaks, the retention time range was set to 2 min to 26 min, peak detection sensitivity was set to two steps away from “exhaustive” to avoid peak detection at noise level.

The result file consists of features that are identified through MS/MS spectral matching with library at various confidence levels, and features with only MS spectrum alone. The features that were identified with high confidence were exported from SCIEX OS software. Further, these identified features were copied into Excel where a vlookup function was used to consolidate the data into a single list of unique compounds with peak areas from each sample associated.

After consolidation of the data, this data was imported as text into MarkerView software for multivariate statistical analysis and generation of 2D graphical plots.

Comprehensive data collection using SWATH DIA

The study of the metabolome is key to understanding the basic biological mechanisms that regulate plant growth and stress responses, however the metabolome is highly complex and comprehensive interrogation is difficult using traditional data

dependent approaches. The SWATH DIA workflow was used here because it can acquire high-resolution MS and MS/MS data on all detectable precursors in the sample.² The method is sensitive and reproducible and provides a comprehensive view of the metabolome using a relatively short 30-minute run. Both the high-resolution MS and MS/MS spectra are used for metabolite identification by library search (Figure 1), and quantitative peak areas can be obtained from both the MS and MS/MS data.

Identifying metabolites using library matching

Feature identification through library matching was performed using Analytics in SCIEX OS software, using the NIST 2017, Natural product and Metabolite_HR-MS-MS libraries. Overall, 508 metabolites were identified across the dataset with high confidence.

Hyperoside (Quercetin 3-galactoside/Hyperin), m/z of 463.0891 eluted at chromatographic retention time of 5.87 min, was one of several features that was identified with very high confidence using NIST library MS-MS spectrum that showed significant difference apical leaves to flowers (Figure 3).

Sample Name	Component Name	Library Hit	Component Group Name	Actual Conc.	Area	Retention Time	Abund. Change	Precursor Mass	Mass Error	Library Conf.	Form. Conf.	Mass Error (ppm)
1.DM_SWATH_H02 (1)	311.0412 / 5.68	Gallic acid (NIST) [Smart Confirmation]	281.0225 / 5.68	NA	6.41763	5.67	NA-H+	211.041	0.0	✓	✓	16.6
1.DM_SWATH_H02 (2)	179.0211 / 5.64	Gallic acid (NIST) [Smart Confirmation]	308.1181 / 5.72	NA	6.22864	5.67	NA-H+	179.021	0.0	✓	✓	16.6
1.DM_SWATH_H02 (3)	161.0241 / 5.60	Gallic acid (NIST) [Smart Confirmation]	308.1181 / 5.72	NA	6.22863	5.67	NA-HCO-H+	161.024	0.0	✓	✓	16.6
1.DM_SWATH_H02 (4)	473.0718 / 5.65	鞣花酸 Gallic acid (NIST) [Smart Confirmation]	478.0718 / 5.65	NA	5.15845	5.67	NA-H+	473.072	0.0	✓	✓	16.6
1.DM_SWATH_H02 (5)	167.0219 / 5.71	3-Hydroxy-5-methylsalicylic acid (NIST) [Smart Confirmation]	167.0219 / 5.71	NA	4.43867	5.70	NA-H+	167.022	0.0	✓	✓	16.6
1.DM_SWATH_H02 (6)	168.0216 / 5.77	3-Hydroxybenzoic acid (NIST) [Smart Confirmation]	168.0216 / 5.77	NA	4.77564	5.76	NA-H+	168.022	0.0	✓	✓	16.6
1.DM_SWATH_H02 (7)	463.0891 / 5.87	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	463.0891 / 5.87	NA	6.70567	5.87	NA-H+	463.089	0.0	✓	✓	16.6
1.DM_SWATH_H02 (8)	461.0772 / 5.80	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	461.0772 / 5.80	NA	6.70567	5.87	NA-H+	461.077	0.0	✓	✓	16.6
1.DM_SWATH_H02 (9)	441.0574 / 5.97	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	441.0574 / 5.97	NA	6.70567	5.87	NA-HCO-H+	441.057	0.0	✓	✓	16.6
1.DM_SWATH_H02 (10)	296.0264 / 6.08	(+)-Cyanocinnamic acid (NIST) [Smart Confirmation]	193.2029 / 6.19	NA	2.07764	6.08	NA-H+	296.026	0.0	✓	✓	16.6
1.DM_SWATH_H02 (11)	171.0263 / 6.09	(+)-Cyanocinnamic acid (NIST) [Smart Confirmation]	171.0263 / 6.09	NA	2.07764	6.08	NA-H+	171.026	0.0	✓	✓	16.6
1.DM_SWATH_H02 (12)	377.0247 / 6.12	(+)-Cyanocinnamic acid (NIST) [Smart Confirmation]	377.0247 / 6.12	NA	4.48145	6.12	NA-H+	377.025	0.0	✓	✓	16.6
1.DM_SWATH_H02 (13)	518.1067 / 6.13	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	518.1067 / 6.13	NA	3.22864	6.13	NA-H+	518.106	0.0	✓	✓	16.6
1.DM_SWATH_H02 (14)	543.1016 / 6.22	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	543.1016 / 6.22	NA	4.33364	6.22	NA-HCO-H+	543.102	0.0	✓	✓	16.6
1.DM_SWATH_H02 (15)	555.1012 / 6.12	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	555.1012 / 6.12	NA	1.37664	6.12	NA-H+	555.101	0.0	✓	✓	16.6
1.DM_SWATH_H02 (16)	395.0214 / 6.11	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	395.0214 / 6.11	NA	2.22764	6.11	NA-HCO-H+	395.022	0.0	✓	✓	16.6
1.DM_SWATH_H02 (17)	168.0424 / 6.49	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	168.0424 / 6.49	NA	2.08864	6.50	NA-H+	168.043	0.0	✓	✓	16.6
1.DM_SWATH_H02 (18)	451.0288 / 6.38	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	451.0288 / 6.38	NA	1.22264	6.37	NA-H+	451.028	0.0	✓	✓	16.6
1.DM_SWATH_H02 (19)	443.0271 / 6.40	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	443.0271 / 6.40	NA	2.08864	6.42	NA-HCO-H+	443.027	0.0	✓	✓	16.6
1.DM_SWATH_H02 (20)	160.0212 / 6.47	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	176.0217 / 6.57	NA	4.11163	6.48	NA-HCO-H+	160.021	0.0	✓	✓	16.6
1.DM_SWATH_H02 (21)	464.1119 / 6.40	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	464.1119 / 6.40	NA	1.16664	6.40	NA-HCO-H+	464.112	0.0	✓	✓	16.6
1.DM_SWATH_H02 (22)	338.0278 / 6.43	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	338.0278 / 6.43	NA	2.22264	6.44	NA-H+	338.028	0.0	✓	✓	16.6
1.DM_SWATH_H02 (23)	161.0212 / 6.49	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	161.0212 / 6.49	NA	6.70264	6.49	NA-H+	161.021	0.0	✓	✓	16.6
1.DM_SWATH_H02 (24)	208.0219 / 6.49	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	187.0219 / 6.51	NA	1.15164	6.49	NA-HCO-H+	208.022	0.0	✓	✓	16.6
1.DM_SWATH_H02 (25)	302.0247 / 6.43	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	302.0247 / 6.43	NA	1.33864	6.43	NA-H+	302.025	0.0	✓	✓	16.6
1.DM_SWATH_H02 (26)	171.0268 / 7.17	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	171.0268 / 7.17	NA	1.02164	7.19	NA-H+	171.027	0.0	✓	✓	16.6
1.DM_SWATH_H02 (27)	160.0213 / 7.18	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	160.0213 / 7.18	NA	1.43164	7.18	NA-H+	160.021	0.0	✓	✓	16.6
1.DM_SWATH_H02 (28)	302.0272 / 7.44	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	302.0272 / 7.44	NA	1.12864	7.44	NA-H+	302.027	0.0	✓	✓	16.6
1.DM_SWATH_H02 (29)	549.1424 / 7.44	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	549.1424 / 7.44	NA	4.72864	7.44	NA-HCO-H+	549.143	0.0	✓	✓	16.6
1.DM_SWATH_H02 (30)	302.0274 / 8.44	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	302.0274 / 8.44	NA	2.44864	8.44	NA-H+	302.028	0.0	✓	✓	16.6
1.DM_SWATH_H02 (31)	549.1421 / 8.44	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	549.1421 / 8.44	NA	4.72864	8.44	NA-HCO-H+	549.142	0.0	✓	✓	16.6
1.DM_SWATH_H02 (32)	267.0249 / 8.64	鞣花糖苷 Hyperoside (NIST) [Smart Confirmation]	267.0249 / 8.64	NA	6.15264	8.64	NA-H+	267.025	0.0	✓	✓	16.6

Figure 3. Results from library search in Analytics module of SCIEX OS software. Result file consists of sample name, component name (m/z and retention time information), library hit, compound ID, the data used to search the feature and library confidence and formula confidence. SCIEX OS software provides flexibility of sorting the compounds with high library confidence and formula confidence.

RUO-MKT-02-14667-A

Statistical analysis of identified features through MarkerView software

The peak areas of metabolites identified with high confidence from the library search performed in SCIEX OS software were formatted in excel and imported into MarkerView software for multivariate statistical analysis.

Figure 4 (left) displays the Scores plot (left) for Partial Least Squares Discriminant Analysis (PLS-DA) and Loadings plot (right) for all the biological replicates studied to understand compounds that significantly differentiate sample groups and replicates. Features with higher D1 and D2 loading scores are more differentiating. On a closer examination, flower samples exposed to white or blue/red light have (Figure 4, left) clustered together in second (top right) and fourth (bottom right) quadrants and clearly separated from the leaf samples exposed to white of blue/red light.

All the metabolites that have contributed to sample variation are displayed on the Loadings plot (Figure 4, right) by compound name, library from which it was identified, and the identification algorithm. For instance, hyperoside had higher D1 and D2 scores making it one of the key compounds differentiating flowers exposed to different wavelengths as opposed to leaves.


Figure 4. PCA analysis of identified metabolites from SWATH DIA data. (Left) Scores plot. PLS-DA showed the distinguishable segregation of flower (F) and apical leaf (AL) samples exposed to white and blue/red right wavelengths. (Right) Loadings plot of supervised multivariate statistical analysis (dimension reduction), shows all the compounds identified with high confidence that were used in performing PLS-DA. From the PLS-DA plots, the compounds that are highlighted in red color seems to have significantly impacted the groups that are observed in the score plots.

Several other compounds, luteolin 7- glucoside, betulinic acid, salvianolic acid and cichoric acid seem be more important compounds that are separating flower and leaf samples apart.

T-test analysis

Next, the pairwise *t*-test was performed to find specific differences between individual groups of interest. Figure 5 (left) shows the volcano plot with log fold change on x axis and *p*-value on y axis constructed from a pairwise *t*-test for 2 groups from batch 2, specifically the apical leaf samples grown under white light (AP L_WL) and the flower samples grown under white light (F_B/R). Compounds with a low *p*-value and high fold change are ones that differ in abundance between the groups with high confidence. Hyperoside and salvianolic acid with *p*-value ≤ 0.05 and fold change ≥ 3 , are compounds of interest for future experiments. Furthermore, about 85% of all the compounds in the Figure 5 (left) in red text were comparable to those identified in the previous work with the alternative workflow^{4,5} suggesting robustness of both data processing workflows.

Peak area response plots comparing peltatoside (Quecetin-3-O-araboglucoside) and luteolin 7- glucoside for apical leaves and


Figure 5. T-test analysis of leaf vs. flower grown under white light. The volcano plot (left) provides both statistical significance (p-value) on the Y axis and fold change for compounds on x-axis used in t-test. The compounds names marked in red are same compounds identified from the interest list generated from MarkerView software in previous technical note², suggesting the robustness of the of both workflows. Area response plots (right) from various samples depict higher levels of peltatoside and luteolin 7- glucoside in flowers compared to apical leaves.

flowers exposed to white light showed higher levels of these compounds in flowers as opposed to leaves (Figure 5, right).

In this study, popular flavonoids such as hyperoside, peltatoside, naringenin, quercetin, luteolin glucosides are upregulated in flowers. It can be also established that flowers either exposed to white or blue/red light have shown higher levels of flavonoids than leaves. Flavonoids are well known for their UV protection and their functional role ROS scavenging under high light stress conditions.⁶ Additionally, these compounds are largely responsible for colors and aroma of the flower.⁷ Further, it will be interesting to investigate the role of Eriodictyol-7-O-glucoside in flower tissues due to its relative abundance in flowers exposed to white light compared to blue/red light.

Conclusions

In this study, *O. basilicum* plants grown under different light wavelengths and SWATH DIA was used to examine the metabolic changes that occurred in different plant parts. Here, an alternative data processing approach was used for metabolite identification through library searching in SCIEX OS software.² Then the identified metabolites were analyzed using MarkerView

software to find the differentially regulated plant metabolites. Both workflows are supported within SCIEX OS software.

- Rich dataset for untargeted plant metabolomics was collected using data independent acquisition (SWATH DIA), providing high-resolution MS and MS/MS data for all detectable precursors in a sample in a single run
- Metabolites were identified using a library searching that leverages the MS and MS/MS from SWATH DIA data
- Statistically significant metabolite differences and trends were easily visualized using MarkerView software, using the various multivariate statistical analysis and visualization tools (PCA, PLS-DA, t-test) to distinguish biological samples. With this workflow, metabolite names are displayed rather than m/z / RT features with the previously described workflow.²

References

1. d'Aquino L. *et al.* (2022). *Growth and metabolism of basil grown in a new-concept microcosm under different lighting conditions*. [Scientia Horticulturae 299:111035](#).
2. Global metabolite profiling of leaves and flowers of *Ocimum basilicum* L. grown in a microcosm under different light wavelengths. [SCIEX technical note, RUO-MKT-02-14157-A](#).
3. K Patel, *et al.* (2018) New insights into the medicinal importance, physiological functions and bioanalytical aspects of an important bioactive compound of foods 'Hyperin': Health benefits of the past, the present, the future” [Beni-Suef University Journal of Basic and Applied Sciences, 7\(1\), 31-42](#).
4. MarkerView software - [SCIEX community posts](#).
5. [Natural Products HR-MS/MS Spectral Library 1.1](#).
6. [BMC Plant Biology](#) volume 17, Article number: 64 (2017)
7. Torawane SD *et al.* (2020) Controlled release of functional bioactive compounds from plants. [Enc. Active Mol. Del System 103-110](#).

Analysis of untargeted metabolomics data from an untargeted Zeno data-dependent acquisition (DDA) workflow using SCIEX OS software

A metabolomics data processing pipeline for DDA data acquired on the ZenoTOF 7600 system

Rebekah Sayers¹, Robert Proos² and Paul RS Baker²
¹SCIEX UK; ²SCIEX US

This technical note describes the best practices for analyzing untargeted metabolomics data acquired using the ZenoTOF 7600 system. Previous work has provided a step-by-step overview of instrument parameter settings to obtain optimal untargeted data¹. This work explains the settings and the subtle changes to the processing method in SCIEX OS software that are needed to maximize coverage of the metabolome.

Metabolites are chemically and structurally diverse compounds that are found endogenously at a wide dynamic concentration range, which makes the metabolome challenging to characterize.

Mass spectrometry is typically used for metabolomics analysis in complex biological samples. DDA is a commonly used analytical technique for untargeted metabolomics to detect and potentially quantify all metabolites in a sample (Figure 1). In this experiment, a TOF MS scan detects all precursor ions within a specified mass range, and MS/MS scans are triggered for those precursors that pass pre-defined criteria. Based on adjustments to the criteria, an ideal experiment should provide as much coverage of the metabolites as possible, while maintaining good MS/MS spectral quality. A complementary technical note detailing the optimal instrument parameter settings for DDA analysis using the ZenoTOF 7600 system is available online¹.

Many metabolites are present at low endogenous concentrations due to the nature of the metabolome, which can lead to low-quality spectra and poor coverage. However, the Zeno trap provides significant gains in MS/MS sensitivity and improves MS/MS spectral quality. The use of the Zeno trap can therefore result in higher confidence spectral matching and greater metabolite coverage. Spectral matching is currently reliant on libraries obtained from instruments without a Zeno trap. Therefore, the processing methods must be optimized to realize the true potential and impact of the data.


Figure 1. Targeted and untargeted metabolite screening workflows on the ZenoTOF 7600 system.

Key features of processing Zeno DDA data using SCIEX OS software

- An untargeted metabolomics workflow with the Zeno trap enabled provides high-quality MS/MS spectra
- Current MS/MS spectral libraries contain fragmentation data acquired on platforms without a Zeno trap, which differ in the ratio/intensity of fragment ions
- To maximize the impact of metabolomics data acquired on the ZenoTOF 7600 system, specific parameter settings within the data processing method should be used

Methods

Sample preparation: NIST SRM 1950 plasma was extracted using 4 volumes of ice-cold methanol. After centrifugation to separate the precipitated proteins, the supernatant was directly analyzed by liquid chromatography electrospray ionization tandem mass spectrometry (HPLC-ESI-MS/MS).

Chromatography: Samples were analyzed using an Exion LC system with a Kinetex F5 column (2.1 × 150 mm, 2.6 μm, Phenomenex). A simple linear gradient from 0 to 95% B was used with standard reversed phase mobile phases at a flow rate of 200 μL/min. Mobile phase A was 0.1% formic acid in water and mobile phase B was 0.1% formic acid in acetonitrile. A 1 μL injection volume was used and the column temperature was maintained at 30°C throughout the analysis. The total runtime was 20 min.

Mass spectrometry: The extracted sample was analyzed on the ZenoTOF 7600 system equipped with the OptiFlow Turbo V ion source. Data were collected using a top-40 DDA method, with dynamic background subtraction (DBS) and exclusion for 6 s after 3 occurrences. The method used a TOF MS accumulation time of 100 ms, collision energy (CE) of 30 V and a TOF MS/MS accumulation time of 5 ms. For a detailed explanation of parameter settings used for DDA analysis on the ZenoTOF 7600 system, see reference 1.

Data processing: All data were analyzed using the Analytics module in SCIEX OS software. The MS/MS spectra were compared to the library spectra using the SCIEX All-in-One HR-MS/MS library version 2.0, SCIEX Accurate Mass Metabolite Spectral Library version 2.0 and NIST 2017 MS/MS library. Statistical analysis was performed using MarkerView software.

SCIEX OS software: Project Default Settings

Quantitative processing parameters (peak integration)

SCIEX OS software can be used to process all data types acquired on both triple quadrupole and QTOF instruments, and processing parameters can be adjusted depending on the type of analysis and data processing used (Fig. 2). SCIEX OS software has 3 integration algorithms that can be used to integrate peaks. These options are listed in the Project Default Settings in the Analytics module of SCIEX OS software (Fig. 3).

MQ4: This algorithm selects a low, but not the lowest, concentration standard or quality control sample as the representative sample for the analytical run. This algorithm

requires some user input to define the integration parameters but can enable faster searching. This algorithm is therefore preferentially used for targeted and untargeted identification workflows.

AutoPeak: This algorithm selects a high, but not saturated, concentration standard or quality control sample as the representative sample of the analytical run. Chromatograms within the batch being processed are evaluated to determine which sample is the best peak model for each transition. The peak model is constructed based on a model of 3 Gaussian peaks. This algorithm requires little user input because few parameters are adjustable. This is the optimal choice for targeted quantitative analysis.

Summation: This algorithm does not perform a normal peak search and instead assumes that a peak is present close to the expected retention time. This algorithm works well with flow-injection analysis (FIA) or infusion data.

When performing quantitative mass spectrometry data processing, it is important to determine whether a given peak is significant and exceeds background noise. There are 3 signal-to-noise algorithms available in the software, including relative noise, standard deviation and peak-to-peak. The latter 2 require the noise region and the peak of interest to be selected on the


Figure 2. Process for defining processing parameters in SCIEX OS software, according to workflow



Figure 3. Project default settings. Top: Peak integration and signal-to-noise algorithms and Bottom: library searching algorithm, filtering and scoring settings.

chromatogram and are therefore not suitable for untargeted data analysis. The relative noise algorithm provides a less subjective

means of determining signal-to-noise calculations and allows a background to be calculated even when no peak-free part of the chromatogram is available.

For targeted and untargeted metabolite screening, the MQ4 integration and relative noise algorithms should be selected.

Qualitative processing method parameters

The quantitation and targeted identification workflow can be used for both the quantitative and qualitative analyses of known analytes and uses library searching to confirm the identity of analytes in the sample. There are 3 types of library searches that can be performed.

Candidate search finds the best spectral match from all compounds in the selected libraries. This search will take the longest time to process.

Confirmation search uses the name in the components list to match against compound names in the libraries. The name in the XIC list must exactly match the name in the library or no hit will be found. This should only be used in a targeted identification search, where the target list has been generated from the library compound list.

Smart confirmation search first searches for matching names, and then searches using the spectra if no match is found. Smart confirmation search is recommended for this workflow as the processing time will be the same as the confirmation search if the compound name exists in the library. This search strategy comes with the added flexibility to perform candidate searching when needed. This is the most used search algorithm.

Once the results have been obtained from the library search, it is important to assess the confidence of each of the identifications. There are 3 scores can be computed to sort the results.

Purity score uses all peaks from both spectra when calculating the score. A high value indicates a high probability that the unknown spectrum has been correctly identified and does not contain significant amounts of peaks from additional compounds. A lower value indicates that the match is less certain or additional fragment ion peaks from another compound are present in the unknown spectrum. This is the most used method to sort results.

Fit score is calculated based on the library spectrum and it ignores peaks that are present only in the unknown spectrum. The fit score describes how well the library spectrum is represented in the unknown spectrum. A high fit score with a low purity score indicates that the unknown spectrum is likely impure but contains library compounds.

Reverse fit score is calculated based on peaks that are present in the unknown spectrum and it ignores peaks present only in the library spectrum. The reverse fit score describes how well the unknown spectrum is represented in the library spectrum.

Default settings for quantitative processing and peak integration should be set for each project and cannot be altered within the processing method. Qualitative processing method parameters can be set for the project and can also be modified in the processing method.

SCIEX OS software: Targeted and non-targeted metabolite screening

When performing untargeted data analysis, select either the quantitation and targeted identification or non-targeted screening workflow. The selected project default parameters are used for peak integration and are saved in the processing method file.

For targeted identification, a component list must be entered. This can be inputted manually, imported from a text file or generated from a library database. Inputting the chemical formula and selecting the adduct/charge will automatically

calculate the precursor mass (Da). The peak integration will then be performed using an example data file. This can be reviewed before proceeding or found later in the results table.

If non-targeted screening is used, a component list can be used if desired but the non-targeted search parameters must be defined. If the processing method contains the targeted analytes, then the customized integration parameters for the targeted components will not affect the non-targeted peak integration. To change the project default parameters, the user must create a new non-targeted method. If the parameters are changed in an existing method, the changed parameters will not be implemented.

Library search parameters

Library search parameters have a big impact on the coverage and the quality of spectral matching (Fig. 4). Precursor mass tolerance, collision energy, retention time and polarity settings can be used to filter the data. Precursor mass tolerance should be set according to data type. For accurate mass data acquired on a high-resolution instrument, the precursor mass tolerance should be set to 0.02 Da. Unchecking the Collision Energy and Collision Energy Spread boxes will allow more spectra to be

Figure 4. Library search parameters for DDA data acquired on the ZenoTOF 7600 system

returned. Using CE and CES is recommended if the experimental conditions are the same as those used to build the library, as this will return the most reliable results. Retention time (RT) should only be selected if the search is against a library that includes RTs and if the same chromatographic separation conditions are used. Metabolomics data are typically acquired in both polarities, so selecting this filter will allow library components to be matched correctly to data acquired in the same polarity.

These settings will determine which spectra in the library will be evaluated against the acquired spectra. Good filtering will improve processing time and reduce false positives. Incorrect filtering will cause false negatives. The remaining parameters will affect how the results are scored. Fragment mass tolerance, intensity threshold, minimal purity and intensity are important to consider when comparing data to a library acquired without using the Zenotrapp because of differences in the intensities of fragments.

Intensity threshold sets which fragment ions are used in the search. The default 0.05 value corresponds to 5% and any peak below 5% of the base peak will be ignored. This value can be used to remove small noise peaks from the spectrum and improve purity scores. Decreasing this will allow more of the smaller peaks to be considered during spectral searching, which

might be important to identify a compound. For processing our Zenotrapp DDA data, this number is set to 0.025%.

Minimal purity indicates how well the sample and library spectra match. All peaks from both spectra are used. High values indicate a higher likelihood that the unknown spectrum has been correctly identified and that it does not contain significant amounts of peaks from additional compounds. Lower values indicate that the match is less certain or that fragment ion peaks from another compound are present in the unknown spectrum. Increasing the value of this setting reduces the frequency of lower quality hits.

Intensity factor compensates for differences in peak height between the unknown spectrum and the library spectra. The larger the intensity factor, the greater the difference allowed between the unknown spectrum and the library while still getting a high score. Increasing this value will remove the significance of relative intensity. When comparing spectra generated with the Zenotrapp, the intensity factor should be set to 20 to reduce the importance of variation in the ratio of fragment ions between the unknown and the library spectra. Increasing the intensity factor also increases the purity and these settings change the score calculation. A low intensity factor can yield low scores for true hits, whereas a large intensity factor might introduce false positives.


Figure 5. Confidence level settings. Found under the Flagging Rules tab, these settings enable easy identification of compounds in the results table that may require further investigation before assigning identity

Confidence level settings – Flagging rules

Confidence level settings establish the acceptable boundaries by which data is matched to the spectral libraries and are used to filter any proposed metabolite identifications and characterize the quality of the match. These levels can be set in the Flagging Rules tab (Fig. 5). A score is calculated for the sample metabolite for all criteria selected. A combined score can be generated based on weighting all parameters. The weighing of each qualitative rule can be adjusted according to the number of rules but should total 100%. The results table can be filtered using this traffic light system in which metabolites that fall within the defined confidence intervals will be highlighted in green for acceptable, yellow for marginal and red for unacceptable.

Qualitative rules include mass error (ppm), fragment mass error (ppm), error in retention time, % difference in isotope ratio, library hit score and formula finder score. When performing targeted identification, all but the fragment mass error qualitative rule can be selected. To include a fragment mass error, it must be entered in the components list.

Formula finder

The formula finder algorithm tries to predict the possible chemical formula based on the MS and MS/MS spectra, precursor mass accuracy, isotopic pattern and MS/MS fragmentation.

Formula predictions are made from a defined list of elements and mass tolerance that can be found under the Advanced Tab (Fig. 6). First, select from either naturally occurring compounds or synthetic compounds. In either case, the type of element and/or number of elements to consider must be entered in the Limits section. Proposed formulas are scored based on precursor mass accuracy and average MS/MS mass accuracy of matching fragments. The MS spectrum contributes 67% to the final formula finding score and the MS/MS spectrum contributes 33%. As a result, the ability of the formula to predict the MS mass is the primary influence on the score. However, the matching of the MS/MS fragments also influences the score. The isotope pattern is used to generate the list of found formulas but it is not used to generate the final score. Therefore, a formula with the wrong isotope pattern will probably not be included in the list. A list of possible formulas is determined using precursor mass accuracy, isotopic pattern and MS/MS fragmentation. A high formula finding score does not guarantee that the sample compound is the one identified by the formula finding algorithm because several formulas often match within the mass error. Care must be taken and other confirmatory testing must be done before a compound is identified using formula finding.

SCIEX OS software: Results

After defining the required algorithms and processing parameters discussed above, a results table will be generated (Fig. 7). Extracted ion chromatograms (XICs) are generated for all


Fig 6. Formula Finder settings. An empirical formula can be predicted based on the TOF MS and TOF MS/MS spectra acquired for DDA experiments. These predictions are influenced by the limits on the elemental composition and mass tolerance settings


Figure 7. Results table generated from a targeted identification workflow in SCIEX OS software. Results have been filtered to show only those meeting the highest confidence for each criterion.

metabolites in the library based on thresholds set by the user, such as formula and expected retention time of all target analytes. The MS and MS/MS information is automatically evaluated if the detected XIC signal exceeds the user-defined intensity threshold or signal-to-noise (Figure). Data processing results are ranked based on 4 selectivity criteria to provide a high degree of confidence in assigning compound identifications to detected compounds:

- Retention time matching
- Mass accuracy
- Isotope pattern fit
- MS/MS library searching

In addition, the peak intensity of the analyte can be compared to that of a standard sample of known concentration to obtain quantitative information.

Formula finder results are shown in Figure 7. Clicking on Formula Finder Results will show additional potential candidates. The chemical structure of the selected formula finder results is also shown in the table if the compound has been updated from ChemSpider. By clicking on the ChemSpider icon in the results table, a ChemSpider session is initiated, and a list of all suggested compounds that match the selected formula will be

generated. Selecting a compound in the results list in the app will display the chemical structure, acquired vs matched spectra and a table of matching fragments, together with their mass errors and elemental composition

MarkerView software: Statistical analysis

The MarkerView software is designed to allow the data from several samples to be compared so that differences can be identified. The program uses multivariate analysis (MVA) techniques to compare the samples and provides both supervised and unsupervised methods (Figs. 8 and 9). Supervised methods use prior knowledge about the sample groups (for example, healthy vs diseased) to determine the variables that distinguish the groups. In contrast, unsupervised methods allow the structure within the data to be determined and visualized. The 2 approaches can be combined, for example, as unsupervised methods can be used to determine the groups and then supervised methods can be used to confirm the important variables.

The first stage of data analysis is to perform an unbiased review of the detected features. Using the wizard, the raw data files (*.wiff or *.wiff2) are imported into the software without any prior processing. Data can be assessed by performing a PCA


Figure 8. Unsupervised PCA Analysis. **Panel A:** Two groups of data, for example treated Vs. control samples can be processed via PCA analysis to identify outlier metabolites in the treated group from the loadings plot. These data can be visualized by highlighting the data point(s) from the loadings plot and to generate either a profile plot or a box and whiskers plot as shown in **Panel B**.

analysis. The resulting scores and loading plots enable us to check data reproducibility and quickly identify any outliers. A t-test between groups of interest can show differences between individual features across samples. A volcano plot of p value vs fold-change highlights any significant changes across the samples. These metabolites can be added to an interest list.

These processing steps are used to identify putative features of interest. Data must then be processed in SCIEX OS software using a spectral library to identify the selected features.

Once the features have been identified, a second statistical analysis using the peak list of identified metabolites can be

performed. When the peak list is imported into MarkerView software, the names of the metabolites are shown instead of the feature's m/z and RT. PCA analysis and t-tests on our sample sets. The resulting volcano plot shows which metabolites have

the largest fold-change and which are significantly different between the samples. This approach can enable fast identification of features of interest, which may reveal important biological insights.


Figure 9. Volcano plot data visualization from MarkerView software. Data can be visualized via a volcano plot (Panel A), and as in figure 8, highlighted data can be visualized quantitatively as a profile plot or a box and whiskers plot as shown in (Panel B).

RUO-MKT-02-15619-A

Conclusions

- High-resolution, accurate mass MS/MS data acquired on the ZenoTOF 7600 system at a fixed collision energy may not match as confidently to spectral libraries acquired without the Zeno trap. Herein, optimized parameter settings are provided
- SCIEX OS software allows targeted identification and quantitation and untargeted screening of metabolites. This high degree of flexibility enables the user to optimize processing parameters for a given workflow and/or data type.
- Algorithm and processing parameter choices will impact how many metabolites are confidently identified when comparing Zeno DDA data with conventional spectral libraries
- MarkerView software together with SCIEX OS software provides a complete solution to a metabolomics workflow

References

1. Baker, PRS and Proos R. Untargeted data-dependent acquisition (DDA) metabolomics analysis using the ZenoTOF 7600 system. [SCIEX technical note, RUO-MKT-02-15367-A](#)

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MExplorer Ultimate: Next-Generation Non-Targeted Metabolomics Data Analysis Software Empowers Gout Disease Research

满卓^a, 郝亚芳^b, 司丹丹^a, 龙志敏^a, 刘冰洁^a, 陈爱明^b
Zhuo Man^a, Yafang Hao^b, Dandan Si^a, Zhimin Long^a, Bingjie Liu^a, Aiming Chen^b

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The screenshot displays a software interface with a table of data. The table has multiple columns, likely representing different parameters or data points for the analysis. The interface includes a menu bar at the top and a sidebar on the left.

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The screenshot shows a report generated by MExplorer Ultimate. On the left, there is a header with the SCIEX logo and the text 'MExplorer Ultimate 液-质数据分析报告'. Below this, there are fields for '数据名称: M_010', 'MSE.LS: 100', and '数据日期: 2023-03-10'. There are also two QR codes and the text '大连达瑞 呵护您的每个样本，每个数据' and 'We treasure Your Samples and Data'. On the right, there is a section titled '主要分析概述' followed by a detailed table with columns for '名称', '分子量', '电荷数', '保留时间', '峰面积', '峰高', '峰宽', '峰形', '峰位', '峰速', '峰宽', '峰位', '峰速'. The table contains several rows of data.

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MS-DIAL software parameters for processing untargeted metabolomics data acquired on the ZenoTOF 7600 system

Cagakan Ozbalci¹, Paul RS Baker², and Rebekah Sayers¹

¹SCIEX, UK; ²SCIEX, USA

This technical note demonstrates the importance of parameter settings in MS-DIAL version 4.92 to process untargeted metabolomic data acquired using the ZenoTOF 7600 system. Previous technical notes¹⁻² have shown the importance of instrumental settings for data acquisition and best practices for data processing using SCIEX OS software to produce the highest quality data with broad metabolome coverage. However, different software programs have different algorithms that process data with uniquely defined parameters and settings (**Figure 1**). Furthermore, parameter settings can be instrument-dependent, which requires comprehensive testing to determine the best parameter settings for a particular data set. Herein, MS-DIAL 4.92 software³, a widely used processing tool for metabolomics and lipidomics data analysis, was used to interpret untargeted metabolomics data. Parameter settings were sequentially adjusted through iterative data processing to reveal how setting changes affect metabolomics results. Optimal MS-DIAL software parameter setting values are presented for data acquired using the ZenoTOF 7600 system.

Key features of optimizing MS-DIAL software parameters for processing metabolomics data acquired using the ZenoTOF 7600 system

- MS-DIAL software rapidly processes untargeted metabolomics data acquired on the ZenoTOF 7600 system
- The high sensitivity of the ZenoTOF 7600 system enables low threshold settings within MS-DIAL software that improves the overall coverage of the metabolomics experiment
- MS-DIAL software is compatible with most metabolomic compound libraries


Figure 1. A major bottleneck for untargeted metabolomics workflows is time spent data processing. Understanding the software-specific parameter settings and how they should be adjusted to accommodate specific instrument data is essential to maximize coverage while maintaining high confidence in a timely manner. Herein optimal parameter settings for MS-DIAL 4.92 software are described for untargeted metabolomics data acquired on the ZenoTOF 7600 system.

Introduction

Untargeted metabolomics aims to detect and quantify all observable small biomolecules within a sample to define the metabolic state of an organism and potentially identify biomarkers of disease⁴. High-resolution mass spectrometry (HRMS) analysis using a data-dependent acquisition (DDA) mode is the primary tool for untargeted metabolomics experiments. To increase compound detection and identification (i.e., coverage), experiments are typically run in both positive and negative ion modes. Acquired data are generally processed using software that matches MS/MS spectra to small molecule databases for compound identification.

The breadth of coverage of metabolic compounds by mass spectrometry depends on several factors. First, instrument performance and hardware parameter settings affect metabolomics data quality. Optimal parameter settings for the ZenoTOF 7600 system have been previously reported¹ that leverage the speed and sensitivity of the instrument to maximize coverage from a DDA experiment. Second, data processing parameters significantly impact the identification of metabolites. These parameters can be related to library match score, minimum intensity threshold, etc. Incorrect parameter settings may result in misidentification or missing a compound altogether. Software parameter settings for SCIEX OS software for untargeted metabolomics data acquired using the ZenoTOF 7600 system have been previously determined². MS DIAL software can also process these data; however, its parameter settings are unique from SCIEX OS software. Systematic adjustment of parameter settings and data review are required to find the optimal software parameter settings.

In this technical note, untargeted metabolomics data were acquired from NIST SRM plasma samples using the ZenoTOF 7600 system. Data were processed by MS DIAL 4.92 software using iterative permutations of different parameter settings for *minimum threshold*, *mass slice width*, and *identification cut-off* tabs in the processing workflow. From these processed data, optimal software processing parameters were identified and are reported herein.

Methods

Sample preparation: NIST SRM 1950 samples were extracted by a one-phase liquid extraction. Four volumes of ice-cold ethanol were added to 1 sample volume and vortexed for 30 seconds. Extraction mixtures were centrifuged to separate the precipitated protein debris, and the supernatant was used directly for metabolomics analysis. The supernatant can be stored at -20 °C for future analysis.

Chromatography: Extracted metabolites were resolved using an Exion UHPLC instrument equipped with a Kinetex F5 column (2.1 × 150 mm, 2.6 μm; Phenomenex). The column oven temperature was

40°C with a constant flow rate of 0.2 mL/min. Gradient details are shown in **Table 1**.

Table 1: Chromatographic gradient (flow rate = 0.20 mL/min)

Time (min)	Mobile phase A (%)	Mobile Phase B (%)
0.0	100	0
2.1	100	0
14	5	95
16	5	95
16.1	0	0
20	0	0

Mass spectrometry: Extracted samples were analyzed using a ZenoTOF 7600 system with an OptiFlow Turbo V ion source. A top-40 DDA method with dynamic background subtraction (DBS) and exclusion for 6 s after 3 occurrences was employed for both positive and negative ion modes. Electrospray ionization voltages were set to 5500 V and -4500 V for positive and negative ion modes, respectively. The collision energy (CE) was set to 30 V for positive and -25 V for negative ion mode, and the accumulation time was set to 5ms. The automated calibrant delivery system (CDS) performed automated calibration every 9 samples. A summary of the MS instrument parameters is presented in **Table 2**.

Table 2: ZenoTOF 7600 system source and gas parameters

Parameter	Setting (Pos)	Setting (Neg)
Curtain gas (CUR)	35	35
Ion source gas 1 (GS1)	50	50
Ion source gas 2 (GS2)	70	70
CAD Gas (CAD)	7	7
Source temperature (TEM)	500 °C	500 °C
Ion spray voltage (IS)	5500 V	-4500 V
TOF MS mass range	60-1000 Da	60-1000 Da
Declustering Potential (DP)	40	-40
Time bins to sum	6	6
Accumulation time	5 ms	5 ms
IDA max candidate ions	40	40
Dynamic background subtraction (DBS)	Yes	Yes
Collision energy (CE)	30	-25
Zeno pulsing	Yes	Yes
Zeno threshold	20000	20000

Data processing: All data were processed using MS-DIAL 4.92 open-source software. For identification, two different libraries (*ExpBioInsilico_NEG_VS17.msp* and *ExpBioInsilico_Pos_VS17.msp*) were employed. These libraries can be downloaded from (<http://prime.psc.riken.jp/compms/msdial/main.html>) and can easily be selected in the "Identification" tab of the software.

MKT-30320-A

Project creation and analysis of parameter settings

Project creation

Initial “new project” creation is illustrated in **Figure 2**. Selecting the project file path where the raw files are stored is critical. WIFF and WIFF2 files can be submitted to MS-DIAL without conversion. Click “Soft ionization,” “Chromatography,” and “Conventional LC/MS or dependent MS/MS” for the first three tabs. Other noteworthy options to consider are “Data type” (MS1 and MS2). Choose “Centroid data” for both sections. Lastly, select “Metabolomics” for the target omics and “negative” or “positive” for the ion mode.


Figure 2. New project creation page

Analysis parameter settings

Data collection tab:

In this section, MS-DIAL software parameter settings critically affect coverage and the time needed to process the data files. Significantly, these settings also affect the number of false positive results. In the “Mass accuracy” settings, the default values are recommended. Next, click the “advanced” button, as shown in **Figure 3**. To limit the data processing mass range and potentially reduce the time required for processing, enter the appropriate MS1 and MS2 mass ranges. It is also recommended to define when retention time starts and ends to remove unnecessary data points, such as the washing step, from the analysis. The other settings can be kept as default values. If a powerful workstation computer is used, the “Number of threads” setting can be set to more than 2.


Figure 3. Data collection tab with advanced settings

Peak detection tab:

Under the “Peak detection” tab, there are two parameters to be adjusted for optimal peak detection: “Minimum peak height” and “Mass slice width” (**Figure 4**). One of the goals of this technical note was to investigate how varying these parameters affects the data quality and the overall processing time. For “minimum peak height,” 50, 250, 500, and 1000 amplitude were sequentially selected, and the mass slice was set to 0.05 or 0.1 Da. If the smoothing method needs to be changed, click Advanced in this section for detailed smoothing settings.


Figure 4. Peak detection settings tab with smoothing method details

Identification tab:

In the identification section, an MSP file, which can be downloaded from MS-DIAL's website, should be selected, and the identification cut-off value should be set. The default value for this parameter is 80%, but in this study, values of 70% and 60% were also used to observe the effect on the final results. The remaining parameter settings can be set as shown in Figure 5. If a custom library is to be used, click on the advanced button and select the library. (The library should be stored as a *.txt file.)


Figure 5. Identification tab

Adducts tab:

Metabolite adducts can vary depending on the modifiers added to solvents. For example, in the positive ion mode, a compound can appear as a protonated ion or as a sodium, potassium, or ammonium adduct. The recommended adducts to be selected in both positive and negative ion modes are presented in Figure 6.

Adduct ion setting				User-defined adduct
Molecular species	Charge	Accurate mass [Da]	Included	
[M+H] ⁺	1	1.007276	<input checked="" type="checkbox"/>	
[M+NH4] ⁺	1	18.033823	<input checked="" type="checkbox"/>	
[M+Na] ⁺	1	22.989718	<input checked="" type="checkbox"/>	
[M+CH3OH+H] ⁺	1	33.033489	<input type="checkbox"/>	
[M+K] ⁺	1	38.963158	<input checked="" type="checkbox"/>	

Adduct ion setting				User-defined adduct
Molecular species	Charge	Accurate mass [Da]	Included	
[M-H] ⁻	1	-1.007276	<input checked="" type="checkbox"/>	
[M+2O-H] ⁻	1	-19.01839	<input type="checkbox"/>	
[M+Na-2H] ⁻	1	20.974666	<input type="checkbox"/>	
[M-C] ⁻	1	34.969402	<input checked="" type="checkbox"/>	
[M-E-2H] ⁻	1	36.948606	<input type="checkbox"/>	
[M+FA-H] ⁻	1	44.998201	<input checked="" type="checkbox"/>	
[M+Hac-H] ⁻	1	59.013851	<input checked="" type="checkbox"/>	

Figure 6. Adduct ion setting tab

Alignment tab:

This is the final section to be completed before processing a sample or a sample batch (Figure 7). A reference sample file should be chosen from a pooled QC sample or from the sample with the highest metabolite concentration. The remaining parameters (and those under the Advanced tab) were set to the default values since they are compatible with the ZenoTOF 7600 system data acquisition. Once all parameters have been set, click "Finish" and wait until the results screen appears.


Figure 7. Alignment tab

Results and discussion

Metabolites extracted from NIST SRM 1950 plasma samples were analyzed in positive and negative ion modes, with two different injection volumes for each polarity. For the positive ion mode, 0.4 μL and 2 μL of the extract were analyzed, whereas 1 μL and 5 μL were injected in the negative ion mode. Tables 3 and 4 show how adjusting the different parameter settings affects the numbers of identified metabolites in the positive and negative ion modes. Only the higher volume injection results are displayed in these tables. For both polarities, the higher injection volume resulted in approximately 1.4 more annotated metabolites than the lower volume injections (data not shown).

Polarity	Volume	Min Threshold	Mass slice width	Identification cut off	Number of hits	Time for analysis
Positive	2 μ L	50	0.1	60%	573	Very long
Positive	2 μ L	50	0.05	60%	651	Very long
Positive	2 μ L	250	0.1	60%	505	Long
Positive	2 μ L	250	0.05	60%	556	Long
Positive	2 μ L	500	0.1	60%	432	Normal
Positive	2 μ L	500	0.05	60%	470	Normal
Positive	2 μ L	1000	0.1	60%	342	Short
Positive	2 μ L	1000	0.05	60%	367	Short
Positive	2 μ L	50	0.1	70%	257	Very long
Positive	2 μ L	50	0.05	70%	276	Very long
Positive	2 μ L	250	0.1	70%	237	Long
Positive	2 μ L	250	0.05	70%	258	Long
Positive	2 μ L	500	0.1	70%	213	Normal
Positive	2 μ L	500	0.05	70%	224	Normal
Positive	2 μ L	1000	0.05	70%	176	Short
Positive	2 μ L	1000	0.1	70%	173	Short
Positive	2 μ L	50	0.1	80%	135	Very long
Positive	2 μ L	50	0.05	80%	141	Very long
Positive	2 μ L	250	0.1	80%	129	Long
Positive	2 μ L	250	0.05	80%	134	Long
Positive	2 μ L	500	0.1	80%	119	Normal
Positive	2 μ L	500	0.05	80%	119	Normal
Positive	2 μ L	1000	0.1	80%	100	Short
Positive	2 μ L	1000	0.05	80%	100	Short

Polarity	Volume	Min Threshold	Mass slice width	Identification cut off	Number of hits	Time for analysis
Negative	5 μ L	50	0.05	60%	343	Very long
Negative	5 μ L	50	0.1	60%	307	Very long
Negative	5 μ L	250	0.05	60%	295	Long
Negative	5 μ L	250	0.1	60%	273	Long
Negative	5 μ L	500	0.05	60%	239	Normal
Negative	5 μ L	500	0.1	60%	231	Normal
Negative	5 μ L	1000	0.05	60%	177	Short
Negative	5 μ L	1000	0.1	60%	172	Short
Negative	5 μ L	50	0.05	70%	217	Very long
Negative	5 μ L	50	0.1	70%	198	Very long
Negative	5 μ L	250	0.05	70%	202	Long
Negative	5 μ L	250	0.1	70%	189	Long
Negative	5 μ L	500	0.05	70%	174	Normal
Negative	5 μ L	500	0.1	70%	168	Normal
Negative	5 μ L	1000	0.05	70%	134	Short
Negative	5 μ L	1000	0.1	70%	132	Short
Negative	5 μ L	50	0.05	80%	95	Very long
Negative	5 μ L	50	0.1	80%	94	Very long
Negative	5 μ L	250	0.05	80%	91	Long
Negative	5 μ L	250	0.1	80%	91	Long
Negative	5 μ L	500	0.05	80%	87	Normal
Negative	5 μ L	500	0.1	80%	87	Normal
Negative	5 μ L	1000	0.05	80%	76	Short
Negative	5 μ L	1000	0.1	80%	75	Short

Table 3. Coverage, time for analysis, and confidence under different parameter settings in the positive ion mode.

Each row in **Tables 3 and 4** is color-coded according to their identification confidence. Dark orange and light orange coded rows denote high confidence, while blue coded rows denote lower confidence and may contain more false positive results. Grey-coded rows represent little increase in the number of hits but indicate significantly increased analysis time. For context, the time for analysis of a single sample is presented in **Table 5**. Changes in the “Min Threshold” parameter can dramatically affect the analysis time but have less affect the number of hits. The identification cut-off percentage has a significant impact on the number of hits. And the “mass slice width” affects the number of hits only for data with low confidence.

Table 5. Typical times for analysis of 1 sample at different threshold values

Threshold	Analysis time (min)	Descriptor
50	12	Very long
250	1.7	Long
500	1.4	Normal
1000	1.2	Short

Table 4. Coverage, time for analysis, and confidence under different parameter settings in the negative ion mode

As shown in **Tables 3 and 4**, a balance must be struck between the analysis time, coverage, and confidence. In the positive ion mode, a minimum threshold of >50 and <250 appears to give the best coverage. Due to the high sensitivity and the high signal-to-noise ratio of the data acquired on the ZenoTOF 7600 system, a value of 250 cps was chosen to give the best coverage using an “identification cutoff” value of 70%. As seen in the tables, however, small changes can significantly affect the overall data, so it is recommended that


Figure 8. Peak spot navigator. MS2 acquired and Ref. matched are selected to filter the metabolites with the most confident annotation.

these software processing parameters be a starting point from which adjustments can be made to generate the best data for a given data set. For example, for data processed from a large cohort study. In that case, increasing the “minimum threshold” to 500-1000 cps may be reasonable to decrease the time needed for analysis. As expected, the data also indicate that a higher sample load (i.e., higher analyte concentration) generates higher-quality data.

The confidence level in the processed data can be adjusted in the “Peak spot navigator.” Clicking on the “Ref. matched” display option will only show compounds with a matching MS/MS reference spectrum in the compound library (**Figure 8**).

Although MS-DIAL’s MS/MS fragment library algorithm matches acquired data with a high degree of confidence, it is recommended to manually validate the metabolite

annotations from the right bottom panel (Identification) and change the annotation, if necessary, as shown in **Figure 9**.


Figure 9. Identification panel. Annotated metabolites can be reviewed from this panel. If the annotation is not correct, it can be easily changed using the MS/MS look-up button (red arrow)

Conclusions

- The speed and sensitivity of the ZenoTOF 7600 system enable the generation of reproducible and high-quality data for untargeted metabolomic analysis
- In MS-DIAL 4.92 software, the “Minimum threshold” value of 250 cps and the “Identification cut off” value of 70% are good initial settings for untargeted metabolomics data acquired on the ZenoTOF 7600 system
- For extensive cohort studies, or in situations where it is necessary to minimize the overall analysis time, the “Minimum threshold” parameter setting value can be increased up to 1000 cps
- The number of analytes identified in a sample is proportional to the concentration of the sample; higher sample volume and concentration are recommended

References

1. Baker, PRS and Proos R. Untargeted data-dependent acquisition (DDA) metabolomics analysis using the ZenoTOF 7600 system. SCIEX technical note, [RUO MKT-02-15367-A](#)
2. Sayers R, Proos R, and PRS Baker. Analysis of untargeted metabolomics data from an untargeted Zeno data-dependent acquisition (DDA) workflow using SCIEX OS software. SCIEX technical note, [RUO-MKT-02-15619-A](#)
3. Tsugawa H, Cajka T, Kind T, Ma Y, Higgins B, Ikeda K, Kanazawa M, VanderGheynst J, Fiehn O, Arita M. Nat Methods. 2015; 12(6):523-6. Doi: 10.1038/nmeth.3393
4. Alseikh, S et al. Mass spectrometry-based metabolomics: a guide for annotation, quantification and best reporting practices. Nat Methods. 2021; 18(7):747-756. Doi 10.1038/s41592-021-01197-1

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Global Natural Products Social Molecular Networking system (GNPS) is combined with SCIEX High Resolution Mass Spectrometry system to study natural products

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Key words: Global Natural Products Social Molecular Networking system, SCIEX High Resolution Mass Spectrometry system, Natural Products

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1. Qualitative flexibility combined with quantitative power -Using the SCIEX ZenoTOF 7600 LC-MS/MS system, powered by SCIEX OS software. SCIEX technical note, RUO-MKT-02-13053-A.
2. Rapid analysis and interpretation of metabolomics SWATH acquisition data using a cloud-based processing pipeline. SCIEX technical note RUO-MKT-02-13056-A.
3. What is principal component analysis and how does it work?SCIEX community post RUO-MKT-18-12137-A.
4. What is principal component variable grouping (PCVG) and how do I use it? SCIEX community post RUO-MKT-18- 12137-A.

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Development of a method on research metabolites in Zea mays L. by SCIEX ZenoTOF™ 7600

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- [1] Shen S, Zhan C, Yang C, Fernie AR, Luo J. Metabolomics-centered mining of plant metabolic diversity and function: Past decade and future perspectives. *Mol Plant*. 2023, 16(1): 43-63.
- [2] Shen J, Liu Y, Wang X, et al. A Comprehensive Review of Health-Benefiting Components in Rapeseed Oil. *Nutrients*. 2023, 15(4): 999.
- [3] Yang R, Deng L, Zhang L, et al. Comparative Metabolomic Analysis of Rapeseeds from Three Countries. *Metabolites*. 2019, 9(8): 161.

Large scale targeted metabolomics assay for quantitative plasma profiling

Using the SCIEX Triple Quad™ or QTRAP® Systems

Kranthi Chebrolu, Robert Proos, Christie Hunter
SCIEX, USA

Untargeted serum metabolomics studies have detected, identified and documented over 75% of endogenous compounds.¹ However, only a handful of robust and sensitive targeted methods are in existence to quantify large panels of endogenous metabolites in a single assay. Methods for the quantification of plasma metabolites have great utility for many applications, including disease research, food omics and pharmaceutical development, among others. However, the wide structural diversity of metabolites, and hence the broad chromatographic and ionization properties of these compounds, have posed an immense challenge in developing high-throughput, large scale targeted LC-MS methods for semi-quantitative purposes.

Here, a large, high-throughput serum metabolite method consisting of 336 serum metabolites was developed and tested using SRM 1950 plasma. The chromatographic separations were performed on a Kinetex F5 (Phenomenex, CA) column using standard reversed-phase mobile phases for reliable separations of a majority of serum analytes. Development of the MRM method was done both from standards and from testing in

matrix, as well as leveraging the CAS numbers, accurate mass information, compound nomenclature and groups from Human Metabolite Data Base (HMDB). Both polarity switching and the Scheduled MRM™ Algorithm were used to ensure the robust analysis of this large panel of analytes in a single injection. The assay was tested on both the SCIEX Triple Quad 7500 LC-MS/MS System – QTRAP Ready and the QTRAP 6500+ LC-MS/MS System.³

Key features of the targeted plasma metabolomics method

- Full solution for targeted quantitative metabolomics for easy adoption on any SCIEX Triple Quad or QTRAP System
- MRM-based assay for reproducible quantification of a broad range of metabolites in biological samples
 - ~85% of MRM transitions with <20% CV in 10 replicates for 734 MRMs in single injection
- Detailed standard operating protocol describing sample preparation through data processing
- sMRM Pro Builder tool for streamlining retention time windows, dwell time determination and assay optimization
- Targeted analysis for simplified data processing in SCIEX OS software
- High data reproducibility demonstrates robustness of the implementation on both the SCIEX 7500 System and the QTRAP 6500+ System


Figure 1. Broad compound class coverage. Targeted MRM assay covers a broad range of metabolites typically found in plasma/serum. The method covers 11 metabolite classes as defined in HMDB.

RUO-MKT-02-13259-A

Methods

Sample preparation: NIST SRM 1950 plasma was extracted using 8 volumes of methanol then spiked with internal standard QReSS kit (Cambridge Isotopes). After centrifugation to precipitate the proteins, the supernatant was directly analyzed. More details on the sample preparation can be found in the SCIEX How method.⁴

Chromatography: Samples were analyzed using the ExionLC™ System with a Kinetex F5 column (2.1 x 150 mm, 2.6 μm, Phenomenex). A simple linear gradient from 0 to 95% B was used with standard reversed phase mobile phases (A = 0.1% formic acid in water and B = 0.1% formic acid in acetonitrile) with a flow rate of 200 μL/min. A 1 μL injection volume was used and the column temperature was maintained at 30 °C throughout the analysis. More details on the chromatography can be found in the SCIEX How method.⁴

Mass spectrometry: The sample extract was analyzed on both the QTRAP 6500+ System and the SCIEX 7500 System. Data were collected using methods built with the Scheduled MRM™ Algorithm as well as unscheduled methods. Finally, this high throughput method is used to analyze 336 serum metabolites (a total of 784 MRM transitions) in a 20 min run period. More details on the MS settings and the MRM transition list can be downloaded from the SCIEX How method.⁴

Data processing: All data was analyzed using the Analytics module in SCIEX OS Software 2.0.


Figure 2. Robust MRM assay with multiple MRMs per metabolite. To increase assay robustness and specificity, over 65% of the 336 metabolites in the assay have more than 1 MRM transition per metabolite.

RUO-MKT-02-13259-A


Figure 3. Reversed phase chromatography approach. The Kinetex F5 column provided good retention of the polar metabolites while also providing separation across the gradient of the many other compound classes. The top graph shows the separation of metabolites detected in positive mode and the bottom graph shows the negative mode metabolites.

Chromatography

Chromatography of metabolites can be difficult due to the broad range of properties across the compounds of interest. Balancing the retention of polar metabolites as well as separation across the gradient of many other classes of metabolites is key to the development of a generic global profiling method. A reverse-phase method with a Kinetex F5 column (Phenomenex) was used to retain both nonpolar and polar analytes in a single separation for this broad targeted metabolite assay. Figure 3 highlights the separation achieved for the metabolites detected in positive ion mode and negative ion mode, color coded by metabolite class. The retention time reproducibility was highly consistent across the 10 replicates, with an RSD of 0.38% in plasma matrix.

It should be noted that the method described here is intended for use as a starter method to cover a broad range of metabolites, metabolite classes and associated pathways. Optimization of the chromatography, compounds monitored, and the specific MRM transitions and parameters could quickly be performed depending on the specific matrix, pathways, and analytes of interest.

Importance of time scheduling MRMs

Injections were performed without time scheduling in positive (blue) and negative (green) ion polarities to detect metabolites and determine retention times. Replicate injections using Scheduled MRM™ Algorithm were then performed in positive and negative ion mode to further refine retention times. The increased dwell times afforded by the use of the Scheduled MRM Algorithm provide improved peak area reproducibility across 10 replicates, as seen in Figure 4.

After removal of a small number of transitions that were poorly detected in this particular plasma sample, 264 MRM transitions were monitored in positive mode for the final dataset. A large improvement in peak area reproducibility was observed with time scheduled (blue solid) and without time scheduled (blue dotted). In the unscheduled method, a dwell time of 5 msec was set, giving a cycle time of 1.8 sec. After time-scheduling, the average dwell time was 64 msec and the target cycle time was set at 1 sec.

A more dramatic effect was observed in negative mode where a much larger number of MRM transitions (470 final MRM


Figure 4. Improvement in data quality with time scheduling of MRM transitions. Ten replicate injections were performed without time scheduling (dotted line) to detect metabolites and to determine retention times. Then using Scheduled MRM Algorithm, 10 more replicates were injected to assess reproducibility gains due to time scheduling. In both positive ion (blue) and negative ion (green) mode, large increases in data reproducibility were observed.

RUO-MKT-02-13259-A


Figure 5. Fast polarity switching for broadest compound coverage in single injection assay. High data quality is maintained when both the positive and negative ion MRM transitions are combined into a single method with fast polarity switching. No decline in data quality was observed when compared to the individual single polarity methods: 80% of MRM transitions analyzed over 10 replicates has %CV < 15%, but with over 700 MRM transitions in the assay.

transitions) was monitored (Figure 4, green lines). Again, a large increase in dwell time is obtained upon using the Scheduled MRM Algorithm for scheduling such a large numbers of MRM transitions, hence the increase in data quality and reproducibility. In this case, with many more analytes, a dwell time of 5 msec was used, giving a cycle time of 3.3 sec for the unscheduled method (dotted green lines). After time-scheduling, the average dwell time was 20 msec and the target cycle time was set at 1 sec (solid green lines).

Fast polarity switching

The next step was to combine both the positive mode and negative mode analytes together into a single, final method. This was possible due to the fast polarity switching available on the SCIEX 7500 System. A total of 734 MRM transitions were combined into the final single, time-scheduled MRM assay with polarity switching, to get wide metabolite coverage in a single assay. Ten replicate injections were performed and the peak area reproducibility was compared to the previous single polarity methods (Figure 5, red line). High data quality was maintained when running the complete assay. There is no decay in the cumulative reproducibility curve from what was seen in negative mode only (470 MRMs) when compared to the full polarity switching method (734 MRMs).

In this final dataset, ~70% of analytes were reproducibly measured across 10 replicates with %CV < 10%, and ~85% were measured (or in other words, about 600 MRM transitions)


Figure 6. Reproducibility across the dynamic range. Wide dynamic range was observed in the detected metabolites using this assay, with over 4 orders of magnitude (top). As expected, the average observed peak area reproducibility across the 10 replicate injections track with peak area, with only the very lowest analytes having higher %CVs (bottom).

with %CV < 20%. The final average dwell time was 14 msec and the target cycle time was set at 1 sec.

Dynamic range and coverage

The ability to precisely quantify both the high abundant metabolites as well as the low abundant analytes in a large multiplex panel is essential for researchers. This assay covers a wide range of metabolites in plasma spread across a broad range of abundance (Figure 6, top). On both the QTRAP 6500+ System and the SCIEX 7500 System, the instrument dynamic range is very high, up to 6 orders of linear dynamic range, enabling this large assay to be performed in a single run using simple sample preparation. The reproducibility observed for the higher abundant analytes was below 5% CV in the final assay. As expected, reproducibility decreases slightly as the abundance decreases, but still remains good down to the very low abundant, difficult to detect metabolites (Figure 6, bottom).

Finally, the method encompasses compounds from 69 pathways including the TCA cycle, glycolysis, purine pyrimidine metabolism and amino acid metabolism (Figure 7).

Conclusions

The large scale polarity switching method using Scheduled MRM Algorithm has been developed to quantify ~336 serum metabolites in plasma and can be used for large scale targeted metabolomics studies.

- The sMRM Pro Builder Tool was developed to streamline optimization of the sMRM assay and provide insights on MRM assay data quality
- Time scheduling of MRMs provided significant improvements over unscheduled methods for this large of a number of MRM transitions, with 95% of positive mode data and 80% of negative mode data having peak area %CV $\leq 15\%$
- Using polarity switching to include both positive and negative ion mode metabolites in a single sMRM method also provided high quality data, with no loss in peak area reproducibility compared to individual single polarity sMRM method
- Improved sensitivity of the SCIEX 7500 System showed up to a 3-5 fold increase in signal to noise ratio over the QTRAP 6500+ System, which provides detection and quantification of ~30% more metabolites in plasma.³


Figure 7. Targeted assay contains metabolites from 69 different metabolic pathways. The compounds present in the method were linked to the pathways they belong to using MetaboAnalyst 5.0 online tool. The compound name list was submitted in the Pathway analysis tool under Annotated features. Further, option *Homo sapiens* (SMPDB) was chosen for library search and match to run the pathway analysis tool. Later, the Pathway_results.csv file was downloaded to obtain all the pathways indicated for the list of uploaded compounds.

References

- Adam MG et al., (2021) Identification and validation of multivariable prediction model based on blood plasma and serum metabolomics for the distinction of chronic pancreatitis subjects from non-pancreas disease control subjects. *Gut* 0: 1–9.
- Yuan M, Breitkopf SB, Yang X, Asara JM. (2012) A positive/negative ion-switching, targeted mass spectrometry-based metabolomics platform for bodily fluids, cells, and fresh and fixed tissue. *Nat Protocols*. 7(5), 872-81.
- Significant sensitivity increases provides 30% more polar metabolites quantified in plasma. [SCIEX technical note RUO-MKT-02-12701-A](#).
- SCIEX How solution for Targeted profiling method for metabolomics. SCIEX How, RUO-MKT-02-12750-A.

Significant sensitivity increases provides 30% more polar metabolites quantified in plasma

With the SCIEX Triple Quad™ 7500 LC-MS/MS System – QTRAP® Ready

Robert Proos, Darren Dumlao, Kranthi Chebrolu, Christie Hunter
SCIEX, USA

Studying the impact of biological perturbations on the metabolome requires high sensitivity and specificity, as well as the ability to multiplex large numbers of known metabolites in a single assay. The utility of LC-MS/MS analysis for the quantification of metabolites from biological samples has become widely accepted as the key analytical solution due to the aforementioned requirements.¹

In addition to very good sensitivity and specificity, targeted analysis using multiple reaction monitoring (MRM) on triple quadrupole MS systems has the advantage of throughput and simplified data processing. By focusing on known important metabolites that are key to biology, quantitative results provide key biological insights. These platforms are capable of fast polarity switching, allowing a broad range of analytes to be interrogated in a single injection. Finally, because of broad dynamic range typically required for biological samples, high sensitivity is a critical attribute of the LC-MS system for the analytes with needing lower detection limits that were previously undetectable and for smaller sample sizes.

A targeted LC-MS/MS assay with MRM transitions for over 450 metabolites was used to assess the impact of added sensitivity on the detection and quantification of metabolites in plasma. Using a generic reversed-phase chromatography method and operating the MS systems in both positive and negative polarity modes to obtain broad coverage, a similar method was run on three separate QTRAP 6500+ systems and as well as three separate SCIEX 7500 systems under independent lab settings.


Figure 1. Increased sensitivity on the SCIEX 7500 System improved the number of metabolites detected. A similar analysis was run on both the QTRAP® 6500+ System and the SCIEX 7500 System was performed. The SCIEX 7500 System showed a 25% increase in metabolite detected was observed in positive ion mode, and a 41% increase was observed in negative mode. Note some metabolites were detected in both polarities and thus the overall gain in detected metabolites was 30%.

Key features of the SCIEX 7500 System for quantification of polar metabolites

- Provides significant gains in sensitivity through increased generation, capture, and transmission of ions with the:
 - OptiFlow® Pro Ion Source with E Lens™ Technology and the D Jet™ Ion Guide²
- NIST SRM 1950 plasma was used as reference material to compare platforms
- Higher sensitivity and S/N on the SCIEX 7500 System yielded large improvements in detected, quantifiable of metabolites—up to 30% more vs the QTRAP 6500+ System (Figure 1)
- Targeted assay is provided, complete with a standard operating protocol, an MRM list, and tools to facilitate method development (sMRM Builder template) can also integrate internal standard kits if required.³

RUO-MKT-02-12701-A

Methods

Sample preparation: NIST SRM 1950 plasma was extracted using 8 volumes of methanol then spiked with internal standard QReSS kit (Cambridge Isotopes). After centrifugation to precipitate the proteins, the supernatant was directly analyzed.

Chromatography: Samples were analyzed using the ExionLC™ System with a Kinetex® F5 column (2.1 x 150 mm, 2.6 μm, Phenomenex). Simple linear gradient from 0 to 95% B was used with standard reverse phase mobile phases (A = 0.1% formic acid in water and B = 0.1% formic acid in acetonitrile) with a flow rate of 200 μL/min. A 1 μL injection volume was used and the column temperature was maintained at 30 °C throughout the analysis.

Mass spectrometry: The same extracted sample was analyzed on both the QTRAP® 6500+ System and the SCIEX 7500 System. Data was collected using methods built with the Scheduled MRM™ Algorithm as well as unscheduled methods. Variations of the method were repeated on two additional instrument sets to confirm the observed sensitivity gains. Detailed method information is available in the supplementary method.³ Finally, this high throughput method is used to analyze over 400 serum metabolites in a 20 min run period.

Data processing: All data was analyzed using the Analytics module in SCIEX OS Software.


Figure 2. Signal to noise improvement for negative ion mode metabolites. Example shown here for Xanthine, created from guanine during purine degradation. Here the improved area and S/N observed when running the same assay on the SCIEX 7500 System as compared to the QTRAP 6500+ System is highlighted.

RUO-MKT-02-12701-A

Arginosuccinate 291 → 70


Figure 3. Signal to noise improvement for positive ion mode metabolites. Example shown here of arginosuccinate, an essential amino acid formed from citrulline in the uric acid biosynthetic pathway, highlights the improved area and S/N observed when running the same assay on the SCIEX 7500 System as compared to the QTRAP 6500+ System.

Sensitivity improvements on the SCIEX 7500 System

A number of innovations on the SCIEX 7500 System have resulted in increased sensitivity for MRM analysis. Here, the impact of these sensitivity gains on a targeted metabolomics assay was evaluated by comparing the signals observed after running a similar assay on six instruments (3 QTRAP 6500+ systems and 3 SCIEX 7500 systems). Example comparisons are shown in Figures 2 and 3 for positive and negative ion mode, highlighting the sensitivity gains on the raw signal. Xanthine analyzed in negative ion mode had a small detectable signal on the QTRAP 6500+ System but had a much improved signal on the SCIEX 7500 System (30x area and 8x S/N gains). An example for positive ion mode is arginosuccinate, showing a 13x gain in area and a 4x gain in S/N.

Perhaps more compelling was the ability to detect metabolites on the SCIEX 7500 System that were undetectable (below the noise) on the QTRAP 6500+ System. Two examples are showing in Figure 4, for both positive ion mode (retinal) and negative ion mode (pyridoxic acid). Retinal is a small molecule metabolized through oxidation of β-carotene in the human body and is essential, especially for night vision. Pyridoxic acid is a catabolic product of vitamin B6 metabolism.


Figure 4. Improved S/N leads to quantification of additional metabolites. Compounds that were not detectable in the plasma sample on the QTRAP 6500+ System are now easily quantified on the SCIEX 7500 System. Retinal was easily quantified in positive ion mode (left) and pyridoxic acid in negative ion mode (right) also showed a large increase in signal and S/N on the SCIEX 7500 System.

Increased numbers of metabolites detected

Next, the results across the assay were compared looking at the peak areas, the peak heights and the S/N between the two instruments and by computing the fold increases for each metabolite. This allowed some general conclusions to be drawn between the two platforms (Table 1). In positive ion mode, the average peak area gain was 17x providing an overall signal / noise gain of 3x. In negative ion mode, the sensitivity increases were a bit more pronounced with a peak area gain of 40x and a S/N gain of 5x. The increased S/N and peak areas led to an increase in the number of metabolites that could be reliably detected and quantified in the assay. Peak detection used a minimum S/N value of 3.

Table 1. Average gains observed on the SCIEX 7500 System vs the QTRAP 6500+ System. The average peak area, peak height and S/N were compared between the two platforms to characterize the average observed gains.

	Average fold increase	
	Positive ion mode	Negative ion mode
Peak area	17	40
Peak height	18	43
Signal / noise	3	5

The resulting gains in metabolites detected and quantified is shown in Figure 1. In positive ion mode, 173 metabolites were detected on the QTRAP 6500+ System and 216 were detected on the SCIEX 7500 System, an increase in detected metabolites of 25%. The results were even more significant in negative mode, with 151 metabolites detected on the QTRAP 6500+ System and 213 detected on the SCIEX 7500 System - a gain of 41%. As some metabolites were measured and detected in both polarities, the overall gain in detected metabolites was 30%.

Confirmation of results on multiple instrument platforms

In order to ensure a robust conclusion, sample analysis was repeated across multiple instrument sets: a total of three QTRAP 6500+ systems and three SCIEX 7500 systems. The S/N for each metabolite detected was measured and then compared between both instruments to compute a fold increase in signal to noise. This was done across the dataset such that a normal distribution for each instrument set could be plotted (Figure 5). Due to the wide range of analytes covered in these global targeted methods, global source conditions were used on both instruments. This is likely the source of some of the larger, or smaller, than expected sensitivity differences between platforms. Nevertheless, the data demonstrated the general increase that can be expected upon translating this global metabolomics method from the QTRAP 6500+ System to the SCIEX 7500 System.


Figure 5. Repeating results on multiple instruments. To confirm the observed sensitivity gains from instrument set 1, similar methods were run on two more instrument sets. As shown in the distributions for the three instrument sets, similar sensitivity gains were observed in both polarities. On average, signal to noise gains of ~3.5 were observed for the SCIEX 7500 System for each instrument set tested. Note general source conditions must be used when monitoring such wide ranging analytes, likely resulting in the broader than expected distribution of S/N gains.

Conclusions

- A targeted MRM assay was developed and run on extracted NIST SRM 1950 plasma in order to compare the performance between the QTRAP 6500+ System and the SCIEX 7500 System
- The higher sensitivity of the SCIEX 7500 System allowed for the detection and quantification of ~30% more metabolites in plasma. Specifically, there was a 25% increase in positive mode, and a 41% increase in negative mode
- The assay was repeated on two other instrument sets to confirm the observation. Using the increase in S/N as the measure of improvement, all three instrument sets showed significant improvements in S/N for a large number of metabolites.

References

1. Yuan M, Breitkopf SB, Yang X, Asara JM. (2012) A positive/negative ion-switching, targeted mass spectrometry-based metabolomics platform for bodily fluids, cells, and fresh and fixed tissue. [Nat Protocols. 7\(5\), 872-81.](#)
2. Enabling new levels of quantification. SCIEX technical note [RUO-MKT-02-11886-A.](#)
3. Targeted profiling method for metabolomics.

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ME/TED: New Generation Software for Targeted-metabolomics Data Analysis

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Targeted metabolomics in rat urine using Zeno MS/MS

Analysis on the SCIEX ZenoTOF 7600 system

Kranthi Chebrolu¹, Jason Causon², Robert Di Lorenzo², Christie Hunter¹
¹SCIEX, USA; ²SCIEX, Canada

Often researchers in academic labs and those performing clinical trials are only given very small amounts of sample, and need to get the most information from the modest volume. This requirement drives the need for targeted assays with high sensitivity and specificity, to obtain high-quality, quantitative results in complex biological samples and to deliver biological insight. High-resolution accurate mass (HRMS) systems have excellent utility in untargeted metabolomics and compound identification in complex matrices because of their high specificity and full-scan MS/MS quality. However, HRMS instruments typically have not been the chosen instrument to perform accurate quantification of a targeted panel of analytes because of their lower sensitivity. In addition, some high resolution platforms cannot keep high resolution at the scan speeds required for fast analyses.

With the introduction of the SCIEX ZenoTOF 7600 system, scientists can now achieve quantitative results at high speeds with high mass accuracy. The core innovation on the ZenoTOF 7600 system is the Zeno trap that when activated, provides significant improvements in duty cycle due to the optimization of ion transmission from the collision cell into the accelerator.¹ This duty cycle improvement provides a substantial increase in MS/MS sensitivity and thus enables targeted high resolution workflows. Here, a small, targeted assay for metabolites in urine from a diabetic rat model was used to characterize the MRM^{HR}

workflow and the impact of Zeno MS/MS on the quality of quantitative results.

Urine samples from a Zucker diabetic rat model were obtained and a targeted assay for 13 metabolites was developed. Sensitivity between the Zeno trap on and Zeno trap off methods was compared and a significant gain in MS/MS signal was observed (for example: cAMP fragment in Figure 1 shows a ~12-fold improvement). In this technical note, the complete workflow from data acquisition through to statistical analysis of results is described.


Figure 1. Significant sensitivity gains in MS/MS. Comparison of extraction ion chromatograms (XICs) for cAMP fragments obtained from MS/MS collect with Zeno trap on (blue) vs. Zeno trap off (pink). Signal/noise ratio improved ~12.5 fold when using the Zeno trap. This XIC is from a 2 μ L injection of diluted urine.

Key features of the ZenoTOF 7600 system for targeted metabolomics

- Significant increases in MS/MS sensitivity due to Zeno trap, which delivers $\geq 90\%$ ¹ duty cycle on MS/MS fragments while still maintaining fast acquisition rates (up to 133 Hz)
- MS/MS XIC peak area gains of ~13 fold with Zeno trap activated for high-sensitivity quantification of detected metabolites
- Flexibility of a QTOF instrument for additional workflow
 - SWATH acquisition, data dependent acquisition, electron activated dissociation (EAD)
- Powerful data processing tools in SCIEX OS software for accurate quantification, followed by multivariate statistical analysis using MarkerView software

RUO-MKT-02-13297-B

Methods

Sample preparation: Urine samples were collected from four distinct rat groups: Zucker diabetic fatty (ZDF) rats, male and female; Sprague Dawley (SD) rats, male and female. Urine was collected from N=5 rats per group. 20 μ L of urine was aliquoted and diluted 10-fold with mobile phase A prior to LC-MS/MS analysis.

Chromatography: An ExionLC AD HPLC system (SCIEX) with a Phenomenex Luna Omega Polar C18, 3 μ m 150 x 2.1 mm (00F-4760-AN) was used for sample separation. A simple linear gradient from 0 to 95% B was used with standard reverse phase mobile phases (A = 0.1% formic acid in water and B = 0.1% formic acid in acetonitrile) with a flow rate of 300 μ L/min. Either a 0.2 or 2 μ L injection volume was used and the column temperature was maintained at 40 $^{\circ}$ C throughout the analysis. The total run time was 13.1 min including 2 min of equilibration.

Mass spectrometry: MRM^{HR} data was acquired on the SCIEX ZenoTOF 7600 system in positive ESI mode using SCIEX OS software. The ion source conditions were as follows: CUR 35, GS1 55, GS2 55, ISVF 5500, TEM 600 $^{\circ}$ C. High resolution MS/MS was collected for each metabolite using an accumulation time of 10 msec. A collision energy (CE) of 30 was used for each MS/MS. Methods were built with the Zeno trap both activated and deactivated to enable the sensitivity comparisons. Three replicates were collected on each sample with each method.

Data processing: MS/MS interpretation, peak integration, and quantitative analysis were conducted in SCIEX OS software, then results were imported into MarkerView software for multivariate statistical analysis (Figure 2). To build a processing method for MRM^{HR} data, the MS/MS spectrum was first examined in the Explorer module to select the best fragment ion. This was also compared to the library spectrum from LibraryView software using the SCIEX Accurate Mass Metabolite Library

Figure 2. Workflow diagram. Data was both acquired and processed using SCIEX OS software. MS/MS was interpreted using both Explorer and Analytics, library searching was performed using the Library View and ChemSpider.

RUO-MKT-02-13297-B


Figure 3. Explorer in SCIEX OS software for metabolite identification. The top pane shows the extracted ion chromatogram (XIC) of phenylalanine, whereas the middle pane shows the experimental accurate mass MS/MS which can be used to select the fragment mass of interest from LibraryView software. The bottom pane shows the structure obtained from ChemSpider that provided in silico (theoretical) fragmentation and accurate mass information that was later used in the calculation of mass error.

(AMMSL 2.0). Structural information from ChemSpider was also used to confirm the identity of the fragment and obtain the theoretical m/z of fragment of interest to be used. This fragment accurate mass information obtained in Explorer mode (Figure 3) was then used to build a final processing method in the Analytics module of SCIEX OS software. Peak areas of the fragment ions were then imported into MarkerView software for statistical analysis.

Zeno MRM^{HR} workflow for targeted quantification

When activated, the Zeno trap provides a significant increase in MS/MS signal on the ZenoTOF 7600 system, while maintaining very high acquisition rates, and not sacrificing mass resolution. Using a targeted MRM^{HR} assay for 13 metabolites in urine, the sensitivity gains due to the activation of the Zeno trap was explored. Extracted ion chromatograms (XICs) were compared from the data collected with the Zeno trap on and off to determine gains in sensitivity. As shown in Figure 1, cyclic AMP produces a dominant fragment ion at m/z 136.0618 with significant signal gains of over 10 fold observed.

To address the concerns of limited sample volume, a comparative experiment was performed using a 0.2 μ L injection volume with the Zeno trap on, and a 2 μ L injection volume with Zeno trap off (Figure 4). Even with ten-fold less sample injected on column, the peak area for the XIC of the m/z 136.0618


Figure 4. Better sensitivity in MS/MS with 10x less sample. Comparison of extraction ion chromatograms for cAMP fragments obtained from MS/MS collect with Zeno trap on (0.2 µL injection, blue) vs. Zeno trap off (2 µL injection, pink). Signal/noise ratio improved ~1.5 fold when using the Zeno trap. MS/MS acquisition rates are very fast (10 msec accumulation) providing 15 data points across the peak at base.

fragment mass with the Zeno trap on ~1.5 fold higher. With higher dilution factors or lower loading of complex matrices on column, matrix effects are reduced improving data quality. Also the ability to analyze much lower sample volumes can be an advantages for some researchers with very precious samples.

The comparison between Zeno trap off vs. on was done for each of the 13 metabolites analyzed and the results are summarized in Table 1. The data quality specifications such as fragment mass error, library match and area gains with the Zeno trap on are presented in Table 1. The mass error for the quantifiable fragment ions were <3ppm for 12 out of 13 compounds analyzed in this study. On average, the significant area gains obtained with the Zeno trap on for MS/MS is 14-fold higher compared to the Zeno trap off. It is important to note that the MS/MS acquisition rate was very high (10 msec accumulation time per MS/MS).

MarkerView software for statistical analysis

Here, a small sample set was explored to test the workflow from quantification to statistical analysis (Figure 5, Table 2).^{3,4} Note these metabolites were selected based on previous results from a SWATH acquisition study on the same sample set.² Metabolites of interest were selected and included in this targeted assay study. Unsupervised principle component analysis (PCA) was used to generate the two-dimensional score plots in MarkerView software (MV). The four categories of mouse models, ZDF-male and female, SD-male and female have clustered differentially, are clearly separated, and 97% of the variance was explained by the PC1 and PC2 (Figure 5, top).

Table 1. Increased quality of MS/MS spectra due to Zeno trap.

Activation of the Zeno trap provided significant MS/MS signal increase and therefore large increases in fragment ion XICs (average of 13.6 fold). High mass accuracy and very good library hits were observed for the resulting Zeno MS/MS spectra.

Metabolite	Fragment ion (m/z)	Library match	MS/MS fragment mass error (ppm)	Area gain with Zeno trap on (on/off)
Acetylglutamate	84.0444	✓	-0.14	12.51
Arginine	70.0651	✓	3.73	13.18
Carnitine	103.0401	✓	-4.52	11.12
Creatine	43.0291	✓	3.87	18.11
Cyclic AMP	136.0618	✓	2.56	10.00
Glutamine	84.0444	✓	1.72	18.08
Histidine	110.0713	✓	0.24	26.72*
Leucine	69.0699	✓	0.56	15.83
Methyladenosine	150.0778	✓	1.28	10.28
Phenylalanine	120.0808	✓	2.37	10.38
Tryptophan	118.0651	✓	0.26	11.93
Tyrosine	119.0495	✓	-3.81	8.89
Uric acid	141.0407	✓	4.24	10.29
Average area gain:				13.64

*Zeno trap off peak area very low, hard to measure

The loading plots showed four PCV groups (data not shown). Metabolites showing large changes on the loadings plot were selected and displayed as box and whisker plots across the samples (Figure 5 bottom). cAMP and methyladenosine had similar pattern of difference across the samples, while creatine showed a different pattern across the ZDF and SD urine samples.

SWATH acquisition to MRM^{HR} workflow

Here, the ability to create a SWATH acquisition to MRM^{HR} workflow for metabolomics was also demonstrated, as SWATH acquisition was performed on the same sample set on the same instrument. Metabolites that showed differences in abundance between experiment groups from the SWATH acquisition data were selected along with a few additional metabolites, and used to build a targeted MRM^{HR} method. Good correlation was observed in the fold change results between the different diabetic mice vs. the SD male sample for the eight metabolites measured in both datasets ($r^2 \geq 0.92$ for all the group comparisons, Figure 6). This highlights the feasibility of


Figure 5. Principal component analysis (PCA) highlights clear differences between samples. Peak areas from MRM^{HR} data were imported into MarkerView software and PCA-PCVG was performed. (Top) Scores plot shows both good reproducibility between replicate injections and clear separation of the different rat urine samples, even with just these 13 targeted metabolites. (Bottom) Box and whisker plots for selected analytes highlights the differences seen across the samples.

performing the non-targeted screening workflow as well as a targeted quantification assay on a single HRAM system.

With the SWATH acquisition workflow, a large number of metabolites can be quantified from a single run and provide preliminary quantitative results to find differences between samples. When an MRM^{HR} assay is next developed for the same analytes, a much more narrow Q1 isolation window is used

RUO-MKT-02-13297-B

Table 2. Differences in metabolite abundances between the rat urine samples. Using the male SD rats for comparison, the peak area fold changes (log₂ peak area ratio) were computed.

Metabolite	Log ₂ peak area ratios		
	SD female / SD male	ZDF female / SD male	ZDF male / SD male
Acetylglutamate	-0.94	0.40	0.36
Arginine	-1.32	-4.42	-1.44
Carnitine	1.77	3.59	2.90
Creatine	-0.38	-1.53	-1.60
Cyclic AMP	0.56	1.97	2.39
Glutamine	3.62	6.34	4.92
Histidine	-0.09	-1.08	-2.33
Leucine	-1.35	-2.06	-2.47
Methyladenosine	-1.11	1.78	2.39
Phenylalanine	-0.76	-2.04	-3.48
Tryptophan	-0.51	-3.50	-3.49
Tyrosine	-0.83	-3.01	-3.94
Uric acid	0.17	1.71	0.73

providing higher specificity of detection, and thus providing an addition confirmation of the screening results. And with the ZenoTOF 7600 system, this can be performed on the same system.


Figure 6. Correlation of quantitative results from SWATH acquisition and MRM^{HR} workflow. The log₂ fold change results were computed for eight compounds across the three groups of comparisons: SD female/SD male, ZDF female/SD male and ZDF male/SD male. Very good correlation was obtained between the two different quantitative techniques.

Conclusions

Here, the targeted MRM^{HR} workflow on the ZenoTOF 7600 system has been explored for use in quantification of metabolites in biological samples.

- Zeno MS/MS provided a 13-fold average increase in MS/MS sensitivity, thus providing both high-quality, full-scan MS/MS data for each metabolite for confident compound identification as well as large increases in fragment ion XIC areas for higher sensitivity quantification
- The sensitivity improvements with the Zeno trap provides the user with more workflow options; for instance, greater sample dilution to reduce matrix effects, or to perform small injection volumes for the analysis of volume-limited samples
- Raw data processing in SCIEX OS software and multivariate statistical analysis visualization using MarkerView software delivers the complete workflow from identification to quantification
- In addition, the ability to transition from non-targeted SWATH acquisition studies to targeted MRM^{HR} workflow on a single MS instrument was demonstrated.

References

1. Qualitative flexibility combined with quantitative power - Using the SCIEX ZenoTOF 7600 LC-MS/MS system, powered by SCIEX OS software. SCIEX technical note, RUO-MKT-02-13053-A.
2. Rapid analysis and interpretation of metabolomics SWATH acquisition data using a cloud-based processing pipeline. SCIEX technical note RUO-MKT-02-13056-A.
3. What is principal component analysis and how does it work? [SCIEX community post RUO-MKT-18-12137-A](#).
4. What is principal component variable grouping (PCVG) and how do I use it? [SCIEX community post RUO-MKT-18-12137-A](#).

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A Underivatized Method for the Determination of 20 Common Amino Acids by Liquid Chromatography Tandem Mass Spectrometry

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SCIEX, China

Key words: Amino Acids, Liquid Chromatography Tandem Mass Spectrometry, Underivatized

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Analytes	Intra_day (n=5, CV%)			Inter_day (n=5, CV%)		
	C _H	C _M	C _L	C _H	C _M	C _L
3-Hydroxyanthranilic acid	1.7	0.8	4.7	1.4	1.2	2.3
3-Hydroxy-L-kynurenine	0.7	0.6	1.3	1.5	1.4	3.5
5-Hydroxyindole acetic acid	1.0	1.5	3.2	1.1	0.9	5.1
5-hydroxytryptophan	0.9	1.2	3.2	1.3	1.2	5.3
Cinnabarinic acid	0.4	0.5	9.6	0.9	1.9	9.8
indole	0.9	1.0	14.1	0.6	1.0	10.2
indole acetic acid	1.2	1.6	6.7	1.3	2.2	6.1
Indole-3-acetamide	0.9	0.4	3.6	0.8	1.2	5.3
Indole-3-carboxaldehyde	0.1	0.8	4.9	1.2	0.9	5.8
kynurenine	0.7	0.9	3.3	0.8	1.2	4.7
Serotonin	0.9	1.7	1.9	1.5	2.0	6.3
Tryptamine	0.8	0.8	8.4	1.8	1.7	5.5
tryptophan	1.0	0.7	2.8	1.8	1.7	3.3
Tryptophol	1.5	1.7	7.7	1.5	1.9	8.1
xanthurenic acid	0.9	1.3	9.5	1.1	1.6	4.2

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1. Yuebiao Feng, Yeqing Peng, Xiumei Song, Han Wen, Yanpeng An, Huiru Tang and Jingwen Wang. Anopheline mosquitoes are protected against parasite infection by tryptophan catabolism in gut microbiota, *Nature Microbiology*, <https://doi.org/10.1038/s41564-022-01099-8>.

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Determination of two histidine methylation by LC-MS/MS

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Haihong Zha, Chenchun Zhong, Zhimin Long, Lihai Guo

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Keywords: histidine methylation; SCIEX Triple Quad™ System;

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- [1] Lv M, Cao D, Zhang L, et al. METTL9 regulates N1-histidine methylation of zinc transporters to promote tumor growth[J]. Cold Spring Harbor Laboratory, 2021. DOI:10.1101/2021.04.20.440582.
- [2] Cao R, Zhang X, Liu X, et al. Molecular basis for histidine N1 position-specific methylation by CARNMT1[J]. Cell Research, 2018, 28(4). DOI:10.1038/s41422-018-0003-0.

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The Determination of Acetyl-CoA and Lactoyl-CoA in Yeast Cells Using LC-MS/MS

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Key Words: Acetyl-CoA; Lactoyl-CoA; Cells

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- [1]. Wang S, Wang Z, Zhou L, et al. Comprehensive Analysis of Short-, Medium-, and Long-Chain Acyl-Coenzyme A by Online Two-Dimensional Liquid Chromatography/Mass Spectrometry[J]. Analytical Chemistry, 2017:12902.

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Determination of 8 Metabolites on main Pathways of Dopamine Biosynthesis and Metabolism by SCIEX Triple Quad™ System

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Lei Min, Long Zhimin, Guo Lihai

SCIEX, 中国

Key Words: Dopamine Biosynthesis and Metabolism Pathway, Target Metabolomics, Dopamine and other metabolites, Biomarker

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- [1] Lister, R. et al. Human DNA methylomes at base resolution show widespread epigenomic differences. *Nature* 462, 315-322[J]. PubMed, 2009. DOI:10.1038/nature08514.
- [2] C.-X. Song, C. Yi, C. He, Mapping recently identified nucleotide variants in the genome and transcriptome, *Nat. Biotechnol.* 2012, 30, 1107-1116.

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Targeted Analysis of Energy Metabolites in Cells by Using SCIEX Liquid Chromatography-Tandem Mass Spectrometry

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Key Words: Energy metabolites; Target Analysis; Cells

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- [1] Han J, Liu Y, Wang R, et al. Metabolic Profiling of Bile Acids in Human and Mouse Blood by LC-MS/MS in Combination with Phospholipid-Depletion Solid-Phase Extraction[J]. Analytical Chemistry, 2015, 87(2):1127-1136.
- [2] Wang X, Xie G, Zhao A, et al. Serum Bile Acids Are Associated with Pathological Progression of Hepatitis B-induced Cirrhosis[J]. Journal of Proteome Research, 2015:150512060716008.
- [3] Xie G, Wang Y, Wang X, et al. Profiling of Serum Bile Acids in a Healthy Chinese Population Using UPLC-MS/MS[J]. Journal of Proteome Research, 2015, 14(2):850-859.

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Determination of 90 Amine- and Phenol-Containing Metabolites in Biological Samples on SCIEX Triple Quad™ 4500 system

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SCIEX, China

Keywords: SCIEX Triple Quad™ 4500 System, Amino Acids, Metabolomics, Biological Samples

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Compound	Q1	Q3	Dwell Time	DP	CE	Compound	Q1	Q3	Dwell Time	DP	CE
Glucosamine	413.14	170.1	8	60	35	Argininosuccinic acid	759.25	170.1	8	60	45
2-Amino adipic acid	395.13	170.1	8	60	30	β -alanine	323.11	170.1	8	60	35
4-Aminobutyric acid	337.12	170.1	8	60	30	L-Homoarginine	422.19	170.1	8	60	40
4-Hydroxyproline	365.12	170.1	8	60	30	L-Homocitrulline	423.17	170.1	8	60	40
5-Glutamylcysteine	484.12	170.1	8	60	40	DL-3-Aminoisobutyric acid	337.12	170.1	8	60	30
5-Oxoproline	363.10	170.1	8	60	30	2-Aminoisobutyric acid	337.13	170.1	8	60	25
Alanine	323.11	170.1	8	60	30	5-Aminovaleric acid	351.13	170.1	8	60	30
Alanyl-glutamine	451.16	170.1	8	60	40	L-2-Aminobutyric acid	337.13	170.1	8	60	30
Arginine	204.60	156.1	8	60	30	DL-2,6-Diaminopimelic acid	424.15	170.1	8	60	40
Asparagine	366.11	170.1	8	60	30	6-Aminocaproic acid	365.15	170.1	8	60	30
Aspartic acid	367.10	170.1	8	60	35	L-Norvaline	351.14	170.1	8	60	30
Citrulline	409.15	170.1	8	60	35	D-(-)- α -Phenylglycine	385.12	170.1	8	60	35
Cystathionine	345.09	170.1	8	60	35	L-Norleucine	365.15	170.1	8	60	30
Cysteine	355.08	170.1	8	60	35	Serotonin	643.20	170.1	8	80	35
Cystine	354.07	170.1	8	60	40	Tryptamine	394.16	170.1	8	60	30
Glutamic acid	381.11	170.1	8	60	30	4-aminohippuric acid	428.13	170.1	8	60	40
Glutamine	380.13	170.1	8	60	30	5-Hydroxyindoleacetic acid	425.12	170.1	8	80	35
Glycine	309.09	170.1	8	60	30	3-Aminosalicylic acid	387.11	170.1	8	60	30
Histidine	622.18	170.1	8	60	45	Taurine	359.07	170.1	8	60	25
Isoleucine	365.15	170.1	8	60	30	Hypotaurine	343.08	170.1	8	60	40
DL-Kynurenine	442.14	170.1	8	60	30	Djenkolic Acid	488.10	170.1	8	60	40
Leucine	365.15	170.1	8	60	30	DL-Ethionine	397.13	170.1	8	60	30
Lysine	307.10	170.1	8	60	30	DL-Homocysteine	369.09	170.1	8	60	30
Methionine	383.11	170.1	8	60	30	DL-Lanthionine	675.16	170.1	8	60	40
Methionine sulfoxide	399.10	170.1	8	60	35	L-Carnosine	460.16	170.1	8	60	30
N-Acetylaspartic acid	409.11	170.1	8	60	35	L-Anserine	354.12	170.1	8	60	30
N-Acetylcysteine	397.09	170.1	8	60	35	Ala-leu	436.19	170.1	8	60	40
Ornithine	599.20	170.1	8	60	45	Sarcosine	323.11	170.1	8	60	25
Phenylalanine	399.14	170.1	8	60	35	3-Methyl-L-histidine	403.14	170.1	8	60	35
Pipecolic acid	363.14	170.1	8	60	30	N ϵ ,N ϵ ,N ϵ -Trimethyllysine	422.21	170.1	8	60	40
Proline	349.12	170.1	8	60	30	Asymmetric dimethylarginine	436.20	170.1	8	60	40
Serine	339.10	170.1	8	60	35	O-acetyl-L-serine	381.11	170.1	8	60	30
Threonine	353.12	170.1	8	60	30	N α -Acetyl-L-lysine	422.17	170.1	8	60	40
Tryptophan	438.15	170.1	8	60	30	DL-5-Hydroxylysine	315.11	170.1	8	60	30
Tyrosine	208.10	156.1	8	60	35	5-Hydroxy-L-tryptophan	687.19	170.1	8	60	40
Valine	351.14	170.1	8	60	30	4-Hydroxy-L-isoleucine	381.15	170.1	8	60	30
4-Aminobenzoic acid	371.11	156.1	8	80	45	Prolinamide	348.14	170.1	8	60	30
Guanine	385.11	170.1	8	60	30	2-Amino-2-methyl-1-propanol	323.14	156.1	8	60	40
2-Aminoethanol	295.11	170.1	8	60	30	Methylguanidine	307.12	170.1	8	60	30
4-Hydroxyphenyllactic acid	416.12	170.1	8	60	30	Homocarnosine	354.12	156.1	8	60	30
Histamine	578.19	170.1	8	60	30	Aminolevulinic acid	365.12	170.1	8	60	30
O-Phosphoethanolamine	375.08	170.1	8	60	30	Metanephrine	431.16	170.1	8	60	40
Boc-D-Tyr-OH	515.18	170.1	8	60	45	Salicylic acid	429.11	170.1	8	60	35
D-Homoserine	353.12	170.1	8	60	30	3-Hydroxyanthranilic acid	387.10	170.1	8	60	30
Saccharopine	510.19	170.1	8	60	45	3-Methyl-L-histidamine	359.15	170.1	8	60	35

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1. Casero, R. A., Marton, L. J. Targeting polyamine metabolism and function in cancer and other hyperproliferative diseases. Nat. Rev. Drug Discovery, 2007, 6, 373-390;
2. Khuhawar, M. Y., Qureshi, G. A. Polyamines as cancer markers : applicable separation methods. J. Chromatogr., B 2001, 764, 385-407;
3. Kevin Guo, Liang Li. Differential 12C-/13C-Isotope Dansylation Labeling and Fast Liquid Chromatography/Mass Spectrometry for Absolute and Relative Quantification of the Metabolome. Anal. Chem, 2009, 81, 3919-3932;

4. Khamis M.M. et al. Comparison of accuracy and precision between multipoint calibration, single point calibration, and relative quantification for targeted metabolomic analysis. *Analytical and Bioanalytical Chemistry*, 2018.

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LC-MS/MS Method for Rapid Determination of Volatile Aldehydes in Cells

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Wang Yang

SCIEX China

Key Words: API 4000™ LC-MS/MS; Tumor; Reactive aldehyde species; Derivatization

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High sensitivity quantification of acylcarnitines using the SCIEX 7500 system

Mackenzie Pearson
SCIEX, USA

Carnitine and acylcarnitines play an important role in the trafficking of fatty acids from the cytosol to the mitochondrial matrix with the additional responsibility of keeping the balance of CoA and acyl-CoA levels within the cell. Acylcarnitines are acetylated forms of L-carnitine produced from carnitine acyltransferases. These acyltransferases cover the entire range of fatty acid lengths, and depending on their chain length specificities, can be found in different tissue types.¹

In many laboratories across the world, acylcarnitines and related molecules from a variety of sample types are quantified to investigate metabolic conditions, such as organic acidemias and fatty acid oxidation defects.² They have also been investigated as potential biomarkers for Type 2 Diabetes and other related conditions.³ Whole blood from dried blood spots (DBS) is most commonly analyzed during screening, however, plasma, serum and urine are other common matrices used for analysis.

Carnitines and acylcarnitines vary in their hydrophobicity, ranging from small, polar carnitine to acylcarnitines with very long acyl chain lengths. This variability poses a unique challenge for chromatographic method development. To overcome this, chemical derivatization techniques are often needed to aid in the retention of these compounds on a reverse phase column.

The detection and quantification of acylcarnitine compounds can be challenging since many species occur as isomers and are


Figure 1. Extracted ion chromatogram (XIC) for 8 unlabeled acylcarnitine compounds.

found at very low concentrations in biological samples. Here, several acylcarnitine compounds were quantified in matrix, using diluted standards and the SCIEX 7500 system. Linearity, lower limits of quantification (LLOQ) and dynamic range were characterized and the sensitivity and reproducibility of the SCIEX 7500 system were evaluated. These experiments were performed without derivatization and with minimal sample preparation to quantify these low-level analytes in a biological matrix.

Key features of the SCIEX 7500 system for the quantification of acylcarnitines

- A high-throughput LC-MS method able to quantify acylcarnitine species
- A highly sensitive method that does not require derivatization of samples before analysis
- A robust method with high sensitivity and reproducibility and broad dynamic range (LDR)

RUO-MKT-02-13839-A

Methods

Sample preparation: Labeled and unlabeled standards were purchased from Cambridge Isotopes (p/n: NSK-B-1, NSK-B-US-1). The unlabeled standards were serially diluted from stock solutions to concentrations ranging from 1:50 to 1:1,000,000. Labeled standards were spiked in at a 1:100 dilution across all samples.

Ten μL of Wistar rat plasma was extracted using 490 μL of 80:20 methanol:water and then vortexed. The samples were spun for 10 minutes at 15000 rpm. An insert vial was loaded with 150 μL of supernatant and 1 μL was injected. Triplicate injections were analyzed.

Chromatography: An ExionLC system with a Phenomenex Kinetex column (C18, 2.6 μm , 150x4.6 mm; p/n 00F-4462-E0) was used to perform the separation. The total run time was 15.5 min. A 1.0 μL injection was used. Table 1 shows the chromatographic gradient used.

Table 1. Chromatographic gradient.

Time (min)	Flow ($\mu\text{L}/\text{min}$)	B Conc (%)	B Curve
3.0	600	15.0	0
3.5	600	75.0	0
11.0	600	99.0	0
12.5	600	99.0	0
13.50	600	15.0	0

Mobile phase A: 99:1 water:methanol with 0.1% Formic Acid
Mobile phase B: 1:99 water:methanol with 0.1% Formic Acid

Mass spectrometry: Samples were analyzed on the triple quadrupole version of the SCIEX 7500 system using the Scheduled MRM algorithm in the positive ion mode. Table 2 lists the source conditions for analysis.

Table 2. Source conditions for the SCIEX 7500 system.

Parameter	Value
CUR	45
CAD	10 (medium)
Temp	475
IS	1400
GS1	45
GS2	70

Data processing: All data were acquired and processed using SCIEX OS software 2.1.0.

Quantification results

Figure 1 shows extracted ion chromatograms for 8 unlabeled acylcarnitine compounds in diluted 1:50 neat solution. Compounds ranging in size from L-carnitine to short-chain acetylcarnitine, to medium-chain octanoylcarnitine, to long-chain palmitoylcarnitine in neat solution at the 1:50 dilution, to highlight the optimized separation.

Next, the standards were diluted into extracted rat plasma across a broad concentration range to determine the sensitivity of the method. Lower limits of quantification (LLOQs) were calculated for all 8 analytes. Table 3 lists all LLOQs for each standard in the assay. The MRM signals observed at the lowest concentration analyzed for 4 of the analytes are shown in Figure 2. These are integrated peaks for the 1:1,000,000 dilution of the unlabeled L-carnitine, acetylcarnitine (short-chain), octanoylcarnitine (medium-chain) and myristoylcarnitine (long-chain) compounds. While the 1:1,000,000 dilution was the final point in the calibration curve, the response from some analytes in the assay suggest that quantification at even lower concentrations might be possible.

In Figure 3, 5 analytes that exhibited 4.5 orders of linear dynamic range (LDR) are shown. The LDR value for each analyte is shown in Table 3. These values were calculated for each analyte based on 3 replicates at each dilution, ranging from 1:50 to 1:1,000,000, of the stock solutions of the commercially available standards.


Figure 2. MRM chromatograms for the 1:1,000,000 dilution point in the concentration curve. Examples of the signals obtained at the lowest concentration measured. Signal was generated at the following concentrations: 152 pM L-carnitine, 38 pM acetylcarnitine, 7.6 pM octanoylcarnitine and 7.6 pM myristoylcarnitine.


Figure 3. Linear dynamic range for several acylcarnitine compounds. Selected compounds exhibited 4.5 orders of linear dynamic range. Note that wider LDR might be possible for some compounds, as detection at lower concentrations is possible based on observed signals.

To assess reproducibility and robustness, 3 replicates of carnitine and acylcarnitines at their respective LLOQ concentrations in matrix were analyzed. A summary of the accuracy and %CV results is shown in Table 3. The reproducibility observed at the LLOQ was <6% across all the analytes, indicating strong reproducibility. The LDR observed for some analytes was 4.5 orders.

Table 3. Calibration curve statistics. The observed LLOQ, accuracy and %CV at the LLOQ and LDR of the calibration for acylcarnitine compounds in matrix is reported.

Name	LLOQ* (pM)	Accuracy	%CV	LDR
<i>L-carnitine</i>	152	106.51	5.88	4.5
<i>acetylcarnitine</i>	38	107.13	5.67	4.5
<i>propionylcarnitine</i>	15.2	113.43	3.81	4
<i>butyrylcarnitine</i>	76	119.45	0.94	3.5
<i>isovalerylcarnitine</i>	15.2	92.98	0.44	4
<i>octanoylcarnitine</i>	7.6	119.31	2.2	4.5
<i>myristoylcarnitine</i>	14.13	91.46	3.87	4.5
<i>palmitoylcarnitine</i>	15.2	99.8	4.45	4.5

*calculated concentration from calibration curve

Next, the endogenous levels of carnitine and acylcarnitines in 10 μL of extracted Wistar rat plasma were assessed. In brief, this plasma was diluted 1:20 utilizing a simple sample preparation, without derivatization or SPE extraction was employed, and 1 μL of plasma was injected on the column. Table 4 illustrates the sensitivity of the SCIEX 7500 system using this minimal extraction protocol. The concentrations of each acylcarnitine compound range from low μM to high nM and provide excellent reproducibility in a biological matrix.

Table 4. Concentration and %CV for endogenous acylcarnitine compounds in extracted Wistar rat plasma.

Name	Concentration (μM)	%CV
<i>L-carnitine</i>	0.8085	2.11
<i>acetylcarnitine</i>	0.7183	4.61
<i>propionylcarnitine</i>	0.0956	5.91
<i>butyrylcarnitine</i>	0.0879	8.89
<i>isovalerylcarnitine</i>	0.0067	1.52
<i>octanoylcarnitine</i>	0.0007	2.12
<i>myristoylcarnitine</i>	0.0033	2.44
<i>palmitoylcarnitine</i>	0.0131	6.32

Conclusions

- The SCIEX 7500 system can quantify acylcarnitine species using LC-separation, yielding a high-throughput, robust and reproducible method that does not require derivatization of samples before analysis
- The SCIEX 7500 system provides unparalleled sensitivity and could be utilized to quantify other low-level acylcarnitines that are present in biological matrices without the need for extensive sample preparation or derivatization
- Since the SCIEX 7500 provides greater depths of sensitivity using plasma, researchers could further dilute the sample or require less sample for extraction, allowing consideration of alternative sample types for analysis.

References

1. Meierhofer D. (2019) Acylcarnitine profiling by low-resolution LC-MS. [*PLoS ONE*, **14**\(8\), 1-11.](#)
2. Han J, Higgins R, Lim MD, *et al.* (2018) Isotope-labeling derivatization with 3-nitrophenylhydrazine for LC/multiple-reaction monitoring-mass-spectrometry-based quantitation of carnitines in dried blood spots. [*Analytica Chimica Acta*, **1037**: 177-187.](#)
3. Shuo Z, Xiao-Fei, F, Ting H, *et al.* (2020) The Association Between Acylcarnitine Metabolites and Cardiovascular Disease in Chinese Patients With Type 2 Diabetes Mellitus. [*Frontiers in Endocrinology*, **11**, 212.](#)

Separation of hexose phosphate isomers using differential mobility spectrometry (DMS)

Separation of five hexose phosphate isomers using the SelexION device on the QTRAP 6500+ system

Catherine S. Lane¹, Denise Mehl², Arno van Rooij³, Alain J. van Gool³, Richard J.T. Rodenburg³, Dirk J. Lefeber³, Marek J. Noga³

¹SCIEX, United Kingdom; ²SCIEX, Germany; ³Translational Metabolic Laboratory, Radboud Consortium for Glycoscience, Radboudumc, Nijmegen, The Netherlands

Sugar phosphates are important metabolites in many biosynthesis pathways including glycolysis, glycogenolysis and biosynthesis of nucleotide sugar building blocks for glycosylation. As such, the determination of their abundances is of key importance for quantitative metabolic profiling experiments and investigation of some metabolic disorders such as congenital disorders of glycosylation (CDG) and defects in glycolysis and gluconeogenesis.^{1,2} However, analytical measurement of these compounds is far from trivial due in part to the presence of structural isomers with similar MS/MS fragmentation patterns. In addition, their highly hydrophilic nature makes their separation by conventional reversed-phase HPLC highly problematic.³ Here, differential mobility spectrometry (DMS) was employed for the

separation of five isomeric hexose phosphate isomers: glucose 1-phosphate (Glc1P), mannose 1-phosphate (Man1P), fructose 6-phosphate (Fru6P), glucose 6-phosphate (Glc6P) and galactose 6-phosphate (Gal6P).

Key features of DMS for distinguishing hexose phosphate isomers

- The analysis of hexose phosphates is challenging due to their isomeric and highly hydrophilic nature
- Differential mobility spectrometry using the SelexION device allows near-baseline separation of Glc1P, Man1P and Fru6P; and Glc6P and Gal6P without the use of chemical modifier
- When all 5 isomers are analysed together separation into three groups is observed: Glc1P/Glc6P, Man1P/Gal6P, Fru6P
- The introduction of chemical modifiers into the DMS cell maximizes the possibilities for isomer separation
- 1-Propanol chemical modifier enhances separation of Man1P from other isomers at the expense of reduced resolution for Glc1P and Gal6P


Figure 1. LC-DMS-MS/MS analysis of mixtures of hexose phosphate isomers (200 µg/mL) in the absence of chemical modifier, using the generic MRM 259->79 with 33 settings of CoV spanning the range 6.2 to 12.6 V. SV = 4100 V, resolution gas = Medium. Each colored trace shows the MRM 259->79 extracted ion chromatogram at a different value of CoV. (A) Glc1P, Man1P, Fru6P; (B) Glc6P, Gal6P.

RUO-MKT-02-14374-A

Methods

Sample preparation: Isomeric hexose phosphates were purchased from Sigma-Aldrich as disodium salts. They were diluted in 10 mM ammonium bicarbonate, 20% acetonitrile in water or 10 nM ammonium bicarbonate in water to concentrations of between 2 and 200 µg/mL for infusion or LC-DMS-MS/MS analysis. Compounds were analyzed both individually and as mixtures.

SelexION device conditions: Experiments were performed using a QTRAP 6500+ system equipped with a SelexION device (Figure 2), the fundamental properties of which have been described elsewhere.⁴⁻⁶ The DMS cell was mounted between the sampling orifice of the mass spectrometer and the ion source. The temperature of the DMS cell was maintained at 150°C, with nitrogen curtain gas operated at 30 psi. Resolving gas (nitrogen) was employed to enhance the separation of the isomers. The DMS cell was used both with and without chemical modifier (1-propanol) added into the curtain gas flow at 1.5% (mole ratio). For DMS infusion experiments, samples were infused at 7 µL/min in 10 mM ammonium bicarbonate, 20% acetonitrile in water. Separation Voltage (SV) was stepped between values of 0 and up to 4200 V. At each value of SV, compensation voltage (CoV) was scanned over values in the range -80 V to +20 V. At each value of CoV, full scan linear ion trap CID MS/MS data were acquired for the sugar isomers. The resulting plots reveal the optimal CoV at which a particular ion is transmitted through the DMS cell at a fixed value of SV. After determination of optimal CoV values for DMS-separated isomers, CoV values can be fixed to allow the selective transmission of individual isomers in a mixture.

LC conditions: Multiple reaction monitoring (MRM) transitions were optimized targeting fragment ions either common to all 5 hexose phosphate isomers, or unique to/predominantly found in one or a subset of the isomers. These MRMs were used to explore DMS separation on an LC timescale. A SCIEX ExionLC AD system with a Phenomenex Kinetex F5 column (2.6 µm, 100 x 2.1 mm) at 30 °C with a gradient of 2% to 20% mobile phase B in 1.8 min (total analysis time 6 min) was used at a flow rate of 300 µL/min. Mobile phase A was 10 mM ammonium bicarbonate in water and mobile phase B was 10 mM ammonium bicarbonate in 80% acetonitrile, 20% water. The injection volume was set to 10 µL.

Figure 2. Configuration of the SelexION device. The DMS cell is mounted in the atmospheric region between the ion source and sampling orifice of the mass spectrometer, allowing the user to easily install the DMS device in minutes.

DMS separation of hexose phosphate isomers by infusion

DMS-MS was employed for the analysis of deprotonated molecules of Glc1P, Man1P, Fru6P, Glc6P and Gal6P, with and without the presence of 1-propanol chemical modifier. The controlled addition of 1-propanol vapor in the DMS cell induces different shifts in optimal CoV for individual isomers, resulting in separation in CoV space that can be orthogonal to the separation obtained when no chemical modifier is employed. Isomers were analysed both individually and as mixtures.

Without modifier, baseline separation of Glc1P/Glc6P, Man1P/Gal6P and Fru6P was observed. Only partial separation of Glc1P from Glc6P was observed, and Gal6P behaved similarly to Man1P (Figure 2). With 1-propanol modifier, separation of Glc1P/Gal6P, Man1P, Glc6P and Fru6P was observed, but Gal6P behaved similarly to Glc1P (Figure 3).



Figure 3 Separation of hexose phosphate isomers in the absence of chemical modifier. Isomers were infused individually at 2 µg/mL (A) or as mixtures at 5 µg/mL (B & C) using SV = 4100 V; resolving gas set to High.


Figure 4 Separation of hexose phosphate isomers in the presence of 1-propanol chemical modifier. Isomers were infused individually at 5 µg/mL (A) or as mixtures (B & C) at 5 µg/mL using SV = 3750 V; resolving gas set to Medium.

RUO-MKT-02-14374-A

LC-DMS-MS/MS analysis of hexose phosphate isomers

LC-DMS-MS/MS was used for the analysis of hexose phosphate isomers individually and as mixtures. Experiments were performed in the absence of chemical modifier; SV was set to 4100 V; resolving gas was set to Medium. For initial experiments, the MRM transition 259->79 (common to all isomers and corresponding to the parent ion to PO₃- fragment ion transition) was monitored using 33 different values of CoV spanning the range 6.2 V to 12.6 V with a step size of 0.2 V. The hexose phosphate isomers, chromatographically indistinguishable under these conditions, could be separated in CoV space (Figures 1 and 5).

The specificity of the LC-DMS-MS/MS separation for Glc1P, Man1P and Fru6P was calculated by observing the amount of non-specific, or breakthrough, signal for the generic 259->79 MRM transition at optimal CoV values for each of the three isomers during single-isomer analysis. Non-specific signal is reported as the area of any breakthrough peak as a percentage of the MRM peak area at the isomer's optimal CoV (Table 1).

Table 1. The specificity of the LC-DMS-MS/MS separation of Glc1P, Man1P and Fru6P. Non-specific signal observed for the MRM 259->79 as peak area percentage of signal observed at optimized CoV for LC-DMS-MS/MS analysis of individual isomers. SV = 4100 V, resolving gas = Medium.

	Glc1P CoV 7 V	Man1P CoV 9.2 V	Fru6P CoV 11 V
Glc1P 200 µg/mL	-	0.11%	0.07%
Man1P 200 µg/mL	0.29%	-	0.02%
Fru6P 200 µg/mL	0.08%	0.16%	-


Figure 5. Determination of optimal CoV values for LC-DMS-MS/MS analysis. LC-DMS-MS/MS analysis of individual hexose phosphate isomers (200 µg/mL) in the absence of chemical modifier, using the generic MRM 259->79 with 33 settings of CoV spanning the range 6.2 to 12.6 V. SV = 4100 V, resolving gas = Medium. Results for the individual experiments are shown as follows: (A) Glc1P, (B) Man1P, (C) Fru6P, (D) Glc6P, (E) Gal6P.

RUO-MKT-02-14374-A

Conclusions

Here, differential mobility spectrometry (DMS) has been used for the analysis of five hexose phosphate isomers. Baseline separation of a mixture of three isomers, Glc1P, Man1P and Fru6P, was observed both in the absence and presence of 1-propanol chemical modifier. Likewise, when Glc6P and Gal6P were analysed together, they could be separated under both sets of conditions. The DMS separations were demonstrated both by infusion and by LC-DMS-MS/MS analysis.

When all five isomers were analysed as a mixture in the absence of chemical modifier, baseline separation of Glc1P/Glc6P, Man1P/Gal6P and Fru6P was observed. With 1-propanol modifier, separation of Glc1P/Gal6P, Man1P, Glc6P and Fru6P was observed. Future work will focus on the evaluation of a wider range of chemical modifiers to establish whether a set of conditions can be found under which Glc1P and Gal6P can be separated from the other isomeric forms.

References

1. Orvisky E *et al.* (2003) Phosphomannomutase activity in congenital disorders of glycosylation type Ia determined by direct analysis of the interconversion of mannose-1-phosphate to mannose-6-phosphate by high-pH anion-exchange chromatography with pulsed amperometric detection. *Anal. Biochem.* **317**(1), 12-18.
2. Tegtmeyer LC *et al.* (2014) Multiple Phenotypes in Phosphoglucomutase 1 Deficiency. *N. Engl. J. Med.* **370**(6), 533-542.
3. Hinterwirth H *et al.* (2010) Selectivity issues in targeted metabolomics: Separation of phosphorylated carbohydrate isomers by mixed-mode hydrophilic interaction/weak anion exchange chromatography. *J. Sep. Sci.* **33**, 3273-3282.
4. Shvartsburg AA, (2008) Differential Ion Mobility Spectrometry. CRC Press, Boca Raton
<https://doi.org/10.1201/9781420051070>.
5. Krylov EV, Nazarov EG, Miller RA. (2007) Differential mobility spectrometer: Model of operation. *Int. J. Mass Spectrom.* **226**, 76-85.
6. Schneider BB, *et al.* (2017) Maximizing Ion Transmission in Differential Mobility Spectrometry. *J. Am. Soc. Mass Spectrom.* **28**, 2151-2159. [RUO-MKT-02-3251-C](https://doi.org/10.1002/jms.3951).

Simultaneous quantification of trimethylamine oxide and its precursors from gut microbial metabolic pathway

Using the QTRAP 4500MD system

Cui Jingwen¹, Zhang Yuanyuan¹ and Li Guoqing¹
¹SCIEX, China

Trimethylamine N-oxide (TMAO), a gut microbiota-derived metabolite, is found in higher concentrations in plasma after ingestion of L-carnitine and phosphatidylcholine. Ingestion of choline, L-carnitine, betaine and foods such as red meat cause significant increase in trimethyl amine (TMA) due to the enzymatic activity of gut microbiota. TMA is absorbed into the bloodstream and metabolized by flavin-containing mono oxygenase (FMO) enzymes to trimethylamine oxide (TMAO) in the liver. TMAO has been implicated in the pathology of multiple diseases, including heart failure.¹

Therefore, establishing a rapid and accurate quantitative method for routine monitoring of metabolites involved in TMAO production is needed to support gut microbiome research and screening of hundreds samples from a large research cohorts.

Here, an LC-MS/MS method for the simultaneous quantification of acetyl-carnitine, carnitine, choline, TMAO and betaine in plasma has been developed. All target compounds were extracted using a simple sample extraction procedure then diluted before analysis using LC-MS/MS.²

Key features of QTRAP 4500MD system for gut microbiome studies

- Robust LC-MS/MS system for routine quantification studies with the Turbo V ion source and Curtain Gas interface
- Simplified plasma sample preparation, using a one-step extraction and dilution
- Fast, reversed-phase chromatography provides good separation of the target analytes
- High accuracy and precision values for all analytes monitored, meeting bioanalytical requirements.


Figure 1. Fast chromatography. A chromatographic peak of trimethylamine N-oxide Q1/Q3 76.1/58.1 using a rapid 5.5 min method for quantification.

RUO-MKT-02-13913-A

Methods

Sample preparation: A 20 μL plasma sample was extracted with 100 μL of acetonitrile and further spiked with 100 μL of a 5 ng/mL deuterated internal standard mixture. The sample mixture was centrifuged at 14000 rpm for 5 min, then 5 μL of the supernatant was injected on the column for analysis. The intra-day and inter-day variability over 5 days of the assay was evaluated.

Chromatography: A Jasper HPLC system (SCIEX) with a Phenomenex Kinetex C18, 2.6 μm , 2.1x100 mm column (00D-4462-AN) was used for sample separation using gradient conditions (Table 1) and a flow rate of 0.6 mL/min. A 5 μL injection volume was used and the column temperature was maintained at 40°C throughout the analysis. The total run time was 5.5 min.

Mass spectrometry: A QTRAP 4500MD system with a Turbo V ion source was used in positive ionization mode. All source and compound parameters were optimized and are outlined in the SCIEX How method.³ The MS method included optimized MRM transitions for multiple analytes.

Data processing: The method development and data acquisition was carried out using Analyst MD software. Data processing was performed using MultiQuant MD software.

Table 1. Fast gradient conditions.

Time	Mobile phase A	Mobile phase B
0	10	90
3.5	50	50
3.6	10	90
5.5	10	90

Mobile phase A – water with 0.1% formic acid in 10 mM ammonium formate.

Mobile phase B - 90% acetonitrile with 0.1% formic acid in 10 mM ammonium formate

Rapid analysis of TMAO and related species

First, phosphate buffered saline (PBS) was used to dilute the standards to make a standard calibration curve and obtain a first rough assessment of sensitivity. The concentration range evaluated and the linearity of the calibration curves is outlined in Table 2 along with the observed LLOQ concentrations for TMAO, betaine, L-carnitine, acetyl-carnitine, and choline. All analytes showed good gaussian peak shape (Figure 2).

Triplicate injections of TMAO, betaine, L-carnitine, acetyl-carnitine and choline were each analyzed at different concentrations and across 5 days, as outlined in Table 3. The intra-day reproducibility was analyzed across triplicate injections from a single day. The accuracy for each metabolite ranged from 89.87% to 112.03% and the %RSD ranged from 1.81% to

Table 2. Concentration curves. Compounds were diluted with phosphate buffered saline (PBS) to create standard concentration (ng/mL) curves for each analyte. Good linearity was achieved for each of the standards. Observed lower limits of quantification in buffer are highlighted in bold, using standard bioanalytical requirements of <20% CV and \pm 20% accuracy.

Series	Acetyl-carnitine	Carnitine	Betaine	Choline	TMAO
STD1	0.03	0.05	0.05	0.05	0.13
STD2	0.05	0.1	0.1	0.1	0.25
STD3	0.1	0.2	0.2	0.2	0.5
STD4	0.25	0.5	0.5	0.5	1.25
STD5	1.25	2.5	2.5	2.5	6.25
STD6	2.5	5	5	5	12.5
STD7	5	10	10	10	25
STD8	10	20	20	20	20
Linear equation	$Y=1.54803x-0.009$	$Y=1.43472x+0.04162$	$Y=0.63080X+0.00622$	$Y=0.93725X-0.00116$	$Y=9.55724x-0.54357$
Correlation coefficient	0.9999	0.99118	0.99713	0.99987	0.99792

RUO-MKT-02-13913-A


Figure 2. Lower limits of quantification (LLOQ) for acetyl-carnitine, L-carnitine, choline, TMAO and betaine. Good peak shape and signal/noise were observed at the LLOQ for each analyte in PBS matrix.

8.87%. The inter-day reproducibility was measured by comparing results from triplicate injections across days. The inter-day accuracy ranged from 92.93% to 107.15% and the %RSD ranged from 1.77% to 8.01%. The LLOQs observed for each analyte are marked in bold in Table 2. Standard bioanalytical requirements were followed.

Next, this rapid LC-MS method was tested in plasma. A simple extraction protocol was used, which provided good detection and good signal/noise of TMAO and its precursors analytes in plasma (Figure 3).

Table 3. Quantification results. Inter-day and intra-day results for the analysis of acetyl-carnitine, carnitine, choline, TMAO and betaine in PBS. Precision and accuracy results for a few of the low concentrations points are shown.

Compound Name	Theoretical concentration ng/mL	Inter-day			Intra-day		
		Mean (ng/mL)	RSD%	Accuracy%	Mean (ng/mL)	RSD%	Accuracy%
Acetyl-carnitine	0.08	0.08	4.88%	95.97%	0.078	4.29	97.32
	0.15	0.15	6.44%	96.96%	0.141	4.91	94.33
	7.5	7.36	7.01%	98.14%	6.74	3.16	89.87
Carnitine	0.15	0.16	8.01%	103.33%	0.143	7.89	95.57
	0.3	0.31	6.78%	103.22%	0.291	8.28	97.03
	15	13.94	1.77%	92.93%	13.971	1.97	93.14
Choline	0.15	0.16	4.53%	107.15%	0.166	3.45	110.55
	0.3	0.32	5.07%	105.63%	0.336	2.04	112.03
	15	14.91	2.88%	99.40%	14.509	2.3	96.73
TMAO	0.38	0.38	5.18%	100.35%	0.402	2.93	105.9
	0.75	0.71	3.53%	94.71%	0.716	3.75	95.52
	37.5	37.06	6.91%	98.83%	35.023	1.81	93.4
Betaine	0.15	0.15	7.89%	101.26%	0.15	8.87	100.04
	0.3	0.32	5.45%	106.87%	0.327	2.27	108.99
	15	14.14	3.95%	94.26%	13.788	4.11	91.92

Conclusions

Here, a method for the rapid quantification of TMAO, betaine, L-carnitine, acetyl-carnitine and choline in plasma has been developed on the QTRAP 4500MD system. This method has the advantages of high specificity, linearity and high accuracy. As shown, this method is suitable to support rapid monitoring of TMAO metabolites:

- Rapid 5.5-minute run time provided good peak separation for detection in a shorter run time than several recently published methods for faster screening.⁴
- Consistent quantification results were obtained for all standards tested, with correlation coefficients ranging from 0.991 to 0.999
- Inter- and intra-day reproducibility and accuracy meet standard bioanalytical requirements of %CV less than $\leq 15\%$ in for all analytes analyzed in plasma
- Method was evaluated in plasma matrix using a simple extraction protocol and good signal/noise was observed.


Figure 3. Detection in plasma. The extracted ion chromatograms (XIC) show the separation of TMAO and its precursors in 3.5 min run time from extracted from a plasma sample.

References

1. The gut, the heart, and TMAO. [Cleveland clinic website](#).
2. Velasquez MT, et al. (2016) Trimethylamine N-oxide: The good, the bad and the unknown. [Toxins 8, 326](#).
3. SCIEX How link
4. *Katrina et al.*, (2021) A simplified LC-MS/MS method for the quantification of the cardiovascular disease biomarker trimethylamine-N-oxide and its precursors. [J. Pharma. Anal. 11, 523-528](#).

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SWATHtoMRM: Development of High-Coverage Targeted Metabolomics Method Using SWATH[®] Technology for Biomarker Discovery

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Haihong Zha¹, Yuping Cai^{2,3}, Yandong Yin^{2,3}, Zhuozhong Wang^{2,4}, Kang Li⁴, Zheng-Jiang Zhu^{2,3}, Lihai Guo¹

¹SCIEX, China; ²Interdisciplinary Research Center on Biology and Chemistry, Shanghai Institute of Organic Chemistry, Chinese Academy of Sciences, China

³University of Chinese Academy of Sciences, China; ⁴Harbin Medical University, China

Key Words: Targeted metabolomics; multiple-reaction monitoring (MRM); high coverage; SWATH-MS; biomarker discovery; colorectal cancer (CRC)

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- [1] Nicholson JK, Lindon JC, Holmes E., "Metabonomics": understanding the metabolic responses of living systems to pathophysiological stimuli via multivariate statistical analysis of biological NMR spectroscopic data," *Xenobiotica*, 1999, 29, 1181-9
- [2] Rinschen MM, Ivanisevic J, Giera M, Siuzdak G., "Identification of bioactive metabolites using activity metabolomics," *Nat Rev Mol Cell Biol.*, 2019, 20, 353-367
- [3] Wishart DS, Tzur D, Knox C, Eisner R, Guo AC, et al, "HMDB: the Human Metabolome Database," *Nucleic Acids Res*, 2007, 35, D521-6
- [4] Tsugawa H, Rai A, Saito K, Nakabayashi R., "Metabolomics and complementary techniques to investigate the plant phytochemical cosmos," *Nat Prod Rep*, 2021, 38, 1729-1759
- [5] Kusano M, Tabuchi M, Fukushima A, Funayama K, Diaz C, Kobayashi M, Hayashi N, Tsuchiya YN, Takahashi H, Kamata A, et al, "Metabolomics data reveal a crucial role of cytosolic glutamine synthetase 1;1 in coordinating metabolic balance in rice," *Plant J.*, 2011, 66, 456-66

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Development of a High Coverage metabolomics method of Ipomoea batatas (L.) Lam by SCIEX QTRAP® System

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Lei Min, Long Zhimin, Guo Lihai

SCIEX, 中国

Key word: Ipomoea batatas (L.) Lam, High coverage metabolomics method, Plant metabolomics, Biomarker

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- [1] F Zheng, X Zhao, Z Zeng, L Wang, G Xu, Development of a plasma pseudotargeted metabolomics method based on ultra-high-performance liquid chromatography-mass spectrometry. *Nature Protocols*. 2020.
- [2] Zhou S, Kremling KA, et al. Metabolome-Scale Genome-Wide Association Studies Reveal Chemical Diversity and Genetic Control of Maize Specialized Metabolites. *Plant Cell*. 2019,31(5):937-955.

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Two different activated dissociation techniques provide complementary and information rich for structure elucidation based on SCIEX ZenoTOF™ 7600 system

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Quantitative analysis and structural characterization of bile acids using the ZenoTOF 7600 system

Paul RS Baker¹, Robert Proos¹, Maxim D Seferovic² and Thomas D. Horvath^{3,4}

¹SCIEX, USA; ²Dept of Obstetrics and Gynecology, Baylor College of Medicine, Houston, TX; ³Dept of Pathology, Texas Children's Hospital, Houston, TX; ⁴Dept of Pathology & Immunology, Baylor College of Medicine, Houston, TX

This technical note demonstrates the power and capability of the ZenoTOF 7600 hybrid time-of-flight mass spectrometry (TOFMS) system to sensitively detect, quantify, and structurally characterize the bile acid content of biological specimens. Analysis of bile acids by nominal mass instruments, such as triple-quadrupole MS (TQMS) systems, is challenging because of the high chemical background found in several precursor ions to precursor ion-based multiple-reaction monitoring (MRM) transitions used in current state-of-the-art assays (1-4). High-resolution mass spectrometry (HRMS) generates a full product ion spectrum for each targeted bile acid, and extracting fragment ions with a narrow mass-to-charge (m/z) window can reduce background chemical interferences and improve the signal-to-noise (S/N) of the assay. The detection of individual bile acid isomers currently depends on chromatographic resolution; collision-induced dissociation (CID)-based fragmentation cannot distinguish these isomeric metabolites. Electron-activated dissociation (EAD; 5-7) is a complementary fragmentation mode to

CID that can provide structurally diagnostic fragment ions to differentiate bile acid isomers (Figure 1).

Herein, the ZenoTOF 7600 system was used to quantify the bile acid content of human plasma sample extracts with a sensitivity comparable to that of high-end TQMS systems. Furthermore, the EAD-based fragmentation generated several structurally diagnostic fragment ions that are capable of distinguishing between the bile acid isomers. These data demonstrate the capability of the ZenoTOF 7600 system to perform highly sensitive and selective analysis of endogenous bile acid content in human plasma.

Key features of bile acid analysis on the ZenoTOF 7600 system

- The ZenoTOF 7600 system can detect and quantify bile acids with a sensitivity comparable to that of high-end TQMS instruments
- The narrow fragment ion extraction window possible with HRMS improves assay sensitivity by reducing background chemical interferences and, consequently, increases the S/N of the assay
- EAD can generate diagnostic fragment ions to enable the distinction between bile acid isomers, which improves the overall specificity of the analysis


Figure 1. Bile acid isomer structural characterization by electron-activated dissociation (EAD). Bile acids have multiple structural isomers per each subclass. Deoxycholic acid (DCA), chenodeoxycholic acid (CDCA), and ursodeoxycholic acid (UDCA) have the same empirical formula and exact mass, but they differ structurally with the position and/or stereochemistry of their hydroxyl functional groups. (A) The TOFMS total ion chromatogram (TIC) for m/z 415.3 (DCA sodium adduct) shows the importance of chromatography to resolve bile acid isomers. (B) The EAD MS/MS spectra for DCA+Na and CDCA+Na are presented in an inverted overlay. Underlined in red are EAD-based fragments unique to DCA+Na; the proposed fragmentation of the DCA+Na fragment with m/z 255.2107 is shown on the DCA+Na structure. (C) The DCA+Na-specific fragment at m/z 255.2107 shows no interference in the XIC for chromatographically resolved CDCA+Na. These data suggest EAD fragmentation may be used instead of chromatographic systems with long gradients to distinguish bile acid isomers.

Introduction

Primary bile acids are cholesterol-derived molecules that are synthesized in the liver and are collected and pre-concentrated in the bile contained in the gall bladder (8). During gastric emptying, bile is secreted into the duodenum via the common bile duct, where bile acids play a crucial role in emulsifying and absorbing dietary fats (8,9). The two primary bile acids, cholic acid (CA) and chenodeoxycholic acid (CDCA) are bioconverted into chemically distinct, bacterially produced secondary bile acids and an expanding list of microbially-conjugated bile acids (MCBAs) by host microbiota (10-16). Of considerable recent interest is the potential role of these bacterial-derived secondary bile acids in human physiology and their contributions to pathologies such as inflammatory bowel disease and cancer. Consequently, there is great interest in identifying and quantifying these compounds in diverse samples.

Multiple methods for bile acid analysis have been reported in the literature (1-4). These methods are typically TQMS-based, and their effectiveness highly depends on sample prep, instrument sensitivity, and chromatography to deliver high-quality bioanalytical results. One common drawback of these published methods is that they often use the negative ion mode. At the mass range of interest (~400-500 Da), there is considerable isobaric chemical interference from the matrix and the mobile phase solvents, depending on their source. This is especially problematic for low-resolution TQMS instruments ($R_s \sim 100-150$ equating to ± 0.1 Da mass error). HRMS offers a potential solution to this problem due to its high mass accuracy and high resolution. Fragment ions can be extracted from a high-resolution MS/MS spectrum using a narrow m/z tolerance window (~20 mDa), which may exclude many isobaric interferences that cause the high background in TQMS instruments. The ZenoTOF 7600 system can acquire data in this mass range with a mass accuracy of ± 2 ppm and a resolution of $>35K$.

Bile acids are cholesterol-derived amphipathic molecules of saturated hydroxylated C-24 sterols. As shown in **Figure 2**, the location and stereochemistry of the hydroxyl functional groups define the different isomers within each sub-class. CID-based fragmentation of bile acids primarily yields product ions related to the bile acid head group, and relatively few fragments are generated from the ring structure; hence, CID alone does not generate diagnostically useful fragment ions for identification or quantification purposes. Using chromatography to resolve bile acid isomers requires relatively long chromatographic gradients, which are not conducive to high-throughput analysis. An alternative to lengthy chromatographic gradients and CID-based fragmentation is EAD. EAD-based fragmentation provides ample fragments derived from the sterol ring structure (**Figure 1A & 1B**), and it can distinguish the isomers of each bile acid subclass without extensive chromatographic method development.


Figure 2. The deoxycholic acid sub-class of bile acids. Bile acids are cholesterol-derived amphipathic molecules of saturated hydroxylated C-24 sterols. DCA has a hydroxyl moiety at carbon 12 (C-12) whereas UDCA and CDCA are hydroxylated at carbon-7 (C-7). The latter two isomers are characterized by a different stereochemistry at the number 7 carbon.

The ZenoTOF 7600 system has an EAD cell that provides complimentary fragmentation that can be activated and scheduled within a typical quantitative bioanalytical method to provide structural characterization of bile acids. In this technical note, the speed and sensitivity of the ZenoTOF 7600 system were leveraged to quantify bile acids in prepared plasma samples. These results were compared to those acquired on the QTRAP 7500 instrument. To demonstrate the power of EAD for structural characterization, the DCA sub-class was structurally characterized to identify EAD-specific fragments for each isomer.

Methods

Materials: Bile acid standard mixtures were purchased from Cambridge Isotopes Laboratories. The stable isotope-labeled bile acid mixes (unconjugated, cat# MSK-BA1; conjugated, MSK-BA2) were used as internal standards, and unlabeled bile acid mixes (unconjugated, cat# MSK-BA1-US; conjugated, MSK-BA2-US) were used as primary reference standards. All solvents were LCMS grade and obtained from Burdick and Jackson and Fisher Scientific.

Sample preparation: Human blood (research use only) was collected into tubes containing potassium ethylenediamine tetraacetic acid (K_2EDTA) anti-coagulant following the procedure outlined in the institutional review board (IRB) protocol # _____. The

Table 1: Chromatographic gradient (flow rate = 0.30 mL/min)

Time (min)	Mobile phase A (%)	Mobile Phase B (%)
0.0	75	25
1.5	75	25
13	20	80
14.5	20	80
14.6	75	25
17	75	25

plasma was isolated from the red blood cell fraction by centrifugation. A 50 μ L volume of each plasma sample was pipetted into individual 0.6 mL Eppendorf tubes, and a 200 μ L volume of methanol was added to precipitate plasma proteins (primary DF=5-fold). All samples were vortex-mixed for 30 seconds and centrifuged at 17,000g for 5 minutes, and the supernatant for each clarified plasma extract sample was transferred into fresh vials.

A 50 μ L volume of each clarified plasma extract sample was transferred into fresh glass vials. A 450 μ L volume of an internal standard solution (10 nM deuterated BAs) was added to each sample vial (DF=50-fold overall), and the vials were capped and

Table 2: ZenoTOF 7600 and QTRAP 7500 systems source and method parameter settings

Parameter	ZenoTOF 7600 system	QTRAP 7500 instrument
Curtain gas (CUR)	35	40
Ion source gas 1 (GS1)	60	35
Ion source gas 2 (GS2)	80	70
CAD Gas (CAD)	7	9
Source temperature (TEM)	250 °C	500 °C
Ion spray voltage (IS)	-4500 V/+5000V (CID/EAD)	-2000 V
Declustering Potential (DP)	-80 V/+80V (CID/EAD)	N/A
CID Accumulation/dwell time	Variable (sMRM ^{HR})	Variable (sMRM)
CID Collision energy (CE)	Molecule-dependent	Molecule-dependent
Q1/Q3 Mass Resolution	N/A	Unit/Unit
EAD Accumulation/dwell time	40 ms	N/A
CE for TOF MS	-10 V/+10V (CID/EAD)	N/A
TOF MS mass range	70-1000 Da	N/A
TOF MS/MS mass range	40-400 Da	N/A
Time bins to sum	6	N/A
EAD Collision energy (CE)	12 V	N/A
Electron kinetic energy (KE)	13 eV	N/A
Zeno pulsing for MS/MS	Yes	N/A
Zeno threshold	20000	N/A

MKT-30747-B

vortex-mixed for 30 seconds. The final volume of each sample was split into two identically labeled glass autosampler vials (w/ inserts), with each set being injected on the QTRAP 7500 and the ZenoTOF 7600, respectively.

Chromatography: Extracted metabolites were resolved using two different HPLC systems. For the ZenoTOF 7600-based system, an Exion UHPLC instrument equipped with a Phenomenex Kinetex XB-C18 column (100 x 2.1 mm; 2.6 μ m particle size) was used. For the QTRAP 7500 TQMS-based system, a Nexera 40 Series UHPLC system equipped with a Restek Raptor C18 column (100 x 2.1 mm; 5.0 μ m particle size) and an Ultra C18 guard column (10 x 2.1 mm i.d., 5 μ m; Restek) was used. The two different columns produced similar but not identical chromatographic elution profiles. Notably, both columns entirely resolved the multiple bile acid isomers. The autosampler sample bay was maintained at 10 °C for both systems, and a 15 μ L sample volume was injected on column. The column oven temperature was kept at 50°C with a constant mobile phase flow rate of 0.3 mL/min for a total run time of 17 min. The mobile phase compositions were (A) 10 mM ammonium formate in water and (B) pure acetonitrile. Gradient details are shown in **Table 1**

Mass spectrometry: Extracted sample extracts were analyzed using the ZenoTOF 7600 system with an OptiFlow Turbo V ion source and a QTrap 7500 system with an OptiFlow Pro Turbo V ion source using

Table 3: ZenoTOF 7600 and QTRAP 7500 systems compound parameter settings

Bile acid	ZenoTOF 7600 System			QTRAP 7500 System	
	RT (min)	DP	CE	RT (min)	CE
LCA	11.9	-80	-5	10.7	-23
UDCA	7.4	-80	-5	6.2	-23
CDCA	9.0	-50	-5	7.9	-24
DCA	9.3	-40	-49	8.1	-47
CA	7.1	-40	-49	5.9	-47
GLCA	9.6	-80	-48	6.7	-46
GUDCA	5.6	-80	-82	3.6	-76
GCDCA	7.15	-80	-82	5.1	-80
GDCA	7.4	-80	-82	5.4	-80
GCA	5.7	-150	-48	3.9	-46
TLCA	7.4	-80	-137	6.6	-135
TUDCA	4.5	-80	-142	3.8	-145
TCDCa	5.8	-80	-142	5.2	-145
TDCA	6	-80	-142	5.4	-136
TCA	4.7	-190	-66	4.1	-64

Table 4. Bile acid standards. Mixtures of labeled and unlabelled bile acids were obtained from Cambridge Isotopes Laboratories. Bile acids were analyzed in the negative ion mode, targeting the $[M - H]^-$ ions for quantitation. For qualitative analysis of bile acids using EAD, the $[M + Na]^+$ ion was used. In cases where the precursor and the product ion were the same in the sMRM^{HR} transition, the TOFMS spectrum was used for quantitation.

<i>Bile Acid</i>	<i>Molecular formula</i>	<i>Abbreviation</i>	<i>[M - H]⁻</i>	<i>[M - H]⁻ Product ion</i>	<i>[M + Na]⁺</i>
Lithocholic Acid	C ₂₄ H ₄₀ O ₃	LCA	375.290	375.290	399.288
D4-Lithocholic Acid	C ₂₄ H ₃₆ [2H] ₄ O ₃	D4-LCA	379.316	379.316	403.313
Ursodeoxycholic Acid	C ₂₄ H ₄₀ O ₄	UDCA	391.285	391.285	415.283
D4-Ursodeoxycholic Acid	C ₂₄ H ₃₆ [2H] ₄ O ₄	D4-UDCA	395.310	395.310	419.308
Chenodeoxycholic Acid	C ₂₄ H ₄₀ O ₄	CDCA	391.285	391.285	415.283
D4-Chenodeoxycholic Acid	C ₂₄ H ₃₆ [2H] ₄ O ₄	D4-CDCA	395.310	395.310	419.308
Deoxycholic Acid	C ₂₄ H ₄₀ O ₄	DCA	391.285	391.285	415.283
D4-Deoxycholic Acid	C ₂₄ H ₃₆ [2H] ₄ O ₄	D4-DCA	395.310	395.310	419.308
Cholic Acid	C ₂₄ H ₄₀ O ₅	CA	407.280	343.266	431.278
D4-Cholic Acid	C ₂₄ H ₃₆ [2H] ₄ O ₅	D4-CA	411.305	411.305	435.303
Glycolithocholic Acid	C ₂₆ H ₄₃ N ₀₄	GLCA	432.312	388.324	456.309
D4-Glycolithocholic Acid	C ₂₆ H ₃₉ [2H] ₄ N ₀₄	D4-GLCA	436.337	392.349	460.335
Glycoursodeoxycholic Acid	C ₂₆ H ₄₃ N ₀₅	GUDCA	448.307	74.025	472.304
D4-Glycoursodeoxycholic Acid	C ₂₆ H ₃₉ [2H] ₄ N ₀₅	D4-GUDCA	452.332	74.025	476.329
Glycochenodeoxycholic Acid	C ₂₆ H ₄₃ N ₀₅	GCDCA	448.307	74.025	472.304
D4-Glycochenodeoxycholic Acid	C ₂₆ H ₃₉ [2H] ₄ N ₀₅	D4-GCDCA	452.332	74.025	476.329
Glycodeoxycholic Acid	C ₂₆ H ₄₃ N ₀₅	GDCA	448.307	74.025	472.304
D4-Glycodeoxycholic Acid	C ₂₆ H ₃₉ [2H] ₄ N ₀₅	D4-GDCA	452.332	74.025	476.329
Glycocholic Acid	C ₂₆ H ₄₃ N ₀₆	GCA	464.302	402.302	488.299
D4-Glycocholic Acid	C ₂₆ H ₃₉ [2H] ₄ N ₀₆	D4-GCA	468.327	406.331	492.324
Taurolithocholic Acid	C ₂₆ H ₄₅ N _{05S}	TLCA	482.295	79.958	506.292
D4-Taurolithocholic Acid	C ₂₆ H ₄₁ [2H] ₄ N _{05S}	D4-TLCA	486.320	79.958	510.317
Tauroursodeoxycholic Acid	C ₂₆ H ₄₅ N _{06S}	TUDCA	498.289	79.958	522.287
D4-Tauroursodeoxycholic Acid	C ₂₆ H ₄₁ [2H] ₄ N _{06S}	D4-TUDCA	502.315	79.958	526.312
Taurochenodeoxycholic Acid	C ₂₆ H ₄₅ N _{06S}	TCDC	498.289	79.958	522.287
D4-Taurochenodeoxycholic Acid	C ₂₆ H ₄₁ [2H] ₄ N _{06S}	D4-TCDC	502.315	79.958	526.312
Taurodeoxycholic Acid	C ₂₆ H ₄₅ N _{06S}	TDCA	498.289	79.958	522.287
D4-Taurodeoxycholic Acid	C ₂₆ H ₄₁ [2H] ₄ N _{06S}	D4-TDCA	502.315	79.958	526.312
Taurocholic Acid	C ₂₆ H ₄₅ N _{07S}	TCA	514.284	124.008	538.282
D4-Taurocholic Acid	C ₂₆ H ₄₁ [2H] ₄ N _{07S}	D4-TCA	518.310	124.007	542.307

A scheduled, high-resolution multiple reaction monitoring (sMRM^{HR}) or a scheduled multiple reaction monitoring (sMRM) scan mode, respectively. In some cases, a pseudo-MRM transition was used (i.e., precursor ion to precursor ion). On the QTRAP 7500 system, these molecules were quantitated at the MS/MS level with low CE. However, for the data acquired for these bile acids on the ZenoTOF 7600 system, the TOFMS data were used.

To structurally characterize bile acids, an MRM-triggered data-dependent acquisition (DDA) scan mode was used to acquire EAD-based product ion spectra in the positive ion mode (standards only). Although EAD-based fragmentation is possible in the negative ion mode, it is relatively inefficient compared to the positive ion mode; consequently, the EAD-based analysis was performed in the positive ion mode. Because bile acids are anionic, the sodium adduct for each

MKT-30747-B


Figure 3. Chromatographic resolution of bile acids. Bile acid standards were analyzed using the ZenoTOF 7600 instrument. Chromatography was designed to separate bile acid isomers (e.g., DCA, UDCA, and CDCA) with baseline resolution.


Figure 4. Standard curve for cholic acid (CA). The detection of CA was linear over a 4 orders of magnitude concentration range, with a LOQ calculated to be 0.613 nM in a neat internal standard curve preparation.

was monitored and isolated for fragmentation. A summary of the MS instrument parameters is presented in **Tables 2 and 3**, and a comprehensive list of the compounds studied is shown in **Table 4**.

Data processing: All data were processed using SCIEX OS software. The Analytics module was used for quantitation, and the Explore module was used to structurally characterize bile acids from EAD-based MS/MS spectra. Limits of detection (LOD) and of quantitation (LOQ) estimates were calculated using the standard deviation of the y-intercepts and the mean slope of the 3 (ZenoTOF 7600 system) or 4 (QTRAP 7500 instrument) calibration curves calculated for each analyte as the previously described in the ICH harmonized tripartite guideline (17). The internal standard concentration (9 nM) was included in the linear regression calculations performed using SCIEX

OS; consequently, the LOQ and LOD equations were modified to the following formats:

$$LOD = [BA](nM_{BA}) = [D_x - BA](nM_{D_x - BA}) \frac{3.3 * \sigma_{y-intercept}}{\bar{m}}$$

$$LOQ = [BA](nM_{BA}) = [D_x - BA](nM_{D_x - BA}) \frac{10 * \sigma_{y-intercept}}{\bar{m}}$$

Where [BA]=the bile acid concentration with units of nM_{BA}; [D_x-BA]= IS concentration with units of nM_{D_x-BA}. This correction is necessary to offset the LOQ and LOD when the [IS] is used, which otherwise would generate artifactual values when [IS] ≠ 1.


Figure 5. Standard curve chromatograms for DCA. On the ZenoTOF 7600-based LC/HRMS system, DCA eluted at 9.3 min; its isobar, CDCA, eluted at 9.1 min. Presented here are the individual runs for each standard (in triplicate, from left to right) at each concentration of the standard curve. Data were smoothed once for both presentation and integration purposes. Data for 0.01 and 0.05 nM are not shown for DCA and were not included in the standard curve calculations.

Results and Discussion

A bile acid standard curve was generated neat in solvent and analyzed on the QTRAP 7500 and the ZenoTOF 7600 systems. The parameter settings were optimized for each instrument, as detailed in Tables 2 and 3, and the chromatographic gradient conditions resulted in the separation of each isomeric bile acid with baseline resolution (**Figure 3**). The concentrations within the internal standard curve ranged from 0.01 to 1000 nM, covering a dynamic range > 5 orders of magnitude. Due to the sensitivity of the QTRAP 7500 instrument, the 1000 nM concentration standards caused detector saturation; this was not an issue on the ZenoTOF 7600 system. Consequently, to better match the conditions on each instrument, the final standard curve calculations and subsequent *in vivo* measurements were truncated to a linear dynamic range of 0.1 to 100 nM. Due to this consideration, it is important to note that the results presented here do not necessarily reflect the maximum extent of the dynamic range capability of the instrument for the measurement of bile acids. The range tested, however, was sufficient to measure bile acids in rat plasma. An example curve acquired for cholic acid on the ZenoTOF 7600 system and its associated ion statistics are presented in **Figure 4**.

The individual chromatograms for each injection of the standard curve from 0.1 to 100 nM DCA are presented in **Figure 5**. Its isomer, CDCA, appears in the chromatogram at ~9.1 min, but the chromatographic resolution is sufficient so that there is no interference in the measurement of these two analytes. The measurement of DCA was performed at the TOFMS level rather than

MKT-30747-B


Figure 6 Benefits of a narrow XIC window with HRMS. (Top) The TOFMS mass spectrum of DCA is zoomed around its mass of 391.2854. The blue box overlaid on the spectrum represents the resolution window of a typical QTMS system. The three peaks in the spectrum would all contribute to the overall calculated concentration if analyzed at the MS1 level on the QTMS. (Bottom) Using HRMS, only the peak associated with DCA (Blue) is considered, which reduces interference and chemical background noise (Red 1 and 2). For presentation purposes, the spectral and chromatographic peaks have been smoothed.

the sMRM^{HR} level because even at very low CE values (-5 V), the precursor ion peak intensity was significantly reduced. Some collision energy is required, even in precursor ion to precursor ion transitions, because this voltage helps push ions through Q2 to the Zeno trap and the pulser (for a better understanding of the ZenoTOF 7600 system ion path, see reference 18). The power of HRMS at the MS1 level of analysis is shown in **Figure 6**. In the top spectrum, the

DCA peak is 391.2831 Da. Two potentially interfering peaks are very close to DCA, located at 391.0722 Da and 391.1560 Da, respectively. These two contaminating peaks would be included in any MS1 level analysis on a TQMS system because the isolation width of a nominal mass instrument is ~ 0.7 Da, as indicated by the transparent blue box overlaid on the spectrum. These additions to DCA's total peak area calculation would affect its measurement accuracy. However, because it is possible with HRMS to extract a narrow window around the ion of interest (XIC width of 0.020 Da), these contaminants, as shown in red in the bottom chromatogram, can be excluded, and the correct peak area for DCA can be determined, which increases the overall accuracy of the assay. TQMS systems rely on analysis at the MS2 level to reduce this type of interference.

A standard curve for each analyte for each injection (triplicate injections) was generated, and the data were used to calculate the LOD and LOQ of each bile acid measured with the ZenoTOF 7600 system or the QTRAP 7500 instrument (Table 5). These values were computed using the standard deviation of the y-intercepts and the mean slope for each calibration curve (17) to estimate best the instruments' relative sensitivity for bile acid analysis. With the exceptions of CA, GLCA and GUDCA, the sensitivity of the ZenoTOF 7600 system is comparable to that of the QTRAP 7500 instrument. The QTRAP 7500 had an approximately 10-fold better sensitivity for the three noted exceptions. Still, the LOQs for all compounds on both systems indicate that either instrument can measure bile acids in rat

Table 5. LOD and LOQ calculations for bile acids measured on the ZenoTOF 7600 system and the QTRAP 7500 instrument

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MKT-30747-B


Figure 7. XIC chromatogram of bile acids in human plasma. Eleven of 15 targeted bile acids were detected in a single rat plasma sample. In the figure, the major peaks are labeled, and the concentrations (μM) of each bile acid \pm the standard deviation are shown in the inset table.


Figure 8. Pearson correlation plot. A Pearson correlation calculation was performed on the quantitative plasma data acquired using the ZenoTOF 7600 system and the QTRAP 7500 instrument. Show here are plots, for multiple primary conjugated and unconjugated bile acids. All plots show a strong linear correlation between the two differently acquired data sets indicating that these two instruments are equally capable of measuring endogenous bile acids in plasma.


Figure 9. Product ion analysis of deoxycholic acid (DCA). A DCA standard (100 nM) was analyzed using the ZenoTOF 7600 system by CID- and EAD-based fragmentation (left panels). The middle panels show the .mol file for DCA, and the panels on the right show the theoretical fragments for CID (i.e., heterolytic fragmentation) and EAD (radical-based fragmentation). The highlighted fragment in blue corresponds to structure highlighted in bold on the middle figures and is indicated by a red arrow in the MS/MS spectra on the left.

plasma. **Figure 7** shows the XIC chromatogram for a single rat plasma extract. The concentrations of each targeted bile acid are displayed to the right of the peaks. Only the major peaks are labeled, but 11 of the 15 bile acids were detected *in vivo* using the ZenoTOF 7600 system. Similar data were obtained using the QTRAP 7500 instrument. **Figure 8** shows a Pearson correlation plot of the bile acid concentration data acquired from the plasma samples using the two different instrument types. The plots show a linear correlation between the two data sets for multiple primary conjugated and unconjugated bile acids in plasma. The *r* values of >0.99 and the *p* values <0.0001 indicate the data sets are equivalent and that both instruments are equally capable of measuring bile acids in plasma.

The EAD-based fragmentation functionality of the ZenoTOF 7600 system can provide significantly more structural details than that generated by CID-based fragmentation. **Figure 9** shows the product ion spectrum for DCA. This experiment analyzed a DCA standard (100 nM) using the ZenoTOF 7600 system by CID- and EAD-based fragmentation. The CID spectrum in the top left panel has relatively few fragments compared to the EAD-based spectrum in the bottom left panel. The middle panels show the .mol files for DCA and indicate the highlighted fragment in the right panels in bold. This feature of SCIEX OS, called fragment pane, enables the structural characterization of compounds by matching the actual MS/MS data to *in silico*-generated fragments based on the .mol file. The predicted CID fragments are shown in red in the top left panel. Note that most

of these fragments are not present in the spectrum. However, all the predicted fragments appear for DCA fragmented via EAD. This ability to generate fine structural detail was extended to determine whether EAD-based fragmentation could distinguish bile acid isomers.

Figure 2 shows that DCA, CDCA, and UDCA have the same empirical formula and exact mass but differ structurally in the position and stereochemistry of their hydroxyl functional groups. **Figure 1A** shows the TOFMS total ion chromatogram (TIC) for *m/z* 415.3 (DCA sodium adduct) and the importance of chromatography in resolving bile acid isomers. Using EAD-based fragmentation, unique fragment ions were identified for DCA and CDCA, the two isomers that elute closely (**Figure 1B**). These fragments are underlined in red on the spectra. Notably, the DCA-specific fragment at *m/z* 255.2107 shows no interference in the XIC for chromatographically resolved CDCA. Thus, the *m/z* fragment at 255.2017 Da can be used as a DCA-specific fragment for sMRM^{HR} analysis and may enable a shorter gradient for higher throughput. EAD-based fragmentation may enable the complete structural characterization of these molecules without requiring purification or other structural diagnostic tools. More importantly, this type of fragmentation may help identify novel conjugated bile acids (10-15). For example, Quinn et al. recently reported the discovery of several new MCBAs (specifically phenylalanochoic acid (Phe-CA), tyroschoic acid (Tyr-CA), and leucochoic acid (Leu-CA)) that were discovered in the intestinal contents collected from mice (11). Another study found evidence of

the amino acids glutamate, glutamine, aspartate, asparagine, methionine, histidine, lysine, serine, tryptophan, valine, alanine, and arginine forming conjugates with CDCA, DCA, or CA (12). Considering the emerging importance of MCBAAs in human disease (19-20) and the critical role of mass spectrometry in their characterization (21), the quantitative sensitivity and qualitative

capabilities of the ZenoTOF 7600 may prove invaluable to the growth of bile acid metabolism research.

Conclusions

- **The specific method details are provided for the analysis of endogenous bile acids using the ZenoTOF 7600 system and the QTRAP 7500 instrument**
- **The ZenoTOF 7600 system is capable of measuring endogenous bile acids with comparable sensitivity to the QTRAP 7500 system**
- **The narrow XIC window possible on HRMS systems enables the MS1-based collection of data without the inclusion of chemical background interferences**
- **EAD-based fragmentation provides structural characterization of molecules and may abrogate the need for extensive chromatographic separation of isomers to ensure specificity**

References

1. Liu D, Li G, Liu D, Shi W, Wang H, Zhang Q, Shen M, Huang X, Lin H. *Quantitative Detection of 15 Serum Bile Acid Metabolic Products by LC/MS/MS in the Diagnosis of Primary Biliary Cholangitis*. *Chem Biodivers*. 2023; Mar;20(3):e202200720. doi: 10.1002/cbdv.202200720
2. Gao T, Hu S, Xu W, Wang Z, Guo T, Chen F, Ma Y, Zhu L, Chen F, Wang X, Zhou J, Lv Z, Lu L. *Targeted LC-MS/MS profiling of bile acids reveals primary/secondary bile acid ratio as a novel biomarker for necrotizing enterocolitis*. *Anal Bioanal Chem*. 2024; Jan;416(1):287-297. doi: 10.1007/s00216-023-05017-7
3. Yang J, Lin J, Wang A, Yang X, Wang Y, Zhang Y, Dong H, Tian Y, Zhang Z, Wang M, Song R. *Study on the effect of calibration standards prepared with different matrices on the accuracy of bile acid quantification using LC-MS/MS*. *J Pharm Biomed Anal*. 2024; 237:115785. doi: 10.1016/j.jpba.2023.115785.
4. Engevik MA, Thapa S, Lillie IM, Yacyshyn MB, Yacyshyn B, Percy AJ, Chace D, Horvath TD. *Repurposing dried blood spot device technology to examine bile acid profiles in human dried fecal spot samples*. *Am J Physiol Gastrointest Liver Physiol*. 2024; 326:G95-G106. Doi: 10.1152/ajpgi.00188.2023
5. Baba T, Campbell JL, Le Blanc JCY, Baker PRS, Hager JW, Thomson BA. *Development of a Branched Radio-Frequency Ion Trap for Electron Based Dissociation and Related Applications*. *Mass Spectrom (Tokyo)*. 2017;6(1):A0058. doi: 10.5702/massspectrometry.A0058.
6. Che P, Davidson JT, Kool J, Kohler I. *Electron activated dissociation - a complementary fragmentation technique to collision-induced dissociation for metabolite identification of synthetic cathinone positional isomers*. *Anal Chim Acta*. 2023; 1283:341962. doi: 10.1016/j.aca.2023.341962
7. Baba T, Campbell JL, Le Blanc JCY, Baker PRS. *Distinguishing Cis and Trans Isomers in Intact Complex Lipids Using Electron Impact Excitation of Ions from Organics Mass Spectrometry*. *Anal Chem*. 2017 Jul 18;89(14):7307-7315. doi: 10.1021/acs.analchem.6b04734. Epub 2017 Jun 28. PMID: 28613874.
8. Russell DW. *The enzymes, regulation, and genetics of bile acid synthesis*. *Annu Rev Biochem*. 2003; 72:137-74. doi: 10.1146/annurev.biochem.72.121801.161712
9. Ridlon JM, Kang DJ, Hylemon PB. *Bile salt biotransformations by human intestinal bacteria*. *J Lipid Res*. 2006 ;47(2):241-59. doi: 10.1194/jlr.R500013-JLR200
10. Guzior DV, Quinn RA. *Review: microbial transformations of human bile acids*. *Microbiome*. 2021; 9(1):140. doi: 10.1186/s40168-021-01101-1
11. Quinn RA, Melnik AV, Vrbanc A, Fu T, Patras KA, Christy MP, Bodai Z, Belda-Ferre P, Tripathi A, Chung LK, Downes M, Welch RD, Quinn M, Humphrey G, Panitchpakdi M, Weldon KC, Aksenov A, da Silva R, Avila-Pacheco J, Clish C, Bae S, Mallick H, Franzosa EA, Lloyd-Price J, Bussell R, Thron T, Nelson AT, Wang M, Leszczynski E, Vargas F, Gauglitz JM, Meehan MJ, Gentry E, Arthur TD, Komor AC, Poulsen O, Boland BS, Chang JT, Sandborn WJ, Lim M, Garg N, Lumeng JC, Xavier RJ, Kazmierczak BI, Jain R, Egan M, Rhee KE, Ferguson D, Raffatellu M, Vlamakis H, Haddad GG, Siegel D, Huttenhower C, Mazmanian SK, Evans RM, Nizet V, Knight R, Dorrestein PC. *Global chemical effects of*

- the microbiome include new bile-acid conjugations.* *Nature.* 2020; 579(7797):123-129. doi: 10.1038/s41586-020-2047-9
12. Lucas LN, Barrett K, Kerby RL, Zhang Q, Cattaneo LE, Stevenson D, Rey FE, Amador-Noguez D. *Dominant Bacterial Phyla from the Human Gut Show Widespread Ability To Transform and Conjugate Bile Acids.* *mSystems.* 2021; e0080521. doi: 10.1128/mSystems.00805-21
 13. Foley MH, Walker ME, Stewart AK, O'Flaherty S, Gentry EC, Patel S, Beaty VV, Allen G, Pan M, Simpson JB, Perkins C, Vanhoy ME, Dougherty MK, McGill SK, Gulati AS, Dorrestein PC, Baker ES, Redinbo MR, Barrangou R, Theriot CM. *Bile salt hydrolases shape the bile acid landscape and restrict Clostridioides difficile growth in the murine gut.* *Nat Microbiol.* 2023; 8(4):611-628. doi: 10.1038/s41564-023-01337-7
 14. Gentry EC, Collins SL, Panitchpakdi M, Belda-Ferre P, Stewart AK, Carrillo Terrazas M, Lu HH, Zuffa S, Yan T, Avila-Pacheco J, Plichta DR, Aron AT, Wang M, Jarmusch AK, Hao F, Syrkin-Nikolau M, Vlamakis H, Ananthakrishnan AN, Boland BS, Hemperly A, Vande Castele N, Gonzalez FJ, Clish CB, Xavier RJ, Chu H, Baker ES, Patterson AD, Knight R, Siegel D, Dorrestein PC. *Reverse metabolomics for the discovery of chemical structures from humans.* *Nature.* 2024; 626(7998):419-426. doi: 10.1038/s41586-023-06906-8
 15. Guzior DV, Okros M, Shivel M, Armwald B, Bridges C, Fu Y, Martin C, Schillmiller AL, Miller WM, Ziegler KM, Sims MD, Maddens ME, Graham SF, Hausinger RP, Quinn RA. *Bile salt hydrolase acyltransferase activity expands bile acid diversity.* *Nature.* 2024; Online Ahead of Print. doi: 10.1038/s41586-024-07017-8.
 16. Horvath TD, Haidacher SJ, Hoch KM, Auchtung JM, Haag AM. *A high-throughput LC-MS/MS method for the measurement of the bile acid/salt content in microbiome-derived sample sets.* *MethodsX.* 2020; 7:100951. doi: 10.1016/j.mex.2020.100951
 17. *ICH harmonized tripartite guideline: validation of analytical procedures: text and methodology Q2(R1)*, International Conference of harmonization of technical requirements for registration of pharmaceuticals for human use, 2005
 18. Baker, PRS and Proos R. *Untargeted data-dependent acquisition (DDA) metabolomics analysis using the ZenoTOF 7600 system.* SCIEX technical note, [RUO MKT-02-15367-A](#)
 19. McDonald D, Hyde E, Debelius JW, Morton JT, Gonzalez A, Ackermann G, Aksenov AA, Behsaz B, Brennan C, Chen Y, DeRight Goldasich L, Dorrestein PC, Dunn RR, Fahimipour AK, Gaffney J, Gilbert JA, Gogul G, Green JL, Hugenholtz P, Humphrey G, Huttenhower C, Jackson MA, Janssen S, Jeste DV, Jiang L, Kelley ST, Knights D, Kosciulek T, Ladau J, Leach J, Marotz C, Meleshko D, Melnik AV, Metcalf JL, Mohimani H, Montassier E, Navas-Molina J, Nguyen TT, Peddada S, Pevzner P, Pollard KS, Rahnavaard G, Robbins-Pianka A, Sangwan N, Shorestein J, Smarr L, Song SJ, Spector T, Swafford AD, Thackray VG, Thompson LR, Tripathi A, Vázquez-Baeza Y, Vrbancac A, Wischmeyer P, Wolfe E, Zhu Q; *The American Gut Consortium; Knight R. American Gut: an Open Platform for Citizen Science Microbiome Research.* *mSystems.* 2018; 3(3):e00031-18. doi: 10.1128/mSystems.00031-18
 20. Lloyd-Price J, Arze C, Ananthakrishnan AN, Schirmer M, Avila-Pacheco J, Poon TW, Andrews E, Ajami NJ, Bonham KS, Brislawn CJ, Casero D, Courtney H, Gonzalez A, Graeber TG, Hall AB, Lake K, Landers CJ, Mallick H, Plichta DR, Prasad M, Rahnavaard G, Sauk J, Shungin D, Vázquez-Baeza Y, White RA 3rd; IBDMDB Investigators; Braun J, Denson LA, Jansson JK, Knight R, Kugathasan S, McGovern DPB, Petrosino JF, Stappenbeck TS, Winter HS, Clish CB, Franzosa EA, Vlamakis H, Xavier RJ, Huttenhower C. *Multi-omics of the gut microbial ecosystem in inflammatory bowel diseases.* *Nature.* 2019; 569(7758):655-662. doi: 10.1038/s41586-019-1237-9
 21. Wang M, Jarmusch AK, Vargas F, Aksenov AA, Gauglitz JM, Weldon K, Petras D, da Silva R, Quinn R, Melnik AV, van der Hooft JJJ, Caraballo-Rodríguez AM, Nothias LF, Aceves CM, Panitchpakdi M, Brown E, Di Ottavio F, Sikora N, Elijah EO, Labarta-Bajo L, Gentry EC, Shalappour S, Kyle KE, Puckett SP, Watrous JD, Carpenter CS, Bouslimani A, Ernst M, Swafford AD, Zúñiga EI, Balunas MJ, Klassen JL, Loomba R, Knight R, Bandeira N, Dorrestein PC. *Mass spectrometry searches using MASST.* *Nat Biotechnol.* 2020; 38(1):23-26. doi: 10.1038/s41587-019-0375-9

Quantitative and qualitative bile acid analysis on the ZenoTOF 8600 system using EAD

Paul RS Baker¹, Remco van Soest¹, Santosh Kapil¹, Robert Proos¹, Maxim D Seferovic² and Thomas D. Horvath^{3,4}

¹SCIEX, USA; ²Dept of Obstetrics and Gynecology, Baylor College of Medicine, Houston, TX; ³Dept of Pathology, Texas Children's Hospital, Houston, TX; ⁴Dept of Pathology & Immunology, Baylor College of Medicine, Houston, TX

This technical note demonstrates the power and capability of the ZenoTOF 8600 system to sensitively detect, quantify, and structurally characterize the bile acid content in human plasma.

The analysis of bile acids in human samples has become increasingly important as insights into the prominent biological role of this class of molecules are discovered [1,2]. Analysis of bile acids by triple-quadrupole MS (TQMS) systems is challenging because of the high chemical background found in several precursor ions to precursor ion-based multiple-reaction monitoring (MRM) transitions used in current state-of-the-art assays [1-4]. High-resolution mass spectrometry (HRMS) enables the extraction of fragment ions with a narrow mass-to-charge (m/z) window, which can reduce background chemical interferences and improve the signal-to-noise (S/N) of the assay as previously demonstrated [3].

The detection of individual bile acid isomers currently depends on chromatographic resolution because collision-induced dissociation (CID)-based fragmentation cannot distinguish these isomeric

metabolites. However, a novel fragmentation method, electron-activated dissociation (EAD) can generate unique fragments, even among isomers that may help reduce the time required for analysis.

Here, the ZenoTOF 8600 system was used to quantify the bile acid content of human plasma sample extracts using both CID- and EAD-based fragmentation methods. EAD-based fragmentation was strategically used when isomeric pairs overlapped to provide specificity and enable faster analysis times.

Key features of bile acid analysis on the ZenoTOF 8600 system

- The ZenoTOF 8600 system can detect and quantify bile acids with an approximately 12-fold better sensitivity than the ZenoTOF 7600 system
- The narrow fragment ion extraction window possible with HRMS improves assay sensitivity by reducing background chemical interferences and, consequently, increases the S/N of the assay
- EAD can generate diagnostic fragment ions to enable the distinction between bile acid isomers, which improves the overall specificity and speed of the analysis


Figure 1. Detection and quantitation of bile acid isomers in human plasma using EAD. Using a fast, 10-minute gradient, bile acids were analyzed on the ZenoTOF 8600 system. Isomers such as GCDCA and GDCA are indistinguishable using CID-based fragmentation. Using EAD, unique, diagnostic fragments were generated, which enabled their accurate quantitation despite the chromatographic overlap. The upper panels show the two isomers overlapping when using a CID-based fragmentation method; the lower panels show the ability to isolate each isomer separately when a unique, EAD-based fragment is used.

RUO-MKT-35189-A

Introduction

Bile acids are amphipathic steroidal molecules synthesized in the liver from cholesterol and play a crucial role in the digestion and absorption of dietary lipids [1]. Primary bile acids, such as cholic acid and chenodeoxycholic acid, are conjugated with glycine or taurine before being secreted into bile. Once released into the intestinal lumen, these acids undergo microbial transformation to form secondary bile acids like deoxycholic acid and lithocholic acid that are collectively termed microbially conjugated bile acids (MCBAs). Enterohepatic circulation maintains a dynamic balance of bile acid synthesis, secretion, reabsorption, and excretion. Beyond their digestive functions, bile acids are now recognized as signaling molecules involved in metabolic regulation and immune responses, acting through nuclear receptors such as FXR and membrane-bound G-protein coupled receptors [2].

Mass spectrometry (MS) has become the analytical method of choice for the quantitation and profiling of bile acids due to its sensitivity, specificity, and ability to differentiate structurally similar compounds. [3-5] The use of HPLC ESI-MS/MS allows for the simultaneous detection of multiple bile acid species, including conjugated and unconjugated forms, in biological matrices such as plasma, bile, and feces. Targeted MS approaches using multiple reaction monitoring (MRM) provide high-throughput and quantitative data essential for clinical and research applications in liver disease, metabolic disorders, and gut microbiome studies.

Recently, we reported the analysis of bile acids using the ZenoTOF 7600 system, which is a hybrid quadrupole time of flight instrument (QTOF) that measures compounds with high resolution (~35K) and mass accuracy (<2 ppm) [TN ref]. In the negative ion mode, bile acid analysis by nominal mass instruments is challenging due to solvent-based high chemical background noise. With accurate mass systems, this issue can be mitigated using narrow product ion extraction windows (XIC), which minimize noise and enable a similar calculated LOQ for the ZenoTOF 7600 system as with the highly sensitive triple quadrupole SCIEX 7500 system. Using the ZenoTOF 8600 system, these studies were extended to capitalize on the instrument's improved sensitivity.

The ZenoTOF 8600 system has unique design features that make it particularly suited for targeted bile acid analysis. First, the instrument has a complementary fragmentation mode, EAD, that provides unparalleled structural characterization of small molecules [6-10]. Second, improved hardware, such as the front end with the Optiflow Pro source, a larger instrument orifice, and an advanced optical detector, enables higher sensitivity compared to the ZenoTOF 7600 system. Third, the instrument's front end is equipped with a DJET+ with Mass Guard, which has been shown to reduce instrument contamination on the SCIX 7500+ system [11], an important consideration for lipid and metabolite analysis, especially in the


Figure 2. The glucodeoxycholic acid subclass of bile acids. Bile acids are cholesterol-derived amphipathic molecules of saturated hydroxylated C-24 sterols. GDCA has a hydroxyl moiety at carbon 12 [C-12], whereas GUDCA and GCDCA are hydroxylated at carbon 7 [C-7]. The latter two isomers are characterized by a different stereochemistry at the number 7 carbon. Quantitative specificity for this class of molecules is highly dependent on chromatographic resolution when analyzed using CID-based fragmentation.

complex matrices common to bile acid analysis, including plasma, feces, and bile.

Bile acids are cholesterol-derived amphipathic molecules of saturated hydroxylated C-24 sterols. As shown in **Figure 2**, the location and stereochemistry of the hydroxyl functional groups define the different isomers within each sub-class. CID-based fragmentation of bile acids primarily yields product ions related to the bile acid head group, and relatively few fragments are generated from the ring structure; hence, CID alone does not generate diagnostically useful fragment ions for identification or quantitation purposes. Using chromatography to resolve bile acid isomers requires relatively long chromatographic gradients, which are not conducive to high-throughput analysis. An alternative to lengthy chromatographic gradients and CID-based fragmentation is EAD. EAD-based fragmentation provides ample fragments derived from the sterol ring structure, and it can distinguish the isomers of each bile acid subclass without extensive chromatographic method development.

In this technical note, the speed and sensitivity of the ZenoTOF 8600 system were leveraged to quantify bile acids in prepared plasma samples. To demonstrate the power of EAD to improve structural specificity during quantitative analysis, unique EAD-based fragments were identified for potentially co-eluting isomer pairs and used strategically in conjunction with CID-based fragmentation to enable a

faster assay while maintaining quantitative specificity, accuracy, and precision.

Methods

Materials: Bile acid standard mixtures were purchased from Cambridge Isotopes Laboratories. The stable isotope-labeled bile acid mixes (unconjugated, cat# MSK-BA1; conjugated, MSK-BA2) were used as internal standards, and unlabeled bile acid mixes (unconjugated, cat# MSK-BA1-US; conjugated, MSK-BA2-US) were used as primary reference standards. All solvents were of LC-MS grade and obtained from Burdick and Jackson and Fisher Scientific.

Sample preparation: Human blood (research use only) was collected into tubes containing potassium ethylenediamine tetraacetic acid (K₂EDTA) anti-coagulant following the procedure approved by the institutional review board (IRB) of BCM. The plasma was isolated from the red blood cell fraction by centrifugation. A 50 µL volume of each plasma sample was pipetted into individual 0.6 mL Eppendorf tubes, and a 200 µL volume of methanol was added to precipitate plasma proteins (primary DF=5-fold). All samples were vortex-mixed for 30 seconds and centrifuged at 17,000g for 5 minutes, and the supernatant for each clarified plasma extract sample was transferred into fresh vials. A 50 µL volume of each clarified plasma extract sample was transferred into fresh glass vials. A 450 µL volume of an internal standard solution (10 nM deuterated BAs) was added to each sample vial (DF=50-fold overall), and the vials were capped and vortex-mixed for 30 seconds. And prepared for analysis by HPLC ESI-MS/MS.

Chromatography: Extracted metabolites were resolved using a Nexera 40 Series UHPLC system equipped with a Restek Raptor C18 column (100 x 2.1 mm; 5.0 µm particle size). The autosampler sample bay was maintained at 10 °C, and a 5 µL sample volume was injected on column. The column oven temperature was kept at 50°C with a constant mobile phase flow rate of 0.3 mL/min for a total run time of 17 min. The mobile phase compositions were (A) 10 mM ammonium formate in water and (B) pure acetonitrile. Gradient details are shown in **Table 1**. Alternatively, the gradient was altered to enable a 10 min total run time (**Table 2**) in which EAD-based fragmentation was selectively employed to maintain compound specificity for co-eluting

Table 1. Chromatographic gradient conditions

Time [min]	Mobile phase A [%]	Mobile Phase B [%]
0	75	25
1.5	75	25
13	20	80
14.5	20	80
14.6	75	25
17	75	25

Table 2. 10-minute chromatographic gradient conditions

Time [min]	Mobile phase A [%]	Mobile Phase B [%]
0	75	25
1.5	75	25
6	20	80
7.5	20	80
7.6	75	25
10	75	25

isomers that are not distinguishable using CID-based fragmentation, an example of which is shown in **Figure 1**.

Mass spectrometry: Sample extracts were analyzed using the ZenoTOF 8600 system with an OptiFlow ProTurbo V ion source and a scheduled, high-resolution multiple reaction monitoring (sMRMHR) scan mode. CID-based fragmentation was used for the initial assessment of the system's performance using the 17-minute gradient. This experiment was repeated using EAD-based fragmentation. To speed up the experiment and reduce the total time to 10 min, an sMRMHR experiment was performed using mixed-mode fragmentation. For isomers that overlapped when the gradient was shortened, an EAD-based fragment that is diagnostic for the overlapping isomers was used. Otherwise, CID was used for analysis due to its greater sensitivity and shorter accumulation times. MS instrument parameters are presented in **Table 3**, and a comprehensive list of the compounds studied is shown in **Table 4**.

Table 3. Instrument parameter settings

Parameter	ZenoTOF 8600 system
Curtain gas (CUR)	35
Ion source gas 1 (GS1)	60
Ion source gas 2 (GS2)	80
CAD Gas (CAD)	7
Source temperature (TEM)	250 °C
Ion spray voltage (IS)	-4500 V/+5000V (CID/EAD)
Declustering Potential (DP)	-80 V/+80V (CID/EAD)
CID Accumulation/dwell time	Variable (sMRM ^{HR})
CID Collision energy (CE)	Molecule- dependent
Q1/Q3 Mass Resolution	N/A
EAD Accumulation/dwell time	40 ms
CE for TOF MS	-10 V/+10V (CID/EAD)
TOF MS mass range	70-1000 Da
TOF MS/MS mass range	40-400 Da
Time bins to sum	6
EAD Collision energy (CE)	12 V
Electron kinetic energy (KE)	13 eV
Zeno pulsing for MS/MS	Yes
Zeno threshold	20000

RUO-MKT-35189-A

Table 4. Bile acid compound parameters

Bile Acid	Abbreviation	[M - H] ⁻	[M - H] ⁻ Product ion	[M + Na] ⁺	Collision energy [V]
Lithocholic Acid	LCA	375.29	375.29	399.288	-5
D4-Lithocholic Acid	D4-LCA	379.316	379.316	403.313	-5
Ursodeoxycholic Acid	UDCA	391.285	391.285	415.283	-5
D4-Ursodeoxycholic Acid	D4-UDCA	395.31	395.31	419.308	-5
Chenodeoxycholic Acid	CDCA	391.285	391.285	415.283	-5
D4-Chenodeoxycholic Acid	D4-CDCA	395.31	395.31	419.308	-5
Deoxycholic Acid	DCA	391.285	391.285	415.283	-49
D4-Deoxycholic Acid	D4-DCA	395.31	395.31	419.308	-49
Cholic Acid	CA	407.28	343.266	431.278	-49
D4-Cholic Acid	D4-CA	411.305	411.305	435.303	-49
Glycolithocholic Acid	GLCA	432.312	388.324	456.309	-48
D4-Glycolithocholic Acid	D4-GLCA	436.337	392.349	460.335	-48
Glycoursodeoxycholic Acid	GUDCA	448.307	74.025	472.304	-82
D4-Glycoursodeoxycholic Acid	D4-GUDCA	452.332	74.025	476.329	-82
Glycochenodeoxycholic Acid	GCDCA	448.307	74.025	472.304	-82
D4-Glycochenodeoxycholic Acid	D4-GCDCA	452.332	74.025	476.329	-82
Glycodeoxycholic Acid	GDCA	448.307	74.025	472.304	-82
D4-Glycodeoxycholic Acid	D4-GDCA	452.332	74.025	476.329	-82
Glycocholic Acid	GCA	464.302	402.302	488.299	-48
D4-Glycocholic Acid	D4-GCA	468.327	406.331	492.324	-48
Taurolithocholic Acid	TLCA	482.295	79.958	506.292	-137
D4-Taurolithocholic Acid	D4-TLCA	486.32	79.958	510.317	-137
Tauroursodeoxycholic Acid	TUDCA	498.289	79.958	522.287	-142
D4-Tauroursodeoxycholic Acid	D4-TUDCA	502.315	79.958	526.312	-142
Taurochenodeoxycholic Acid	TCDC	498.289	79.958	522.287	-142
D4-Taurochenodeoxycholic Acid	D4-TCDC	502.315	79.958	526.312	-142
Taurodeoxycholic Acid	TDCA	498.289	79.958	522.287	-142
D4-Taurodeoxycholic Acid	D4-TDCA	502.315	79.958	526.312	-142
Taurocholic Acid	TCA	514.284	124.008	538.282	-66
D4-Taurocholic Acid	D4-TCA	518.31	124.007	542.307	-66

During sMRMhr analysis, a full product ion spectrum is acquired. To identify diagnostic, EAD-based fragments for co-eluting bile acid isomers such as GCDCA and GDCA,

The Analytics module within SCIEX OS software was used for quantitation, and the Explore module was used to structurally characterize bile acids from EAD-based MS/MS spectra. The limit of quantitation (LOQ) was determined from the lowest standard injected with the %CV value < 15 and an accuracy within 85% and 115% of the expected values for the internal standard curve.

Results and Discussion

A bile acid standard curve was generated neat in solvent and analyzed on the ZenoTOF 8600 system. The parameter settings were optimized, as detailed in **Tables 3 and 4**, and the chromatographic gradient conditions resulted in the separation of each isomeric bile acid with

baseline resolution (**Figure 3**). The concentrations within the internal standard curve ranged from 0.011 to 100 nM, covering a dynamic range of > 4 orders of magnitude. An example curve acquired for glycoursodeoxycholic acid (GUDCA) is presented in **Figure 4 (panel A)**, and the individual chromatograms for each injection of the standard curve from 0.011 to 50 nM GUDCA are shown (**panel B**). The associated quantitative statistics for GUDCA quantitation are shown in **Table 5**. The data derived from the internal standard curves were used to calculate the LOQ of each bile acid measured with the ZenoTOF 8600 system (**Table 6**). These values were determined as the lowest concentration standard injected that gave a %CV value <15% and an accuracy within the range of 85% to 115% of the expected value.

The LOQ values for the ZenoTOF 7600 system were calculated in a different way [BA TN] than they were in this report [see above]. Additionally, on the ZenoTOF 7600 system, a 15 µL sample was injected, whereas only 5 µL was injected onto the ZenoTOF 8600

RUO-MKT-35189-A


Figure 3. Elution profile of targeted bile acids using CID-based fragmentation in the negative ion mode. The 50 nM standard was injected and analyzed using the sMRM scan mode. Total run time was 17 min.


Figure 4. Example internal standard curve with representative peaks for GUDCA. **[A]** Internal standard curve for GUDCA from 0.011 to 50 nM. **[B]** Example chromatograms at each concentration level of the curve; blank region highlighted in grey.

Table 5. Quantitative statistics for GUDCA

Standard	Standard concentration (nM)	Calculated concentration (nM)	Std dev	%CV	Average accuracy (n = 5)
GUDCA	0.011	0.011	0.0014	12.9	100.4
GUDCA	0.033	0.033	0.0017	5.13	101.0
GUDCA	0.098	0.093	0.0074	7.90	95.0
GUDCA	0.391	0.385	0.0087	2.27	98.5
GUDCA	1.56	1.57	0.0873	5.58	100.3
GUDCA	6.25	6.39	0.1008	1.58	102.3
GUDCA	25	26	0.6598	2.59	102.0
GUDCA	50	50	1.3860	2.76	100.5

RUO-MKT-35189-A

system. To compare the relative sensitivity of the two different systems, the on-column injection amount at the LOQ was used as a reference (**Table 7**). The data show that the ZenoTOF 8600 system is an average of ~20-fold more sensitive than the ZenoTOF 7600 system [range = 1.1 to 61-fold improvement]. The wide range may be, in part, due to the difference in the sources used with the two instruments. The ZenoTOF 7600 system is equipped with the Turbo V twin spray source, which optimizes differently and may have a different ionization capacity, or a differential ionization efficiency compared to the Optiflow Pro Turbo V ion source on the ZenoTOF 8600 system. The ion path of the ZenoTOF 8600 system is also different from the previous version, which may also slightly affect transmission in a compound-dependent manner.

The assay developed using neat standards was used to measure endogenous bile acids in human plasma. **Figure 5** shows the combined XIC data for all targeted bile acids measured with the sMRM workflow in human plasma extract. Except for lithocholic acid (LCA), all bile acids were detected and quantified using the ZenoTOF 8600 system. The data shown here were collected using the CID-based fragmentation mode, and the LC gradient extended over 17 minutes. The concentrations of bile acids detected in two representative plasma samples are shown in **Table 8**. The same experiment was performed using EAD-based fragmentation for all compounds (not shown), and those data were used to identify diagnostic fragment ions for bile acid isomers that co-elute using the faster, 10 min gradient.

Table 6. LOQ values for bile acids on the ZenoTOF 8600 system

Bile acid	LOQ (nM)
LCA	0.098
UDCA	0.098
CDCA	0.033
DCA	0.033
CA	0.033
GLCA	0.033
GUDCA	0.011
GCDCA	0.011
GDC	0.011
GCA	0.098
TLCA	0.033
TUDCA	0.011
TCDC	0.033
TDCA	0.033
TCA	0.033

Table 7. Comparison of on-column injection load at LOD between the ZenoTOF 7600 and 8600 systems

Bile Acid	fmol on-column injection at LOQ		Fold-increase in sensitivity
	ZenoTOF 7600 system	ZenoTOF 8600 system	
LCA	6.075	0.490	12.4
UDCA	2.160	0.490	4.4
CDCA	0.446	0.165	2.7
DCA	0.188	0.165	1.1
CA	9.195	0.165	55.7
GLCA	10.035	0.165	60.8
GUDCA	2.025	0.055	36.8
GCDCA	1.277	0.055	23.2
GDC	0.777	0.055	14.1
GCA	3.150	0.490	6.4
TLCA	3.165	0.165	19.2
TUDCA	1.115	0.055	20.3
TCDC	4.635	0.165	28.1
TDCA	1.515	0.165	9.2
TCA	3.810	0.165	23.1


Figure 5. Bile acid detection and quantitation in human plasma using sMRMhr analysis on the ZenoTOF 8600 system

Table 8. Example quantitative results from 2 different human plasma samples

Bile Acid	Abbreviation	Sample 1		Sample 2	
		Concentration (nM)	%CV (n=3)	Concentration (nM)	%CV (n=3)
Lithocholic Acid	LCA	ND	N/A	ND	N/A
Ursodeoxycholic Acid	UDCA	0.714	4.9	0.490	4.3
Chenodeoxycholic Acid	CDCA	12.4	1.7	0.612	1.1
Deoxycholic Acid	DCA	12.3	2.4	0.708	6.0
Cholic Acid	CA	6.03	2.4	5.32	7.2
Glycolithocholic Acid	GLCA	0.395	6.0	0.041	14
Glycoursodeoxycholic Acid	GUDCA	0.181	1.0	2.48	1.0
Glycochenodeoxycholic Acid	GCDCA	797	1.4	112	3.7
Glycodeoxycholic Acid	GDC	146	3.8	0.863	6.1
Glycocholic Acid	GA	1570	2.7	206	0.9
Taurolithocholic Acid	TLCA	1.13	0.89	0.093	3.7
Tauroursodeoxycholic Acid	TUDCA	2.41	4.3	1.53	4.1
Taurochenodeoxycholic Acid	TCDC	1510	0.37	344	6.0
Taurodeoxycholic Acid	TDCA	246	2.5	254	0.30
Taurocholic Acid	TCA	3780	1.2	313	3.7

The EAD-based fragmentation functionality of the ZenoTOF 7600 system can provide significantly more structural details than that

generated by CID-based fragmentation. Previously, we reported the identification of unique fragments for the bile acid isomers CDCA and

RUO-MKT-35189-A


Figure 6. Comparison of EAD-based fragmentation of GDCA and GCDCA. A portion of the product ion spectra for GCDCA and GDCA is highlighted, with the spectrum for GDCA overlaid and inverted. Multiple fragments were identified for each bile acid isomer that appear to be unique and could serve as a selective fragment ion during analysis (red arrows). The fragment ions with m/z 262.1403 and 316.1888 were chosen for GDCA and GCDCA, respectively, due to their relatively high abundances.


Figure 7. Resolution of bile acid isomers using EAD-based fragmentation. GCDCA and GDCA are bile acid isomers that have identical CID-based MS/MS spectra. Using EAD, unique, diagnostic fragments were observed for each isomer (A and B) in a standard mixture (6.56 nM). Using the fragments pane in SCIEX OS, structures for each fragment, m/z 316.1886 and 262.1403, were proposed for GUDCA and GDCA, respectively (C). Using a 10 mDa XIC window, panel D shows the improved specificity of the analysis for the two standards.

DCA [TN BA]. These isomers elute very closely and are responsible, in part, for the relatively long analysis time so that they remain chromatographically resolved. It was hypothesized that using EAD-

based fragmentation for these isomer pairs would allow for a truncated chromatographic gradient and enable faster sample analysis. Here, we extended these findings using the ZenoTOF 8600

RUO-MKT-35189-A

system and leveraged the ~10-fold increase in sensitivity to strategically use EAD-based fragmentation for specific bile acids during sMRMhr analysis. **Figure 6** shows a region of the EAD-derived product ion spectrum for GCDCA and GDCA acquired as neat standards; the spectrum for GCDCA is inverted so the differences in the spectra are more apparent. Fragment ions at m/z 262.1403 and 316.1888 were identified as unique, diagnostic ions for GDCA and GCDCA, respectively. Other apparently unique fragment ions for each bile acid isomer are indicated with red arrows. Having multiple unique fragment ions may help with quantitation, as fragment intensities can be added during quantitative processing in the Analytics module of SCIEX OS software.

Using the faster, 10-min gradient, the bile acid standards were analyzed using EAD-based fragmentation. **Figure 7, panel A**, shows the elution profile of GCDCA and GDCA using neat standards, with the peak-to-peak retention time difference of ~12 s. The product ion spectra are shown in **panel B**, with the diagnostic ion for each apparent in its respective spectrum. The fragment panel app within the SCIEX OS Explore module proposed 2 possible structures of the fragments, as shown in **Panel C**. In **Panel D**, the XICs for GCDCA and GDCA show minimal to no interference when the EAD-specific fragments are used.

To test whether this technique can be used to measure bile acids in plasma, the same 10-minute, EAD-based method used for **Figure 7** was applied to the analysis of human plasma extracts (**Figure 1**). In contrast to the results with neat standards, in the plasma matrix, the isomers GCDCA and GDCA are not baseline resolved when using a common fragment—in this case, the precursor-to-precursor MRMhr transition. However, when the data were processed using the unique EAD-based fragments, the individual bile acids were completely resolved, allowing for accurate quantitation.

In the data reported here, the ZenoTOF 8600 system was able to measure bile acids in human plasma with ~10-fold greater sensitivity than the ZenoTOF 7600 system using CID-based fragmentation. Of note, for this particular assay, a choice has to be made whether to maximize sensitivity and run the samples using the longer gradient in the negative ion mode with CID-based fragmentation, or to maximize speed and specificity by running the assay with the fast gradient using EAD-based fragmentation and sacrifice ~100-fold sensitivity compared to the longer assay. For the isomers in question here, GDCA, GCDCA, and GUDCA, the EAD-based method was appropriate for biological analysis. Currently, a study is underway to measure amino acid conjugates of bile acids. These molecules are measured in the positive ion mode and are amenable to a combined analysis of CID- and EAD-based analysis, using EAD only when needed to improve specificity.

RUO-MKT-35189-A

Conclusions

- The ZenoTOF 8600 system is ~20-fold more sensitive for the analysis of bile acids compared to the ZenoTOF 7600 system narrow XIC window possible on HRMS systems enables the MS1-based collection of data without the inclusion of chemical background interferences
- EAD-based fragmentation provides structural characterization of molecules and may abrogate the need for extensive chromatographic separation of isomers to ensure specificity

References

1. Collins, S. L., Stine, J. G., Bisanz, J. E., Okafor, C. D., & Patterson, A. D. [2023]. Bile acids and the gut microbiota: metabolic interactions and impacts on disease. *Nature reviews. Microbiology*, 21(4), 236–247. <https://doi.org/10.1038/s41579-022-00805-x>
2. Lin, S., Wang, S., Wang, P., Tang, C., Wang, Z., Chen, L., Luo, G., Chen, H., Liu, Y., Feng, B., Wu, D., Burrin, D. G., & Fang, Z. [2023]. Bile acids and their receptors in regulation of gut health and diseases. *Progress in lipid research*, 89, 101210. <https://doi.org/10.1016/j.plipres.2022.101210>
3. Baker, PRS., Proos, R., Sefevoric, MD., and Horvath, TD. Quantitative analysis and structural characterization of bile acids using the ZenoTOF 7600 system. SCIEX technical note MKT-30747-B. <https://sciex.com/tech-notes/life-science-research/metabolomics/quantitative-analysis-and-structural-characterization-of-bile-ac>
4. Liu, D., Li, G., Liu, D., Shi, W., Wang, H., Zhang, Q., Shen, M., Huang, X., & Lin, H. [2023]. Quantitative Detection of 15 Serum Bile Acid Metabolic Products by LC/MS/MS in the Diagnosis of Primary Biliary Cholangitis. *Chemistry & biodiversity*, 20(3), e202200720. <https://doi.org/10.1002/cbdv.202200720>
5. Gao, T., Hu, S., Xu, W., Wang, Z., Guo, T., Chen, F., Ma, Y., Zhu, L., Chen, F., Wang, X., Zhou, J., Lv, Z., & Lu, L. [2024]. Targeted LC-MS/MS profiling of bile acids reveals primary/secondary bile acid ratio as a novel biomarker for necrotizing enterocolitis. *Analytical and bioanalytical chemistry*, 416(1), 287–297. <https://doi.org/10.1007/s00216-023-05017-7>
6. Yang, J., Lin, J., Wang, A., Yang, X., Wang, Y., Zhang, Y., Dong, H., Tian, Y., Zhang, Z., Wang, M., & Song, R. [2024]. Study on the effect of calibration standards prepared with different matrix on the accuracy of bile acid quantification using LC-MS/MS. *Journal of*

pharmaceutical and biomedical analysis, 237, 115785.

<https://doi.org/10.1016/j.jpba.2023.115785>

7. Baba T, Campbell JL, Le Blanc JCY, Baker PRS, Hager JW, Thomson BA. *Development of a Branched Radio-Frequency Ion Trap for Electron Based Dissociation and Related Applications*. *Mass Spectrom [Tokyo]*. 2017;6(1):A0058. doi: 10.5702/massspectrometry.A0058
8. Baba, T., Campbell, J. L., Le Blanc, J. C., & Baker, P. R. [2016]. In-depth sphingomyelin characterization using electron impact excitation of ions from organics and mass spectrometry. *Journal of lipid research*, 57(5), 858–867. <https://doi.org/10.1194/jlr.M067199>
9. Baba, T., Campbell, J. L., Le Blanc, J. C., & Baker, P. R. [2016]. Structural identification of triacylglycerol isomers using electron impact excitation of ions from organics [EIEIO]. *Journal of lipid research*, 57(11), 2015–2027. <https://doi.org/10.1194/jlr.M070177>
10. Baba, T., Campbell, J. L., Le Blanc, J. C. Y., Baker, P. R. S., & Ikeda, K. [2018]. Quantitative structural multiclass lipidomics using differential mobility: electron impact excitation of ions from organics [EIEIO] mass spectrometry. *Journal of lipid research*, 59(5), 910–919. <https://doi.org/10.1194/jlr.D083261>
11. Baba T, Campbell JL, Le Blanc JCY, Baker PRS. Distinguishing Cis and Trans Isomers in Intact Complex Lipids Using Electron Impact Excitation of Ions from Organics Mass Spectrometry. *Anal Chem*. 2017 Jul 18;89(14):7307–7315. <https://doi:10.1021/acs.analchem.6b04734>
12. Selen, E., Baghla, R., Moore, I., Galermo, E., Zhang, Z., and Jones, E. Redefine bioanalysis with enhanced robustness on the SCIEX 7500+ system. SCIEX technical note MKT-31350-A. <https://sciex.com/tech-notes/pharma/bioanalysis-pk/redefine-bioanalysis-with-enhanced-robustness-on-the-sciex-7500-plus-system>

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High-coverage Targeted Lipidomics Analysis of Trace Cells based on Laser Microdissection-LC-MS/MS Method

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Key words: lipidomics; single cell; laser microdissection; Triple Quad 7500; MRM

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序号	代谢物	母离子	子离子	碰撞能量(V)	序号	代谢物	母离子	子离子	碰撞能量(V)
193	Carbamaoyl Asparatate	175.0	132.0	-12	228	Deoxyribose-Phosphate	213.0	79.0	-35
194	Maleic Acid	115.0	71.0	-13	229	UDP	403.0	305.0	-15
195	Parahydroxybenzoate	137.0	93.0	-21	230	Uridine	243.0	200.0	-21
196	Acetylphospahte	139.0	79.0	-22	231	dTDP	411.0	159.0	-29
197	2-Keto-Isovalerate	115.1	71.1	-13	232	GDP	442.0	159.0	-19
198	Methylmalonic Acid	117.0	73.1	-13	233	GTP	522.0	424.0	-23
199	Pyroglutamic Acid	128.0	82.1	-19	234	dGTP	506.1	159.0	-30
200	Citraconic Acid	129.0	85.1	-13	235	CDP	402.0	159.0	-29
201	2-Ketohaxanoic Acid	129.0	101.3	-13	236	CDP-Choline	487.0	428.0	-21
202	N-Acetyl-L-Alanine	130.0	88.0	-14	237	Adenosine 5-Phosphosulfate	426.0	346.0	-22
203	Hydroxyisocaproic Acid	131.0	85.1	-16	238	ADP	426.0	134.1	-24
204	p-Aminobenzoate	136.1	92.0	-18	239	X5P	365.0	97.0	-11
205	p-Hydroxybenzoate	137.0	93.0	-23	240	Cyclic-AMP	328.0	134.1	-33
206	Acetylphosphate	139.0	79.0	-24	241	3-Hydroxybutyryl-Coa	852.2	408.0	-43
207	Phenylpropionic Acid	145.0	101.0	-18	242	D-Glucose	179.1	89.1	-12
208	2-Hydroxy-2-Methylbutanedioic Acid	147.0	85.1	-17	243	Hexose Phosphate	259.0	79.0	-40
209	3-Methylphenylacetic Acid	149.0	105.0	-12	244	Glucose-6-Phosphate	259.0	97.1	-15
210	Hydroxyphenylacetic Acid	151.0	107.0	-18	245	F6P	259.0	97.1	-17
211	2,3-Dihydroxybenzoic Acid	153.0	109.0	-19	246	Fructose-6-Phosphate	259.0	169.0	-16
212	Allantoin	157.1	114.0	-17	247	Fructose-1,6-Bisphosphate	339.0	97.0	-30
213	Indole-3-Carboxylic Acid	160.0	116.0	-20	248	D-Glyceraldehyde-3-Phosphate	169.1	97.0	-10
214	Phenylpyruvate	163.0	91.0	-13	249	Dihydroxy-Acetone-Phosphate	169.0	79.0	-40
215	Atrolactic Acid	165.0	119.0	-21	250	1,3-Bisphopshoglycerate	265.0	79.1	-20
216	Phenyllactic Acid	165.0	103.1	-21	251	1,3-Diphopshateglycerate	265.0	79.0	-37
217	Allantoate	175.0	132.0	-14	252	3-Phosphoglycerate	185.0	97.5	-35
218	2-Isopropylmalic Acid	175.0	115.0	-19	253	PGA	184.9	97.1	-20
219	Pyrophosphate	176.8	158.6	-16	254	2-Phosphoglycerate	185.5	97.1	-35
220	Hydroxyphenylpyruvate	179.1	107.0	-13	255	Phosphoenolpyruvate	167.0	79.0	-41
221	4-Pyridoxic Acid	182.0	138.0	-18	256	Pyruvate	87.0	43.0	-20
222	Indoleacrylic Acid	186.0	142.0	-20	257	Lactate	89.0	43.2	-16
223	Xanthurenic Acid	204.0	160.0	-19	258	Gluc-D-Lactone	177.0	129.0	-13
224	Lipoate	205.0	171.0	-13	259	6-Phospho-D-Gluco-1,5-Lactone	257.0	97.0	-22
225	Orotate	155.0	111.0	-15	260	6-Phospho-D-Gluconate	275.0	97.0	-13
226	Orotidine-5-Phosphate	367.0	323.0	-18	261	Ribulose-5-Phosphate	229.0	79.4	-26
227	Dihydroorotate	157.0	113.0	-14	262	Ribose-Phosphate	229.0	79.0	-42

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Rapid Analysis of Lipids Components in Biological Sample by SCIEX Echo® MS+ System

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1. 一级质量范围TOF MS, m/z 300-1100, accumulation time, 0.2 s

- TOF MS was followed by 100 product ion scan, with accumulation time of 15 ms
- 二级质量范围MS/MS (CID); m/z 100 – 1100
- 源气参数: IS电压, 5000 V; 气帘气CUR, 35 psi; 雾化气GS1, 90 psi; 辅助气GS2, 45 psi; 离子源温度为350 °C; 碰撞气, 8

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- [1]. Pengwei Guan, Yuting Wang, Tiantian Chen, Jun Yang, Xiaolin Wang, Guowang Xu, Xinyu Liu; Novel Method for Simultaneously Untargeted Metabolome and Targeted Exposome Analysis in One Injection. *Analytical Chemistry*. 2025, 97(7): 3996-4004. <https://doi.org/10.1021/acs.analchem.4c05565>

Targeted flux analysis through the Mass Isotopomer Multi-Ordinate Spectral Analysis (MIMOSA) Workflow

Using the QTRAP 6500+ LC-MS/MS System with SelexION® Technology and Elucidata Software

Darren Dumlao¹, Tiago Alves² and Richard Kibbey²
¹SCIEX, USA, ²Yale University, New Haven, USA

Fluxomics is a relatively recent application of mass spectrometry-based techniques that employs the use of stable labeled isotopes to trace how specific metabolites and molecules are metabolized through biological pathways. Commonly, ¹³C or ¹⁵N labeled precursor molecules are fed into a biological system and allowed to be metabolized for a given amount before the experiment is quenched. The experimental design can include a single time point once the system reaches a steady-state or a series of time points for a kinetic understanding of how the metabolites are changing over time. Unlike isotope labeling experiments, which mainly focus on the incorporation of labeled atoms into certain metabolites, the output of a flux experiment is a rate. A flux rate reveals details about how a biological system changes between 2 or more cohorts. For a metabolite with an increased enrichment of label, a flux rate can distinguish if this is due to the pathway speeding up and forming more product, or accumulation due to subsequent steps in the pathway slowing down. While fluxomics is a powerful technique in understanding the dynamic changes in metabolism, these experiments require

an expert level of understanding of the pathways and are mathematically rigorous.

The Kibbey lab has developed a targeted fluxomics technique, designated M.I.M.O.S.A. (Mass Isotopomer Multi-Ordinate Spectral Analysis), which is focused on energy metabolism.¹ The MIMOSA technique utilizes a targeted MRM-based method to understand each isotopologue as well as the position of the labeled carbon atoms. For separation of the glycolytic and TCA cycle intermediates, the method employs two 5-min LC gradients in conjunction with differential mobility separation. MIMOSA has been fully incorporated into the Polly software (Elucidata) with their LC-MS/MS labeled workflow. This automated workflow has removed the barriers and improved the rigor for fluxomics research.

Key features of MIMOSA fluxomics

- Uses the QTRAP® 6500+ System with SelexION® Technology for high sensitivity and specificity
- Describes fluxes for glycolytic and TCA cycle intermediates for energy metabolism studies
- Deconvolves the labeling patterns at the MS/MS level to measure positional enrichment
- Experimental output calculates normalized rates, pathway directionality, and dynamic metabolic changes between cohorts
- Polly software fully automates the data processing (natural abundance correction, absolute % enrichment, background correction, flux calculations) and data visualization which dramatically reduces the sample-to-graph time


Figure 1. Fluxes were mapped to central carbon pathways to summarize the metabolic changes that occur in response to T-cell activation. The pathway map was generated in the Polly LC-MS/MS labeled workflow (Elucidata). Red indicates increases in flux and green indicates decreases in activated T-cells.

RUO-MKT-02-12698-A

- Polly is flexible and customizable and is capable of handling other stable labeled isotopes

Methods

Sample preparation: A protocol detailing the isolation, culture, and activation of primary human CD4+ T-cells is described in Hiemer *et al.*⁴ Briefly, cells were unstimulated or stimulated with ImmunoCult. Cells were harvested and incubated with fresh cold media without glucose or glutamine (MP Biomedicals) that had been supplemented with sodium bicarbonate (3.7 g/L), 10 mM HEPES, 0.2% essentially FFA free BSA (Sigma), 0.05 mM sodium pyruvate, 4 mM glutamine, and 0.45 mM L-lactate for 1 h. Cells were washed with PBS and resuspended in ¹³C-glucose labeled media that consisted of 5.5 mM [U-¹³C₆]-D-glucose (CLM-2001; Cambridge Isotopes). Cells were incubated in glucose labeling media for 4 h. Cells were subsequently washed with ice-cold PBS, and cellular metabolism was quenched with ice-cold 20% methanol in 0.1% formic acid supplemented with 10 μM d₄-taurine (Cambridge Isotopes), 3 mM sodium fluoride (Sigma), 1 mM phenylalanine (Sigma), and 100 μM EDTA. Cells were frozen in V-bottom plates (Bioexpress, Cat #T-3025-8B) and lyophilized overnight to dryness. Samples were resuspended in water and 2 μL were injected on column.

Chromatography: The M.I.M.O.S.A chromatography has been described in Alves *et al.*¹ In short, glycolytic intermediates were analyzed using an Amino Acid 50 x 2.1 mm, 3 μm column (Imtakt). The TCA cycle intermediates were analyzed using a Hypercarb 100 x 2.1 mm column (ThermoFisher). The gradient conditions for both are described in Table 1.

Mass spectrometry: MS analysis was performed on the QTRAP 6500+ System equipped with SelexION Technology and an IonDrive™ Turbo V Ion Source. The details of the MS

Table 1. Chromatography gradients used.

Glycolytic intermediates		TCA cycle intermediates	
Time (mins)	%B	Time (mins)	%B
Flow rate = 600 μL/min		Flow rate = 400 μL/min	
0	90	0	0
0.75	90	1	0
2	70	1.5	20
3.6	70	3.2	20
4	90	3.3	0
5	90	4	0
Mobile phase A – 100 mM ammonium acetate, pH 6.9 Mobile phase B – acetonitrile		Mobile phase A – 15 mM ammonium formate, 20 μL of 0.5 mM EDTA pH 6.9 Mobile phase B – acetonitrile	

RUO-MKT-02-12698-A

experiments have been described in Alves *et al.*¹ Briefly, the instrument parameters for the glycolytic intermediates experiment were the following: CAD 9, Curtain gas 30, IonSpray Voltage -3000, Temp 450, GS1 60, GS2 60, DMS modifier flow 250 μL/min, DMS temp low, SV 2200, DR off, DMS Modifier 2-propanol, DMS Modifier Composition low, Dwell 10 msec, Pause 3 msec. The instrument parameters used for the TCA cycle intermediates were: CAD 9, Curtain gas 30, IonSpray Voltage -3000, Temp 450, GS1 55, GS2 55, DMS modifier flow 325 μL/min, DMS temp medium, SV 2600, DR open, DMS Modifier 2-propanol, DMS Modifier Composition low, Dwell 10 msec, Pause 3 msec.

Data processing: The output .wiff data files were converted to mzML using MSConvert (version 3.0.19172, Proteowizard). Peak alignment and integration were performed using EI-Maven v0.12.0. Natural abundance correction, fractional enrichment, background correction, flux calculations, and graphical data visualization were calculated using Polly LC-MS/MS labeled workflow (Elucidata).

Flux analysis on T-cells

T-cells are known to employ a different metabolic program depending upon their activation state.² While in their resting state, T-cells utilize β-oxidation of fatty acids as a major energy source. Once stimulated, T-cells exhibit the Warburg effect and switch to aerobic glycolysis producing increased amounts of lactate. The uptake of glucose, glutamine, and fatty acids from the exogenous environment is required for T-cells to rapidly divide and mount an immune response.³ Here, the steady state metabolic changes that occur in human T-cells 24 hours after activation was studied using a targeted MIMOSA fluxomic approach with SCIEX QTRAP 6500+ System equipped with SelexION Technology and the Polly platform for data processing and visualization.

Data processing using EI-Maven

Flux data acquired on QTRAP 6500+ System with SelexION Technology were converted to mzML to process in EI-Maven for peak alignment and peak integration. Alternatively, MultiQuant™ Software (SCIEX) could have been used for peak integration, then uploaded into Polly as a .csv file with the proper format. In Figure 2 (top), M+2 glutamate (148/43) is used to showcase how the data are visually represented in a single plot with multiple overlaid traces. The bar graph on the right side denotes the peak area for each sample allowing for collective assessment of the data. Viewing the data in this manner allows for fast peak review, instead of the hundreds or thousands of mouse clicks performed when reviewing each scan independently. Importantly, manual integrations of the data can


Figure 2. The output from Polly LC-MS/MS labeled workflow after natural abundance correction. Additional tabs include fractional enrichment data, background correction, and a merged table. All data can be easily exported as a .csv file.

be performed as needed. EI-Maven outputs processed data as a table that contains detailed information about MRM transitions, peak area, retention time, S/N, and number of peaks detected as how many are considered good peaks (Figure 2, bottom). This table is used to QC the data where the user can denote a detected metabolite as either good or bad with a checkmark or an X, respectively. The data can be automatically transferred to Polly using EI-Maven's upload to the cloud function which requires an active Polly account.

Polly LC-MS/MS labeled workflow

The required files for workflow include a peak areas file (EI-Maven output), a metadata file detailing the metabolites detected (molecular formula, fragment, isotope label), and a sample metadata file that denotes which cohort each sample belongs to for comparison. After the files are uploaded, the user can access the NA correction tool to perform the natural abundance

correction. Natural abundance correction involves numerous calculations that considers the natural occurrence of the isotope ($^{13}\text{C} = 1.1\%$), the number of carbons in the parent molecular and resultant fragment with all the possible labeling combinations. Setting up these calculations is laborious, time consuming, and the introduction of user error at this step is quite common. In contrast, with this automated workflow, natural abundance is calculated with a single mouse click (Figure 3).

Label	Sample	NA Corrected	Name	Formula	metab_name	Pool_Total	Truth Value	NA Corrected with zero	
1309	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glucose labeled overnight 1 (sample 36)	248900.184854	Lactate 89/89	C3-000	Lactate	3273297.814206	false	249900.184854
1310	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glucose labeled overnight 2 (sample 39)	246001.69169	Lactate 89/89	C3-103	Lactate	4290077.72	false	246001.69169
1311	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glucose labeled overnight 3 (sample 40)	257657.622876	Lactate 89/89	C3-203	Lactate	3944329.856972	false	257657.622876
1312	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glucose labeled overnight 4 (sample 141)	234629.981196	Lactate 89/89	C3-303	Lactate	3004944.402940	false	234629.981196
1313	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glucose labeled overnight 5 (sample 144)	236201.642099	Lactate 89/89	C3-403	Lactate	3003436.482499	false	236201.642099
1314	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glucose labeled overnight 6 (sample 120)	192260.066528	Lactate 89/89	C3-503	Lactate	3315050.156528	false	192260.066528
1315	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glutamine labeled hr-1 (sample 125)	6121162.93993	Lactate 89/89	C3-003	Lactate	8121142.93993	false	6121142.93993
1316	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glutamine labeled hr-2 (sample 126)	4567933.9705	Lactate 89/89	C3-103	Lactate	4370915.268611	false	4567933.9705
1317	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glutamine labeled hr-3 (sample 62)	4572183.62205	Lactate 89/89	C3-203	Lactate	4599356.343209	false	4572183.62205
1318	C13_89_0_89_0	20170320_TcellTimeCourse-4L-Activated Glutamine labeled hr-4 (sample 116)	5617023.78285	Lactate 89/89	C3-303	Lactate	5570505.663035	false	5617023.78285

Figure 3. The output of stable isotope enrichment results calculated using Polly LC-MS/MS labeled workflow Phi (Beta).

In addition, the workflow also calculates the fractional enrichment of metabolite species and performs a background correction. The embedded search function allows the table to be filtered to only display the desired metabolite. The data can be represented graphically using the built-in visualization tools, Figure 4 shows the natural abundance corrected data for lactate. The natural abundance corrected data (represented as peak area, Figure 4A) for lactate and its isotopologues show that all isotopologue species have increased in activated T-cells (blue) compared to rested T-cells (orange). The fractional enrichment data (represented as the percentage of total lactate pool) show similar amounts of M+3 species between activated and resting T-cells (Figure 4B). An increased relative flux from glucose through lactate is suggested in activated T-cells by the lower dilutions from other non-labeled pathways. Further examination of the fractional enrichment profile of TCA intermediates show complex labeling patterns (Figure 4C and D). Citrate has an increased amount of ^{13}C incorporated into the M+2 species in general but there is more dilution between citrate and malate in activated T-cells indicating. While these data are suggestive of different utilization patterns of citrate, a more definitive interpretation requires the interpretation of positional enrichment differences between precursor and product metabolites (see below). All the data (tables and graphs) can be easily exported in a .csv or .png picture file, respectively. This fully automated workflow removes

the time consuming rigorous manual labor involved in processing these types of experiments.

Phi – Flux calculations

Once the data is corrected for naturally occurring isotopes, the individual flux rates (Phi or ϕ) can be calculated. This is accomplished in the Polly Phi (Beta) tab that contains three subtabs: Calculate Phi, Results: Glucose Labeled, and Results: Generic Labeled. The Calculate Phi tab is used to define the labeled isotopes used and to display the different cohorts. Also, Polly offers customizable options for metabolites with different labeling schemes and fragmentation patterns. The two results tabs display the phi values in table and graph form for glucose and generic (used for other labeled precursors like glutamine), respectively.

In this fluxomics experiment, only uniformly ^{13}C labeled glucose ($[\text{U-}^{13}\text{C}_6]\text{-D-glucose}$) was used. Figure 5 depicts the results for the glucose labeled output with the different phi values (rows) for each sample (columns) listed. Here 45 individual phi values and positional isotopomers were calculated when using $[\text{U-}^{13}\text{C}_6]\text{-D-glucose}$ as the tracer and these can be expanded to include others by user. Polly has built-in visualization tools to graphically display the flux data. Representative graphs with the expression formulas are shown to highlight the changes in metabolic flux between activated and rested cells (Figure 6). Activated T-cells


Figure 4. Visualization of natural abundance correction (NA) and fractional enrichment (FE). Representative graphs using Polly's built-in visualization tools for natural abundance correction (NA) and fractional enrichment (FE) in activated (blue bars) and rested (orange bars) T-cells. (A) depicts the natural abundance corrected isotope profile for lactate. The data is reported as average peak area. (B, C, and D) are the fractional enrichment plots for the isotope profiles of lactate, citrate, and malate, respectively. The data is reported as % of total metabolite pool. Each cohort had an n value = 6.

were observed to have a 2.5-fold increased enrichment in the pyruvate dehydrogenase (PDH) derived carbons of citrate (Figure 6A) despite only a slightly greater enrichment in lactate (Figure 4B). This identifies a greater contribution of carbohydrate oxidation to overall TCA flux in the activated cells. In contrast, there is a 70% greater loss of ^{13}C enrichment between Acetyl CoA (carbons C1,2) and those same carbons in glutamate (C4,5, Figure 6C). The same observations made by an independent calculation quantitatively validated this dramatic finding (Figure 6B). To represent the data in a biological context, metabolic flux changes that occurred in response to T-cell activation were mapped to central carbon pathways (Figure 1). Representing the data in this fashion can be useful for understanding metabolic changes, nutrient utilization, and pathway directionality.

Discussion

T-cells undergo a metabolic reprogramming when activated by a stimulus event. Using $[\text{U-}^{13}\text{C}_6]\text{-D-glucose}$, these metabolic changes were assessed after 24 hours of activation. An elevated flux of glucose through glycolysis (compared to other pyruvate generating fluxes) was observed in activated T-cells with the accumulation of ^{13}C -carbon in lactate. The dependence on glycolysis, which can rapidly generate energy, is thought to be in response to the cell's accelerated growth in mounting an immune response. Glucose was also utilized to generating citrate via PDH by a greater extent, but reduced label incorporations were seen in other TCA cycle intermediates and metabolites (malate and glutamate). One explanation of this observation is that a dilution of ^{13}C -carbon is occurring after citrate. Consistent with this, it has been reported that activated T-cells exhibit increased uptake of exogenous glutamine. Unlabeled glutamine could be

responsible for the observed dilution entering the TCA via glutamate to α -ketoglutarate exchange. A flux experiment using labeled glutamine (for instance, [1,2- ^{13}C]-L-glutamine) would be needed to provide an orthogonal confirmation.

This experiment was used to showcase the utility of a targeted fluxomics approach to better understand the dynamic metabolic changes that occur during T-cell activation. The sensitivity and speed of the QTRAP 6500+ System enabled the highest quality data to be collected for the 100+ MRMs within the method. The orthogonal separation of the SelexION Technology was instrumental in separating many of the isobaric compounds in glycolysis and the TCA cycle.

The Polly LC-MS/MS labeled workflow removed many of the barriers and bottlenecks associated with processing fluxomic data. This was accomplished by streamlining the processing from peak integration to natural abundance correction. In addition to reducing data processing time, the automated workflow eliminated user derived errors commonly associated with multi-tab spreadsheets. Another positive feature of Polly is the visualization tools. Comparisons between different cohorts can quickly be selected and graphed. The data and graphs are easily exported for use in manuscripts, posters, or presentations. The reducing the time from sample to graph allows the user to spend more time in understanding the biological impact of their experiment.


Figure 6. Representative graphs using Polly's built-in visualization tools for calculated fluxes (Phi values) for activated (blue) and rested (red) T-cells. (A) A greater enrichment in the pyruvate dehydrogenase (PDH) derived carbons of citrate was observed in activated T-cells. The next segment of the TCA cycle showed the opposite, namely there was a greater entry of glutamate carbon (via GDH or transamination pathways) in activated T-Cells as demonstrated by two independent calculations from citrate to glutamate **(B)** Phi_Cit_Glut and from Acetyl-CoA to citrate to glutamate **(C)** Phi_AcCitG, respectively. Each cohort had an n = 6.

RUO-MKT-02-12698-A

Phi Results for Glucose Labelled (for TCA Cycle)

Calculated Phi Actual Phi Identifier Expression Visualization

Steps 1 > 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45

MetID	20170320_TcellTimeCourse-Act-Activated Glucose labeledWeight-1 (sample 20)	20170320_TcellTimeCourse-Act-Activated Glucose labeledWeight-2 (sample 49)	20170320_TcellTimeCourse-Act-Activated Glucose labeledWeight-3 (sample 98)	20170320_TcellTimeCourse-Act-Activated Glucose labeledWeight-4 (sample 147)	20170323_TcellTimeCourse-Act-Activated Glucose labeledWeight-5 (sample 14)	20170326_TcellTimeCourse-Act-Activated Glucose labeledWeight-6 (sample 120)	20170326_TcellTimeCourse-Act-Activated Glucose labeledWeight-7 (sample 8)
1 OAA_PIC	0.0322992395456	0.2342627962387	0.008962769059031	0.0274297047932	0.00744855430687618	0.030547289492229	0.00688749195021193
2 OAA_TCA	0.013407073300272	0.134010649442052	0.02391112296172676	0.0193612162042409	0.0347268421852872	0.0214921172328441	0.0756629771732311
3 PHU_Oxif	0.1891601600279	0.17639842893271	1.19086271790639	0.005490304980387	0.145900786514654	0.27050580505886	1.3079214689156
4 OLCIC	0.2113441574	0.00574544383	0.005770032189129	0.015457983005553	0.0018721272	0.014700949189228	0.0010481846
5 PHU_AcEtS	0.127149104072381	0.1903621320782	0.14444314958213	0.146912814991030	0.129200677172141	0.148320262160227	0.219852191867580
6 PHU_Oxif	0.17194275448791	0.37122310250276	0.406915296162143	0.49140364481122	0.5440044754113	0.494565361571356	0.3470370791306
7 OAA_M2	0.10418880291216	0.114248348778646	0.12374274988217	0.12770561982119	0.12024866515781	0.108056267677726	0.1991936793145
8 OAA_M3	0.049693881576291	0.0322892881443	0.02228948633317	0.049401186257208	0.042132074831823	0.056248426123959	0.083464587713421
9 OAA_M4	0.002401084928157	0.012615349357791	0.037351105418844	0.004290075722956	0.0379479683641227	0.006994589133317	0.0254407881181250
10 PHU_M1	0.3423787144279	0.6221781348153	0.5026780148236	0.5423211165763	0.666882217956	0.767888818837	0.7320963227414

Showing 1 to 10 of 45 entries

Previous 1 2 3 4 5 Next

Figure 5. Phi results for TCA cycle calculated from MIMOSA deconvolved isotopomers from $[U-^{13}C_6]$ -D-glucose tracer labeling of T-Cells in Polly Phi (beta). Shown are 10 of the 45 Phi calculations (identifier column) for the different T-Cell samples. Individual mathematical descriptions of the different Phis shown in a different tab (see Figure 6). Data from individual samples (identifier descriptions provided) are shown as fractional values (multiply x 100 for percent).

References

- Alves, TC *et al.* (2015) Integrated, Step-Wise, Mass-Isotopomeric Flux Analysis of the TCA Cycle. *Cell Metab.* **22**, 936-947.
- Geltink, RIK, Kyle, RL, & Pearce. (2018) Unraveling the Complex Interplay Between T Cell Metabolism and Function. *Annu. Rev. Immunol.* **36**, 461-488.
- van der Windt, GJW & Pearce, EL (2012) Metabolic switching and fuel choice during T-cell differentiation and memory development. *Immunol. Rev.* **249**, 27-42.
- Heimer S *et al.* (2019) Integrated Metabolomic and Transcriptomic Profiling Reveals Novel Activation-Induced Metabolic Network in Human T cells. [Pre-publication](#).

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Rapid analysis of 11 energy-rich phosphate compounds by SCIEX Echo™ MS System

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Keywords: SCIEX Echo™ MS System, Energy Rich Phosphate Compounds, Nucleotide Compounds, dATP, dCTP, dTTP.

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- [1]. Henryk S, Markus N, Marco H; Enhanced nucleotide analysis enables the quantification of deoxynucleotides in plants revealing connections between nucleoside and deoxynucleoside metabolism. THE PLANT CELL 2021; 33: 270–289.
- [2]. Rapid Analysis of 11 Energy Rich Phosphate Compounds by SCIEX Triple Quad Mass Spectrometer. SCIEX technical note RUO-MKT-02-14446-ZH-A.

The main body of the page is a large, empty white space, likely intended for text or content.

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© 2025 DH Tech. Dev. Pte. Ltd. MKT-37001-A