



BSA Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments

Protocol

1 Product Description

Use the BSA Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments to test instrument function, optimize instrument parameters, and calibrate the mass scale using a known standard. The kit is designed for use in the Peptide Mass Standards Kit for Calibration of AB SCIEX MALDI-TOF Instruments with the Voyager™ Biospectrometry™ Workstation or with any AB SCIEX MALDI-TOF (matrix-assisted laser desorption ionization time-of-flight) mass spectrometry system.

The calibration standard in the kit requires minimal preparation and covers a mass range of 60,000 to 70,000 daltons (Da). The standard was chosen for this kit on the basis of high purity and consistent results when analyzed using MALDI-TOF mass spectrometry.

Applications

Using the BSA Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments, you can monitor results in a mass range from 60,000 to 70,000 Da for mass assignment, resolution, and sensitivity.

2 Materials

Materials Provided

The BSA Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments includes —lyophilized bovine serum albumin, **BSA** (1 vial, 2 mg/vial)

Matrix and Diluent

- **Matrix B:Sinapinic Acid**—3,5-Dimethoxy-4-hydroxycinnamic acid, 1 vial (7 to 10 mg/vial)
- **Matrix B Diluent**—30% acetonitrile in 0.3% TFA, 1 vial (1 ml/vial)

3 Preparing Matrix and Calibration Standards

See Section 7, Storing the Kit, for storage and stability conditions of prepared reagents.

▲ WARNING CHEMICAL HAZARD. Diluent (with acetonitrile) is a flammable liquid and vapor. It may cause eye, skin, and respiratory tract irritation, central nervous system depression, and heart, liver, and kidney damage. Please read the MSDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

▲ CAUTION CHEMICAL HAZARD. Sinapinic acid may cause eye, skin, and respiratory tract irritation. Please read the MSDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Reconstitute BSA

Store the BSA Calibration Standard at –20° C prior to use. Sinapinic Acid Matrix from a Sequazyme kit is required to use this standard.

1. Add 0.5 ml HPLC-grade water to the **BSA** vial.
2. Recap and vortex for 1 minute.

You can use the standard immediately.

Sinapinic Acid Matrix

Use sinapinic acid matrix B from the Peptide Mass Standards Kit for Calibration of AB SCIEX MALDI-TOF Instruments.

1. Add the entire contents of the **Matrix B Diluent** vial to the **Matrix B:Sinapinic Acid** vial.
2. Vortex for 1 minute.
3. Warm the solution, if necessary, to completely dissolve the solid. Use the supernatant for sample preparation.

Section	Page
1 Product Description.....	1
2 Materials.....	1
3 Preparing Matrix and Calibration Standards	1
4 Loading Matrix:Standard on Sample Plates	2
5 Analyzing Calibration Standards	2
6 BSA Calibration Standard Spectrum and Masses.....	2
7 Storing the Kit.....	2
8 Accessories, Spare Parts, and Ordering Information	3
9 Technical Support.....	3

Standard

Final concentration of the calibration standard after mixing with matrix is 4 pmol/ μ l.

1. Pipette 1 μ l of the reconstituted BSA solution into a microcentrifuge tube.
2. Add 29 μ l of the sinapinic acid matrix.
3. Vortex briefly at low speed.

4 Loading Matrix:Standard on Sample Plates

Each matrix:standard preparation yields 30 μ l, enough to load approximately 30 sample positions.

1. Pipette 1.0 μ l of sample on the MALDI plate.
2. Allow the mixture to air dry until all solvent is evaporated, usually less than 5 minutes.

Analyze within one day for best results.

5 Analyzing Calibration Standards

Standard methods on your system

Your Voyager™ workstation includes standard methods created during factory testing of your instrument. Standard methods contain typical parameters for a given mass or mass range. Standard methods are useful as starting points when you create methods to analyze the standards in the Peptide Mass Standards Kit for Calibration of AB SCIEX MALDI-TOF Instruments.


Analyzing standards

Analyze the standard according to the needs of your application.

Refer to Section 6, BSA Calibration Standard Spectrum and Masses, for representative spectra and mass assignments for standards and matrices in the kit.

6 BSA Calibration Standard Spectrum and Masses

Figure 1 shows a representative spectrum for the test mixture. Masses are included in the spectra for peak identification only. Use the precise masses listed in Table 1 for calibration.



Standard	Charge (n)	(M+nH) ⁿ⁺ Average
BSA (Bovine Serum Albumin)	+1	66,431
	+2	33,216

7 Storing the Kit

Store the BSA Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments and components of the kit under the following conditions. Avoid prolonged exposure to light.

Kit Component	Storage Temperature	Stability
Unopened kit	-20°C	1 year from date of shipment
Reconstituted standards	-20°C	6 months
Reconstituted matrix	4°C	1 week

8 Accessories, Spare Parts, and Ordering Information

Description	Quantity (Kits)	Part Number
Peptide Mass Install Kit for AB SCIEX MALDI-TOF Instruments includes: <ul style="list-style-type: none">• Peptide Mass Standards Kit for Calibration of AB SCIEX MALDI-TOF Instruments• BSA Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments• IgG Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments	1	4316866
Peptide Mass Standards Kit for Calibration of AB SCIEX MALDI-TOF Instruments	1	P2-3143-00
BSA Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments	1	2-2158-00
IgG Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments	1	GEN602151

9 Technical Support

AB Sciex Pte. Ltd. is committed to meeting the needs of your research through enabling technologies like the BSA Calibration Standard Kit for AB SCIEX MALDI-TOF Instruments. Our dedicated support staff is available to answer questions about using this product to the fullest extent possible.

For further details or for answers to questions related to other products contact AB Sciex Pte. Ltd. Refer to the back page for contact information.

© 2012 AB SCIEX. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

Information in this document is subject to change without notice. AB Sciex Pte. Ltd. assumes no responsibility for any errors that may appear in this document. This document is believed to be complete and accurate at the time of publication. In no event shall AB Sciex Pte. Ltd. be liable for incidental, special, multiple, or consequential damages in connection with or arising from the use of this document.

AB Sciex Pte. Ltd.

500 Old Connecticut Path
Framingham, MA 01701 USA
Phone: +1 877.740.2129
Fax: +1 650.627.2803

Technical Resources and Support

For the latest technical resources and support information for all locations, please refer to our Web site at:

www.sciex.com

www.sciex.com

Printed in the USA, 12/2012
Part Number 4317219 Rev. D