# : Post-translational modification (PTM) profiling on fusion protein aflibercept using a novel fragmentation technique

### Xuezhi Bi<sup>1</sup>, Zoe Zhang<sup>2</sup> <sup>1</sup>Bioprocessing Technology Institute, A\*STAR Research Entities, Singapore; <sup>2</sup>SCIEX, Redwood City, USA

14000

13000-

12000

11000-

10000-

9000-

5000-

4000-

1500

1000-

## INTRODUCTION

Aflibercept is an Fc fusion protein used to treat wet macular degeneration and metastatic colorectal cancers. It contains five different N-linked glycosylation sites with extensive glycan distribution, which exhibit additional levels of complexity in terms of characterization. Glycosylation including glycan structures and glycan occupancy, plays a vital role in stability and efficacy of a biotherapeutic. Therefore, the glycan profile is often considered a critical quality attribute (CQA). Bottom-up workflows are frequently used for glycopeptide characterization as protein sequence information and glycan profiling can be obtained, simultaneously. However, traditionally used collision-induced dissociation (CID) struggles to provide accurate information on glycosylated peptides. Alternative fragmentation on the other hand is challenged by providing such information easily and reproducibly in order to be used routinely within the biopharma industry. Here within, the ability of electron activated dissociation (EAD) for peptide mapping MS/MS analyses is examined.

EAD, a newly developed fragmentation technique unique to the ZenoTOF 7600 system (SCIEX), allows for optimized peptide fragmentation with a platform method.<sup>1,2</sup> This allows for accurate localization of the linked glycans alongside confident identification of the peptide through high MS/MS sequence coverage of the peptide backbone.<sup>3,4</sup> In addition, Zeno EAD enables fast and highly sensitive data-dependent acquisition (DDA). This approach further overcomes long existing challenges of alternative fragmentation, such as long reaction times, limited charge state fragmentation, and low sensitivity. EAD allowed for analysis of a wide range of peptide types derived from digested aflibercept in one single DDA. Singularly charged and large sized peptides, PTM localization, and differentiation of amino acid isomers were studied.

# **MATERIALS AND METHODS**

#### Sample preparation:

Aflibercept was denaturated with 7.2 M guanidine hydrochloride (HCI) in a buffer of 50 mM Tris-HCI, pH 7.0, and free cysteines were capped with 5 mM iodoacetamide. Digestion was performed with either trypsin/Lys-C or trypsin/Lys-C/Asp-N at 30° C for 16 h. The reaction was stopped with 1% formic acid.

#### **HPLC conditions:**

The digested peptides were separated with a CSH C18 column (1.7 µm particle size, 130 Å, 2.1 × 100 mm, Waters) using an ExionLC system. The mobile phase A consisted of water with 0.1% formic acid, while the organic phase B was acetonitrile 0.1% formic acid. A gradient profile was used at a flow rate of 300 µL/min. The column temperature was maintained at 50°C.

#### MS and MS/MS conditions:

Data were acquired with an information dependent acquisition (IDA) method using the SCIEX ZenoTOF 7600

#### Data analysis:

Data was processed using the new Biologics Explorer software (SCIEX) with the following nodes:





![](_page_0_Figure_17.jpeg)

![](_page_0_Picture_18.jpeg)

Figure 2. 3D visualization for identified peptides at RT = 27-37 min. X-axis shows RT, y-axis m/z, z-axis intensity.

Protein	Location	Modification	EAD_1 [%]	EAD_2 [%]	EAD_3 [%]
Protein	N36	A2G2F1S2	48.7	48.9	49.0
Protein	N36	A2G2F1S1	30.1	29.9	30.7
Protein	N36	A3F1S3	7.3	6.9	6.3
Protein	N36	A3F1S2	6.3	6.2	6.0
Protein	N36	A2G2F1	5.4	5.8	6.0
Protein	N36	A2G1F1S1	0.7	0.7	0.7
Protein	N36	A4F1S3	0.4	0.4	0.4
Protein	N36	A4F1S4	0.2	0.3	0.0
Protein	N36	A3F1	0.2	0.2	0.2
Protein	N36	M4A2G2S1	0.2	0.2	0.2
Protein	N36	Deamidated	0.1	0.1	0.1

#### Figure 3. Glycans identified at site N36. Three replicates

demonstrate the reproducibility of the identification independently of abundance of the glycoform.

![](_page_0_Figure_23.jpeg)

peptide with two oxidations.

![](_page_0_Figure_25.jpeg)

Figure 4. Data for unmodified and oxidized peptide SDTGRPFVEMYSEIPEIIHMTEGR. Rich fragmentation leading to high fragment coverage and confident identification was obtained with the EAD platform method for all peptides shown. Mirror plots of MS/MS EAD spectra for (A) non-modified peptide vs. peptide with one oxidation at M10 and (B) nonmodified peptide vs. peptide with two oxidations at M10 and M20. (C) 3D plot for identified peptides visualizes all relevant parameters easily: RT and *m*/*z* shift for oxidized peptides in comparison to un-modified peptide and intensity difference can quickly be reviewed. Blue circle: unmodified peptide. Red circles: peptide with one oxidation. Green circle:

Figure 5. Data for unmodified and glycosylated peptide GFIISNATYK. MS/MS EAD spectra of (A) the +1 charged state of the aglycosylated peptide and (B) of the peptide carrying 0.5% G3S3F at N68. Despite the low charge state of +1, rich fragmentation was obtained leading to 89% fragment coverage with EAD for the unglycosylated peptide. EAD is an alternative fragmentation mechanism, which can fragment all charge states of peptides with an easy-to-use platform DDA method. Despite the low abundance of the glycosylated peptide of ~0.5%, fragments with intact glycosylation attached were obtained, providing unambiguous localization information. In addition, peptide backbone information was achieved. The MS/MS data was achieved using an automatically operating ion trap (Zeno trap) prior to the TOF region, which significantly enhances the duty cycle up to 95 %.

Protein	Location	Modification	EAD_1 [%] <sup>1</sup>	EAD_2 [%]	EAD_3 [%]	High
Protein	N196	A2G2S1	37.8	39.1	39.5	
Protein	N196	A2G2	21.7	18.0	19.5	
Protein	N196	A2G2S2	15.8	16.2	16.2	
Protein	N196	M5	8.3	10.6	8.8	
Protein	N196	A2G1	3.0	2.5	2.7	
Protein	N196	M4A1G1	2.4	1.6	1.5	
Protein	N196	M4A1G1S1	2.2	2.4	2.3	
Protein	N196	A2G1S1	2.1	2.4	2.5	
Protein	N196	A2G2F1	1.9	1.8	2.0	
Protein	N196	M5A1G1S1	1.8	1.9	1.9	
Protein	N196	M5A1G1	1.4	1.8	1.5	
Protein	N196	M4	0.6	0.4	0.3	
Protein	N196	A2G2F1S2	0.3	0.3	0.3	
Protein	N196	M4A2G2S1	0.3	0.3	0.3	
Protein	N196	M4A2G2	0.3	0.2	0.2	
Protein	N196	A2G2Sg1	0.1	0.1	0.1	
Protein	N196	A2G1B	0.1	0.1	0.0	
Protein	N196	A1	0.1	0.1	0.0	
Protein	N196	A2G1F1	0.0	0.1	0.1	
Protein	N196	A2G2B	0.0	0.2	0.2	Low

Figure 6. Glycans identified on site N196. 20 different glycans were identified on the glycopeptide. Three replicates demonstrate great reproducibility of the identification and relative quantification.

## CONCLUSIONS

- quality, even for low-abundant species.
- powered by SCIEX OS software.

## **REFERENCES**

- spectrometer. Submitted to JASMS. MKT-02-12639-A.
- technical note. RUO-MKT-02-12980-A.

# **TRADEMARKS/LICENSING**

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-10-13941-A.

![](_page_0_Picture_44.jpeg)

MS/MS fragment detection can be significantly enhanced compared to traditional accurate mass MS/MS analysis when using the Zeno trap. This enables confident identification due to excellent MS/MS raw data

• The robust, reproducible and easy-to-use alternative fragmentation technique of EAD enables users to directly answer challenging analytical questions with a platform method on the ZenoTOF 7600 system

• Automatic data analysis with Biologics Explorer software provides an accurate data analysis tool and powerful visualization feature, which enables complex glycan profiling.

Baba T et al. (2015) Electron capture dissociation in a branched radio-frequency ion trap, Anal Chem, 87, 785–792. Baba T et al. (2020) Dissociation of biomolecules by an intense low-energy electron beam in a high sensitivity time-of-flight mass

Comprehensive peptide mapping of biopharmaceuticals utilizing electron activated dissociation (EAD). SCIEX technical note, RUO-

A new electron activated dissociation (EAD) approach for comprehensive glycopeptide analysis of therapeutic proteins, SCIEX