

Shane Needham<sup>1</sup>, Eshani Nandita<sup>2</sup>, Lei Xiong<sup>2</sup>, Elliott Jones<sup>2</sup>, Zoe Zhang<sup>2</sup>, Kerstin Pohl<sup>3</sup> <sup>1</sup>Alturas Analytics Inc., Moscow, Idaho; <sup>2</sup>SCIEX, Redwood City, California; <sup>3</sup>SCIEX, Framingham, Massachusetts

# **ABSTRACT**

Given the increased focus on producing new protein and peptide therapeutics, there is a resulting demand for highly sensitive and robust quantitative bioanalytical techniques to ensure proper testing for their safety and efficacy. Bioanalysis of peptide therapeutics is often faced with analytical challenges such as inadequate sensitivity and complexity of the matrix resulting in poor selectivity. High-resolution accurate mass spectrometry (HRAMS) has been increasingly adopted in bioanalytical workflows as it provides high selectivity with narrow mass extraction windows. As part of this work, evaluation of enhanced duty cycle by a novel ion beam to timeof-flight (TOF) pulser efficiency was performed. The new quantitative enhancements were evaluated for peptide quantification in a complex biological matrix.

# INTRODUCTION

- Traditional workflows for quantitative bioanalyses, such as immunological assays, have been displaced by LC-MS/MS analysis on a triple quadrupole mass spectrometer. Immunoassays often lack selectivity, specificity, and have a limited linear dynamic range
- While the triple quadrupole platform provides excellent sensitivity and quantitative performance, there can be some limitations with background interference based on the lower resolution considering the type of mass analyzer. Background interference is a common issue for workflows where analytes are present in a highly complex matrix. High-resolution accurate mass spectrometry (HRAMS) has been increasingly adopted for quantitative bioanalysis<sup>1,2</sup>
- The ZenoTOF 7600 system offers an exceptional combination of mass resolution, sensitivity, and acquisition speed for quantitative analysis. It also aids in more accurate and automated integration, the potential for less ion path tuning, the ability to change measured fragments post-acquisition, and improved reproducibility and LDR when interferences are mitigated
- These attributes complement the excellent sensitivity of a nominal mass triple quadrupole system for a biopharmaceutical lab which requires a full range of capabilities

# MATERIALS AND METHODS

## Sample preparation:

Digested universal proteomics standard (Sigma Aldrich, 48 human proteins) was spiked into rat plasma and serially diluted to generate a calibration curve. Samples were denatured using N-octyl-glucoside followed by reduction using dithiothreitol, and alkylation using iodoacetamide. Protease digestion was performed using Trypsin/Lys-C at 37°C overnight, at an enzyme:protein ratio of 1:25. Digestion was stopped using formic acid. Supernatant was subjected to LC-MS/MS analysis.

## LC conditions:

An ExionLC system was used for analyte separation. A volume of 20 µL was injected for analysis. Mobile phase A consisted of water with 0.1% FA in water, while organic phase B was composed of 0.1% FA in acetonitrile. For analyte separation, the operating flow rate was set to 0.5 mL/min using a Phenomenex Kinetex C18 column (3 x 50 mm, 2.6 µm, 100 Å). The column oven temperature was set to 40°C.

## **MS** conditions:

Samples were analyzed in triplicate. Method details such as source and gas parameters and MS conditions are summarized in Table 2. Sample analysis was performed using scheduled Zeno MRM<sup>HR</sup> on the ZenoTOF 7600 system. The ZenoTOF 7600 system provides a scan speed of 133 Hz.

## Data processing:

MRM data were processed using SCIEX OS software 2.0. Integration was performed using the MQ4 algorithm. Linear regression with 1/x weighting was used for quantification of all peptides. The XIC peak width was set to 0.05 Da for both MS/MS and MS1 quantification.

# RESULTS

- With traditional time-of-flight MS/MS, fragment ions arriving from the collision cell are often lost in transmission between TOF pulses due to differences in velocity
- The Zeno trap ensures greater ion transmission by controlling the ion beam from the collision cell into the TOF accelerator (Figure 1)









Figure 1. The Zeno trap enables ion beam control from the collision cell before entering the TOF accelerator.



### Figure 2. Greater sensitivity for peptide quantification was observed for Zeno MRM<sup>HR</sup>.

- Significantly lower LLOQs were achieved using the Zeno trap
- The level of improvement in sensitivity varied across the peptide of interest (Figure 2)

Figure 3. Factor of improvement in LLOQ for peptide quantification with Zeno MRM<sup>HR</sup> compared to MRM<sup>HR</sup>.

- The Zeno trap enabled significant improvements in MS/MS sensitivity
- Out of 48 peptides measured, 75% of peptides showed ≥3-fold improvement in LLOQ
- On average, a 5-fold improvement in LLOQ for peptide quantification was observed with Zeno MRM<sup>HR</sup> in comparison to standard MRM<sup>HR</sup> (Figure 3)



- The accessibility of TOF MS/MS data can be advantageous as postacquisition data decisions can be made on which measured fragments can be utilized for MRM<sup>HR</sup>
- A 3-fold improvement in LLOQ was achieved with quantification using summation of multiple dominant fragment ions (Figure 5)

# Enhanced selectivity and sensitivity for peptide quantification in a complex matrix using high-resolution LC-MS/MS workflow

Figure 4. ZenoTOF 7600 system offered improved mass resolution between the target peptide and matrix-related components.



Figure 5. Summation of multiple fragment ions enhances assay sensitivity.



Figure 6. Strong linearity was achieved with Zeno MRM<sup>HR</sup>.

## **CONCLUSIONS**

- enhances the duty cycle through the accumulation of ions during each TOF pulse
- accumulation of ions during each TOF pulse for enhanced duty cycle
- ZenoTOF 7600 system
- resulting in a 3-fold enhancement in LLOQ
- overall data integrity

## REFERENCES

# **TRADEMARKS/LICENSING**

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <u>www.sciex.com/diagnostics</u>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks) © 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-10-13893-A. AB SCIEX™ is being used under license.



P05413_LILTLTHGTAVC[CAM]TR y12			P02741_ESDTSYVSLK y6			P63165_FLFEGQR y5		
Conc. (ng/mL)	CV (%)	Accuracy (%)	Conc. (ng/mL)	CV (%)	Accuracy (%)	Conc. (ng/mL)	CV (%)	Accuracy (%)
1177.33	1.78	100.43	1842.46	1.66	97.90	3103.01	5.27	99.24
412.06	3.20	102.60	644.86	5.54	106.65	1086.05	2.87	101.85
144.22	5.22	91.90	225.70	3.61	99.43	380.12	3.60	104.36
50.48	11.05	91.71	78.99	3.92	96.48	133.04	12.48	92.80
17.67	0.70	98.30	27.65	8.66	99.54	46.56	8.27	87.47
6.18	9.38	102.62	N/A	N/A	N/A	16.30	14.75	114.28
2.16	13.00	112.44	N/A	N/A	N/A	N/A	N/A	N/A

- Strong linearity was achieved for all peptides analyzed (Figure 6)
- Excellent accuracy and precision was achieved using the Zeno trap (Table 1)
- Accuracy at the LLOQ was within 80%-120%, while for all other non-zero calibrators, accuracy was within 85%-115% of the nominal concentration
- Overall, precision was <15%, demonstrating high reproducibility</li>

• An average of 5-fold in LLOQ improvement was achieved for peptide quantification in this sample set using the Zeno trap, which

A highly accurate and reproducible quantitative workflow for peptides was demonstrated using the ZenoTOF 7600 system

Higher sensitivity, based on LLOQ levels, was achieved for peptide quantification using the Zeno trap by improving duty cycle through the

• Greater selectivity was reached between target peptides and matrix-related components with the higher mass resolution offered by the

Improved LLOQs were reached by summing of multiple highly abundant fragment ions along with the availability of TOF MS/MS data,

• Automated and accurate peak integration was easily attainable on the ZenoTOF 7600 system, with greater mass resolution ensuring

1. Mike-Qingtao Huang, Zhongping (John) Lin, Naidong Weng (2013). Applications of high-resolution MS in bioanalysis. Bioanalysis 5(10):1269-1276.

2. Yuan-Qing Xia, Jim Lau, Timothy Olah, Mohammed Jemal (2011). Targeted quantitative bioanalysis in plasma using liquid chromatography/high-resolution accurate mass spectrometry: an evaluation of global selectivity as a function of mass resolving power and extraction window, with comparison of centroid and profile modes. Rapid Communications in Mass Spectrometry **25**(19):2863-2878.