# Automated mass spectra comparison algorithm for high-throughput compound QC applications

**SCIEX**: Chang Liu; Gordana Ivosev; Thomas R. Covey **Pfizer**: Alandra Quinn; Jeff Chin; Hui Zhang

ASMS 2021, Philadelphia, PA, November 1, 2021



#### Compound management and distribution

# **//// Echo**<sup>®</sup> MS System



The Pou

# Compound QC challenge

# **Echo** MS System



© 2021 DH Tech. Dev. Pte. Ltd.

### AEMS system







# AEMS for high-throughput analysis

# **//// Echo**<sup>®</sup> MS System



#### TOF based AEMS system prototype





OPP MS	Found	Not Found	Review	total cpds	
Plate 1	304	10	6	320	
Plate 2	299	20	1	320	
Plate 3	308	10	2	320	
Plate 4	287	31	2	320	
Plate 5	285	34	1	320	
Plate 6	285	30	5	320	
Plate 7	281	28	11	320	
Plate 8	288	31	1	320	
Plate 9	297	16	7	320	
Plate 10	291	14	15	320	
Plate 11	317	3	0	320	
Plate 12	307	11	2	320	
Plate 13	304	15	1	320	
Plate 14	308	12	0	320	
Plate 15	311	9	0	320	
Plate 16	309	10	1	320	
Plate 17	303	17	0	320	
Plate 18	289	26	5	320	
Plate 19	298	16	6	320	
Plate 20	73	3	4	80	
Plate 21	310	9	1	320	
Plate 22	300	16	4	320	
Total cpds	6354	371	75	6800	
	% mass match	<mark>93%</mark>			
	% not Found	5%			
	% review	1%			

22 plates by AEMS within <1 day

Same samples would take > 1 month with current LC/MS based QC



### TOF based AEMS data

# Echo<sup>®</sup> MS System



SCIEX The Power of Precision

© 2021 DH Tech. Dev. Pte. Ltd.

#### **TOF based AEMS data**

# Echo<sup>®</sup> MS System





### Automated data processing overview

- Data splitting: correlate well position with each signal peak
- Integration and visualization
- Background identification and subtraction
- Spectra comparison



Echo<sup>®</sup> MS System

### Data splitting

# Echo<sup>®</sup> MS System

#### CORRELATE WELL POSITION WITH EACH SIGNAL PEAK





### Integration and report

**//// Echo** MS System

- Input information
  - Split data (MS signal correlated with well position)
  - Compound information table
  - Integration settings (MS window width, etc.)
- Report table
  - Intensity, mass accuracy, isotope pattern matching score, S/N, etc.

sample	formula_group_1	charge_agent_group_1	formula_group_2	charge_agent_group_2	formula_group_3	charge_agent_group_
A1	C5H9NO2	H+	C5H9NO3	H+	C5H9NO4	H+
A2	C6H13NO	H+	C5H9NO3	H+	C5H9NO4	H+
A3	C6H11NO2	H+	C6H11NO3	H+	C6H11NO4	H+
A4	C6H11NO2	H+	C6H11NO3	H+	C6H11NO4	H+
A5	C5H9NO2	H+	C5H9NO3	H+	C5H9NO4	H+
A6	C5H9NO2	H+	C5H9NO3	H+	C5H9NO4	H+
A7	C5H9N	H+	C5H9NO	H+	C5H9NO2	H+
A8	C5H9N	H+	C5H9NO	H+	C5H9NO2	H+
A9	C7H13N	H+	C7H13NO	H+	C7H13NO2	H+
A10	C7H13N	H+	C7H13NO	H+	C7H13NO2	H+
A11	C6H13N	H+	C6H13NO	H+	C6H13NO2	H+
A12	C6H13N	H+	C6H13NO	H+	C6H13NO2	H+
A13	C6H13NO	H+	C6H13NO2	H+	C6H13NO3	H+
A14	C6H13NO	H+	C6H13NO2	H+	C6H13NO3	H+
A15	C6H9NO2	H+	C6H9NO3	H+	C6H9NO4	H+
A16	C6H9NO2	H+	C6H9NO3	H+	C6H9NO4	H+
A17	C7H11NO2	H+	C7H11NO3	H+	C7H11NO4	H+
A18	C7H11NO2	H+	C7H11NO3	H+	C7H11NO4	H+
A19	C7H11NO2	H+	C7H11NO3	H+	C7H11NO4	H+
A20	C7H11NO2	H+	C7H11NO3	H+	C7H11NO4	H+
A21	C6H9NO2	H+	C6H9NO3	H+	C6H9NO4	H+

0,7	Sample 🔍 📎	MZ (1) 0 🖉 📎	Status (1) 🖓 🖉	Intensity (1) 💷 📉	5/N (1) 0 🖉 📎	Abs Mz Error (1) 0 / 🔍	Error Type (1) 0 🖉 🖄	AVG Ratio Diffs (1)
)	A17	470.1145	9K 👻	379850.1	17067.9	1	Positive	0.0364375691427
)	A18	465.1356	08 🗸 👻	680719.1	10490.3	1	Positive	0.041457953915
)	A19	474.1247	00 👻	705874.2	32371.8	0		0.0409172184878
)	A20	465.1356	08 🔍 🕶	218335.2	3208.2	2	Positive	0.0361845186635
)	A21	471.1713	08 🔍 🕶 -	23219.5	413.1	1	Positive	0.0785033672792
)	A22	465.0913	0K 🚽 👻	1670878.0	63492.6	0	-	0.0485186411713
)	473	471 1713		70104 0	1617 3	٩	Positiva	A A37AA75385586



### Data review and visualization

Echo<sup>®</sup> MS System

#### HEAT MAP



#### MASS SPECTRA AND XIC VISUALIZATION





XIC for Sample I7, Group #1 M/Z: 472.1553, +/- 0.02 Da





#### Validation rules

# Echo<sup>®</sup> MS System



- Validation rules
  - Mass accuracy
  - Isotope pattern
  - S/N



Spectrum for Sample O23, Group #1 around the MZ: 458.1509 (+/- 8 Da)







Pte. Ltd.

### Background mass spectra pemoval

Background identification and subtraction (with well-specific scaling factor)



© 2021 DH Tech. Dev. Pte. Ltd.

Echo<sup>®</sup> MS System

### Background mass spectra removal

Echo<sup>®</sup> MS System

Background identification and subtraction (with well-specific scaling factor)





© 2021 DH Tech. Dev. Pte. Ltd.

#### Spectra comparison

**Echo** MS System

- Auto pairing from two sets of data (reference and test sets)
- Calculating the spectra similarity score

L1_20201120_OPP_JChinTackle_split	formula1	ChargeAgent1	ind1	L2_20201117_OPP_ChinHipsplit	formula2	ChargeAgent2	ind2	SumOfSquares_orig	SumOfSquares_log	SumOf_ABS_orig	DotProd_orig	DotProd_log	chebychev_orig	chebychev_log	hamming
A1	C11H13N3O	H+	1	К7	C11H13N3O	H+	77	0.959456172	0.955820259	0.645442795	0.959456172	0.95582026	0.796880899	0.849142716	0.7540124
A2	C21H25N3S	H+	2	L1	C21H25N3S	H+	78	0.012769056	0.035536743	0.024073165	0.012769056	0.03553674	0	0	0.63627819
A3	C16H14N4O4S	H+	3	P5	C16H14N4O4S	H+	110	0.376369085	0.521629264	0.431661871	0.376369085	0.52162926	0	0	0.57399973
A4	C15H16CINO3	H+	4	L7	C15H16CINO3	H+	84	0.575450241	0.798650943	0.59424423	0.575450241	0.79865094	0	0	0.57098524
A5	C20H30N2O2	H+	5	E2	C20H30N2O2	H+	30	0.246982531	0.514444398	0.360448312	0.246982531	0.5144444	0	0	0.57540668
A6	C22H28N2O2	H+	6	B5	C22H28N2O2	H+	12	0.934506833	0.934983036	0.673043026	0.934506833	0.93498304	0.710370764	0.815177024	0.69882028
A7	C19H24N2O	H+	7	J5	C19H24N2O	H+	68	0.975984359	0.975293725	0.717011186	0.975984359	0.97529373	0.807628839	0.876087665	0.67489348
B1	C20H20N2OS	H+	8	M4	C20H20N2OS	H+	88	0.942930052	0.923079481	0.579466376	0.942930052	0.92307948	0.760249987	0.803224657	0.61001602
B2	C18H20O2	H+	9	A5	C18H20O2	H+	5	0.629982577	0.806256346	0.571039906	0.629982577	0.80625635	0	0	0.58906072
B3	C19H26N2O	H+	10	A3	C19H26N2O	H+	3	0.97830618	0.969743791	0.668430451	0.97830618	0.96974379	0.865575518	0.880806797	0.61427126
B4	C28H31N7O	H+	11	N6	C28H31N7O	H+	97	0.455716805	0.55414278	0.388890379	0.455716805	0.55414278	0	0.057437431	0.56637181
B5	C19H24N4O3	H+	12	M3	C19H24N4O3	H+	87	0.720422371	0.737543874	0.595474874	0.720422371	0.73754387	0.02396553	0.325962775	0.58712939
B6	C22H18N4O3	H+	13	N2	C22H18N4O3	H+	93	0.784752158	0.801449044	0.491496274	0.784752158	0.80144904	0.285153155	0.57117613	0.58073362
B7	C29H31N7O2	H+	14	G1	C29H31N7O2	H+	43	0.513106546	0.689543154	0.458382234	0.513106546	0.68954315	0	0	0.57883667
C1	C16H23NO2	H+	15	O2	C16H23NO2	H+	100	0.600978976	0.677178713	0.508897909	0.600978976	0.67717871	0	0.308279677	0.57361002
C2	C17H18N4	H+	16	D4	C17H18N4	H+	25	0.953061801	0.942710098	0.601272225	0.953061801	0.9427101	0.803503805	0.86625866	0.61351478
C3	C20H24N2	H+	17	M7	C20H24N2	H+	91	0.961009919	0.965121483	0.700107091	0.961009919	0.96512148	0.782007216	0.863901431	0.70986386
C4	C21H23N3OS2	H+	18	H2	C21H23N3OS2	H+	51	5.56523E-05	0.000763311	0.000416035	5.56523E-05	0.00076331	0	0	0.58622389
C5	C14H18N2O	H+	19	P6	C14H18N2O	H+	111	0.800553535	0.772565728	0.518654281	0.800553535	0.77256573	0.424909407	0.548170368	0.60166026
C6	C16H14N2O3	H+	20	K5	C16H14N2O3	H+	75	0.665656906	0.768897344	0.604907656	0.665656906	0.76889734	0	0	0.58993183
C7	C18H23N3O2	H+	21	L6	C18H23N3O2	H+	83	0.639113904	0.758830027	0.499500931	0.639113904	0.75883003	0	0	0.56209077
D1	C12H14N2O3	H+	22	J2	C12H14N2O3	H+	65	0.667574896	0.774567625	0.550112697	0.667574896	0.77456763	0	0	0.57328049
D2	C12H12F3N3S	H+	23	17	C12H12F3N3S	H+	63	0.74729223	0.729210257	0.449139073	0.74729223	0.72921026	0	0	0.60106138
D3	C25H28N2O3	H+	24	O3	C25H28N2O3	H+	101	0.008131822	0.03555204	0.019647503	0.008131822	0.03555204	0	0	0.62086475
D4	C13H16N4O	H+	25	H7	C13H16N4O	H+	56	0.948149344	0.929838929	0.61596809	0.948149344	0.92983893	0.759258155	0.800922094	0.63014892
D5	C23H34N2O2	H+	26	H4	C23H34N2O2	H+	53	0.959585554	0.950042239	0.681816781	0.959585554	0.95004224	0.812018621	0.828515384	0.61561804
D6	C21H23NO2	H+	27	J4	C21H23NO2	H+	67	0.955406317	0.941799395	0.591993253	0.955406317	0.9417994	0.829509871	0.855498614	0.62816314
E1	C16H10N2OS	H+	28	F2	C16H10N2OS	H+	37	0.794237262	0.775753339	0.471712523	0.794237262	0.77575334	0.41179044	0.583293401	0.58171648

© 2021 DH Tech. Dev. Pte. Ltd.

#### Spectra comparison



- Background removing feature enabled
- UI for data review
- PCA analysis included for direct visualization





**Echo** MS System

- AEMS enables the high-throughput compound QC
- The automated data processing workflow enabling
  - Target ion integration
  - Data visualization
  - Background processing
  - Spectra comparison across data set



The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Beckman Coulter® is being used under license.

Echo and Echo MS are trademarks or registered trademarks of Labcyte, Inc. in the United States and other countries, and are being used under license.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see <a href="https://www.sciex.com/trademarks">www.sciex.com/trademarks</a>).

© 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-11-13981-A

