

//// Acoustic Ejection Mass Spectrometry (AEMS) / Differential Ion Mobility for Ultra-High Throughput Mass Spectrometry

Leigh Bedford; Brad Schneider; Eva Duchoslav; Aaron Stella; Subhasish Purkayastha; Chang Liu; Thomas R Covey SCIEX, 71 Four Valley Drive, Concord, ON, L4K 4V8 Canada

INTRODUCTION

Acoustic Ejection Mass Spectrometry (AEMS) is a form of flow injection analysis that achieves sampling rates of greater than 1 Hz and is applicable for quantitative high throughput assays. HPLC is by-passed but the need for chemically based separations to distinguish isobars remains. Low field ion mobility has been suggested but most require pulsed ion beams which introduces a duty cycle loss when coupled to MRM acquisition for quantitation. Mobility measurements must be in the low milliseconds to keep pace. Differential mobility spectrometry (DMS) is unique among forms of mobility as it has all three of the key attributes, millisecond flight times, continuous ion beam operation, and chemically based separations.

MATERIALS AND METHODS

Experiments conducted on a triple quadrupole (QqQ) with the DMS analyzer at atmospheric pressure between the pneumatically-assisted ESI and vacuum aperture as shown in Figure 1. An open port interface (OPI) was positioned over the microtiter plate sample wells and acoustic transducer that dispensed 2.5 nL samples into the OPI. Isobaric opioid and benzodiazepine drugs were analyzed, including morphine, hydromorphone, norhydrocodone, noroxycodone, oxymorphone, dihydrocodeine, mirtazapine, desmethyldoxepin, 7aminoflunitrazepam, diazepam, 7-aminoclonazepam, oxazepam, chlordiazepoxide, temazepam, olanzapine, desmethylclozapine, flunitrazepam, amoxapine, clonazepam, midazolam, clozapine, fentanyl, and norfentanyl Compounds were spiked in human urine and analyzed directly or after desalting. DMS separations were optimized by adding isopropanol, acetonitrile, or ethylacetate chemical modifiers to the transport gas to induce the cluster-decluster mechanism which separates molecules based on their unique gas phase ion chemistry.





RESULTS

The first step in a high-throughput workflow that utilizes modifier based DMS separations is developing a means to determine optimal separation parameters for the analytes in an automated fashion. For each drug a tuning signal could be generated using an acoustic droplet injection frequency of 20 Hz which consumed sample at 1 µL/min. Figure 2 shows the difference between running the system at a frequency of 10Hz (left side of A) which results in peak widths on the order of 1s for pseudo continuous mode, and 400Hz (right side of A) which results in ejection of a series of 6 droplets over a time period of approximately 2.5 ms to generate a discrete peak comprising 15 nL of analyte diluted in the transport flow stream.



Figure 2. A) Switching between pseudo-continuous signal mode (left) and pulsed mode (right). The pseudocontinuous mode provides a stable signal for tuning purposes such as ramping the compensation voltage to determine optimal separation parameters (B). Switching to pulsed mode enables running samples at rates up approximately 1 Hz to generate. C) Pulsed data acqu using the optimal CoV determined from (B).







DMS TRANSMISSION

Figure 3. AEMS data for reserpine ions taken with (black traces) and without DMS installed (red traces). The bottom pane shows a zoom of the baseline.

to Jired		2.40 2	2.45 2.50	2.55
Ana	lyte	Count Rate with DMS Installed (cps)	Count Rate with DMS Removed (cps)	Ratio of Area with DMS Off/On
Reser	rpine	70,200	90,200	1.3
Minc	oxidil	43,000	125,000	2.9
Pro	line	5,680	18,900	3.3

-20 -18 -16 -14 -12 -10 -8

0.8 -

0.6 -

0.4 -

0.2 -

0.8

0.6 -

0.4 –

0.2 -

-Fluorouracil

Faurocholic Acid

 Table 1.
 Table of transmission data for 5 different
compounds taken with the DMS installed on the MS system and with the DMS removed from the MS system. Typical signal reductions ranged from 1.3 – 3.3X for these samples.

3,100,000

3,530,000

5,090,000

1.6

The AEMS process samples directly from wells of a sample plate and therefore does not include provisions for LC separation. As a result, the general applicability of this approach is limited by the presence of chemical interferences. DMS separation can substitute for a traditional LC separation on a time-scale that is compatible with AEMS analysis (\sim 1 sec per sample).



Figure 4A. Analytical signal for measured in desmethyldoxepin MRM channel (DMS off, bottom pane).



The elimination of isobaric interferences for these species enables accurate quantitation for either species in the presence of the other. Figure 6A and B shows examples of quantitation of mirtazapine (A) and desmethyldoxepin (B) from a desalted urine matrix without an internal standard, using an injection volume of 50 nL.



Figure 6A. Calibration curve for mirtazapine in desalted urine at an injection volume of 50 nL.

Linear calibration curves were generated from a total of 10 replicates for the blanks and desalted urine samples. High quality data were generated with LOQs of 0.6 ng/mL and 0.4 ng/mL for mirtazapine & desmethyldoxepin, respectively.





Figure 7. Group of 5 compounds with m/z 313-316 with interferences.



Figure 4B. Analytical signal for nirtazapine (top pane) and interference mirtazapine (top pane) and interference measured in desmethyldoxepin MRM channel (DMS on, bottom pane)

> desmethyldoxepin (bottom pane) mirtazapine MRM channel (DMS

Mirtazapine (blue) and Desmethyldoxepin (red).





Y = 605.9X + 1299 R^{*} = 0.9994 Concentration (ng/mL)

Figure 6B. Calibration curve for desmethyldoxepin in desalted urine at an injection volume of 50 nL.







Figure 9A. AEMS data for injection of olanzapine standard (top pane) in the absence of DMS or LC separation.

CONCLUSIONS

High throughput screening approaches that target analysis of tens of thousands of samples per day with throughput on the order of 1 Hz or greater are not presently compatible with liquid chromatography separations. The elimination of LC significantly reduces analytical selectivity, substantially increasing the likelihood of isobaric interferences in a tandem mass spectrometry measurement from species endogenous or exogenous to the sample matrix. While its peak capacity is generally lower than LC approaches, DMS presents an alternative method to augment selectivity for high throughput screening. Limits of quantitation in the ng/mL range, CVs < 10%, and three orders of magnitude calibration curve linearity were achieved for all compounds including those that would otherwise interfere with each other's MRM transitions.

TRADEMARKS/LICENSING

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-10-13996-A



Figure 9B. AEMS data for olanzapine for the group of five with a DMS installed for separation.