

# Comprehensive site-specific glycan profiling of a protein-based flu vaccine using electron activated dissociation (EAD)

## Featuring alternative fragmentation and streamlined data analysis software

Andy Mahan<sup>1</sup>, Hirsh Nanda<sup>1</sup> and Zoe Zhang<sup>2</sup>

<sup>1</sup>Johnson & Johnson Pharmaceutical Research & Development, Spring House, PA, USA; <sup>2</sup>SCIEX, Redwood City, CA, USA

### ABSTRACT

In this work a comprehensive glycopeptide analysis of a glycoprotein vaccine is presented, using novel electron activation dissociation (EAD)<sup>1,2</sup> as part of an automated data-dependent acquisition (DDA) workflow. Data reduction was accomplished with Protein Metrics Inc. software providing extensive glycosylation libraries and highly accurate results.

### Introduction

Glycoproteins are an essential modality in vaccine development, with many approved applications to date. Glycosylation, including glycan structures and occupancy, plays a vital role in the stability and efficacy of protein-based vaccines. A peptide mapping workflow is generally used for glycopeptide characterization, as it provides protein sequence information and glycan profiling, simultaneously. Traditional collision-induced dissociation (CID) preferentially dissociates glycans and often fails to provide accurate confirmation of peptide sequences and localization information. Alternative fragmentation can provide more descriptive fragmentation, however, long reaction times and the need for peptide-based tuning has limited its usage as a routine characterization method. Here, the ability of latest ECD technology to characterize a vaccine protein with multiple glycosylation sites was examined.

### MATERIALS AND METHOD

**Sample preparation:** A sample of glycoprotein-based vaccine, expressed in a CHO cell line of 25 kDa containing 6 N-linked glycosylation sites was denaturated with 7.2 M guanidine hydrochloride, 100 mM Tris buffer pH 7.2, followed by reduction with 10 mM DL-dithiothreitol and alkylation with 30 mM iodoacetamide. Digestion was performed with trypsin/Lys-C enzyme at 37°C for 16 h.

**Chromatography:** A total of 3 µL (4 µg) of the Lys-C digest were separated with a Waters ACQUITY CSH C18 column (2.1 × 100 mm, 1.7 µm, 130 Å) using an ExionLC AD system (SCIEX). The mobile phase consisted of water with 0.1% formic acid, while the organic phase was acetonitrile with 0.1% formic acid. A 60 min gradient profile was used at a flow rate of 250 µL/min. The column temperature was maintained at 50°C.

**Mass spectrometry:** Data were acquired using DDA with the ZenoTOF 7600 system (Figure 1). The electron energy for the alternative fragmentation in the EAD cell was set to a value of 7 eV.

**Data processing:** Data were processed in Byos software (Protein Metrics Inc.). To achieve side by side comparison, the standard PTM workflow was modified to include two MS/MS identification using the Byonic processing nodes, one for CID data processing, one for EAD data processing. Peptide identification and fragment mass tolerance were set at 6 ppm and 20 ppm, respectively.






Figure 1. The ZenoTOF 7600 system from SCIEX.

### RESULTS & DISCUSSION



**Figure 2. Comparison of MS/MS spectra of a glycopeptide using Zeno EAD and Zeno CID.** MS/MS spectra are shown for a glycopeptide carrying G0 at position N38. Blue and red hash marks depict fragment ion coverage. In addition to higher fragment coverage with better signal-to-noise ratio (S/N), EAD provided unambiguous information on the localization of the glycosylation through comprehensive ion series with intact glycosylation (such as, c11++ onwards). CID resulted in high-abundance oxonium ions and the loss of the fragile modification, limiting confident localization of the modification.



**Figure 3. Comparison of MS/MS spectra of glycosylated and aglycosylated peptide using Zeno EAD.** MS/MS spectra are shown for a glycopeptide carrying HexNac(5)Hex(5)NeuAc(1) and its aglycosylated counterpart at position N38. The red circle highlight exemplary c-ions that differ between the two peptides based on the presence or absence of the modification. Blue and red hash marks depict fragment ion coverage with all N-terminal fragment ions in blue and the C-terminal ions in red.

