

ABSTRACT

INTRODUCTION

Gain more useful MS/MS information in every experiment, particularly on lower abundance species, with a 5-20x gain in MS/MS sensitivity coupled with either EAD or CID fragmentation

New

Zeno trap

New

Wide dynamic range
5GHz, 10bit ADC with 40GHz TOF timing with 25 psec detection rate
High speed pulse counting to maintain resolution and mass accuracy up to 133kHz and over 6 orders LOD

New

New Q0 design

Improved ion transmission and maintenance


New

Electron activated dissociation
Fine tune the fragmentation energies specific to your molecules of interest and acquire key MS/MS features for large molecules, peptides and lipids and small molecules.

SCIEX
The Power of Precision

© 2021 SCIEX Tech. Div. Pte. Ltd.

- Increases duty cycle to 90% at MS/MS level
 - Rich MS/MS spectra



MATERIALS AND METHODS

Data processing: SWATH DIA and Zeno SWATH DIA data were processed using DIA-NN software² (v. 1.8.1 beta 11) using (a) an in-house generated pH fractionated spectral library of HeLa and K562 cell lines created with ProteinPilot in OneOmics suite (cloud, searched against SwissProt-Human canonical+isoform FASTA file) using the 'Per-File' search strategy followed by fraction level retention time alignment and merging using the Extractor application in OneOmics suite (11269 protein groups and 169395 peptide precursors) or (b) a SwissProt-Human canonical+isoform (Jan. 2021) FASTA database for library-free searches. The pg.matrix.tsv and pr.matrix.tsv reports were used for reporting protein groups and precursors and for calculating identifications at %CV thresholds. Zeno DDA files were processed with the ProteinPilot application search engine in OneOmics suite (cloud) against a SwissProt-Human canonical+isoform FASTA database. Iodoacetamide was selected as a fixed modification and results were reported at 1% FDR.

RESULTS

Protein groups identified (black bars)
Protein groups 20% CV (pink bars)
Protein groups 10% CV (teal bars)

Annotations:

- 3.3x more ID**
5x more quantified at 20% CV with Zeno SWATH DIA
- 2.9x more ID**
4x more quantified at 20% CV with Zeno SWATH DIA
- 1.7x more ID**
2x more quantified at 20% CV with Zeno SWATH DIA

95% at CV < 20%
82% at CV < 10%

Concentration	Protein groups identified	Protein groups 20% CV	Protein groups 10% CV
Blank	0	0	0
0.25 ng	1003	428	176
0.5 ng SWATH DIA	449	138	55
0.5 ng	1481	690	344
1 ng SWATH DIA	795	322	158
1 ng	2341	1282	661
1 ng FASTA	1689	1000	591
5 ng SWATH DIA	2409	1298	690
5 ng	4215	2708	1692
5 ng FASTA	3282	2262	1430
10 ng	5021	3595	2389
25 ng	6138	5105	3785
50 ng	6344	5682	4486
200 ng	7079	6690	5774
200 ng FASTA	6909	6628	5838

Figure 1 consists of two bar charts, (a) and (b), showing the number of protein groups and precursors identified in the 1000 most abundant protein groups across three different conditions: 0.25 ng, 0.5 ng, and 1 ng. The y-axis for both charts is labeled 'protein groups/precursors' and ranges from 0 to 14000. The x-axis for both charts is labeled with the sample amounts: 0.25 ng, 0.5 ng, and 1 ng. The legend indicates that blue bars represent protein groups and red bars represent precursors. In chart (a), the number of protein groups is approximately 1000 for 0.25 ng, 1500 for 0.5 ng, and 2000 for 1 ng. In chart (b), the number of precursors is approximately 3000 for 0.25 ng, 5500 for 0.5 ng, and 10500 for 1 ng. Error bars are shown for each bar.

Sample Amount	Protein Groups (a)	Precursors (b)
0.25 ng	~1000	~3000
0.5 ng	~1500	~5500
1 ng	~2000	~10500

Venn diagram illustrating the overlap of protein identifications between DIA-FASTA and DIA-SpecLib. The diagram shows two overlapping circles: a blue circle for DIA-FASTA and a red circle for DIA-SpecLib. The numbers within the circles are: 142 (DIA-FASTA only), 1425 (Intersection), and 703 (DIA-SpecLib only).

Figure 1: Proteomic analysis of the 200 ng and 400 ng fractions.

Left Panel: Protein CV Distribution

Violin plots showing the distribution of protein CV for two fractions: 250 pg K562 (green) and 1000 pg K562 (blue). The y-axis represents Percent CV (0 to 130). The median CV is 19.9% for the 250 pg K562 fraction and 16.3% for the 1000 pg K562 fraction.

Right Panel: Protein Groups and Peptides Identified

Bar charts showing the number of protein groups and peptides identified for the 200 ng and 400 ng fractions. The left y-axis represents Protein groups (0 to 6000), and the right y-axis represents peptides (0 to 70000).

Fraction	Protein groups identified	Peptides identified
200 ng	4580	43015
400 ng	5134	55949

Figure 4. Protein and peptide identifications with Zeno DDA for 200 and 400 ng K562, using <1% FDR filter.

CONCLUSIONS

- Zeno trap improves duty cycle at MS/MS level to over 90%, thereby improving MS/MS sensitivity 5-10x
- At single-cell level protein loads, Zeno SWATH DIA increases protein detection over 3x and increases number of proteins quantified by 4-6x relative to SWATH DIA (Figures 1a, 1b). Consistent detection of low loads was achieved across 30 days:
 - 0.25 ng load: 900-1100 protein groups (3 sets over 1 month, Figure 1c)
 - 0.5 ng load: 1400-1500 protein groups (3 sets over 1 month, Figure 1c)
 - 1.0 ng load: 2100-2300 protein groups (3 sets over 1 month, Figure 1c)
 - 0.25 and 1 ng K562 loads had median CVs of 20% and 16%, respectively (Figure 3)
- For a 200 ng protein load, Zeno SWATH DIA enabled detection of over 7000 protein groups with 95% of detections having a CV < 20% (Figure 1a) and over 66,000 precursors (Figure 1b)
- Library-free (FASTA database) search performed well and eliminated the need to generate spectral libraries for different biological cohorts (Figures 1a, 1d)
- Excellent peptide spectral matching at 0.25 ng and corresponding fragment peptide and protein intensity up to 20 ng (Figure 2a) were observed with median CVs of 20% and 16% for 0.25 and 1.0 ng loads, respectively (Figure 3)
- In Zeno DDA, over 4500 protein groups with 43,000 peptides were identified for a 200 ng K562 load and over 5100 protein groups with 56,000 peptides were identified for a 400 ng load (Figure 4)

REFERENCES

- ## TRADEMARKS/LICENSING

