

Analysis of post-translational modifications using fast electron-activated dissociation (EAD)



Alexandra Antonoplis¹, Nick Morrice¹, and Christie Hunter²
¹SCIEX, USA; ²SCIEX, UK

INTRODUCTION

Characterization of human cellular proteomes enables identification of disease biomarkers and new therapeutic targets. Mass spectrometry-based techniques, such as bottom-up proteomics, are widely used to analyze cellular proteomes. For analysis of protein digests for protein identification, data-dependent acquisition (DDA or IDA) is often used. The ZenoTOF 7600 system, a QTOF platform, uses the Zeno trap to increase duty cycle to ≥90% across the entire fragment ion mass range, resulting in signal gains of 4- to 25-fold in MS/MS mode.¹ These gains enable up to a 40% improvements in proteins identified using DDA workflows.²

For proteomics experiments, there are two fragmentation options available on the ZenoTOF 7600 system: collision induced dissociation (CID) and electron activated dissociation (EAD). Unlike CID, which relies on the collision of ions with nitrogen gas to fragment molecules at their weakest bonds, EAD involves the capture of electrons by molecular ions to form a radical state that then fragments. CID and EAD result in different fragments for the same peptide, with CID yielding y and b ions and EAD yielding additional c' and z+ ions. EAD fragmentation provides complementary sequence information to CID and often preserves post-translational modifications (PTMs) that undergo neutral loss in a CID experiment.^{3,4} This work focuses on DDA method development and analysis of HeLa digest using EAD, with an emphasis on site-specific localization of phosphorylation sites. In addition to phosphorylated peptides, glycosylated human peptides were also investigated using EAD.

MATERIALS AND METHODS

Sample preparation: A 100 µg sample of digested HeLa cell lysate was fractionated into 44 samples using high pH RP-HPLC, as previously described.⁵ Then, 10% of each fraction was injected for LC-MS/MS analysis. Human plasma was used as the source of human serum albumin (HSA) for glycation modification analysis.

Chromatography: HeLa peptide fractions were separated using a Waters ACQUITY UPLC M-class system with a flow rate of 6 µL/min in direct inject mode. A Phenomenex Kinetex 2.6 µm XB-C18 LC column (100 Å, 150 x 0.3 mm) was used with a rapid linear gradient from 5 – 30% B over 21 minutes to interrogate the fractions. Mobile phase A was water with 0.1% formic acid and mobile phase B was acetonitrile with 0.1% formic acid.

Mass spectrometry: A SCIEX ZenoTOF 7600 system equipped with the OptiFlow Turbo V ion source⁶ using a low microflow probe and electrode was used for all data acquisition. DDA parameters that were implemented for all experiments included a TOF MS accumulation time of 250 ms and an exclusion time of 6 seconds. Only precursors with charge states in the range 2-5 with intensities greater than 100 cps were selected for fragmentation. For CID acquisition, the maximum number of candidate ions per cycle was 45 and the accumulation time was 20 ms. For EAD analysis of HeLa digest, the maximum number of candidate ions was 10 per cycle and an accumulation time of 50 ms was used. The reaction time for EAD was 20 ms was used, with an electron beam current of 3000 nA with a KE of 0 eV. For analysis of glycosylated peptides, 20 candidate ions were used per cycle with a 20 ms reaction time and 25 ms total acquisition time.

Data processing: MASCOT software was used for all data processing. HeLa digest files were searched against the Swissprot_Human database using the SCIEX EAD algorithm for EAD and the ESI-QUAD-TOF algorithm for CID. For both phosphorylated and glycosylated peptides, high quality figures were generated using Bio Tool Kit in SCIEX OS software.⁷



To receive a copy of this poster:

- Scan the code with your phone camera
- Complete the form



Proteomic coverage of HeLa digest obtained using EAD

In order to select the acquisition parameters for Zeno EAD DDA analysis of HeLa digest, a variety of parameters were tested, including electron beam current, electron kinetic energy, reaction time and accumulation time. For peptides, it was determined that using an electron KE of 0 eV and a beam current of 3000 nA promoted fragmentation and yielded high fragment ion intensities. To evaluate the performance of EAD for protein identification, the same fractionated HeLa cell lysate was analyzed with data-dependent analysis using both CID and EAD fragmentation for MS/MS. The analysis identified 52,905 peptides at 1% FDR (Figure 1) using EAD. Fewer total numbers of peptides were identified using CID due to the slightly slower acquisition speed (2.5x slower). However, the orthogonal fragmentation mode provided identifications of peptides not found using CID, yielding 11% more peptides. This highlights that while slower, the use of EAD is easily fast enough to be compatible with data-dependent acquisition. An example peptide is shown in which nearly complete c and z+ (z+1) ion series were detected, with CID fragmentation yielding only partial y and b ion series (Figure 2).




Figure 1. Comparison of peptide identifications at 1% FDR as a function of fragmentation mode on the ZenoTOF 7600 system. Due to the higher frequency of sampling with CID, peptide identifications were approximately 2-fold higher than with the EAD fragmentation mode, but EAD provided over 10,000 new peptide identifications.





Figure 2. Comparing CID and EAD fragmentation of a non-phosphorylated peptide. (Top) CID fragmentation yielded incomplete y and b ion series for the LHLASGADVTSANAKAK peptide. (Bottom) Nearly complete z+ (z+1) in table) and c ion series were detected for the same peptide when fragmented by EAD. Detected target fragment ion masses (+1 charge state) are shown in bold red in the tables below each associated spectrum.

Phosphorylation site localization using Zeno EAD MS/MS spectra

Overall, 52,905 peptides were identified at 1% FDR from the EAD data, including 158 phosphopeptides. The sample preparation for the digest was not optimized for preservation of modifications and no enrichment step was performed. EAD spectra from phosphopeptides were used to identify the sequence and site of phosphorylation in peptides, including peptides with multiple serines (Figure 3). Other challenging phosphorylation sites were located with EAD spectra, including on a threonine residue near the peptide N-terminus (Figure 4) and on a peptide with 17 amino acids and a molecular weight of 1889.88 Da (Figure 5). The associated c and z+ ion series clearly detailed the location of the phosphorylation sites on the peptides. EAD fragmentation provided a means to identify challenging phosphorylation sites otherwise unattainable with CID.

