



## Improved compound identification in non-target screening

**David Cox, Craig Butt, Janna Anichina, Adrian Taylor | SCIEX**

ASMS 2022

RUO-MKT-11-14749-A © 2022 DH Tech. Dev. Pte. Ltd.



**Searching chemical databases** by m/z or molecular formula often generates a **very long list** of possible compounds. The **correct answer** is often the one that is **easiest to read**.

# Compound ID difficulty: targeted, suspect, unknown

## Targeted

Compound is in library, has been measured by this lab

- Mass error
- Retention time
- Ion ratio
- Fragment mass error
- Isotope match
- Spectral library match

## Suspect

Compound is in library, has **never** been measured by this lab before

- Mass error
- Fragment mass error
- Isotope match
- Spectral library match

## Unknown


Compound is in PubChem database, no spectrum ever measured

- Mass error
- Isotope match
- Predict formula



## General Unknown

Not in PubChem database

## Finding “real” features




## Identifying a compound



### Propamocarb


Molecular Formula	C <sub>9</sub> H <sub>20</sub> N <sub>2</sub> O <sub>2</sub>
Average mass	188.267 Da
Monoisotopic mass	188.152481 Da
ChemSpider ID	30114

# Feature finding



# What is a “good” peak?


I can't describe it, but I know it when I see it ...




# Using intensity to sort “good” from “bad” peaks

While features with a high intensity are typically real, and features with a low intensity are often junk, where do you draw a cut-off?


1st



8 971





13 394



# What is a “good” peak?

We can use other measures.

But how do we decide on cut-offs? How to keep track of many measures at once?




## Intensity

S:N, Area Height, Region Height, Quality, Retention Time Delta min\_, Total Width, Width at 50 , Baseline Delta Height, Width at 5, Width at 10, Slope of Baseline, Tailing Factor, Asymmetry Factor, Points Across Baseline, Points Across Half Height, AutoPeak Asymmetry, AutoPeak Candidate Model Quality, AutoPeak Group Confidence, AutoPeak Integration Quality, AutoPeak Num Peaks, AutoPeak Peak Width Confidence



# Enabling machine learning with Analytics prototype in SCIEX OS software

- Modified version of processing software
- Used Microsoft's open source ML.NET framework
- Implemented a classifier that uses the Stochastic Dual Coordinate Ascent (SDCA) method
  - Trained using a custom column in the results table, where a user enters classifications ("good" or "bad")
  - Results of classification are output to new columns (predicted classification and score)



<https://livebook.manning.com/book/grokking-machine-learning/appendix-b/v-15/29>

- For more details on ML.NET and machine learning in general:

<https://dotnet.microsoft.com/en-us/apps/machinelearning-ai/ml-dotnet>


<https://rubikscore.net/2021/02/01/machine-learning-with-ml-net-ultimate-guide-to-classification/>

<https://www.kdnuggets.com/2020/05/5-concepts-gradient-descent-cost-function.html>

<https://www.jmlr.org/papers/v14/gonen13a.html>


# Review a few peaks to train the model

With minimal keystrokes ( $\uparrow \downarrow$  g b), classification of a **few hundred** peaks is fast and easy. This set of data is used to train the model so it can be used to automatically classify the remaining ~13 000 features.




# Machine learning classification finds more “real” features

1st




Deciding what to focus our attention on is much easier now. We can still find true features, even at low intensities.


8 971



13 394




# How well does machine learning perform compared to other scoring?




A **randomly** generated number has **no predictive** power. It generates as many true positives as false positives.

A (mythical) perfect scoring algorithm would obtain all true positives before any false positives occur.




# Intensity sorting is better than a random score $\bar{\_}(\text{ツ})\_/\bar{\_}$



For example, classifying peaks above 0.004% relative intensity obtained **90% of the 559 true peaks** in the test set, but generated **222 false positives**.




# Machine learning scoring performs better than intensity scoring





ML score above 50% obtained **90% of the 559 true peaks** in the test set, and only generated **56 false positives**.

## Finding “real” features



## Identifying a compound



### Propamocarb

Molecular Formula	C <sub>9</sub> H <sub>20</sub> N <sub>2</sub> O <sub>2</sub>
Average mass	188.267 Da
Monoisotopic mass	188.152481 Da
ChemSpider ID	30114


# Library searching is common for identification

This only works if the MS/MS for the compound has been seen before, and has been entered into the spectral library you have.






# Searching a formula in online databases



## Sorted by Reference Count

- 1 N-(2-Acetamidoethyl)-4-[5-(chloromethyl)-1,2,4-oxadiazol-3-yl]benzamide
- 2 (4-Chloro-1-methyl-1H-pyrazol-5-yl)[4-(2-furoyl)-1-piperazinyl]methanone
- 3 4-Chloro-N-[(1-ethyl-3-methyl-1H-pyrazol-4-yl)methyl]-3-nitrobenzamide
- 4 N-(3-Chloro-2-methylphenyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)acetamide
- 5 2-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-4-nitrobenzamide
- 6 N-(6-Amino-1-benzyl-2,4-dioxo-1,2,3,4-tetrahydro-5-pyrimidinyl)-2-chloro-N-methylacetamide
- 7 N-(2-Chlorobenzyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)acetamide
- 8 5-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-2-nitrobenzamide
- 9 4-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-3-nitrobenzamide
- 10 4-Chloro-N-[(1-ethyl-3-methyl-1H-pyrazol-4-yl)methyl]-2-nitrobenzamide
- ...



# Searching a formula in online databases



## Sorted by Reference Count

- 1 N-(2-Acetamidoethyl)-4-[5-(chloromethyl)-1,2,4-oxadiazol-3-yl]benzamide
- 2 (4-Chloro-1-methyl-1H-pyrazol-5-yl)[4-(2-furoyl)-1-piperazinyl]methanone
- 3 4-Chloro-N-[(1-ethyl-3-methyl-1H-pyrazol-4-yl)methyl]-3-nitrobenzamide
- 4 N-(3-Chloro-2-methylphenyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)acetamide
- 5 2-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-4-nitrobenzamide
- 6 N-(6-Amino-1-benzyl-2,4-dioxo-1,2,3,4-tetrahydro-5-pyrimidinyl)-2-chloro-N-methylacetamide
- 7 N-(2-Chlorobenzyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)acetamide
- 8 5-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-2-nitrobenzamide
- 9 4-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-3-nitrobenzamide
- 10 4-Chloro-N-[(1-ethyl-3-methyl-1H-pyrazol-4-yl)methyl]-2-nitrobenzamide
- ...
- ...
- 301 4-Chloro-N'-[(Z)-(3,4-dimethoxyphenyl)methylene]-1-methyl-1H-pyrazole-5-carbohydrazide
- 302 2-Oxo-2-[(1,3,5-trimethyl-1H-pyrazol-4-yl)amino]ethyl 6-chloronicotinate
- 303 2-Oxo-2-[(1,3,5-trimethyl-1H-pyrazol-4-yl)amino]ethyl 2-chloronicotinate
- 304 **Cycloxaprid**
- 305 N-(4-Chlorophenyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)propanamide
- 306 N-(2-Chloro-5-nitrobenzyl)-3-ethyl-N-methyl-1H-pyrazole-5-carboxamide
- 307 Ethyl 2-[(5-chloro-3-methyl-1-phenyl-1H-pyrazol-4-yl)carbonyl]hydrazinecarboxylate
- 308 N-[3-Chloro-2-(dimethylamino)phenyl]-2-(2,4-dioxo-3,4-dihydro-1(2H)-pyrimidinyl)acetamide
- 309 N-[3-Chloro-4-(1H-pyrazol-1-yl)phenyl]-N'-(1-hydroxy-2-propanyl)ethanediamide
- 310 4-Chloro-2-[2-[(1,3-dimethyl-1H-pyrazol-5-yl)amino]-2-oxoethoxy]benzamide

# Searching a formula in online databases



Randall Munroe:  
<https://xkcd.com/1334/>

I HATE FEELING DESPERATE ENOUGH TO VISIT  
THE SECOND PAGE OF ~~GOOGLE~~ RESULTS.


*chemical*

## Sorted by Reference Count

- 1 N-(2-Acetamidoethyl)-4-[5-(chloromethyl)-1,2,4-oxadiazol-3-yl]benzamide
- 2 (4-Chloro-1-methyl-1H-pyrazol-5-yl)[4-(2-furoyl)-1-piperazinyl]methanone
- 3 4-Chloro-N-[(1-ethyl-3-methyl-1H-pyrazol-4-yl)methyl]-3-nitrobenzamide
- 4 N-(3-Chloro-2-methylphenyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)acetamide
- 5 2-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-4-nitrobenzamide
- 6 N-(6-Amino-1-benzyl-2,4-dioxo-1,2,3,4-tetrahydro-5-pyrimidinyl)-2-chloro-N-methylacetamide
- 7 N-(2-Chlorobenzyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)acetamide
- 8 5-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-2-nitrobenzamide
- 9 4-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-3-nitrobenzamide
- 10 4-Chloro-N-[(1-ethyl-3-methyl-1H-pyrazol-4-yl)methyl]-2-nitrobenzamide
- ...
- ...
- 301 4-Chloro-N'-[(Z)-(3,4-dimethoxyphenyl)methylene]-1-methyl-1H-pyrazole-5-carbohydrazide
- 302 2-Oxo-2-[(1,3,5-trimethyl-1H-pyrazol-4-yl)amino]ethyl 6-chloronicotinate
- 303 2-Oxo-2-[(1,3,5-trimethyl-1H-pyrazol-4-yl)amino]ethyl 2-chloronicotinate
- 304 **Cycloxaprid**
- 305 N-(4-Chlorophenyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)propanamide
- 306 N-(2-Chloro-5-nitrobenzyl)-3-ethyl-N-methyl-1H-pyrazole-5-carboxamide
- 307 Ethyl 2-[(5-chloro-3-methyl-1-phenyl-1H-pyrazol-4-yl)carbonyl]hydrazinecarboxylate
- 308 N-[3-Chloro-2-(dimethylamino)phenyl]-2-(2,4-dioxo-3,4-dihydro-1(2H)-pyrimidinyl)acetamide
- 309 N-[3-Chloro-4-(1H-pyrazol-1-yl)phenyl]-N'-(1-hydroxy-2-propanyl)ethanediamide
- 310 4-Chloro-2-[2-[(1,3-dimethyl-1H-pyrazol-5-yl)amino]-2-oxoethoxy]benzamide
- ...
- ...

... tell me when you see something  
you can **pronounce**.

# The correct answer is usually the one that is easiest to read



## Sorted by Readability

1	<b>Cycloxaprid</b>
2	(5s,8r)-cycloxaprid
3	avadomide hydrochloride (usan)
4	uracil, 1-[p-[3-(2-chloroethyl)ureido]benzyl]-
5	5-chloro-n-(oxan-2-yloxy)-2-(1h-1,2,4-triazol-1-yl)benzamide
6	8-chloro-4-(2-methyl-1,4-oxazepan-4-yl)-6-nitroquinazoline
7	n-(1-tert-butyl-1h-pyrazol-3-yl)-4-chloro-2-nitrobenzamide
8	1-(5-chloro-2-methoxyphenyl)-3-(6-ethoxyimidin-4-yl)urea
9	1-(5-chloro-2-methoxyphenyl)-3-(2-(6-oxopyridazin-1(6h)-yl)ethyl)urea
10	2-(7-chloro-5-nitro-1h-indazol-1-yl)-1-(piperidin-1-yl)ethan-1-one
...	
...	

## Sorted by Reference Count

1	N-(2-Acetamidoethyl)-4-[5-(chloromethyl)-1,2,4-oxadiazol-3-yl]benzamide
2	(4-Chloro-1-methyl-1H-pyrazol-5-yl)[4-(2-furoyl)-1-piperazinyl]methanone
3	4-Chloro-N-[(1-ethyl-3-methyl-1H-pyrazol-4-yl)methyl]-3-nitrobenzamide
4	N-(3-Chloro-2-methylphenyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)acetamide
5	2-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-4-nitrobenzamide
6	N-(6-Amino-1-benzyl-2,4-dioxo-1,2,3,4-tetrahydro-5-pyrimidinyl)-2-chloro-N-methylacetamide
7	N-(2-Chlorobenzyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)acetamide
8	5-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-2-nitrobenzamide
9	4-Chloro-N-[(1-ethyl-5-methyl-1H-pyrazol-4-yl)methyl]-3-nitrobenzamide
10	4-Chloro-N-[(1-ethyl-3-methyl-1H-pyrazol-4-yl)methyl]-2-nitrobenzamide
...	
...	
301	4-Chloro-N'-[(Z)-(3,4-dimethoxyphenyl)methylene]1-methyl-1H-pyrazole-5-carbohydrazide
302	2-Oxo-2-[(1,3,5-trimethyl-1H-pyrazol-4-yl)amino]ethyl 6chloronicotinate
303	2-Oxo-2-[(1,3,5-trimethyl-1H-pyrazol-4-yl)amino]ethyl 2chloronicotinate
304	<b>Cycloxaprid</b>
305	N-(4-Chlorophenyl)-2-(3,5-dimethyl-4-nitro-1H-pyrazol-1-yl)propanamide
306	N-(2-Chloro-5-nitrobenzyl)-3-ethyl-N-methyl-1H-pyrazole-5-carboxamide
307	Ethyl 2-[(5-chloro-3-methyl-1-phenyl-1H-pyrazol-4-yl)carbonyl]hydrazinecarboxylate
308	N-[3-Chloro-2-(dimethylamino)phenyl]-2-(2,4-dioxo-3,4-dihydro-1(2H)-pyrimidinyl)acetamide
309	N-[3-Chloro-4-(1H-pyrazol-1-yl)phenyl]-N'-(1-hydroxy-2-propanyl)ethanediamide
310	4-Chloro-2-[2-[(1,3-dimethyl-1H-pyrazol-5-yl)amino]-2-oxoethoxy]benzamide
...	
...	


When compounds are sorted by how easy their name is to read, using a common children's book readability score, the correct answer is often in the first few rows.

While sorting by how many references a compound has will often reveal the correct answer, for newly registered pesticides, this does not work. In this example it is row 304, and it is on page 16 of a chemical compound search.


# How do we measure readability?

- Simple algorithms (some dating back to 1923) for estimating how difficult it is to read the text in a book


## Flesch-Kincaid Grade Level



[https://en.wikipedia.org/wiki/Green\\_Eggs\\_and\\_Ham#/media/File:Green\\_Eggs\\_and\\_Ham.jpg](https://en.wikipedia.org/wiki/Green_Eggs_and_Ham#/media/File:Green_Eggs_and_Ham.jpg)



<https://readable.com/readability/flesch-reading-ease-flesch-kincaid-grade-level/>



# Using readability scores on a larger data set



## Customer sample

- 242 compounds relevant to food safety
- 53 of these not in available spectral library



Can the compounds that are **not** in the spectral library be identified?

# Automate the workflow using Python

---

- Pre-calculate
  - Use textstat to score readability of all compound name synonyms from PubChem database
- For a given m/z or predicted chemical formula
  - Find all possible PubChem matches using PubChemPy
- Lookup the readability of the synonyms for all matches
  - Keep anything that is readable (from infant up to my thesis supervisor)
- Search anything remaining in Google along with search terms relevant to the sample
  - Using Google custom search API








# Using reference counting, Isofetamid is found on page 5

## Isofetamid

On the fifth page of results, it finally shows up, buried amongst other structures with similar reference counts

Found 2513 results

Search term: **MF = 'C\_{20}H\_{25}NO\_3S'**

1 2 3 4 5 6 7 8 9 10 ...					
ID	Structure	Molecular Formula	Molecular Weight	# of Data Sources	# of References ▼
<a href="#">13137952</a> 0/1 defined		C <sub>20</sub> H <sub>25</sub> NO <sub>3</sub> S	359.4824	9	10
<a href="#">13140800</a>		C <sub>20</sub> H <sub>25</sub> NO <sub>3</sub> S	359.4824	9	10
<a href="#">13147899</a> 0/1 defined		C <sub>20</sub> H <sub>25</sub> NO <sub>3</sub> S	359.4824	9	10
<b>27473807</b>		C <sub>20</sub> H <sub>25</sub> NO <sub>3</sub> S	359.4824	9	10
					

# Using readability, Isofetamid is successfully identified


## Isofetamid

	A	B	C	D	E	F
1	▼	cid ▼	name ▼	googleHits ▼↑	readability ▼	msmsMatch1 ▼
2	0	71657865	isofetamid	8350	12.6	0.980003843
3	1	58742385		0	300	-1
4	2	117979745		0	300	-1
5	3	25192457	1-(4-ethylphenylsulfonyl)-4-phenylazepan-4-ol	0	146	-1
6	4	25164282		0	300	-1
7	5	25129861		0	300	-1

Using readability on its own, or combined with Google searching or theoretical MS/MS matching, enables successful identification of the correct compound


# More unknowns identified using readability

Identified using **readability** and number of **Google page hits**



Scoring the number of **Google page hits** for only those compound names that are **readable** gave the **correct** identification within the **top 2 compounds**.

Identified using **reference counts**



While using **reference counts** does work for most unknowns, this technique did **not** work for several compounds. For 3 of these examples, the correct identification was **not** in the **top 25 compounds**. Purchasing this many compounds to confirm a **single** compound would be prohibitively **expensive**.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to [www.sciex.com/diagnostics](http://www.sciex.com/diagnostics). All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see [www.sciex.com/trademarks](http://www.sciex.com/trademarks)).

© 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-11-14749-A

**Searching chemical databases** by m/z or molecular formula often generates a **very long list** of possible compounds. The **correct answer** is often the one that is **easiest to read**.



## Questions and answers






Does **readability** have anything to do with **spectral library** matching?



# Readability scoring is almost as effective as library searching



The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to [www.sciex.com/diagnostics](http://www.sciex.com/diagnostics). All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see [www.sciex.com/trademarks](http://www.sciex.com/trademarks)).

© 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-11-14749-A