

Determination of 18 Steroid Hormones in Human Serum Using Rapid Protein Precipitation Method Coupled With Liquid Chromatography-Tandem Mass Spectrometry

Dan Liu; Fengmei Hu; Chao Huang; Xianglong Zhao; Lihai Guo Shanghai AB Sciex Analytical Instrument Trading Co., Ltd., Shanghai, PR China

ABSTRACT

Steroid hormones play an important role in maintaining normal endocrine, regulating sexual function, and immune regulation. The traditional biochemical immunoassay may be inaccurate as its false-positive result. Liquid chromatography-tandem mass spectrometry (LC-MS/MS) has obvious advantages like strong specificity and high sensitivity. There are many methods for the detection of various steroid hormones in serum samples by LC-MS/MS. Liquid-liquid extraction (LLE) and solid-phase extraction (SPE) are the mainly reported pretreatment methods^[1-2], which are complicated, time-consuming, and expensive. Therefore, a rapid assay for the simultaneous and accurate detection of 18 serum steroid hormones was established by using rapid protein precipitation pretreatment method.

INTRODUCTION

Based on SCIEX Triple Quad[™] 6500+ system, a liquid chromatograph-tandem mass spectrometry method was developed for the accurate qualitative and quantitative detection of 18 steroid hormones. Standards of 18 compounds generate standard curves at more than 7 different concentration points. Each compound showed good linear correlation within the range of standard curves. The linearity of the assay was demonstrated by all the R values greater than 0.995. All 18 compounds had low LOQ while meeting the requirements of clinical physiological concentration level. For example, LOQ concentrations were 10 pg/mL for testosterone, aldosterone and androstenedione, and 20 pg/mL for estriol, estrone and dihydrotestosterone.

In order to overcome the disadvantage of using SPE method and LLE method in the pretreatment process, the protein precipitation pretreatment method was developed as it is more convenient and easier to operate and automate.

In the process of method development, the matrix effects and recovery of analytes of different types of protein precipitators were studied. Zinc sulfate($ZnSO_4$) is a cationic precipitator, and its effect on protein precipitation is better than traditional organic reagents. Considering the precipitation efficiency and recovery, the co-treatment method of cationic precipitator ($ZnSO_4$) and organic precipitator (acetonitrile) was confirmed.

This method uses less samples and reagent, which ensures good removal efficiency of interference in complex matrix and obtains better sensitivity at the same time. The most important point is that this method is simple and fast, which only takes 5~10 mins to process a sample, saving at least 3~5 folds of time compared with traditional SPE method.

MATERIALS AND METHODS

Sample Preparation:

Only 50 μ L of serum was needed, protein precipitator(ZnSO₄ and acetonitrile) and internal standard were added, and 18 steroid hormones could be extracted from the sample at the same time through a simple vortex and centrifugation process. (Figure 1)

HPLC Conditions:

A Kinetex C18, 100x3.0 mm, 2.6 µm column at 40° C with a gradient of eluent A 0.5 mM NH₄F aqueous solution and eluent B ethanol was used at a flow rate of 450 µL/min. The injection volume was set to 15 µL. The liquid phase gradient time is 12 minutes.

MS/MS Conditions:

On the SCIEX Triple Quad[™] 6500+ LC-MS/MS system with Turbo V[™] source and Electrospray Ionization (ESI) probe was used. 18 steroid hormones were detected using 2 MRM transitions per compound to ensure quantitation and identification. The positive and negative ions were monitored simultaneously and *Scheduled* MRM[™] mode was used.

RESULTS

Under the developed method, each compound obtains good peak shape and good response. The multi-pair isomers can also achieve baseline separation (Figure 2). The chromatogram of some compounds in human serum samples is shown in Figure 3, with good signal-to-noise ratio and good peak shape.



The results of standard curve and linear range of each substance to be measured are shown in Table1. All compounds have good linearity within their linear range, and the R value is greater than 0.995, which met the relevant requirements. The recovery of all compounds was in the range of 86.8%-113.8% when different concentrations of solution were added to serum from different human (Table 2). The results show the accuracy of this method.



Figure 2. Detection of 18 hormones by LC-MS/MS

Table 1. Standard curve and linear range of 18 hormones				
Name	linear range	linear	R	
17α-Hydroxyprogesterone	0.05 - 50	Y=0.22269x+0.00533	0.9951	
Corticosterone	0.1-100	Y=0.04669x+0.00644	0.9955	
Dihydrotestosterone;DHT	0.02-20	Y=0.09402x+0.00029	0.998	
Pregnenolone	0.02-20	Y=0.04924x+0.00106	0.9987	
Progesterone	0.05-100	Y=0.87480x+0.02038	0.9975	
11-Doxycortisol	0.01-20	Y=0.64005x+0.00247	0.9974	
21-Doxycortisol	0.02-20	Y=3.81310x+0.02050	0.9973	
21-Hydroxyprogesterone	0.02-20	Y=1.52072x+0.01628	0.9969	
Androstenedione	0.01-20	Y=1.71114x+0.00899	0.9968	
Cortisone	0.05-100	Y=3.33320x+0.12860	0.9981	
Testosterone	0.01-20	Y=1.51109x+0.00509	0.9967	
Cortisol	0.05-100	Y=0.44718x+0.01197	0.9971	
DHEA	0.25-100	Y=1.51278x-0.23593	0.9954	
Aldosterone	0.01-20	Y=0.40852x-0.00197	0.9971	
Estriol	0.02-20	Y=0.75346x+0.00389	0.9970	
Estrone	0.02-20	Y=0.22695x+0.00025	0.9985	
DHEAS	2.5-1000	Y=0.19739x+0.18901	0.9974	
Estradiol	0.05-100	Y=0.87181x+0.00495	0.9952	



Figure 3. Chromatogram of hormones in Human Serum Sample (Part)

Table 2. Recovery of 18 hormones were added to different human serum				
Analyte	R-L	R-M	R-H	
17α-Hydroxyprogesterone	105.60%	102.70%	86.80%	
Corticosterone	112.30%	96.90%	87.60%	
Dihydrotestosterone	100.00%	104.70%	100.50%	
Pregnenolone	100.00%	94.00%	105.00%	
Progesterone	109.00%	95.80%	102.60%	
11-Doxycortisol	106.30%	102.30%	100.40%	
21-Doxycortisol	100.00%	104.50%	110.80%	
21-Hydroxyprogesterone	94.70%	87.30%	86.80%	
Androstenedione	104.10%	100.60%	98.30%	
Cortisone	113.20%	88.50%	110.40%	
Testosterone	100.60%	101.00%	97.70%	
Cortisol	105.40%	105.80%	105.50%	
DHEA	91.50%	106.80%	100.80%	
Aldosterone	101.50%	94.10%	93.10%	
Estriol	110.00%	90.00%	102.50%	
Estrone	100.00%	103.40%	112.80%	
DHEAS	113.80%	87.50%	107.70%	
Estradiol	104.80%	90.50%	89.10%	

CONCLUSIONS

A liquid chromatography-tandem mass spectrometry method for the accurate qualitative and quantitative detection of 18 kinds of hormones was established. The method has high sensitivity, good reproducibility and can meet the needs of clinical analysis. Protein precipitation pretreatment method is used for sample purification. Mass spectrometry detection parameters are optimized by SCIEX LC-MS/MS. In summary, this method has the following advantages

- The lower limit of quantitation is pg/mL:
- Pretreatment protein precipitation method is simple and fast, reducing the time of sample pretreatment:
- No need to purchase SPE consumables and equipment, saving cost
- Positive and negative ion switching, Scheduled MRM[™] test mode, simultaneous determination of multiple steroid hormones.

REFERENCES

TRADEMARKS/LICENSING

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures. Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks). © 2023 DH Tech. Dev. Pte. Ltd. MKT-27973-A.



Small sample size, only 50 µL sample can meet the analysis requirements;

Steve S., Paul W.N. L., Andrew L, et al. 2008: Simultaneous Measurement of Serum Testosterone and Dihydrotestosterone by Liquid Chromatography–Tandem Mass Spectrometry. Clin Chem 54, 1855-1863 2 Yongjuan J., Xingli L., Li X., et al. 2021: Liquid chromatography-tandem mass spectrometry measurement of 26 steroid hormones in human serum and plasma samples. J Sep Sci 44, 2358-2370.



To receive a copy of this poster:

- Scan the code with your phone camera
- Complete the form