



Confident characterization and identification of glucuronide metabolites using diagnostic fragments from orthogonal MS/MS data

Rahul Baghla | Manager Scientific Marketing

SCIEX



© 2023 DH Tech. Dev. Pte. Ltd.

Presentation outline

- Challenges with small molecule metabolite identification
- Introduction to instrumentation and electron activated dissociation (EAD) technology
- Diagnostic fragments and confident structure assignment using EAD
- Conclusions
- Acknowledgments

Challenges with identification of small molecule metabolites

Extensive fragmentation coverage for small molecules

- Low abundance of relevant circulating metabolites
- Ionization efficiency in complex matrices
- Qualitative information and a lack of fragmentation coverage with CID
- Phase II conjugations with non-determinative CID fragments to confirm
 - I. Type of conjugation (relies on accurate mass and observed neutral losses)
 - II. Glucuronide position (for example, N-, O-, acyl-)

System hardware advancements




Zeno trap

- Improved MS/MS duty cycle gain of $\geq 90\%$



New Q0 design

- Improved ion transmission and maintenance



Wide dynamic range

10-bit ADC with 40 GHz TDC timing and a 25 psec detection rate

- 5 GHz
- High-speed pulse counting to maintain resolution and mass accuracy up to 133 Hz and a linear dynamic range (LDR) of over 5 orders of magnitude



Complementary fragmentation

- Increased sensitivity using the EAD cell

What is duty cycle?


... And why is duty cycle important?

- What is duty cycle?

- % of ions injected into the TOF
- Typically, ~5-25%
 - Dependent on
 - Fragment mass
 - Scan range upper limit

- Why is duty cycle not 100%?


- Ion losses occur when combining:
 - Pulsed measurement technique
 - TOF
 - Continuous ion beam
 - Quadrupole




Zeno trap

FOR SENSITIVITY GAINS IN MS/MS



- The Zeno trap provides control of the ion beam from the collision cell into the TOF accelerator
- Ions exit the Zeno trap in an ordered release based on potential energy
 - Ions are generally released from a high m/z to a low m/z
 - All ions now arrive in the TOF accelerator at the same time and location




EAD technology



- Free electrons are captured by ions and form a radical state that then fragments
- Electrons introduced with different energies will induce fragmentation in different molecule types
- EAD cell enables you to perform ECD, hot ECD and EIEIO in a single instrument





EAD leads to more confident
structure assignment

Experimental details

Verapamil, buspirone, darunavir and nefazodone metabolism

- **Sample preparation:**

- Verapamil, buspirone, darunavir and nefazodone were incubated at 37°C in rat hepatocytes at a starting concentration of 1 µM

- **ExionLC AD system:**

- Column: Phenomenex Luna Omega Polar C18 column (2.1 x 100 mm, 3 µm, 100 Å)
- Mobile phase A: 0.1% (v/v) formic acid in water
- Mobile phase B: 0.1% (v/v) formic acid in acetonitrile

- **ZenoTOF 7600 system:**

- The samples were analyzed in data dependent acquisition (DDA) mode using Zeno CID DDA and Zeno EAD DDA
- CID conditions: CE of 35 V with CE spread of 15 V
- EAD conditions for verapamil, buspirone and nefazodone:
 - Electron kinetic energy: 10 eV
 - Electron beam current: 5000 nA
- EAD conditions for darunavir:
 - Electron kinetic energy: 11 eV
 - Electron beam current: 8000 nA

Chromatographic gradient.

Time (min)	Mobile phase A (%)	Mobile phase B (%)
0.00	95	5
0.50	95	5
1.50	85	15
3.50	50	50
4.75	5	95
5.75	5	95
5.80	95	5
6.50	95	5


Injection volume: 5 µL

Darunavir-o-glucuronide fragmentation comparison


Alternative fragmentation provides more structural information

Darunavir was incubated at 37°C in rat hepatocytes at a starting concentration of 1µM

Darunavir_t120_Zeno_CID_IDA_5uL - Precursor: 724.3 Da, +1, CE: 35.0



Darunavir_t120_Zeno_EAD_IDA_5uL_KE10 - Precursor: 724.3 Da, +1



Fragments which do not contain glucuronide





Fragments which do contain glucuronide

What does EAD offer?

COMPLEMENTARY AND INFORMATION-RICH STRUCTURE ELUCIDATION WITH EAD

- **CID is a soft, thermal fragmentation technique**
 - Often leads to cleavage of most labile sites
 - Results in few diagnostic fragments
 - Insufficient cleavage without protonation sites
- **EAD offers complementary fragmentation information**
 - Radical dissociation mechanism
 - Can maintain labile modifications
 - Potential to result in many diagnostic fragments



Parent Structure	Structure Candidates
Rank by: EAD	0 % CID 100 %
EAD	100 % CID 0 %

Interpretation pane from Molecule Profiler software

Features of Molecule Profiler software

Potential Metabolites: 12 of 12 Peaks

Report	Peak ID	Name	Formula	Assigned	Neutral Mass	Average Mass	m/z	Charge	ppm	R.T. (min)	Peak Area	% Area	% Score	MS/MS Spectra
1	M207	Glucuronidation [M-H] ⁻	C33H45N3O13S	✓	723.27	723.53	724.2740	1	-0.9	4.04	3.38E+05	0.90	89.1	2
2	M196	Glucuronidation [M-H] ⁻	C33H45N3O13S	✓	723.27	723.58	724.2742	1	-0.5	3.92	5.22E+04	0.14	83.1	2

Interpretation Develop Prepare Options Generate Apply Remove More Selected neutral formula: C33H45N3O13S

EAD and CID MS/MS spectra displaying unique fragments

Assigned: 81 of 90 peaks, score for 81 proposed assignments in total: 1980.0

Use	Mass (m/z)	Ion Formula	Error (ppm)	Intensity (cps)	RDB	Proposed Structures	Score	Fragmentation Type
<input checked="" type="checkbox"/>	576.2363	C28H38N3O8S	-1.9	93.0	12.0	3	16.5	EAD
<input checked="" type="checkbox"/>	612.2197	C27H38N3O...	-3.9	152.0	11.0	13	30.5	EAD
<input checked="" type="checkbox"/>	723.2677	C33H45N3O...	1.3	309.0	13.5	1	31.0	EAD
<input checked="" type="checkbox"/>	724.2780	C33H45N3O...	4.7	8001.0	13.0	1	32.0	EAD

Flexibility to edit structures and assign modified structures

Composition: C33H45N3O13S, Mass: 723.2673, Selected: C33H45N3O13S, Mass: 723.2673

Structure 1 of 3, rank = 1

Proposed formula and structure information for matched fragments

Use	Mass (m/z)	Ion Formula	Error (ppm)	Intensity (cps)	RDB	Proposed Structures	Score	Fragmentation Type
<input checked="" type="checkbox"/>	576.2363	C28H38N3O8S	-1.9	93.0	12.0	3	16.5	EAD
<input checked="" type="checkbox"/>	612.2197	C27H38N3O...	-3.9	152.0	11.0	13	30.5	EAD
<input checked="" type="checkbox"/>	723.2677	C33H45N3O...	1.3	309.0	13.5	1	31.0	EAD
<input checked="" type="checkbox"/>	724.2780	C33H45N3O...	4.7	8001.0	13.0	1	32.0	EAD

Predicted structure ranking with EAD and CID spectra weighting

Rank	Relative Evidence	Apply to Results
1	No structure	<input type="checkbox"/>
2		<input checked="" type="checkbox"/>
3		<input type="checkbox"/>
4		<input type="checkbox"/>

Neutral loss information


Use	Mass	Formula
<input checked="" type="checkbox"/>	112.0582	C6H8O2

- The interpretation pane displays detailed information on both spectra and unique fragments, along with user options for modifying and reassigning metabolites
- EAD spectra show unique fragments (m/z 476.1695, 463.2411, 411.2115 and 525.2081), supporting the identification of o-glucuronide darunavir metabolite

Identification of midazolam N-glucuronide

CID FRAGMENTS ORIGINATED FROM THE PRIMARY MIDAZOLAM STRUCTURE



- The 240-min incubation sample showed a significant peak for midazolam N-glucuronide at a retention time of 4.21-min. Zeno CID DDA did not indicate any specific fragments for midazolam N-glucuronide.
- Instead, fragments from CID originated from the primary midazolam structure. The Mass-Metastite software predicted 2 possible sites of metabolism with CID MS/MS spectra.



Characterization of N-glucuronide midazolam

UNIQUE FRAGMENTS USING EAD ENABLED CONFIDENT IDENTIFICATION

- EAD showed unique fragments at m/z 309.0589 and m/z 354.0804 and confirmed N-glucuronide conjugation on the imidazole ring




N-glucuronide midazolam



Aromatic/aliphatic hydroxylation and o-glucuronide conjugation

EAD ENABLES SITE-SPECIFIC IDENTIFICATION

- Due to the absence of any glucuronide-specific fragments with CID, 4 possible sites of metabolism on the benzene ring of midazolam were predicted
- Site-specific information from EAD enabled the identification of the peak as 1-hydroxymidazolam o-glucuronide




CID



EAD


1-hydroxymidazolam o-glucuronide




Identification of an N-dealkylated midazolam N-glucuronide metabolite

UNIQUE FRAGMENTS USING EAD ENABLED CONFIDENT IDENTIFICATION

- EAD provided rich MS/MS spectra, enabling the identification of an N-dealkylated midazolam N-glucuronide metabolite
- EAD spectra included all fragments generated using CID along with a glucuronide-specific fragment at m/z 386.0702, confirming the site of conjugation




Glucuronidation (tertiary amine) – N Dealkylation



Correlation analysis

Study of drug metabolism over time

- A correlation analysis was performed for drugs and metabolites using peak area based on precursor ion peaks from a TOF MS experiment
- The relative quantitation of drugs and metabolites were correlated and revealed a decrease in drug concentration relative to an increase in different metabolite concentrations over time
- The high mass accuracy enabled the confident prediction of metabolites in an in vitro metabolism study



Metabolite and fragment identification were performed with <10 ppm mass error

Conclusion

- Electron activated dissociation combined with the Zeno trap shows great promise as a new fragmentation technique for delivering information that is complementary to CID
- Using examples such as darunavir and midazolam, n- and o-glucuronide structures along with dealkylation and aromatic/aliphatic hydroxylation glucuronide metabolites were identified with good mass accuracy (MS and MS/MS level)
 - Structures were successfully assigned with unique fragmentation pathways using EAD
- Quick analytical methodology with the software-aided approach for metabolite identification was demonstrated to support high throughput workflows for drug discovery process.

Acknowledgments

- SCIEX team
 - Jason Causon
 - Eshani Nandita
 - Disha Thakkar
 - Eva Duchoslav
- Merck team
 - Kevin Bateman

Thank You

