

# The Application of Research Grade MetabolitePilot<sup>™</sup> Software for the Determination of Exenatide Catabolites using HRAM with SWATH Acquisition

<sup>1</sup> Altasciences Clinical Research, Laval, QC; <sup>2</sup> SCIEX, Concord, ON

## Overview

#### Purpose

To implement research grade MetabolitePilot<sup>™</sup> software for the determination of *in vitro* catabolites of exenatide, a 4.2 kDa GLP-1 agonist.

#### Method

- Incubation of exenatide in rat whole blood at 37°C
- Extraction by protein precipitation in ethanol:acetonitrile (7:3)
- HRAM measurements of SWATH data acquired using a SCIEX TripleTOF5600+

#### Results

The MS/MS spectra of chromatographically unresolved exenatide(3-39) and (4-39) catabolites whose precursor masses were transmitted through the same SWATH window were successfully identified using the advanced spectral deconvolution algorithm in research grade MetabolitePilot<sup>™</sup>.

#### Introduction

The stability of peptide/protein biotherapeutics directly impacts their pharmacokinetic profile, efficacy, and safety, making it essential to characterize potential metabolic soft spots. To facilitate an accurate mass workflow for the confirmation of peptide biotransformations and their profiling across a time course, a research grade version of MetabolitePilot<sup>TM</sup> software has been  $y_4$  ion at m/z 396.2245) were also incorporated into the peak finding algorithm (Figure 1a, b). engineered with an expanded peak finding strategy that (i) supports higher charge states, (ii) generates putative catabolic products by cleavage of the amide backbone and disulfide bonds, (iii) considers isotopic distribution of catabolic products, and (iv) interprets MS/MS data using conventional peptide fragmentation patterns. In the current investigation, MetabolitePilot™ was implemented to determine the *in vitro* catabolites of exenatide, obtained using HRAM in SWATH acquisition mode (Sequential Window Acquisition of All Theoretical MS).



#### Methods

Rat whole blood fortified with exenatide (1  $\mu$ g/mL) and incubated at 37°C was sampled at t<sub>0</sub>, t<sub>30</sub>,  $t_{60}$ ,  $t_{120}$  and  $t_{240}$  and subsequently precipitated with ethanol:acetonitrile (7:3). Extracts were chromatographed on a Halo Peptide ES-C18 column (2.1 x 150 mm, 5 µm) under gradient conditions at 50°C with an acidified water/acetonitrile mobile phase, ramped from 98% aqueous to 70% organic over 30 minutes. Data was acquired with a SCIEX TripleTOF 5600<sup>+</sup> operated in SWATH mode using an accumulation time of 100 ms per experiment, where each experiment coincided to fixed 50 Da wide MS/MS windows at 45 eV collision energy for precursor masses 200 - 1250 Da. A research grade version of SCIEX MetabolitePilot<sup>™</sup> Software was used for post-acquisition processing.

### **Results and Discussion**

The post-acquisition processing workflow in research grade MetabolitePilot<sup>™</sup> initially involved defining potential bio-transformations. In addition to the provided default modifications in the Biologics set, custom bio-transformations can be entered, and for exenatide included deamidation (N/Q), oxidative deamination to alcohol (K), demethylation (A/T) and demethylation + oxidation (A/T). The MetabolitePilot<sup>™</sup> processing algorithm always considers hydrolytic cleavages, and for exenatide, generated 702 theoretical catabolites.

The next stage in the MetabolitePilot<sup>™</sup> workflow defined the Peak Finding Strategy, which leverages the TOF-MS data derived from the first experiment in the SWATH acquisition. Within the Processing Parameters dialog, the peptide sequence of exenatide was entered along with an MS/MS reference spectrum from which two diagnostic product ions ( $y_3$  ion at m/z 299.1717 and

\Exenatide_Hydrolytic.xml - Processing Paran	neters				×
New Open Save	Save As	Delete	Method type: Peptic	les	How Do I?
Compound Information Select	ct From Library	Sequence			
Compound name: Exenatide		1 HC	EGTFTSDL SKQMEEEAVR LFIEWLKNGG	PSSGAPPPS-[Ami]	*
Chemical formula: C184H282N50O60S					
Polarity:   Positive  Negative					
Charge state: From: 2 To: 5					
Ion type: [M+5H] <sup>5+</sup>					
m/z: 837.8127					
					<b>v</b>
	Generic Paramet	ers Compound-Spec	cific Parameters		
Peak Finding Strategy	Catabolites Iso	tope Pattern Produc	t lons and Neutral Losses		
~ TOF MS					
Predicted metabolites	Potential	Hydrolytic Cleavag	88		
Generic peak finding	Max. pept	de bonds to break:	2 Max. cross-links to break: 0	Min. AA count: 3	
Apply mass defect filter	Catabolite	s selected: 702			
Apply charge state filter					
Mass defect		AA Index	Name	Neutral Formula	Neutral Mass
Isotope pattern		1:25-35	WLKNGGPSSGA	C47H72N14O15	1072.530
TOF MSMS	V	1:24-33	EWLKNGGPSS	C47H71N13O16	1073.514
Find characteristic product ions		1:19-26	VRLFIEWL	C54H82N12O11	1074.622
All specified ions		1:26-37	LKNGGPSSGAPP	C46H76N14O16	1080.556
At least 2 ions		1:4-13	GTFTSDLSKQ	C46H74N12O18	1082.524
Find characteristic neutral losses		1:3-12	EGTFTSDLSK	C46H73N11O19	1083.508
<ul> <li>All specified losses</li> </ul>	<b>V</b>	1:23-32	IEWLKNGGPS	C50H77N13O15	1099.566
At least     Insses		1:13-21	QMEEEAVRL	C45H77N13O17S	1103.528
Consider internal neutral losses		1:20-27	RLFIEWLK	C55H85N13O11	1103.649
Isotope pattern (SWATH® Only)					
Save Default Settings Restore Defaults					Save and Close
Save Delauri Settings					Save and Close

Figure 1a. Processing Parameters workflow in research grade MetabolitePilot™. Both TOF-MS and TOF-MS/MS data acquired in SWATH mode were incorporated in the peak finding strategy for exenatide catabolites. For a minimum of three amino acid residues, 702 putative catabolites were proposed based on hydrolytic cleavages alone.

Jeff Plomley<sup>1</sup>, Yi Zhang<sup>2</sup>, Eva Duchoslav<sup>2</sup>, Daniel Villeneuve<sup>1</sup>, Kevork Mekhssian<sup>1</sup> and Anahita Keyhani<sup>1</sup>

# **Results and Discussion (Continued)**

Generic Parameters Compound	-Specific Parameters			
Catabolites Isotope Pattern Pr	roduct lons and Neutra	l Losses		
Reference MS/MS Spectru     wiff file Browse	m © Compound Library	r [		
100 50 0	397.2283 379.1980 299.1717			83
200	400	600		80
	Assign Fragments		_	
Filters	, todigit i raginionito	m/z	Z 1	v3
m/z: From: 100.0 T	o: 1500.0	379.1980	1	, y4-ľ
Charge From: 1 T	īo: 5 🔽	396.2245	1	y4
Only show product 5	%	837.8156	5	
Mass accuracy 10.00	) ppm 🚹	•	_	
		Add prod	luct id	ons,

#### Figure 1b. Compound specific parameters indicating the use of y<sub>3</sub> and y<sub>4</sub> product ions incorporated into the peak finding algorithm.

With the Bio-transformation and Processing Parameters established, incubated samples were interrogated against control samples. Proposed catabolites presented in the Results workspace (Figure 2) were considered only if the measured parent mass was within 10 ppm of theoretical and the response was three-fold greater than that observed in control samples. This same mass accuracy was applied to the assignment of product ions measured against theoretical b- and yion masses for proposed catabolite sequences (Figure 2).



assignment.



Figure 2. The Results workspace of research grade MetabolitePilot<sup>™</sup> highlights proposed catabolite sequences for exenatide including elemental composition, measured *m/z*, parent mass error and chromatographic retention time. The Interpretation pane plots the catabolite MS/MS spectrum correlating fragment ion mass with sequence information, here for exenatide(3-39). The measured mass error for each fragment ion is also reported. In the Results pane, the XIC from a TOF-MS scan can be plotted against control sample and both the isotope pattern and MS/MS spectrum can be displayed and are incorporated into a scoring

# **Results and Discussion (Continued)**

Results from each time-point were compiled in the Correlation workspace and potential Regardless, each minor putative catabolite generated the diagnostic  $y_3$ - and  $y_4$ - ions of the catabolites plotted (Figure 3). In the case of exenatide, only one catabolite demonstrated exenatide reference spectrum, and each parent mass was measured within 5 ppm of theoretical. increased response with incubation time, and coincided to the N-terminal HG clipping Of particular note, exenatide(3-39) and (4-39) were chromatographically unresolved and their biotransformation product exenatide(3-39), whose chromatographic profile and MS/MS spectrum [M+5H]<sup>+5</sup> precursor masses were simultaneously transmitted through the same SWATH window are compared to exenatide in Figure 4. Exenatide(4-39), (5-39), and (7-39) catabolites were also (i.e. *m*/z 749 – 800), thereby generating a mixed MS/MS spectrum. While a formidable challenge detected at incubation times  $\geq t_{120}$ , and therefore could not be fully correlated within the timeto deconvolute this complex scenario in applications such as PeakView, the advanced algorithm used in MetabolitePilot<sup>™</sup> successfully re-constructed the MS/MS spectrum derived from each of frame of the experiment (Figure 3). exenatide(3-39) and (4-39), thereby aligning the TOF-MS XICs with confirmatory b- and yfragment ions, as outlined in Figure 5.



Figure 3. Incubation time profiles for proposed catabolites are generated in the **Correlation workspace.** The major catabolite for exenatide resulted from N-terminal HG clipping whilst exenatide(4-39), (5-39), and (7-39) were only expressed after t<sub>120</sub>. XICs from the TOF-MS scan of the SWATH acquisition can be plotted at each incubation time point as illustrated for exenatide(3-39).



Figure 4. Comparison of chromatographic retention time and MS/MS SWATH spectra for exenatide and its major catabolite, exenatide(3-39) indicating the highly abundant diagnostic  $y_3/y_4$  ions.

# **Results and Discussion (Continued)**



Figure 5. The XIC in PeakView for the  $y_4$  diagnostic ion (*m*/z 396.223) derived from the SWATH MS/MS experiment for m/z 749 – 800 (a) demonstrates a chromatographic profile suggestive of co-eluting catabolites (b). Since the SWATH MS/MS spectrum is derived from a 50 Da mass window, it is difficult to correlate specific fragment ions with precursor mass. However, the advanced MetabolitePilot<sup>™</sup> algorithm could identify the putative catabolites as exenatide(3-39) and exenatide(4-39) despite their co-elution, and moreover could re-construct the MS/MS spectrum derived from each of the parent masses 798.9973 (c), and 773.1894 (d), respectively. With the MS/MS spectra properly de-convoluted, confirmation of the proposed catabolites was possible.

#### Conclusions

The research grade version of MetabolitePilot<sup>™</sup> software, when combined with accurate mass measurements from a TripleTOF platform operated in SWATH acquisition mode, represents a formidable platform for protein biotherapeutic catabolism studies. The advanced capabilities of the software algorithm were exemplified in properly assigning co-eluting catabolites of exenatide whose parent masses were simultaneously transmitted through the same SWATH window. A commercial version of MetabolitePilot<sup>™</sup> (v 2.0) is now available which incorporates the peptide sequencing capabilities presented herein.