

Protein disulfide bond characterization with DMS and ECD



J.C.Yves Le Blanc, Brendon Seale, Lyle Burton and Takashi Baba ; SCIEX, 71 Four Valley Drive, Concord, ON, L4K 4V8 Canada

ABSTRACT

Evaluating a collision cell with dual CID and ECD capability to characterize disulfide bond location in digested proteins.

INTRODUCTION

Identification and characterization of disulfide linkage yields critical information pertaining to biologic compound integrity throughout the entire development. Many of the approaches proposed have been relying on differential MS mapping of reduced and non-reduced digested sample. This is mainly due to the complexity of the MSMS fragmentation by CID, and the lack of automated tools for interpretation. It has been shown that ECD and ETD could yield simplified fragmentation that would improve the confidence in the identification and linkage of peptides [1-3]. When combined with differential mobility (DMS), additional selectivity could be obtained over LC separation [3] for sequence analogues. Here we explored the roles that both of these techniques (ECD and DMS) could play in the characterization of biologic compounds.


MATERIALS AND METHODS

Sample Preparation: Bovine serum albumin (BSA) was digested with trypsin after alkylation but without the reduction step. BSA was used a simple model to generate mixture of intra- and cross-linked peptides. Table 1 give a list of the cross-linked peptides observed and detected with BioPharmaview Software(Sciex).

HPLC Conditions: A Shimadzu Prominence-XR LC system with an Aries XB C18, 100x2.1mm, 2.6µm column (Phenomenex) at 45°C with a gradient of eluent A water/acetonitrile (98/2) + 0.1% formic acid and eluent B water/acetonitrile (5/95) + 0.1% formic acid was used at a flow rate of 300µL/min. Peptide elution was performed with a gradient to 45% of eluent B in 24 minute.

MS/MS Conditions: An Sciex Turbo V™ source and Electrospray Ionization (ESI) probe was used. The mass spectrometer was a research grade hybrid quadrupole time-of-flight system that share many of the ion optic element typically used on a TripleTOF® 5600 LC/MS system. The CID/ECD cell operational details were reported previously [4]. The ECD-enabled quadrupole time-of-flight (Q-ToF) mass spectrometer that can operated in 2 modes: 1) flow thought for CID and ECD and 2) simultaneous trapping of precursor ions with continuous beam of electron. This is achieved with a unique RF ion guide (Figure 1), whose branched structure enables the inter-crossing of an adjustable energy electron beam with the conventional analytical ion beam of the mass spectrometer.

Figure 1. Details of the Chimera-ECD device in terms of location on the ion optics of the QqTOF system



RESULTS

When digested under non-reducing conditions, disulfide linked peptide yield larger charge state peptides are easily separated from others ones with DMS. Using BSA as a model protein, most of the disulfide linked peptides generated have charge state ranging from +3 to +7 (Table 1). Though all peptide originate from a single protein chain, a wide range of intra- and inter- linked disulfide peptides are generated and represent a good pool of peptides to evaluate both ECD and DMS as tools for characterization of proteins. Our attention focused on a particular case where it was not possible to distinguish between 2 possible sequence combination; neither high mass accuracy or CID fragmentation. The isobaric nature of the cases requires alternative tool to distinguish which combination is present in the sample and detected (Figure 1). Since ECD fragmentation yields dominant fragment ions that originate from the cleavage of the S-S bond, complementary sequence information can be used to identify the proper peptide.




Table 1. List of detected peptides containing disulfide links originating from BSA when digestion is performed with alkylation and without reduction (allowed single missed cleavage). A mix of intra- and inter-sequence peptide are generated to evaluate CID and ECD fragmentation.

Theoretical Mono m/z	RT	Observed Mono m/z	Error (PPM)	Score	Charge	XIC Area	Sequence	Disulfide Bonds	Peptide	Chains
445.48431	2.43476	445.84852	0.46627	-	1	112923.277	TCVADESHAGCEK	T7@2(513)-T7@1(162)	T7	1
1347.53038	2.43534	1347.53392	2.62945	-	1	1793.24401	TCVADESHAGCEK	T7@2(513)-T7@1(162)	T7	1
674.26883	2.43456	674.26898	0.22907	-	2	6345.8183	TCVADESHAGCEK	T7@2(513)-T7@1(162)	T7	1
855.6136	5.15146	855.59897	-17.09619	-	5	842.08478	LCVLHEK.CCTESLVNR.RPCFSALTPDETVPKAFD EK	T66@1(475)-T67-68@3(486), T66@2(476)-T63@2(460)	T63;T66;T67-68	1;1;1
921.69767	4.58202	921.69925	1.71198	-	4	1389.98883	LCVLHEK.CCTESLVNR.RPCFSALTPDETVPK	T66@1(475)-T63@2(460)	T63;T66;T67	1;1;1
737.55959	4.57615	737.55869	-1.22106	-	5	8218.68246	LCVLHEK.CCTESLVNR.RPCFSALTPDETVPK	T66@1(475)-T63@2(460)	T63;T66;T67	1;1;1
614.80087	4.57856	614.80025	-1.02101	-	6	11833.4599	LCVLHEK.CCTESLVNR.RPCFSALTPDETVPK	T66@1(475)-T63@2(460)	T63;T66;T67	1;1;1
527.11607	4.57603	527.11559	-0.9268	-	7	3132.89757	LCVLHEK.CCTESLVNR.RPCFSALTPDETVPK	T66@1(475)-T63@2(460)	T63;T66;T67	1;1;1
1228.59447	4.57087	1228.59426	-0.1667	-	3	253.49691	LCVLHEK.CCTESLVNR.RPCFSALTPDETVPK	T66@1(475)-T63@2(460)	T63;T66;T67	1;1;1
954.45561	3.63965	954.45932	3.88686	-	4	408.14896	YICDNQDTISSK.LKECCDKPPLLEK.SHCIAEVEK	T42@3(288)-T40-41@2(277), T39@3(278)-T39-40@3(264)	T39;T40-41;T42	1;1;1
636.6395	3.63616	636.63994	0.69633	-	6	5727.28223	YICDNQDTISSK.LKECCDKPPLLEK.SHCIAEVEK	T41@2(277)-T42@3(288), T41@3(278)-T39-40@3(264)	T39-40;T41;T42	1;1;1
545.8349	3.63562	545.83483	-0.12884	-	7	3232.45392	YICDNQDTISSK.LKECCDKPPLLEK.SHCIAEVEK	T41@2(277)-T42@3(288), T41@3(278)-T39-40@3(264)	T39-40;T41;T42	1;1;1
763.76594	3.63604	763.76601	0.0889	-	5	1882.21442	YICDNQDTISSK.LKECCDKPPLLEK.SHCIAEVEK	T41@2(277)-T42@3(288), T41@3(278)-T39-40@3(264)	T39-40;T41;T42	1;1;1
954.45561	3.63965	954.45932	3.88686	-	4	408.14896	YICDNQDTISSK.LKECCDKPPLLEK.SHCIAEVEK	T41@2(277)-T42@3(288), T41@3(278)-T39-40@3(264)	T39-40;T41;T42	1;1;1
763.76594	3.63604	763.76601	0.0889	-	5	1882.21442	YICDNQDTISSK.LKECCDKPPLLEK.SHCIAEVEK	T41@2(277)-T42@3(288), T41@3(278)-T39-40@3(264)	T39-40;T41;T42	1;1;1
636.6395	3.63616	636.63994	0.69633	-	6	5727.28223	YICDNQDTISSK.LKECCDKPPLLEK.SHCIAEVEK	T41@2(277)-T42@3(288), T41@3(278)-T39-40@3(264)	T39-40;T41;T42	1;1;1
545.8349	3.63562	545.83483	-0.12884	-	7	3232.45392	YICDNQDTISSK.LKECCDKPPLLEK.SHCIAEVEK	T41@2(277)-T42@3(288), T41@3(278)-T39-40@3(264)	T39-40;T41;T42	1;1;1
445.38449	3.96898	445.38407	-0.93754	-	5	7768.94516	CASIQK.EECHGDLLECCADDR	T28@1(199)-T36-37@6(245), T37@2(244)-T37@1(162)	T28;T37	1;1
741.63596	3.96587	741.63639	0.01288	-	3	4276.59977	CASIQK.EECHGDLLECCADDR	T28@1(199)-T36-37@6(245), T37@2(244)-T37@1(162)	T28;T37	1;1
1111.95031	3.96394	1111.95306	2.48089	-	2	513.70399	CASIQK.EECHGDLLECCADDR	T28@1(199)-T36-37@6(245), T37@2(244)-T37@1(162)	T28;T37	1;1
556.47879	3.96729	556.47904	0.45441	-	4	23455.1979	CASIQK.EECHGDLLECCADDR	T28@1(199)-T36-37@6(245), T37@2(244)-T37@1(162)	T28;T37	1;1
863.04339	3.53919	863.04717	4.37561	-	3	3379.29984	CASIQK.VHKECHGDLLECCADDR	T28@1(199)-T36-37@6(245), T37@2(244)-T37@1(162)	T28;T36-37	1;1
647.53436	3.54143	647.53437	0.01288	-	4	14364.5793	CASIQK.VHKECHGDLLECCADDR	T28@1(199)-T36-37@6(245), T37@2(244)-T37@1(162)	T28;T36-37	1;1
518.22895	3.54251	518.22931	0.69678	4.27	5	41674.827	CASIQK.VHKECHGDLLECCADDR	T28@1(199)-T36-37@6(245), T37@2(244)-T37@1(162)	T28;T36-37	1;1
432.02533	3.5418	432.02508	-0.38654	-	6	7641.97721	CASIQK.VHKECHGDLLECCADDR	T28@1(199)-T36-37@6(245), T37@2(244)-T37@1(162)	T28;T36-37	1;1
639.29827	5.48553	639.28913	-14.30031	-	6	876.87475	LKPDNPTL.CDFE.K.VNGVQECCEQAEKDGACLLPK	T12-22@8(167)-T21-22@9(168), T14@9(123)-T21-22@9(168)	T14;T21-22	1;1
958.44377	5.30873	958.4458	2.11756	-	4	507.43848	LKPDNPTL.CDFE.K.VNGVQECCEQAEKDGACLLPK	T12-22@8(167)-T21-22@9(168), T14@9(123)-T21-22@9(168)	T14;T21-22	1;1
766.95647	5.30838	766.95731	1.08323	-	5	1845.59535	LKPDNPTL.CDFE.K.VNGVQECCEQAEKDGACLLPK	T12-22@8(167)-T21-22@9(168), T14@9(123)-T21-22@9(168)	T14;T21-22	1;1




Figure 3. CID and ECD spectra for peptide depicted in Figure 1. The precursor charge state ranged from z=+4 to +7. ECD conditions were kept constant; 3eV with 10ms fill and 30ms reaction time.




Figure 4. ECD spectra of the T7 peptide which contains a single intra-chain disulfide link. The reduction of the S-S bond is detected by cleavage of the backbone of the peptide. Similar to the cross-linked peptide, the higher charge states tend to provide more efficient fragmentation.

Figure 2 shows that ECD fragmentation yields fragmentation information that confirm detection of the T39-40 and T41 peptide in both the M(SH) and M(S) form as doubly and singly charged species. As show in Table 2, the higher charge states (Z=7+ and 6+) yielded fragment ions that confirmed the proper peptide. Additional fragmentation allowed confirmation that the proper peptide detected. From the data, it is possible to extract additional fragmentation information that confirmed the proper linkage S-S linkage (depicted in Figure 2)

	SH for [M+H] ⁺	SH for [M+2H] ²⁺	S for [M+H] ²⁺	2xS for [M] ²⁺	Observed for Z
YICDNQDTISSK	1627.79959	814.40343	813.89952	n.a.	7 & 6
ECCDKPPLLEK	1177.55915	589.28321	589.28321	588.27539	7 & 6
SHCIAEVEK	1015.4877	508.24749	507.74357	n.a.	7 & 6
YICDNQDTISSK	1386.62057	693.81392	693.31001	n.a.	n.d.
LKECCDKPPLLEK	1418.73818	709.87273	709.36881	708.8649	n.d.
SHCIAEVEK	1015.4877	508.24749	507.74357	n.a.	7 & 6

Table 2. Expect list of masses predicted and observed to confirm which of the peptide is present (green highlight). Of note, the higher charge states tend to generates the diagnostic fragment ions that confirm the T39-40;T41;T42.

CONCLUSIONS

ECD fragmentation offers additional complementary information over CID in the case of disulfide containing peptides. When the link is intra-peptide, ring opening is observed and minimal fragmentation information is obtained. When the disulfide link is across peptide originating from different section of the sequence, cleavage of the bond generates individual chains which can yield information about the individual chains as well as some sequence information for each. This provides improved confidence in the location of the disulfide bonds.

Future efforts will concentration of development of comprehensive IDA logic to ensure data collection in automated way to optimize sequence information from the ECD and CID perspective. This future work would also include integration of DMS into the work flow. DMS has demonstrated that multiply charged species can be segregated ahead of mass selection, thus further simplifying the data reduction.

REFERENCES

- Huilin Li, Peter O'Connor, J.Am. Soc. Mass Spectrom. (2012) 23, 2001-2010
- Barbara Ganis, Kathrinn Breuker, Chemistry OPEN (2012) 1, 260-268
- Shiaw-Lin Wu, et al. Anal. Chem. (2009) 81, 112-122
- Takashi Baba, Larry Campbell, J.C.Yves Le Blanc, Bruce Thomson, Anal. Chem. (2014) 87, 785-792

TRADEMARKS/LICENSING

AB Sciex is doing business as SCIEX.
© 2017 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO- MKT-10- 5831