# Towards a Complete Non-Linear Peptide MS/MS Characterization

## Eva Duchoslav; Xu Guo; Tanmaykumar Desai SCIEX, 71 Four Valley Drive, Concord, ON, L4K 4V8 Canada

# **ABSTRACT**

This presentation highlights a generic workflow for characterization of non-linear peptides using high resolution MS/MS data. Non-linear peptides are represented as strings of amino acid residues with links outlining the attachment chemistry; this allows for the prediction of both single-chain and linked fragments. Comprehensive MS/MS annotation is possible by accommodating for side-chain losses as well as multiplycharged fragments and their isotope patterns in MS/MS. With this enhanced approach peptide sequence coverage and percent annotated MS/MS TIC were similar for data from IDA and DIA methods. The annotation coverage increased with the best outcome for purely cyclic peptides, where the improvements of up to 100% from the previous approaches<sup>1,2</sup> were noted.

# INTRODUCTION

Biotherapeutic drugs are typically non-linear, cyclic peptides. These peptides and their catabolites need to be characterized within the drug development process.

While identification and annotation of linear peptides from MS/MS spectra has been well understood, the analysis of large bio-therapeutic peptides, such as cyclic ones, has presented challenges since their fragmentation does not follow conventional fragmentation pathways<sup>3,4,5</sup>. In addition, large bio-therapeutic peptides tend to generate larger, often multiply charged fragments and confident detection of the monoisotopic ion can be difficult due to its potential absence or poor signal to noise. Meanwhile, the traditional naming logic for peptide fragments needs to be adjusted to fit non-linear sequences. This presentation will demonstrate the improvements of peptide annotation strategy addressing the above issues.

# MATERIALS AND METHODS

### Sample Preparation:

Several non-linear peptide standards, oxytocin, cyclo(-GRGDSP), cyclosporine A and insulin (AnaSpec, 34801 Campus Drive, Fremont, CA 94555, USA) were prepared in 0.1% formic acid buffer and analyzed by LC/MS.

## HPLC Conditions

A Shimadzu Prominence LC system with an Agilent ZORBAX Eclipse Plus C18, 50x2.1mm, 1.8µm column at 40° C with a gradient of eluent A water + 0.1% Formic acid and eluent B acetonitrile+ 0.1% Formic acid was used at a flow rate of  $200\mu$ L/min. The injection volume was set to  $5\mu$ L.

### **MS Conditions:**

Data was collected on a SCIEX X500B QTOF System with SCIEX OS 1.2. TOF MS acquisition covered m/z 100 to 1500 with 250ms accumulation time. IDA data were acquired using the peptide workflow using the charge state filter with dynamic background subtraction, intensity threshold 100 cps and maximum candidate ion 5. Data independent acquisition was performed using SWATH<sup>®</sup> Acquisition with 22 fixed fifty Da wide SWATH Acquisition windows. The MS/MS scan from each SWATH Acquisition covered m/z 100 to 1500 with a 100ms accumulation time. Total scan time for the SWATH Acquisition method was 2.5sec. Calibration delivery system (CDS) was used to calibrate the system between injections.

### Data processing:

LC/MS data were processed with a research version of MetabolitePilot<sup>™</sup> software. The sequences of nonlinear peptides were provided in a string representation with the side-chain attachment descriptors. The non-natural amino acids were described in terms of their substructure building blocks, that were handled at the fragment level as variable modifications. In addition, characteristic fragments and neutral losses for amino acids from a peptide dictionary were utilized.

The processing workflow is outlined in Figure 1.

	Pre
_	
	Prepare
	Μ
$\bigcap$	

# RESULTS

Abu: [ BMe]: [1Me]:

+••	FIOF MS/	VIS OT U
Intensity, cps	1200 1000 800 600 400 200 0	100.11

1551	gnea. Z	ie or o r i pea	aks, MolMo Feak Alea Assigned. 60.9%, Sequence Co	verage. If of it consecutive a	amino aci	us											
Fragments: 219 of 996 Proposed Formulae			<b>6</b> 7	m/z	Fragment sequence	type	charge	error (ppm)	intensity								
	Use	Mass (m/z)	Sequence	Ion	Charge	Error (ppm)	Intensity (cps)		397.2807	ALLV	0 b4	1	-0.6	3494.6			
97	1	269.1856	T[_BMe]Abu	4 b2	1	-1.4	671.7		411.2962	AL[1Me]LV	0 b4	1	-0.9	733.4			
8	<b>V</b>	298.2124	ALL	0 b3	1	-0.3	532.3		425.312	AL[1Me]L[1Me]V	0 b4	1	-0.6	2826.9			
0	1	425.3123	ALL[1Me]V[1Me]	0 b4	1	0.3	510.6				- 1 -						
3		128.1066	L[1Me]	1 b1	1	-3.1	462.3		439.3275	AL[1Me]L[1Me]V[1Me]	0 b4	1	-0.9	258.8			
3	<b>V</b>	284.1969	AbuG[1Me]L[1Me]	5 b3	1	0.1	387.6										
7	1	693.4914	ALL[1Me]V[1Me]T[_BMe]Abu	0 b6	<b>≺</b>				CV	clic pentide fragr	nont	type	e indica	ha			
94	<b>V</b>	166.1224	T[_BMe]	4 b1 (- H2O)	1	-1.7	347.2		Cyclic peptide nagment types indicate								
1	-	114.0911	L	1 b1	1	-2.4	311.5	<u> </u>	the	- the AA index of the ring opening followed							
7	<b>V</b>	171.1489	L[1Me]A	9 a2	1	-1.5	298.9		by	by a art b fragmant of the requilting linear							
7	-	545.3854	ALLVT[_BMe]AbuG[1Me]L[1Me]VL[1Me]	0 b10	2	-1.0	281.4		by a or b fragment of the resulting linear								
10	1	933.6766	L[1Me]L[1Me]V[1Me]T[_BMe]AbuG[1Me]L[1Me]V	1 b8	1	2.0	247.9		sec	quence.							
52	<b>V</b>	323.2327	LLV[1Me]	1 b3 (- NH3)	1	-0.7	230.7		000	1							
-																	





Figure 1: Workflow for the peptide MS/MS fragment annotation

Peptide: Cyclosporine A, A[\*1] L[1Me] L[1Me] V[1Me] T[\_BMe] Abu G[1Me] L[1Me] V L[1Me] A[\*1,(OH)-1]

> 2-Aminobutyric acid butenyl-methylmethyl

[\*x, (OH)-1]: attachment of index x, loss of OH group from C-terminal



Assigned: 140 of 495 peaks, MSMS Peak Area Assigned: 71.8%, Sequence Coverage: 11 of 11 consecutive amino acids

Requested types of theoretical fragments are generated from a given sequence.

ragmentation Settings	
ragment Types:	🔽 a 🔽 b 🔽 y 🔲 I 📄 Backbone Loss
laximum bonds to break:	2 🔽
reak linkages:	
lax Side Chain Losses:	3 🔻
ollapse Isotopic Signal:	

Experimentally observed neutral losses are validated against peptide structure and respective amino acids.

In case of isomeric fragments, the "winner" fragment requires the least bonds to break as well as evidence of a predecessor.



Comprehensive fragment annotation considers any combination of losses of peptide side chains. As such, the AL[1Me]L[1Me]V[1Me] fragment was found in 4 forms.



Peptide: Insulin, M=5729.6008 +TOF MS at 27.17 min 968.7818 969.1166 967.0 967.5 968.0 968.5 969.0 969.5 970.0 970.5 971.0

> Since series of isotope peaks are selected in Q1 and fragmented in parallel, peptide fragment isotopes in MS/MS spectra can be used in confident annotation.

As biotherapeutic peptides tend to generate the large fragments, the monoisotopic ions are usually weak or, sometimes, absent. We investigated a number of approaches to handle this challenge by using the information from fragment isotope peaks; by collapsing fragment isotopic signal to signal of fragment monoisotopic peak, we improve S/N of minor large fragments and enable their contribution to sequence confirmation. For example, in case of Insulin, the contribution of doubly charged fragments to the overall assignment rose from 3.6% to 6.9%.

Immonium ions (primary and secondary) are provided for the annotation from the amino acid dictionary. Cascading losses of water and ammonia are validated against the fragment type and





Figure 2: IDA (panel A) and SWATH acquisition (panel B) of cyclosporine MS/MS, with fragment annotation and distribution of m/z error, processed with a research MetabolitePilot<sup>™</sup> software .





Figure 3: MS/MS of peptide Insulin 6+ (MS/MS triggered on the first isotope).



Figure 4: MS/MS of Oxytocin triggered on doubly charged mono-isotopic peak. Predicted fragments include breaking of the disulphide bond and cascading losses of ammonia from both side chains and C-terminus.

243.1098 QN

226.0829 QN

v6 | b5

v6 | b5 (-NH3 \* 1

Peptide	Category	Sequence	Peaks (A)	Peaks (B)	Peaks (C)	TIC assigned (C)	Settings
Evonatido	Linoar	HGEGTFTSDL SKQMEEEAVR LFIEWLKNGG PSSGAPPPS-[Ami]	25	57	174	43.6%	S/N >10
Exenative	Linear						error< 5ppm
Overtagin	Lasso	C[*1]YIQNC[*1]PLG-[NH2]	52	81	81	61%	S/N >10
Oxytocin							error< 5ppm
Cyclo	Cualia		31	88	91	51.7%	S/N >10
(-GRGDSP)		G[*1]KGDSP[*1,(OH)-1]					error< 10ppm
Cyclosporin	Cyclic with A[*1]L[1Me]L[1Me]V[1Me]T[_BMe]AbuG[1Me]L[1Me]VL[1		100*	225	254	CE 70/	S/N >10
А	side chains	A[*1,(OH)-1]	160*	225	251	65.7%	Error< 5ppm
Inculin		GIVEQC[*1]C[*2]ASV C[*1]SLYQLENYC[*3] N/ FVNQHLC[*2]GSH		100	220		S/N >1000
insuin	iviuiti-chain	LVEALYLVC[*3]G ERGFFYTPKA	157	102	320	20.5%**	error< 5ppm

**Table 1.** Summary of peptide interpretation results with different in-silico annotation strategies Peaks (A) column indicates number of peptide MS/MS fragments having unique m/z, found considering a, b, y, yla and ylb fragments; *Peaks (B) column* indicates number of peptide MS/MS fragments considering Peaks A and up to 4 side-chain losses; *Peaks (C) column* indicates number of peptide MS/MS fragments considering Peaks B and collapsing fragment isotope peaks to one peak. \*Variable side-chain methyl losses were considered.

\*\*Fragments with up to 3 broken bonds were proposed in-silico.

Using the expanded algorithms, the number of annotated peptide fragments increased. The percentage of peak area assignment rose by a factor of 2.1 for peptide Oxytocin, and a factor of 2.0 for peptide Cyclo(-GRGDSP) compared to the previous annotation strategy.

# CONCLUSIONS

- sequences and their parts.
- fragment identification and sequence coverage.
- from larger multiply-charged fragments.
- localization and characterization of modifications.

# REFERENCES

Error (mDa) intensity

1 132.6

0.6 217.3

- 3. Wei-Ting Liu, W.-T. et al., Anal. Chem. 2009, 81, 4200–4209
- 4. Novak, J., Lemr, K., Schug, K.A., Havlicek, V., JASMS, 2015 (10), 1780-1786

# **TRADEMARKS/LICENSING**

AB Sciex is doing business as SCIEX. © 2018 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX<sup>™</sup> is being used under license.

Document number: RUO-MKT-10-7848-A



■ The workflow involving the SCIEX X500B QTOF System with a research version of MetabolitePilot<sup>™</sup> software effectively utilizes complex peptide fragments leading to the characterization of non-linear peptides. The enhanced non-linear peptide representation can effectively describe various types of non-linear peptide

With the implementation of non-natural amino acids and custom side-chain modifications into the amino-acid dictionary, the software exploits core of the cyclic structure as well as any side chain losses, thus improving

The enhanced algorithm recognizes fragment isotopic peaks and is well suited for leveraging information

The proposed workflow can be used not only for peptide sequence validation, but also for the detection,

Moore, I., Patel J., Rapid Peptide Catabolite ID using the SCIEX Routine Biotransform Solution, 2017, https://sciex.com/Documents/tech%20notes/Rapid-Peptide-Catabolite-ID-using-the-SCIEX-Routine.pdf 2. Patel J., Dindyal-Popescu A., Guo X., Moore I., Metabolite Profiling and ID of Cyclic Peptides Using Automated Software Processing and High Resolution Mass Spectrometry, ASMS 2018, ThP 143

5. Niedermeyer, T.H.J., Strohalm, M., PLOS One, September 2012, Volume 7, Issue 9, e44913