

Master Your Data to Get to the Answer, View Results with Ease and Confidence

Using TripleTOF® Technology and MasterView™ Software to Identify and Quantify Food Residues

André Schreiber¹, David Cox¹, Nick Zhu², and Cheng Yuan Cai²

¹AB SCIEX Concord, Ontario (Canada); ²AB SCIEX Shanghai, (China)

Overview


Here we present results of using a novel data processing approach of comparative multi-target screening to rapidly and accurately identify and quantify chemical residues in food. Fruit and vegetable samples were extracted using a QuEChERS procedure and analyzed with reversed phase LC. High resolution and accurate mass MS and MS/MS information was collected in a single injection on AB SCIEX TripleTOF® 4600 and 5600 systems. Data were processed using the new MasterView™ software.

Introduction

LC-MS/MS is a powerful analytical tool for the analysis of food residues and contaminants. In particular, triple quadrupole based mass analyzers are popular for targeted quantitation of hundreds of food contaminants in a single analysis because of their extra degree of selectivity and sensitivity when operated in Multiple Reaction Monitoring (MRM) mode.

Advancements in LC-MS/MS technology, including hybrid systems like triple quadrupole linear ion trap (QTRAP®) and quadrupole-quadrupole Time-of-Flight (TripleTOF®), now provide the ability to perform both targeted and non-targeted screening in food samples on a routine basis. However, full scan chromatograms are very rich in information and contain easily thousands of ions from both chemical compounds present in the sample as well as from the sample matrix itself. Thus, powerful software tools are needed to explore the high resolution and accurate mass data generated to get answers and results from these complex data.

TOF-MS and MS/MS information was used to identify and quantify targeted food residues. Quantitation was achieved by performing single concentration standard addition at the level of the Maximum Residue Limit (MRL). Identification was based on retention time, accurate mass of the quasi-molecular ion, isotopic pattern and MS/MS library searching. The molecular fingerprint saved into MS/MS spectra allow us to differentiate isomeric



compounds and greatly reduce the number of potential false positive results.


The new MasterView™ software allows quick processing and easy results review and reporting capabilities.

Experimental


- Fruit and vegetables samples from a local supermarket and samples from a EU proficiency tests for pesticides and fruits and vegetables
- Quantitation using the Sciex iD Quant™ standards kit for pesticide analysis
- QuEChERS extraction following guideline EN 15662/2007 using commercial QuEChERS kits
- 5-20x dilution of sample extracts to minimize possible matrix effects
- UHPLC using a Shimadzu UFLC_{XR} system with Restek Ultra Aqueous C18 (100 x 2.1 mm) 3 µm column
- Gradient of water and methanol with 10 mM ammonium formate at a flow rate of 0.5 mL/min and injection volume of 10 µL

- AB SCIEX TripleTOF[®] 4600 and 5600 system with DuoSpray[™] source operated in ESI mode
- Continuous recalibration between injections using the Calibrant Delivery System (CDS)
- Information Dependent Acquisition (IDA) using a TOF-MS survey scan 100-1000 Da (100 ms) and up to 10 or 20 dependent TOF-MS/MS scans 50-1000 Da (100 ms) using Collision Energy (CE) of 35 V with Collision Energy Spread (CES) of ± 15 V
- Dynamic background subtraction (DBS) was activated for best IDA coverage, no inclusion list was used to allow retrospective unknown identification without the need for a second injection to acquire MS/MS data

3. Review results using 'traffic lights' and display of MS and MS/MS spectra




4. Filter results based on identification criteria



5. Select results for:

- Reporting
- Export into a LIMS system
- Further processing in MultiQuant™ software
- Further processing in LibraryView™ software



Quantitation in MasterView™ Software

The quantitative comparison in MasterView™ software allows users to quickly identify compounds above a concentration of interest, such as the Maximum Residue Limit (MRL).

Load samples into MasterView™ software and specify the control sample. Load the target XIC table and click the 'Process' button to start data processing. Compounds identified with high confidence are indicated using green check symbols. Sort the table by the 'traffic light' column and/or filter by 'identification criteria' for easy result review and select compounds for report.

The batch of 20 samples shown in Figure 2 (containing one standard and 19 food extracts) was processed on a Windows 7 64-bit computer (8 core, 16 GB memory) in only 1:30 minutes.




Figure 2. Quantitative comparison in MasterView™ software: four pesticides were identified in the broccoli with high confidence and three of them were present at a concentration above 10 µg/kg in the sample (metalaxyl, imidacloprid, and spirotetramat)

Results of Samples of the EU Proficiency Test

Two samples from the EU proficiency test for pesticides and fruits and vegetables were extracted and analyzed for pesticides (Figures 3 and 4).

Retention times errors were between 0 and 2.4% and mass errors were between -1.7 and 2.1 ppm (positive polarity) and -0.7 and 2.4 ppm (negative polarity) and were well below the required 5 ppm (SANCO/12495/2011). The results summary is listed in Tables 1 and 2.

Table 1. Result summary of the leek sample of the EU proficiency test

In Sample	Identified	RT Error (%)	Mass Error (mDa)	Isotope Ration Difference (%)	MS/MS PUR (%)	Comment
<i>Carbofuran</i>	Carbofuran	0.1	0.2	1.7	99.2	
<i>Cyromazine</i>	Cyromazine	1.1	-0.1	0.3	99.4	
<i>DEET</i>	DEET	0.1	0.1	0.2	100	
<i>Diuron</i>	Diuron	0.1	-0.2	2.8	100	also identified in negative polarity
<i>Ethoxyquin</i>	Ethoxyquin	2.4	-0.1	16.2	45.7	low QuEChERS recovery – low MS/MS fit
<i>Fenpropidin</i>	Fenpropidin	0.1	0.4	0.4	99.1	
<i>Fenpyroximate</i>	Fenpyroximate	0.1	-0.5	8.4	96.8	
<i>Furathiocarb</i>	Furathiocarb	0.2	-0.1	2.3	98.4	
<i>Isofenphos-methyl</i>	Isofenphos-methyl	0.1	-0.3	1.9	100	
<i>Isoproc carb</i>	Isoproc carb	0.1	0.2	0.4	100	
<i>Mecarbam</i>	Mecarbam	0.2	-0.2	2.8	92.4	
<i>Metolachlor</i>	Metolachlor	0.3	0.1	1.1	97.8	
<i>Metribuzin</i>	Metribuzin	0.3	-0.2	5.2	97.7	
<i>Mevinphos</i>	Mevinphos	0.3	-0.1	0.5	100	
<i>Phorate</i>	Phorate	0.3	-0.5	14.7	77.4	also Phorate-sulfoxide identified
<i>Picolinafen</i>	Picolinafen	0.1	-0.3	0.3	99.5	also identified in negative polarity
<i>Prometryn</i>	Prometryn	0.3	0.0	0.8	99.7	isomer Terbutryn with different MS/MS
<i>Propazine</i>	Propazine	0.3	-0.2	2.4	100	isomers sec- and tert-Butylazine with different MS/MS
<i>Propoxur</i>	Propoxur	0.1	0.1	0.5	99.7	
<i>Pyraclostrobin</i>	Pyraclostrobin	0.3	0.2	2.4	98.4	
<i>Quinoclamine</i>	Quinoclamine	0.1	0.1	1.9	100	also identified in negative polarity
<i>Simazine</i>	Simazine	0.1	0.3	3.3	99.8	
<i>minimum</i>		0.1	-0.5	0.2	45.7	
<i>maximum</i>		2.4	0.4	16.2	100.0	
<i>average</i>		0.3	-0.1	3.2	95.5	

Table 2. Result summary of the pear sample of the EU proficiency test

In Sample	Identified	RT Error (%)	Mass Error (mDa)	Isotope Ration Difference (%)	MS/MS PUR (%)	Comment
<i>Bromacil</i>	Bromacil	0.0	0.1	10.3	100.0	also identified in negative polarity
<i>Bromoxinyl</i>	Bromoxinyl (negative)	0.1	0.3	13.7	99.0	
<i>Cadusafos</i>	Cadusafos	0.2	0.5	1.1	99.8	
<i>Diflubenzuron</i>	Diflubenzuron	0.3	0.4	2.7	99.9	also identified in negative polarity
<i>Diniconazole</i>	Diniconazole	0.3	0.7	6.5	100.0	
<i>Ethoxyquin</i>						low QuEChERS recovery – not detected
<i>Etrinfos</i>	Etrinfos	0.1	0.4	1.2	99.3	
<i>Fenpropidin</i>	Fenpropidin	0.1	0.2	0.1	99.8	
<i>Flufenacet</i>	Flufenacet	0.1	0.3	2.0	100.0	
<i>Flutolanil</i>	Flutolanil	0.3	0.3	0.5	98.0	also identified in negative polarity
<i>Imazapyr</i>	Imazapyr	1.3	0.3	0.5	99.5	also identified in negative polarity
<i>Ioxynil</i>	Ioxynil (negative)	0.2	0.2	0.5	99.7	
<i>Isoproturon</i>	Isoproturon	0.3	0.1	0.3	98.6	
<i>Metazachlor</i>	Metazachlor	0.2	0.2	16.9	100.0	interference on 13C isotope
<i>Napropamide</i>	Napropamide	0.1	0.2	1.5	100.0	
<i>Prometryn</i>	Prometryn	0.2	0.2	1.7	100.0	isomer Terbutryn with different MS/MS
<i>Propaquizafop</i>	Propaquizafop	0.3	-0.3	1.8	99.4	
<i>PyrifenoX</i>	PyrifenoX	0.1	0.1	7.3	100.0	
<i>Spinosad</i>	Spinosyn A + D	0.2	-0.9	1.8	100.0	
<i>Terbacil</i>	Terbacil (negative)	0.3	0.5	3.8	99.2	
<i>Terbumeton</i>	Terbumeton, Secbumeton	0.0	0.4	1.0	97.3, 99.6	no LC separation of Terbumeton and Secbumeton, isomer Prometon with different MS/MS
<i>Vamidothion</i>	Vamidothion	0.3	0.4	0.5	99.8	
<i>minimum</i>		0.0	-0.9	0.1	97.3	
<i>maximum</i>		1.3	0.7	16.9	100.0	
<i>average</i>		0.2	0.2	3.6	99.5	




Figure 3. 22 (out of 22) pesticides were identified in the leek sample



References


- ¹ EU Commission Decision 'concerning the performance of analytical methods and the interpretation of results' #2002/657/EC
- ² SANCO Document: 'Method Validation and Quality Control Procedures for Pesticide Residues Analysis in Food and Feed' #SANCO/12495/2011
- ³ A. Schreiber and C. Seto: 'Target and Non-Target Screening for Chemical Residues in Food Samples using the AB SCIEX TripleTOF® 4600 System and Intuitive Data Processing Tools' Application Note AB SCIEX (2012) #5680212-01
- ⁴ A. Schreiber and C. Borton: 'Target and Non-Target Screening for Pesticide Residues in Food Samples using the AB SCIEX TripleTOF® 5600 System' Application Note AB SCIEX (2010) #0460110-02
- ⁵ http://www.crl-pesticides.eu/docs/public/tmpl_article.asp?CntID=878&LabID=200&Lang=EN

The Complete Series of MasterView™ Software Workflow Demos:

Video 1:

Screen for **targeted** compounds in your unknown samples


Click www.absciex.com/MasterVideo1 to watch the video.



Video 2:

Easy set-up of data processing and data review parameters

Click www.absciex.com/MasterVideo2 to watch the video.



Video 3:

Screen for **non-targeted** compounds in your unknown samples

Click www.absciex.com/MasterVideo3 to watch the video.



For Research Use Only. Not for use in diagnostic procedures.

© 2013 AB SCIEX. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Publication number: 7830413-01