

# Confident Identification For Comprehensive Screening

## Using the Wiley Registry of Tandem Mass Spectral Data: MSforID

Pierre Negri<sup>1</sup>, Adrian M. Taylor<sup>2</sup>


<sup>1</sup>SCIEX, US; <sup>2</sup>SCIEX, CA

High resolution mass spectrometry (HRMS) using the SCIEX TripleTOF® or X500R QTOF systems provides forensic investigators a powerful tool for the detection and identification of toxic compounds or novel drugs of abuse. High MS and MS/MS acquisition rates ensure the acquisition of accurate mass spectral data on large numbers of analytes, even at low analyte concentrations. The comprehensive analyte-specific MS/MS fragment spectra in the compound identification workflow increases the selectivity, specificity and confidence in the results.

**Combining this data with the Wiley Registry of Tandem Mass Spectral Data: MSforID** enables accurate compound detection and identification through library spectral matching. Automated processing matches the MS data to the analyte mass and the MS/MS is matched to the analyte fragmentation pattern in the library, providing a “fit score” for assessment of confidence of compound identification. As shown in the example below, multiple points of high-quality data are used to detect and identify codeine and its main metabolite, norcodeine, from a postmortem case sample with excellent fit scores of 100% for both analytes.

## Key Advantages of the Wiley High Mass-Accuracy LC-MS/MS Spectral Library

- The Wiley Registry of Tandem Mass Spectral Data: MSforID spectral library includes a variety of compounds including pesticides, pharmaceutical compounds, illicit drugs and other small bioorganic molecules, useful for researchers in metabolomics, pharmaceutical research, forensics and toxicology
- Includes a collection of 10,000 positive and negative mode high resolution MS/MS spectra for over 1,200 compounds, all carefully measured in a series of controlled conditions to enable accurate, reliable, and reproducible search results in a variety of settings
- Library is compatible with SCIEX OS Software to process data from the TripleTOF®, X-Series QTOF and QTRAP® Systems with SCIEX OS Software
- Provides the ability to create customized libraries specific to workflows by simply selecting the analytes of interest in LibraryView™ Software



**Figure 1: Achieve Wide Analyte Coverage and Confident Compound Identification Using the Wiley Registry of Tandem Mass Spectral Data: MSforID.** XIC, TOF MS and MS/MS spectra obtained showing confident and accurate identification of (A) codeine and (B) its main metabolite, norcodeine, from a postmortem case sample. The sample was analyzed using IDA on the SCIEX X500R QTOF System and the acquired data was searched against the Wiley High Mass-Accuracy LC-MS/MS Spectral Library.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries.

© 2020 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-10085-A. AB SCIEX™ is being used under license.