



Enhanced protein therapeutic characterization utilizing SCIEX CE, icIEF-UV/MS and EAD MS/MS for superior PQA identification

Marcia Santos, Rashmi Madda, Haichuan Liu and Zoe Zhang
SCIEX, USA

This technical note demonstrates the utility of orthogonal technologies, such as capillary isoelectric focusing (cIEF), imaged cIEF (icIEF)-UV/MS, and LC-MS based peptide mapping and middle-down workflows in the realm of biopharmaceuticals, to enable an understanding of critical quality attributes (CQAs). These streamlined workflows can be leveraged separately or jointly for charge variant analysis at all stages of biotherapeutic discovery and development, empowering scientists to make informed decisions at critical stages to ensure the highest level of drug quality and safety.

Charge heterogeneity is an important CQA, which not only needs to be fully characterized, but also closely monitored to ensure the quality, safety and efficacy of a biotherapeutic product.^{1,2} The presence of various post-translation modifications (PTMs), such as deamidation, glycosylation, glycation and C-terminal clipping, leads to a complex charge profile of biotherapeutics.^{1,2} cIEF is the gold standard for charge variant analysis in biotherapeutics development; by accurately measuring migration time and determining the percent compositions of the charge variants, biopharma

scientists can quickly evaluate biotherapeutics' quality attributes, ensuring their safety and efficacy. However, there are technical and resource challenges that limit the ability to identify proteoforms in charge variant peaks in a timeframe that keeps pace with process development.

In this technical note, we address the challenge of proteoform identification and characterization in cIEF results, by harnessing the power of orthogonal workflows- cIEF, icIEF-UV/MS, peptide mapping, and middle-down assays to gain a comprehensive understanding of the attributes contributing to the complex charge profile of infliximab. (Figure 1).

Key features of full characterization toolset from SCIEX

- **High throughput:** Enable rapid charge heterogeneity assessment for high-throughput screening, monitoring and method optimization.
- **Confident identification:** icIEF-UV/MS and LC-MS workflows provide sensitive detection and confident identification of proteoforms in each charge variant.
- **Power of EAD:** EAD-based peptide mapping and middle-down workflows provide accurate localization of labile




Figure 1. Comprehensive characterization of infliximab and its charge variants using orthogonal analytical workflows. The multi-capillary cIEF technology offered by the BioPhase 8800 system provides reproducible charge variant separation in a high throughput manner. The icIEF-UV/MS capability of the Intabio ZT system coupled to ZenoTOF 7600 system offers icIEF separation, UV quantitation, and mass measurement of the charge variants in a single assay. LC-MS based peptide mapping and middle-down workflows using the ZenoTOF 7600 system enable confident sequence confirmation, PTM localization and isomer differentiation on the peptide or subunit level.


PTMs, rapid disulfide bond mapping, differentiation of amino acid isomers and high sequence coverage.

- **Platform assays:** The streamlined cIEF, icIEF-UV/MS and EAD-based peptide mapping workflows can be easily adopted for routine biotherapeutic characterization for the entire development pipeline.


Introduction

Antibody-based therapeutics are highly heterogeneous due to the presence of degradation products, aggregates and charge variants.^{1,2} The charge variants are predominantly formed from different modifications of the main species, such as C-terminal Lys pyroglutamate and amidation for the basic species and deamidation and glycation for the acidic species.² Many of these species are critical quality attributes whose levels must be closely monitored and controlled to ensure product quality.¹ While SCIEX cIEF approaches^{3,4} offer high-resolution and reproducible separation of the charge variant profile in a high throughput manner, it is challenging to identify the proteoforms in each charge variant with just this system. The Intabio ZT system enables an integrated icIEF-UV/MS workflow and effectively bridges the gap, enabling confident and precise characterization of various proteoforms within biotherapeutic charge variants.⁵⁻⁷ Further, to provide detailed information for biotherapeutics - such as sequence confirmation, PTM localization, isomer differentiation and disulfide bond mapping¹. LC-MS based methods, such as peptide mapping and middle-down, enable comprehensive characterization. The implementation of electron activated dissociation (EAD) for these LC-MS workflows significantly expands their capabilities for biotherapeutic characterization.⁸⁻¹⁰


In a previous cIEF study of intact infliximab, an intriguing charge profile was observed with abundant basic species¹¹. A different study using icIEF-UV/MS suggested that the C-terminal Lys residue primarily contributes to these species.⁵ Therefore, employing multiple orthogonal technologies together to elucidate the charge variant profile of biotherapeutics comprehensively is important. This work leverages cIEF, icIEF-UV/MS and LC-MS approaches to achieve a complete characterization of the charge profile of infliximab, as shown in Figure 2.




cIEF workflow observes abnormally high basic variant



icIEF/UV-MS workflow achieves the same separation as cIEF and identify the proteoforms contributing to the basic variants



Middle down confirm the existence of high level of lysine at Fc region



Peptide map confirms the existence of high level of C-terminal lysine

Figure 2. Flow chart to identify high basic variant peak observed in cIEF assays.

Methods

Detailed information about sample preparation, data acquisition and data analysis for the cIEF, icIEF-UV/MS and LC-MS workflows can be found in the previous technical notes.³⁻¹⁰ Below is a brief description of these steps for each workflow.

Samples preparation:

cIEF: All buffers and reagents were prepared following instructions in the Capillary Isoelectric Focusing (cIEF) Kit for the BioPhase 8800 system application guide.^{3,4} The infliximab sample was thoroughly mixed with the master solution comprised of 4M urea-cIEF gel, cathodic stabilizer, anodic stabilizer, pharmalyte 3-10 and pI markers 10.0 and 4.1. An aliquot of the mixture was transferred to a sample plate, which was spun at 20 *g* for 4 minutes.

icIEF-UV/MS: Infliximab was desalted and then mixed with the master mix solution containing arginine, pharmalyte 5 to 8 (Cytiva), pharmalyte 8 to 10.5 and 6.0 µg/mL peptide pI markers. The solution was vortexed and degassed by centrifugation.

LC-MS: Enzymatic digestion of infliximab was performed using the Trypsin/Lys-C (Promega) and FabBRICATOR (IdeS, Genovis) proteases for peptide mapping for middle-down analyses, respectively.

Data acquisition:

cIEF: The cIEF data was acquired using a BioPhase 8800 system (SCIEX) equipped with a UV detector. The cIEF separation method was executed at 20°C, with a focusing step of 25 kV for 15 min and the mobilization step at 30 kV for 30 min.

icIEF-UV/MS: During icIEF-UV/MS analysis, the infliximab sample was separated with an Intabio cartridge (SCIEX) installed on the Intabio ZT system (SCIEX). UV absorbance measurements were collected at 1 Hz during the focusing and mobilization steps. The samples were introduced into the ZenoTOF 7600 system by a metered 3 µL/min flow of chemical mobilizer, and the data was acquired

LC-MS: In the peptide mapping experiment, tryptic peptides were separated using an ACQUITY BEH C18 column (2.1 × 150 mm, 1.7 µm, 130 Å, Waters). A flow rate of 0.25 mL/min was used for the chromatographic separation. The column was kept at 60°C in the column oven of an ExionLC AD system (SCIEX). Mobile phase A was 0.1% formic acid (FA) in water and mobile phase B was 0.1% FA in acetonitrile. The peptide mapping data was acquired using a joint CID/EAD data-dependent acquisition (DDA) method described previously.⁸

The LC, Fc/2 and Fd subunits from IdeS treatment of infliximab were separated at a flow rate of 0.3 mL/min using an ACQUITY UPLC

Protein BEH C4 column (2.1 × 50 mm, 1.7 µm, 300 Å, Waters). The accurate mass of each subunit was measured using a TOF MS method. Middle-down analysis of the subunits, including LC, Fd and the Fc/2 subunit with or without the C-terminal Lys, was performed using an MRM^{HR} EAD method. Three charge states per subunit were selected for EAD fragmentation with an electron kinetic energy of 1 eV.

Data analysis:

cIEF: The cIEF data were analyzed using the BioPhase 8800 system software, version 1.0. The integrated peak areas of the main species and charge variants were used for percent composition calculations.

icIEF-UV/MS: UV traces and mass spectra from icIEF-UV/MS analysis were interpreted using Biologics Explorer software (SCIEX). Each peak in the charge profile was integrated using Intabio software to determine its peak area and percentage composition. Intact masses of the main peak and charge variants were determined from the deconvolution of the corresponding mass spectra.

LC-MS: Intact mass analysis of infliximab subunits was conducted using an intact workflow template in Biologics Explorer software. CID/EAD DDA and EAD MRM^{HR} data were interpreted using the peptide mapping and middle-down workflow templates, respectively.

Separation of charge variants using cIEF

cIEF is a powerful analytical separation technology for the characterization of protein therapeutics. cIEF is sensitive to small changes in the overall charge of the protein. For example, the retaining of the C-terminal Lys, pyroglutamate, amidation and succinimide formation yield basic variants, whereas deamidation, sialylation and glycation lead to acidic variants.² The percent compositions of the main, basic and acidic species are a key matrix for quality assessment of a biopharmaceutical because of the impact of charge variants on the safety and efficacy of the therapeutic product. The multi-capillary platform offered by the BioPhase 8800 system enables reproducible separation of the charge variants in a high throughput manner, facilitating rapid assessment of the product quality.

Figure 3 shows the cIEF charge profile of intact infliximab and 2 pI markers. A clear separation of the main, basic and acidic species was achieved. The percent compositions of the 3 species were calculated from their integrated peak areas using BioPhase 8800 system software. The result shows that the total composition of the basic variants (46.3%) is slightly higher than that of the main species (39.6%), while the acidic variants account for the remaining 14.1%. The high percentage of basic variants poses a question of what attributes contributed to the abundance of the basic peak. Further investigation is needed to determine the underlying modifications responsible for this observation, as C-terminal Lys, pyroglutamate, amidation and succinimide can all contribute to basic variants. Therefore, icIEF-UV/MS is employed to provide proteoform




Figure 3. cIEF separation of intact infliximab. cIEF provided excellent separation of the basic (grey shade) and acidic (green shade) variants from the main species (blue shade). The numbers in the parenthesis are the percent compositions (%) calculated from the integrated peak areas of 3 species.

identification in each charge variant. The relative peak abundance from cIEF separation agree well with those obtained from icIEF-UV/MS analysis in the previous and present studies, as will be described below.




Figure 4. icIEF-UV/MS analysis of intact infliximab. A) The icIEF-UV profile showcases the high-resolution separation of the main species from the basic and acidic charge variants using the Intabio ZT system. The numbers on top of each peak are the pI values. The integrated peak areas were used to obtain the percent compositions (%) shown in the parenthesis. The major basic and acidic peaks are B1, B2 and A1, respectively. B) The MS BPE acquired using the ZenoTOF 7600 system shows that peak separation was well maintained after the mobilization step. C) An overlay of deconvoluted mass spectra for the major species labelled, identified by Biologics Explorer in Figures A and B. Charge separation and accurate mass measurement led to a confident assignment of the C-terminal Lys species for the basic variants B1 and B2 and sialylated glycoforms for the acidic variant A1.

Separation and identification of charge variants by icIEF-UV/MS

Simultaneous separation, quantitation and identification of biotherapeutic charge variants in a single assay is challenging. The Intabio ZT system addresses this challenge by coupling charge separation and UV detection offered by icIEF-UV with accurate mass measurement using the high-resolution ZenoTOF 7600 system.⁵⁻⁷ The separation efficiency and resolution of icIEF are maintained after the mobilization step, enabling confident identification of different charge variants by MS. It was demonstrated previously that this powerful workflow provides reproducible separation and detection of intact mAbs and their charge variants.⁵⁻⁷

Figure 4 shows the icIEF-UV/MS results of intact infliximab obtained using the Intabio ZT system. The icIEF-UV profile reveals a high resolution separation of the main species and charge variants (Figure 4A). This result matches well with that obtained using cIEF, except that the basic and acidic species appear on the opposite sides of the main peak in 2 profiles, as the data was collected using icIEF separation (Figures 3 and 4A). The MS Base Peak Electropherogram (BPE) shows all the peaks observed on icIEF-UV profile were also detected in icIEF-MS profile without compromising separation resolution (Figure 4B), demonstrating that the separation power was maintained after chemical mobilization. From the icIEF-UV profile, the percent compositions of the main, acidic and basic species were calculated to be 50.6%, 12.6% and 36.8%, respectively.

These percentage numbers correlate with the values from the cIEF separation in this work (Figure 3) and the previous icIEF-UV/MS study⁵, highlighting the reproducibility and transferability of this technology.

The icIEF-UV/MS workflow enables confident assignment of each charge variant due to the ability to perform accurate mass measurement using the ZenoTOF 7600 system. Biologics Explorer software offers streamlined automatic identification of proteoforms from each charge variant peak, and scientists can select glycan libraries for different cell lines from the software. Figure 4C shows an overlay of the deconvoluted mass spectra for the major species labelled in Figures 4A and 4B, including the main peak, 2 basic species (B1 and B2), and 1 acidic species (A1). The main peak consists of 3 major glycoforms, including G0F/G0F, G0F/G1F and G1F/G1F, without the C-terminal Lys (green trace in Figure 4C). The 2 high abundant basic species (B1 and B2) are confirmed to be glycoforms with 1 or 2 C-terminal Lys residues (red and blue traces in Figure 4C). The acidic variant A1 is a mixture of the species carrying deamidation and sialylation (NeuAc) based on pI shift and the molecular mass detected in the mass spectra (purple trace in Figure 4C). Low-abundant potential glycosylated and sialylated species were also observed for the main peak. These results demonstrate the power of icIEF-UV/MS for elucidating the charge profile of biotherapeutics despite high sample heterogeneity. With this information we are able to determine that the high abundant basic variants observed during high throughput screening on the BioPhase 8800 system can be attributed to a potential C-terminal




Figure 5. Peptide mapping result confirmed the high abundant C-terminal Lys variant of infliximab. The joint CID/EAD DDA method combines the advantages of CID and EAD to achieve high sequence coverage, PTM localization, isomer differentiation and disulfide bond mapping.⁸ In this work, CID/EAD DDA provided a sequence coverage of >95% for light chain (LC) and heavy chain (HC) in a single injection (A). The peptides with or without the C-terminal Lys were confidently identified from MS/MS of the singly charged precursors by highly sensitive CID DDA (B and C).

lysine. This identification was accomplished with a significant reduction in time to result from weeks to approximately 1 hour for routine samples, compared to ion exchange chromatography (IEX) with fraction collection⁸.

Peptide mapping and middle-down analyses

To confirm the identification in icIEF-UV/MS workflow, a combination of techniques was employed, including peptide mapping and EAD based middle-down analysis. This comprehensive approach allows for a thorough characterization of infliximab, providing a deeper understanding of quality attributes. Previous studies have demonstrated the capability, versatility and flexibility of the ZenoTOF 7600 system equipped with the Zeno trap and EAD cell to achieve unparalleled peptide mapping and middle-down results for comprehensive biotherapeutic characterization in a single injection.⁹⁻¹¹ A 5- to 10-fold sensitivity increase in the detection of MS/MS fragments offered by the Zeno trap leads to high sequence coverage of biotherapeutics on the peptide or subunit levels. The ability of EAD to preserve labile PTMs and differentiate amino acid isomers greatly enhances the capability of peptide mapping and middle-down workflows for biotherapeutic characterization. In addition, the flexibility to create DDA or MRM^{HR} methods using CID, EAD or the combination of 2 fragmentation techniques in SCIEX OS software ensures that different needs for biotherapeutic development can be met.

Figure 5 shows selected peptide mapping results obtained using a joint CID/EAD DDA method. The combined CID and EAD result led to >95% sequence coverage of LC and HC in a single injection (Figure 4A). The HC peptides with or without the C-terminal Lys residues were confidently identified by the highly sensitive CID DDA (Figures 5B and 5C). In addition to achieving high sequence coverage, the joint CID/EAD DDA method provided the benefits of 2 fragmentation techniques. While CID provides effective fragmentation of the singly charged short peptides, as was the case for the C-terminal peptide without Lys (Figure 5B). EAD enables accurate localization of glycosylation and unambiguous differentiation of amino acid isomers. Figure 6 shows complementary CID and EAD MS/MS spectra of the G0F-containing glycopeptide EEQYNSTYR from the same data file. EAD produced a series of sequence fragments (e.g. *c*₅₋₇ and *z*₅₋₇) containing the G0F moiety (Figure 6B), enabling the accurate localization of the glycan.

Middle-down MS approach has proven to be a powerful tool for biotherapeutic characterization.⁹ It provides confident sequence confirmation, accurate PTM localization and rapid mapping of intra-chain disulfide bonds while benefiting from low modification artefacts due to simple sample preparation.⁹ Figure 7 shows the deconvoluted mass spectrum of the Fc/2 subunit from IdeS treatment of infliximab. The accurate mass measurement provided by the ZenoTOF 7600 system led to confident assignment of the




Figure 6. Complementary CID and EAD fragmentation of the G0F-containing glycopeptide EEQYNSTYR. While CID led to abundant oxonium ions at low *m/z* by cleaving the glycan (A), EAD preserved this labile PTM in the fragments (e.g. *c*₅₋₇ and *z*₅₋₇) for its accurate localization (B).




Figure 7. Deconvoluted mass spectrum of the Fc/2 subunit from infliximab. Accurate mass measurement of infliximab subunits enabled a confident assignment of the glycoforms detected for the Fc/2 subunit. The result confirmed the presence of abundant Lys variants in infliximab. The measured mass errors of the major glycoforms were shown in the parenthesis.

glycoforms detected for this subunit. The deconvoluted data revealed the presence of abundant Lys variants (e.g. G0F+K and G1F+K) in the sample, further supporting the data from cIEF and icIEF-UV/MS analyses (Figures 3 and 4). The Fc/2 G0F subunit with




Figure 8. Middle-down results confirmed the C-terminal sequences of infliximab and its Lys variant. EAD provided rich MS/MS spectra and high sequence coverage ($\geq 75\%$) of the Fc/2 G0F subunit (A) and its Lys variant (B) in a single injection. The C-terminal sequences of these 2 species can be confidently confirmed based on the detection of z- and y-series fragments. The insets show an expanded view of the fragments y_2 (A) for the main species and y_3 for the Lys variant.

or without the C-terminal Lys was selected for EAD-based middle-down analysis in an MRM^{HR} experiment. The middle-down data provided high sequence coverage of these 2 species ($\geq 75\%$) in a single injection (Figure 8). The detection of rich z- and y-series fragments enabled rapid confirmation of the C-terminal sequences. For example, the y_2 and y_3 fragments shown in the insets of Figure 8 correspond to the C-terminal sequences of PG and PGK for the main species and its Lys

variant, respectively. The detailed peptide mapping and middle-down results of infliximab will be described in a separate technical note.

In summary, the data presented in this technical note demonstrates the power of orthogonal cIEF, icIEF-UV/MS and LC-MS workflows for a complete characterization of biotherapeutics and their charge variants. SCIEX characterization portfolio offers high-throughput cIEF assays for monitoring throughout the development pipeline. icIEF provides proteoform identification for each charge variant, while middle-down and peptide mapping confidently confirm the

results. Enabling decisions to be made earlier in the development pipeline while utilizing less resource intensive assays.

The multi-capillary cIEF workflow provides reproducible separation of the charge variants in a high throughput manner. Integrated MS detection in the icIEF-UV/MS workflow leads to simultaneous charge variant separation, accurate quantitation and identification, providing an extra layer of depth in biotherapeutic characterization by identifying the cause of charge heterogeneity. LC-MS based peptide mapping and middle-down workflows offer an added confidence in sequence confirmation on the peptide or subunit level. All the workflows described in this technical note are streamlined and easy to implement separately or jointly for routine charge variant analysis of biotherapeutics. This can help biopharma scientists make informed decisions earlier in the development process.

Conclusions

- SCIEX provides a wide range of analytical tools for biotherapeutic characterization.
 - Orthogonal cIEF, icIEF-UV/MS and LC-MS workflows offer an in-depth characterization of biotherapeutics and their charge variants.
 - The multi-capillary platform offered by the cIEF workflow provides charge variant separation in a high throughput manner suitable for product monitoring.
 - The icIEF-UV/MS workflow provides a 30-min sample analysis which is significantly faster than conventional cIEF and IEX workflows requiring fractionation for the following identification.
 - LC-MS based peptide mapping and middle-down workflows offer the advantages of sequence confirmation and PTM localization on the peptide or subunit level.
 - All the workflows described in this technical note are streamlined and can be implemented separately or jointly for routine charge variant analysis of biotherapeutics in the realm of biopharmaceuticals.
8. Ostrowski M. Rapid multi-attribute characterization of intact bispecific antibodies by a microfluidic chip-based integrated icIEF-MS technology. [Electrophoresis. 2022 Oct;1-9](#)
 9. A single-injection workflow for enhanced peptide mapping using collision-induced dissociation (CID) and electron-activated dissociation (EAD). [SCIEX technical note, MKT-28039-A](#).
 10. An enhanced single-injection middle-down workflow to achieve high sequence coverage and disulfide mapping of antibody subunits. [SCIEX technical note, MKT-30575-A](#).
 11. Differentiation of leucine and isoleucine for enhanced sequence variant analysis using electron activated dissociation. [SCIEX technical note, MKT-30799-A](#).
 12. Native State Charge Variant Analysis of Commercialized Monoclonal Antibodies in Minutes. [SCIEX technical note, RUO-MKT-02-10446-A](#)

References

1. Anna Robotham and John Kelly. (2020) LC-MS characterization of antibody-based therapeutics: recent highlights and future prospects. [Approaches to the Purification, Analysis and Characterization of Antibody-Based Therapeutics. Chapter 1: 1-33](#).
2. Yi Du *et al.* (2012) Chromatographic analysis of the acidic and basic species of recombinant monoclonal antibodies. [MABS. 4\(5\):578-585](#).
3. Purity assessment and charge heterogeneity analysis of intact and subunit of etanercept on a single CE platform. [SCIEX technical note, MKT-29756-A](#).
4. N-linked glycan analysis of NISTmAb, etanercept and trastuzumab. [SCIEX technical note, MKT-27527-A](#).
5. Direct and rapid multi-attribute monitoring of multiple intact monoclonal antibodies with a wide isoelectric-point range of 7.3 to 9.1. [SCIEX technical note, MKT-26996-A](#).
6. A comprehensive workflow to characterize deglycosylated NISTmAb using imaged capillary isoelectric focusing (icIEF)-UV/MS. [SCIEX technical note, MKT-27945-A](#).
7. Charge variant analysis of antibody-drug conjugates using an icIEF-UV/MS workflow. [SCIEX technical note, MKT-29837-A](#).

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2023 DH Tech. Dev. Pte. Ltd. MKT-31116-A



Headquarters

500 Old Connecticut Path | Framingham, MA 01701
USA
Phone 508-383-7700
sciex.com

International Sales

For our office locations please call the division headquarters or refer to our website at sciex.com/offices