

Oligonucleotide sequencing using a streamlined ion-pair reverse-phase LC-MS/MS workflow with automated data analysis

Haichuan Liu, Jingwen Ding and Zoe Zhang
SCIEX, USA

This technical note highlights a streamlined oligonucleotide (oligo) sequencing workflow using ion-pair reverse-phase liquid chromatography-tandem mass spectrometry (IP-RP LC-MS/MS) and automatic data analysis. This workflow combines excellent chromatographic separation and high-quality MS/MS data of oligos with intuitive and automated software tools, leading to rapid oligo sequencing with improved user experience.

Oligo therapeutics are promising drugs targeting a variety of diseases through various mechanisms of action.¹ Synthetic oligos often carry various modifications and product-related impurities. IP-RP LC-MS and MS/MS are commonly used to characterize oligos and their impurities. A streamlined workflow—from chromatographic separation to data acquisition and analysis—is essential to enable rapid mass measurement, sequence confirmation, and impurity assessment, accelerating the development and decision-making process of therapeutic oligos.

Key features of the streamlined IP-RP LC-MS/MS workflow for oligo analysis

- **Streamlined:** Oligo analysis is streamlined from chromatographic separation to MS/MS data acquisition and interpretation
- **Flexibility:** SCIEX OS software allows the user to build data-dependent CID methods using a fixed collision energy (CE) with spread or custom dynamic CEs, providing flexibility to obtain optimal MS/MS fragmentation of oligos in different lengths or charge states
- **Simplicity:** Biologics Explorer software provides easy-to-use data analysis workflows for oligo sequence creation, intact deconvolution, mass mapping, and MS/MS sequencing
- **Consistent recovery:** Biozen Oligo LC column combines core-shell particle technology and bioinert hardware, significantly increasing the chromatographic efficiency




Figure 1. A streamlined IP-RP LC-MS/MS workflow for oligo sequence analysis. In this workflow, oligos are separated based on IP-RP chromatography using a Phenomenex Biozen Oligo LC column, followed by LC-MS and MS/MS data acquisition using a ZenoTOF 7600+ system or a ZenoTOF 8600 system. Oligo deconvolution and sequencing are performed using intuitive data analysis workflows within Biologics Explorer software, v8.0.3. This powerful software provides ease and flexibility in defining oligo sequences, creating custom-modifications, performing intact deconvolution based on monoisotopic or average mass, and conducting oligo sequencing using MS and MS/MS data.

Methods

Sample preparation: A synthetic oligo 18mer with full phosphorothioate (PS) linkages was purchased from Integrated DNA Technologies (IDT). An antisense oligo (ASO) 20mer with a maleimide C2 linker (ASO-MC) on the 3' end was ordered from VectorBuilder. A stock solution of 1 mg/mL was prepared for 2 oligos. This solution was diluted to 10 ng/ μ L, from which 5 μ L (50 ng) was injected for LC-MS analysis.

IP-RP LC separation: Oligos were separated based on IR-RP chromatography using a [Phenomenex Biozen Oligo LC column](#) [2.6 μ m, 50 x 2.1 mm] installed on a Nexera XS inert HPLC system [Shimadzu]. The separation was performed using the gradient shown in Table 1 with a column temperature of 60°C and a flow rate of 0.3 mL/min. Mobile phases A and B are water and methanol/water (1:1 v/v), respectively, with 50 mM 1,1,3,3,3-hexafluoroisopropanol [Sigma-Aldrich] and 15 mM diisopropylethylamine [Sigma-Aldrich].

Table 1. Gradient for IP RPLC separation.

Time [Min]	Mobile phase A [%]	Mobile Phase B [%]
Initial	80	20
1	80	20
10	30	70
11	5	95
13	5	95
13.1	80	20
15	80	20

Mass spectrometry: MS data were acquired using a data-dependent acquisition (DDA) method with CID in the negative mode on a [ZenoTOF 7600+ system](#) [SCIEX] or a [ZenoTOF 8600 system](#) [SCIEX]. The key instrument parameters are listed in Table 2.

Table 2. Source and CID DDA parameters.

Parameter	MS	MS/MS [CID]
Workflow	Intact proteins	
Curtain gas		35
CAD gas		9
Ion source gas 1		50 psi
Ion source gas 2		50 psi
Temperature		400°C
Spray voltage		-4,500 V or -3,000 V ⁱ
Time bins to sum		8
Decustering potential		-80 V or -40 V ⁱ
Start mass	500 Da	150 Da
Stop mass	3,000 Da	2,000 Da
Accumulation time	0.2 s	0.08 s
Collision energy	-10 V	-(30-50) V or custom ⁱⁱ
CE spread	-	10 or 0
Max. candidate ions	-	8
Charge states	-	2-12
Q1 resolution	-	Low

ⁱSettings specific to the ZenoTOF 8600 system

ⁱⁱA CE equation of $CE = 0.0273 \times [m/z] + 0.8961$ was used for all charge states with the dynamic CE option enabled.²

Data analysis: Oligo sequence creation and intact deconvolution of TOF MS data were described previously.³ CID DDA data were processed using an oligo sequencing workflow (Figure 2) within [Biologics Explorer software](#), Version 8.0.3 [SCIEX].




Figure 2. Oligo sequencing workflow from Biologics Explorer software. Biologics Explorer software offers a prebuilt, intuitive data analysis workflow for oligo sequencing using MS/MS data. This workflow consists of 3 main sections, including data selection and processing, sequence candidate generation, and mass and MS/MS mapping. The activities marked with a star contain MS/MS specific settings.

Oligo sequence creation and deconvolution

Built upon powerful algorithms from Genedata Expressionist®, Biologics Explorer software offers intuitive data analysis workflows for comprehensive characterization of protein- or oligo-based therapeutics. These streamlined workflows require minimal user's input and hence are quick to learn and easy to run. For advanced users, the workflows offer the flexibility to fine-tune various settings or algorithms.

Oligo sequences are created using 1-letter abbreviations for the building blocks, including sugars, linkers and bases.³ For example, an "s" and "o" represent PS and phosphate (PO) linkers, respectively. Custom building blocks can be added using a Nucleotide Editor within Biologics Explorer software. In this work, 5-methyl-cytosine and 2'-methoxyethylribose [2'-MOE] were defined as Q and e, respectively, within the Nucleotide Editor. With the building blocks defined, the sequences of the 18mer and ASO-MC analyzed in this study were defined as the following:

18mer

eT-s-eQ-s-eA-s-eQ-s-eT-s-eT-s-eT-s-eQ-s-eA-s-eT-s-eA-s-eA-s-eT-s-eG-s-eQ-s-eT-s-eG-s-eG

ASO-MC

eQ-s-eT-s-eA-s-eG-s-eA-s-dT-s-dG-s-dC-s-dA-s-dA-s-dT-s-dG-s-dT-s-dT-s-dG-s-eQ-s-eT-s-eQ-s-eA-s-eT

The PS to PO modification was defined for 2 oligos in the Unlocalized tab of the Modifications activity within the sequence candidate generation section of the data analysis workflow (Figure 2). The ASO-MC analyzed in this work was modified with an MC C2 linker on its 3' end. The gains of the

formula of the linker [C₁₄H₂₁N₂O₇P] and its hydrated form [C₁₄H₂₃N₂O₈P] were specified in the Localized tab of the Modification activity (Figure 3).

Modifications [1] - Settings

Unlocalized Localized Display

Limit to: 2 Modifications per Sequence

Modifications: Table

Name	Gain	Loss	On Characters	On Locations (Backbone)
- MC-C2	C ₁₄ H ₂₁ N ₂ O ₇ P			3'
- MC-C2-H2O	C ₁₄ H ₂₃ N ₂ O ₈ P			3'

Figure 3. Add oligo modifications. Biologics Explorer software allows ease and flexibility in generating candidate sequences for MS or MS/MS mapping. For ASO-MC, the MC linker and its hydrated form were specified for the 3' end in the Localized tab of the Modification activity.




Figure 4. Intact deconvolution of the 18mer. Intact deconvolution using the monoisotopic mass led to accurate mass measurement of the 18mer and its PO species. The values in parentheses are the measured mass errors in ppm.




Figure 5. Intact deconvolution of the ASO-MC. The intact form of the ASO-MC and its hydrated species (+H₂O) were chromatographically separated using the Phenomenex Biozen Oligo LC column (A). Intact deconvolution using the monoisotopic mass led to accurate mass determination of these 2 forms and their PO species, with the measured mass errors of <5 ppm (B and C). The values in parentheses are the measured mass errors in ppm.

Intact oligo deconvolution

Similar high-quality TOF MS and CID DDA data of oligos were obtained using the ZenoTOF 7600+ system and the ZenoTOF 8600 system. The data from the ZenoTOF 7600+ system were used here to demonstrate the streamlined intact deconvolution and oligo sequencing workflows of Biologics Explorer software.

The intact deconvolution workflow for oligos was described in a previous technical note.³ The deconvolution can be performed using the monoisotopic or average masses. In this work, the monoisotopic masses of oligo candidate sequences were used to map the isotopically resolved spectra of the 18mer and ASO-MC.

Figures 4 and 5 show the deconvolution spectra of the 18mer and ASO-MC, respectively. High mass accuracies (<5 ppm) were obtained for the intact forms of these 2 oligos and their PO species. The MC linker of the ASO-MC is prone to hydrolysis,

leading to the formation of a hydrated (+H₂O) species. The Phenomenex Biozen Oligo LC column provided good separation of the ASO-MC and its hydrated species (Figure 5A) for their accurate mass determinations (Figures 5B and 5C).

Oligo sequencing

The CID DDA data of the oligos were acquired using a fixed CE with spread or dynamic CEs with a custom equation of $CE = 0.0273 \times [m/z] + 0.8961$ described in a previous study². The oligo sequencing workflow within Biologics Explorer software allows a quick comparison of the data acquired using different settings, accelerating method optimization for oligo analysis.

Figures 6 shows the CID DDA results of the 18mer obtained using 3 different CE settings. The MS/MS spectra of the charge states -3 to -10 were combined and deisotoped prior to MS/MS




Figure 6. 18mer sequencing using CID DDA methods with different CE settings. Biologics Explorer software enabled a quick comparison and visualization of the CID DDA results of the 18mer using different CE settings, such as a fixed CE of 30V (A) and 40V (B) with a spread of 10V and dynamic CEs with a custom equation of $CE = 0.0273 \times [m/z] + 0.8961$ (C). While all CE settings provided a complete sequence coverage of the 18mer (A-C), the custom dynamic CE setting led to the most balanced fragmentation of the oligo (C). The sequence coverage map of the 18mer obtained using this optimized setting is shown in D. The MS/MS spectra shown in A-C were deisotoped from the combined MS/MS spectra of charge states -3 to -10.




Figure 7. MS/MS sequencing of the ASO-MC. CID DDA with the custom dynamic CE setting resulted in excellent fragmentation of the ASO-MC (A) and its hydrated form (B), leading to complete sequence coverage of these 2 oligo species (C and D). The peaks marked with asterisks correspond to the *wb* or *wb* - H₂O species.

mapping against the oligo candidate sequences (Figures 6A-6C). While all CE settings provided a complete coverage of the 18mer, the dynamic CEs with a custom CE equation resulted in the most balanced fragmentation of the oligo (Figure 6C). The complete sequence coverage map of the 18mer obtained using the custom CE setting is shown in Figure 6D. The CID DDA method with the custom dynamic CE setting also provided high-quality MS/MS data and complete sequence coverage of the ASO-MC and its hydrated species, confidently confirming the oligo sequence and the location of the MC modification on the 3' end (Figure 7).

In summary, this technical note highlights a streamlined IP-RP LC-MS/MS workflow for accurate mass measurement and sequence analysis of oligos. The oligo sequencing workflow offered by Biologics Explorer software allows rapid method optimization and sequence confirmation, accelerating the development of oligo therapeutics.

Conclusions

- A streamlined IP-RP LC-MS/MS workflow was developed for comprehensive oligo sequence analysis, leveraging the

optimized DDA method with custom dynamic CEs using the ZenoTOF 7600+ system or the ZenoTOF 8600 system and the intuitive oligo sequencing workflow from Biologics Explorer software.

- Rapid data processing using the oligo sequencing workflow allows quick comparison of DDA data acquired using different settings, accelerating method optimization.
- Compared to fixed CEs with spread, the custom dynamic CEs provided balanced fragmentation of oligos in different charge states, leading to high-quality MS/MS spectra and complete sequence coverage of the oligos.
- Biozen Oligo LC column provided good separation of the ASO-MC and its hydrated form for confident characterization of these 2 species.

References

1. Pourshahian S. [2021] Therapeutic oligonucleotides, impurities, degradants, and their characterization by mass spectrometry. [Mass Spectrom. Rev. 40:75-109.](#)
2. Karasawa K et al. Method development for metabolite and impurity profiling of oligonucleotide therapeutics. ASMS *Conference on Mass Spectrometry and Allied Topics*, June 2-6, 2019, Atlanta, Georgia.
3. A streamlined LC-MS workflow for intact oligonucleotide analysis with enhanced MS sensitivity. [SCIEX Technical Note, MKT-36686-A.](#)

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product[s] not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries [see www.sciex.com/trademarks].

© 2026 DH Tech. Dev. Pte. Ltd. MKT-37808-A



Headquarters
250 Forest Street | Marlborough,
MA 01752 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division
headquarters or refer to our website at
sciex.com/offices