

Simple and fast-tracking of DAR distribution using intact multiple attribute methodology (Intact MAM)

DAR determination for cystine-linked ADC (Vorsetuzumab) featuring SCIEX OS software

Kazuko Inoue¹, Ryo Yokoyama², Chie Inagaki², Greg Roman³, and Zoe Zhang³

¹Eisai Co.,Ltd., Japan ²SCIEX Japan, ³SCIEX USA

This technical note highlights the utility of SCIEX's flexible solution for intact MAM within SCIEX OS software. This example shows the drug antibody ratio (DAR) monitoring for antibody-drug conjugates (ADCs). The streamlined MAM data flow from data acquisition to data processing in compliance-ready SCIEX OS software will also be demonstrated.

ADCs are biotherapeutics composed of small molecule drugs conjugated to an antibody scaffold. An ADC molecule consists of three components: monoclonal antibody, linker, and cytotoxic drug. As a result, ADCs have a very complex molecular structure with macromolecules and small molecule characteristics. The information about the average DAR is essential for assessing ADC qualities. For instance, a low DAR could reduce the efficacy of a drug. Therefore, DAR values are one of the critical quality attributes (CQAs) for the therapeutic index of ADCs.

In this technical note, we demonstrate SCIEX's intact MAM workflow in DAR determination of a cystine-linked ADC from

monkey and mouse plasma. A time course study was performed to showcase the streamlined LC-MS based workflow, combining SCIEX high-resolution QTOF mass spectrometry with user-friendly data analysis software.

Key features of SCIEX's intact MAM solution for DAR monitoring

- **High throughput:** DAR monitoring for a large sample set
- **Intuitive software:** Powerful product attribute definition, tracking, and quantification with flexible custom calculations for CQAs based on specific user needs.
- **Streamlined solution:** Complete software solution for acquisition and quantitative tracking of DAR changes.
- **Compliance-ready:** DAR quantification and monitoring are performed in compliance-ready SCIEX OS software.



Streamlined Sciex flexible solution for MAM




Figure 1. A streamlined workflow for DAR monitoring using SCIEX intact MAM

Methods

Sample preparation:

Vorsetuzumab Mc-VC-PAB-MMAE ADC was purchased from Creative Biolab Inc., (Shirley, NY). Human, monkey, mouse plasma, and PBS samples were incubated with vorsetuzumab MMAE ADC (30 μ L, 0.1 mg/mL) for time intervals of 0, 24, 48, and 72 hrs. These samples were diluted with HBS-EP buffer (260 μ L), followed by the addition of biotinylated human CD27 ligand as antigen of vorsetuzumab (10 μ L of 0.2 mg/mL aqueous solution). The mixture was incubated at room temperature for 15 minutes to generate conjugates with biotinylated human CD27 and vorsetuzumab MMAE ADC. Subsequently, streptavidin immobilized magnetic beads (Dynabeads M-280 streptavidin) were suspended into the resulting mixture and incubated at room temperature for 15 minutes with continuous mixing. After incubation, the resulting magnetic beads collected from the mixture were washed with HBS-EP buffer twice.

Table 2. Chromatography for intact and subunit analysis.

Time [min]	Mobile Phase A [%]	Mobile Phase B [%]
<i>Initial</i>	80	20
8	50	50
9	10	90
12	10	90
12.1	80	20
17	80	20

The biotinylated vorsetuzumab MMAE ADC was released from the streptavidin magnetic beads using a pre-mixed PNGaseF reaction solution (1 μ L of rapid PNGaseF diluted with 37 $^{\circ}$ C prewarmed MilliQ water). Then it was incubated at 37 $^{\circ}$ C for 60 minutes to conduct on-bead deglycosylation. After collecting the magnetic beads from the rapid PNGaseF reaction mixture, the beads were washed with HBS-EP buffer. The vorsetuzumab MMAE ADC captured on the beads was eluted with 50 μ L of acetonitrile/MilliQ water/formic acid (10/90/0.01 V/V/V). The eluates were treated with 10 mM TCEP at a final concentration for 30 minutes at 37 $^{\circ}$ C to reduce the intra-chain disulfide bond of the ADC. The reduced ADC sample was injected into a TripleTOF 6600+ system with a UPLC system.

Chromatography:

Subunits were separated using an ACQUITY BEH protein C4 1.7 μ m, 2.1 x 50 mm analytical column (Waters), which was kept at 80 $^{\circ}$ C in the column oven of a Nexera 30A (Shimadzu) UHPLC system. Table 1 shows the LC gradient used for subunit separation at a flow rate of 0.3 mL/min with mobile phases A and

B consisting of 0.1% formic acid in water and 0.1% FA in acetonitrile, respectively.

Mass spectrometry:

LC-MS data were acquired using the TripleTOF 6600+ system. The key TOF parameters are listed in Table 2.

Table 2. MS parameters for subunit and intact mass analysis.

Parameter	Setting
<i>Scan mode</i>	TOF-MS
<i>Polarity</i>	Positive
<i>Intact protein mode</i>	ON
<i>Gas 1</i>	60 psi
<i>Gas 2</i>	60 psi
<i>Curtain gas</i>	30 psi
<i>Temperature</i>	450 $^{\circ}$ C
<i>Ion spray voltage</i>	5500 V
<i>CAD gas</i>	6
<i>Time bins to sum</i>	80
<i>Accumulation time</i>	0.5 s
<i>Start mass</i>	500 m/z
<i>Stop mass</i>	4,000 m/z
<i>Declustering potential</i>	120 V
<i>Collision energy</i>	10 V

Data processing:

Reconstruction and average DAR calculation of cystine-linked ADCs were performed in compliance-ready SCIEX OS software, version 2.2.




Figure 2. Reconstructed spectrum from Cys-linked ADC incubated in human plasma. Top: Heavy chain. Bottom: light chain. Peaks corresponding to the molecule carrying a higher number of drug payloads decreased throughout incubation time.

Data collection for DAR calculations

Compared to traditional antibodies, ADCs go through more potential biotransformations with a payload *in vivo* and *in vitro*, which can change DAR distribution significantly.¹ The change of *in vivo* DAR distribution of trastuzumab emtansine was recently

reported in literature.^{2,3} DAR is one of the CQAs that must be closely monitored during ADC manufacture and storage as it can affect drug efficacy and safety. In this work, a time course stability study was designed to investigate the changes of DAR distribution from incubation of vorsetuzumab in different plasma. Data were collected, evaluated, and reconstructed using SCIEX




Figure 3. Parameter definition for each subunit of vorsetuzumab MMAE ADC in SCIEX OS software. SCIEX OS software allows for convenient deconvolution of the respective subunits related to the ADC.

OS software over the time course study for vorsetuzumab MMAE ADC in human plasma, as shown in Figure 2. The ADC peaks with a different number of drug payloads were labeled in the reconstructed spectra. The dominant species observed for the heavy chain carried 0-3 drug molecules, while the light chain was mainly conjugated with 0 or 1 drug payload. The DAR distributions of both chains were shifted toward the lower payloads throughout incubation.

Method generation for DAR calculation

SCIEX OS software offers a simple and streamlined intact MAM solution⁴. Attributes can be easily defined, and integration parameters for individual components can be adjusted to achieve accurate quantitation in SCIEX OS software, as shown in the previous technical note.⁴ Figure 3 displays reconstruction and peak integration parameters for developing intact MAM. Using the advanced reconstruction parameters, the user can define the number of spectra to average, S/N threshold, resolution, and mass range for the assay. The attribute level can be calculated within each attribute group using flexible custom calculations. In this technical note, a method was built for

automated DAR calculation, as shown in Figure 4. The percentage of peak area for each payload from each subunit was calculated using the formula displayed in Figure 4A. As vorsetuzumab MMAE ADC has two heavy chains and two light chains, the total ratio of each payload needs to be multiplied by 2, as shown in Figure 4B. Finally, the average DAR was calculated by adding all the ratios from each payload (Figure 4C).

Intact MAM for automated DAR monitoring

Vorsetuzumab MMAE ADC samples incubated with plasma were purified through immunocapture. On-bead reduction and deglycosylation were performed to reduce sample complexity, and a TripleTOF 6600+ system was employed to measure the ADC's heavy chain and light chain. Molecular mass and quantitative information of each subunit were obtained using SCIEX's intact MAM solution. The average DAR of the entire ADC can be calculated by building a custom formula within SCIEX OS software, as illustrated in Figure 4. Detailed result tables can be obtained and customized for accelerating data review (Figure 5). The table includes all the quantification results from intact MAM workflow, such as the ADC's absolute and

Workflow (A) ← Accept changes and return to Calculated Columns × Discard

Use the calculator to create a new formula.

Formula name: % Area of HC Drug_1

COUNT	MAX	STDEV	Clear
SUM	MIN	MEDIAN	(
MEAN	ABS	IF)
GET	GETGROUP	GETSTAT	+
/	*	-	=

[HC Drug_1 Area] / [HC Total Area]

Note: The "Original text" option is recommended for formulas that contain functions, such as the IF function, that compare non-numeric values to numeric values.

Workflow (B) ← Accept changes and return to Calculated Columns × Discard

Use the calculator to create a new formula.

Formula name: DAR of HC Drug_1

COUNT	MAX	STDEV	Clear
SUM	MIN	MEDIAN	(
MEAN	ABS	IF)
GET	GETGROUP	GETSTAT	+
/	*	-	=

2 * [0 * [% Area of HC Drug_0]]

2 * [1 * [% Area of HC Drug_1]]

2 * [2 * [% Area of HC Drug_2]]

2 * [3 * [% Area of HC Drug_3]]

Note: The "Original text" option is recommended for formulas that contain functions, such as the IF function, that compare non-numeric values to numeric values.

Workflow (C) ← Accept changes and return to Calculated Columns × Discard

Use the calculator to create a new formula.

Formula name: Average DAR

COUNT	MAX	STDEV	Clear
SUM	MIN	MEDIAN	(
MEAN	ABS	IF)
GET	GETGROUP	GETSTAT	+
/	*	-	=

[DAR of HC Drug_0] + [DAR of HC Drug_1] + [DAR of HC Drug_2] + [DAR of HC Drug_3] + [DAR of LC Drug_0] + [DAR of LC Drug_1]

Note: The "Original text" option is recommended for formulas that contain functions, such as the IF function, that compare non-numeric values to numeric values.

Figure 4. Custom formula for DAR calculation. (A) The calculation for the % area of each payload. Drug_1 is shown as an example. (B) The calculation for the DAR of each payload. Drug_0, Drug_1, Drug_2, and Drug_3 on the heavy chain are shown as an example. (C) The calculation for the average DAR of the entire ADC.

Index	Sample Name	Average DAR	Component Name	Area	Expected MW	MW	Retent. Time	*HC Total...	*LC Total...	*% Area o...	*% Area o...	*% Area o...	*% Area o...
1	Human_0hr	4.005	HC Drug_0	217708.684	49031.46	49034.06	3.28	1100783...	7632423...	0.198	0.422	0.216	0.164
7	Human_24hr	3.320	HC Drug_0	310455.958	49031.46	49034.06	3.28	1100783...	7632423...	0.380	0.281	0.210	0.129
13	Human_48hr	3.093	HC Drug_0	675670.600	49031.46	49034.06	3.28	1100783...	7632423...	0.383	0.348	0.192	0.077
19	Human_72hr	2.732	HC Drug_0	720611.281	49031.46	49034.06	3.28	1100783...	7632423...	0.451	0.302	0.193	0.053
25	Monkey_0hr	3.963	HC Drug_0	224700.232	49031.46	49033.46	3.26	1101100...	7290826...	0.204	0.414	0.219	0.164
31	Monkey_24hr	3.203	HC Drug_0	520477.272	49031.46	49032.91	3.29	1388537...	1114825...	0.375	0.327	0.185	0.113
37	Monkey_48hr	2.662	HC Drug_0	778331.693	49031.46	49032.91	3.14	1825259...	1494614...	0.426	0.325	0.171	0.077
43	Monkey_72hr	2.586	HC Drug_0	942245.913	49031.46	49032.87	2.95	2263696...	1442856...	0.416	0.419	0.121	0.044
49	Mouse_0hr	4.157	HC Drug_0	116379.661	49031.46	49033.48	3.29	6988035...	3793062...	0.167	0.420	0.218	0.196
55	Mouse_24hr	3.670	HC Drug_0	298439.038	49031.46	49033.56	3.28	1015596...	8005060...	0.294	0.347	0.196	0.163
61	Mouse_48hr	3.384	HC Drug_0	476676.363	49031.46	49033.11	2.95	1515097...	1069604...	0.315	0.361	0.211	0.114
67	Mouse_72hr	3.224	HC Drug_0	427207.493	49031.46	49033.29	2.95	1230837...	8286971...	0.347	0.357	0.199	0.097
73	PBS_0hr	3.942	HC Drug_0	111708.068	49031.46	49034.89	3.28	5382101...	3731481...	0.208	0.448	0.192	0.153
79	PBS_24hr	3.886	HC Drug_0	92896.534	49031.46	49035.07	3.28	4222307...	2686647...	0.220	0.443	0.182	0.155
85	PBS_48hr	3.911	HC Drug_0	76542.650	49031.46	49035.34	3.29	3259591...	1981727...	0.235	0.409	0.190	0.166
91	PBS_72hr	3.723	HC Drug_0	119146.311	49031.46	49034.98	3.28	4762897...	2909359...	0.250	0.416	0.181	0.153

Manual Integration

Retention Time (RT)

Expected RT: 3.985 min

XIC Gaussian Smooth Width: 0.0 points

RT Half Window: 30.0 sec


Mass Peak Selection

Expected MW: 49031.46 Da

Mass Half Window: 5 Da

Figure 5. Result table and metric plot of DAR in SCIEX OS software. (A) Result table showing reconstructed peak area and DAR calculated customized formula. TIC, MS, and reconstructed spectrum were displayed on the same page. (B) Metric plots allow for easy visualization of DAR distribution changes.

percentage peak areas with a different number of drugs and



1; Human_0hr, 2; Human_24hr 3; Human_48hr 4; Human_72hr
5; Monkey_0hr, 6; Monkey_24hr, 7; Monkey_48hr, 8; Monkey_72hr
9; Mouse_0hr, 10; Mouse_24hr, 11; Mouse_48hr, 12; Mouse_72hr
13; PBS_0hr, 14; PBS_24hr, 15; PBS_48hr, 16; PBS_72hr

Figure 6. Average DAR of Intact ADC and DAR of each payload results. A; Average DAR monitoring from each sample. B; DAR monitoring of each payload result.

average DAR (Figure 5A). The results can be sorted by sample, targeted attribute, or modification event.

A metric plot can be created to visualize the change of each attribute (Figure 5B), offering a quick way of detecting changes in the molecule. The metric plot in Figure 5B clearly shows that the DAR distribution decreased throughout incubation of vorsetuzumab MMAE ADC in different plasma (human, monkey, and mouse). By comparison, the control samples incubated with PBS did not show this trend, as the DAR remained consistent over 3 days of incubation. The change in DAR distribution can also be visualized using a bar graph, as displayed in Figure 6. These results suggested that vorsetuzumab MMAE ADC underwent a potential biotransformation process of the payload in plasma, leading to a decrease in the average DAR value.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-15063

Conclusions

- SCIEX's flexible solution for intact MAM within SCIEX OS software offers a powerful workflow for ADC analysis by providing a streamlined and compliant software package, from data acquisition through data analysis
- The combination of SCIEX high resolution mass spectrometry and streamlined software presents a cutting-edge solution for attribute monitoring in process development, enabling faster decision making
- The TripleTOF 6600+ system and SCIEX OS software provide an excellent DAR monitoring tool.

References

1. Simone, chadt et al. (2019) Are Biotransformation Studies of therapeutic proteins needed? Scientific considerations and technical challenges. [Drug Metab Dispos. 2019 Dec;47\(12\):1443-1456.](#)
2. Keyang, Xu et al (2013) Characterization of the drug-to-antibody ratio distribution for antibody–drug conjugates in plasma/serum. *Bioanalysis*. 2013 May;5(9):1057-71.
3. Kaur, Surinder (2013) Bioanalytical assay strategies for the development of antibody–drug conjugate biotherapeutics [Bioanalysis. 2013 Jan;5\(2\):201-26](#)
4. Compliant Attribute Monitoring for Biopharmaceutical Product Quality Attributes Employing Intact Mass Analysis. SCIEX technical note, RUO-MKT-02-11314-A.



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division headquarters or refer to our website at sciex.com/offices