

Sensitive multiplexed quantification of protein biomarkers for early drug discovery and development

Increasing sensitivity for quantification of surfactant proteins in matrix using the SCIEX 7500 system

Bo An¹, Timothy Sikorski¹, John Mehl¹, Eshani Nandita², Remco van Soes² and Elliott Jones²

¹Glaxosmithkline, Upper Providence, Pennsylvania; ²SCIEX, Redwood City, California

This technical note demonstrates a highly sensitive quantification workflow for surfactant proteins at a signature peptide level using a high-end triple quadrupole mass spectrometer. Enhanced sensitivity for quantification of surfactant proteins was achieved with the improved front-end technology to generate, capture and transmit more ions on the SCIEX 7500 system. Low-ng/mL lower limits of quantification (LLOQs) were achieved for the 4 surfactant proteins analyzed in human plasma.

Despite the recent developments in LC-MS based protein biomarker quantification, sensitive multiplexed quantification of proteins remains a challenging step in the early stages of drug discovery and development. Biomarker targets such as surfactant proteins (SFTPs) are essential for facilitating and modulating the inflammatory response in the respiratory system. SFTP quantification in plasma is highly informative, as it enables the monitoring of lung health in an easily accessible matrix.

Quantification of protein biomarkers such as SFTPs is commonly performed using LC-MS techniques as it allows for the ability to multiplex. However, advanced sample pretreatment is required to extract protein biomarkers from matrix and sensitive analytical platforms are necessary to achieve quantification for low abundant targets. Target enrichment such as immunoprecipitation is often applied for quantification. However, it can be challenging and time consuming to develop an effective antibody particularly in the case of early drug discovery and development. Additionally, interferences such as anti-drug antibodies and low drug tolerance can further hinder immunocapture assay sensitivity.

Here, a novel antibody-free platform (ABFP) was applied for the enrichment of peptides in complex biological matrices.¹ ABFP employs a highly efficient chromatographic fractionation technique for signature peptide enrichment. Peptide analysis was performed using a microflow LC with trap-and-elute setup on the SCIEX 7500 system. Greater sampling efficiency of the microflow LC paired with the improved front-end technology on the SCIEX 7500 system substantially increased the overall sensitivity of the assay for SFTP biomarker quantification.

Key features of the protein biomarker quantification assay on the SCIEX 7500 system

- Accomplish sensitive multiplexed quantification of surfactant proteins in plasma at low-ng/mL level using the SCIEX 7500 system
- Achieve 2- to 5-fold improvement in LLOQ using the improved front-end technology on the SCIEX 7500 system compared with previous generation instrumentation
- Reach outstanding accuracy, precision and linearity for quantification of multiple protein targets in complex matrices




Figure 1. Schematic of the workflow for SFTP analysis. Samples were first fractionated and monitored online. Simultaneously, a set of fractions were collected and analyzed using the M5 microLC system in trap-and-elute mode coupled to the SCIEX 7500 system.

Methods

Sample preparation: Human SFTPs (hSFTPs) including hSFTP-A, hSFTP-B, hSFTP-C and hSFTP-D were enzymatically digested in human plasma. Details of the sample preparation protocol are provided in the manuscript.¹

Chromatography: Sample fractionation was on BioZen NX-C18 Column (150 x 2.1 mm), using high pH (pH=9.5) mobile phase, Separation was performed on a M5 MicroLC system operated in a trap-and-elute mode. A Phenomenex Luna Omega polar trap column (0.5 x 10 mm, 5 µm, 100 Å) was used for analyte trapping at a flow rate of 50 µL/min. The trap column was operated at room temperature. Sample separation was performed on a Phenomenex Luna C18 (2) column (0.3 x 100 mm, 3 µm, 100 Å) at a flow rate of 10 µL/min. The column temperature was set to 50°C. Mobile phase A was 0.1% formic acid in water and mobile phase B was 0.1% acetic acid in acetonitrile for both analyte trapping and separation.

Table 1. Gradient conditions for analyte trapping.

Time (min)	Mobile phase A (%)	Mobile phase B (%)
0.0	95	5
1.0	95	5

The gradient conditions are summarized in Tables 1 and 2 for analyte trapping and separation, respectively. A 50 µL aliquot of the sample was injected for LC-MS/MS analysis.

Table 2. Gradient conditions for analyte separation.

Time (min)	Mobile phase A (%)	Mobile phase B (%)
0.0	90	10
1.0	90	10
9.0	65	35
9.3	10	90
10.8	10	90
11.0	90	10
12.0	90	10

Mass spectrometry: A SCIEX 7500 system with an OptiFlow pro source with E Lens probe and micro electrode was used in positive MRM mode. The source conditions and MRM transitions used are summarized in Tables 3 and 4, respectively.

Table 3. Source conditions for the SCIEX 7500 system.

Parameter	Settings
Polarity	Positive
Curtain gas	35 psi
Ion source gas 1	20 psi
Ion source gas 2	65 psi
CAD gas	11
Ion spray voltage	3750 V
Source temperature	200°C

Data processing: The MRM data were processed using the Analytics module in the SCIEX OS software, version 2.0 using the MQ4 integration algorithm. A weighting of $1/x^2$ was used for quantification.

Table 4. MRM transitions for hSFTP-A, hSFTP-B, hSFTP-C and hSFTP-D.

Protein	Sequence	Precursor ion (m/z)	Product ion (m/z)	CE (V)	CXP (V)
hSFTP-A	NPEENEAIASFVK	724.357	1107.568	36	15
hSFTP-B	SPTGEWLPR	521.769	858.447	21	15
hSFTP-C	IAPESIPSLEALNR	755.417	663.357	27	15
hSFTP-D	SAAENAALQLVVAK	756.923	969.609	33	15

Overview of the workflow for surfactant protein quantification assay

A total of 4 surfactant proteins with different physiochemical properties were analyzed in plasma. This included hSFTP-A, hSFTP-B, hSFTP-C and hSFTP-B. The ABFP technique was applied to enrich the signature peptides. A concentration range between 10 ng/mL and 10,000 ng/mL was measured. A total of 2 replicates were assessed for each calibration point. Analysis was performed on the SCIEX Triple Quad 6500+ system and the SCIEX 7500 system to compare sensitivity for analysis of SFTPs.

Quality control (QC) samples at low, mid and high concentrations were also assessed. QC low, mid and high were prepared at 8000 ng/mL, 800 ng/mL and 100 ng/mL, respectively. A total of 4 replicates were analyzed for each QC level.

ABFP provided an effective enrichment method for a multitude of signature peptides in human plasma. Separation analysis using microflow LC with trap-and-elute setup enabled increased sampling efficiency enabling sensitive analysis of critical protein biomarkers. Furthermore, coupling microflow LC analysis to the SCIEX 7500 system with front-end improvements such as increase in ion generation, capture and transmission enabled quantification of low-level targets at ng/mL levels.

Quantification results

ABFP employs a solid-phase extraction (SPE) cleanup along with an orthogonal fractionation step to separate the peptide targets. ABFP was compared with direct analysis to assess the level of improvement in sensitivity for hSFTP-C (Figure 2). ABFP achieved a 5-fold S/N improvement at 50 ng/mL over traditional direct analysis. This demonstrates that the SPE fractionation




Figure 2. Extracted ion chromatograms (XICs) were used for the quantification of hSFTP-C. A S/N improvement of 5-fold was observed when comparing direct analysis and ABFP of hSFTP-C at 50 ng/mL using the SCIEX Triple Quad 6500+ system (A-B). With the application of the SCIEX 7500 system, a 9-fold S/N enhancement was achieved for the quantification of hSFTP-C at 10 ng/mL (C-D). A representative XIC of the IS intensity is shown on the bottom panel per analysis.

method substantially improves the sensitivity for quantification of signature peptides through targeted enrichment.

In addition, a comparative analysis was performed between the SCIEX 7500 system and the SCIEX Triple Quad 6500+ system using hSFTP-C at a concentration of 10 ng/mL. The application of the SCIEX 7500 system resulted in a S/N enhancement of up

to 9-fold compared with the previous generation instrument (Figure 2). As a result, the LLOQs were improved by a factor of 2x to 5x for the hSFTPs using the SCIEX 7500 system (Table 5).

Figure 3 shows the representative XICs at the LLOQ levels of hSFTP-A and hSFTP-D. An LLOQ of 5 ng/mL and 10 ng/mL were achieved for hSFTP-A and hSFTP-D, respectively.




Figure 3. XICs were used for the quantification of hSFTP-A and hSFTP-D. The LLOQs for hSFTP-A and hSFTP-D were 5 ng/mL and 10 ng/mL, respectively. Additionally, XICs showing endogenous levels of hSFTP-A and hSFTP-D in a healthy subject demonstrated that the assay can accurately quantify hSFTPs in disease and healthy states.

Moreover, the endogenous levels in healthy subjects were within the assay limits enabling the presented method to successfully quantify differences between healthy and disease states (Figure 3).

A highly sensitive, reproducible and efficient analysis of hSFTPs was demonstrated using ABFP enrichment coupled to a microflow LC and SCIEX 7500 system supporting quantification of low-level targets in early drug discovery and development.

Table 5. Summary of the LLOQs reached using direct analysis, ABFP and comparison of the SCIEX Triple Quad 6500+ system and SCIEX 7500 system.

Target protein	Direct analysis SCIEX Triple Quad 6500+ system	ABFP SCIEX Triple Quad 6500+ system	ABFP SCIEX 7500 system	LLOQ improvement between SCIEX Triple Quad 6500+ system and SCIEX 7500 system
<i>hSFTP-A</i>	1000 ng/mL	10 ng/mL	5 ng/mL	2x
<i>hSFTP-B</i>	200 ng/mL	25 ng/mL	10 ng/mL	2.5x
<i>hSFTP-C</i>	250 ng/mL	25 ng/mL	10 ng/mL	2.5x
<i>hSFTP-D</i>	200 ng/mL	50 ng/mL	10 ng/mL	5x

Table 6. Quantitative performance summarizing assay accuracy, precision and linearity at QC low, mid and high concentrations.

Target protein	Linear range (ng/mL)	Correlation coefficient (r^2)	Accuracy \pm CV (%)		
			QC Low	QC Mid	QC High
<i>hSFTP-A</i>	10-10,000	0.998	83.8 \pm 12.4	92.0 \pm 13.5	96.1 \pm 7.8
<i>hSFTP-B</i>	25-10,000	0.983	112 \pm 16.4	80.5 \pm 6.7	98.6 \pm 5.9
<i>hSFTP-C</i>	25-10,000	0.939	114 \pm 10.2	93.4 \pm 9.8	98.1 \pm 4.5
<i>hSFTP-D</i>	50-10,000	0.999	101 \pm 11.5	97.1 \pm 14.5	103 \pm 12.3

QC low, mid and high were evaluated for all surfactant proteins as shown in Table 6. For *hSFTP-A*, QC accuracies were within 16.2% of the nominal concentration with %CV less than 13.5%. QC accuracies for *hSFTP-B* were within 19.5% of the nominal concentration with %CV less than 16.4%. For *hSFTP-C*, QC accuracies were within 14% of the nominal concentration with a %CV less than 10.2%. Finally, for *hSFTP-D*, QC accuracies were within 3.0% of the nominal concentration with %CV less than 14.5%. Overall, excellent accuracy and %CV were reached for the quantification of hSFTPs. Additionally, strong linearity was achieved with correlation coefficients (r^2) greater than 0.983 (Table 6).

Conclusions

- Low-ng/mL LLOQs were achieved for quantification of hSFTPs in human plasma using the SCIEX 7500 system with improved front-end technology for ion generation, capture and transmission
- A 2x to 5x improvement in LLOQ were reached using the SCIEX 7500 system compared with the SCIEX Triple Quad 6500+ system
- Excellent linearity, accuracy and %CV were achieved for the quantification of hSFTPs using the SCIEX 7500 system

References

1. Bo An, Timothy W Sikorski, John F Kellie, Zhuo Chen, Nicole A Schneck, John Mehl, Huaping Tang, Jun Qu, Tujin Shi, Yuqian Gao, Jon M Jacobs, Eshani Nandita, Remco van Soest and Elliott Jones (2022). An antibody-free platform for multiplexed, sensitive quantification of protein biomarkers in complex biomatrices. [*Journal of Chromatography A*, 1676: 463261](#).

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-15128-A



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division headquarters or refer to our website at sciex.com/offices