



# Automated CE-based mRNA-LNP encapsulation efficiency and stability analyses with sample and reagent preparation on a Biomek i7 Hybrid Automated workstation

Jane Luo<sup>1</sup>, Sarah Gomes de Oliveira<sup>2</sup>, Tingting Li<sup>1</sup>, Mervin Gutierrez<sup>1</sup>, Zachary Smith<sup>2</sup>, Partha Banerjee<sup>2</sup>, and Sahana Mollah<sup>1</sup>

1. SCIEX, USA, 2. Beckman Coulter Life Sciences, USA.


This technical note describes an automated workflow for mRNA-LNP encapsulation efficiency [EE%] determination by capillary electrophoresis [CE]. The workflow integrates automated sample and reagent preparation on the Biomek i7 Hybrid Automated workstation with separation and analysis on the BioPhase 8800 system using the RNA 9000 Purity & Integrity kit.

The rapid expansion of mRNA-based vaccines and therapeutics has increased demand for efficient, reproducible workflows for EE% determination and mRNA purity and integrity assessment. The automated workflow described here enables streamlined, reproducible analysis while minimizing manual intervention.


## Key features

- **Automated EE% analysis with integrated RNA purity assessment:** Enables analysis of EE% and RNA purity and integrity with size estimation in a single CE workflow
- **Minimized manual intervention:** Significantly reduces hands-on time and mitigates the risk of RNase contamination
- **Reproducible, automated, RNA purity, integrity, and EE% analysis:** %RSD values of  $\leq 2.5\%$  for purity and  $< 0.5\%$  for EE%
- **Stability-indicating capability:** Detects temperature- and time-dependent mRNA degradation in stressed mRNA-LNPs


### Calibration curve generation



### mRNA quantity determination



### EE% calculation



**Figure 1: Automated workflow to determine the EE% of mRNA-LNPs.** The sample and reagent preparation were performed on the Biomek i7 Hybrid Automated workstation [A], CE analysis was carried out with the RNA 9000 Purity & Integrity kit [B] on the BioPhase 8800 system [C]. CPA: corrected peak area.

## Introduction

The success of COVID-19 vaccines has demonstrated the effectiveness of LNP-encapsulated mRNA and highlighted its broader potential for both vaccines and therapeutics. Beyond the physicochemical properties of LNPs, mRNA EE% is a critical quality attribute (CQA) for evaluating and screening vaccine formulation and assembly processes of mRNA-LNP products.<sup>1-2</sup> Rapid expansion of mRNA-LNP-based therapeutics has increased the demand for highly reproducible analytical workflows for EE% determination.<sup>3</sup> In addition, mRNA integrity analysis is essential for ensuring that mRNA drug substances remain intact and functional throughout development and manufacturing. However, manual sample and reagent preparation can limit efficiency and introduce variability. Furthermore, manual procedures inherently pose a risk of RNase contamination, which can lead to RNA degradation.

To address these challenges, a fully automated workflow integrating liquid handling automation with CE-based analysis has been developed. In this study, the workflow is implemented using the Biomek i7 Hybrid Automated workstation in combination with the BioPhase 8800 system and the RNA 9000 Purity & Integrity kit. This approach enables streamlined, reproducible EE% determination and assessment of mRNA purity and integrity.

While demonstrated using the Biomek i7 Hybrid Automated workstation, the workflow can be adapted to other automated liquid handling systems.

## Methods

**Materials:** [The RNA 9000 Purity & Integrity kit](#) [P/N: C48231], sample loading solution [SLS, P/N: 608082], and the [BioPhase BFS capillary cartridge - 8 x 30 cm](#) [P/N: 5080121] were from SCIEX (Marlborough, MA). Rainin LTS filter tips were from Mettler Toledo (Oakland, CA). Nuclease-free water [NFW] [P/N: AM9932] Surfact-Amps X-100 [10% (v/v) [Triton X-100, P/N: 28314], RNA storage solution [P/N: AM7000, 1mM sodium citrate], 10X Tris EDTA [TE] buffer [P/N: 75834] and 10X phosphate buffered saline [PBS, P/N: AM9624] were obtained from Thermo Fisher Scientific (Waltham, MA). Sucrose [P/N: S9378] was from Sigma-Aldrich (St. Louis, MO). The firefly luciferase [FLuc] mRNA [P/N: L-7602] was manufactured by TriLink BioTechnologies (San Diego, CA). The FLuc mRNA contains a 5'-cap and a 3'-poly(A) tail, mimicking an mRNA drug substance. The FLuc mRNA and heat-stressed

FLuc mRNA-loaded mRNA-LNP samples were provided by Acuitas Therapeutics Inc. (Vancouver, Canada). The 0.45  $\mu$ m syringe filter [P/N: 4654] was from PALL (Port Washington, NY). Part numbers for pipette tips, tubes, reservoirs, and plates used in the sample and reagent preparation workflow are listed in Table 1. A pair of custom-made outlet plate adapters was provided by Beckman Coulter Life Sciences (Indianapolis, IN).

**Table 1. List of labware and consumables required in this sample and reagent preparation workflow for one 96-well plate of sample and standards.**

Description	Vendor	Part number	Quantity
190 $\mu$ L pipette tips, sterile	Beckman Coulter Life Sciences	C41863	3
1025 $\mu$ L pipette tips, sterile	Beckman Coulter Life Sciences	B85955	3
Tube rack holder with 11 mm inserts	Beckman Coulter Life Sciences	373661 373696	1 10
Nunc microplate lids	Thermo Fisher Scientific	250003	2
Reservoir	Agilent	201244-100	2
Modular reservoir quarter module	Beckman Coulter Life Sciences	372790	1
96 PCR well plate	Axygen	P-96-450V-C	2
RNase-free microfuge Tubes, 2 mL	Thermo Fisher Scientific	AM12480	1
BioPhase sample and reagent plates [4,4,8]	SCIEX	5080311	1

**Liquid handler for sample and reagent preparation:** The Biomek i7 Hybrid Automated workstation [P/N C02613, Figure 1, panel A], equipped with a multichannel head, a Span-8 head, 2 ColdPlates, an orbital shaker, and a tip washing station, were from Beckman Coulter Life Sciences. The Biomek software [P/N: B87585] version was 5.1. Detailed lists of the required part numbers for automation hardware and software are provided in the Instructions for Use [IFU].<sup>4</sup>


**Storage of the RNA 9000 Purity & Integrity kit:** The Acid Wash/Regenerating Solution and CE Grade water were stored at room temperature upon receipt. The Nucleic Acid Extended Range Gel and the LIF Performance Test Mix were refrigerated at 2-8°C. The ssRNA Ladder and the SYBR™ Green II RNA Gel Stain\* were kept at -35°C to -15°C.<sup>5</sup>

**Preparation of the gel buffer:** The Nucleic Acid Extended Range Gel was pre-warmed to room temperature and filtered through a 0.45 µm syringe filter. The SYBR™ Green II RNA Gel Stain was thawed at room temperature in the dark and then added to the filtered gel at a 500-fold dilution.

**Preparation of 0.4% Triton X-100 solution:** The 0.4% Triton X-100 solution [Surfact-Amps X-100] was prepared by adding 1.6 mL of 10% Triton X-100 to 38.4 mL of NFW in a 50-mL conical centrifuge tube. The solution was mixed well by carefully inverting the tube at least 20 times and avoiding bubbles.

**Automated preparation of the diluent buffers for making the 8-point calibration curves:** The TE diluent buffer #1 [1 mL]

containing 1X TE, 1/5X 300mM sucrose in PBS, and 1/5X 1mM sodium citrate [SC] was prepared on the Biomek i7 Hybrid Automated workstation by mixing 100 µL of 10X TE buffer, 200 µL of 300mM sucrose in PBS, 200 µL of 1mM SC, and 500 µL of NFW in a 2 mL nuclease-free tube. The Triton diluent buffer #1 [1 mL] containing 1/25X 300mM sucrose in PBS, 1/25X 1mM SC, and 0.2% Triton X-100 was made by mixing 40 µL of 300mM sucrose in PBS, 40 µL of 1mM SC, 500 µL of 0.4% Triton solution, and 420 µL of NFW in a 2 mL nuclease-free tube. The components in the diluent buffers were mixed well by pipetting up and down 10 times.



**Figure 2. Automated workflows for free and total mRNA quantitation to determine mRNA-LNP EE%.** Panel A: Overall workflow for the EE% analysis with the SCIEX RNA 9000 Purity & Integrity kit and the BioPhase 8800 system, with sample and reagent preparation automated on the Biomek i7 Hybrid Automated workstation. Panel B: Automated parallel workflows for preparation of mRNA standards and mRNA-LNP samples enable determination of total and free mRNA using dual calibration curves.

### Automated preparation of the diluent buffers for mRNA-LNP

**sample treatments:** The TE diluent buffer #2 [7 mL] was prepared on the Biomek i7 Hybrid Automated workstation by mixing 0.875 mL of 10X TE buffer, 1.75 mL of 1mM SC, and 4.375 mL of NFW in a 4-divider reservoir. The Triton diluent buffer #2 [8 mL] was prepared by mixing 0.336 mL of 1mM SC, 4.168 mL of 0.4% Triton solution, and 3.496 mL of NFW in a 4-divider reservoir. The components in the diluent buffers were mixed well by pipetting up and down 10 times.


### Preparation of the sample plates for processing on the Biomek i7 Hybrid Automated workstation:

Two Axygen plates were used. For the analysis of total mRNA, 2 µL of the mRNA-LNP samples were transferred from the stock vial into each well of columns 2-6 on an Axygen plate. For the analysis of free mRNA, 16 µL of the mRNA-LNP samples were transferred from the stock vial into each well of columns 8-12 on another Axygen plate.

### Automated sample and reagent preparation for the determination of encapsulation efficiency of the mRNA-LNP:

To demonstrate the Biomek i7 Hybrid Automated workstation's capability for supporting mRNA-LNP encapsulation efficiency determination, a full plate was processed [40 samples in addition to an 8-point calibration curve, each in TE buffer or in Triton buffer]. The entire automated workflow is shown in Figure 2. After launching the method, selections were made for 40 samples at a final analysis volume of 55 µL, an 8-point free mRNA calibration curve in TE buffer, an 8-point total mRNA calibration curve in Triton buffer, sample outlet plates, reagent inlet plates, and reagent outlet plates on the MOS, as shown in Figure 3. The MOS provides a structured, step-by-step approach to experiment setup, complemented by the Guided Labware Setup [Figure 4] for accurate deck configuration and reagent volume calculations. Users can choose from 1 to 40 samples, at a final volume of 50 to 100 µL. For calibration curves, users can choose between 3 and 8 dilutions. This automated workflow is flexible and supports customizable sample processing and throughput. Once the deck setup was completed, diluent buffers for calibration curves or sample dilution were prepared in TE or Triton X-100 using the Biomek i7 Hybrid Automated workstation. Then, the Triton diluent buffer #2 was added to the samples at 48 µL per well to release mRNA from the LNPs. The calibration curve in the Triton diluent buffer #1 was generated by serial dilution of the FLuc mRNA standard from 1 mg/mL to 7.8 µg/mL. The diluted standards [column 1] and samples [columns 2-6] were incubated at room temperature for 20 minutes on a shaker


at 200 rpm. After that, equal volumes of formamide were added to each well, and the samples were incubated at 70°C for 5 minutes. Then, the plate was transferred to a cold plate for snap-chill. Next, FLuc mRNA standards were serially diluted in TE diluent buffer #1 to generate the calibration curve in column 7, followed by the addition of TE diluent buffer #2 [64 µL/well] to mRNA-LNP samples [columns 8-12] for free mRNA analysis. Subsequently, all preparation for the sample inlet plate, reagent inlet plate, sample outlet plate, and reagent outlet plate was also performed on the Biomek i7 Hybrid Automated workstation. Finally, sample separation and data processing were performed on the BioPhase 8800 system.



**Figure 3: Method Options Selector (MOS).** The automated MOS enables users to select the number and volume of samples and choose to run each section individually or from start to finish through a single user interface.

### Automated release of mRNA payload from mRNA-LNP samples in a stability study:

A master mix of a deformation solution containing Triton X-100 was prepared on the Biomek i7 Hybrid Automated workstation, dispensed into the mRNA-LNP samples, and thoroughly mixed by pipetting up and down. The released mRNA samples were diluted with formamide, heated at 70°C for 5 minutes, and snap-cooled at 4°C on the deck. Sample inlet, reagent inlet, sample outlet, and reagent outlet plates were prepared on the Biomek i7 Hybrid Automated workstation before the mRNA sample separation by CE on the BioPhase 8800 system.




**Figure 4: Deck set up for automated sample and reagent plates preparation for EE% analysis.** #1 & #9: Nunc lid; #2: Full reservoir with CE Grade water; #3: Sample outlet plate on an adapter; #4: Reagent outlet plate on an adapter; #5: Quarter modular reservoir for sample diluent preparation; #6, #12, and #13: 1025F tips; #7: Tube rack with 2 mL tubes for holding the acid conditioning solution, reagents for making diluents for calibration curves; #8, #14, and #15: 190F tips; #10: Full reservoir with the separation gel and dye; #11: Reagent inlet plate. BS1: Orbital shaker for deformulation incubation; CP1 and CP2: cold plates for heating or cooling the samples.

**The CE Instrument and software:** The BioPhase 8800 system [Figure 1, panel C] with UV/LIF [P/N: 5089278]—equipped with a laser-induced fluorescence [LIF] detector utilizing an excitation wavelength of 488 nm and an emission wavelength of 520 nm—was from SCIEX [Marlborough, MA]. Automated data acquisition and analysis were performed using BioPhase software version 1.4 from SCIEX.

## Results and discussions

**Fully automated mRNA-LNP EE% analysis:** The RNA 9000 Purity & Integrity kit is a well-established analytical tool for evaluating RNA integrity and purity.<sup>6-7</sup> It was also used for manual EE% analysis previously.<sup>8</sup> In this technical note, automated EE% analysis is performed with FLuc mRNA and FLuc mRNA-loaded mRNA-LNP samples as examples of the mRNA drug substance and the mRNA-LNP drug product. As shown in Figure 1 and described in the Methods section, the automation workflow for EE% analysis with the RNA 9000 Purity & Integrity kit on the Biomek i7 Hybrid Automated workstation is flexible, supporting the processing of 1–40 samples, accommodating calibration curves with 8 or fewer data points, and allowing certain steps to be performed off deck to save reagents. As a practical demonstration of its robustness, EE% was performed on 40 replicates of the same mRNA-LNP sample. Serial dilutions of the FLuc mRNA standard in either TE or Triton buffer were performed by the Biomek i7 hybrid Automated workstation. The standards

were then analyzed on the BioPhase 8800 system. Figures 5C and 5D show overlays of electropherograms obtained with serially diluted FLuc mRNA standards in TE or Triton buffer. CPA values of the FLuc mRNA peaks were plotted against their corresponding concentrations to generate the calibration curves in TE buffer [Figure 5A] and Triton buffer [Figure 5B].  $R^2$  values of 0.9958 and 0.9970 demonstrated excellent linearity. Stacked electropherograms of the mRNA-LNP samples analyzed for the free mRNA in TE buffer or the total mRNA in the Triton buffer are shown in Figures 6A and 6B. Consistent peak profiles were observed across replicates. No results were excluded, confirming the robustness of sample and reagent preparation by the Biomek i7 Hybrid Automated workstation. The average CPA value for free mRNA was 2.60, corresponding to a concentration of 20  $\mu\text{g/mL}$  based on the calibration curve generated in TE buffer. The average CPA value for total mRNA was 2.63, corresponding to a concentration of 461  $\mu\text{g/mL}$  based on the calibration curve generated in Triton buffer. The EE% for each mRNA-LNP replicate was calculated and ranged from 94.6% to 96.4%, with an average EE% of 95.6% and a %RSD of <0.5%. These EE% results were comparable to the 95% EE% determined manually for the same mRNA-LNP sample in a previous study.<sup>8</sup> Together, these results demonstrate that the fully automated EE% analysis workflow using the Biomek i7 Hybrid Automated workstation, the RNA 9000 Purity & Integrity kit, and the BioPhase 8800 system supports automated EE% determination





**Figure 5. Calibration curves in TE buffer and Triton buffer with excellent  $R^2$  values.** Serial dilutions of the FLuc mRNA standard in either TE or Triton buffer were performed by the Biomek i7 Hybrid Automated workstation. The standards were then analyzed on the BioPhase 8800 system. Panels C and D show overlays of electropherograms obtained with serially diluted FLuc mRNA standards in TE or Triton buffer. CPA values of the FLuc mRNA peaks were plotted against their concentrations to generate the calibration curves in TE buffer [A] and Triton buffer [B]. HMW: high molecular weight.

of mRNA-LNP samples. In addition, the purity of the total mRNA was also assessed in these 40 mRNA-LNP samples. The average purity was 97% with a %RSD of <0.3%. These results demonstrate that the same automated workflow enables both EE% determination and assessments of mRNA purity and integrity.

**Automated purity and size analysis of stressed mRNA and mRNA-LNP samples:** The FLuc mRNA-loaded mRNA-LNP


samples were subjected to thermal stress at 37°C for 2 and 5 days, or at 80°C for 20 minutes off deck. Further sample and reagent preparation were performed using the Biomek i7 Hybrid Automated workstation as described in the Methods section. Results in Figure 7 and Table 2 showed reduced full-length FLuc mRNA CPA% values, from 91.5% in controls to 80.4% and 68.2% in 2- and 5-day stressed samples, indicating temperature- and time-dependent RNA degradation. Two distinct degradant peaks



**Figure 6. Stacked electropherograms of the mRNA-LNP samples analyzed for the free mRNA [A] in TE buffer or the total mRNA [B] in the Triton buffer.** Panel A shows results obtained under non-denaturing conditions, which account for the additional peaks corresponding to RNA species with higher-order structure or HMW.

were observed in samples incubated at 37°C, with their CPA% increasing from day 2 to day 5, similar to previous results obtained manually with the same samples.<sup>8</sup> Their sizes were determined as 0.9 and 1.1 kb using a calibration curve generated with the ssRNA ladder from the RNA 9000 Purity & Integrity kit.

Additionally, LNP-encapsulated FLuc mRNA exhibited a smaller loss of the intact FLuc mRNA CPA% [12%] after stress at 80°C for 20 minutes compared with unencapsulated FLuc mRNA (a 28% loss), highlighting the protective effect of LNP formulation.



**Figure 7. Purity and size analysis of mRNA-LNP samples in an accelerated stability study.** The FLuc mRNA-loaded mRNA-LNPs were stressed at 37°C for 0, 2, and 5 days, or at 80°C for 20 minutes. The control sample showed a main FLuc mRNA peak at 1.9 kb, a FLuc mRNA peak with higher-order structure, a tailless peak at 1.8 kb, and a small degradation smear in front of the tailless peak. Two additional impurity peaks (0.9 and 1.1 kb) emerged after 2 days of stress at 37°C and increased further by day 5, possibly through temperature-induced mRNA degradation. Correspondingly, the CPA% of the FLuc mRNA peaks decreased from 91.5% [control] to 80.4% and 68.2%, while impurity levels increased from 8.5% to 19.6% and 31.8%, respectively. Treatment of FLuc mRNA-loaded mRNA-LNP samples at 80°C for 20 minutes similarly reduced the CPA% of the FLuc mRNA peaks to 79.3%, indicative of degradation. The presence of a small amount of FLuc mRNA with higher-order structure is attributed to mRNA refolding during sample handling.

**Table 2. Purity analysis of heat-stressed FLuc mRNA and FLuc mRNA-loaded mRNA-LNP samples in an accelerated stability study.**

FLuc mRNA	Control	80°C, 20 min		
Average CPA%_FLuc mRNA	90.2%	62.6%		
%RSD (n=8)	1.0%	0.7%		
FLuc mRNA-loaded mRNA-LNP	Control	80°C, 20 min	37°C, 2 Days	37°C, 5 Days
Average CPA%_FLuc mRNA	91.5%	79.3%	80.4%	68.2%
%RSD (n=8)	2.5%	2.1%	2.3%	2.5%

## Conclusions

- **Single CE-based workflow for comprehensive mRNA-LNP analysis:** Enables simultaneous assessment of EE%, mRNA purity, and integrity.
- **Accurate, reproducible, and robust automated EE% workflow:** Delivers EE% values [95.6% average] consistent with manual results [95%] across 40 replicates, with strong linearity ( $R^2 \geq 0.995$ ) and reproducibility [%RSD of <0.5%].
- **Scalable, efficient EE% analysis without compromising data quality:** Processes up to 40 samples within 9 hours with consistent performance, while reducing hands-on time and minimizing the risk of RNase contamination.
- **LNP encapsulation enhances mRNA thermal stability and is detectable by the automated assay:** Reveals temperature- and time-dependent degradation, with reduced full-length mRNA loss for encapsulated versus unencapsulated samples under thermal stress.

Acknowledgment: We sincerely thank Acuitas Therapeutics Inc. for providing the FLuc mRNA and FLuc mRNA-loaded mRNA-LNP samples.

## References

1. Makkar, S.K. Advances in RNA-based therapeutics: current breakthroughs, clinical translation, and future perspectives. [Front. Genet. 2025; 16:1675209.](#)
2. Camperi J. et al. Current Analytical Strategies for mRNA-Based Therapeutics. [Molecules 2025, 30\(7\):1629.](#)
3. Sadaqa E. et al. The evolving landscape of mRNA therapeutics. [Biointerface research in applied chemistry. 2026;16\(1\):81.](#)
4. Biomek i5 and Biomek i7 Automated Workstations Instructions for Use. [Beckman Coulter, Inc. P/N: B54473AD \[March 2025\]](#)
5. RNA 9000 Purity & Integrity kit for the BioPhase 8800 system, application guide, [RU0-IDV-05-13438-D.](#)
6. Analytical procedures for the quality of mRNA vaccines and therapeutics [Draft guidelines: third edition]. [2024. USP.](#)
7. Mantri P. et al. Comparison of capillary electrophoresis-based methods for the analytical characterization of purity and stability of in vitro transcribed mRNA. [J Pharm Biomed Anal. 2024 Oct 15:249:116352.](#)
8. Analysis of quality attributes for mRNA vaccines and therapeutics by capillary electrophoresis. SCIEX white paper, [MKT-36904-A.](#)

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product[s] not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries [see [www.sciex.com/trademarks](http://www.sciex.com/trademarks)].

\* SYBR™ is a trademark of the Life Technologies Corporation. SYBR™ Green II RNA gel stain is not available for resale.

Biomek i-Series Automated workstations are not intended or validated for use in the diagnosis of disease or other conditions.

© 2026 DH Tech. Dev. Pte. Ltd. MKT-38275-A.



**Headquarters**  
250 Forest Street, Marlborough,  
MA 01752 USA  
Phone 508-383-7700  
[sciex.com](http://sciex.com)

**International Sales**  
For our office locations please call the division  
headquarters or refer to our website at  
[sciex.com/offices](http://sciex.com/offices)