

Differentiation of aspartic and isoaspartic acid using electron activated dissociation (EAD)

Featuring the SCIEX ZenoTOF 7600 system using EAD and Protein Metrics Inc. software

Zoe Zhang¹, Pavel Ryumin², Takashi Baba², Bill Loyd², Jason Causon², Kerstin Pohl¹

¹SCIEX, USA; ²SCIEX, Canada

The data presented here show the unambiguous identification and localization of aspartic (Asp) and isoaspartic acid (isoAsp) of a biotherapeutic monoclonal antibody (mAb). The differentiation of Asp and isoAsp was performed using rapid data dependent acquisition (DDA) and automated data interpretation with Protein Metrics Inc. software as part of a standard peptide mapping analysis. With this workflow, regular and advanced characterization leveraging electron activated dissociation (EAD)^{1,2} is achievable in one injection, enabling a streamlined characterization accessible to all user levels.

Ensuring drug safety and efficacy is essential for biotherapeutics, which drives the need for in-depth characterization during their development. This includes the identification and the localization of post-translational modifications (PTMs), such as deamidations and the differentiation of Asp and isoAsp. Isomerizations are one of the most difficult PTMs to analyze due to the ambiguity of fragments in the MS/MS spectra. Bottom-up approaches are the method of choice, enabling the simultaneous identification and localization of modifications. Traditionally, isoAsp and Asp are differentiated using chromatographic retention time.³ However, retention times and elution orders can vary depending on the peptides and chromatographic setup being used. With increasingly more complex protein therapeutics in development,

the need to confidently distinguish between Asp and isoAsp has dramatically increased. Furthermore, it is likely that more analytical questions will need to be answered which require more effective and efficient MS tools. Previous techniques of alternative fragmentation suffer from long reaction times, low sensitivity and irreproducibility.

A new fragmentation type based on EAD^{1,2} was used to identify Asp and isoAsp in multiple ambiguous conditions from a commercial mAb. The data were acquired using an untargeted 10 Hz DDA method and interpreted with Protein Metrics Inc. software. This workflow demonstrates a streamlined characterization with a new level of data quality in a routine manner.

Key features of the SCIEX ZenoTOF 7600 system

- **New depths of peptide mapping analysis:** EAD with fast DDA enables alternative fragmentation for routine, in-depth analysis of next generation protein therapeutics and standard mAbs
- **Higher levels of structural information:** Changing the mechanism of fragmentation by tuning the electron energy may provide a higher level of structural information, particularly for isomer differentiation such as Asp and isoAsp differentiation
- **Higher MS/MS sensitivity:** Increased detection of fragments (5 to 10 fold) using the Zeno trap enables higher confidence in data assignment
- **High reproducibility:** Reproducible fragmentation with EAD for singly, doubly, and multiply charged ions enables analysis of more precursors than other alternative and low reproducibility fragmentation techniques
- **Streamlined and easy-to-use:** Fully automated data acquisition in DDA mode using EAD with SCIEX OS software, and automated data interpretation with Byos software (Protein Metrics Inc.) simplifies the entire user experience



Figure 1. The SCIEX ZenoTOF 7600 system.

Methods

Sample preparation: The mAb sample (adalimumab) was denatured with 7.2M guanidine hydrochloride, 100mM Tris buffer pH 7.2, followed by reduction with 10mM DL-dithiothreitol and alkylation with 30mM iodoacetamide. Digestion was performed with trypsin/Lys-C enzyme at 37°C for 16 h.

Chromatography: 10 µL (4 µg) of the trypsin/Lys-C digest were separated with a CSH C18 column (2.1×100 mm, 1.7 µm, 130 Å, Waters) using an ExionLC system. The mobile phase A consisted of water with 0.1% formic acid, while the organic phase B was acetonitrile 0.1% formic acid. A gradient profile was used at a flow rate of 300 µL/min (Table 1). The column temperature was maintained at 50°C.

Table 1. Chromatography for peptide mapping analysis.

Time [min]	Mobile phase A [%]	Mobile phase B [%]
<i>Initial</i>	98	2.0
5.0	98	2.0
6.0	90	10
40	55	45
44	10	90
46	10	90
47	98	2.0
50	98	2.0
51	10	90
54	10	90
55	98	2.0
60	98	2.0

Mass spectrometry: Data were acquired in positive ionization mode with an information dependent acquisition (IDA) method using the SCIEX ZenoTOF 7600 system. The electron energy for the alternative fragmentation in the EAD cell was set to a value of 7 eV. Detailed method parameters are summarized in Table 2.

Table 2. MS parameters.

Parameter	MS	MS/MS
<i>Scan mode</i>	TOF-MS	IDA dependent
<i>Gas 1</i>		50 psi
<i>Gas 2</i>		50 psi
<i>Curtain gas</i>		35 psi
<i>Source temperature</i>		550°C
<i>Ion spray voltage</i>		5500 V
<i>Declustering potential</i>		80 V
<i>Collision energy</i>		12 V
<i>CAD gas</i>		7
<i>Maximum candidate ion</i>		10
<i>Intensity threshold</i>		125 cps
<i>Charge states</i>		2 to 10
<i>Exclusion time</i>		6 s after 2 occurrences
<i>Start mass</i>	200 m/z	100 m/z
<i>Stop mass</i>	2,000 m/z	3,000 m/z
<i>Electron KE</i>	NA	7 eV
<i>Electron beam current</i>	NA	4500 nA
<i>ETC</i>	NA	100
<i>Zeno trap</i>	NA	ON
<i>Accumulation time</i>	0.25 s	0.10 s
<i>Time bins to sum</i>	8	12

Data processing: Data were processed using Byos software (Protein Metrics Inc.).

The what, why and how

The formation of isoAsp has two routes: one is the deamidation of asparagine (Asn) and the other is the isomerization of Asp (Figure 2).⁴ The formation of isoAsp is a non-enzymatic process through a succinimide intermediate, usually increasing with prolonged storage time. Since the isomerization of Asp introduces an additional methylene group in the protein backbone, it can lead to a dramatic impact on the protein folding⁵, which can affect the drug's efficacy and safety. It was reported that isomerization of Asp in the complementary determining regions (CDRs) of antibodies decreased binding efficiency.⁶ PTMs with such an effect on the product are classified as critical quality attributes (CQA) and require comprehensive characterization and monitoring. Bottom-up approaches are the method of choice for characterization of product quality attributes, enabling the simultaneous identification and localization of modifications. Peptide mapping analyses with collision-induced dissociation (CID) are commonly utilized to investigate deamidation sites, as each deamidation will introduce a 0.98 Da mass shift. However, it is a challenge to distinguish isoAsp and Asp as they are isomeric, resulting in the same *b*- and *y*-fragments derived from CID. Therefore, their differentiation often leveraged chromatographic retention times.³ Since the elution order of Asp and isoAsp can change depending on the chromatographic conditions⁶, this indirect approach lacks reliability.

In contrast, alternative fragmentation techniques can produce *c*- and *z*-fragments indicative of isoAsp. After an electron is attached to an isoAsp, electron rearrangements can lead to a bond breakage between the α -carbon and the methylene group in the peptide backbone, producing *c*+57 and *z*-57 fragments (Figure 3).⁷ These ions are not present when fragmenting Asp due to the lack of a methylene group in the peptide backbone. Although alternative fragmentation technology has been used for




Figure 2. Scheme of Asp and isoAsp formation. Asparagine can undergo deamidation resulting in a succinimide intermediate, which can either react further to form Asp or isoAsp through the addition of water.




Figure 3. Formation of *c*+57 ion and *z*-57 ion from isoAsp by EAD.

this type of identification in the past, the application suffered from overall low sensitivity, with the usage of nano-flow regimes and selective reaction monitoring instead of automated DDA workflows.^{4,7} This largely limited the adoption in the biopharmaceutical industry.

With the SCIEX ZenoTOF 7600 system, an alternative fragmentation is introduced, enabling scientists to get an in-depth picture of their samples by using analytical flow liquid chromatography (LC) separation in combination with a fast scanning DDA method and automated processing using Protein Metrics Inc. software. This breakthrough technology realizes the dream of high-throughput fragmentation in the biopharmaceutical industry capable of answering complex questions in a routine manner.

Identification and localization of deamidations

The study focused on the characterization of a commercialized mAb: adalimumab. Multiple potential deamidation sites are present in the sequence of adalimumab and previous high pH stress studies on the intact protein suggest that several sites are highly accessible for deamidation events and isomerization of Asp.⁸ However, the intact mass analysis approach could not reveal the exact localization of these events.

Here, a DDA approach in combination with the patented Zeno EAD was chosen. With this approach, routine peptide mapping analyses can be performed, while the EAD enables advanced characterization in the same, single analysis. Furthermore, the detection of fragment ions and thus the correct identification of low abundant deamidations or Asp isomerizations is enhanced by the Zeno EAD. This approach allowed for the straight-forward data interpretation using Byos software (Protein Metrics Inc.).

Deamidation sites could be accurately identified as shown for the peptide NSLYLQMNSLR (Figure 4). The high MS/MS sequence coverage of >80% with excellent signal-to-noise level provided great confidence in the assignment of both, non-modified and the deamidated species. The modification site could be clearly localized by the mass shift of z4 and c8 (and greater) of 0.98 amu compared to the unmodified peptide.

Differentiation of Asp and isoAsp

The identification and localization of deamidations and associated isomerizations can be further complicated if multiple potential sites are present in the same peptide: Peptide SGTASVVCLLN~~N~~FYPR contains two potential deamidation sites (asparagine; N) right next to each other. The main form without any modification and three deamidated forms were found (Figure 5). Interestingly, two modified forms eluted before the main peak and one thereafter (in Figure 5A and 5B). The Zeno EAD spectra revealed that in two cases the site of modification was found to be the N-11, clearly identified by a mass shift of z6, but not z5 ions, compared to the unmodified peptide. This excludes the N-12 position as the site of modification. The difference of the two modified peptides in Figure 5A and 5B is localized to the side chain of the asparagine, which results in the presence or absence of c10+57 and the z6-57. For the peak eluting prior to the main species, the diagnostic ions for isoAsp (c10+57 and z6-57) were found, confirming that the peptide eluting first contains the isoAsp (Figure 5A). In addition, the deamidated version of the peptide at position N-12 was found to be eluting prior to the main species (Figure 5C). Again, this peptide and its modification was clearly identified by characteristic ions. In this case, the elution time is consistent with previous predictions that isoAsp has an earlier elution time compared to the non-modified peptide, and the deamidated version.¹ However, the elution of isoaspartic acid does not always follow this rule.³ The peptide VVSVLTVLHQDWLNGK in adalimumab shows an example for a deamidation and associated isoAsp formation for which the mentioned elution order is partly reversed (Figure 6). Two deamidation peaks were observed for this peptide, which elute after the non-deamidated peak (Figure 6A and 6B). The automated assignment of the EAD spectra proves that N-14 represents the modification site for both peaks as z-ions greater than z-4 all showed a 0.98 amu mass shift (Figure 6A and 6B). However, the diagnostic ions z3-57 and c13+57 were only detected in the peak eluting right after the main form, confirming isoAsp formation (Figure 5A).

Unambiguous identification, localization and differentiation of deamidations and isoAsp formations in one single DDA run with automatized data interpretation was achieved with the SCIEX ZenoTOF 7600 system using Zeno EAD. This is an example of how distinguishing isomers previously thought a challenge by LC-MS/MS can be simplified in a reproducible manner by EAD in combination with the Zeno trap (Zeno EAD). This strategy can also be applied for further PTM analysis as well as to the differentiation of leucine and isoleucine for full confidence in the correct sequence of biopharmaceuticals in development.^{9,10}




Figure 4. Zeno EAD spectra for peptide NSLYLQMNSLR using Zeno EAD. (A) Deamidated, (B) Non-deamidated.




Figure 5. Data for peptide SGTASVVCLLNFFYPR using Zeno EAD. Extracted ion chromatograms with dots indicating the time point for MS/MS data acquisition on the left and associated EAD spectra on the right hand side for (A) isoaspartic acid formation for position N11, (B) aspartic acid for position N11, (C) aspartic acid for position N12, (D) non-deamidated main form. For (A) diagnostic fragment ions z6-57 and c10+57 were observed (encircled).




Figure 6. Data for peptide VVSVLTVLHQDWLNGK using Zeno EAD. Extracted ion chromatograms with red dots indicating the time point for MS/MS data acquisition on the left and associated EAD spectra on the right hand side for (A) isoaspartic acid (B) aspartic acid (C) non-deamidated main form. For (A) diagnostic fragment ions z3-57 and c13+57 (encircled) were observed.

Conclusions

- The exact localization and differentiation of Asp and isoAsp in peptides with multiple potential sites was achieved with EAD, a novel fragmentation technique
- MS/MS fragment detection was significantly enhanced compared to traditional high resolution MS/MS analysis enabling great data quality for confident fragment assignment even for precursors with medium or very low intensities such as modified peptides utilizing the Zeno trap
- The robust, reproducible and easy-to-use alternative fragmentation enables users to answer challenging analytical questions more streamlined with the SCIEX ZenoTOF 7600 system controlled by SCIEX OS software
- Automatic data processing enables the routine and advanced characterization of complex biotherapeutics and standard mAbs in a reproducible manner using Protein Metrics Inc. software

References

1. Baba T *et al.* (2015) Electron capture dissociation in a branched radio-frequency ion trap, [Anal Chem, 87, 785-792.](#)
2. Baba T *et al.* (2021) Dissociation of biomolecules by an intense low-energy electron beam in a high sensitivity time-of-flight mass spectrometer. *Accepted by JASMS.*
3. Geiger Terrence *et al.* (1987) Deamidation, isomerization, and racemization at asparaginyl and aspartyl residues in peptides. [J Biological Chemistry, 262\(2\):785-797.](#)
4. Ni, Wenqin *et al.* (2011) Analysis of isoaspartic acid by selective proteolysis with Asp-N and electron transfer dissociation mass spectrometry. [Anal Chem, 82\(17\):7485-7491.](#)
5. O'Connor Peter *et al.* (2006) Differentiation of aspartic and isoaspartic acids using electron transfer dissociation. [JASMS, 17\(1\):15-19.](#)
6. Harris Reed *et al.* (2001) Identification of multiple sources of charge heterogeneity in a recombinant antibody. [J Chromatography B, 752\(2001\):233-245.](#)
7. DeGraan-Weber *et al.* (2016) Distinguishing aspartic and isoaspartic acids in peptides by several mass spectrometric fragmentation methods. [J Am Soc Mass Spectrom, 27\(12\):2041-2053.](#)
8. Elucidation and tracking of charge variants in antibodies: combining ease of use with a gain in information. [RUO-MKT-02-12073-B.](#)
9. Comprehensive peptide mapping of biopharmaceuticals utilizing electron activation dissociation (EAD). SCIEX technical note RUO-MKT-02-12639-B.
10. Differentiation of leucine and isoleucine using electron activated dissociation (EAD). SCIEX technical note RUO-MKT-02-12605-B.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-12550-B