

# High-resolution LC-MS/MS solution for improved quantification of peptides in a complex matrix

#### Featuring the ZenoTOF 7600 LC-MS/MS system

Shane Needham<sup>1</sup>, Eshani Nandita<sup>2</sup>, Lei Xiong<sup>2</sup>, Elliott Jones<sup>2</sup>, Zoe Zhang<sup>2</sup>, Kerstin Pohl<sup>2</sup> <sup>1</sup>Alturas Analytics Inc, Moscow, Idaho; <sup>2</sup>SCIEX, USA

Improved selectivity and greater sensitivity were achieved for peptide quantitification in a complex matrix using the ZenoTOF 7600 system featuring the Zeno trap. The greater mass resolution provides selectivity between target peptide versus matrix-derived interferences and general high background effects. This enhanced selectivity aids in superior quantitative results in terms of LOQ, precise data integration, improved reproducibility, and linear dynamic range. In addition, enhanced sensitivity for peptide quantification was achievable using the Zeno trap to improve the duty cycle versus traditional time-offlight systems.

Traditional workflows for quantitative bioanalyses, such as immunological assays, have been displaced by LC-MS/MS analysis on a triple quadrupole mass spectrometer. Immunoassays often lack selectivity, specificity, and have a limited linear dynamic range. While the triple quadrupole platform provides excellent sensitivity and quantitative performance, there can be some limitations with background interference based on the lower resolution considering the type of mass analyzer. Background interference is a common issue for workflows where analytes are present in a highly complex matrix. High-resolution



Figure 1. ZenoTOF 7600 system offers greater selectivity for peptide quantification compared with a triple quadrupole system. XIC of the matrix blank from the triple quadrupole system shows significant matrix interference at retention time of analyte resulting in poor quantitative performance. Overall, the ZenoTOF 7600 system reduces time for method development and optimization.

accurate mass spectrometry (HRAMS) has been increasingly adopted for quantitative bioanalysis.<sup>1,2</sup> The ZenoTOF 7600 system offers an exceptional combination of mass resolution, sensitivity, and acquisition speed for quantitative analysis. It also aids in more accurate and automated integration, the potential for less ion path tuning, the ability to change measured fragments post-acquisition, and improved reproducibility and LDR when interferences are mitigated. These attributes complement the excellent sensitivity of a nominal mass triple quadrupole system such as the SCIEX Triple Quad 7500 LC-MS/MS system – QTRAP Ready or SCIEX Triple Quad 6500+ LC-MS/MS system for a biopharmaceutical lab which requires a full range of capabilities.

The ZenoTOF 7600 system offers a high-resolution MS/MSbased acquisition mode for peptide quantification with Zeno MRM<sup>HR</sup>, along with improved sensitivity using the Zeno trap. The implementation of the Zeno trap allows for improvement of the duty cycle by  $\geq$ 90 %, thereby improving overall MS/MS sensitivity.

## Key features of the ZenoTOF 7600 system for peptide quantification

- Gain higher sensitivity for peptide quantification using the Zeno trap, by enhancing duty cycle through an accumulation of ions during each TOF pulse
- Achieve greater selectivity between target peptides and matrix-related components with the higher mass resolution offered by the ZenoTOF 7600 system
- Ensure exceptional accuracy and precision for quantitative workflows using the ZenoTOF 7600 system
- Achieve superior sensitivity for peptide quantification using HRAMS MS/MS in comparison with single MS mode
- Perform highly automated and accurate peak integration given the improvements in selectivity with higher resolution and ensure overall data integrity
- Easily acquire, process, and manage data on a single platform using the SCIEX OS software

### Methods

**Samples and reagents:** Universal Proteomics Standard (UPS) was purchased from Sigma-Aldrich. Rat plasma (Sprague Dawley, K2 EDTA) was purchased from BioIVT.

**Sample preparation:** The calibration curve was prepared by spiking digested UPS into rat plasma digest followed by serial dilution.

Samples were denatured by incubating with N-octyl-glucoside (OGS), followed by reduction with dithiothreitol (DTT) and alkylation with iodoacetamide (IAM). A trypsin/Lys-C digestion was performed at 37 °C overnight, with an enzyme-protein ratio of 1:25. Formic acid was spiked into the samples to abort digestion. The samples were centrifuged at a speed of 12,000 g and the supernatant was then injected for LC-MS analysis.

Proteins used for this study had limited starting concentrations. Therefore, the final LDRs were narrow for the peptides analyzed.

**Chromatography:** An ExionLC system was used for analyte separation. A volume of 20  $\mu$ L was injected for analysis. Mobile phase A consisted of water with 0.1% FA in water, while organic phase B was composed of 0.1% FA in acetonitrile. For analyte separation, the operating flow rate was set to 0.5 mL/min using a Phenomenex Kinetex C18 column (3 x 50 mm, 2.6  $\mu$ m, 100 Å). The column oven temperature was set to 40 °C.

Chromatographic conditions for analyte separation are shown in Table 1.

#### Table 1. Chromatographic conditions for analyte separation.

Time (min)	Mobile phase A (%)	Mobile phase B (%)		
0.50	98.0	2.0		
1.00	88.0	12.0		
6.00	68.0	32.0		
6.10	10.0	90.0		
7.00	10.0	90.0		
7.10	98.0	2.0		
8.00	98.0	2.0		

*Mass spectrometry:* Samples were analyzed in triplicate. Method details such as source and gas parameters and MS conditions are summarized in Table 2. Sample analysis was performed using scheduled Zeno MRM<sup>HR</sup> on the ZenoTOF 7600 system. The ZenoTOF 7600 system provides a scan speed of 133 Hz. **Data processing:** MRM data were processed using SCIEX OS 2.0 software. Integration was performed using the MQ4 algorithm. Linear regression with 1/x weighting was used for quantification of all peptides. The XIC peak width was set to 0.05 Da for both MS/MS and MS1 quantification.

#### Table 2. Source and MS conditions.

Parameter	Value	Parameter	Value		
Curtain gas	30 psi	Source temperature	550 °C		
lon source gas1	65 psi	lon source gas 2	65 psi		
CAD gas	12	lon spray voltage	5500 V		
MS accumulation time	e 40 ms	MS/MS accumulation time 20 ms			
TOF MS start mass	350 Da	TOF MS stop mass	1500 Da		
TOF MS/MS start mass	300 Da	TOF MS/MS stop mass	≥1000* Da		
ZOD threshold	20,000 cps				

\*TOF MS/MS stop mass depends on peptide analyzed

### Zeno trap provides greater sensitivity

In traditional time-of-flight MS/MS acquisition, ions are often lost in between TOF pulses due to differences in velocity, with typical duty cycle values range approximately between 5-25%.



Figure 2. Zeno trap enables ion beam control from the collision cell prior to entrance into the TOF accelerator. Gains in ion transmission improve overall MS/MS sensitivity.



Due to loss in ion transmission, there is a subsequent decline in overall sensitivity as fewer ions arrive at the detector. The Zeno trap gains back the ion transmission by providing control of the ion beam from the collision cell into the TOF accelerator (Figure 2). Ions exit the Zeno trap in an ordered release based on potential energy.

With the design of the Zeno trap, significant enhancements in MS/MS sensitivity were observed. As a result, considerable improvements in LLOQ levels were achieved using Zeno MRM<sup>HR</sup> in comparison with standard MRM<sup>HR</sup> (Figure 3). From the three peptide examples shown in Figure 3, greater than 8 times improvement in LLOQ sensitivity was observed.



**Figure 3. Improved sensitivity was observed with Zeno trap.** Significantly better LLOQs were achieved with enhancements in duty cycle with Zeno MRM<sup>HR</sup> (right).

# Selectivity enhancement for peptide quantification

Selectivity issues often arise in quantification workflows in the form of matrix interference and high background. Matrix interference is the presence of a defined peak in the matrix blank at the retention time of the analyte where the precursor and fragment ions are the same as that of the target analyte. High background arises as a result of fragment ions that elute at the retention time of the analyte, which often limits the ability to reach trace levels of quantification.

Higher resolution mass spectrometers can maximize selectivity by offering a greater mass resolution between target analyte and any matrix-related components or background ions. With the high resolution capabilities of the ZenoTOF 7600 system, common bioanalysis challenges with matrix interference and high background can be mitigated in comparison with a nominal mass system (Figure 4).

The ZenoTOF 7600 system also achieved strong linearity, accuracy and precision, demonstrating GLP-level quantitative performance (Figure 4, Table 3). Accuracy at the LLOQ was within 80%-120%, while accuracy for all other non-zero calibrators were within 85%-115% of the nominal concentration. Overall, precision was <11%, demonstrating high reproducibility.

## Comparison between MS/MS and MS1 based quantification

In addition to the MS/MS information for MRM<sup>HR</sup> analysis, quantification can also be performed using the precursor ion associated with the target peptide. Implementation of the precursor ion for quantification does not include additional confirmation from any fragment ions to verify the measured peptide structure. Mainly, when the peptide is present in a challenging matrix, precursors might be present that do not correspond to the target peptide, which may hinder quantification by introducing interferences. As a result, the selectivity of the target peptide can be compromised.

In comparison with MS/MS quantification, MS1 based quantification lacks additional selectivity, causing the S/N to often be lower (Figure 6). When MS/MS data is applied for quantification, S/N is typically considerably enhanced with the fragment ion data providing greater selectivity of the target peptide. Therefore, for peptides present in a complex matrix, MS/MS-based quantification offers a highly selective and specific approach.



Figure 4. ZenoTOF 7600 system offered improved mass resolution between the target peptide and matrix-related components. Enhanced selectivity was observed when comparing the XICs of the matrix blank and non-zero calibrator between a triple quadrupole system and the ZenoTOF 7600 system.





Figure 5. Calibration curves show strong linearity on the ZenoTOF 7600 system.

#### Automated and accurate peak integration

Efficiency of integration can often be driven by the presence or absence of matrix interference or high background ions. With a high mass resolution system, target peptide and any matrix interferences and background ions can be easily resolved by mass. With lack of interferences, an XIC with greater selectivity was generated, allowing for automated and accurate peak integration. As shown in Figure 7, co-eluting interferences are separated by greater mass resolution on the ZenoTOF 7600 system when compared to a standard triple quadrupole system.



Figure 6. An example of quantification using fragment ion XIC in comparison to precursor ion XIC. Significantly better S/Ns were observed using fragment ion XICs (left).

For bioanalysis workflows, this can be of great advantage, as it allows for more convenient and efficient peak integration. Overall, data integrity can also be enhanced because there is less of a need for manual integration given the greater selectivity provided by high mass resolution.



Table 3. Accuracy and precision values per concentration level. Excellent accuracy and precision was achieved on the ZenoTOF 7600 system.

O76070_EGVVGAVEK   y5			P02144_VEADIPGHGQEVLIR   y10			P41159_VTGGLDFIPGLHIPLTLSK   y11		
Concentration (ng/mL)	Percent CV (%)	Accuracy (%)	Concentration (ng/mL)	Percent CV (%)	Accuracy (%)	Concentration (ng/mL)	Percent CV (%)	Accuracy (%)
1228.17	3.42	99.93	1363.28	1.46	96.07	1291.73	2.01	97.97
429.86	4.97	100.97	477.15	4.65	108.50	452.10	6.21	105.51
150.45	2.66	96.38	167.00	1.49	109.40	158.24	1.30	103.50
52.66	10.17	104.64	58.45	2.14	98.47	55.38	4.48	93.83
18.43	10.49	98.07	20.46	4.90	92.94	19.38	9.69	94.32
N/A	N/A	N/A	7.16	5.95	94.61	6.78	8.84	104.87



Figure 7. Ease of automated peak integration was observed on the ZenoTOF 7600 system. XICs of the matrix blank and at the LLOQ level are displayed for both systems. Improved mass resolution between target peptide and co-eluting components allowed for a more accurate and automated peak integration.

#### Conclusions

- Higher sensitivity, based on LLOQ levels, was achieved for peptide quantification using the Zeno trap by improving duty cycle through the accumulation of ions during each TOF pulse for enhanced duty cycle
- Greater selectivity was reached between target peptides and matrix-related components with the higher mass resolution offered by the ZenoTOF 7600 system
- GLP compliant accuracy and precision was observed on the ZenoTOF 7600 system, demonstrating GLP-level quantitative performance
- Significantly better S/N was observed for quantification at the MS/MS level in comparison to the MS1 level

 Automated and accurate peak integration was easily attainable on the ZenoTOF 7600 system, with greater mass resolution ensuring overall data integrity



#### References

- Mike-Qingtao Huang, Zhongping (John) Lin, Naidong Weng (2013). <u>Applications of high-resolution MS in bioanalysis.</u> <u>Bioanalysis 5(10):1269-1276</u>.
- Yuan-Qing Xia, Jim Lau, Timothy Olah, Mohammed Jemal (2011). Targeted quantitative bioanalysis in plasma using liquid chromatography/high-resolution accurate mass spectrometry: an evaluation of global selectivity as a function of mass resolving power and extraction window, with comparison of centroid and profile modes. <u>Rapid</u> <u>Communications in Mass Spectrometry **25**(19):2863-2878.</u>

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-13323-A



Headquarters 500 Old Connecticut Path | Framingham, MA 01701 USA Phone 508-383-7700 sciex.com

International Sales For our office locations please call the division headquarters or refer to our website at sciex.com/offices