



Sensitive quantitation of vitamin A and E in human serum

Pierre Negri¹, Fabian Reijn², and Judith Boom²

¹ SCIEX, USA, ² Diagnostix, The Netherlands

This technical note describes a rapid protein precipitation sample preparation procedure and a robust and sensitive LC-MS/MS method using the QTRAP 6500+ system for the quantitation of vitamin A and E in human serum. Excellent linearity was observed across clinically relevant concentrations. Low- $\mu\text{mol/L}$ level sensitivity was achieved at the lowest calibrator with signal-to-noise ratios [S/N] of 53:1 for vitamin A at 0.41 $\mu\text{mol/L}$ and 753:1 for vitamin E at 5.04 $\mu\text{mol/L}$. In addition, the method showed excellent precision and accuracy at low-level concentrations, demonstrating the quantitative performance of the assay.

Key benefits of vitamin A and E analysis from human serum using the QTRAP 6500+ system

Rapid sample preparation: The two fat-soluble vitamins were extracted from human serum samples using a protein precipitation procedure

Low- $\mu\text{mol/L}$ level sensitivity: excellent sensitivity at the lowest calibrator with S/N of 53:1 for vitamin A [0.41 $\mu\text{mol/L}$] and 753:1 for vitamin E [5.04 $\mu\text{mol/L}$]

Excellent linearity: calibration curves for vitamin A and vitamin E showed r^2 values above 0.99 across the calibration range

Excellent quantitative performance: sensitive quantitation of vitamin A and vitamin E was achieved with excellent precision [2.7% for vitamin A and 3.4% for vitamin E] and accuracy 105.4% for vitamin A and 102.9% for vitamin E at the lowest calibrator levels [0.41 $\mu\text{mol/L}$ for vitamin A and 5.04 $\mu\text{mol/L}$ for vitamin E]




Figure 1. Chromatogram of vitamin A (blue) and vitamin E (pink) extracted from serum matrix. Chromatogram of calibration standards in matrix for vitamin A [retinol] at 0.41 $\mu\text{mol/L}$ and vitamin E [alpha-tocopherol] at 5.04 $\mu\text{mol/L}$ shows a S/N of 53.0:1 for vitamin A and 753.3:1 for vitamin E, respectively, based on a peak-to-peak algorithm.

Introduction

Vitamins A and E are fat-soluble vitamins that play an essential role in various bodily functions such as vision, antioxidant protection, and immune system support. Deficiencies in these vitamins are typically caused by malnutrition or underlying health conditions. As a result, the ability to accurately quantify these vitamins is critical to diagnose deficiencies.

Methods

Sample preparation: Vitamin A and vitamin E were extracted from human serum using a protein precipitation procedure. Briefly, 25 μL of blank serum was spiked with vitamin A and vitamin E at six concentration levels and added to a centrifuge tube. Then, 975 μL of an organic deproteinization solution was spiked with deuterated vitamin A and vitamin E internal standards. The tube was vortexed and centrifuged for 5 minutes at 10,000 \times g. The supernatant was transferred to a UPLC vial for analysis.

Liquid chromatography: Chromatographic separation was achieved using a [Phenomenex Kinetex Luna Omega Polar C18 column](#) (50 \times 2.1 mm, 1.6 μm , 00B-4748-AN). Mobile phase A was water, and mobile phase B was methanol. The total run time was 6.5 min, and the flow rate was 600 $\mu\text{L}/\text{min}$. The injection volume was 15 μL . The LC gradient program is presented in **Table 1**.

Table 1: Chromatographic gradient for the analysis of vitamin A and E in human serum.

Time [min]	Mobile phase A [%]	Mobile phase B [%]
0.0	30	70
1.6	5	95
2.5	5	95
2.51	0	100
4.6	0	100
4.61	30	70
6.5	30	70

Mass spectrometry: Data was collected using a [QTRAP 6500+ system](#) with an IonDrive Turbo V ion source and operated in electrospray ionization (ESI) positive mode. The Scheduled MRM algorithm was used in [SCIEX OS software](#) (version 3.1.6) to

collect sufficient data points to produce robust, quantifiable data. Source and gas conditions are presented in **Table 2**. Compound-dependent parameters were optimized by infusion.

Table 2: Source and gas parameters for the analysis of vitamin A and E in human serum using the QTRAP 6500+ system.

Parameter	Value
Polarity	Positive
Ion source gas 1	60 psi
Ion source gas 2	70 psi
Curtain gas	30 psi
Source temperature	500 $^{\circ}\text{C}$
Ion spray voltage	5250 V
CAD gas	High

Data processing: Data processing was performed using [SCIEX OS software](#) (version 3.1.6). Peak integration was achieved using the MQ4 algorithm. Quantitative analysis was conducted in the Analytics module of SCIEX OS, where calibration curves, concentration calculations, assay precision, and accuracy statistics were automatically generated.

Results and discussion

Figure 1 shows the chromatographic separation of vitamin A and E in a control human serum sample at a final concentration of 0.41 $\mu\text{mol}/\text{L}$ for vitamin A and 5.04 $\mu\text{mol}/\text{L}$ for vitamin E, respectively. The extracted ion chromatograms showed a S/N of 53:1 for vitamin A and 753:1 for vitamin E, at the lowest matrix calibrator measured (0.41 $\mu\text{mol}/\text{L}$ for vitamin A and 5.04 $\mu\text{mol}/\text{L}$ for vitamin E), as calculated using the peak-to-peak algorithm.

Figure 2 shows the representative extracted ion chromatograms [XICs] for A) vitamin A and B) vitamin E across their respective concentration ranges (0.41-4.51 $\mu\text{mol}/\text{L}$ for vitamin A and 5.04-49.50 $\mu\text{mol}/\text{L}$ for vitamin E). The signals shown for the lowest matrix calibrators measured (0.41 $\mu\text{mol}/\text{L}$ for vitamin A and 5.04 $\mu\text{mol}/\text{L}$ for vitamin E) are well above the blank signal for both analytes.

The quantitative performance of the method was investigated by injecting a series of calibrator samples spiked at concentrations ranging from 0.41-4.51 $\mu\text{mol}/\text{L}$ for vitamin A and 5.04-49.50 $\mu\text{mol}/\text{L}$ for vitamin E, respectively. Linearity, accuracy, and precision were assessed across the calibration ranges for the two analytes. **Figure 3** shows the calibration




Figure 2. Extracted ion chromatograms for the two analytes targeted in this method. The XICs show the signal for A) vitamin A and B) vitamin E across their respective concentration ranges [0.41-4.51 µmol/L for vitamin A and 5.04-49.50 µmol/L for vitamin E].

curves for vitamin A [left] and vitamin E [right] over the analytes' respective calibration ranges. The plots show excellent linear responses across the calibration series, with r^2 values greater than 0.99 for both analytes.

The accuracy and precision values were calculated using 6 replicates of the lowest matrix calibrators analyzed [0.41 µmol/L for vitamin A and 5.04 µmol/L for vitamin E]. The accuracy was 105.4% for vitamin A and 102.9% for vitamin E, respectively. The precision [%CV] was 2.7% for vitamin A and 3.4% for vitamin E, respectively.

Conclusions

A fast and sensitive LC-MS/MS method for detecting vitamin A and E extracted from human serum samples was developed using the QTRAP 6500+ system. The method demonstrated:

- Fast sample preparation, which consisted of a simple protein deproteination, requiring 25 µL of human serum sample
- Excellent sensitivity at the lowest calibrator level, resulting in S/N of 53:1 for vitamin A at 0.41 µmol/L and 753:1 for vitamin E at 5.04 µmol/L
- Excellent linear responses across the calibration series consisting of 5 calibrators, with r^2 values greater than 0.99 for both analytes
- High quantitation performance of the method, resulting in excellent precision [2.7% for vitamin A and 3.4% for vitamin E] and accuracy 105.4% for vitamin A and 102.9% for vitamin E at the lowest calibrator levels [0.41 µmol/L for vitamin A and 5.04 µmol/L for vitamin E]




Figure 3. Linear calibration curves for vitamin A (left) and vitamin E (right) extracted from the serum matrix. The calibration curves were run across the following concentration ranges [0.41-4.51 $\mu\text{mol/L}$ for vitamin A and 5.04-49.50 $\mu\text{mol/L}$ for vitamin E]. The curves were generated using linear regression and 1/x weighting for vitamin A and vitamin E in serum, resulting in r^2 values of 0.9999 for vitamin A and 0.998 for vitamin E, respectively.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries [see www.sciex.com/trademarks].

© 2025 DH Tech. Dev. Pte. Ltd. MKT-34904-A



Headquarters
 250 Forest Street,
 Marlborough,
 MA 01752 USA
 Phone 508-383-
 7700
sciex.com

International Sales
 For our office locations please call the division
 headquarters or refer to our website at
sciex.com/offices