



Quantitative analysis of 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ in human serum

This technical note demonstrates the sensitive and accurate quantitation of 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ in human serum using a liquid-liquid extraction procedure with analysis using the SCIEX QTRAP 4500 system. Excellent linearity was observed across clinically relevant concentrations with r^2 values above 0.99 across the calibration range [1-100 ng/mL]. The method showed excellent accuracy (93.5% for 25-OH-vitamin D₃, 82.7% for 3-epi-25-OH-vitamin D₃ and 119.3% for 25-OH-vitamin D₂) at the lowest calibrator level (1 ng/mL), highlighting the assay's strong quantitative performance.

Key benefits of 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ analysis from human serum using the QTRAP 4500 system

Chromatographic separation: Optimized LC conditions enabled fast, 6 min chromatographic separation of 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂

Excellent linearity. Calibration curves for 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ showed r^2 values above 0.99 across the calibration range [1-100 ng/mL]

Low-ng/mL level sensitivity and excellent quantitative performance. Sensitive quantitation of 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ was achieved with excellent accuracy (93.5% for 25-OH-vitamin D₃, 82.7% for 3-epi-25-OH-vitamin D₃ and 119.3% for 25-OH-vitamin D₂) at the lowest calibrator levels [1 ng/mL] in human serum




Figure 1. Overlaid extracted ion chromatograms (XICs) of 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-Vitamin D₂ extracted from serum matrix and analyzed using the SCIEX Triple Quad 4500 system. The 6-min gradient in combination with the column selection and the choice of mobile phase composition resulted in baseline separation of the 25-OH-vitamin D₃ and 3-epi-25-OH-vitamin D₃ epimers.

Introduction

5-OH-vitamin D₃ and 25-OH-vitamin D₂ are the main circulating forms of vitamin D used to assess vitamin D status, with 25-OH-D₃ being the predominant and more biologically active form; 3-epi-25-OH-vitamin D₃ is a less active epimer that may interfere with accurate measurement in some assays. Deficiencies in these metabolites can lead to impaired calcium absorption and bone demineralization, and, as a result, the ability to accurately quantify vitamin D is critical.

Methods

Sample preparation: 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ were extracted from human serum using a liquid-liquid extraction procedure. Briefly, 200 µL of blank serum spiked with 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ at seven concentration levels were vortexed and then extracted by a liquid-liquid extraction (LLE) method using methyl-tertiary-butyl ether (MTBE). Following mixing and centrifugation, the organic layer was separated by snap-freezing and evaporated to dryness. It was then reconstituted in 175 µL of methanol/water [80:20, v/v].

Liquid chromatography: Chromatographic separation was achieved using a [Phenomenex Kinetex PFP analytical column](#) [100 x 3.0 mm, 2.6 µm, 00D-4477-Y0]. Mobile phase A was formic acid in water and mobile phase B was formic acid in methanol. The total run time was 6 minutes at a flow rate of 400 µL/min.

Mass spectrometry: Data was collected using a [QTRAP 4500 system](#) with an IonDrive Turbo V source and operated in electrospray ionization (ESI) positive mode. The scheduled MRM algorithm was used in [SCIEX OS software](#) [version 3.1.6] to collect 10-12 data points to produce robust, quantifiable data. Compound-dependent parameters were optimized by infusion.

Data processing: Data processing was performed using [SCIEX OS software](#) [version 3.1.6]. Peak integration was achieved using the MQ4 algorithm. Quantitative analysis was conducted in the Analytics module of SCIEX OS, where calibration curves, concentration calculations, assay accuracy statistics were automatically generated.

Results and discussion

Figure 1 shows the chromatographic separation of 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ in a control human serum sample. The combination of the 6 min gradient, selected column, and optimized mobile phase composition enabled baseline separation of the 25-OH-vitamin D₃ and 3-epi-25-OH-vitamin D₃ epimers.

The quantitative performance of the method was investigated by injecting a series of serum calibrator samples spiked at concentrations ranging from 1-100 ng/mL. Linearity and accuracy were assessed across the calibration ranges for each of the three analytes. **Figure 2** shows the calibration curves for A) 25-OH-vitamin D₃, B) 3-epi-25-OH-vitamin D₃ and C) 25-OH-vitamin D₂. The plots show excellent linear responses across the calibration series, with r² values greater than 0.99 for all three analytes.

The mean accuracy was evaluated in 3 matrix spike replicates at the lowest matrix calibrators measured [1 ng/mL]. The accuracy was 93.5% for 25-OH-vitamin D₃, 82.7% for 3-epi-25-OH-vitamin D₃ and 119.3% for 25-OH-vitamin D₂ at the lowest calibrator levels [1 ng/mL].

Conclusion

A fast and sensitive LC-MS/MS method for the detection of 25-OH-vitamin D₃, 3-epi-25-OH-vitamin D₃ and 25-OH-vitamin D₂ extracted from human serum samples was developed using the 4500 system. The method demonstrated:

- Chromatographic separation of 1,25-dihydroxyvitamin D₂ and D₃
- Excellent linear responses across the calibration series consisting of 7 calibrators, with r² values greater than 0.99 for all three analytes
- High quantitation performance of the method, resulting in excellent accuracy [93.5% for 25-OH-vitamin D₃, 82.7% for 3-epi-25-OH-vitamin D₃ and 119.3% for 25-OH-vitamin D₂] at the lowest calibrator levels [1 ng/mL]




Figure 2. Linear calibration curves for A) 25-OH-vitamin D₃, B) 3-epi-25-OH-vitamin D₃, and C) 25-OH-Vitamin D₂ extracted from serum matrix. The calibration curves using the matrix calibrators were run across the concentration range [1-100 ng/mL]. The curves were generated using linear regression and 1/x weighting, resulting in a r² value of 0.9966 for 25-OH-vitamin D₃, 0.9939 for 3-epi-25-OH-vitamin D₃ and 0.9958 for 25-OH-vitamin D₂.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries [see www.sciex.com/trademarks].

© 2025 DH Tech. Dev. Pte. Ltd. MKT-36051-A



Headquarters
250 Forest Street, Marlborough,
MA 01752 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division
headquarters or refer to our website at
sciex.com/offices