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¹ CE & Biopharma, Sciex, China; ² Beijing BioCEart Technology Institute, China; ³ Division of Recombinant Biological Products, National Institute of Food and Drug Control, Beijing, China

Key Word: Erythropoietin; Capillary zone electrophoresis; Rapid; Variants

Abstract: Erythropoietin (EPO) is a highly glycosylated glycoprotein hormone secreted by the kidney and liver, whose main function is to promote the production of red blood cells. EPO has a high content of glycosylation modification, leading to the microheterogeneity of EPO molecules and becoming different variants (isoforms) carrying different glycosylated chain. The type and content of variants are important parameters for quality control of EPO products. Capillary zone electrophoresis (CZE) is a basic separation model in capillary electrophoresis, which separates different analytes on the basis of the difference of charge/volume. Owing to the differences of glycosylation sites, number and glycoform of EPO variants, they can be isolated in CZE. Currently, the European pharmacopoeia uses CZE (Called EP-CZE) to characterize EPO variants. However, the capillary condition and separation will take a long time, and the repeatability of migration time is poor in traditional EP-CZE method. In this study, we developed a rapid CZE method (Called Rapid-CZE), which can realize rapid analysis of EPO variants and provides better repeatability. Compared with the EP-CZE method, the number and relative content of EPO variants are consistent. The resolution of the 5th and 6th variants is 1.39, which meet the requirements of European pharmacopoeia. The Rapid-CZE method has good repeatability. The RSD of migration time of each variant is about 1%, and the RSD of the corrected peak area% of the main variants are less than 5%. The Rapid-CZE method has faster separation time and capillary condition time than the EP-CZE method, which is very appropriate for the rapid analysis of EPO variants.

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© 2019. AB Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.
AB Sciex is operating as SCIEX.