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## 参考文献（reference）：

- 1 Tamika T. J. Williams and H. Perreault , "Selective detection of nitrated polycyclic aromatic hydrocarbons by electrospray ionization mass spectrometry and constant neutral loss scanning" Rapid Commun. Mass Spectrom. 14, 1474-1481 (2000)
- 2 R.P. Barreto, F.C. Albuquerque, Annibal D. Pereira Netto, "Optimization of an improved analytical method for the determination of 1-nitropyrene in milligram diesel soot samples by high-performance liquid chromatography-mass spectrometry", Journal of Chromatography A, 1163 (2007) 219-227
- 3 Mahmoud Bataineh, Urte Lübcke-von Varel, Heiko Hayen, and Werner Bracka, "HPLC/APCI-FTICR-MS as a Tool for Identification of Partial Polar Mutagenic Compounds in Effect-Directed Analysis", J Am Soc Mass Spectrom 2010, 21, 1016-1027
- 4 Zoltán Nyiri, Márton Novák, Zsolt Bodai, Bálint Sámuel Szabó, Zsuzsanna Eke, Gyula Záray, Tamás Szigeti, "Determination of particulate phase polycyclic aromatic hydrocarbons and their nitrated and oxygenated derivatives using gas chromatography - mass spectrometry and liquid chromatography - tandem mass spectrometry", J Chromatogr A. 2016 Nov 11;1472:88-98
- 5 Mahmoud Bataineh, Alya A. Arabi, Jibrán Iqbal, Fares M. Howari, and Werner Brack, "Method Development for Selective and Nontargeted Identification of Nitro Compounds in Diesel Particulate Matter", Energy Fuels, 2017, 31 (11), pp 11615-11626

- 6 María Guiñez , Cristian Bazan , Luis D. Martinez , Soledad Cerutti, "Determination of nitrated and oxygenated polycyclic aromatic hydrocarbons in water samples by a liquid-liquid phase microextraction procedure based on the solidification of a floating organic drop followed by solvent assisted back-extraction and liquid chromatography-tandem mass spectrometry", *Microchemical Journal* 139 (2018) 164-173
- 7 Christian Schauer · Reinhard Niessner · Ulrich Pöschl, "Analysis of nitrated polycyclic aromatic hydrocarbons by liquid chromatography with fluorescence and mass spectrometry detection: air particulate matter, soot, and reaction product studies", *Anal Bioanal Chem* (2004) 378 : 725-736
- 8 Crystal Cheely, Tracy Hefner etc. 卷烟主流烟气中硝基-多环芳香烃的分析, 第58届烟草科学研究会议论文集

For Research Use Only. Not for use in Diagnostics Procedures.

AB Sciex is operating as SCIEX.

© 2019. AB Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

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