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- [1] GB/T 5749-2006《生活饮用水卫生标准》[S]. 中华人民共和国卫生部  
GB/T 5749-2006 Standard for drinking water quality[S] Ministry of Public Health, PRC
- [2] GB/T 14848-2017《地下水环境质量标准》[S]. 全国国土资源标准化技术委员会  
GB/T 14848-2017 Standard for groundwater quality[S] National technical committee for standardization of land and resources
- [3] 王敏, 李广益, 高慧, 等. 高效液相色谱-柱后衍生法测定水中草甘膦[J]. 食品科技, 2017, 42(02):295-297.  
WANG M, LI GY, GAO H, et al. Determination of glyphosate in water by high performance liquid chromatography with post column derivatization[J]. Food Science and Technology, 2017, 42(02):295-297 (in Chinese).
- [4] IOS 21458, 2008. Water quality-Determination of glyphosate and AMPA-Method using high performance liquid chromatography (HPLC) and fluorometric detection[s].
- [5] 王天玉, 武秀停, 赵静, 等. 固相萃取-气相色谱法测定薏苡仁、白茯苓和山药中的草甘膦残留[J]. 分析科学学报, 2014, 30(01):63-66.  
WANG TY, WU XT, ZHAO J, et al. Determination of Glyphosate Residue in Semen Coicis, Poria Cocos and Rhizoma Dioscoreae by Solid Phase Extraction-Gas Chromatography[J]. Journal of Analytical Science, 2014, 30(01):63-66 (in Chinese).
- [6] Xie Peijin, Li Gang, Yang Hongda, et al. Identification of Glyphosate Herbicide with Ion Chromatographic Column Switching Method[J]. Plant Diseases and Pests, 2017, 8(06):41-44.
- [7] 韩宛汝. 离子色谱法测定水中草甘膦[J]. 供水技术, 2018, 12(01): 56-58.  
HAN WR. Determination of glyphosate in water by ion chromatography[J]. Water Technology, 2018, 12(01):56-58 (in Chinese).
- [8] Santilio Angela, Pompili Chiara, Giambenedetti Arianna, et al. Determination of glyphosate residue in maize and rice using a fast and easy method involving liquid chromatography-mass spectrometry (LC-MS/MS). Journal of Environmental Science and Health, Part B, 2019, 54 (3), 205-210.
- [9] Luisina Delma Demonte, Nicolás Michlig, Monica Gaggiotti, et al. Determination of glyphosate, AMPA and glufosinate in dairy farm water from Argentina using a simplified UHPLC-MS/MS method. Science of the Total Environment, 2018, 645, 34-43.

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© 2019 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-10719-ZH-A