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SCIEX Tandem Mass Spectrometry Method for Quantification of Chondroitin Sulfate and Dermatan Sulfate by Methanolysis Derivatization

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Jing Lianyan, Li Guoqing

Keywords: Chondroitin Sulfate, Dermatan Sulfate, Methanolysis Derivatization, LC-MS/MS

Abstract

Here we report on an alternative high-performance liquid chromatography tandem mass spectrometry (HPLC-MS/MS) method for measurement of chondroitin sulfate and dermatan sulfate. This method describes the quantification of derivatized dimers chemically derived from DS and CS. Dried urine samples are derivatized using methanolic HCl. These products are dried and mixed with internal standards that are deuterio-labeled dimers prepared by deuteriomethanolysis of CS standards. The CS and DS-derived dimers are individually quantified by comparing the signals of the GAG-specific dimer products with the internal standard, and conversion to concentration by means of calibration curves. The accuracy and precision values for all QC samples of DS and CS-derived dimers met the acceptance criteria. This validated HPLC-MS/MS method was successfully applied to the quantitation of DS and CS.

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- [1] E. Neufeld, J. Muenzer, The mucopolysaccharidoses, in: C.R. Scriver, A.L. Beaudet, W.S. Sly, D. Valle (Eds.), The Metabolic and Molecular Bases of Inherited Disease.
- [2] J. Muenzer, The mucopolysaccharidoses: a heterogeneous group of disorders with variable pediatric presentations.
- [3] H. Zhou, P. Fernhoff, R.F. Vogt, Newborn bloodspot screening for lysosomal storage disorders.
- [4] H Zhang, S P Young, D S Millington. Quantification of glycosaminoglycans in urine by isotope-dilution liquid chromatography-electrospray ionization tandem mass spectrometry. [J]. Current Protocols in Human Genetics.

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