



# SCIEX OS calculated columns and custom flagging for PFAS applications: EPA methods 533 and 1633

Craig M. Butt<sup>1</sup>, Megumi Shimizu<sup>1</sup>, Karl Oetjen<sup>1</sup>  
SCIEX, USA

This technical note demonstrates the formulation of several quality control (QC) calculations and flagging rules directly within the SCIEX OS software for EPA PFAS Methods 533 and 1633. Using the Calculated Columns feature in the processing method, SCIEX OS eliminates the need to export the data to a third-party software platform, minimizing data processing and review time [Figure 1]. In addition, the development of Custom Flagging rules within SCIEX OS allows for the quick review of data outside of the EPA criteria. The flexibility of SCIEX OS is shown, highlighting the ability to adapt to changing QC requirements.

## Key benefits of using SCIEX OS calculated columns and custom flagging for EPA PFAS methods 533 and 1633

- **Flexibility to build QC calculations for EPA PFAS methods.** Calculated Columns feature in SCIEX OS was used to build extensive QC equations for EPA methods 533 and 1633; calculation output directly within Results Table
- **Rapidly identify out-of-range QC samples.** Custom Flagging rules were developed in SCIEX OS to flag QC samples that were outside of acceptance criteria range
- **Data processing and review time saved.** Ability to formulate QC calculations and flagging rules directly in SCIEX OS reduces the time spent on data processing and review

The figure illustrates the application of SCIEX OS for PFAS analysis. It includes chemical structures of PFAS compounds, a workflow diagram, and two screenshots of the software interface.

**[MQ4] Modify Method - Calculated Columns:** This screenshot shows the 'Calculated Columns' tab. The 'Formula name' is 'EPA QC calculation'. The 'Description' is 'Type a description for the formula'. The 'Formula' field contains the following table:

COUNT	MAX	STDEV	Clear
SUM	MIN	MEDIAN	(
MEAN	ABS	IF	)
LEFT	RIGHT	ISNUMBER	+

The 'Formula Details' section shows 'I:0' and a search bar for columns.

**[MQ4] Modify Method - Flagging Rules:** This screenshot shows the 'Flagging Rules' tab. The 'Rule name' is 'EPA QC flagging rule'. The 'Flag a results column' is 'Select a results column'. The 'Flagging criteria' is 'Range'. The 'Step 1: Define the values for the flagging criteria' section shows 'Value for all components' selected, with 'Lower limit' and 'Upper limit' fields. The 'Step 2: Apply the values to the following sample types' section shows 'Unknowns' and 'Standards' selected, with 'Only if the sample name contains...' options.

Figure 1. Screenshots of the Calculated Columns and Flagging Rules features in SCIEX OS

## Introduction

The US EPA has published several LC-MS/MS methods for PFAS analysis in environmental samples, including EPA Method 533 for drinking water and EPA Method 1633 for non-potable water, solids (soil, biosolids, sediment), and tissue samples.<sup>1,2</sup> Both methods contain extensive criteria for initial and ongoing calibration verification, qualitative peak identification and quantitative determination. Performing the required quality control calculations through exporting the raw data into third-party software packages is time-consuming. Further, flagging out-of-range values is labor-intensive and potentially susceptible to error. This technical note highlights the capability of the SCIEX OS software to build custom calculations and flagging rules directly within the data results table, achieving the criteria requirements of EPA Methods 533 and 1633 while reducing the burden of data processing for the analyst.

## Methods

**Data acquisition.** Samples were prepared and analyzed according to EPA Methods 533<sup>1</sup> and 1633<sup>2</sup>, and full details are provided in their individual technical notes<sup>3,4</sup>. Extracts for both methods were analyzed using the SCIEX 5500+ system.

**Data processing:** Data processing was performed in SCIEX OS software, version 3.4.

## EPA 1633: Ion abundance ratio

**Ion abundance ratio.** EPA Method 1633 specifies that the ion abundance ratio (IAR) must be between 50% to 150% of the IAR in either the mid-point calibration standard or the calibration verification (CV) standard injected at the beginning of the analytical batch [Sections 14.3.5 and 15.1.3]. The IAR requirement does not apply to PFBA, PFPeA, NMeFOSE, NEtFOSE, PFMPA and PFMBPA since these compounds exhibit only one stable transition. Since, the default SCIEX OS calculation uses all samples defined as “standards” for the ion ratio equation, a custom calculation is necessary to ensure that the EPA 1633 calculation is only applied to the specific samples.

**Step 1.** The “Group” column must be completed in the Components pane (**Figure 2**). In the SCIEX OS software, the first row in a group represents the quantifying transition, and the second represents the qualifying transition.

**Step 2.** In the Calculated Columns pane, the equation is entered to calculate the sample IAR relative to the mid-point calibration standard [specified as “Cal 4” in this example] or the CV standard from the beginning of the batch [defined as “CCV”] (**Figure 3**).

**Step 3.** In the Flagging Rules pane, the custom rule is created to flag IARs outside of +/- 50% acceptance criteria (**Figure 4**).

Row	IS	Group	Name	Precursor (Q1) Mass (Da)	Fragment (Q3) Mass (Da)	XIC Width (Da)	Retention Time Mode	Retention Time (min)	IS Name	Experiment Index
1	<input type="checkbox"/>	PFBA	PFBA	213	169	0.02	RT value	3.200	EIS-PFBA-13C4	1 MRM (115 transitions)
2	<input type="checkbox"/>	PFPeA	PFPeA_1	263	219	0.02	RT value	4.346	EIS-PFPeA-13C5	1 MRM (115 transitions)
3	<input type="checkbox"/>	PFPeA	PFPeA_2	263	68.9	0.02	RT value	4.392	EIS-PFPeA-13C5	1 MRM (115 transitions)
4	<input type="checkbox"/>	PFHxA	PFHxA_1	313	269	0.02	RT value	5.031	EIS-PFHxA-13C5	1 MRM (115 transitions)
5	<input type="checkbox"/>	PFHxA	PFHxA_2	313	118.9	0.02	RT value	5.031	EIS-PFHxA-13C5	1 MRM (115 transitions)
6	<input type="checkbox"/>	PFHpA	PFHpA_1	363.1	319	0.02	RT value	6.067	EIS-PFHpA-13C4	1 MRM (115 transitions)
7	<input type="checkbox"/>	PFHpA	PFHpA_2	363.1	169	0.02	RT value	6.060	EIS-PFHpA-13C4	1 MRM (115 transitions)
8	<input type="checkbox"/>	PFOA	PFOA_1	413	369	0.02	RT value	7.837	EIS-PFOA-13C8	1 MRM (115 transitions)
9	<input type="checkbox"/>	PFOA	PFOA_2	413	169	0.02	RT value	7.837	EIS-PFOA-13C8	1 MRM (115 transitions)
10	<input type="checkbox"/>	PFNA	PFNA_1	463	419	0.02	RT value	10.132	EIS-PFNA-13C9	1 MRM (115 transitions)
11	<input type="checkbox"/>	PFNA	PFNA_2	463	219	0.02	RT value	10.125	EIS-PFNA-13C9	1 MRM (115 transitions)
12	<input type="checkbox"/>	PFDA	PFDA_1	512.9	469	0.02	RT value	10.394	EIS-PFDA-13C6	1 MRM (115 transitions)
13	<input type="checkbox"/>	PFDA	PFDA_2	512.9	219	0.02	RT value	10.394	EIS-PFDA-13C6	1 MRM (115 transitions)
14	<input type="checkbox"/>	PFUnA	PFUnA_1	563.1	519	0.02	RT value	10.517	EIS-PFUnA-13C7	1 MRM (115 transitions)

Figure 2. Components pane of the processing method in SCIEX OS Analytics showing the “Group” column [labelled as “Component Group Name” in the Results table]. The Group column must be completed for several of the EPA 1633 calculations to be performed.

[MQ4] Modify Method

Workflow Components Integration Library Search Calculated Columns Flagging Rules Advanced Formula Finder Non-targeted Peaks

← Accept changes and return to Calculated Columns X Discard

Use the calculator to create a new formula.

Formula name: ion ratio / cal 4 ion ratio

Description: Type a description for the formula

Formula: 
$$\text{ABS}(1 - ([\text{Ion Ratio}] / \text{MEAN}([\text{Ion Ratio}]))) * 100$$

Note: The "Original text" option is recommended for formulas that contain functions, such as the IF function, that compare non-numeric values to numeric values.

Treat resulting text values as: Error (N/A)

Select appropriate sample type:

☒ Standards  
☒ Only if the sample name contains... Cal4

☒ QCs  
☒ Only if the sample name contains... CCV

Formula Details

Columns: Search

Regression parameters

MEAN value will be calculated using the following sample types:

☐ Unknowns  
☐ Only if the sample name contains... Type comma separated text

☐ Standards  
☐ Only if the sample name contains... Cal4

☐ QCs  
☐ Only if the sample name contains... CCV

☐ Blanks  
☐ Solvent ☐ Blank ☐ Double blank  
☐ Only if the sample name contains... Type comma separated text

Process & Close Print Close Help

Figure 3. The Calculated Columns pane of the SCIEX OS processing method showing the IAR calculation. Either the "Standards" or "QCs" sample type is selected, depending on the choice of reference sample used. The user also needs to check the "Only if the sample name contains ..." box to input the specific sample name used for the selected reference sample.

[MQ4] Modify Method

Workflow Components Integration Library Search Calculated Columns Flagging Rules Advanced Formula Finder Non-targeted Peaks

← Accept changes and return to Flagging Rules X Discard

Rule name: Ion Ratio 50%

Flag a results column: ion ratio / cal 4 ion ratio

Flagging criteria: Upper limit

Step 1: Define the values for the flagging criteria

☒ Value for all components ☐ Values per component type

Upper limit: 50

By component

Step 2: Apply the values to the following sample types

☒ Unknowns  
☐ Only if the sample name contains... Type comma separated text

☒ Standards  
☐ Only if the sample name contains... Type comma separated text

☒ QCs  
☐ Only if the sample name contains... Type comma separated text

☐ Blanks  
☐ Solvent ☐ Blank ☐ Double blank  
☐ Only if the sample name contains... Type comma separated text

Treat resulting text values as: Error

Process & Close Print Close Help

Figure 4. The Flagging Rules pane of the SCIEX OS processing method showing the development of the custom rule for flagging IARs outside of +/- 50% acceptance criteria. The flag is applied to the custom-built IAR calculated column in the previous step and is applied to unknowns, standards and QCs.

## EPA 1633: Retention time verification

EPA Method 1633 specifies that the retention time [RT] for all target analytes with exact corresponding stable-isotope analogues must be within  $\pm 0.1$  min of the associated extracted internal standard [EIS]. This criterion is applicable to 24 out of the 40 target PFAS compounds.

**Step 1.** A conditional lookup table is created under the Calculated Columns pane (Add > Conditional lookup). The “Component Group Name” and “Equals” are specified under the “Column” and “Condition” headings, respectively. The PFAS compounds with exact EIS matches are manually entered. The “Output” column is entered as “true”, “Default output” is “false” [or left empty]. This table creates a new results table

column (“ISD EIS Match”) which identifies PFAS analytes with exact EIS matches as “true” and is used in Step 2 for the RT verification calculation [Figure 5].

**Step 2.** In the Calculated Columns pane, the equation is entered to calculate the analyte RT and the EIS RT difference, if the EIS is an exact match (indicated by “true” in the “ISD EIS Match” column). The output column is “RT verification” [Figure 5].

**Step 3.** In the Flagging Rules pane, a custom rule is added (“RT verification”) to flag instances when the analyte and exact match EIS RTs exceed  $\pm 0.1$  min. Specifically, the “Flagging criteria” is set to “Upper limit”, the “Value for all components” is set to “Upper limit = 0.1” and all sample types are selected [Figure 6].

**[MQ4] Modify Method**

Workflow | Components | Integration | Library Search | **Calculated Columns** | Flagging Rules | Advanced | Formula Finder | Non-targeted Peaks

Accept changes and return to Calculated Columns | Discard

Name: ISD EIS match  
Description: specifies which analyte groups have matching EIS

Table Properties: Column count 1 | Row count 25

Select All Rows | Delete Selected Rows

Column	Condition	Output
Component Group Name	Equals	
PFBA		true
PFPeA		true
PFHxA		true
PFHpA		true
PFOA		true
PFNA		true

Default output: false

Process & Close | Print | Close | Help

### Conditional lookup table



### Calculated column

**[MQ4] Modify Method**

Workflow | Components | Integration | Library Search | **Calculated Columns** | Flagging Rules | Advanced | Formula Finder | Non-targeted Peaks

Accept changes and return to Calculated Columns | Discard

Use the calculator to create a new formula.

Formula name: RT verification  
Description: Type a description for the formula

COUNT	MAX	STDEV	Clear
SUM	MIN	MEDIAN	(
MEAN	ABS	IF	)
LEFT	RIGHT	ISNUMBER	+
/	*	-	=

Formula editor: IF([ISD EIS match]=true; ABS([Retention Time]-[IS Retention Time]); 'false')

I:18

Formula Details

Columns: Search | Regression paramet

IS Width at 50%

**Figure 5. Conditional lookup table and Calculated Columns equation for the development of the EPA 1633 retention time verification calculation.** The lookup table creates a new Result Table column which designates PFAS analytes with exact EIS matches as “true”. The Calculated column determines the analyte RT and the EIS RT difference, if the EIS is an exact match.

Figure 6. The custom flagging rule to flag instances where the analyte retention time and the associated EIS retention times are >0.1 min. The flag is applied to the custom-built calculation in the previous step.

## EPA 533: Laboratory fortified sample matrix [LFSM] recovery

**Laboratory fortified sample matrix [LFSM] recovery.** EPA Method 533 specifies that at least one LFSM sample must be prepared with each extraction batch [section 9.2.6]. For spikes at concentrations  $\leq 2\times$  of the minimum reporting level [MRL], the LFSM recovery must be within 50-150% of the true value, whereas, spikes at higher levels must be within 70-130%. An important note is that the fortified samples must be corrected for the PFAS levels in the unfortified samples. The LFSM recovery is calculated by the equation:

$$\%R = \frac{(A - B)}{C} \times 100$$

Where,

A = measured concentration in the fortified sample,  
B = measured concentration in the unfortified sample, and  
C = fortification concentration.

*Step 1.* In the Calculated Columns pane, the equation is entered for the LFSM recovery using the IF function to apply the calculation to LFSM samples only, and the GETSAMPLE function which pulls in the analyte concentration in the unfortified sample only. LFSM samples are designated as QC samples and the spike concentration is set as the actual concentration. The unfortified sample is treated as an unknown. Since the LFSM recovery criteria differ depending on the spiking level, the LFSM sample name needs to be specified, such as "LFSM-MDL" and "LFSM" to distinguish between the LFSM spiked at near-MDL and higher levels [Figure 7].

*Step 2.* In the Flagging Rules pane, custom rules are added to flag instances where the LFSM recovery is outside of the 50-150% criteria for LFSM samples  $\leq 2\times$  MRL and outside of the 70-130% criteria for higher level spikes [Figure 7].



## Step 1. Calculated column

1. Rename Formula name and "Sample Name" if the LFSM samples are  $\leq 2x$  MRL (for example, "LFSM-MDL")

[MQ4] Modify Method

Workflow

Components

Integration

Library Search

Calculated Columns

Flagging Rules

Advanced

Formula Finder

Non-targeted Peaks

Use the calculator to create a new formula.

Formula name: LFSM Recovery

Description: Type a description for the formula

IF([Sample Name]='LFSM' | [Sample Name]='LFSMD'; ([Calculated Concentration]-GETSAMPLE([Calculated Concentration]; 'Unfortified matrix'))/[Actual Concentration] \*100; '' )

## Step 2. Flagging rule

1. Rename Rule name and "QCs" name if the LFSM samples are  $\leq 2x$  MRL (for example, "LFSM-MDL")

[MQ4] Modify Method

Workflow

Components

Integration

Library Search

Calculated Columns

Flagging Rules

Advanced

Formula Finder

Non-targeted Peaks

Rule name: LFSM Recovery

Flag a results column: LFSM Recovery

Flagging criteria: Range

Step 1: Define the values for the flagging criteria

☒ Value for all components ☐ Values per component type

Lower limit: 70

Upper limit: 130

Step 2: Apply the values to the following sample types

☐ Unknowns

☐ Standards

☒ QCs

☐ Blanks

Only if the sample name contains... LFSM

Treat resulting text values as: Error

2. Change limits to 50 (lower limit) and 150 (upper limit) if the LFSM samples are  $\leq 2x$  MRL

Figure 7. Calculated column and custom flagging rule for the determination of the laboratory fortified sample matrix (LFSM) recovery in EPA 533. The calculated column determines the recovery of the LFSM sample, after subtracting the unfortified sample concentration. The flagging rule identifies instances where the LFSM recovery is outside of the acceptable 70-130% range. For LFSM samples with spiking levels  $\leq 2x$  MDL, a separate calculated column and flagging rule is created since the acceptable range is 50-150% for these QC samples.

## EPA 533: Continuing calibration check (CCC) flagging

EPA Method 533 specifies that CCC standards must be analyzed at the beginning and end of each analysis batch and after every tenth field sample [section 10.4]. The purpose of the CCC QC samples is to verify the existing calibration accuracy. The CCC at the beginning of the batch must be at, or below, the MRL for each analyte, and the accuracy must be within 50-150% of the true value. The remaining CCCs in the batch may be between the mid and high calibration levels, and the accuracy must be within 70-130%.

Step 1. Since the quality control criteria are different for the beginning CCC ["CCC-low"] and remaining CCCs, unique flagging rules must be created [Figure 8].

Step 2. In the Flagging Rules pane, custom rules are created to flag instances where the CCC accuracies are outside of the specified criteria [Figure 8]. "Accuracy" is selected as the "Flag a results column", and the "Flagging criteria" is "Range". Depending on the specific CCC type, the lower and upper ranges are entered. Finally, the "QCs" box is checked under the sample type option and the appropriate CCC sample names are entered [either "CCC-Low" or "CCC-Mid" and "CCC-High"].

### 1. Flagging rule for CCC samples run at beginning of batch ["CCC-low"]

The screenshot shows the 'Flagging Rules' configuration window for a rule named 'CCC-Low'. The 'Flag a results column' is set to 'Accuracy' and the 'Flagging criteria' is set to 'Range'. The 'Values per component type' table is highlighted with a red circle, showing lower and upper limits for various components.

Component	Lower limit	Upper limit
PFBA	50	150
PFPeA	50	150
PFHxA	50	150
PFHxA_2	50	150
PFHpA	50	150
PFHpA_2	50	150
PFOA	50	150
PFOA_2	50	150
PFNA	50	150

Step 2: Apply the values to the following sample types. The 'QCs' checkbox is checked, and the sample name 'CCC-Low' is entered in the 'Only if the sample name contains...' field.

### 2. Flagging rule for CCC samples run during the mid- and end of batch ["CCC"]

The screenshot shows the 'Flagging Rules' configuration window for a rule named 'CCC'. The 'Flag a results column' is set to 'Accuracy' and the 'Flagging criteria' is set to 'Range'. The 'Values per component type' table is highlighted with a red circle, showing lower and upper limits for various components.

Component	Lower limit	Upper limit
PFBA	70	130
PFPeA	70	130
PFHxA	70	130
PFHxA_2	70	130
PFHpA	70	130
PFHpA_2	70	130
PFOA	70	130
PFOA_2	70	130
PFNA	70	130

Step 2: Apply the values to the following sample types. The 'QCs' checkbox is checked, and the sample names 'CCC-Mid, CCC-High' are entered in the 'Only if the sample name contains...' field.

Figure 8. The custom flagging rule for the continuing calibration check (CCC) sample accuracy in EPA 533. The flagging rule identifies instances where the CCC accuracy is outside of the acceptable 50-150% range [CCC samples at the batch start] or 70-130% [remaining batch CCC samples].

## Conclusions

- QC calculations for the EPA methods 533 and 1633 can be built directly within SCIEX OS using the Calculated Columns feature in the processing method
- Ability to perform custom calculations within the SCIEX OS software remove the need to validate calculations performed in third-party software
- Custom Flagging rules allow the user to rapidly identify samples that are outside of the EPA criteria
- Time for data processing and review saved through formulating QC calculations and flagging rules directly within the SCIEX OS software

## References

1. Method 533: Determination of per- and polyfluoroalkyl substances in drinking water by isotope dilution anion exchange solid phase extraction and liquid chromatography/tandem mass spectrometry. U.S. Environmental Protection Agency, Office of Water [MS-140], November 2019. [EPA 815-B-19-020](#).
2. Method 1633. Analysis of per- and polyfluoroalkyl substances (PFAS) in aqueous, solid, biosolids, and tissue samples by LC-MS/MS. U.S. Environmental Protection Agency, Office of Water [4304T], Office of Science and Technology, Engineering and Analysis Division, Washington, DC, January 2024. [EPA 821-R-24-001](#).
3. EPA method 533 for PFAS analysis in drinking water at low parts-per-trillion level. SCIEX technical note, **2021**, [RUO-MKT-02-12843-A](#).
4. Analysis of per- and polyfluoroalkyl substances (PFAS) in aqueous, solid, biosolid and tissue samples following EPA method 1633. SCIEX technical note, **2023**, [MKT-29278-A](#).

1. Method 533: Determination of per- and polyfluoroalkyl substances in drinking water by isotope dilution anion exchange solid phase extraction and liquid chromatography/tandem mass spectrometry. U.S.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries [see [www.sciex.com/trademarks](http://www.sciex.com/trademarks)].

© 2024 DH Tech. Dev. Pte. Ltd. MKT-33079-A



### Headquarters

500 Old Connecticut Path | Framingham, MA 01701  
USA  
Phone 508-383-7700  
[sciex.com](https://sciex.com)

### International Sales

For our office locations please all the division headquarters or refer to our website at [sciex.com/offices](https://sciex.com/offices)