

Trace level detection of glyphosate in water and beer samples

Using LC-MS/MS analysis on the QTRAP® 6500+ System

André Schreiber¹, Wen Jin¹, and Paul Winkler²

¹SCIEX, Canada; ²SCIEX, USA

Glyphosate (N-(phosphonomethyl)glycine) is a widely-used, broad-spectrum, systemic herbicide and crop desiccant. Generally, glyphosate is considered as safe and not toxic to humans.¹⁻³ However, glyphosate is a topic with an extraordinary degree of public attention and concern since the International Agency for Research on Cancer (IARC), a branch of the World Health Organization, classified glyphosate as a probable human carcinogen.⁴ Traces of glyphosate have been found in surface water, many foods (such as bread, breakfast cereals, dairy, and beer) and also in human urine and breast milk.⁵⁻⁹

Glyphosate can be analyzed using an enzyme-linked immunosorbent assay (ELISA). Although relatively quick and simple to perform, ELISA tests are limited in selectivity and are susceptible to cross-reactivity, which can lead to false positive or false negative results. When analyzed using LC, glyphosate is derivatized with FMOC to improve its retention, as it is very polar. This derivatization step complicates the analysis and there is a growing need for a method that can detect glyphosate and AMPA in their underivatized forms. Anion exchange, HILIC, porous graphitized carbon and mixed-mode columns were used with LC-MS/MS to determine underivatized polar pesticides with limited success.^{7, 10-12}



Here, an LC method was developed using a mixed-mode column and a low pH mobile phase (pH 2.9). LOQs as low as 100 ng/L in water and 200 ng/L in beer were achieved by utilizing large volume injections (50 μ L) and high sensitivity detection with the SCIEX QTRAP 6500+ System. The method was successfully applied to the analysis of 40 different beers, yielding comparable results to previously reported results where available.⁹

Key features of the method

- Enables the identification and quantification of underivatized glyphosate and its metabolite AMPA in water and beer samples
- High sensitivity achieved using the SCIEX QTRAP 6500+ LC-MS/MS System
- Very simple sample preparation using large volume injections of either water or diluted beer
- High confidence in identification was achieved by monitoring 4 MRM transitions per compound
- Limits of quantification (LOQ) for both compounds of 100 ng/L in water samples and 200 ng/L in beer samples
- Excellent repeatability and linearity was observed




Figure 1. Separation of AMPA and glyphosate. MRM chromatograms of 10 ng/mL of AMPA and glyphosate using a 10 μ L (bottom) and 50 μ L injection volume (top).

Experimental

Sample preparation: Tap water was obtained from the laboratories in the SCIEX office in Concord, Ontario (Canada). Store-bought samples were obtained from the Liquor Control Board of Ontario stores (LCBO). One home-made ale brewed with Toronto tap water was obtained in addition to a commercial barley malted beer. All samples were degassed and diluted 2x with LC grade water.

Chromatography: The ExionLC™ AD System was used to perform the separation, using an Acclaim Trinity Q1 column (100 x 3 mm, 3 µm). Mobile phase A was water with 50 mM ammonium formate/formic acid (pH = 2.9) and mobile phase B was acetonitrile. Injection volume was 50 µL.

Mass spectrometry: Mass analysis was performed on a SCIEX QTRAP 6500+ System equipped with an IonDrive™ Turbo V Ion Source and an electrospray ionization (ESI) probe. Both analytes were detected using negative polarity using the MRM transitions outlined in Table 1. Data acquisition was done using Analyst® Software 1.6.3.

Data processing: Data processing was performed using MultiQuant™ Software 3.0.2.

Table 1. MRM transitions.

Compound	Q1	Q3	DP (V)	CE (V)
Glyphosate	168	63	-30	-26
	168	150	-30	-14
	168	124	-30	-16
	168	81	-30	-20
AMPA	110	63	-15	-26
	110	79	-15	-36
	110	81	-15	-16
	110	80	-15	-24

Results

Large volume injection was used to achieve the desired LOQ of 100 ng/L. Figure 1 shows MRM chromatograms of AMPA and glyphosate using 10 and 50 µL injection volumes. It can be seen that the larger injection volume increases the glyphosate signal by a factor of 5, but also results in peak broadening of the earlier eluting AMPA.




Figure 2. Tap water results. (Top) Glyphosate was spiked into tap water and the signal for the 4 different MRM transitions at the LOQ of 100 ng/mL is shown. Replicate injections had %CV of 3.32% (n=5). (Bottom) Similar data is shown for AMPA. Here, the %CV was 11.4% at the LOQ of 100 ng/mL.

The LOQ was evaluated by repeat analysis of low level standards spiked into tap water (which was tested previously to not contain glyphosate and AMPA). Figures 2a and 2b show the 4 MRM transitions of both compounds at a concentration of 100 ng/L. After 5 injections the coefficient of variation (%CV) was 3.32% for glyphosate and 11.4% for AMPA, respectively.

Linearity for quantification was evaluated over a range from 100 ng/L to 100 µg/L. Linearity was excellent, with coefficients of regression better than 0.999 using linear fit with 1/x weighting (Figure 3). Accuracies were all between 80 and 120% at all concentration levels.

After initial verification, the new LC-MS/MS method was applied to the analysis of glyphosate and AMPA in commercial and homemade beers. Glyphosate was frequently detected. Example chromatograms are shown in Figure 4. AMPA was not detected in any samples.




Figure 3. Linear dynamic range explored. The linearity for glyphosate (left) and AMPA (right) across the concentration range of 100 ng/mL to 100 µg/mL is shown, using a linear fit and a 1/x weighting.




Figure 4. Glyphosate findings in different beers. Glyphosate levels were measured in the various beverage samples: German pilsner (GP), 21.6 µg/L; American light beer (AL), 3.8 µg/L; Irish stout (IS), 16.2 µg/L; Canadian craft India pale ale (IPA), 9.5 µg/L; German weissbier (GW), 0.2 µg/L; and home-made ale (HA), 0.7 µg/L.

High confidence in identification was achieved by the detection of 4 MRM transitions and calculation of quantifier-qualifier ratios. An MRM ratio tolerance of 30% was applied to identification. Due to interferences, some beer samples had failing MRM transition ratios, but in all cases except the German weissbeer, at least 2 transitions were present at the correct ratio for all of the beers tested. Quantitative results are listed in Table 2. These results correlate well with previously reported data from the Environmental Institute (Umweltinstitut München, Germany), which is surprising, considering that samples were purchased in different stores at different times.

The glyphosate concentrations in beer analyzed by LC-MS/MS ranged from 0.22 to 23.78 µg/L. There is no obvious correlation between the concentration of glyphosate and the origin or style of the beer. However, beers brewed with adjuncts such as rice (typical for American light beers) or wheat (German weissbier) tend to have a lower concentration of glyphosate.

These results support the hypothesis that glyphosate originated from the malted barley used for brewing and not from other ingredients, such as water, hops, and yeast.

Table 1. Glyphosate concentrations measured in 40 different beers. Where comparison was possible, results correlated well with published values.

Beer	Glyphosate (µg/L)	Previously reported ⁹
American IPA	2.72	
American Light (1)	3.56	
American Light (2)	1.55	
American Light (3)	1.13	
American Light (4)	0.69	
American Light (5)	0.22	
Canadian Ale (1)	9.99	
Canadian Ale (2)	7.54	
Canadian Bock (1)	6.52	
Canadian Bock (2)	4.21	
Canadian IPA (1)	15.71	
Canadian IPA (2)	14.93	
Canadian IPA (3)	13.97	
Canadian IPA (4)	10.63	
Canadian IPA (5)	10.48	
Canadian IPA (6)	9.48	
Canadian IPA (7)	7.10	
Canadian IPA (8)	6.97	
Canadian IPA (9)	6.83	
Canadian IPA (10)	5.61	
Canadian IPA (11)	5.14	
Canadian IPA (12)	5.09	
Canadian IPA (13)	3.06	
Canadian IPA (14)	2.31	
Canadian Stout	9.84	
Czech Pilsner (1)	13.96	
Czech Pilsner (2)	6.18	
Czech Pilsner (3)	6.15	
Czech Pilsner (4)	3.95	
German Pilsner (1)		29.74
German Pilsner (2)	23.78	23.04
German Pilsner (3)	7.21	20.73

Table 1 continued. Glyphosate concentrations measured in 40 different beers. Where comparison was possible, results correlated well with published values.

Beer	Glyphosate (µg/L)	Previously reported ⁹
German Pilsner (4)	6.77	12.01
German Pilsner (5)	4.98	
German Pilsner (6)	3.41	0.50
German Pilsner (7)	2.78	
German Pilsner (8)	0.87	2.99
German Pilsner (9)	0.76	0.55
German Pilsner (9)	0.27	
German Pilsner (9)		5.78
German Pilsner (9)		3.86
German Pilsner (9)		3.35
German Weissbier (1)	0.75	
German Weissbier (2)		2.92
German Weissbier (3)		0.66
German Weissbier (4)		0.49
Home-made Ale	18.65	
Irish Stout	13.96	

Summary

Here, the analytical results for underivatized glyphosate and its metabolite, AMPA, in water and beer samples using LC-MS/MS was investigated. The method, using a SCIEX QTRAP 6500+ System and direct injection of 50 µL of liquid, provided excellent sensitivity, repeatability, and linearity. Water samples were injected directly resulting in an LOQ of 100 ng/L and beer samples were injected after degassing and 1/1 dilution with water resulting in an LOQ of 200 ng/L. High confidence in identification was achieved by monitoring 4 MRM transitions per compound. 40 beers samples were analyzed with glyphosate findings between 0.22 to 23.78 µg/L. Results correlate well with previous reported data.⁹

References

- European Food Safety Authority (EFSA): Conclusion on the peer review of the pesticide risk assessment of the active substance glyphosate. [EFSA Nov 2015](#).
- Environmental Protection Agency (EPA): Glyphosate. [EPA Nov 2015](#).
- International Agency for Research on Cancer (IARC): Evaluation of five organophosphate insecticides and herbicides. [IARC Monographs 112 \(2015\)](#).
- W. A. Battaglin *et al.* (2002) Glyphosate, other herbicides, and transformation products in Midwestern streams, [Journal of the American Water Resources Association \(JAWRA\) 4, 323-332](#).
- M. Krüger *et al.* (2014) Detection of Glyphosate Residues in Animals and Humans. [J Environ Anal Toxicol 4, 1-5](#).
- N. Chamkasem *et al.* (2015) Direct Determination of Glyphosate, Glufosinate, and AMPA in milk by Liquid chromatography/tandem mass spectrometry. [Journal of Regulatory Science 02, 20-26](#).
- Glyphosate testing full report: findings in American mothers breast milk, urine and water. [Moms Across America 2014](#).
- <http://www.umweltinstitut.org/aktuelle-meldungen/meldungen/umweltinstitut-findet-glyphosat-in-deutschem-bier.html>
- J. Dahlmann *et al.* (2006) Direct injection detection using LC/MS/MS: Analysis of dissociated organo-phosphorus pesticides. [CLB Chemie in Labor und Biotechnik 57, 356-359](#).
- A. Vass *et al.* (2016) Study of different HILIC, mixed-mode, and other aqueous normal-phase approaches for the liquid chromatography mass spectrometry-based determination of challenging polar pesticides. [Anal. Bioanal. Chem. 408, 4857-4869](#).
- M. Anastassiades *et al.* (2016) Quick Method for the Analysis of numerous Highly Polar Pesticides in Foods of Plant Origin via LC-MS/MS involving Simultaneous Extraction with Methanol. [EURL-SRM v9](#).

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries. © 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-4543-A. AB SCIEX™ is being used under license.



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division headquarters or refer to our website at sciex.com/offices