

# Achieve sensitive quantitation of signature peptides using microflow LC and accurate mass spectrometry

*Simplifying quantitation methods using the ZenoTOF 7600 system, powered by SCIEX OS software*

Xin Zhang, Eshani Nandita, Lei Xiong, Zoe Zhang and Rahul Baghla  
SCIEX, USA

This technical note describes a sensitive quantitation method for signature peptides in rat plasma using a microflow LC coupled to the ZenoTOF 7600 system. Streamlined method development was performed by collecting fragment ion information over a wide MS/MS range for each peptide. Multiple highly abundant fragment ions were summed, resulting in a 2- to 5-fold improvement in the lower limit of quantitation (LLOQ, Figure 1). The LLOQs observed ranged from 2.64 pg/mL to 16.8 pg/mL. The broad linear dynamic range (LDR) spanned more than orders of magnitude.

Quantitation of peptide and protein therapeutics in biological matrices is important for the development of biotherapeutics. Serving as an orthogonal technology to the traditional ligand binding assays (LBAs), LC-MS has been routinely adopted for quantitative measurement of protein levels in bioanalytical laboratories. The triple quadrupole platform offers great sensitivity and quantitative performance and has been a key driver for most bioanalytical methods. However, accurate mass spectrometry has increasingly been adopted for quantitative bioanalysis.<sup>1,2</sup> With the inherent advantage of greater selectivity, improved mass resolution and the flexibility of TOF MS/MS data analysis, the ZenoTOF 7600 system provides excellent quantitative performance in multiple dimensions.<sup>3</sup> Accurate mass spectrometry platforms, such as traditional time-of-flight (TOF) systems, often lack sensitivity due to a loss of ion transmission between TOF pulses. The Zeno trap controls the ion beam from the collision cell to facilitate greater ion transmission to the TOF accelerator. The Zeno trap boosts the duty cycle to  $\geq 90\%$ , which enhances overall MS/MS sampling efficiency. In this technical note, 3 signature peptides were quantified using a microflow LC method paired with the ZenoTOF 7600 system.

## Key features of peptide quantitation using a microflow LC with the ZenoTOF 7600 system

- **Low-level quantitation:** Achieve sensitive quantitation (LLOQs between 2.64 pg/mL and 16.8 pg/mL) of signature peptides in complex matrices using the Zeno MRM<sup>HR</sup> method
- **Simplify quantitation methods:** Streamline quantitation methods by collecting data over a wide MS/MS range for flexible post-acquisition data processing
- **Improve quantitative sensitivity:** Reach enhanced sensitivity by summing multiple highly abundant fragment ions using TOF MS/MS data and the Zeno trap
- **Quantitative performance:** Ensure accurate and highly reproducible (%CV <15%) quantitative methods with a LDR spanning  $\geq 4$  orders of magnitude using the ZenoTOF 7600 system
- **Streamline data management:** Easily acquire and process data on a single platform using SCIEX OS software




Figure 1. Extracted ion chromatograms (XICs) at the LLOQ using a single fragment ion (left) or multiple summed fragment ions (right). A 2- to 5-fold improvement in LLOQ was observed when multiple highly abundant fragment ions were summed for 3 peptides.

## Methods

**Sample preparation:** Plasma proteins were precipitated with cold methanol. After centrifugation, the supernatant was discarded. The pellet was solubilized in 200mM ammonium bicarbonate in 10:90 (v/v), methanol/water. Digestion was performed using trypsin. The solution was heated for 1 hour at 60°C and then acidified by the addition of formic acid.<sup>4</sup> The digested plasma was diluted by 200x using a solution containing 5% acetonitrile, 1% formic acid and 94% water by volume. Synthesized peptides (Table 1) were spiked into the digested plasma solution and the resulting solution was then serially diluted in matrix. The final injection volume was 10  $\mu$ L.

**Table 1. List of peptide targets.**

Peptide sequence	Description
FNWYVDGVEVHNAK	Conserved sequence in human immunoglobulin G (IgG)
AGLIVAEGVTK*	Synthetic peptide with C terminal K heavy isotope labeled (C <sup>13</sup> N <sup>15</sup> )
LGLDFDSFR*	Synthetic peptide with C terminal R heavy isotope labeled (C <sup>13</sup> N <sup>15</sup> )

**Chromatography:** The separation was performed using a Waters M Class system. A YMC-Triart C16 column (0.3 x 50 mm, 3  $\mu$ m, 120 Å) was used for separation with a YMC-Triart C16 1/16 capillary guard (0.3 x 5 mm, 3  $\mu$ m, 120 Å). The flow rate was 40  $\mu$ L/min. The column oven temperature was set to 40°C. Mobile phase A was 0.1% formic acid in water and mobile phase B was 0.1% formic acid in acetonitrile. The gradient conditions used are summarized in Table 2. A volume of 10  $\mu$ L was injected for analysis. All samples were analyzed in triplicate.

**Table 2. Analyte separation conditions.**

Time (min)	Mobile phase A (%)	Mobile phase B (%)
0.0	95	5
1.0	95	5
2.5	60	40
3.0	60	40
3.5	5	95
8.0	5	95
8.1	95	5
10.0	95	5

**Mass spectrometry:** Data were acquired in positive mode using Zeno MRM<sup>HR</sup> on a ZenoTOF 7600 system. The source was operated in positive ion mode. Collision energy (CE), source and MS parameters were optimized for all the signature peptides. The source and MS parameters and the Zeno trap settings are summarized in Table 3. The MRM<sup>HR</sup> parameters and fragments used to quantify each of the signature peptides are summarized in Table 4. Unit Q1 resolution was used for analysis.

**Table 3. Source and MS conditions.**

Parameter	Value	Parameter	Value
Curtain gas	25 psi	Source temperature	100°C
Ion source gas 1	20 psi	Ion source gas 2	20 psi
CAD gas	7	Ion spray voltage	5500 V
MS accumulation time	80 ms	MS/MS accumulation time	10 ms
TOF MS start mass	400	TOF MS stop mass	800
TOF MS/MS start mass	100	TOF MS/MS stop mass	1200
Zeno threshold	20,000 cps		

**Data processing:** Data were processed using the Analytics module in SCIEX OS software, version 3.0 with the MQ4 integration algorithm. A 1/x<sup>2</sup> weighting was used for quantitation.

**Table 4. MRM<sup>HR</sup> parameters and fragments used for quantitation.**

Peptide	Q1 mass (m/z)	Fragment mass (m/z)	DP (V)	CE (V)
AGLIVAEGVTK* (y6)	533.32	612.344	80	32
AGLIVAEGVTK* (y7)	533.32	711.388	80	32
AGLIVAEGVTK* (b3)	533.32	242.125	80	32
LGLDFDSFR* (y5)	540.27	681.322	80	34
LGLDFDSFR* (y6)	540.27	796.351	80	34
LGLDFDSFR* (y8)	540.27	966.457	80	34
FNWYVDGVEVHNAK (y4)	560.27	469.252	80	30
FNWYVDGVEVHNAK (y6)	560.27	697.363	80	30

## Summation of multiple fragment ions

FNWYVDGVEVHNAK, AGLIVAEGVTK\* and LGLDFDSFR\* were quantified using a microflow LC method on the ZenoTOF 7600 system. In this method, the quantitation of signature peptides was performed using Zeno MRM<sup>HR</sup>.

The accessibility of TOF MS/MS data can be advantageous as post-acquisition data decisions can be made on which measured fragments can be utilized for MRM<sup>HR</sup>. For MRM<sup>HR</sup>, quantitation can be performed using a single fragment ion or by summing multiple dominant fragment ions. When multiple highly abundant fragment ions are generated from the target peptide, summed XICs can further enhance assay sensitivity.

For quantitation using multiple fragments, the most abundant ions were summed for optimal assay sensitivity. Fragment ions  $y_5$ ,  $y_6$  and  $y_8$  were summed for the peptide LGLDFDSFR\*, while fragment ions  $y_4$  and  $y_6$  were summed for the peptide FNWYVDGVEVHNAK. For peptide AGLIVAEGVTK\*, fragment ions  $y_6$ ,  $y_7$  and  $b_3$  were summed for quantitation.


Figure 2 shows the application of SCIEX OS software for the summation of fragment ions. The MS/MS spectrum of  $m/z$  533.3 for peptide AGLIVAEGVTK\* was generated using the Explorer

module in SCIEX OS software. Highly abundant fragment ions including  $b_3$ ,  $y_6$  and  $y_7$  were selected for summation. The product ion information was input into the Analytics module in SCIEX OS software and the feature to sum multiple ions was applied to generate a total XIC for quantitation of the peptide AGLIVAEGVTK\*.


## Quantitative performance for the analysis of signature peptides

The LLOQ was determined based on the requirements that the %CV of the average concentration must be below 20% and the accuracy must be between 80% and 120%. For concentrations above the LLOQ, the %CV of the mean calculated concentration was required to be below 15% and the accuracy was required to be between 85% and 115%.

Summed fragment ions provided LLOQs of 2.64 pg/mL, 10.7 pg/mL and 16.8 pg/mL for peptides AGLIVAEGVTK\*, LGLDFDSFR\* and FNWYVDGVEVHNAK, respectively (Figure 3).




**Figure 2. Summation of multiple fragment ions in SCIEX OS software.** Peptide AGLIVAEGVTK\* was selected as a representative example. Using the Explorer module in SCIEX OS software, the MS/MS spectrum of  $m/z$  533.3 was generated (A). Fragment ions  $b_3$ ,  $y_6$  and  $y_7$  were selected and input into the Analytics module in SCIEX OS software (B). The software enables users to easily select and sum multiple ions to generate a total XIC.



**Figure 3.** XICs of the matrix blank and at the LLOQ using the summation of multiple fragment ions. The LLOQs observed using multiple summed fragment ions for AGLIVAEGVTK ( $y_4$  and  $y_6$ ), LGLDFDSFR ( $y_6$ ,  $y_7$  and  $b_3$ ) and FNWYVDGVEVHNAK ( $y_5$ ,  $y_6$  and  $y_8$ ) were 2.64 pg/mL, 10.7 pg/mL and 16.8 pg/mL, respectively. Minor interference (<20% of the LLOQ) was observed in the matrix blank of the AGLIVAEGVTK\* peptide.

Minor interference (<20% of the LLOQ) was observed in the matrix blank of the AGLIVAEGVTK\* peptide. Figure 1 shows that a 2- to 5-fold improvement in the LLOQ was observed when using the summation of multiple highly abundant fragment ions compared to using the single fragment ion for quantitation.

All calibration points were measured in triplicate. An LDR spanning  $\geq 4$  orders of magnitude was achieved for all peptides studied with a coefficient of determination ( $r^2$ ) of  $>0.99$  (Figure 4). The calculated concentrations of each calibration point were within  $\pm 13\%$  of the nominal value and the overall %CV was  $<14\%$  (Table 5). Accuracy was within  $\pm 6\%$  of the nominal concentration at the LLOQ.



**Figure 4.** Calibration curves for the quantitation of signature peptides using multiple summed fragment ions. An LDR spanning  $\geq 4$  orders of magnitude was achieved with an  $r^2 > 0.99$ .

This assay demonstrates the ability of a microflow LC method to produce greater sensitivity levels for peptide quantitation compared to a previous study using a high-flow LC solution.<sup>4</sup> The LLOQs were improved 2.5- to 5-fold with the application of a microflow LC providing enhanced sampling efficiency for peptide analysis.

**Table 5. Calculated concentration, precision and accuracy for quantitation using summed fragment ions.**

AGLIVAEGVTK*			LGLDFDSFR*			FNWYVDGVEVHNAK		
Concentration (pg/mL)	Accuracy (%)	CV (%)	Concentration (pg/mL)	Accuracy (%)	CV (%)	Concentration (pg/mL)	Accuracy (%)	CV (%)
2.64	94.8	8.48	2.67	N/A	N/A	4.19	N/A	N/A
5.29	109	7.37	5.35	N/A	N/A	8.39	N/A	N/A
10.6	98.9	4.89	10.7	101	2.45	16.8	99.8	13.6
21.1	106	0.64	21.4	101	6.45	33.6	97.5	7.85
52.9	103	10.5	53.5	92.3	10.4	83.9	104	6.37
106	107	5.57	107	97.6	6.49	168	101	4.80
211	107	3.13	214	100	9.04	335	107	1.69
529	99.0	6.27	534	100	6.72	839	105	5.30
1057	98.7	2.85	1069	104	1.44	1678	102	2.37
2114	98.2	0.71	2138	101	2.64	3356	96.6	0.49
5286	96.7	1.44	5346	86.9	1.12	8389	98.6	3.33
10572	97.4	3.27	10692	89.4	1.36	16778	98.2	3.53
21145	96.0	3.44	21384	95.2	0.99	33556	100	1.32
52862	88.3	4.14	53459	104	0.50	83892	97.5	3.81
105725	N/A	N/A	106918	113	1.90	167784	91.7	1.96
211451	N/A	N/A	213836	113	2.92	335567	N/A	N/A

## Conclusions

- Low levels of quantitation (LLOQs between 2.64 pg/mL and 16.8 pg/mL) for signature peptides were achieved using a Zeno MRM<sup>HR</sup> method on the ZenoTOF 7600 system
- Easy method development for peptide quantitation allowed the collection of fragment ion information over a wide MS/MS range for more efficient post-acquisition data processing
- Summation of multiple fragment ions provided a 2- to 5-fold enhancement in LLOQ with the availability of TOF MS/MS data and analysis using the Zeno trap
- Streamlined data acquisition and processing were performed using SCIEX OS software, in which users can efficiently integrate and evaluate quantitative data using features such as the summation of multiple fragment ions
- Microflow LC provided greater sensitivity for peptide quantitation compared to a previously reported high-flow method.<sup>4</sup> A 2.5- to 5-fold improvement in LLOQ was observed given the greater gas phase ion production that microflow LCs afford.

## References

1. Mike-Qingtao Huang, Zhongping (John) Lin, Naidong Weng (2013). [Applications of high-resolution MS in bioanalysis. \*Bioanalysis\* 5\(10\):1269-1276.](#)
2. Yuan-Qing Xia, Jim Lau, Timothy Olah, Mohammed Jemal (2011). Targeted quantitative bioanalysis in plasma using liquid chromatography/high-resolution accurate mass spectrometry: an evaluation of global selectivity as a function of mass resolving power and extraction window, with comparison of centroid and profile modes. [Rapid Communications in Mass Spectrometry 25\(19\):2863-2878.](#)
3. Enhanced sensitivity for peptide quantitation in a complex matrix using high-resolution LC-MS/MS. [SCIEX technical note, RUO-MKT-02-13324-A.](#)
4. Sensitive signature peptide quantification in a complex matrix using accurate mass spectrometry. [SCIEX technical note, RUO-MKT-02-14193-A.](#)

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see [www.sciex.com/trademarks](http://www.sciex.com/trademarks)).

© 2023 DH Tech. Dev. Pte. Ltd. MKT-27700-A



**Headquarters**  
500 Old Connecticut Path | Framingham, MA 01701 USA  
Phone 508-383-7700  
[sciex.com](http://sciex.com)

**International Sales**  
For our office locations please call the division  
headquarters or refer to our website at  
[sciex.com/offices](http://sciex.com/offices)