

Characterization and identification of therapeutic peptide catabolites

Finding, identifying and confirming peptide catabolites using the ZenoTOF 7600 system

Rahul Baghla, Eshani Nandita and Zoe Zhang
SCIEX, USA

This technical note demonstrates the rapid characterization and identification of therapeutic peptide catabolites from rat liver s9 fractions using the ZenoTOF 7600 system. Collision-induced dissociation (CID) using the Zeno trap generated fragment ions for the confident and sensitive identification of peptide catabolites (Figure 1). Catabolite identification, substrate disappearance and metabolite formation were monitored at different time points using Molecule Profiler software. A streamlined workflow was developed to efficiently characterize and identify therapeutic peptide catabolites using a 10-minute run.

There has been a significant increase in peptide therapeutics approved by the Food and Drug Administration (FDA) recently.¹ Peptide therapeutics have many advantages over small

molecule drugs, including higher affinity and lower toxicity. Since peptide-like drugs are generally biotransformed into smaller peptides and amino acids, catabolite identification assays are crucial during drug development. Such assays are important to confidently identify structures containing modified amino acids and metabolic soft spots that facilitate their modification to more stable sequences. In vitro models can provide valuable insight into peptide drug catabolism. Efficient identification and verification of catabolites are crucial at multiple drug discovery and development stages.¹ Therapeutic peptide catabolism studies present various challenges, including the possibility of unspecific hydrolysis with multiple charge states for both substrate and potential catabolites. Given a large panel of potential catabolites, more difficulties with identification are presented than with small molecule metabolites.¹ A streamlined workflow for peptide catabolite characterization and identification is demonstrated here using the ZenoTOF 7600 system and Molecule Profiler software.




Figure 1. Comparison of Zeno CID DDA with CID DDA using representative MS/MS spectra for the PyrHWSYLLR catabolite. Zeno CID DDA provided a 20-fold gain in overall MS/MS intensity over CID DDA. Confident identification of lower abundant catabolites can be performed using the Zeno trap. Catabolite and fragment identification were performed with <10 ppm mass error.

Key features of catabolite identification using the ZenoTOF 7600 system

- Comprehensive characterization and identification:** Achieve comprehensive characterization and identification of therapeutic peptide catabolites from rat liver s9 incubations of leuprolide using the ZenoTOF 7600 system
- Detection of low-level catabolites:** Identify all critical catabolites present in peptide catabolism studies with enhanced MS/MS sensitivity provided by the Zeno trap
- Site-specific identification:** Utilize a dedicated peptide processing workflow to easily identify the site of catabolism using Molecule Profiler software
- Integrated correlation function:** Perform catabolism comparison across multiple samples for time course studies or inter-species comparison, using both MS/MS and analog data
- Streamlined data processing:** Develop confident structure-metabolic stability analyses for drug products with a quick, easy-to-use methodology from acquisition to analysis using SCIEX OS software

Methods

Sample preparation: Leuprolide was incubated at 37°C in rat liver s9 fractions at a starting concentration of 10 µM. Incubations were performed with NADPH. Samples were removed from incubation and quenched with acetonitrile at 0-, 30-, 60-, 90-, 150- and 270-minute intervals. Each sample was centrifuged at 9400 rcf for 10 min at room temperature and was diluted 1:1 (v/v) with water for analysis.

Chromatography: LC separation was performed on an ExionLC system with a [Phenomenex Kinetex Polar C18 column \(2.1 x 100 mm, 2.6 µm, 100 Å\)](#) at a flow rate of 0.5 mL/min. The operating column temperature was 40°C. Mobile phase A was 0.1% (v/v) formic acid in water and mobile phase B was 0.1% (v/v) formic acid in acetonitrile. An injection volume of 2 µL was used for analysis.

The chromatographic gradient conditions used are summarized in Table 1.

Table 1. Chromatographic gradient.

Time (min)	Mobile phase A (%)	Mobile phase B (%)
0.0	95	5
1.0	95	5
5.0	5	50
7.0	5	95
9.0	5	95
9.1	95	5
10.0	95	5

Mass spectrometry: The samples were analyzed using data-dependent acquisition (DDA) with Zeno CID on the ZenoTOF 7600 system. The method conditions used are summarized in Table 2. The source and gas conditions used are summarized in Table 3.

Table 2. Zeno CID DDA parameters.

Parameter	Setting
Method duration	10 min
TOF MS start–stop mass	100–1500 Da
Maximum candidate ions	8
Accumulation time (TOF MS)	0.1 s
TOF MS/MS start–stop mass	50–1500 Da
Accumulation time (TOF MS/MS)	0.06 s
Collision energy (CID)	35 V
Collision energy spread (CID)	15 V

Table 3. Source and gas conditions.

Parameter	Setting
Curtain gas	40 psi
Ion source gas 1	50 psi
Ion source gas 2	60 psi
CAD gas	7
Ion spray voltage	5500 V
Source temperature	500°C

Data processing: SCIEX OS software, version 3.0 was used for data acquisition. Molecule Profiler software was used to predict catabolism sites using Zeno CID DDA data.

Identification of peptide catabolites

The Zeno trap in the ZenoTOF 7600 system significantly boosts the sensitivity of MS/MS fragments. Therefore, Zeno DDA data enables good MS/MS coverage for TOF MS peaks of interest in CID acquisitions. For data processing, Molecule Profiler software automatically predicts the catabolites based on MS1 data and performs catabolite identification using the peptide precursor and catabolite-specific fragment ions.

Leuprolide (pyrHWSYLLRP-[C₂H₆N]) was used as an experimental model to demonstrate catabolite identification using Molecule Profiler software. Catabolite and fragment identification were performed with <10 ppm mass error. Fragment ions with a signal-to-noise ratio (S/N) ≥5 were

selected for peptide sequence assignment to minimize false positives.

Zeno CID DDA generated abundant diagnostic fragments to identify sites of catabolism from leuprolide incubations. Molecule Profiler software predicted potential catabolites from the parent peptide sequence entered. With software features including candidate ranking and a/y/b ion labeling, Molecule Profiler software facilitates automatic catabolite confirmation and identification.

The 270-minute incubation sample used in this technical note yielded 13 catabolites, including hydrolytic catabolites (Figure 2). The raw data was presented in extracted ion chromatograms (XICs), MS and MS/MS spectra.²




Figure 2. Results view in Molecule Profiler software. The software shows the number of potential catabolites with detailed information about the structure and scoring, options to view the XICs of selected catabolites and TOF MS spectra. Molecule Profiler software also displays the mass defect, isotope pattern, CID and EAD spectra of the selected catabolite.

The software identified PyrHWSYLLR (hydrolysis at R8) as a positive catabolite. The interpretation pane for pyrHWSYLLR in Figure 3 shows the MS/MS spectra for the catabolite and the fragment pane with proposed formulas. In addition to MS/MS matching, Molecule Profiler software uses the overall intensity

of the fragments and the peak area assignments to confidently determine the catabolites. Complete sequence coverage for 8 amino acids was achieved with only 1 sequence candidate for the selected catabolite. Here, 86.5% of the fragments from pyrHWSYLLR were assigned based on the




Figure 3. Interpretation pane for catabolite PyrHWSYLLR. The MS/MS spectra for the catabolite PyrHWSYLLR are shown with common product ions. The parent peptide is highlighted in yellow. Complete sequence coverage was achieved for all 8 amino acids and only 1 sequence candidate was predicted, demonstrating comprehensive MS/MS coverage for confident catabolite identification.

MS/MS spectrum peak area using a S/N ≥5. Another hydrolysis catabolite, YLLRP-[C₂H₆N], underwent hydrolysis at R4 and P5, had a retention time of 3.39 minutes and a charge state of +1. YLLRP-[C₂H₆N] also showed complete sequence

coverage for 5 amino acids with 87.5% of fragments assigned based on MS/MS spectrum peak area with a S/N ≥5. Confident assignment was performed with only 1 possible candidate based on Zeno CID DDA MS/MS fragments (Figure 4).




Figure 4. Interpretation pane for the YLLRP catabolite. The MS/MS spectra of a hydrolysis catabolite at m/z 688.5 provided complete sequence coverage, capturing all 5 amino acids in the sequence. One sequence candidate was predicted based on the Zeno CID DDA results, demonstrating high confidence and comprehensive catabolite identification.

Correlation analysis

The correlation workspace of Molecule Profiler software was used to confirm the presence of the significantly identified catabolites across several time points. The correlation workspace has an interactive graph window in which results can be displayed as an x-y plot, bar graph or table. In addition to the correlation workspace, overlaid XICs, TOF MS and TOF MS/MS spectra can be displayed for the chosen catabolites or metabolites. Figure 5 shows the correlation plot of leuprolide and 2 major catabolites discussed in this study.

insights into the initial steps of peptide degradation and elimination. This study was therefore designed to evaluate the ability to measure the disappearance of the substrate peptide and to identify catabolites in the liver model (liver s9 fraction). The results generated by combining Molecule Profiler software with the ZenoTOF 7600 system support the streamlined workflow to efficiently characterize and identify therapeutic peptide catabolites.¹




Figure 5. Correlation analysis using Molecule Profiler software.

Correlation analysis showed increased catabolite concentration with time and decreased parent concentration in incubated samples.

The metabolism of peptide drugs is not solely dependent on the activity of cytochrome P450 (CYP) enzymes found in the liver. As a result, using liver-derived models alone is insufficient to accurately predict the peptide therapeutic catabolites. Other tissues and blood can also contribute to the metabolic pathways of the peptide drugs. However, despite their inadequacy, hepatic models, blood/plasma and kidneys are frequently used and are irreplaceable systems to study peptide catabolism. Together, these models provide valuable

Conclusions

- Characterization and identification of therapeutic peptide catabolites from rat liver s9 incubations of leuprolide were demonstrated on the ZenoTOF 7600 system
- A sensitive workflow enables the detection of low-level catabolites and can be easily adapted for in vivo catabolism studies with the enhanced sensitivity provided by the Zeno trap
- Thirteen catabolites from leuprolide were characterized and identified using Molecule Profiler software
- Diagnostic fragment ions were used to identify the site of catabolism for several hydrolysis catabolites using Zeno CID DDA
- A streamlined data processing method was utilized for ease of data reduction and development of confident structure-metabolic stability relationships for drugs

References

1. Jyrkas, J. et al. (2021). Hepatic in vitro metabolism of peptides; Comparison of human liverS9, hepatocytes and Upcyte hepatocytes with cyclosporine A, leuprorelin, desmopressin and cetrorelix as model compounds. [Journal of Pharmaceutical and Biomedical Analysis. 196\(2021\): 113921.](#)
2. Therapeutic Peptide Catabolite Identification using MetabolitePilot™ 2.0 Software and the TripleTOF® 6600 System. SCIEX technical note, [RUO-MKT-02-5477](#)

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2023 DH Tech. Dev. Pte. Ltd. MKT-27886-A



Headquarters

500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales

For our office locations please call the division headquarters or refer to our website at sciex.com/offices