



Sensitive quantitation of an antibody-drug conjugate [ADC] using LC-MRM^{HR} in human serum

Lakshmanan Deenadayalan¹, Sashank Pillai¹, Rahul Baghla² and Eshani Galermo²

¹SCIEX, India and ²SCIEX, USA

This technical note demonstrates a sensitive method for quantitation of trastuzumab deruxtecan (TDx) in human serum using high-resolution accurate mass spectrometry. A lower limit of quantitation (LLOQ) of 0.05 µg/mL (0.3 ng on column) was achieved in extracted human serum samples (Figure 1).

TDx was recently approved by FDA [2019] for treating breast, gastric and gastroesophageal where it functions by binding to HER2 of the malignant cells causing targeted DNA damage in cancer cells. It is an anti-HER2 ADC that comprises a HER2 monoclonal AB conjugate via internal cysteine residues.¹ TDx has been known to express antitumor activity including in low HER2-expressing cancers.^{2,3} As a result, it is critical to accurately measure levels of ADCs such as TDx in biological matrices for toxicokinetic and pharmacokinetic profiles to meet safety and efficacy requirements.

Key benefits for analysis of trastuzumab deruxtecan using the ZenoTOF 7600 system

- **Sensitive quantitation of ADC:** Achieve 0.05 µg/mL LLOQ for quantitation of TDx in human serum
- **Low serum consumption:** Low-level quantitation was achieved using 50 µL human serum with increased MS/MS sampling efficiency using Zeno MRM^{HR}
- **Effortlessly meet critical quantitative performance criteria:** Achieve accurate quantitative performance with %CV <13% at all concentration levels across a linear dynamic range (LDR) of 3.7 orders of magnitude
- **Streamlined data management:** SCIEX OS software, a 21 CFR Part 11-compliant platform, simplifies data acquisition and processing




Figure 1: Representative extracted ion chromatograms (XICs) of matrix blank, 0.05 µg/mL (LLOQ) and 5 µg/mL for TDx in extracted human serum.

Introduction

ADCs have emerged as the most promising type of targeted cancer treatment with high target specificity.⁴ The structure of ADC comprises a monoclonal antibody (mAb) attached to a cytotoxic drug through a chemical linker. The mAb specifically targets the tumor antigen while the payload can be released to the site for its cytotoxic effects.⁵ Due to the specificity of ADCs towards cancer cells, such structures present high efficacy with minimal side effects. TDx is an anti-HER2 ADC with a specific target for tumor cells. Thus, it is very efficient for treating patients with breast and gastric cancers.

Due to its high potency in treating cancers, highly sensitive assays are necessary to ensure precise and accurate detection and quantitation when assessing pharmacokinetic and pharmacodynamic effects.

Methods

Standard preparation: 1 mg of TDx stock was procured from Medchem Express and dissolved in water. The dilutions were made in PBS for further processing.

Sample preparation: Dilutions were performed in human serum to prepare spiked sample concentrations ranging from 0.05 µg/mL to 250 µg/mL. Sample preparation was performed using 50 µL of spiked human serum, 200 µL of protein A beads and 200 µL of PBS. Protein A beads were washed three times with PBS before use. Samples were gently shaken for 45 minutes at room temperature, followed by two PBS washes. The beads were resuspended in 150 µL of digestion buffer containing 150mM ammonium bicarbonate and 1mM calcium chloride before denaturation at 95°C for 10 minutes. After allowing samples to cool to room temperature, 10 µg of trypsin was added to each sample, followed by on-bead digestion for 2 hours at 50°C. Digestion was stopped by adding 3 µL of formic acid. Samples were separated from the beads and transferred to vials for LC-MS/MS analysis.

Chromatography: Analytical separation was performed on the ExionLC AE system using a [Phenomenex Aries peptide XB C18 \(2.1 × 100 mm, 2.6 µm\)](#) column at a 0.5 mL/min flow rate. Mobile phase A was 0.1% [v/v] formic acid in water and mobile phase B was 0.1% [v/v] formic acid in acetonitrile. The column temperature was set to 40°C. The gradient conditions used are

summarized in Table 1. A 20 µL sample was used for LC-MS/MS analysis.

Table 1: LC gradient conditions for TDx.

Time [min]	Mobile phase A [%]	Mobile phase B [%]
0.0	95	5
5.0	30	30
5.1	5	95
5.6	5	95
5.7	95	5
8.0	95	5

Table 2: Source, gas and ZenoTOF 7600 system conditions.

Parameter	MS	MS/MS
Scan mode	TOF MS	MRM ^{HR}
Polarity		Positive
Gas 1		60 psi
Gas 2		60 psi
Curtain gas		35 psi
Source temperature		650° C
Ion spray voltage		5500 V
CAD gas		12
Declustering potential	80 V	See Table 3
Start mass	m/z 100	m/z 100
Stop mass	m/z 1500	m/z 1200
Q1 resolution	NA	Unit
Accumulation time	0.1 s	0.1 s
Collision energy	10 V	See Table 3
CE spread	0 V	0 V
Zeno trap	NA	ON
ZOD threshold [CID]	NA	20,000 cps
Time bins to sum	6	6

Table 3: Zeno MRM^{HR} parameters used for quantitation.

ID	Precursor ion [m/z]	Fragment ion [m/z]	CE [V]	DP [V]
IYPTNGYTR	542.70	404.701	20	80
IYPTNGYTR	542.70	808.394	20	80

Data processing: Analysis was performed using SCIEX OS software, version 3.4.0. Peaks were integrated using the MQ4 algorithm, and a weighting of $1/x^2$ was used for TDx quantitation. An XIC peak width of 0.05 Da was used for quantitation.

Quantitative performance on the ZenoTOF 7600 system

The Zeno MRM^{HR} technique provides superior sensitivity and selectivity for measuring TDx in human serum. For quantitative workflows requiring high sensitivity, the Zeno trap on the ZenoTOF 7600 system enhances the effectiveness of total MS/MS sampling.

The Zeno MRM^{HR} method achieved excellent linearity, accuracy and precision from 0.05 µg/mL to 250 µg/mL for TDx. The on-bead digestion immunoprecipitation workflow resulted in clean samples that enabled a minimal human serum volume of 50 µL for the sensitive quantitation of TDx.

Linearity was achieved across concentrations ranging from 0.05 µg/mL to 250 µg/mL [Figure 2], achieving an LDR of 3.7 orders of magnitude.

Calibration for IYPTNGYTR_01: $y = 2004.80695x \dots = 0.99726$, $r^2 = 0.99453$ (weighting: $1/x^2$)




Figure 2: Calibration curve for quantitation of TDx with weighting factor $1/x^2$.

Analytical performance was evaluated based on the requirement that the accuracy of the calculated mean should be between 80% and 120% at the LLOQ and between 85% and 115% at higher concentrations. The %CV of the calculated mean of the concentration should be below 20% at the LLOQ and below 15% at all higher concentrations.⁴

This assay's accuracy was within $\pm 15\%$ of the nominal concentration and %CV was <13 for the quantitation of TDx in human serum [Figure 3]. Calculated % accuracy and %CV values were within the acceptance criteria at each concentration level.

Row	Component Name	Actual Concentration	Num. Values	Mean	Standard Deviation	Percent CV	Average Accuracy across Replicates
1	IYPTNGYTR_01	0.050	3 of 3	0.050	0.003	6.50	100.
2	IYPTNGYTR_01	0.100	3 of 3	0.100	0.012	12.1	99.6
3	IYPTNGYTR_01	5.000	3 of 3	5.139	0.101	1.97	103.
4	IYPTNGYTR_01	10.000	3 of 3	10.906	0.270	2.47	109.
5	IYPTNGYTR_01	50.000	3 of 3	47.074	1.475	3.13	94.1
6	IYPTNGYTR_01	100.000	3 of 3	96.584	1.228	1.27	96.6
7	IYPTNGYTR_01	250.000	3 of 3	244.158	4.464	1.83	97.7

Figure 3: Quantitative performance for TDx (m/z 542.70 → m/z 404.701) analysis. Reproducibility and accuracy results were determined from the calibration curve standards across 3 replicates at each concentration. Statistical results were summarized using the Analytics module in SCIEX OS software.

Compliance-ready SCIEX OS software

Equivalent SCIEX OS software capabilities for regulated bioanalysis can be executed on the ZenoTOF 7600 system, ensuring high fidelity when performing method transfers while retaining critical compliance features.

SCIEX OS software is a closed system and requires records and signatures to be stored electronically, meeting the regulations outlined by 21 CFR Part 11. SCIEX OS software can open raw data files from any visible storage location within a closed network by using designated processing workstations. Figure 4 illustrates the features of SCIEX OS software that are used to monitor the audit trail, acquire and process data, and configure user access. The audit trail feature enables users to audit critical user actions and locks in data integrity. The Central

Administrator Console [CAC] feature allows users to centralize acquisition and processing using a single platform to maximize efficiency for multi-instrument laboratories, independent of compliance standards. The configuration module allows users to assign roles and access as the administrator, method developer, analyst, and reviewer.



Audit Trail

Easily search and filter for specific high-risk events in audit trail viewer. Built-in data integrity features allow you to tailor each functionality specifically to meet compliance needs and data security requirements.



Filter Audit Trail

Find instances where

Sample Name: [dropdown] Contains: Standard [dropdown]

Reason: [dropdown] Contains: manual integration [dropdown]

<No Filter> [dropdown] Contains: [dropdown]

☐ And where time and date are: From: 5/9/2023 12:00 AM To: 5/9/2023 12:00 AM

OK Cancel



Central Administration

Users can manage groups, role definitions, workstations and projects across all systems using the SCIEX OS software Central Administrator Console (CAC). It supports all regulated and non-regulated compliance standards.



Centralized Management

- User Management
- User Roles and Permissions
- Workstation Management
- Project Management
- Workgroup Management

Manage Workgroups and Assignments

Method Developers

Reviewers

Users

Workstations

Projects

☒ Use Central Settings for Projects

Assign or Unassign Project Root from Project Pool

Project Root Directory

Users who have local drive (local drive) assigned to CAC Demo

CAC Demo

Confirm Change Events

Change Events

Timestamp	Event Name	Event Description
5/9/2023 4:48:40 PM	Project settings enabled/ disabled in a workgroup	Project based central settings enabled for workgroup Reviewers
5/9/2023 4:48:40 PM	Projects assigned/unassigned to a workgroup	Users who have local drive (local drive) assigned to CAC Demo (Default assigned to Workgroup Reviewers)

E-Signature

Timestamp: 5/9/2023 4:48:40 PM


User Name: ARSCEDEV\BUT.Project

Full User Name: BUT.Project

Reason: Project assigned

Password: [password field]

Cancel Change Rollback Change



Configuration

Assign users and access to administrator, method developer, analyst and reviewer roles under the audit trail module. Easily customize the role and specify level of access.



Devices

Projects

User Management

Queue

Licenses

LIMS Communication

General

Software Updates

CAC

About

Users, Roles, Security

User Roles and Permission Categories

Permission	Administrator	Method Developer	Analyst	Reviewer
Batch				
Submit unlocked methods	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Open	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Save as	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Submit	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Save	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Save ion reference table	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Add data sub-folder	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Configure Decision Rules	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Configuration				
General tab	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
General: change regional setting	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
General: full screen mode	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
LIMS communication tab	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Audit maps tab	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Queue tab	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Queue: instrument idle time	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Queue: max. number of acquired samples	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Queue: other queue settings	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Duplicate a User Role

Create a new user role based on an existing role.

Existing user role: Administrator [dropdown]

New role name: [text field]

New role description: [text field]

OK Cancel

Figure 4: Features of SCIEX OS software for monitoring user access and evaluating the audit trail. The audit trail view allows users to filter for high-risk events easily and enables data integrity features to meet compliance requirements. The software features a Central Administrator Console [CAC] to manage users and groups, role definitions, workstations and projects across all systems. The CAC feature supports both regulated and non-regulated compliance standards. The configuration module enables users to quickly set up roles and levels of access for the administrator, method developer, analyst and reviewer levels.

Conclusions

- An LLOQ of 0.05 µg/mL was achieved for the quantitation of TDx in human serum
- Low-level quantitation was achieved using 50 µL human serum given the increased MS/MS sampling efficiency with Zeno MRM^{HR}
- Linearity was achieved at concentrations ranging from 0.05 µg/mL to 250 µg/mL, achieving an LDR of 3.7 orders of magnitude
- Comparable quantitative performance was demonstrated with accurate and highly reproducible (%CV <13%) results on the ZenoTOF 7600 system
- A single platform for streamlined data acquisition, processing, and management with SCIEX OS software was presented
- Retain data management and compliance-readiness (21 CFR Part 11) features using SCIEX OS software to support regulated bioanalysis on the ZenoTOF 7600 system

References

1. FDA approves fam-trastuzumab deruxtecan-nxki for HER2-positive gastric adenocarcinomas. [Jan 2021](#)
2. Nagai Y, Oitate M, Shiozawa H, Ando O. Comprehensive preclinical pharmacokinetic evaluations of trastuzumab deruxtecan [DS-8201a], a HER2-targeting antibody-drug conjugate, in cynomolgus monkeys, *Animal Pharmacokinetics and metabolism*. 49 (2019) 1086-1096. <https://doi.org/10.1080/00498254.2018.1531158>
3. Toshihiko Doi, Hiroji Iwata, Junji Tsurutani, Shunji Takahashi, Haeseong Park, Charles H. Redfern, Kohei Shitara, Chikako Shimizu, Hiroya Taniguchi, Tsutomu Iwasa, Shinichiro Taira, Albert C. Lockhart, Jennifer M. Fisher, Takahiro Jikoh, Yoshihiko Fujisaki, Caleb C. Lee, Antoine Yver, and Kenji Tamura, Safety, pharmacokinetics, and antitumour activity of trastuzumab deruxtecan [DS-8201], a HER2-targeting antibody-drug conjugate, in patients with advanced breast and gastric or gastro-oesophageal tumours: a phase 1 dose-escalation study, *The Lancet oncology*. 18(11) (2017) 1512-1522. [https://doi.org/10.1016/S1470-2045\(17\)30604-6](https://doi.org/10.1016/S1470-2045(17)30604-6)
4. Malik, P.; Phipps, C.; Edginton, A.; Blay, J. Pharmacokinetic considerations for antibody-drug conjugates against cancer. *Pharm. Res.* 34 (2017), 2579–2595. DOI: [10.1007/s11095-017-2259-3](https://doi.org/10.1007/s11095-017-2259-3)
5. Nareshkumar Jain, Sean W. Smith, Sanjeevani Ghone Bruce Tomczuk, *Current ADC linker Chemistry*, 32 (2015) 3526–3540. <https://doi.org/10.1007/s11095-015-1657-7>
6. [Bioanalytical Method Validation, May 2018.](#)

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2024 DH Tech. Dev. Pte. Ltd. MKT-33769-A



Headquarters

500 Old Connecticut Path | Framingham, MA 01701
USA
Phone 508-383-7700
sciex.com

International Sales

For our office locations please all the division headquarters or refer to our website at sciex.com/offices