

Low-pg/mL quantification of TL 13-112, a proteolysis targeting chimera (PROTAC) in rat plasma

Excellent quantitative performance for PROTACs in matrix using the SCIEX 7500 system

Ebru Selen, Rahul Baghla and Eshani Nandita
SCIEX, USA

This technical note demonstrates a highly sensitive, robust and rapid workflow to quantify TL 13-112, a selective ALK degrader, and its inactive control, TL 13-110, in rat plasma using a high-end triple quadrupole mass spectrometer. A lower limit of quantification (LLOQ) of 10 pg/mL was achieved for both TL 13-112 and TL 13-110 using a fast and simple protein precipitation method with a 10-minute LC-MS/MS analysis (Figure 1).

The interest in targeted protein degradation has shifted from academia to industry after the therapeutic potential of a PROTAC was documented in 2001.¹ PROTACs have emerged as a therapeutic modality and several candidates have moved into clinical trials.² The potential of PROTACs is coded in its structure. A linker connects a protein of interest (POI) binding moiety to a ubiquitin E3 ligase recognition moiety (Figure 2A). The heterobifunctional structure enables PROTACs to bring the POI and E3 ligase closer in proximity. This induces the

ubiquitination of the POI, which is then targeted by the disposal machinery of the cell.²

One of the many attractive hallmarks of PROTACs is their high potency in nanomolar drug concentrations.³ While their potential is well-documented,¹ challenges remain for the analysis of PROTACs. Sensitive and selective assays for high-confidence detection and quantification of PROTACs are needed to ensure the safety and efficacy in the drug development pipeline.

Here, a highly sensitive assay for the quantification of PROTACs in a complex matrix was demonstrated. The quantitative performance of the assay was evaluated using the commercially available TL 13-112 (PROTAC) and TL 13-110 (inactive control) structures. Quantification at low-pg/mL levels was achieved for both analytes in rat plasma using the SCIEX 7500 system. The front-end enhancements of the system facilitated greater sensitivity, which improved overall ion generation, capture and transmission.




Figure 1. Representative extracted ion chromatograms (XICs) of the selective ALK degrader (PROTAC) and its inactive control in matrix and at the LLOQ level. Low-pg/mL quantification without matrix interference was achieved using the SCIEX 7500 system equipped with SCIEX OS software.

Key features of the quantification of PROTACs using the SCIEX 7500 system

- **New levels of quantification of low-dose, high-potency drug modalities:** Achieve low-pg/mL level LLOQs for the quantification of PROTACs in rat plasma on the SCIEX 7500 system equipped with an innovative front-end design
- **Robust analytical performance:** Reach exceptional quantitative performance with strong linearity and excellent accuracy and precision for low-level quantification
- **Streamlined data management:** Employ fast, intuitive and integrated data acquisition and processing using SCIEX OS software

Methods

Sample preparation: Commercially available individual PROTAC degrader (TL 13-112) and its inactive control (TL 13-110) were reconstituted in DMSO. PROTACs were spiked into 100 μ L of rat plasma at concentrations ranging from 10 pg/mL to 15000 pg/mL. Protein precipitation was performed with 600 μ L of 1:1 (v/v), acetonitrile/methanol. Samples were vortexed for 30 seconds and then centrifuged at 13000 rpm for 12 minutes at room temperature. The supernatant was transferred to a new Eppendorf tube and dried under nitrogen flow. Samples were reconstituted using 200 μ L of 1:1 (v/v), methanol/acetonitrile prior to analysis.

Chromatography: Sample separation was performed using an ExionLC system at a flow rate of 0.3 mL/min using a [Phenomenex Kinetex XB-C18 \(2.1 x 50 mm, 1.7 \$\mu\$ m, 100 \$\text{\AA}\$ \)](#) column. A 10-minute gradient was used for analysis (Table 1). Mobile phase A was 0.1% formic acid in water and mobile phase B was 0.1% formic acid in acetonitrile. The column temperature was kept at 40°C. An injection volume of 10 μ L was used for analysis. A mixture with equal parts by volume of acetonitrile, methanol and water was used as the needle wash solvent.

Table 1. Chromatography gradient for the ALK degrader and its inactive control.

Time (min)	Mobile phase A (%)	Mobile phase B (%)
0	85	15
0.2	85	15
5	50	50
5.5	5	95
8.5	5	95
8.6	85	15
10	85	15

Mass spectrometry: The MRM transitions used are summarized in Table 2. The optimized analyte-dependent MRM parameters are listed in Table 2. The optimized source and gas parameters are listed in Table 3.

Table 2. MRM parameters used for quantification.

ID	Precursor ion (m/z)	Fragment ion (m/z)	CE (V)	CXP (V)
TL 13-112	1002.2	584.1	70	15
TL 13-110	988.7	542.4	70	20

Table 3. Source and gas parameters.

Parameter	Value
Polarity	Positive
Ion source gas 1	55 psi
Ion source gas 2	65 psi
Curtain gas	45 psi
Source temperature	600°C
Ion spray voltage	3000 V
CAD gas	7

Data processing: Data collection, analysis and quantification were performed using SCIEX OS software, version 3.0. Peaks were automatically integrated using the MQ4 algorithm and a weighting of $1/x^2$ was used for quantification.

Quantitative performance

Given the high potency of PROTACs, sensitive and robust bioanalytical methods are needed for accurate quantification to ensure proper safety and efficacy during pre-clinical evaluation. This technical note demonstrates a low-pg/mL level quantification assay of a PROTAC and its inactive control in rat plasma using the SCIEX 7500 system.

A calibration curve was prepared, as described in the sample preparation section, for concentrations ranging from 10 pg/mL to 15000 pg/mL. Individual concentrations were run in triplicate.

An LLOQ of 10 pg/mL was achieved for both TL 13-112 and TL 13-110 (Figure 2B). No interferences were observed in the matrix blank (rat plasma) for either analyte (Figure 2B). Strong linearity was achieved for both analytes and the linear dynamic range (LDR) spanned 3.2 orders of magnitude (Figure 3).

Analytical performance was evaluated based on the requirement that the accuracy of the calculated mean should be between 80% and 120% at the LLOQ and between 85% and 115% at higher concentrations. The %CV of the calculated mean of the concentration should be below 20% at the LLOQ and below 15% for all higher concentrations.

Accuracy was within $\pm 11\%$ and $\pm 12\%$ of the nominal concentration for TL 13-112 and TL 13-110, respectively (Table 4). The %CV was $< 10\%$ for both analytes, as listed in Table 4. Calculated values for accuracy and %CV were within the acceptance criteria at each concentration level (Table 4).




Figure 3. Calibration curves for TL 13-112 and TL 13-110 in rat plasma. Strong linearity was achieved for TL 13-112 (top panel) and TL 13-110 (bottom panel) in rat plasma, with a correlation coefficient (r^2) of 0.99 for both targets. Each concentration was run in triplicate.




Figure 2. Low-pg/mL level quantification was achieved for the PROTAC and its inactive control. The PROTAC used for this assay was TL 13-112 and the inactive control was TL 13-110. The structure of TL 13-112 contains an additional carbonyl oxygen on the pomalidomide group compared to the structure of TL 13-110 (see the part of structure in bold, A). The blue rectangle highlights the POI binding moiety and the magenta rectangle highlights the E3 ligase binding moiety (A). Representative XICs of the PROTAC and its inactive control are shown (B). An LLOQ of 10 pg/mL was achieved for the PROTAC and its inactive control. No matrix interferences were observed in the plasma blank for either analyte.

Table 4. Accuracy and %CV at each concentration level measured for TL 13-112 and TL 13-110. Each concentration level was run in triplicate.

Concentration (pg/mL)	TL 13-112		TL 13-110	
	Accuracy (%)	CV (%)	Accuracy (%)	CV (%)
10	103	8.76	102	7.34
15	93.6	3.74	108	9.44
30	102	9.33	107	2.02
75	98.3	5.75	108	2.19
375	103	5.94	98.1	1.70
1500	110	1.71	97.7	1.37
7500	100	3.07	95.0	1.26
15000	89.0	2.45	88.2	3.57

Conclusions

- An LLOQ of 10 pg/mL was reached for the quantification of PROTACs in rat plasma with minimal sample preparation
- A highly sensitive assay for the quantification of PROTACs was demonstrated on the SCIEX 7500 system with an improved front-end technology for better ion generation, capture and transmission
- Excellent linearity, accuracy and precision were achieved for the concentrations analyzed, demonstrating exceptional quantitative performance
- Streamlined data acquisition, processing and management were performed using SCIEX OS software

References

1. Kathleen M. Sakamoto, Kyung B. Kim, Akiko Kumagai, Frank Mercurio, Craig M. Crews and Raymond J. Deshaies (2001). Protacs: Chimeric molecules that target proteins to the Skp1–Cullin–F box complex for ubiquitination and degradation. [PNAS 98\(15\): 8554-8559.](#)
2. Miklós Békés, David R. Langley and Craig M. Crews. PROTAC targeted protein degraders: the past is prologue. [Nat Rev Drug Discov 21, 181–200 \(2022\).](#)
3. Ming He, Wenxing Lv and Yu Rao (2021). Opportunities and Challenges of Small Molecule Induced Targeted Protein Degradation. [Frontiers in Cell and Developmental Biology, 9.](#)

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2022 DH Tech. Dev. Pte. Ltd. MKT-26576-A



Headquarters
 500 Old Connecticut Path | Framingham, MA 01701 USA
 Phone 508-383-7700
sciex.com

International Sales
 For our office locations please call the division headquarters or refer to our website at sciex.com/offices