

# Rapid determination of nitrofurantoin metabolite residues in aquatic products

Using the SCIEX Triple Quad™ 3500 LC-MS/MS System

Zhai Nannan, Jia Yanbo, Jin Wenhai  
SCIEX, China

Nitrofurantoin drugs are widely used in the livestock, poultry and aquaculture industries because of their low price and effectiveness in treating enteritis, scabies, red fin disease, and ulcers caused by Escherichia coli or Salmonella. Due to the carcinogenic and teratogenic side effects of nitrofurantoin and their metabolites, individual countries have banned the use of nitrofurantoin in livestock, poultry, and aquatic animal foods, and have strictly enforce residue detection. Per Announcement No. 235 issued by the Ministry of Agriculture of the People’s Republic of China on December 24, 2002 and Announcement No. 560 issued on October 28, 2005, nitrofurantoin should not be detectable in animal foods. Since the announcements were issued, the use of nitrofurantoin drugs in the feeding of animals has become illegal.


There are four common nitrofurantoin drugs: furazolidone, furantone, nitrofurantoin, and furancillin. Because nitrofurantoin prototype drugs are rapidly metabolized in the organism, and their metabolites (AOZ, AMOZ, AHD, SEM) and proteins are quite stable, the detection of metabolites is often used to reflect the residual status of nitrofurantoin drugs.



Here, an accurate quantitative method for nitrofurantoin metabolites was developed on the SCIEX Triple Quad 3500 LC-MS/MS System. This method provides a simple and fast solution to the problem of nitrofurantoin compound residues in animal-derived foods.

## Key features of method

- Highly sensitive method—a quantification limit of 0.05 µg/kg was achieved. This sensitivity is ~10 times better than required by the GB method (GB 21311-2007, Ministry of Agriculture 783-1-2006 and 781-4-2006 have a quantification limit of 0.5 µg/kg).
- High reproducibility of detection—less than 5.42% CV at quantification limit
- Method extraction recovery rate between 70%-120%, ensuring accuracy and reliability of actual sample detection



| Row | Component Name | Actual Concentration | Num. Values | Mean   | Standard Deviation | Percent CV |
|-----|----------------|----------------------|-------------|--------|--------------------|------------|
| 1   | SEM 1          | 0.0500               | 4 of 4      | 0.0500 | 0.0015             | 3.07       |
| 2   | AOZ 1          | 0.0500               | 4 of 4      | 0.0500 | 0.0010             | 2.06       |
| 3   | AHD 1          | 0.0500               | 4 of 4      | 0.0500 | 0.0027             | 5.42       |
| 4   | AMOZ 1         | 0.0500               | 4 of 4      | 0.0500 | 0.0015             | 3.06       |

Figure 1. Typical chromatogram of nitrofurantoin metabolites. (Top) Signal for nitrofurantoin metabolites at a concentration of 0.05 µg/kg. (Bottom) High reproducibility was obtained with %CV < 5.42.

## Methods

**Sample preparation:** Weigh 2g of homogenized sample in a 50 mL centrifuge tube. Add 0.05 mL of mixed isotope internal standard (100 ng/mL) and vortex for 30 seconds. Add 5 mL hydrochloric acid solution (0.2 mol/L) and 0.2 mL 2-nitrobenzaldehyde solution (0.05 mol/L), then vortex for 30 seconds and derivatize in a constant temperature water bath at 37 °C for 16 hours

Next, take the centrifuge tube to cool to room temperature. Add 3-5 mL of dipotassium hydrogen phosphate solution (1.0 mol/L) and adjust the pH to 7.0-7.5. Add 4 mL of ethyl acetate and vortex for 1 minute then centrifuge at 6,000 rpm for 5 minutes. Next, put 2 mL of the supernatant into a clean, 10 mL glass tube, then add 4 mL of ethyl acetate to the residue and vortex for 1 minute. Centrifuge at 6,000 rpm for 5 minutes. Take 3 mL of the supernatant and combine it into the above 10 mL glass tube. Combine the supernatant and blow dry with nitrogen at 40 °C.

Reconstitute the samples in 1 mL acetonitrile/water = 1/9 (v/v), perform LC-MS/MS analysis.

**Chromatography:** Separation was performed using the ExionLC™ System and a Phenomenex Kinetex C18 column (50x3.0 mm, 2.6 µm) and a flow rate of 0.4 mL/min. Column temperature of 40 °C and an injection volume of 20 µL.

**Table 1. Chromatography.**

| Time (mins) | % A | % B |
|-------------|-----|-----|
| 0.0         | 97  | 3   |
| 1.0         | 97  | 3   |
| 6.0         | 49  | 60  |
| 6.5         | 5   | 95  |
| 7.0         | 5   | 95  |
| 7.1         | 97  | 3   |
| 9.0         | 97  | 3   |

Mobile phase A: water (5 mM ammonium formate)

Mobile phase B: acetonitrile

**Mass spectrometry:** MS analysis was performed on the SCIEX Triple Quad 3500 System using an ESI source operated in positive ion mode.


Ion source parameters were:

- IS voltage (ISV): 5500 V
- Curtain gas (CUR): 20 psi

- GS1: 50 psi
- GS2: 60 psi
- Source temperature TEM: 550 °C
- Collision gas CAD: 9 psi

**Table 1. Nitrofuran metabolites and isotope internal standard mass spectrometry parameters.**

| Compound   | Q1    | Q3    | DP | CE |
|--|-------|-------|----|----|
| AOZ  | 236.1 | 133.9 | 80 | 17 |
|  |       | 103.9 | 80 | 31 |
| AOZ-D <sub>4</sub>                                 | 240.0 | 134.0 | 80 | 17 |
| AHD  | 249.2 | 134.1 | 80 | 17 |
|  |       | 104.1 | 80 | 27 |
| AHD- <sup>13</sup> C <sub>3</sub>                  | 252.0 | 134.1 | 80 | 17 |
| SEM  | 209.2 | 166.2 | 80 | 14 |
|  |       | 192.1 | 80 | 16 |
| SEM- <sup>13</sup> C- <sup>15</sup> N <sub>2</sub> | 212.0 | 168.0 | 80 | 14 |
| AMOZ   | 335.2 | 291.1 | 80 | 17 |
|  |       | 262.2 | 90 | 23 |
| AMOZ-D <sub>5</sub>                                | 340.0 | 296.0 | 80 | 17 |



**Figure 2. Calibration curve of nitrofuran metabolites.** Good linearity was obtained across the concentration range of 0.05-5.0 µg/kg ( $r > 0.99$ ).

**Table 3. Extraction and recovery rate of nitrofurantoin metabolites.**

| Compound | Extraction recovery rate | CV% |
|----------|--------------------------|-----|
| AOZ      | 98.9                     | 3.6 |
| AHD      | 101.7                    | 1.9 |
| SEM      | 94.5                     | 1.5 |
| AMOZ     | 101.2                    | 4.2 |

## Summary

Here, a rapid and accurate LC-MS/MS detection method for nitrofurantoin metabolites was established on the SCIEX Triple Quad 3500 System. The sensitivity of this method is 10 times better than the GB method. It fully meets the quantification limit requirements of GB/T 21311-2007, Ministry of Agriculture 783-1-2006 and 781-4-2006. And the method has a high extraction recovery rate, which can ensure the accuracy of real sample measurement.

## Results

Chromatography was developed that enabled the separation of nitrofurantoin metabolites using a 9 min run time (Figure 1). Very good reproducibility of signal was observed at a concentration of 0.05 µg/kg. Concentration curves were generated across a concentration range of 0.05 to 5.0 µg/kg and demonstrated very good linearity (Figure 2) ensuring accurate quantification of samples at different concentration levels. Under the conditions of this method, the extraction recovery rate for the 4 compounds tested was between 94.5 and 101.7% (Table 3).

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries.

© 2020 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-12418-A. AB SCIEX™ is being used under license.



**Headquarters**  
500 Old Connecticut Path | Framingham, MA 01701 USA  
Phone 508-383-7700  
[sciex.com](https://sciex.com)

**International Sales**  
For our office locations please call the division headquarters or refer to our website at [sciex.com/offices](https://sciex.com/offices)