

Rapid determination of sulfonamide and quinolone residues in aquatic products

Quantitative accuracy with the SCIEX Triple Quad™ 3500 LC-MS/MS System

Zhai Nannan, Jia Yanbo, Jin Wenhai
SCIEX, China

Sulfonamides and quinolones are common synthetic antibiotics. Because of their effectiveness and low cost, they are widely used in livestock and aquaculture. The extensive use and large dosage leads to the presence of residues, which is a serious health concern. Ingestion of these antibiotics, by any means, will lead to accumulation in the human body. Sulfonamide residues will destroy the body's hematopoietic functions, and quinolone residues are potentially carcinogenic and genotoxic. This presents a significant health safety hazard for the population.

China has clear requirements for the limits of sulfonamides and quinolones in aquatic products. Therefore, aquatic products must be tested for these two types of drug residues. Many products require testing, and the requirements are rigorous. The basis of existing sulfonamides detection in China is Ministry of Agriculture 1025-23-2008, with a quantification limit of 0.5 µg/kg. The basis of quinolone detection is GB 20366-2006, with a quantification limit of 1.0 µg/kg. Both of these drugs are addressed by Ministry of Agriculture Announcement No. 1077-1-2008, which requires a quantification limit of 2.0 µg/kg.



Here, an accurate, quantitative method for the simultaneous analysis of sulphonamides and quinolones was established on the SCIEX Triple Quad™ 3500 LC-MS/MS System, to meet the specific regulation requirements. The method provides a quick and simple solution to the problem of sulfonamide and quinolone residues in aquatic products.

Key Feature of method

- Highly sensitive: the quantification limit for sulfonamides and quinolones is 0.1 µg/kg, which is 5 times lower than the limit required by Ministry of Agriculture 1025-23-2008, 10 times lower than the required limit for quinolone in GB 20366-2006, and 20 times lower than the limit required by Ministry of Agriculture Announcement No. 1077-1-2008
- The method extraction recovery rate is in the range of 70%-120%, showing the accuracy and reliability of actual sample detection




Figure 1. Representative chromatography of various substances extracted from fish substrate at a concentration of 1 µg/kg

Methods

Sample Preparation: Weigh 2.5 g of homogenized fish sample in a 15 mL centrifuge tube. Add 5 mL of 5% perchloric acid and vortex for 1 minute. Centrifuge at 10,000 rpm for 10 minutes. Purify all the supernatants by SPE. Dry the eluent with nitrogen. Reconstitute the samples in 0.5 mL methanol/water = 1/9 (v/v) before performing LC-MS/MS analysis.

Chromatography: Separation was performed using the ExionLC™ System and a Phenomenex Kinetex C18 column (50×3.0 mm, 2.6 µm) and a flow rate of 0.4 mL/min. Column temperature of 40 °C and an injection volume of 20 µL.

Table 1. Chromatography.

Time (mins)	% A	% B
0.0	97	3
1.1	97	3
3.5	70	30
4.5	70	30
6.4	25	75
6.46	5	95
8.2	5	95
8.21	97	3
9.5	97	3

Mobile phase A: water + 0.1% formic acid

Mobile phase B: methanol + 0.1% formic acid

Mass spectrometry: MS analysis was performed on the SCIEX Triple Quad 3500 System using an ESI source operated in positive ion mode.

Ion source parameters were:

- IS voltage (ISV): 5500 V
- Curtain gas (CUR): 20 psi
- GS1: 50 psi
- GS2: 60 psi
- Source temperature TEM: 550 °C
- Collision gas CAD: 9 psi

Table 1. Nitrofurans metabolites and isotope internal standard mass spectrometry parameters.

Compound	Q1	Q3	DP	CE
Sulfamethazine	279.1	186.1	60	23
		156.0	60	27
Sulfadoxine	311.1	156.1	70	30
		108.2	70	37
Sulphadimethoxine	311.1	156.1	70	28
		218.0	70	28
Sulfamethoxazole	254.1	156.0	65	22
		108.0	65	36
Sulfabenzamide	277.1	156.0	60	19
		108.0	60	32
Sulfaquinoxaline	301.1	156.0	80	24
		108.0	80	36
Trimethoprim	291.1	230.1	95	33
		123.1	95	34
Sulfaphenazole	315.0	156.0	90	27
		108.0	90	40
Lomefloxacin	352.0	265.0	80	33
		308.1	80	28
Nalidixic acid	233.0	215.0	68	18
		187.0	68	34
Flumequine	262.1	244.1	77	23
		202.1	77	42

Results


Chromatography was developed to separate sulfonamides and quinolones using a 9.5 min run time (Figure 1). Concentration curves were generated and each substance demonstrated good linearity between 0.1-10.0 µg/kg ($r > 0.99$), ensuring accurate quantification (Figure 2, top). Statistics for the concentration curves showed very good reproducibility, the method had %CVs of 5.73 at a concentration of 0.1 µg/kg (Figure 2, bottom).

The extraction recovery rate was also evaluated. Under the conditions of this method, the extraction recovery rate of the 11 compounds is 71.2%-116.4% (0.25 µg/kg).

Table 3. Extraction and recovery rate of sulfonamides and quinolones.

Compound	Extraction recovery rate	CV%
<i>Sulfamethazine</i>	71.2	5.0
<i>Sulfadoxine</i>	84.6	3.2
<i>Sulphadimethoxine</i>	92.8	0.9
<i>Sulfamethoxazole</i>	92.8	0.9
<i>Benzosulfonamide</i>	87.2	3.3
<i>Sulfaquinoxaline</i>	91.4	1.5
<i>Trimethoprim</i>	90.1	5.1
<i>Sulfapyrazole</i>	93.8	5.5
<i>Lomefloxacin</i>	116.4	2.5
<i>Nalidixic acid</i>	106.9	1.5
<i>Flumequine</i>	116.2	5.1

• Calibration for Sulfamethazine $1y = 2.7229065x + 3089.6303 (r = 0.9933)$ (weighting $1/x^2$)
 • Calibration for Sulfadoxine $1y = 2.193945x + 10392.11465 (r = 0.9904)$ (weighting $1/x$)
 • Calibration for Sulphadimethoxine $1y = 3.95726x + 7603.9228 (r = 0.9935)$ (weighting $1/x^2$)
 • Calibration for Sulfamethoxazole $1y = 7.703304x + 868.39935 (r = 0.9949)$ (weighting $1/x^2$)
 • Calibration for Benzosulfonamide $1y = 151.0264x + 1039.12152 (r = 0.9919)$ (weighting $1/x^2$)
 • Calibration for Sulfaquinoxaline $1y = 1.02956x + 2959.49918 (r = 0.9919)$ (weighting $1/x^2$)
 • Calibration for Trimethoprim $1y = 2.225246x + 2999.61752 (r = 0.9933)$ (weighting $1/x^2$)
 • Calibration for Sulfapyrazole $1y = 5.616504x + 625.1656 (r = 0.9946)$ (weighting $1/x^2$)
 • Calibration for Lomefloxacin $1y = 3.861594x + 339.27510 (r = 0.9958)$ (weighting $1/x^2$)
 • Calibration for Nalidixic acid $1y = 2.13816x + 4729.0303 (r = 0.9919)$ (weighting $1/x^2$)
 • Calibration for Flumequine $1y = 1.29756x + 1775.64918 (r = 0.9913)$ (weighting $1/x^2$)



Row	Component Name	Actual Concentration	Num. Values	Mean	Standard Deviation	Percent CV
1	Sulfamethazine 1	0.1000	3 of 3	0.1000	0.0051	5.08
2	Sulfadoxine 1	0.1000	3 of 3	0.1000	0.0057	5.73
3	Sulphadimethoxine 1	0.1000	3 of 3	0.1000	0.0045	4.47
4	Sulfamethoxazole 1	0.1000	3 of 3	0.1000	0.0045	4.47
5	Sulfabenzamide 1	0.1000	3 of 3	0.1000	0.0055	5.53
6	Sulfaquinoxaline 1	0.1000	3 of 3	0.1000	0.0032	3.22
7	Trimethoprim 1	0.1000	3 of 3	0.1000	0.0044	4.37
8	Sulfapyrazole 1	0.1000	3 of 3	0.1000	0.0036	3.64
9	Lomefloxacin 1	0.1000	3 of 3	0.1000	0.0048	4.77
10	Nalidixic acid 1	0.1000	3 of 3	0.1000	0.0038	3.78
11	Flumequine 1	0.1000	3 of 3	0.1000	0.0029	2.88

Figure 2. Linear curves of sulfonamides and quinolones. Reproducibility results of sulfonamides and quinolones in fish substrate (0.1 µg/kg) are shown.

Conclusions

Here, a set of LC-MS/MS methods were developed for the detection of 11 sulfonamides and quinolones in aquatic products using the SCIEX Triple Quad 3500 System. This solution fully meets the requirements of the existing GB method detection standards. The method has simple and rapid sample processing, high method sensitivity, and demonstrates accurate quantification, from sample preparation to LC-MS/MS.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries.

© 2020 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-12702-A. AB SCIEX™ is being used under license.



Headquarters
 500 Old Connecticut Path | Framingham, MA 01701 USA
 Phone 508-383-7700
sciex.com

International Sales
 For our office locations please call the division headquarters or refer to our website at sciex.com/offices