

Detection of sesame peptides and 12 other food allergens using the SCIEX vMethod for food allergen testing

Sesame allergen screening using the SCIEX 7500 system

Simon Roberts¹, Karl Oetjen¹ and Jianru Stahl-Zeng²
¹SCIEX, USA, ²SCIEX, Germany

In 2021, the Food Allergy Safety, Treatment, Education, and Research (FASTER) Act was passed in the United States. The FASTER Act requires all foods sold in the United States that contain sesame to declare it as an ingredient or state “Contains: Sesame” immediately after the ingredient list.¹

The SCIEX vMethod application for food allergen testing previously provided a workflow for sample preparation and LC-MS/MS detection of 12 distinct allergens, including egg, milk, almond, Brazil nut, cashew, hazelnut, pine nut, pistachio, pecan, walnut, peanut and soy^{2,3}. After evaluating 24 different sesame peptides⁴, 2 of the most sensitive peptides were selected (Figure 1) and added to the SCIEX vMethod (Table 1).




Figure 1. LC-MS/MS chromatograms for Ses3, Ses4, Ses6 and Ses7 peptides of sesame proteins.

Ten replicate injections of cookie dough spiked with 5 ppm sesame flour (Figure 2) showed good repeatability with a CV of 3.7%. Good linearity of $r^2 = 0.998$ was observed in a curve prepared with 5, 25, 100 and 250 ppm sesame flour. By including the MRM parameters for sesame peptides in the SCIEX vMethod for food allergen testing, the SCIEX 7500 system was able to reach a detection limit of 5 ppm for both sesame peptides tested.




Figure 2. Chromatograms for 2 peptides in sesame-free cookie dough (light) and cookie dough spiked with 5 ppm sesame flour (dark).

Table 1. MRM parameters for 2 sesame peptides.

Peptide	RT (min)	Q1	Q3	CE
ISGAQPSLR_1	3.7	464.8	815.4	22
ISGAQPSLR_2	3.7	464.8	472.3	22
AFYLAGGVPR_1	6.1	525.8	556.3	23
AFYLAGGVPR_2	6.1	525.8	669.4	23

References

- [H.R.1202-](#) FASTER Act of 2021 117th Congress (2021-2022)
- [SCIEX vMethod Application](#) for Food Allergen Testing.
- New, Lee Sun *et al.* “Simultaneous Analysis of Multiple Allergens in Food Products by LC-MS/MS.” *Journal of AOAC International* **101,1** (2018): 132-145.
- Ma, Xiuli *et al.* “Comprehensive quantification of sesame allergens in processed food using liquid chromatography-tandem mass spectrometry.” *Food Control* **107** (2020): 106744.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. Or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks). © 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-14260-A