

Rapid, high-resolution isobaric separation of isotopically labeled metabolites using differential mobility separation

Improved speed and selectivity with the SCIEX TripleTOF® 6600 LC-MS/MS System and SelexION® Differential Mobility Separation Technology

Qiushi Sun¹, Richard M. Williams¹, Richard G. Kibbey¹
¹Yale University, New Haven, CT, USA

Mass spectrometry for metabolomics studies has become an important research tool for identifying and quantifying the vast diversity of small molecules present within biological systems. However, a major limitation of many metabolomics studies is that they provide a static picture of the components present at a specific point in time. In the past, stable isotope labeling has been used to provide a more dynamic picture by enabling the study of metabolic fluxes and the differentiation of metabolic pathways and networks in a variety of biological systems. They often only examine the mass of a metabolite, however, without identifying the location of the label, and they only extract a minimum of kinetic information.¹

Liquid chromatography (LC) is typically coupled with mass spectrometry to provide separation of metabolites prior to mass spectrometry analysis. Due to the enormous diversity and complexity of the metabolome, complete separation of certain metabolites such as isomeric compounds can be challenging. differential mobility spectrometry (DMS) is a method of separating ions based on the difference between ion mobility in high and low electric fields at atmospheric pressure.^{2,3} It adds a new dimension of selectivity and performance, which can reduce background noise and eliminate interferences.

In previous work, a targeted Mass Isotopomer Multi-Ordinate Spectral Analysis (MIMOSA) approach on a QTRAP® 5500 LC-




Figure 1. Separation of 3-phosphoglycerate (3PG) and 2-phosphoglycerate (2PG) using SelexION Technology. A separation voltage (SV) of 2000V and compensation voltages of -12V and -15.5V achieved separation of isotopically labeled 2PG and 3PG, respectively.



MS/MS System equipped with a SelexION Differential Mobility Separation Device was developed to measure metabolic flow of central carbon metabolism of glycolysis and the TCA cycle after isotopic label (¹³C) introduction.^{4,5}

In this study, a TripleTOF 6600 LC-MS/MS System was used with the MIMOSA workflow.⁶ The TripleTOF system has higher resolving power than the QTRAP® 5500 System, aiding in the differentiation of isotopomers. Additionally, the LC-DMS-MRM^{HR} method used on the TripleTOF 6600 System is fast and able to reduce the mass overlap issue (separate isobars) while simultaneously measuring all possible fragments from each metabolite.

Key features of LC-DMS-MRM^{HR} for MIMOSA fluxomics measurements

- Isobaric separation of isotopically labeled metabolites, that differ only in the location of the label, using differential mobility separation
- High-resolution MS/MS data for all fragments of the target isotopomers in one experiment, reducing cycle time
- Increased confidence in the identification of the isotope label position
- High resolution and reduced cycle time supports expanding analyses to a larger number of metabolites

Methods

Sample preparation: Labeled and unlabeled cells were cultured and seeded into 6-well plates. Cells were first washed quickly with 5 mM HEPES (2 mL/well), then harvested with quench buffer (150 μ L/well). Cell lysates were transferred into a 96-well plate and then lyophilized. Samples were prepared by resuspending the cell pellet in 50 μ L D₄-taurine water solution (used as the internal standard to normalize data) and then using 5 μ L supernatant for analysis.

Chromatography: Separation was performed using an ExionLC™ system and a Thermo Scientific Hypercarb Carbon column (100x4.6 mm, 3 μ m). Flow rate was 1 mL/min with a linear gradient as outlined in Table 1.

Mass spectrometry: Samples were analyzed using the TripleTOF 6600 System equipped with a SelexION Device. Data were collected using the MRM^{HR} workflow which consisted of a TOF MS scan (100 msec) and looped high resolution MS/MS scans (15 msec each) with a mass range of 15-550 for all scans. The ion source conditions were: Temp = 600 °C, ISC = -1500V, GS1 = 55, GS2 = 65 and CUR = 40. Data were also collected on a QTRAP 5500 System, also equipped with a SelexION Device, using targeted MRM scans. The ion source conditions were: Temp = 600 °C, ISC = -2000V, GS1 = 55, GS2 = 65 and CUR = 40. Compound dependent parameters for both systems are outlined in Table 2.

Table 1. Chromatographic separation.

Time	% A	%B
0	100	0
0.5	100	0
1	60	40
1.5	60	40
2	100	0
6	100	0

Mobile phase A: 15 mM ammonium formate, 10 μ M EDTA, 0.03% acetyl acetone

Mobile phase B: 60% ACN, 35% IPA, 15 mM ammonium formate

SelexION Device parameters: The DMS separation of all targeted analytes was optimized using infusion and parameters were found to be: DMS temperature (DT) = low and DMS modifier (DM) = isopropanol (IPA) at low compensation for the TripleTOF 6600 System, or IPA at high compensation for QTRAP 5500 System. The optimized separation voltages (SV) and compensation voltages (CoV) are listed in Table 2. Note that values were selected to make sure there were no interferences among isobaric compounds.

Data processing: EI-MAVEN software was used for peak picking and curation. Polly software was used for natural abundance correction and quantification.

Table 2. Parameter settings for the TripleTOF 6600 System and QTRAP 5500 System for all compounds in the MIMOSA experiment.

Compounds	TripleTOF 6600 System parameters					QTRAP 5500 System parameters					
	DP	CE	CES	SV	CoV	DP	EP	CE	CXP	SV	CoV
pyruvate	-30	-3.7	0	1500	-9.5	-30	-6	-6	-7	2400	-30.75
lactate	-40	-4	0	1500	-8.5	-40	-10	-8	-5	2800	-44.5
succinate	-35	-13	0	2500	-39.5	-35	-13	-15	-8	2600	-41.25
malate	-15	-15	5	2500	-19.75	-40	-5	-15	-11	2400	-25.5
glutamine	-30	-20	10	3000	-49.25	-60	-5	-15	-8	3400	-57.5
glutamate	-20	-15.7	0	3000	-35	-40	-10	-15	-12	3400	-51.5
PEP	-20	-15.6	0	2600	-28.25	-40	-10	-14	-6	3000	-42.5
DHAP	-20	-55	0	2600	-37.5	-26	-4	-16	-8	3000	-47.25
glycerol 3P	-30	-22	0	2200	-21.75	-47	-8	-22	-13	2500	-25.25
3-phosphoglycerate (3PG)	-10	-50	0	2000	-15.5	-40	-12	-45	-9	3400	-44.25
2-phosphoglycerate (2PG)	-40	-50	0	2000	-12	-40	-9	-21	-10	3000	-32
citrate	-20	-31	0	3400	-1.3	-40	-10	-17	-10	3400	-8
D ₄ -taurine	-20	-35	0	2200	-31.5	-40	-10	-25	-8	2800	-48.25




Figure 4. Comparison background and peaks for quantitation. (A) Extracted ion chromatogram (XIC) of fragment ion 111 from full scan MS/MS of citrate (m/z 191) from the TripleTOF 6600 System. (B) Background observed in the MS/MS for the citrate precursor for the TripleTOF 6600 System. (C) MRM data of citrate 191/111 from the QTRAP 5500 System. (D) Background observed at the citrate MRM transitions for the QTRAP 5500 System.

Resolve isobaric isomers with DMS

3-Phosphoglycerate (3PG) and 2-Phosphoglycerate (2PG) are the main metabolites in the glycolytic pathway. The only difference in their structures is the position of the phosphate group (Figure 5, top). They have the same precursor/fragment ions (185/79) and cannot be separated by LC which prevents accurate quantification (Figure 5, middle).

SeleXION Technology uses differential mobility spectrometry (DMS) to separate ions based upon their differences in ion mobility in high and low electric fields in gases at or near atmospheric pressure. With DMS, the separation waveform (SV) radially displaces ions towards one or the other electrode, depending upon high and low field mobility characteristics. The compensation voltage (CoV) restores the trajectory for a given ion to allow them to transmit through the DMS device and enter the mass spectrometer.

When DMS was used to separate 3PG and 2PG, a separation voltage of 2000V was used. The compensation voltages were tuned for both 2PG and 3PG and were found to have quite different CoVs, -12 and -15.5, respectively (Figure 5, bottom). This allows complete separation of the two metabolites with DMS for accurate quantification with minimal overlap.




Figure 5. LC Separation of 2PG and 3PG. Incomplete separation using LC alone inhibits accurate quantification for phosphoglycerate isomers.

Conclusions

Resolving mass overlap between metabolites can be challenging using traditional LC-MS techniques. This can be further compounded during MIMOSA studies when ^{13}C isotopes are introduced. The unique selectivity of the SelexION Differential Mobility Separation Device and the high-resolution of the TripleTOF 6600 LC-MS/MS System can provide clear separation of the labeled isobars, resulting in increased confidence in isotope labeling studies. Additionally, because MRM^{HR} provides high resolution data for all fragment ions of an isotopologue simultaneously in a single scan, this reduces cycle time and provides an opportunity to expand analyses to larger numbers of metabolites. This increased efficiency also allows for exploring additional metabolic position information from other fragments without adding more transitions to the method.

The MIMOSA data generated here was consistent with previous ^{13}C -enrichment results for the central carbon metabolites⁴ and provides the following advantages:

- SelexION Technology provides an additional level of separation and selectivity beyond LC
- The high-resolution MS/MS of the TripleTOF 6600 System provides clear separation of labeled isotopes for accurate identification and quantification
- The full scan high resolution MS/MS data generated using MRM^{HR} workflow decreases the total number of targeted experiments required to accurately track and quantify metabolite isotopologues allowing more metabolites and positional information to be tracked in a single experiment.

References

1. Chokkathukalam A *et al.* (2014) Stable isotope-labeling studies in metabolomics: new insights into structure and dynamics of metabolic networks. [Bioanalysis 6\(4\): 511-524.](#)
2. SelexION® Technology: The solution to selectivity challenges in quantitative analysis. [SCIEX technical note, RUO-MKT-02-3251-A.](#)
3. Differential mobility spectrometry resolves isobaric metabolite overlap for metabolic flux analysis. [SCIEX technical note, RUO-MKT-02-6150-A.](#)
4. Alves *et al.*, (2015) Integrated, Step-Wise, Mass-Isotopomeric Flux Analysis of the TCA Cycle. [Cell Metabolism 22, 936-947.](#)
5. Targeted flux analysis through the Mass Isotopomer Multi-Ordinate Spectral Analysis (MIMOSA) Workflow. [SCIEX technical note, RUO-MKT-02-12698-A.](#)
6. A Fast and Robust LC-DMS-MRM^{HR} Method to Increase Isobar Separation Power for Tracking Isotope Labels in Central Carbon Metabolism. Sun *et al.*, Citation ID 302491, MP404, Proceedings of the 68th ASMS Conference on Mass Spectrometry and Allied Topics, 2020.

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries.

© 2021 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-12892-A. AB SCIEX™ is being used under license.



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division headquarters or refer to our website at sciex.com/offices