

Analysis of post-translational modifications using fast electron-activated dissociation (EAD)

Using the ZenoTOF 7600 system with Zeno MS/MS

Nick Morrice¹, Alexandra Antonoplis² and Christie Hunter²
¹SCIEX, UK; ²SCIEX, USA

Characterization of human cellular proteomes enables the identification of disease biomarkers and new therapeutic targets. Mass spectrometry-based techniques, such as bottom-up proteomics, are widely used to analyze cellular proteomes. During bottom-up proteomics experiments, intact proteins are cleaved by enzymes to yield peptides and the resulting peptide sequences are determined by liquid chromatography-mass spectrometry (LC-MS/MS).

Several factors influence peptide MS/MS spectral data quality and depth of sequence coverage, including the data acquisition speed and fragmentation mode implemented. Data-dependent acquisition (DDA) is often used to analyze protein digests for protein identification. MS/MS sensitivity is a crucial performance attribute in DDA workflows. The ZenoTOF 7600 system uses the Zeno trap to increase the duty cycle to $\geq 90\%$ across the entire fragment ion mass range, resulting in average signal gains of 5- to 6-fold in MS/MS mode for peptides.¹⁻² These gains enable significant improvements in proteins identified using DDA workflows.³

For proteomics experiments, 2 fragmentation options are available on the ZenoTOF 7600 system: collision-induced dissociation (CID) and electron-activated dissociation (EAD). Unlike CID, which relies on the collision of ions with nitrogen gas to fragment molecules at their weakest bonds, EAD involves the capture of electrons by molecular ions to form a radical state that



fragments. CID and EAD result in different fragments for the same peptide, with CID yielding b and y ions and EAD yielding additional c' and z• ions. EAD fragmentation provides complementary sequence information to CID and often preserves post-translational modifications (PTMs) that undergo neutral loss in a CID experiment.⁴⁻⁵ Many different modifications can be labile during MS analysis, including phosphorylation, sulfation and O-GlcNAc modifications. This work focuses on DDA method development and analysis of HeLa digest using EAD, with an emphasis on site-specific localization of phosphorylation sites. In addition to phosphorylated peptides, glycosylated human peptides were also investigated using EAD.

Key features of the ZenoTOF 7600 system for Zeno EAD DDA analysis

- The ZenoTOF 7600 system supports collision-induced dissociation (CID) and electron-activated dissociation (EAD), which can provide complementary peptide sequence information for routine proteomic analysis
- The use of the Zeno trap yields a 5- to 6-fold gain in MS/MS sensitivity for peptides, enabling the identification of more high-confidence proteins and peptides.²⁻³
- Zeno EAD DDA methods using microflow chromatography were developed for sensitive analysis of HeLa digest
- EAD complements CID by providing a fast, orthogonal fragmentation strategy that can identify additional peptides
- EAD is easy to implement for DDA workflows in SCIEX OS software with automated EAD tuning and streamlined method building




Figure 1. Comparison of peptide identifications at 1% FDR as a function of fragmentation mode on the ZenoTOF 7600 system. A) Due to the higher frequency of sampling with CID, peptide identifications were approximately 2-fold higher than EAD with this fragmentation mode but EAD provided over 10,000 new peptide identifications. B) Illustration of typical fragment ions produced by CID (marked with red lines) and EAD (marked with green lines).

Methods

Sample preparation: A 100 µg sample of digested HeLa cell lysate was fractionated into 44 samples using high pH RP-HPLC, as previously described.⁶ Then, 10% of each fraction was injected for LC-MS/MS analysis. For analysis of glycosylated peptides, human serum albumin was used from human plasma.

Chromatography: Peptide digests were separated using a Waters ACQUITY UPLC M-class system with a flow rate of 6 µL/min in direct-inject mode. A Phenomenex Kinetex 2.6 µm XB-C18 LC column (100 Å, 150 x 0.3 mm) was used. A rapid linear gradient from 5–30% B over 21 minutes was used to interrogate the fractions. Mobile phase A was water with 0.1% formic acid and mobile phase B was acetonitrile with 0.1% formic acid.

Mass spectrometry: A ZenoTOF 7600 system that was equipped with the OptiFlow Turbo V source⁷ using a low microflow probe and electrode was used for all data acquisition. DDA parameters that were implemented for all experiments included a TOF MS accumulation time of 250 ms and an exclusion time of 6 seconds. Only precursors with charge states in the range of 2–5 with intensities greater than 100 cps were selected for fragmentation. For CID acquisition, the maximum number of candidate ions per cycle was 45 and the accumulation time was 20 ms. For EAD analysis of HeLa digest, the maximum number of candidate ions was 10 per cycle and an accumulation time of 50 ms was used. The reaction time for EAD was 20 ms with an electron beam current of 3000 nA with a KE of 0 eV. For analysis of glycosylated peptides, 20 candidate ions were used per cycle with a 20 ms reaction time and 25 ms total acquisition time.

Data processing: Mascot software (Matrix Science, version 2.6) was used for all data processing. Data files were searched against the Swissprot_Human database using the SCIEX EAD algorithm for EAD and the ESI-QUAD-TOF algorithm for CID. For analysis of phosphorylated peptides in HeLa digest, the following modifications were searched in Mascot software: fixed: carbamidomethyl; variable: phosphorylation (on serine, threonine and tyrosine); methionine oxidation and n-terminal acetylation. Scaffold software (Proteome Software Inc., version 4.8.5) was used to confirm MS/MS-based peptide and protein identifications. For analysis of glycation sites, hexose variable modifications were searched in Mascot software with a fixed carboxymethyl modification. For both phosphorylated and glycosylated peptides, figures were generated using Bio Tool Kit in SCIEX OS software.⁷

Determining Zeno EAD DDA parameters for proteomic analysis of labile PTMs

In order to select the acquisition settings for Zeno EAD DDA analysis of HeLa digest, a variety of parameters were tested, including electron beam current, electron kinetic energy, reaction time and accumulation time. For peptides, it was determined that using an electron KE of 0 eV and a beam current of 3000 nA promoted fragmentation and yielded high fragment ion intensities.




Figure 2. Comparing CID and EAD fragmentation of a non-phosphorylated peptide. (Top) CID fragmentation of the peptide, DSYVGDEAQSQR, yields y and b ions for peptide sequencing, with detected fragment ion masses (+1 charge state) shown in bold red in the table beneath the spectrum. Red italics indicate a charge state other than +1. (Bottom) Nearly complete z* (z+1 in table) and c' ion series were detected for the same peptide when fragmented by EAD. Together, CID and EAD provide complementary information for sequencing peptides.

A reaction time of 20 ms coupled with an accumulation time of 50 ms enabled 2 EAD reactions to be performed for each MS/MS spectrum and these parameters were employed for all experiments. The Zeno trap was turned on to improve fragment ion sensitivity, as Zeno trap activation results in a restoration of duty cycle to $\geq 90\%$ across the entire mass range of the MS/MS spectrum.¹

For analysis of glycosylated peptides, further optimization was performed and it was determined that using a method with 20 MS/MS candidates and a 20 ms EAD reaction time yielded accurate identification of challenging modifications.

Proteomic coverage of HeLa digest obtained using EAD


Forty-four fractions of HeLa digest were separated using a fast 21-minute gradient and analyzed with DDA. Both CID and EAD MS/MS fragmentation were evaluated and results were processed and visualized using Mascot software and Scaffold software. In the CID experiment, 93,866 peptides were identified at 1% FDR using Mascot software as the database search engine (Figure 1). This search was performed using the standard settings for CID fragmentation and included 1 fixed and 5 variable modifications in the search space.

When the fragmentation mode was switched to EAD, 52,905 peptides were identified at 1% FDR (Figure 1). Fewer total numbers of peptides were identified using EAD due to the slightly slower acquisition speed (2.5x slower). However, the orthogonal fragmentation mode provided complementary sequence information (Figure 2) and identifications of peptides not found using CID, yielding 11% more peptides. This highlights that EAD is able to run at a frequency compatible with the LC timescale and with DDA.

Phosphorylation site localization using Zeno EAD MS/MS spectra

Although the HeLa cell digest was not specifically prepared to preserve PTMs, 157 phosphopeptides were identified at 1% FDR with EAD fragmentation. The raw phosphopeptide MS/MS spectra were evaluated in Explorer in SCIEX OS software. Using the Bio Tool Kit, the peptide sequences could be overlaid on the spectra to evaluate the specific fragments and the fragmentation coverage (Figures 3-5).⁷


EAD spectra can be used to identify the sequence and phosphorylation site in peptides, including peptides with multiple serines (Figure 3). Other phosphorylation sites were located using EAD spectra, including a site on a 1889.88 Da peptide with 17 amino acids (Figure 4) and a site on a threonine residue near



Symbol	Res. Mass	# (N)	b	c	y	z	z + 1	z + 2	# (C)
▶ I	113.08406	1	114.09134	131.11789	960.39791	<i>933.27136</i>	934.37918	935.38701	8
S	87.03203	2	201.12337	218.14992	837.31384	<i>820.28729</i>	821.29512	822.30294	7
S	87.03203	3	288.15540	305.18195	750.28181	733.25526	734.26309	735.27091	6
S	87.03203	4	375.18743	392.21398	663.24978	646.22324	647.23106	648.23889	5
S[Pho]	166.99836	5	542.18579	559.21234	576.21776	559.19121	560.19903	561.20686	4
F	147.06841	6	689.25420	706.28075	409.21939	392.19285	393.20067	394.20850	3
S	87.03203	7	<i>776.28623</i>	793.31278	262.15098	245.12443	246.13226	247.14008	2
R	156.10111	8	<i>332.38734</i>	175.11895	158.09240	159.10023	160.10805	1	


Figure 3. Localization of phosphorylation site in a peptide with multiple serines. An EAD spectrum for peptide sequence ISSSS[Pho]FSR contains a nearly complete z^{\bullet} ($z+1$) ion series. The associated fragment ion masses are shown in the table below the spectrum for reference. The modification site was localized with 97.42% probability in Mascot software. Bold red indicates the fragments found with +1 charge state and red italics indicate fragments found with a different charge state.

the peptide N-terminus (Figure 5). The associated c' and z^{\bullet} ($z+1$) ion series clearly detailed the location of the phosphorylation sites on the peptides.



Symbol	Res. Mass	# (N)	b	c	y	z	z + 1	z + 2	# (C)
▶ R	156.10111	1	157.10839	174.13494	1890.88918	<i>1873.86263</i>	1874.87045	1875.87828	17
A	71.03711	2	228.14550	245.17205	1734.78807	<i>1717.76152</i>	1718.76934	1719.77717	16
P	97.05276	3	325.19827	342.22481	1663.75095	1646.72440	1647.73223	1648.74005	15
S[Pho]	166.99836	4	492.19663	509.22318	1566.69819	1549.67164	1550.67946	1551.68729	14
V	99.06841	5	591.26504	608.29159	1399.63983	1382.61328	1383.62110	1384.62893	13
A	71.03711	6	662.30215	679.32870	1300.63141	1283.60486	1284.61269	1285.62051	12
N	114.04293	7	776.34508	793.37163	1229.59430	1212.56775	1213.57558	1214.58340	11
V	99.06841	8	875.41350	892.44004	1115.55137	1098.52482	1099.53265	1100.54047	10
G	57.02146	9	932.43496	949.46151	1016.48296	999.45641	1000.46423	1001.47206	9
S	87.03203	10	<i>1079.88859</i>	1036.49354	959.46149	942.43494	943.44277	944.45059	8
H	137.05891	11	1156.52590	1173.55645	872.42947	855.40292	856.41074	857.41857	7
D	160.03065	12	1316.55655	1333.58710	735.39869	718.34400	719.35183	720.35965	6
D	115.02594	13	1431.58349	1448.61404	575.33990	558.31135	559.31918	560.32701	5
L	113.08406	14	1544.66756	1561.69411	460.31296	443.28641	444.29424	445.30206	4
S	87.03203	15	1631.69958	1648.72613	347.22590	330.20235	331.21017	332.21800	3
L	113.08406	16	1744.73365	1761.75720	260.19637	243.17032	244.17814	245.18597	2
K	128.09496	17	1872.87861	1889.88918	147.11280	130.08626	131.09408	132.10191	1

Figure 4. Identification of phosphorylation site in a 17-residue peptide. The presence of c' and z^{\bullet} ($z+1$) fragment ions enable the identification of a phosphorylated serine in the sequence of a 1889.88 Da peptide. The modification site was located with 100% probability in Mascot software. Detected fragment ions (+1 charge state) are tabulated in bold red, with red italics indicating a charge state other than +1.



Symbol	Res. Mass	# (N)	b	c	y	z	z + 1	z + 2	# (C)
V	99.06841	1	100.07569	117.10224	1645.77742	1628.75087	1629.75070	1630.76652	12
Y	163.06333	2	263.13902	280.16567	1545.70207	1529.68246	1530.69829	1531.68111	11
T[Pho]	181.01401	3	444.15303	461.17958	1383.64568	1366.61913	1367.62696	1368.63478	10
H	137.05891	4	581.21194	598.23849	1202.63167	1185.60512	1186.61296	1187.62077	9
E	129.04259	5	710.25454	727.28108	1065.57276	1048.54621	1049.55403	1050.56186	8
V	99.06841	6	809.32295	826.34960	936.53016	919.50362	920.51144	921.51927	7
V	99.06841	7	908.39136	925.41791	837.46175	820.43520	821.44303	822.45086	6
T	101.04768	8	1009.43904	1026.46569	738.39334	721.36679	722.37461	723.38244	5
L	113.08406	9	1122.52311	1139.54965	637.34566	620.31911	621.32693	622.33476	4
W	186.07931	10	1308.60242	1325.62897	524.26169	507.23504	508.24287	509.25069	3
Y	163.06333	11	1471.66575	1488.69230	338.18228	321.15573	322.16356	323.17138	2
R	156.10111	12	1627.76686	175.11895	158.09240	159.10023	160.10805	1	

Figure 5. Localization of a phosphorylation site near the N-terminus of a peptide with tyrosine and threonine. The c' fragment ion series enables localization of a phosphorylated threonine adjacent to a tyrosine. The c' ion series is shown in the spectrum and the modification site was located with 99.98% probability in Mascot software. Both c' and z• (z+1) ion series are highlighted in the table below the EAD spectrum and the detected fragments (+1 charge state) are highlighted in bold red.


Glycation site localization using Zeno EAD MS/MS

A panel of human peptides was analyzed to evaluate hexose and modifications using Zeno EAD DDA. While a hexose modification could not be localized using CID, EAD enabled unambiguous localization of a hexose modification on a lysine residue on a long, 24 amino acid peptide (Figure 6).

Conclusions

The Zeno EAD DDA workflow developed in this work enabled the robust identification of thousands of unique peptides in HeLa digest. High pH fractionation of HeLa digest followed by LC-MS/MS analysis with the ZenoTOF 7600 system generated a large-scale Zeno EAD DDA dataset for interrogation. Microflow chromatography was used to improve sensitivity and the Zeno trap was activated to both improve MS/MS sensitivity and yield high fragment ion intensities. This work highlights that Zeno EAD is fast and compatible with DDA experiments.

EAD provides an alternative fragmentation option for the study of proteomic samples and enables better characterization of some types of peptides, such as long tryptic peptides and peptides containing labile PTMs, compared to CID. In the EAD spectra evaluated here, intact z• (z+1) and c' ion series were readily detectable for sequencing and localization of phosphorylation and glycation sites.



Symbol	Res. Mass	# (N)	b	c	y	z	z + 1	z + 2	# (C)
L	113.08406	1	114.09134	131.11789	3019.34417	3002.31762	3003.32544	3004.33327	24
V	99.06841	2	213.15975	230.18630	2906.26010	2889.23355	2890.24138	2891.24920	23
T	101.04768	3	314.20743	337.23208	2807.19169	2790.16514	2791.17297	2792.18079	22
D	115.02694	4	429.23438	446.26093	2706.14401	2689.11746	2690.12529	2691.13311	21
L	113.08406	5	542.31844	559.34499	2591.11707	2574.09052	2575.09834	2576.10617	20
T	101.04768	6	643.36612	660.39267	2478.03300	2461.00645	2462.01428	2463.02210	19
K[Hex]	290.14779	7	933.51391	950.54045	2376.98533	2359.95878	2360.96660	2361.97443	18
V	99.06841	8	1032.58232	1049.60887	2086.83754	2069.81099	2070.81881	2071.82664	17
H	137.05891	9	1169.64123	1186.66778	1987.76912	1970.74258	1971.75040	1972.75823	16
T	101.04768	10	1270.68891	1287.71546	1850.71021	1833.68366	1834.69149	1835.69931	15
E	129.04259	11	1399.73150	1416.75805	1749.66253	1732.63598	1733.64381	1734.65164	14
C[CAM]	160.03065	12	1559.76215	1578.78870	1620.61994	1603.59339	1604.60122	1605.60904	13
C[CAM]	160.03065	13	1719.79280	1738.81935	1460.58929	1443.56274	1444.57057	1445.57839	12
H	137.05891	14	1856.85171	1873.87626	1300.55864	1283.53209	1284.53992	1285.54774	11
G	57.02146	15	1913.87318	1930.89973	1163.49973	1146.47318	1147.48101	1148.48883	10
D	115.02694	16	2022.90012	2045.92667	1106.47827	1089.45172	1090.45954	1091.46737	9
L	113.08406	17	2141.98419	2169.01073	891.45182	874.42477	875.43260	876.44042	8
L	113.08406	18	2255.06825	2272.09480	878.36726	861.34071	862.34854	863.35636	7
E	129.04259	19	2384.11084	2401.13739	765.28320	748.25665	749.26447	750.27230	6
C[CAM]	160.03065	20	2544.14149	2561.16804	636.24060	619.21405	620.22188	621.22970	5
A	71.03711	21	2615.17861	2632.20516	476.20996	459.18340	460.19123	461.19906	4
D	115.02694	22	2730.20555	2747.23210	405.17284	388.14629	389.15411	390.16194	3
D	115.02694	23	2845.23249	2862.25904	290.14590	273.11935	274.12717	275.13500	2
R	156.10111	24	3001.33360	175.11895	158.09240	159.10023	160.10806	1	

Figure 6. EAD fragmentation of a glycosylated peptide from human serum albumin protein. A hexose modification was localized on a lysine using a nearly complete c' ion series. Detected fragment ion masses are shown in red in the table below the EAD spectrum. The hexose modification was located with 100% probability in Mascot software. Bold red indicates the fragments found with +1 charge state, while the red italics indicate fragments found with a different charge state.

References

1. Qualitative flexibility combined with quantitative power. [SCIEX technical note, RUO-MKT-02-13053-A.](#)
2. Large-scale, targeted, peptide quantification of 804 peptides with high reproducibility, using Zeno MS/MS. [SCIEX community post, RUO-MKT-18-13346-A.](#)
3. Over 40% More Proteins Identified with Zeno MS/MS. [SCIEX technical note, RUO-MKT-02-13286-B.](#)
4. Tunable Electron-Assisted Dissociation (EAD) MS/MS to preserve particularly labile PTMs. [SCIEX technical note, RUO-MKT-02-13006-B.](#)
5. Peptide fragmentation using electron activated dissociation (EAD). [SCIEX community post, RUO-MKT-18-13296-A](#)
6. Large scale protein identification using microflow chromatography on the ZenoTOF 7600 system. [SCIEX technical note, RUO-MKT-18- 14415 -A.](#)
7. Single source solution for low flow chromatography. [SCIEX technical note, RUO-MKT-02-9701-A.](#)
8. Using Bio Tool Kit for interpreting peptide MS/MS data obtained using electron activated dissociation (EAD). [SCIEX community post, RUO-MKT-18-13296-A.](#)

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks).

© 2022 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-14795-A



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division
headquarters or refer to our website at
sciex.com/offices