

Quantification of 11-nor-9-Carboxy-THC and Panel of 22 Drugs in Hair using a Hybrid Triple Quadrupole Linear Ion Trap Mass Spectrometer



Xiang He¹, Adrian Taylor² and Alexandre Wang¹
¹SCIEX, 1201 Radio Rd, Redwood City, CA 94065, USA; ²SCIEX, 71 Four Valley Drive, Concord, Ontario, L4K 4V8 Canada.

INTRODUCTION

Why Hair Testing?

- Investigate past exposure towards xenobiotics (history of use over time)
 - Hair provides an opportunity to identify the drug, or drugs, many months later; even after single dosage
 - Determining recent past drug use as well as examining long-term drug history through segmental analysis
 - Can be used as indication of abstinence
- No active metabolism or excretion occurs within hair structure to remove drugs once they have been deposited
- Sample collection is non-invasive

Compound Metabolite Detection and Quantification

- Confirming presence of compounds with confidence
 - Monitoring both parent drugs and their metabolites in hair.
- 11-nor-9-Carboxy-Tetrahydrocannabinol (THC-COOH) is the secondary metabolite for Tetrahydrocannabinol with cut off level of 0.2 pg/mg.
 - Requires a unique mass spectrometric analysis (MRM³) that is only reliably achievable on a QTRAP[®] instrument.
 - We present a method on a SCIEX QTRAP[®] 6500+ LC-MS/MS system that reproducibly detects and quantifies 11-nor-9-Carboxy-THC levels in hair down to 0.04 pg/mg.

Goal


- To evaluate MRM³ quantitation method of THC-COOH (in hair) using QTRAP[®] 6500+ LC-MS/MS system.
- Develop a QTRAP[®] workflow that enables simultaneous identification and confirmation of compounds from hair through the use of *Scheduled MRMTM Pro Algorithm* – Information Dependent Acquisition – Enhanced Product Ion (sMRM –IDA-EPI).

Calibrator levels and Analytes Analyzed by the sMRM-IDA-EPI Method


Analyte	Analyte	IS	Level	ng/mg of each compound in hair	Inj. Volume (µL)	Amount on column (pg)
Amphetamine	Cocaine	THC-D3	Cal 1	0.005	1	0.18
Methamphetamine	Methadone	Benzoylgonine-D3	Cal 2	0.05	1	1.8
MDA	THC	Cocaine-D3	Cal 3	0.2	1	7.2
Ecgoninemethyl ester	Buprenorphine	Morphin-D3	Cal 4	1.0	1	36.3
MDEA	Cocaine	Amphetamine-D5	Cal 5	2.5	1	90.9
EDDP	Norocaine	Methamphetamine-D3				
Morphine	Heroin	Methadone-D3				
Benzoylgonine	Ahydroecgonine methyl ester	Buprenorphine-D4				
Cannabidiol	Cannabiol					

Sample preparation

For compounds analysed by sMRM-IDA-EPI Method



THC-COOH Sample Preparation



The procedure outlined in *Journal of Analytical Toxicology*, Vol. 32, June 2008; S. Hegstad, H.Z. Khabrani, L. Kristoffersen, N. Kunee, P.P. Lobaier and A.S. Christophersen was followed.

LC-MS/MS

Phenomenex Kinetex C18 column were used. Mobile phase A (MPA) was acetic acid in water and mobile phase B (MPB) was acetic acid in methanol. The LC flowrate was 0.5 mL/min and the LC runtime was 8.0 minutes for the MRM³ method and 6.5 minutes for the sMRM-IDA-EPI method. Data acquisition was done with Analyst 1.6.3 using *Scheduled MRMTM Pro Algorithm*.

QTRAP[®] 6500+ LC-MS/MS System Conditions used in the Analysis of THC-COOH

- MS Detection
 - Turbo VTM source, **negative** electrospray ionization
- Source/gas
 - IS: -4500
 - Cur: 30 psi
 - TEM: 660°C
 - GS1: 50
 - GS2: 60
 - CAD: 10 (MRM) or High (MRM³)

THCA	MS/MS (MS3)
Scan Type	343.1/299.2/245.2(+1)
Transition	Negative
Polarity	343.10 m/z
1st Precursor	DP -70
DP	28
CE	299.10 Da
2nd Precursor	Resolution Q3 10,000 Da/s
Resolution Q1	Unit LIT
Resolution Q3	Scan Rate
Unit	Yes (245.2 / 2 Da)
Q0 trapping	Yes
LIT fill time	250.00 msec
Fragmentation	Yes
AF2	0.17 V
Excitation Time	20.00 msec

IS used for MRM³ method

MRM³

sMRM-IDA-EPI Settings on QTRAP[®] 6500+ LC-MS/MS System


- MS Detection
 - Turbo VTM source, **positive** electrospray ionization
- Source/gas
 - IS: 4000
 - Cur: 25 psi
 - TEM: 500°C
 - GS1: 90
 - GS2: 70
 - CAD: 10
- Acquisition method
 - Mass Spec 6.505 min
 - Period 6.505 min
 - MRM
 - IDA Criteria
 - EPI
- IDA Criteria
 - Triggers EPI after Dynamic Background Subtraction of MRM Survey Scan

Scheduled MRMTM Pro Algorithm

- 2 MRM transitions per compound for ion ratio confirmation


MRM³ is necessary

0.2 pg/mg THC-COOH (cutoff), 0.4 pg on-column



Accuracy, precision, matrix effect, overall recovery


Row	Component Name	Actual Concentration	Num. Values	Mean	Standard Deviation	Percent CV	Accuracy
1	THC-COOH	0.040	2 of 3	0.040	0.004	10.52	99.18
2	THC-COOH	0.100	3 of 3	0.099	0.008	7.94	97.84
3	THC-COOH	0.200	3 of 3	0.206	0.024	11.43	103.12
4	THC-COOH	0.400	3 of 3	0.400	0.020	5.11	99.94
5	THC-COOH	1.000	3 of 3	0.996	0.059	5.93	99.63



sMRM-IDA-EPI Workflow Results


Elution Profile of 22 Analytes from 1 ng/mg sample (1 µL injection; 36.3 pg on column)

Identification and Quantification - All from a single injection




Examples of detection, quantification and simultaneous confirmation by ion ratio, library matching and RT for: Above right, MDEA calibrator series, Below cocaine metabolite from unknown spiked sample.

Unknown Spiked Sample





Representative Example of Quantitative Performance



CONCLUSIONS

- QTRAP[®] technology provides unique advantages in the ability to maximize selectivity when confirming and quantifying low level metabolites in difficult hair matrices.
 - We developed an MRM³ method that reproducibly detects and quantifies 11-nor-9-Carboxy-THC levels in hair down to 0.04 pg/mg
- QTRAP[®] workflows enable simultaneous identification and confirmation of compounds from hair through the use of *Scheduled MRMTM Pro Algorithm* –IDA – EPI
 - Confirmation by MS/MS library matching (scores >75% for all compounds) and ion ratios (<20%CV for all analytes)
 - Successful quantification of each compound in the panel with high precision.
 - Linearity from 5.0 pg/mg to 2.5 ng/mg for all compounds except heroin, ecgonine methyl ester (5.0 pg/mg to 1.0 ng/mg) with R² values >0.99


MRM³



Information Dependent Acquisition (IDA)

MRM → EPI

Multi-Target Screening and Quantitation with MS/MS Identification



MATERIALS and METHODS

Compound list and spiking solutions:


Calibrator levels for THC-COOH

Calibrator	[THC-COOH] in hair	[THC-COOH-D9] in hair	Hair weight (mg)	Recon. Vol. (µL)	Injection vol. (µL)	THC-COOH on column (pg)
Cal 1	0.02 pg/mg	10 pg/mg	20	200	20	0.04
Cal 2	0.04 pg/mg	10 pg/mg	20	200	20	0.08
Cal 3	0.1 pg/mg	10 pg/mg	20	200	20	0.2
Cal 4	0.2 pg/mg	10 pg/mg	20	200	20	0.4
Cal 5	0.4 pg/mg	10 pg/mg	20	200	20	0.8
Cal 6	1 pg/mg	10 pg/mg	20	200	20	2

RESULTS and DISCUSSION

MRM³ for quantitative workflows eliminate interference and background, provides high sensitivity MS³ fragmentation leading to improved LOQ with reproducible % CVs when dealing with tough interferences compared to an MRM approach; demonstrated by the following figures.

Detection of THC-COOH with MRM compared to MRM³



TRADEMARKS/LICENSING

For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB SCIEX Pte. Ltd. or their respective owners. AB SCIEX[™] is being used under license.