

# Central Administrator Console Software 3.3.1 Release Notes



---

## Introduction

The Central Administrator Console (CAC) software 3.3.1 benefits network administrators in regulated environments where management of large groups of people, projects, and workstations can be costly and time consuming. To help administrators manage resources more effectively, the CAC software supplies the option of to manage users, roles, and projects centrally.

The CAC software is enabled by a license and can be installed on any workstation that supports the SCIEX OS software 3.0 and Windows Server 2019.

The CAC software and the SCIEX OS software (acquisition or processing) cannot be installed on the same workstation.

---

**Note:** The CAC software uses TCP ports 63333 and 44144 for communication. The CAC server must permit inbound traffic on port 63333 and outbound traffic on ports 44144 and 63333. The computer with the SCIEX OS software must permit inbound traffic on ports 63333 and 44144 and outbound traffic on port 63333.

---

---

**Note:** New permissions have been added to the Administrator role. The user can create customized roles with permissions to get access to and administer the CAC software.

---

---

**Note:** The CAC software does not support the Analyst software.

---

## New in Version 3.3.1

This section gives a description of the changes in the Central Administrator Console (CAC) software 3.3.1. To see the enhancements and fixed issues for a previous version of the CAC software, refer to the document: *Release Notes* that came with that version of the software.

### Features and Enhancements in Version 3.3.1

- The following user permissions have been added:
  - Analytics workspace: **Modify Weight/Volume**
  - Queue workspace: **Edit Sample**
- The project audit trail shows additional detail for the following events:
  - Analytics workspace

- **Processing method saved**
- **Project default settings changed**
- Batch workspace
  - **Batch saved**
  - **Batch submitted**
- Explorer workspace
  - **Print**
- MS Method workspace
  - **MS method saved**
- LC Method workspace
  - **LC method saved**
- Queue workspace
  - **Sample acquisition has completed**
  - **Sample Starts to Acquire**
- The following audit events have been added to the project audit trail:
  - Analytics workspace: **Weight/Volume changed**
  - Queue workspace: **Sample edited**
- The *Help* has been updated to an `html` file format.

## Notes on Use and Known Issues

### Notes on Use

Regulated customers: If user management settings are imported after software validation, then we recommend that customers use their internal change control process to document the configuration changes.

- The Central Administrator Console (CAC) software is compatible with the SCIEX OS software 3.0 or later.
- The CAC software only supports the SCIEX OS software. The Analyst Administrator Console (AAC) is used to manage the Analyst software settings.
- The CAC software can be installed on the same computer as the Flexnet Server.
- The CAC software supports cross-domain users. All users do not have to be in same Active Directory domain.

## Known Issues

- The workstation report does not have the workstation **Status** and last **Status Update Time**. (ONYX-24475)

## Contact Us

### Customer Training

- In North America: [NA.CustomerTraining@sciex.com](mailto:NA.CustomerTraining@sciex.com)
- In Europe: [Europe.CustomerTraining@sciex.com](mailto:Europe.CustomerTraining@sciex.com)
- Outside the EU and North America, visit [sciex.com/education](https://sciex.com/education) for contact information.

### Online Learning Center

- [SCIEX Now Learning Hub](#)

## SCIEX Support

SCIEX and its representatives maintain a staff of fully-trained service and technical specialists located throughout the world. They can answer questions about the system or any technical issues that might arise. For more information, visit the SCIEX website at [sciex.com](https://sciex.com) or contact us in one of the following ways:

- [sciex.com/contact-us](https://sciex.com/contact-us)
- [sciex.com/request-support](https://sciex.com/request-support)

## CyberSecurity

For the latest guidance on cybersecurity for SCIEX products, visit [sciex.com/productsecurity](https://sciex.com/productsecurity).

## Documentation

This version of the document supercedes all previous versions of this document.

The latest versions of the documentation are available on the SCIEX website, at [sciex.com/customer-documents](https://sciex.com/customer-documents).

---

**Note:** To request a free, printed version of this document, contact [sciex.com/contact-us](https://sciex.com/contact-us).

---

---

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate such products as those manufacturers' products and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are the sole and exclusive representations, warranties, and obligations of SCIEX. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

(GEN-IDV-09-10816-E)

For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd., or their respective owners, in the United States and/or certain other countries (see [sciex.com/trademarks](http://sciex.com/trademarks)).

AB Sciex™ is being used under license.

© 2023 DH Tech. Dev. Pte. Ltd.



AB Sciex Pte. Ltd.  
Blk33, #04-06 Marsiling Industrial Estate Road 3  
Woodlands Central Industrial Estate, Singapore 739256