

# X500R QTOF

A COMPENDIUM OF  
RELATED RESEARCH  
APPLICATIONS



1 2



# Contents

## Food Testing

Pesticide Analysis in Food

Antibiotic Analysis in Food

Mycotoxin Analysis in Food

Pesticide and Potency Testing for the Cannabis Industry

Using the X500R QTOF System and SCIEX OS Software to Identify and Quantify Food Residues

Using the X500R QTOF System and SCIEX OS Software to Quickly Identify Unknowns in Food Samples

## Water Analysis

Pesticide Analysis in Water

PPCP Analysis in Water

Illicit Drugs Analysis in Water

## Forensic Analysis

Forensic Identification and Quantitation Workflows Delivered on a Revolutionary Designed QTOF and SCIEX OS Software

Forensic Drug Screening Analysis - Urine

Forensic Drug Screening Analysis - Blood


# X500R QTOF FOOD TESTING

# Food Method

## Pesticide analysis in food

### Elevate your food testing with the X500R QTOF System

Method details and access to HR-MS/MS libraries to detect, quantify, and confirm pesticides in food extracts using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software..





### LC Method

Column	Phenomenex Kinetex Biphenyl, 50 x 2.1 mm, 2.6 µm	
Mobile Phase A	5 mM ammonium formate in water	
Mobile Phase B	5 mM ammonium formate in methanol	
Flow rate	0.5 mL/min	
Column temperature	40°C	
Injection volume	2 µL	
Gradient profile	<b>Time (min)</b>	<b>% B</b>
	0	10
	0.5	10
	2.00	30
	9.0	60
	11.0	80
	12.0	95
	15.0	95
	16.0	10
	20.0	10

Suggested sample preparation conditions based on the QuEChERS method (QuEChERS European standard method 15662).






# Food Method

## Antibiotic analysis in food

### Elevate your food testing with the X500R QTOF System

Method details and access to HR-MS/MS libraries to detect, quantify, and confirm antibiotic vet drugs in tissue extracts using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software..






### LC Method

<i>Column</i>	Phenomenex Gemini 3µm C18 110Å column, 50 x 2.0mm	
<i>Mobile Phase A</i>	0.1% formic acid in water	
<i>Mobile Phase B</i>	0.1% formic acid in methanol	
<i>Flow rate</i>	0.5 mL/min	
<i>Column temperature</i>	40°C	
<i>Injection volume</i>	10 µL	
<i>Gradient profile</i>	<b>Time (min)</b>	<b>% B</b>
	0	2
	0.3	2
	7.27	80
	7.37	99
	10.9	99
	11	2
	15	2

Sample prep protocol adopted from:  
Mastovska & Lightfield, J. Chrom. A., 2008, 1202, 118-123

# MS Method



Suggested IDA (Information Dependent Acquisition) conditions for routine food contaminant testing as displayed in SCIEX OS.

# Data Processing



Review your results with utmost efficiency using SCIEX OS for simultaneous quantitation and MS/MS library confirmation.

# X500R High Resolution Libraries

Download the XIC List

Download a FREE Trial of the MS/MS Library

Download a free XIC compound list detailing a full list of antibiotic compounds including molecular formula and accurate mass.

Download a free trial of the antibiotic high resolution MS/MS library, containing 244 compounds.

[Learn more at sciex.com/X-Antibiotics.](http://sciex.com/X-Antibiotics)

AB Sciex is doing business as SCIEX.

© 2015 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-07-2858A

# Food Method

## Mycotoxin analysis in food

*Elevate your food testing with the X500R QTOF System*

Method details and access to HR-MS/MS libraries to detect, quantify, and confirm mycotoxins in food extracts using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software.



❏

### Sample Prep

Step 1

- Add 4mL extraction solvent (79/20/1. V/V/V) Acetonitrile/Water/Acetic Acid per 1g sample.

Step 2

- Extract on a rotary shaker for 90 minutes

Step 3

- Centrifuge.

Step 4

- Dilute 350µL of extract with 350µL 20/79/1 (v/v/v) Acetonitrile/Water/Acetic Acid.

Step 5

- Filter into vial for analysis.

Sample prep protocol based : Sulyok M, Krska R, Schumacher R (2010) Food Chem 119:408-416

❏

### LC Method


<i>Column</i>	Phenomenex Gemini C18, 150 x 4.6 mm, 5 µm	
<i>Mobile Phase A</i>	5 mM ammonium acetate + 1% acetic acid in water	
<i>Mobile Phase B</i>	5 mM ammonium acetate + 1% acetic acid in methanol	
<i>Flow rate</i>	1.0 mL/min	
<i>Column temperature</i>	25°C	
<i>Injection volume</i>	5 µL	
<i>Gradient profile</i>	<b>Time (min)</b>	<b>% B</b>
	0	0
	2	0
	14	100
	18	10
	18.1	0
	20.5	0



SCIEX OS can deliver faster method set-up



# MS Method



Suggested IDA (Information Dependent Acquisition) conditions for routine food contaminant testing as displayed in SCIEX OS.

# Data Processing



Review your results with utmost efficiency using SCIEX OS for simultaneous quantitation and MS/MS library confirmation.

# X500R High Resolution Libraries

Download a free XIC compound list detailing a full list of mycotoxin compounds include molecular formula and accurate mass.

Download a free trial of the mycotoxin high resolution MS/MS library, containing 288 compounds.

Download a free XIC compound list detailing a full list of mycotoxin compounds include molecular formula and accurate mass.

Download a free trial of the mycotoxin high resolution MS/MS library, containing 288 compounds.

[Learn more at sciex.com/X-Mycotoxins.](http://sciex.com/X-Mycotoxins)

AB Sciex is doing business as SCIEX.

© 2015 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-07-2860-A

# Food Method

## Pesticide and Potency Testing for the Cannabis Industry

*Elevate your confidence in cannabis testing with the X500R QTOF System*

Method details and access to HR-MS/MS libraries to detect, quantify, and confirm pesticides, mycotoxins, cannabinoids and terpenes in plant edible samples using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software.



⬢ **Sample Prep**

⬇ **Step 1**

- Weight 0.1g sample

⬇ **Step 2**

- Sonicate in 9mL methanol for 5 minutes

⬇ **Step 3**

- Bring to 10mL final volume

⬇ **Step 4**

- Centrifuge and draw off 1mL for analysis.


⬢ **LC Method**

<i>Column</i>	Restek Raptor ARC-18 50 x 2.1mm, 2.7µm	
<i>Mobile Phase A</i>	0.1% formic acid, 5mM ammonium formate in water	
<i>Mobile Phase B</i>	0.1% formic acid, 5mM ammonium formate in acetonitrile	
<i>Flow rate</i>	0.4 mL/min	
<i>Column temperature</i>	40°C	
<i>Injection volume</i>	5 µL	
<i>Gradient profile</i>	<b>Time (min)</b>	<b>% B</b>
	0	30
	0.5	30
	4	95
	5	95
	5.1	30



SCIEX OS can deliver faster method set-up

# MS Method



Suggested IDA (Information Dependent Acquisition) conditions for routine food contaminant testing as displayed in SCIEX OS.

# Data Processing



Review your results with utmost efficiency using SCIEX OS for simultaneous quantitation and MS/MS library confirmation.

# X500R High Resolution Libraries

**Download the XIC List**

**Download a FREE Trial of the MS/MS Library**

Download a free XIC compound list detailing a full list of applicable compounds including molecular formula and accurate mass.

Download a free trial of our All in One high resolution MS/MS library, containing 2135 compounds.

[Learn more at sciex.com/X-All.](http://sciex.com/X-All)

AB Sciex is doing business as SCIEX.

© 2015 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-07-2859



# Perfect Balance to Elevate your Lab's Performance

*Using the X500R QTOF System and SCIEX OS Software to Identify and Quantify Food Residues*

André Schreiber<sup>1</sup>, April Quinn-Paquet<sup>2</sup>, and Jianru Stahl-Zeng<sup>3</sup>

<sup>1</sup>SCIEX Concord, Ontario (Canada), <sup>2</sup>SCIEX Framingham, Massachusetts (USA), <sup>3</sup>SCIEX Darmstadt (Germany)

## Overview

Here we present results using a new method to identify and quantify pesticide residues in food using the SCIEX X500R QTOF system. Samples were extracted using a QuEChERS method and analyzed by LC-HR-MS/MS. Limits of quantitation of 10 µg/kg were achieved for every compound after 10x dilution of the extract to minimize possible matrix effects.

Target compounds were automatically identified by matching retention time, accurate mass and isotope pattern of the molecular ion and MS/MS library searching using SCIEX OS software. In the same data processing step, compounds were quantified and unknown samples were flagged when a user-defined reporting level was exceeded.

## Introduction

Recent advancements in LC-MS/MS technology, including hybrid systems like quadrupole-quadrupole Time-of-Flight (QTOF), now provide the ability to perform targeted and non-targeted screening in food samples on a routine basis.<sup>3</sup>

The SCIEX X500R QTOF system is a robust, high performance high resolution MS/MS system designed for routine use providing:

- Sensitivity to easily detect compounds at maximum residue levels
- Resolving power to remove interference from complex food matrices
- Linearity to quantify over up to 3 orders of magnitude
- Mass accuracy to identify compounds following regulatory guidelines
- Confident identification using MS/MS spectra and ion ratios
- Industry leading robustness of Turbo V™ source and Curtain Gas™ interface

Full scan chromatograms are very rich in information and easily contain thousands of ions from any residue present in the sample, including the food matrix itself. Powerful software is



needed to explore the high resolution MS/MS spectra generated to get answers and results from these complex data.

The SCIEX OS software is a single platform for MS control, data processing, and reporting and provides:

- Simple software workflows that deliver reliable results
- Simultaneous identification and quantitation
- Quick data review and reporting utilizing customizable flagging and filtering of results

## Experimental

### Standards

A standard mix of 200 pesticides was used to prepare serial dilutions for quantitative analysis.

### Sample preparation

EU proficiency test samples and food samples from a local supermarket were extracted using a QuEChERS procedure following guideline EN 15662/2007. Sample extracts were diluted 10x to minimize possible matrix effects.

### LC Separation

LC separation was performed using a SCIEX ExionLC™ AC system with a Phenomenex Kinetex Biphenyl 2.6u (50 x 2.1mm) column and a fast gradient of water and methanol with 5 mM ammonium formate buffer at a flow rate of 0.5 mL/min (see Table 1 for the gradient profile). The injection volume was 5 µL.



**Table 1.** Gradient conditions used for the separation of pesticides

Step	Time (min)	A (%)	B (%)
0	0.0	90	10
1	0.5	90	10
2	2.0	70	30
3	9.0	40	60
4	11.0	20	80
5	12.0	5	95
6	15.0	5	95
7	16.0	90	10
8	20.0	90	10

### MS/MS Detection

The SCIEX X500R QTOF system with Turbo V™ source and Electrospray Ionization (ESI) was used.

Mass calibration was achieved using the integrated calibrant delivery system (CDS) with the TwinSprayer probe (dual ESI needle).

High resolution data were acquired using an IDA method consisting of a TOF-MS survey (100-1000 Da for 100 msec) and up to 20 dependent MS/MS scans (50-1000 Da for 35 msec). MS/MS fragmentation was achieved using CE of 35 V with a collision energy spread (CES) of ±15 V.

Dynamic background subtraction (DBS) was activated for best MS/MS coverage, and no inclusion list was used to also allow retrospective unknown identification without the need for a second injection to acquire MS/MS data.

### Data Acquisition and Processing

All data were acquired and processed using SCIEX OS software version 1.0, which showcases a thoughtfully designed user interface that is fast to learn and delivers improved lab productivity.


## Results and Discussion

### X500R Performance Characteristics

Resolution > 20,000 (at full width half height) and mass accuracy <5 ppm are often sufficient to separate the analytes of interest from interfering matrices and, thus, are identified as the set requirements for compound identification in various guidelines.<sup>1,2</sup>

The X500R QTOF system utilizes N-optics design to maximize resolution while maintaining benchtop design and a minimized footprint. Its resolving power increases with mass range providing ~30000 to 40000 for the typical molecular weight range of pesticides.

The 4 mm orifice leading into the TOF accelerator delivers resolution without compromise in sensitivity. The sensitivity of the X500R QTOF system is comparable to a SCIEX QTRAP® 5500 system operated in MRM mode, allowing extract dilution to minimize ion suppression while detecting easily at 10 µg/kg levels (Figure 1).




**Figure 1.** Sensitivity and resolution of different pesticides, left: XIC of the molecular ion of each compound ± 5 mDa at 1 ng/mL (Emamectin at 10 ng/mL), right: TOF-MS spectrum of molecular ion with achieved resolution value (average of seven X500R QTOF systems)

The X500R QTOF system achieves stable mass accuracy of less than 2 ppm by using a heated TOF configuration, with 6 heater drones throughout the TOF path to maintain mass accuracy and robustness. In addition, the integrated CDS with the TwinSprayer probe provides an independent calibrant delivery path for reliable auto-calibration. The CDS setup maintains mass accuracy over long periods of time by automatically calibrating in batch mode (it is recommended to infuse a calibrant standard every hour or two).

Furthermore, the X500R QTOF's mass accuracy is supplemented by legendary dynamic transmission control and dynamic background calibration, introduced in 2010 with the TripleTOF® system and optimized over time.

Figure 2 shows an example of mass accuracy for a selected pesticide detected over a wide concentration range. Paclobutrazol was quantified from 0.1 to 1,000 ng/mL with good linearity ( $r^2 = 0.9993$ ). Excellent mass accuracy was achieved (-0.2 to 0.91 ppm) at all levels, even at the highest concentration of 10,000 ng/mL which was above the upper limit of quantitation for this analyte.




**Figure 2.** Detection of Paclobutrazol from 0.1 to 10,000 ng/mL with good linearity (0.1 to 1,000 ng/mL) and mass errors of < 1 ppm even at the highest concentration above the upper limit of quantitation

Despite the high selectivity of high resolution MS detection, there is a risk of false positive findings due to interfering isomers and matrix signals. As a result food testing guidelines require the detection of the “molecular ion” and “at least one fragment ion”, and for “a higher degree of confidence in identification, further evidence may be gained from additional mass spectrometric information. For example, evaluation of full scan spectra, isotope pattern, adduct ions, additional accurate mass fragment ions... (in MS/MS)”<sup>2</sup>

The example shown in Figure 3 highlights the need of fragment ion detection to confidently differentiate between isomers.

The pesticides Prometon and Terbumeton have identical molecular formulae ( $C_{10}H_{19}N_5O$ ) and as a result the identical molecular ion and isotope pattern. The retention time difference of less than 0.1 min, due to highly similar structures, is not sufficient to differentiate both pesticides.



**Figure 3.** Confident identification of isomers Prometon and Terbumeton using characteristic MS/MS fragment ions and MS/MS library searching

However, the two compounds have unique and characteristic fragment ions,  $C_7H_{14}N_5O^+$  and  $C_6H_{12}N_5O^+$ , respectively, which can be used for identification. Molecular and fragment ions have been measured with good mass accuracy of < 5 ppm and less < 1 mDa, respectively.

### Processing Workflow for Targeted Identification and Quantitation in SCIEX OS Software

Extracted Ion Chromatograms (XIC) of all target analytes are generated based on user input (chemical formula and expected retention time). MS and MS/MS information is automatically evaluated if an XIC signal is detected and compounds are identified by matching retention time, accurate mass and isotope pattern of the molecular ion and MS/MS library searching. Qualitative rules are defined in the processing method and can be used for results review and filtering (Figures 4a and b).

Name	Chemical Formula	Adducts	Precursor Mass (Da)	Fragment Mass (Da)	MS Width (Da)	Retention Time (min)	IS Name
2,4-D	C <sub>8</sub> H <sub>10</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	226.07327	152	0.00	16.25	
2,4-DB	C <sub>8</sub> H <sub>10</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	226.07468	152	0.00	16.00	
2,4-DB	C <sub>8</sub> H <sub>10</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	249.09799	152	0.00	15.38	
2,4-DB	C <sub>8</sub> H <sub>10</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	268.10623	152	0.00	15.68	
2-Hydroxyurea	C <sub>4</sub> H <sub>8</sub> N <sub>2</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	226.10719	152	0.00	8.30	
2-Keto-urea	C <sub>4</sub> H <sub>6</sub> N <sub>2</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	204.09173	152	0.00	15.59	
2,5-Dimethylurea	C <sub>6</sub> H <sub>12</sub> N <sub>2</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	246.11796	152	0.00	8.26	
ACEA	C <sub>8</sub> H <sub>10</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	187.02565	152	0.00	10.79	
2-Hydroxy-urea	C <sub>4</sub> H <sub>8</sub> N <sub>2</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	408.12423	152	0.00	8.79	
2-Hydroxy-urea	C <sub>4</sub> H <sub>8</sub> N <sub>2</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	408.12423	152	0.00	8.79	
2-Hydroxy-urea	C <sub>4</sub> H <sub>8</sub> N <sub>2</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	272.09449	152	0.00	14.83	
2-Hydroxy-urea	C <sub>4</sub> H <sub>8</sub> N <sub>2</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	278.10826	152	0.00	11.17	
2-Chloro-urea	C <sub>4</sub> H <sub>7</sub> N <sub>2</sub> O <sub>2</sub>	[M+H] <sup>+</sup>	207.03387	152	0.00	8.66	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	144.02918	152	0.00	8.49	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	109.22736	152	0.00	11.03	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	223.07761	152	0.00	15.2	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	270.22219	152	0.00	12.23	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	270.22219	152	0.00	8.79	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	266.03747	152	0.00	10.29	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	270.22219	152	0.00	10.22	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	270.22219	152	0.00	10.23	
Acetone	C <sub>3</sub> H <sub>6</sub> O	[M+H] <sup>+</sup>	400.12262	152	0.00	15.23	

**Figure 4a.** Method editor in SCIEX OS software, user input for target compounds including chemical formula to calculate precursor ion mass and expected retention time



**Figure 4b.** Method editor in SCIEX OS software, user input for qualitative rules (traffic lights) to enable easy results review and filtering

In the same data processing step standard calibration lines are generated to automatically calculate concentrations in unknown samples (Figure 2).

### Results of EU Proficiency Test Samples

Two samples of an EU proficiency test for pesticides and fruits and vegetables were extracted and analyzed for pesticides. Results are listed in Table 2. Retention time errors were less than 0.1 min and mass errors were between -1.20 and 1.17 ppm and were well below the required 5 ppm (SANTE/11945/2015).

Concentrations were assigned for pesticides present in the SCIEX iDQuant™ standards kit for pesticide analysis.

**Table 2.** Pesticides identified and quantified in two EU proficiency test (EUPT) samples based on matching retention time (RT), accurate mass and isotope pattern and MS/MS library searching

Pesticide	RT error (min)	Mass error (ppm)	Isotope ratio error	MS/MS FIT (%)	Conc. (µg/kg)
<b>EUPT 1</b>					
<i>Acetamiprid</i>	0.00	0.09	2.2	100.0	449
<i>Acrinathrin</i>	0.00	0.61	1.0	98.9	-
<i>Buprofezin</i>	0.01	0.32	1.1	100.0	204
<i>Chlorpyrifos</i>	0.00	-0.78	3.3	95.2	-
<i>Cypermethrin</i>	0.01	-0.27	4.9	99.2	-
<i>Cyprodinil</i>	0.01	-0.17	1.1	100.0	374
<i>Diazinon</i>	0.00	-0.20	1.7	100.0	-
<i>Difenoconazole</i>	0.00	0.22	1.8	100.0	1092
<i>Fenamiphos</i>	0.00	-1.74	1.3	99.9	-
<i>Fenamiphos-sulfone</i>	0.00	-0.26	1.7	100.0	-
<i>Fenamiphos-sulfoxide</i>	0.00	-0.94	1.3	97.1	-
<i>Fenhexamid</i>	0.02	0.16	0.6	100.0	871
<i>Fludioxonil (-)</i>	0.01	-0.69	0.8	99.6	236
<i>lambda-Cyhalothrin</i>	0.00	0.42	2.4	99.0	-
<i>Methoxyfenozide</i>	0.02	0.63	12.2	100.0	94.0
<i>Pirimicarb</i>	0.02	-0.37	0.3	100.0	478
<i>Pyridaben</i>	0.01	0.41	3.1	100.0	1063
<i>Spinosyn A</i>	0.01	-0.24	3.3	100.0	366
<i>Spinosyn D</i>	0.01	1.17	13.3	N/A	57.4
<i>Tetraconazole</i>	0.01	-0.36	9.3	100.0	111



Table 2. cont. (sample 2)

Pesticide	RT error (min)	Mass error (ppm)	Isotope ratio error	MS/MS FIT (%)	Conc. (µg/kg)
<b>EUPT 2</b>					
Atrazine	0.00	0.12	7.3	100.0	
Cadusafos	0.00	-1.20	2.3	99.2	
Carbetamide	0.02	-1.02	16.3	100.0	
Demeton-S-methyl-sulfone	0.00	0.21	0.4	99.7	
Ethoprophos	0.00	-0.47	1.7	98.7	
Fenprovidin	0.00	-0.34	2.2	100.0	
Fipronil (-)	0.00	0.20	7.3	100.0	
Flubendiamide (-)	0.00	0.11	8.9	0.0	
Fluometuron	0.01	-0.03	0.9	99.9	
Fuberidazole	0.02	-0.56	1.3	99.7	
Furathiocarb	0.01	-0.31	2.3	100.0	
Metosulam	0.00	-0.42	1.7	100.0	
Prosulfocarb	0.00	-0.54	1.2	100.0	
Secbumeton	0.00	0.06	1.6	100.0	
Spiromesifen	0.01	-0.84	5.9	99.0	

(-): identified in negative polarity

Figures 5a and 5b show screenshots of the result table used for pesticide identification in proficiency test samples.



Figure 5a. Pesticides identified in proficiency test sample 1 in positive polarity based on matching retention time, accurate mass, isotope pattern and MS/MS library searching (note: Fludioxonil was identified in negative polarity)



Figure 5b. Pesticides identified in proficiency test sample 2 in positive polarity based on matching retention time, accurate mass, isotope pattern and MS/MS library searching (note: Fipronil and Flubendiamide were identified in negative polarity)

No false positive results were reported. MS/MS data and mass spectral library searching were crucial to differentiate and correctly identify structural isomers. Library searching results were reported as FIT and in all cases were above 90%.

The pesticide Flubendiamide was not present in our MS/MS libraries. Here the built-in 'Fragments Tool' of SCIEX OS was used to compare the structure of the suspected compound with the high resolution MS/MS spectrum. All measured fragment ions matched the theoretical fragmentation pathway, resulting in a tentative identification of Flubendiamide.




Figure 6. Tentative identification of Flubendiamide based on a comparison of the HR-MS/MS spectrum with the theoretical fragmentation pathway



## Results of Store-bought Samples

Fruit and vegetable samples obtained from a local supermarket were extracted and tested for pesticide residues. Results above 10 µg/kg are listed in Table 3.

**Table 3.** Pesticides identified and quantified in store-bought fruit and vegetable samples based on matching retention time (RT), accurate mass and isotope pattern and MS/MS library searching

Sample / Pesticide	RT error (min)	Mass error (ppm)	Isotope ratio error	MS/MS FIT (%)	Conc. (µg/kg)
<b>Banana</b>					
<i>Buprofezin</i>	0.01	0.32	3.5	100.0	341
<i>Imazalil</i>	0.02	0.79	15.1	91.5	565
<i>Thiabendazole</i>	0.01	-1.51	13.9	97.6	444
<b>Blueberry</b>			n.d.		
<b>Carrot</b>			n.d.		
<b>Grapes</b>					
<i>Boscalid</i>	0.01	-0.80	8.8	97.2	115
<i>Buprofezin</i>	0.01	0.22	7.3	99.6	17.3
<i>Cyprodinil</i>	0.01	-0.87	3.3	94.8	412
<i>Imidacloprid</i>	0.01	-0.58	14.6	96.1	82.5
<i>Pyraclostrobin</i>	0.00	-1.31	4.8	100.0	46.7
<b>Lemon</b>					
<i>Imazalil</i>	0.02	0.74	7.3	94.7	1080
<i>Pyrimethanil</i>	0.01	-0.77	1.0	99.2	164
<i>Pyriproxyfen</i>	0.01	0.43	11.4	95.3	31.6
<b>Organic banana</b>					
<i>Spinosyn D</i>	0.00	2.33	19.8	100.0	12.6
<b>Organic strawberry</b>					
<i>Spinosyn A</i>	0.01	0.55	9.1	100.0	13.9
<i>Spinosyn D</i>	0.01	1.63	6.0	99.4	33.3
<b>Spinach</b>			n.d.		
<b>Strawberry</b>					
<i>Acetamiprid</i>	0.08	-0.35	6.5	98.7	19.2

**Table 3.** cont.

<i>Boscalid</i>	0.00	-0.49	4.9	99.3	161
<i>Myclobutanil</i>	0.00	-0.31	13.9	100.0	85.0
<i>Pyraclostrobin</i>	0.00	1.33	16.3	99.0	40.5
<i>Pyrimethanil</i>	0.00	0.32	4.7	97.3	391
<b>Tomato (n.d.)</b>			n.d.		

n.d.: no pesticide detected

## Summary

A new method to identify and quantify pesticide residues in food samples was developed using the SCIEX X500R QTOF system. Qualitative and quantitative data processing was performed in SCIEX OS software.

The method was successfully applied to EU proficiency test samples and store-bought fruit and vegetable samples. Samples were extracted using a QuEChERS procedure and analyzed using LC-HR-MS/MS. Limits of quantitation of 10 µg/kg were achieved for all compounds after 10x dilution the extracts to minimize possible matrix effects.

Pesticides were automatically identified by matching retention time, accurate mass and isotope pattern of the molecular ion and MS/MS library searching using SCIEX OS software. In the same data processing step compounds were quantified and unknown samples were flagged when a user-defined reporting level was exceeded.

## Acknowledgement

The authors thank Amadeo Fernandez-Alba (EURL) Almeria, Spain for providing EUPT samples.

## References

- <sup>1</sup> EU Commission Decision 'concerning the performance of analytical methods and the interpretation of results' #2002/657/EC
- <sup>2</sup> EU Commission Guidance Document: 'on analytical quality control and method validation procedures for pesticides residues analysis in food and feed' #SANTE/11945/2015
- <sup>3</sup> André Schreiber et al.: 'Using the X500R QTOF System and SCIEX OS Software to Quickly Identify Unknowns in Food Samples' Application Note SCIEX (2016) # RUO-MKT-02-3761-A

# Perfect Balance to Elevate your Lab's Performance

*Using the X500R QTOF System and SCIEX OS Software to Quickly Identify Unknowns in Food Samples*

André Schreiber<sup>1</sup>, Yuji Aita<sup>2</sup>, and Jianru Stahl-Zeng<sup>3</sup>

<sup>1</sup>SCIEX Concord, Ontario (Canada), <sup>2</sup>SCIEX Tokyo (Japan), <sup>3</sup>SCIEX Darmstadt (Germany)

## Overview

Here we present results using a new method to identify unexpected chemical residues and contaminants in food using the SCIEX X500R QTOF system. Samples were extracted using a QuEChERS method and analyzed by LC-HR-MS/MS.

Unknown compounds were automatically identified by using a non-target peak finding algorithm followed by sample-control-comparison to separate matrix and sample specific signals from true contaminations. TOF-MS and MS/MS data for ions of interest were automatically processed using formula finding and searched against mass spectral libraries and online databases, such as ChemSpider, for identification. The SCIEX OS software offers an easy to use and intuitive workflow to tentatively identify unexpected chemicals in food.

## Introduction

Hybrid LC-MS/MS systems like quadrupole-quadrupole Time-of-Flight (QTOF) provide the ability to perform targeted and non-targeted screening in food samples on a routine basis.

The SCIEX X500R QTOF system is a robust, high performance high resolution MS/MS system designed for routine use providing:

- Sensitivity to easily detect compounds at relevant concentrations
- Resolving power to remove interference from complex food matrices
- Linearity over up to 3 orders of magnitude to identify compounds at different concentration levels
- Mass accuracy to identify compounds following regulatory guidelines
- Confident identification using MS/MS spectra and ion ratios
- Industry leading robustness of Turbo V™ source and Curtain Gas™ interface

Full scan chromatograms are very rich in information and easily contain thousands of ions from any chemical present in the sample, including the food matrix itself. Powerful software is



needed to explore the high resolution MS/MS spectra generated to get answers and results from these complex data.

The SCIEX OS software is a single platform for MS control, data processing and reporting, and provides:

- Simple software workflows that deliver reliable results
- Automated identification of unknowns
- Quick data review and reporting utilizing customizable flagging and filtering of results

## Experimental

### Sample preparation

Food samples from a local supermarket were extracted using a QuEChERS procedure following guideline EN 15662/2007. Sample extracts were diluted 10x to minimize possible matrix effects.

### LC Separation

LC separation was performed using a SCIEX ExionLC™ AC system with a Phenomenex Kinetex Biphenyl 2.6u (50 x 2.1mm) column and a fast gradient of water and methanol with 5 mM ammonium formate buffer at a flow rate of 0.5 mL/min (see Table 1 for the gradient profile).

The injection volume was 5 µL.

**Table 1.** Gradient conditions used for unknown screening

Step	Time (min)	A (%)	B (%)
0	0.0	90	10
1	0.5	90	10
2	2.0	70	30
3	9.0	40	60
4	11.0	20	80
5	12.0	5	95
6	15.0	5	95
7	16.0	90	10
8	20.0	90	10

### MS/MS Detection

The SCIEX X500R QTOF system with Turbo V™ source and Electrospray Ionization (ESI) was used.

Mass calibration was achieved using the integrated calibrant delivery system (CDS) with the TwinSprayer probe (dual ESI needle).

High resolution data were acquired using an IDA method consisting of a TOF-MS survey (100-1000 Da for 100 msec) and up to 20 dependent MS/MS scans (50-1000 Da for 35 msec). MS/MS fragmentation was achieved using CE of 35 V with a collision energy spread (CES) of  $\pm 15$  V.

Dynamic background subtraction (DBS) was activated for best MS/MS coverage, and no inclusion list was used to also allow retrospective unknown identification without the need for a second injection to acquire MS/MS data.

### Data Acquisition and Processing

All data were acquired and processed using SCIEX OS software version 1.0, which showcases a thoughtfully designed user interface that is fast to learn and delivers improved lab productivity.


## Results and Discussion

### X500R Performance Characteristics

Resolution > 20,000 (at full width half height) and mass accuracy <5 ppm are often sufficient to separate the analytes of interest from interfering matrices and, thus, are identified as the set requirements for compound identification in various guidelines.<sup>1,2</sup>

The X500R QTOF system utilizes N-optics design to maximize resolution while maintaining benchtop design and a minimized footprint (Figure 1). Its resolving power increases with mass range providing ~30000 to 40000 for the typical molecular weight range of pesticides.<sup>3</sup>

The 4 mm orifice leading into the TOF accelerator delivers resolution without compromise in sensitivity. The sensitivity of the X500R QTOF system is comparable to a SCIEX QTRAP® 5500 system operated in MRM mode, allowing extract dilution to minimize ion suppression while detecting easily at 10 µg/kg levels.<sup>3</sup>



**Figure 1.** N-optics design of the X500R QTOF system to maximize resolution while maintaining benchtop design and a minimized footprint, 6 heater drones are integrated into the TOF path to maintain mass accuracy and robustness


The X500R QTOF system achieves stable mass accuracy of less than 2 ppm by using a heated TOF configuration, with 6 heater drones throughout the TOF path to maintain mass accuracy and robustness. In addition, the integrated CDS with the TwinSprayer probe provides an independent calibrant

delivery path for reliable auto-calibration. The CDS setup maintains mass accuracy over long periods of time by automatically calibrating in batch mode (it is recommended to infuse a calibrant standard every hour or two).


Furthermore, the X500R QTOF's mass accuracy is supplemented by legendary dynamic transmission control and dynamic background calibration, introduced in 2010 with the TripleTOF<sup>®</sup> system and optimized over time.

While accurate mass measurement of the molecular ion is important for empirical formula finding, this is not the only information available. Combining all available accurate mass MS and MS/MS information is crucial to minimize the list of potential formulae. Figures 2, 3 and Table 2 illustrate that the number of formulae can be reduced from over 200 to a single match by not only using the accurate mass of the molecular ion but also including the isotope pattern and MS/MS matching in the formula-finding algorithm.

Using the combined scoring of MS and MS/MS matches, SCIEX OS lists the most likely chemical formula at the top of results table. Also, SCIEX OS downloads a ChemSpider hit count for each calculated formula which further assists in identifying the correct result (Figure 2).



**Figure 2.** TOF-MS and MS/MS spectra used for empirical formula finding, results are ranked by a combined score using MS and MS/MS information, and when combined with the ChemSpider hit count, can be used to quickly find the correct match



**Figure 3.** Number of matching molecular formulae depending on the information and mass accuracy used for empirical formula finding (elements allowed C<sub>49</sub>H<sub>75</sub>Br<sub>3</sub>Cl<sub>5</sub>F<sub>3</sub>I<sub>3</sub>N<sub>10</sub>O<sub>10</sub>PS<sub>3</sub>)


**Table 2.** Ranking of matching formulae using MS and MS/MS information collected for Trifloxystrobin, the MS rank combines mass accuracy and isotope pattern matching and the MS/MS rank combines mass accuracy and number of ions (n)

Hit	Formula	MS Rank	ppm	MS/MS Rank	ppm (n=11)
1	C <sub>20</sub> H <sub>19</sub> F <sub>3</sub> N <sub>2</sub> O <sub>4</sub>	2	0.3	2	2.0
2	C <sub>21</sub> H <sub>15</sub> F <sub>3</sub> N <sub>6</sub>	9	-2.9	4	3.0
3	C <sub>18</sub> H <sub>16</sub> N <sub>6</sub> O <sub>4</sub>	4	0.9	6	4.8
4	C <sub>15</sub> H <sub>17</sub> FN <sub>8</sub> O <sub>5</sub>	11	-1.9	5	4.8
5	C <sub>16</sub> H <sub>13</sub> FN <sub>12</sub> O	7	-5.2	10	9.0
6	C <sub>14</sub> H <sub>20</sub> F <sub>3</sub> N <sub>6</sub> O <sub>3</sub> P	22	2.8	1	2.0
7	C <sub>16</sub> H <sub>21</sub> N <sub>6</sub> O <sub>5</sub> P	7	-3.1	11	9.4
8	C <sub>23</sub> H <sub>18</sub> F <sub>2</sub> N <sub>2</sub> O <sub>3</sub>	9	3.1	14	9.4
9	C <sub>21</sub> H <sub>23</sub> F <sub>2</sub> O <sub>4</sub> P	1	-0.9	24	22.1
10	C <sub>19</sub> H <sub>21</sub> FN <sub>2</sub> O <sub>7</sub>	16	-8.4	12	9.4

In addition to more efficient formula finding, MS/MS spectra are also needed for structural elucidation. Without MS/MS spectra it is impossible to conclude a correct structure from a molecular formula alone.

The example shown in Figure 4 highlights the need of fragment ion detection to confidently differentiate between isomers Prometon and Terbumeton.





**Figure 4.** Confident identification of isomers Prometon and Terbumeton using characteristic MS/MS fragment ions and MS/MS library searching

### Processing Workflow for Unknown Identification in SCIEX OS Software

Extracted Ion Chromatograms (XIC) are generated using a non-target peak finding algorithm. No masses or retention times are provided to find chromatographic features. Sample-control-comparison is used to separate matrix and sample-specific signals from true contaminations.


High resolution TOF-MS and MS/MS data of ions of interest are automatically processed using:

- MS/MS library searching to identify compounds already present in existing libraries
- Empirical formula finding based on TOF-MS and MS/MS
- ChemSpider searching
- Comparison of structures retrieved from ChemSpider against the acquired HR-MS/MS spectra

The method editor in SCIEX OS software to setup parameters and criteria for unknown identification is shown in Figure 5a-c.



**Figure 5a.** Method editor in SCIEX OS software for unknown identification, selection of sample and control-sample for non-target peak finding



**Figure 5b.** Method editor in SCIEX OS software for unknown identification, configuration of library search parameters



**Figure 5c.** Method editor in SCIEX OS software for unknown identification, configuration of formula finding options


SCIEX offers true HR-MS/MS spectral libraries for over 2500 compounds, including pesticides, veterinary drugs, toxins, fluorochemicals, pharmaceuticals, and illicit drugs.

### Results of Unknown Identification

Two samples of bell pepper, including an organic pepper, were extracted and analyzed using the developed LC-HR-MS/MS method in positive and negative polarity. Both samples were processed using the described non-target workflow.

A total of 2358 (positive polarity) and 1563 (negative polarity) chromatographic features were identified using the non-target peak finding algorithm. Less than 50 features were found to be characteristic for the contaminated bell pepper after sample-control-comparison using an area ratio of 10.

Results can be sorted and filtered for easy data review after performing sample-control-comparison. Library searching and formula finding results and scores are listed in the result table. More details and a visual display of XIC, TOF-MS and MS/MS for both samples can be found in peak review (Figure 6).





**Figure 6.** Results display after non-target screening, library searching and formula finding results are displayed in the table (top) and chromatograms and spectra with result details can be reviewed (bottom)

Formula finding results are displayed below the TOF-MS spectrum in the peak review window. Results are automatically ranked by mass accuracy (MS and MS/MS) and the matching of the isotope pattern. In addition the ChemSpider hit count is listed to quickly identify the correct match. The formulae can be searched against ChemSpider. Structural information from ChemSpider will be automatically compared against the acquired MS/MS spectrum to provide feedback for a quick identification.


Examples of tentatively identified pesticides in the bell pepper sample are shown in Figures 7, 8 and 9.




**Figure 7a.** Review of XIC of m/z 239.150 at RT 5.3 min and spectra with a found formula of  $C_{11}H_{18}N_4O_2$




**Figure 7b.** The ChemSpider search and automatic elucidation of the MS/MS spectrum led to the tentative identification of Pirimicarb (top) and also of its metabolite Desmethyl-pirimicarb (bottom), both compounds were confirmed by MS/MS library searching




**Figure 8a.** Review of XIC of m/z 226.134 at RT 8.2 min and spectra with a found formula of  $C_{14}H_{15}N_3$ , although ranked second based on mass accuracy the high ChemSpider hit count revealed the correct match



**Figure 8b.** The ChemSpider search and automatic elucidation of the MS/MS spectrum led to the tentative identification of Cyprodinil, this compound was confirmed by MS/MS library searching



**Figure 9a.** Results display after non-target screening of the negative polarity data, review of XIC of m/z 367.203 at RT 6.7 min and spectra with a found formula of  $C_{22}H_{28}N_2O_3$



**Figure 9b.** The ChemSpider search and automatic elucidation of the MS/MS spectrum led to the tentative identification of Methoxyfenozide

© 2016 AB Sciex. For Research Use Only. Not for use in diagnostic procedures.

The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Publication number: RUO-MKT-02-3761-A

## Summary

A new method to identify unexpected chemical residues and contaminants in food samples was developed using the SCIEX X500R QTOF system. Store-bought food samples were extracted using a QuEChERS procedure and analyzed by LC-HR-MS/MS.

Data processing was performed in SCIEX OS software. The processing workflow consists of peak finding using a non-target algorithm (no masses or retention times were provided to find chromatographic features). Automatic sample-control-comparison was used to separate matrix and sample specific signals from true contaminations. In a final step, tools such as empirical formula finding, MS/MS library searching and online database searching was used for identification.

The method was successfully applied to tentatively identify pesticide residues in vegetable samples.

## References

- 1 EU Commission Decision 'concerning the performance of analytical methods and the interpretation of results' #2002/657/EC
- 2 EU Commission Guidance Document: 'on analytical quality control and method validation procedures for pesticides residues analysis in food and feed' #SANTE/11945/2015
- 3 André Schreiber et al.: 'Using the X500R QTOF System and SCIEX OS Software to Identify and Quantify Food Residues' Application Note SCIEX (2016) # RUO-MKT-02-3760-A

# X500R QTOF WATER TESTING



# Environmental Method

## Pesticide analysis in water

*Elevate your environmental testing with the X500R QTOF System*

Method details and access to HR-MS/MS libraries to detect, quantify, and confirm pesticides in water samples using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software.




### Sample Prep

- Step 1
- Direct aqueous injection



SCIEX OS can deliver faster method set-up



SCIEX OS delivers enhanced data exploration of your acquired TOF MS and TOF MS/MS data

### LC Method

Column	Phenomenex Kinetex Biphenyl, 100 x 2.1 mm, 2.6 um	
Mobile Phase A	5 mM ammonium formate in water	
Mobile Phase B	5 mM ammonium formate in methanol	
Flow rate	0.5 mL/min	
Column temperature	30°C	
Injection volume	100 uL	
Gradient profile	<b>Time (min)</b>	<b>% B</b>
	0	0
	10	90
	13	90
	13.1	10
	15	10



# Environmental Method

## PPCP analysis in water

*Elevate your environmental testing with the X500R QTOF System*

Method details and access to HR-MS/MS libraries to detect, quantify, and confirm pharmaceuticals and personal care products in water samples using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software.

### Sample Prep

- Step 1 • Direct aqueous injection




SCIEX OS can deliver faster method set-up



### LC Method

Column	Phenomenex Kinetex Biphenyl, 100 x 2.1 mm, 2.6 µm column	
Mobile Phase A	0.1% formic acid in water	
Mobile Phase B	0.1% formic acid in methanol	
Flow rate	0.6 mL/min	
Column temperature	30°C	
Injection volume	100 µL	
Gradient profile	<b>Time (min)</b>	<b>% B</b>
	0	2
	1	2
	7	65
	7.1	100
	9	100
	9.1	2
	12	2



SCIEX OS delivers enhanced data exploration of your acquired TOF MS and TOF MS/MS data



# MS Method

**Method Overview**  
 Device: X500 QTOF  
 Ion Source: TurboSpray

Method duration: 20 min | Total scan time: 0.760605 sec  
 Estimated cycles: 1577

**Source and Gas Parameters**  
 Ion source gas 1: 50 psi | Curtain gas: 30 | Temperature: 450 °C  
 Ion source gas 2: 70 psi | CAD gas: 7

**Experiment** IDA  
 Polarity: Positive V | Spray voltage: 5500 V

**TOF MS**  
 TOF start mass: 100 Da | Declustering potential: 80 V | Collision energy: 10 V  
 TOF stop mass: 1000 Da | DP spread: 0 V | CE spread: 0 V  
 Accumulation time: 0.2 sec

**IDA Criteria** Small molecule  
 Maximum candidate ions: 10 |  Dynamic background subtraction  
 Intensity threshold exceeds: 100 cps |  Exclude former candidate ions  
 For: 1 sec  
 After: 1 occurrences

**TOF MSMS**  
 Precursor ion: 1.27 Da | Declustering potential: 80 V | Collision energy: 35 V  
 TOF start mass: 50 Da | DP spread: 0 V | CE spread: 15 V  
 TOF stop mass: 1000 Da | Accumulation time: 0.05 sec

Suggested IDA (Information Dependent Acquisition) conditions for routine environmental testing as displayed in SCIEX OS.

# Data Processing

Index	Sample Name	Sample Type	Elution	Retention Time	Actual Concentration	Peak Area	Peak Height	Retention Time (min)	Retention Time (min)	Component Name	Calculated Concentration	Accuracy	Precursor Ion	Mass	Retention Time (min)
12	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
13	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
14	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
15	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
16	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
17	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
18	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
19	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
20	10	Standard	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0

Review your results with utmost efficiency using SCIEX OS for simultaneous quantitation and MS/MS library confirmation.

# X500R High Resolution Libraries

Compound	CAS Number	Formula	Elemental Weight
2,2-Dimethylpropane	50143-07-1	C5H12	72.1508
2,2-Dimethylbutane	50143-08-2	C6H14	86.1719
2,2,3,3-Tetramethylbutane	50143-09-3	C10H22	142.2078
2,2,4,4-Tetramethylpentane	50143-10-4	C11H24	156.2289
2,2,4,6,6-Pentamethylheptane	50143-11-5	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-12-6	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-13-7	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-14-8	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-15-9	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-16-0	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-17-1	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-18-2	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-19-3	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-20-4	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-21-5	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-22-6	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-23-7	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-24-8	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-25-9	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-26-0	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-27-1	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-28-2	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-29-3	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-30-4	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-31-5	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-32-6	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-33-7	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-34-8	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-35-9	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-36-0	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-37-1	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-38-2	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-39-3	C15H32	208.3440
2,2,4,6,6-Pentamethylheptane	50143-40-4	C15H32	208.3440

Download the XIC List

Download a FREE Trial of the MS/MS Library

Download a free XIC compound list detailing a full list of PPCP compounds including molecular formula and accurate mass.

Download a free trial of our All in One high resolution MS/MS library, containing 2135 compounds.

Learn more at [sciex.com/X-All in One](http://sciex.com/X-All in One).

AB Sciex is doing business as SCIEX.

© 2015 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-07-2864-A



# Environmental Method

## Illicit drugs analysis in water

*Elevate your environmental testing with the X500R QTOF System*

Method details and access to HR-MS/MS libraries to detect, quantify, and confirm illicit drugs in water samples using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software.




### Sample Prep

- Step 1**
  - Samples acidified to pH 2 using HCL and stored in the dark at 4 C until required
- Step 2**
  - Transfer 1mL of sample into an analytical vial



SCIEX OS can deliver faster method set-up



SCIEX OS delivers enhanced data exploration of your acquired TOF MS and TOF MS/MS data

### LC Method

Column	Phenomenex Kinetex C18, 100 x 4.6 mm, 5 um column	
Mobile Phase A	0.1% formic acid in water + 2mM ammonium formate	
Mobile Phase B	Acetonitrile	
Flow rate	0.9 mL/min	
Column temperature	30°C	
Injection volume	100 uL	
Gradient profile	<b>Time (min)</b>	<b>% B</b>
	0	2
	1	2
	7	65
	7.1	100
	9	100
	9.1	2
	12	2



# X500R QTOF FORENSIC ANALYSIS

# Forensic Identification and Quantification Workflows Delivered on a Revolutionary Designed QTOF and SCIEX OS Software

## *Igniting your routine forensic testing with the new SCIEX X500R QTOF System*

Xiang He<sup>1</sup> and Adrian Taylor,<sup>2</sup>

<sup>1</sup>SCIEX, 1201 Radio Rd, Redwood City, CA 94065, USA; <sup>2</sup>SCIEX, 71 Four Valley Drive, Concord, Ontario, L4K 4V8 Canada.

### Overview

Quadrupole Time-of-Flight (QTOF) mass spectrometry is becoming the desired technology for sensitive and selective screening workflows in a forensic toxicological setting. The technology overcomes many challenges faced when using traditional techniques and more significantly captures all information about the sample in one injection to allow for retrospectively mining the data. Using the accurate mass and mass resolution information from both TOF-MS and TOF-MS/MS acquired data allows for simultaneous highly specific targeted quantitation and non-targeted screening. Here we describe a new benchtop QTOF system with revolutionary N geometry TOF designed flight path and new, intuitive software for easy adoption of accurate mass technology to forensic testing. We demonstrate that the new hardware and software combined allow a high level of confidence for compound identification and quantification from urine samples in one seamless workflow.

### Introduction

Liquid Chromatography coupled to Tandem Mass Spectrometry (LC-MS/MS) is a widely used analytical tool for the screening of compounds and metabolites. Triple quadrupole based mass analyzers operated in Multiple Reaction Monitoring (MRM) mode have become the preferred method to routinely deliver highly selective and sensitive quantitative results, but are limited to targeted screening only.

With an increasing demand for retrospective and non-targeted analyses of forensic toxicological samples, high resolution, accurate mass and full scan mass analyzers are gaining popularity. The adoption of the technology has been restricted by more complicated to use and more expensive instrumentation compared to their nominal mass counterparts. Here we introduce a revolutionary new Quadrupole Time-of-Flight (QTOF) mass spectrometer that contains advances in engineering design to bring the high performance TOF-MS and TOF-MS/MS capabilities into a compact benchtop platform.



**Figure 1: The SCIEX ExionLC™ AC HPLC system (left), the SCIEX X500R QTOF System (middle) and SCIEX OS Software (right).**

The SCIEX X500R QTOF mass spectrometer is part of a complete workflow from the fully integrated SCIEX ExionLC™ Systems to the freshly designed SCIEX OS software; a new user interface for simultaneous identification and quantification workflows (Figure 1.)

### SCIEX X500R QTOF System

The new benchtop SCIEX X500R QTOF System with revolutionary N geometry TOF designed flight path has been engineered for simplicity, service accessibility and minimized footprint. N TOF geometry, versus V geometry, gives the same effective flight path length for ions and therefore resolution, but in a smaller overall foot print. This has been accomplished with an extra mirror in the TOF chamber without a loss in transmission (Figure 2). To maintain stable mass accuracy the system uses a simple heated TOF vacuum chamber design. This consists of 6 discreet heater drones maintaining a constant TOF chamber temperature, insulating against ambient temperature changes (Figure 2).






Figure 2. SCIEX X500R QTOF System and Technology Advances

The system has been designed to maximize robustness and uptime

- Integrated Calibrant Delivery System and Turbo V™ Source with TwinSpray probe (Figure 3), allows seamless mass accuracy auto-calibrations during long runs.
- Service Accessibility
  - Easy QJet® and Turbo pump access for fast and efficient maintenance, increasing system uptime
  - Segmented TOF vacuum chamber – allows easy access to detector while protecting sensitive accelerator.



Figure 3. Integrated Calibrant Delivery System and Turbo V™ Source with TwinSpray probe

Figure 4 shows the mass accuracy stability of the SCIEX X500R QTOF System when analyzing multiple urine samples, spiked with various concentrations of analytes, without auto-calibration, over a ten hour period. The majority of compounds are shown to be within a 1 ppm mass accuracy over this time period.




Figure 4. Mass Accuracy Stability of the SCIEX X500R QTOF System in the Analysis of Urine Samples

Figure 5 shows the resolution for both TOF-MS and TOF-MS/MS masses sampled over a seven day time period on a SCIEX X500R QTOF System.




Figure 5. Resolution of the SCIEX X500R QTOF System Over a Week's Period for Selected  $m/z$ ; both TOF-MS and TOF-MS/MS

Figure 6 shows a representative linear dynamic range of the SCIEX X500R QTOF System showing 4 orders for the Asenapine compound.

### SCIEX OS Software

SCIEX OS Software is a single software platform for LC and MS control, data processing as well as reporting.



Figure 6. Linearity of the SCIEX X500R QTOF System shown for Asenapine (0.5 ng/mL to 1000 ng/mL)

The SCIEX OS software is intuitive and logical, segregated into Acquisition, Processing and Management work spaces (Figure 7). In the Acquisition work space there are separate method editors for the LC and MS parameters as well as batch creation and queue panes. The Processing allows for simultaneous identification and quantification. The Management workspace allows the adjustment of hardware, software and user settings.




Figure 7. Home Page of SCIEX OS Software. Single Software Platform for LC/MS Control, Data Processing and Reporting.

### Acquisition

The SCIEX OS software has a simplified step by step acquisition method setup with only relevant parameters being visible. Figure 8 shows the setup for an Information Dependent Acquisition method for the analysis of small molecules and the intuitive steps that are taken to input the MS parameter values.

For a quick instrument status check, the Manual Tune guides the user through the steps to perform a quick review of the system performance, perform an auto-calibration and report out the test result prior to running a batch (Figure 9).



Figure 8. SCIEX OS Software MS Acquisition Method




Figure 9. SCIEX OS Software MS Tune Allows for Quick Instrument Status Check via Simple Step by Step Instructions

Building a batch is assisted by the smart grid design allowing copy/paste, fill down, auto increment and import/export. Figure 10 shows the batch editor and the link to the auto-calibration setup.

Once the batch has been submitted to the queue the Auto-Cal samples are inserted as shown in the Queue Manager in Figure 11. The SCIEX OS software allows for detailed instrument status including monitoring and recording of LC pressure traces as well as direct control of the individual components of the system (Figure 11).






Figure 10. SCIEX OS Software Batch Editor and Setup for Auto-Calibration



Figure 12. SCIEX OS Software Allows the Simultaneous Review of Qualitative and Quantitative Results

The SCIEX OS Software allows the user to filter the results to only show compounds that pass acceptance criteria and are detected with user defined confidence (Figure 13)



Figure 11. SCIEX OS Software Queue Manager with Inserted Auto-Cal Samples and Detailed Instrument Status Panel

### Processing- Analytics

Once a results table is generated, quantitative and qualitative results can be reviewed in the same panel (Figure 12). A Traffic light system indicates the confidence of the identification based on accurate mass, retention time, isotopic pattern and library matching. Compounds calculated to be above the cutoff concentration in unknown samples are flagged. In the same workspace the peak integration, spectra and calibration lines can be displayed.



Figure 13. SCIEX OS Software Filtering Criteria

Finally results can be reported out using the *Create Report* functionality (Figure 14)




Figure 14. SCIEX OS Software Report Generation

## Acquisition Workflows on the SCIEX X500R QTOF System with SCIEX OS Software

### Information Dependent Acquisition

Information Dependent Acquisition (IDA) is a non-targeted data acquisition (Figure 15). It allows for TOF-MS quantification and provides high confidence in screening with MS/MS information that uses high selectivity through unit Q1 resolution. IDA-MS/MS provides the most interference-free fragmentation information.




Figure 15. Information Dependent Acquisition

When creating an IDA acquisition, MS and MS/MS settings are all contained in a single User Interface. Figure 16 shows the parameters used in the IDA experiments described in this technical note. In this example, one TOF-MS survey scan and up to 16 dependent TOF-MS/MS scans are triggered from the survey scan, in each data cycle.

Figure 16. SCIEX OS Software Information Dependent Acquisition Method Editor

Due to the high scanning speed (up to 100 Hz for single collision energy) on SCIEX X500R QTOF systems, almost all potential

compound targets in the samples can be confirmed with confident MS/MS library matching.

IDA-MS/MS is a non-targeted data acquisition method and the user needs to define the maximum number of candidates in each data cycle. More intense ions take higher priority within any data cycle, so for less abundant species especially in complex sample matrices, the associated MS/MS information might be missed. Therefore, an unbiased MS/MS data acquisition approach that collects MS/MS information for everything at all times (MS/MS<sup>All</sup>) will solve this potential concern.

### SWATH<sup>®</sup> Acquisition

SWATH<sup>®</sup> acquisition (Figure 17) is non-targeted and provides MS/MS information for everything in the sample, all the time. Each scan cycle in SWATH<sup>®</sup> Acquisition starts with a TOF-MS experiment. The acquisition approach therefore allows for screening and quantification from both TOF-MS and TOF-MS/MS acquired data.

Most of the existing MS/MS<sup>All</sup> techniques collect MS and MS/MS information for all ions in an alternating fashion, i.e. MS scan of all precursor ions, followed by MS/MS scan of the fragments of all precursor ions. Without precursor ion selection, such approaches suffer from insufficient sensitivity, selectivity and narrower linear range compared to IDA-MS/MS.




Figure 17. SWATH<sup>®</sup> Acquisition

SWATH<sup>®</sup> acquisition uses either a fixed or a variable Q1 isolation window, as part of a TOF-MS/MS experiment, which is stepped across the mass range of interest. Figure 18 shows the SWATH<sup>®</sup> acquisition method editor in the SCIEX OS Software, with the example of 16 looped TOF-MS/MS experiments, each with a different (variable) Q1 isolation window, that are required to cover the mass range of interest (120 to 500 *m/z*).





Row	IS	Group	Name	Precursor Ma...	Fragment Mass (Da)	XIC Width (Da)	Retention Time (min)	IS--	Experiment Index
1		6-MAM	6-MAM	136.154	211.0747	0.02	2.04		2 TOF MSMS of 136.2 (40 - 500)
2		7-Ameoclonazepam	7-Ameoclonazepam	286.074	222.1025	0.02	2.67		3 TOF MSMS of 286.1 (40 - 500)
3		7-Hydroxytrazepam	7-Hydroxytrazepam	4.9.223	180.0884	0.02	2.86		4 TOF MSMS of 43.2 (40 - 500)
4		Acetyl Fentanyl	Acetyl Fentanyl	313.212	105.07	0.02	3.08		5 TOF MSMS of 313.2 (40 - 500)
5		Alpha-Hydroxyflutazepam	Alpha-Hydroxyflutazepam	325.087	297.0665	0.02	4.19		6 TOF MSMS of 325.1 (40 - 500)
6		Alpha-Hydroxyflutazepam	Alpha-Hydroxyflutazepam	342.08	203.0777	0.02	4.51		7 TOF MSMS of 342.1 (40 - 500)
7		Alpha-Hydroxyflutazepam	Alpha-Hydroxyflutazepam	339.046	331.0272	0.02	3.94		8 TOF MSMS of 356.0 (40 - 500)
8		Alpha-PPP	Alpha-PPP	204.136	105.0699	0.02	2.04		8 TOF MSMS of 204.1 (40 - 500)
9		Alpha-PVP	Alpha-PVP	212.17	161.0963	0.02	2.61		10 TOF MSMS of 212.2 (40 - 500)
10		Alprazolam	Alprazolam	309.09	285.073	0.02	4.48		11 TOF MSMS of 309.1 (40 - 500)
11		AM-2201 4-Cl pentyl	AM-2201 4-Cl pentyl	876.171	155.0489	0.02	6.07		12 TOF MSMS of 376.2 (40 - 500)
12		Amisulpride	Amisulpride	276.18	117.0702	0.02	4.33		13 TOF MSMS of 276.2 (40 - 500)
13		Amphetamine	Amphetamine	116.112	91.0531	0.02	1.83		14 TOF MSMS of 116.1 (40 - 500)
14		Benzylcyclohexane	Benzylcyclohexane	290.139	168.1021	0.02	2.84		15 TOF MSMS of 290.1 (40 - 500)
15		Buprenorphine	Buprenorphine	378.123	131.07	0.02	1.88		16 TOF MSMS of 378.1 (40 - 500)
16		Buprenorphine	Buprenorphine	448.111	434.2938	0.02	7.67		17 TOF MSMS of 448.1 (40 - 500)
17		Carisoprodol	Carisoprodol	243.181	55.0555	0.02	1.68		18 TOF MSMS of 243.2 (40 - 500)
18		Clonazepam	Clonazepam	25.162	86.0959	0.02	5.24		19 TOF MSMS of 215.2 (40 - 500)
19		Cocaine	Cocaine	300.159	215.1139	0.02	1.87		20 TOF MSMS of 300.2 (40 - 500)
20		Cocaine	Cocaine	177.102	80.0496	0.02	1.83		21 TOF MSMS of 177.1 (40 - 500)
21		Cyclobenzaprine	Cyclobenzaprine	276.175	215.0878	0.02	4.28		22 TOF MSMS of 276.2 (40 - 500)
22		Desallylfurazepam	Desallylfurazepam	293.054	140.0284	0.02	4.42		23 TOF MSMS of 289.1 (40 - 500)
23		Desipramine	Desipramine	247.189	72.0923	0.02	4.28		24 TOF MSMS of 247.2 (40 - 500)
24		Demethylloperan	Demethylloperan	264.154	107.0483	0.02	2.60		25 TOF MSMS of 264.2 (40 - 500)
25		Dextromethorphan	Dextromethorphan	272.201	215.1438	0.02	3.84		26 TOF MSMS of 272.2 (40 - 500)
26		Diiazepam	Diiazepam	293.079	154.0424	0.02	5.43		27 TOF MSMS of 285.1 (40 - 500)

Figure 22. Automatically generating the SCIEX OS Software MRM<sup>HR</sup> Quantification Method from the SCIEX OS Software MRM<sup>HR</sup> Acquisition Method

Alternatively, if the fragment masses are not known at the time of the acquisition method creation, then the traditional MRM<sup>HR</sup> setup is still achievable by inputting the TOF start and stop masses (Figure 23).

Figure 23. Scheduled MRM<sup>HR</sup> Method Editor, MS/MS Full Scan

## Materials and Methods

### Compound list and spiking solutions

Table 1 lists the 93 compounds plus internal standards. All were procured from Cerilliant Corporation (Round Rock, TX). Two spiking solutions in methanol were prepared: one for analytes (SA) and the other for internal standards (SIS). Concentrations of all the analytes in the spiking solution SA are listed in Table 1.

Compounds in black font are in the regular panel (72 analytes) and the ones in blue font are the additional 21 analytes in the extended panel (93 analytes). Internal standards are shown in grey background.

### Calibrator preparation

Blank human urine was spiked with solution SA to prepare calibrators. Four levels of calibrators were prepared. Actual concentrations varied for each compound, however the concentration ratio between these calibrators was always (in descending order): 20:6:2:1. For instance, the four different concentrations (in descending order) for fentanyl in the calibrators were: 20, 6, 2 and 1 ng/mL, while those of gabapentin were: 1000, 300, 100 and 50 ng/mL.

### Sample preparation

- 100 µL urine sample was mixed with 25 µL IMCS Rapid Hydrolysis Buffer, 20 µL IMCSzyme and 10 µL SIS. Both IMCS Rapid Hydrolysis Buffer and IMCSzyme were acquired from IMCS (Columbia, SC). Hydrolysis time was typically between 30 and 60 min at 55°C.
- After hydrolysis was complete, 0.2 mL methanol and 0.625 mL water were added to the mixture.
- The mixture was then centrifuged at 21,000 g for 10 min.
- The supernatant was transferred to a glass vial with insert for analysis by LC-MS/MS.

### Liquid Chromatography

Liquid Chromatography analysis was performed on the SCIEX ExionLC<sup>TM</sup> AC HPLC system at 30°C. Separation was achieved using a Phenomenex Kinetex Phenyl-Hexyl column (50 × 2.1 mm, 2.6 µm, 00B-4495-E0), with a Phenomenex SecurityGuard ULTRA UHPLC Phenyl (AJ0-8774) and ULTRA holder (AJ0-9000). Mobile phase A (MPA) was ammonium formate in water. Mobile phase B (MPB) was formic acid in methanol. The LC flow rate was 0.5 mL/min and the LC run-times investigated were 8.0 and 2.0 minutes. Injection volume was 10 µL.

Table 1: List of analytes and internal standards, and their concentrations in spiking solution (for preparation of calibrators)

Compounds	(ng/mL)	Compounds	(ng/mL)	Compounds	(ng/mL)	Compounds	(ng/mL)
6-MAM	1000	Gabapentin	10000	Naloxone	5000	Pentobarbital	10000
7-Aminoclonazepam	5000	Hydrocodone	5000	Naltrexone	5000	Secobarbital	10000
7-Hydroxymitragynine	1000	Hydromorphone	5000	N-desmethyltapentadol	5000	THC-COOH	2000
Acetyl Fentanyl	200	Imipramine	5000	Norbuprenorphine	2000	6-MAM-d3	
Alpha-Hydroxyalprazolam	5000	JWH 122 5-OH pentyl	1000	Norcodeine	5000	Amphetamine-d5	
Alpha-Hydroxymidazolam	5000	JWH 19 6-OH hexyl	1000	Nordiazepam	5000	Benzoyllecgonine-d3	
Alpha-Hydroxytriazolam	5000	JWH 210 5-OH-pentyl	1000	Norfentanyl	200	Buprenorphine-d4	
Alpha-PPP	1000	JWH-018 4-OH pentyl	1000	Norhydrocodone	5000	Carisoprodol-d7	
Alpha-PVP	1000	JWH-018 pentanoic acid	1000	Normeperidine	5000	Codeine-d6	
Alprazolam	5000	JWH-073 3-OH butyl	1000	Noroxycodone	5000	Fentanyl-d5	
AM-2201 4-OH pentyl	1000	JWH-073-butanoic acid	1000	Norpropoxyphene	10000	Hydrocodone-d6	
Amitriptyline	5000	JWH-250-N-4-OH pentyl	1000	Nortriptyline	5000	Hydromorphone-d6	
Amphetamine	10000	JWH-073-butanoic acid	1000	O-Desmethyltramadol	5000	JWH 018 4-OH pentyl-d5	
Benzoyllecgonine	5000	JWH-250-N-4-OH pentyl	1000	Oxazepam	5000	JWH 019 6-OH hexyl-d5	
Buphedrone	1000	Lorazepam	5000	Oxycodone	5000	MDPV-d8	
Buprenorphine	2000	MDA	10000	Oxymorphone	5000	Meperidine-d4	
Carisoprodol	10000	MDEA	10000	PCP	2500	Mephedrone-d3	
Clomipramine	5000	MDMA	10000	Pregabalin	10000	Meprobamate-d7	
Codeine	5000	MDPV	1000	Propoxyphene	10000	Methadone-d3	
Cotinine	5000	Meperidine	5000	Protriptyline	5000	Methamphetamine-d5	
Cyclobenzaprine	5000	Mephedrone	1000	RCS4-4-OH-pentyl	1000	Methylone-d3	
Desalkylflurazepam	5000	Meprobamate	10000	Ritalinic Acid	5000	Mitragynine-d3	
Desipramine	5000	Methadone	10000	Sufentanil	200	Morphine-d6	
Desmethyldoxepin	5000	Methamphetamine	10000	Tapentadol	5000	Nordiazepam-d5	
Dextromethorphan	5000	Methedrone	1000	Temazepam	5000	Nortriptyline-d3	
Diazepam	5000	Methylone	1000	Tramadol	5000	Oxycodone-d6	
Dihydrocodeine	5000	Methylphenidate	5000	Zolpidem	5000	Oxymorphone-d3	
Doxepin	5000	Midazolam	5000	Amobarbital/pentobarbital	10000	THC-COOH-d3	
EDDP	10000	Mitragynine	1000	Butabarbital	10000	Butalbital-d5	
Fentanyl	200	Morphine	5000	Butalbital	10000	Secobarbital-d5	

Grey background: IS



## SCIEX OS Software Processing

### Identification and Quantification Results

Defining the retention time and accurate precursor and fragment mass for each analyte is performed first (Figure 24) followed by setting up the library searching parameters.




Figure 24. Defining the Retention Time, Accurate mass of Precursor and Fragment Ions

Defining the qualifying components includes setting accuracy tolerance levels for calibrants and controls as well as flagging integration discrepancies. Qualifying definitions also includes defining the identification criteria and setting the confidence levels at which mass error, error in retention time, isotope pattern and library matching scores are deemed an acceptable difference, marginal difference or unacceptable difference (Figure 25).




Figure 25. Defining the Identification and Quantification Qualifying Components in the SCIEX OS Software

## Results and Discussion

As part of evaluating the new SCIEX X500R QTOF to perform simultaneous identification and quantification of compounds from forensically related samples routinely, we investigated two LC gradients. We evaluated each methods capabilities to elute all analytes throughout the entire gradient as evenly as possible in

order to maximize triggering IDA MS/MS for all components, reduce the MRM<sup>HR</sup> concurrency for quality of data (*Scheduled MRM<sup>HR</sup>*), resolve isobaric species and alleviate ion suppression caused by co-elution of excessive number of analytes. Figure 26 shows the Extracted Ion Chromatograms (XICs) for the 8.0 minute run and Figure 27 show the XICs for the 2.0 minute run.




Figure 26. Extracted Ion Chromatograms for Analytes from a Urine Analysis using an 8.0 Minute LC Runtime




Figure 27. Extracted Ion Chromatograms for Analytes from a Urine Analysis using a 2.0 Minute LC Runtime

## Information Dependent Acquisition

With the ability to provide the most interference free fragmentation information for library searching in a non-targeted acquisition, the IDA workflow provides the highest confidence screening using MS/MS information. Figure 28 shows the multiple screening criteria that are used for identification purposes in the SCIEX OS Software's easy to understand user interface.



Component Name	Act. Conc.	Exp. RT	Height	Ret. Time	Precursor Mass	Mass Error	RT Conf.	Isotope Conf.	Library Conf.	Found At Mass	Mass Error (ppm)	Retention Time (%)	Library Hit	Library Score	Comb. Score	Isotope Ratio
6-MAM	N/A	2.05	171193	2.05	328.1544	0.0	✓	✓	✓	328.1544	1.0	0.0	6-Monoacetyl-morphine (Smart Confirmation)	2784	1.25	1.25
7-Aminoclonazepam	N/A	2.67	371193	2.67	286.0742	0.0	✓	✓	✓	286.0742	2.3	0.0	7-Aminoclonazepam	93.4	90.9	5.8
7-Hydroxythiazepam	N/A	2.90	23047	2.97	415.2217	0.0	✓	✓	✓	415.2217	0.8	0.4	7-Hydroxythiazepam	67.0	75.1	10.7
Acetyl Fentanyl	N/A	3.01	41934	3.01	323.2118	0.0	✓	✓	✓	323.2118	1.0	0.1	Acetyl Fentanyl	04.0	143.5	0.8
Alpha-Hydroxytyrosine	N/A	4.19	150463	4.19	325.0851	0.0	✓	✓	✓	325.0851	1.0	0.0	Alpha-Hydroxytyrosine	96.3	95.3	2.2
Alpha-Hydroxytryptamine	N/A	4.51	150463	4.51	342.0804	1.2	✓	✓	✓	342.0804	1.2	0.1	Alpha-Hydroxytryptamine	98.4	96.3	1.3
Alpha-Hydroxytryptolite	N/A	3.95	81889	3.95	359.0461	0.0	✓	✓	✓	359.0461	1.0	0.0	alpha-Hydroxytryptolite	91.5	92.2	2.6
Alpha-PPP	N/A	2.04	155541	2.03	204.1363	0.0	✓	✓	✓	204.1363	1.7	0.4	Alpha-PPP	95.4	89.9	1.6
Alpha-PVP	N/A	2.61	217894	2.61	232.1696	0.0	✓	✓	✓	232.1696	1.2	0.0	Alpha-PVP	99.2	98.8	1.6
Alprazolam	N/A	4.49	393119	4.49	309.0902	0.0	✓	✓	✓	309.0902	1.8	0.1	Alprazolam	82.4	79.7	8.3
Am-2201 4-Ch (M)	N/A	6.07	121905	6.07	376.1707	0.0	✓	✓	✓	376.1707	0.4	2.1	Am-2201 4-Ch (M)	100.0	99.1	2.9
Amphetamine	N/A	4.55	64850	4.52	278.1953	0.0	✓	✓	✓	278.1953	1.4	0.3	Amphetamine	98.3	98.5	1.4
Amphetamine	N/A	1.62	28842	1.64	136.1121	1.0	✓	✓	✓	136.1121	1.0	1.4	Amphetamine	100.0	91.4	10.4
Benzylhexgonine	N/A	2.16	86060	2.17	290.1347	0.0	✓	✓	✓	290.1347	1.8	0.2	Benzylhexgonine	99.4	95.7	0.7
Buprenorphine	N/A	1.97	107591	1.98	178.1226	0.0	✓	✓	✓	178.1226	0.7	0.4	Buprenorphine	99.1	97.2	1.2
Buprenorphine	N/A	3.68	111748	3.68	468.3108	0.0	✓	✓	✓	468.3112	0.7	0.1	Buprenorphine	86.3	96.0	0.9
Carisoprodol	N/A	3.68	55995	3.68	261.1809	0.0	✓	✓	✓	261.1811	1.0	0.4	Carisoprodol	99.6	95.3	1.1




Figure 28 Screening and Identification Results from an IDA Experiment

The importance of acquiring quality MS/MS data for identification purposes, and not to solely rely on the accurate mass of the precursor ion, is demonstrated in Figures 29, 30 and 31. Each figure demonstrates how, by acquiring MS/MS data, we can distinguish between structural isobaric compounds. In each example shown, isobaric compounds are barely chromatographically separated and so the presence of either or both the compounds cannot be identified by either accurate mass of the precursor ions or retention time if there is any drift in retention of the compounds. The highest confidence is gained through library MS/MS comparisons.




Figure 29. High Confidence Identification of Naloxone and 6-MAM Isobaric Compounds Gained through Library MS/MS Comparisons




Figure 30. High Confidence Identification of Buprenorphine and Mephedrone Isobaric Compounds Gained through Library MS/MS Comparisons




Figure 31. High Confidence Identification of Methylphenidate and Normeperidine Isobaric Compounds Gained through Library MS/MS Comparisons

Figures 32 and 33 show selected compound examples of XICs from the TOF-MS information acquired as part of the IDA workflow. This information can be used for quantification purposes.




Figure 32. XICs of  $\alpha$ -PVP in Urine from TOF-MS information (Urine was diluted 10-fold; 10  $\mu$ L injection)




Figure 33. XICs of Sufentanil in Urine from TOF-MS information (Urine was diluted 10-fold; 10  $\mu$ L injection)

Figure 34 shows representative calibration curves obtained from the IDA experiment.




Figure 34. Representative Calibration Curves for Selected Compounds Showing that the TOF-MS information can be used for Quantification in an IDA Workflow

## SWATH® Acquisition Results

SWATH® Acquired data can be processed in a similar way to processing IDA data for screening purposes. Again this uses multiple criteria for confidence in identification; most importantly using MS/MS library matching. Figure 35 shows a result of this from the 8.0 minute LC run which resulted in a high true positive rate of 98%.




Figure 35. Processed SWATH® Acquired Data using Multiple Identification Criteria; including MS/MS Library Matching

With traditional IDA-MS/MS, quantitation can only be performed from TOF-MS mode and not from the *in situ* sporadic TOF-MS/MS data points. In contrast, due to the continual and looped MS/MS scan function, quantification from fragment ions is achievable from SWATH® acquisition. Better selectivity from the fragment ion information (Figure 36) relative to parent ion information, allows more sensitive detection in MS/MS mode of lower concentration species in complex matrices.





Figure 36. Gains in Selectivity with the Ability to Extract Out a Specific Fragment Ion From Variable Window SWATH® Acquired Data Compared to Extracted Accurate Mass of the Precursor Ion

Figure 37 shows identification and quantification results for a synthetic drug obtained from SWATH® Acquisition using the 8.0 min LC run time. This compound was not in the original targeted list but retrospective interrogation of the data from this unknown




sample allowed for its identification without having to re-inject the sample again.




**Figure 37. Identification and Quantification Results for n-Ethylcathinone Ephedrine Metabolite Compound Analysed by SWATH® Acquisition**

The n-ethylcathinone ephedrine metabolite compound was identified based on unique fragment ions and their ratios as well as a library searching match (Figure 38). In a SWATH® acquisition experiment, not only can confirmation of the presence of compounds be made through MS/MS library matching and ion ratio calculations but because of the ability to extract out many unique fragment ions from the SWATH® acquired MS/MS data we can also determine the concentration based on quantification of either or both the precursor and fragment ions depending on which has less interferences.




**Figure 38. Extraction of Unique Fragment Ions From SWATH® Acquisition and Using Both Ion Ratio and Library Matching to Confirm Presence of n-Ethylcathinone Ephedrine metabolite in an Unknown Urine Sample**

When investigating using a 2.0 minute LC run time as part of the SWATH® acquisition, we were able to accomplish good quantification results. Figure 39 shows representative calibration curves obtained from the ultra-fast screening experiment.




**Figure 39. Representative Calibration Curves Generated from the SWATH® Acquisition using a 2.0 minute LC Runtime (n=3)**

Sensitivity examples are shown in Figures 40, 41 and 42 for selected compounds.



**Figure 40. XICs of Buprenorphine at Various Concentrations in Urine (Diluted 10-fold, 10 µL injection)**



**Figure 41. XICs of Fentanyl at Various Concentrations in Urine (Diluted 10-fold, 10 µL injection)**






Figure 42. XICs of Amitriptyline at Various Concentrations in Urine (Diluted 10-fold, 10 µL injection)

In the SWATH® Acquisition, MS/MS information is always available and so we can confirm the presence of the compound through MS/MS library matching (Figures 43 and 44) at the same time as determining how much of the compound is present.




Figure 43. Confident Identification of JWH-073 3-Butyl and Nortapentadol from SWATH® Acquisition Through Library Searching




Figure 44. Confident Identification of Amitriptyline and EDDP from SWATH® Acquisition Through Library Searching; Showing LC Separation Between Isomers was Still Achievable with this Fast Method

At the cutoff concentrations, library matching worked well with 80% of compounds yielding greater than 70% hit score (Figure 45).

Sample Name	Component No.	Actual Conc.	Height	Ret. Time	Precursor Mass	Mass Error	RT Conf.	Isotope Conf.	Library Conf.	Found At Mass	Mass Error	Ret. Time	Library Hit	Library Score	Comb. Score	Isotope Ratio Diff.
EPDCCurve_UCa1	Alprazolam	50.00	39322	1.31	325.0815		✓	✓	✓	325.0815	0.0	0.00	Alprazolam	100.0	95.1	18.6
EPDCCurve_UCa2	Alphacetylmethadol	50.00	14742	1.12	342.0804		✓	✓	✓	342.0804	0.0	0.00	Alphacetylmethadol	100.0	99.1	1.8
EPDCCurve_UCa2	Buprenorphine	20.00	20422	1.24	468.1108		✓	✓	✓	468.1107	-0.1	0.01	Buprenorphine	100.0	97.9	2.9
EPDCCurve_UCa2	Desulfinorphine	50.00	19370	1.33	289.0518		✓	✓	✓	289.0519	0.0	0.01	Desulfinorphine	100.0	98.2	1.4
EPDCCurve_UCa2	EDDP	100.00	768376	1.20	278.1938		✓	✓	✓	278.1938	0.0	0.00	EDDP	100.0	96.1	1.3
EPDCCurve_UCa2	Fentanyl	2.00	10599	1.17	317.2274		✓	✓	✓	317.2274	0.0	0.01	Fentanyl	100.0	97.3	0.7
EPDCCurve_UCa2	Hydrocodone	50.00	61526	0.92	300.1594		✓	✓	✓	300.1597	0.0	0.00	Hydrocodone	100.0	94.9	5.3
EPDCCurve_UCa2	JWH-018 4-OH	10.00	29593	1.39	358.1802		✓	✓	✓	358.1801	-0.1	0.02	JWH-018 Hydroxylated	100.0	98.0	4.6
EPDCCurve_UCa2	JWH-018 6-OH	10.00	44558	1.40	372.1938		✓	✓	✓	372.1939	0.0	0.01	JWH-018 6-Hydroxylated	100.0	98.7	4.4
EPDCCurve_UCa2	JWH-073 3-OH	10.00	15459	1.38	344.1645		✓	✓	✓	344.1648	0.0	0.01	JWH-073 3-Hydroxylated	100.0	94.3	5.1
EPDCCurve_UCa2	JWH-122 5-OH	10.00	44558	1.40	372.1938		✓	✓	✓	372.1939	0.0	0.01	JWH-122 5-Hydroxylated	100.0	96.7	4.2
EPDCCurve_UCa2	Mefenorex	100.00	15875	0.92	310.2277		✓	✓	✓	310.2277	-0.4	0.00	Mefenorex	100.0	98.0	3.8
EPDCCurve_UCa2	Mefenorex	100.00	19347	1.01	234.1489		✓	✓	✓	234.1490	0.0	0.00	Mefenorex	100.0	98.4	0.7
EPDCCurve_UCa2	Nicobutylpropylamine	20.00	27223	1.17	414.2610		✓	✓	✓	414.2611	0.0	0.00	Nicobutylpropylamine	100.0	97.5	2.3
EPDCCurve_UCa2	Nisoxetine	50.00	44884	1.35	371.0613		✓	✓	✓	371.0615	0.0	0.00	Nisoxetine	100.0	91.9	3.7
EPDCCurve_UCa2	Nisoxetine	50.00	23823	0.90	302.1387		✓	✓	✓	302.1388	-0.2	0.01	Nisoxetine	100.0	97.9	4.3
EPDCCurve_UCa2	Tamoxifen	50.00	20036	1.09	222.1052		✓	✓	✓	222.1054	0.0	0.00	Tamoxifen	100.0	95.6	1.1
EPDCCurve_UCa2	Tamoxifen	50.00	18735	1.34	301.0738		✓	✓	✓	301.0738	-0.1	0.01	Tamoxifen	100.0	98.1	1.2
EPDCCurve_UCa2	Tamoxifen	50.00	16948	1.06	264.1958		✓	✓	✓	264.1961	1.0	0.00	Tamoxifen	100.0	95.9	0.6
EPDCCurve_UCa2	Zolpidem	50.00	27738	1.13	308.1737		✓	✓	✓	308.1737	0.0	0.00	Zolpidem	100.0	98.1	1.9
EPDCCurve_UCa2	T-Amoxicillin	10.00	78113	1.12	286.0742		✓	✓	✓	286.0744	0.0	0.01	T-Amoxicillin	99.8	95.0	3.2
EPDCCurve_UCa2	Alprazolam	50.00	12384	1.32	309.0902		✓	✓	✓	309.0904	0.0	0.00	Alprazolam	99.8	95.7	4.0
EPDCCurve_UCa2	Clonidine	50.00	11040	1.16	235.0789		✓	✓	✓	235.0791	0.0	0.01	Clonidine	99.7	95.5	3.4
EPDCCurve_UCa2	Flunitrazepam	50.00	16139	0.94	312.1700		✓	✓	✓	312.1698	-0.4	0.00	Flunitrazepam	99.7	97.3	1.0
EPDCCurve_UCa2	Cyproheptadine	50.00	1763	0.41	312.1387		✓	✓	✓	312.1389	0.8	0.02	Cyproheptadine	99.7	94.2	3.5
EPDCCurve_UCa2	Amphetamine	100.00	1840	0.88	116.1123		✓	✓	✓	116.1120	-0.1	0.00	Amphetamine	99.6	97.1	4.1

Figure 45. Library Searching and Identification of Compounds in the 2.0 Minute Method at Cutoff Concentration Levels

MRM<sup>HR</sup>

MRM<sup>HR</sup> is a purely targeted data MS/MS acquisition and can be unscheduled or scheduled. The only non-targeted and therefore retrospective capability is through the TOF-MS experiment which is performed at the beginning of every scan. The power of the workflow however, is its selectivity capabilities through the accurate mass of unique fragment ions for quantification purposes. This is demonstrated in Figure 46 where MRM<sup>HR</sup> is compared to the MRM analysis, extracted at nominal mass, and the extraction of the accurate mass of the precursor ion from a TOF-MS experiment. The compound is not able to be distinguished from the high background and interferences of the nominal mass experiment and not even by the extraction of the accurate mass of the precursor ion from the full scan TOF-MS experiment. It is not until we extract out two unique accurate mass fragment ions from the MRM<sup>HR</sup> experiment that we achieve the selectivity required to detect this compound by removal of the background and interferences and increase the S/N; improving the quantification capabilities. Another example of this selectivity gain over the accurate mass of the precursor ion is demonstrated in Figure 47 where a visible improvement in S/N is gained for the analysis of buprenorphine by the MRM<sup>HR</sup> approach.




Figure 46. Increased Selectivity with MRM<sup>HR</sup>; Avoiding False Negatives. Example given is a Feed Sample Tested Positive for NP Semicarbazide




Figure 39. Scheduled MRM<sup>HR</sup> Selectivity Compared to TOF-MS; Buprenorphine (5ng/mL in urine, 10 fold dilution, 10 µL injection)

Quantification performance of the MRM<sup>HR</sup> is demonstrated in Figure 40 for the 8.0 minute LC-MS/MS method.

Scheduled MRM<sup>HR</sup>, 372.2 → 169.0644 ± 0.0100 m/z




Figure 40. MRM<sup>HR</sup> Quantification Results for JWH-122 5-OH Pentyl in urine (Urine was diluted 10-fold, 10 µL injection)

### Negative Mode Performance

Figures 41 and 42 show a couple of examples of negative mode performance of the SCIEX X500R QTOF System.



Figure 41. Negative Mode Performance of SCIEX X500R QTOF System for Analysis of Amobarbital



Figure 42. Negative Mode Performance of SCIEX X500R QTOF System for Analysis of THC-COOH

### Conclusion

The arrival of the next generation QTOF, with the launch of the SCIEX X500R QTOF System and SCIEX OS Software, brings the powerful performance capabilities of the high resolution accurate mass technology to the routine identification and quantification forensic workflows.

- Hardware
  - SCIEX ExionLC<sup>TM</sup> Systems
    - Fully controlled by SCIEX OS software
    - Improved software integration for better stability
  - SCIEX X500R QTOF System
    - N-geometry design (same effective flight path length for ions and therefore resolution than V-geometry, but in a smaller overall footprint)
    - Heated TOF path for mass accuracy stability



- Minimized footprint, engineered for simplicity and service accessibility
- Software
  - SCIEX OS Software
    - Intuitive and logical single software platform for LC control, MS control, data processing and reporting.
    - New user interface
    - Simultaneous identification and quantitation

runtime with MS/MS information always being available with this MS/MS<sup>ALL</sup> approach.

We have described the screening and quantification workflows of the SCIEX X500R QTOF System. Each workflow is straightforward to setup in the newly designed SCIEX OS Software and depending on the end users requirements we have demonstrated in this technical note the strengths of each workflow. Each provides TOF-MS and TOF-MS/MS analysis, both data being crucial in confidently identifying and quantifying forensic compounds.

- TOF-MS
- TOF-MS/MS
  - IDA
    - Non-targeted data acquisition
    - MS quantitation
    - Highest confidence screening with MS/MS information
  - MRM<sup>HR</sup>
    - Targeted data acquisition for quantitation purpose
    - Can be performed unscheduled or scheduled
  - SWATH<sup>®</sup> Acquisition (with variable windows)
    - Non-targeted data acquisition
    - MS/MS for everything all the time
    - Screening and quantitation (MS/MS)
      - Library Searching and Ion Ratio

We evaluated different LC runtime methods. The longer method aided eluting all analytes throughout the entire gradient as evenly as possible in order to maximize triggering IDA MS/MS for all components and reduce the MRM<sup>HR</sup> concurrency for quality of data (*Scheduled* MRM<sup>HR</sup>). The library searching worked well for the SWATH<sup>®</sup> Acquisition in the 2.0 minute LC

**For Research Use Only. Not for use in diagnostic procedures.**

© 2016 AB SCIEX. AB Sciex is doing business as Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-02-3786-A



**Headquarters**  
500 Old Connecticut Path | Framingham, MA 01701 USA  
Phone 508-383-7700  
[www.sciex.com](http://www.sciex.com)

**International Sales**  
For our office locations please call the division  
headquarters or refer to our website at  
[www.sciex.com/offices](http://www.sciex.com/offices)



# Forensic Method

## Forensic drug screening analysis

*Elevate your forensic testing with the X500R QTOF System*

Method details and access to HR-MS/MS libraries to screen for forensic drugs in urine samples using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software.

### Sample Prep

Step 1


- Dilute urine samples 5 times with the following diluent: 90% Mobile Phase A (5mM ammonium formate in water): 10% Mobile Phase B (Methanol with 0.05% formic acid)



SCIEX OS can deliver faster method set-up


### LC Method

Column	Phenomenex Kinetex Biphenyl, 100 x 3.0 mm, 2.6 um column	
Mobile Phase A	0.1% formic acid in water	
Mobile Phase B	0.1% formic acid in methanol	
Flow rate	0.6 mL/min	
Column temperature	30°C	
Injection volume	10 µL	
Gradient profile	<b>Time (min)</b>	<b>% B</b>
	0	2
	1	2
	7	65
	7.1	100
	9	100
	9.1	2
	12	2



SCIEX OS delivers enhanced data exploration of your acquired TOF MS and TOF MS/MS data

# MS Method



Suggested IDA (Information Dependent Acquisition) conditions for routine forensic drug screening as displayed in SCIEX OS

## Data Processing



Review your results with utmost efficiency using SCIEX OS for simultaneous quantitation and MS/MS library confirmation.

## X500R High Resolution Libraries



Download the XIC List

Download a FREE Trial of the MS/MS Library

Download a free XIC compound list detailing a full list of forensic drug compounds including molecular formula and accurate mass.

Download a free trial of the forensic high resolution MS/MS library, containing 1703 compounds.

[Learn more at sciex.com/X-Forensics.](http://sciex.com/X-Forensics)

AB Sciex is doing business as SCIEX.

© 2015 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RO-MKT-07-2866-A

# Forensic Method

## Forensic drug screening analysis

### Elevate your forensic testing with the X500R QTOF System

Method details and access to HR-MS/MS libraries to screen for forensic drugs in blood extracts using HPLC coupled with the X500R QTOF system, powered by SCIEX OS Software.



### Sample Prep

- Step 1

  - Add 0.2ml of blood into centrifuge tube
- Step 2

  - Add 0.8mL of cold acetonitrile dropwise while vortexing
- Step 3

  - Centrifuge for 10 minutes at 5000rpm
- Step 4

  - Transfer supernatant to clean glass vials and dry down.
- Step 5

  - Reconstitute in 200µL of 20% methanol in water.
- Step 6

  - Vortex for 1 minute and centrifuge for 10 minutes at 10,000xg
- Step 7

  - Transfer supernatant to clean glass HPLC vials

### LC Method

<i>Column</i>	Phenomenex Kinetex Biphenyl, 100 x 3.0 mm, 2.6 um column	
<i>Mobile Phase A</i>	0.1% formic acid in water	
<i>Mobile Phase B</i>	0.1% formic acid in methanol	
<i>Flow rate</i>	0.6 mL/min	
<i>Column temperature</i>	30°C	
<i>Injection volume</i>	10 µL	
<i>Gradient profile</i>	<b>Time (min)</b>	<b>% B</b>
	0	2
	1	2
	7	65
	7.1	100
	9	100
	9.1	2
	12	2



SCIEX OS can deliver faster method set-up



# MS Method

**Method Overview**  
 Device: X500 QTOF  
 Ion Source: Turbospray

Method duration: 6.5 min | Total scan time: 0.41118 sec  
 Estimated cycles: 948

**Source and Gas Parameters**  
 Ion source gas 1: 50 psi | Curtain gas: 30 | Temperature: 600 °C  
 Ion source gas 2: 50 psi | CAD gas: 7

**Experiment** IDA  
 Polarity: Positive V | Spray voltage: 5500 V

**TOF MS**  
 TOF start mass: 100 Da | Declustering potential: 80 V | Collision energy: 10 V  
 TOF stop mass: 1000 Da | DP spread: 0 V | CE spread: 0 V  
 Accumulation time: 0.1 sec

**IDA Criteria** Small molecule  
 Maximum candidate ions: 10 |  Dynamic background subtraction  
 Intensity threshold exceeds: 10 cps |  Exclude former candidate ions  
 For: 1 sec | After: 1 occurrences

**TOF MSMS**  
 Presursor ion: 2 Da | Declustering potential: 80 V | Collision energy: 35 V  
 TOF start mass: 40 Da | DP spread: 0 V | CE spread: 15 V  
 TOF stop mass: 1000 Da | Accumulation time: 0.025 sec

Suggested IDA (Information Dependent Acquisition) conditions for routine forensic drug screening as displayed in SCIEX OS

## Data Processing



Review your results with utmost efficiency using SCIEX OS for simultaneous quantitation and MS/MS library confirmation.

## X500R High Resolution Libraries

**Download the XIC List**

**Download a FREE Trial of the MS/MS Library**

Download a free XIC compound list detailing a full list of forensic drug compounds including molecular formula and accurate mass.

Download a free trial of the forensic high resolution MS/MS library, containing 1703 compounds.

[Learn more at sciex.com/X-Forensics.](http://sciex.com/X-Forensics)

AB Sciex is doing business as SCIEX.

© 2015 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-07-2865-A

# Your Success is Our Success

## We take it personally

As an SCIEX customer you have access to an excellent customer support organization. Wherever you are, we're there with you as a trusted partner to answer questions, provide solutions, and maximize lab productivity.

Our customer support organization has access to the latest product updates, software revisions, methods and repair procedures to make sure that you stay on top of your game.

When you have questions, we have answers.

Learn more at [sciex.com/customersupport](https://sciex.com/customersupport), or locate your local account representative at [sciex.com/contactus](https://sciex.com/contactus)



AB Sciex is doing business as SCIEX.

© 2016 AB Sciex. For research use only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of the AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

RUO-MKT-03-4305-A



#### Headquarters

500 Old Connecticut Path | Framingham, MA 01701 USA  
Phone 508-383-7700  
[sciex.com](https://sciex.com)

#### International Sales

For our office locations please call the division headquarters or refer to our website at [sciex.com/offices](https://sciex.com/offices)