


Fast and comprehensive detection and characterization of cyclic peptide metabolites using software-aided data processing tools

Ming Yao¹, Eva Duchoslav², Li Ma¹, Silvi Chacko¹, Mingshe Zhu^{1,3}

¹Pharmaceutical Candidate Characterization, Bristol-Myers Squibb, Princeton, NJ 08543, USA. ²Sciex, Concord, Ontario, L4K 4V8, CANADA. ³Current Address: MassDefect Technologies, Princeton, NJ 08540 USA

Introduction

- Cyclic peptide (CP) has better biological activity compared to their linear due to counterparts the conformational rigidity.
- Resistance to hydrolysis due to the lack of both amino and carboxyl termini.
- Some cyclic peptides can cross the cell membrane
- 9 Cyclic peptide drugs have been approved by FDA and/or EMA from 2006 to 2014



Challenges in identifying cyclic peptide metabolites by LC/MS

- Commonly employed LC/MS data processing tools for finding metabolites of small molecule drugs are not well suited for analysis CP metabolites
- CPs often contain unnatural amino acids and peptide linkages
- CP metabolites are difficult to generate high quality and interpretable fragments due to their unique structures (multiple disulphide linkages)
- The number of theoretical metabolite products of a CP drug is huge, searching for all of predicted metabolites is impractical


Objective

- Apply and test recently expanded tool, MetabolitePilot™ software, for automated detection and structural characterization of cyclic peptide metabolites using insulin as model cyclic peptide

Experimental conditions

- Insulin was incubated with Trypsin and Chymotrypsin in buffer for 0 to 3h
- Chromatography separation was performed on a Waters Acquity BEH C18 column (1.7 μm, 2.1 x 100 mm) using mobile phases 0.1% formic acid in water and 0.1% formic acid in acetonitrile.
- High-resolution mass spectrometry data were obtained using the TripleTOF® 5600 system (SCIEX) with IDA acquisitions (CE of 50.0). IDA for MS/MS only include the charge state from 2 to 5, and up to 20 MS/MS spectra with dynamic exclusion.
- The data were analyzed with MetabolitePilot™ software 2.0 based on an insulin-tailored processing method using the TOF-MS predicted metabolite peak finder, multiple-charge filter and TOF-MSMS characteristic product ions peak finder.

Workflow for detecting and identifying metabolites of cyclic peptides using MetabolitePilot™ (Sciex)



Peak finding strategy from MS datasets


- Predicted metabolites strategy finds peaks in accurate mass XICs for parent compound and predicted hydrolytic cleavage products with or without additional biotransformation.
- Generic peak finder strategy finds 3-dimensional peaks in the LC/MS data, then removes entries that do not have charge state within the user selected range.
- If any control samples are used in processing, peaks that are present in control samples are filtered from the results.
- Accurate m/z and peak charge state are confirmed in the TOF MS spectra.

Processing method for searching for peptide metabolites

~46000 predicted catabolites

Identification of insulin metabolites in enzymatic incubation (1h)

Identification of M4, a major metabolite of Insulin in enzymatic incubation



Insulin fragments

161 of both singly and multiply charged fragments (~70% TIC) could be associated with insulin sequence. Due to the insulin structure larger fragments had relatively low abundance but they maintained good mass accuracy.

MS/MS Fragment assignment of GFFYTPKT metabolite

51 (~73% TIC) out of 751 peaks with S/N > 5 were assigned as sequence-related fragments.

Summary of insulin metabolites detected and characterized in enzymatic incubation (1h)


Peak ID	Name	Formula	Neutral	net	Charge	R Time	Peak Area	MS/MS assigned
M1	YTPKT	C18H44N10O9	608.32	305.1688	2	5.4	4.28E+05	✓
M2	QLENYC[*]3N/VLC[*]3GER	C17H30N10O25S2	1442.60	720.305	2	5.8	1.71E+04	✓
M3	PTFTK	C18H40N10O9	755.99	377.9959	3	11.8	3.83E+05	✓
M4	QLENYC[*]3JN/LVC[*]3GER	C18H38N10O25S2	1555.69	777.8425	2	12.4	6.04E+05	✓
M5	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	778.3301	2	12.9	5.20E+04	✓
M6	QLENYC[*]3N/VLC[*]3GER	C17H30N10O25S2	1442.60	720.305	2	14.9	1.83E+05	✓
M7 (M4-H2O)	QLENYC[*]3N/VLC[*]3GER-Loss of Water	C18H36N10O25S2	1537.67	768.8425	2	15.5	5.29E+05	✓
M8	NVC[*]3JN/LVC[*]3GER/[*]2H2O	C18H38N10O25S2	1538.66	770.3301	2	16.1	4.93E+05	✓
M9	QLENYC[*]3N/VLC[*]3GER	C17H30N10O25S2	1442.60	720.305	2	16.4	5.38E+04	✓
M10	GFFYTPK	C48H88N10O12	858.43	429.2147	2	17.7	9.70E+05	✓
M11	GFFYTPK	C48H88N10O12	858.43	429.2147	2	18.2	8.05E+05	✓
M12	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	18.4	2.54E+04	✓
M13	QLENYC[*]3N/VLC[*]3GER	C18H38N10O25S2	1556.66	777.8425	2	18.8	8.15E+04	✓
M14	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY-Loss of Water	C18H36N10O25S2	1537.67	768.8425	2	21.0	2.27E+04	✓
M15	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	21.2	1.46E+04	✓
M16	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	22.4	5.88E+04	✓
M17	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	22.9	2.83E+05	✓
M18 (M4-H2O)	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY-Loss of Water	C18H36N10O25S2	1537.67	768.8425	2	25.1	2.55E+04	✓
M19	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	25.0	3.78E+05	✓
M20	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	25.1	2.55E+04	✓
M21	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	25.4	1.07E+04	✓
M22	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	25.6	2.72E+04	✓
M23	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	26.0	3.78E+04	✓
M24	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	26.2	7.03E+04	✓
M25	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	26.2	4.44E+04	✓
M26	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	27.0	1.83E+05	✓
M27	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	27.4	1.18E+05	✓
M28	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	28.0	2.38E+05	✓
M29	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY	C18H34N10O25S2	1556.66	777.8425	2	28.0	2.78E+05	✓
M30 (M4-H2O)	GVGQC[*]1CT/[*]2TSIC/[*]1JSLY-Loss of Water	C18H36N10O25S2	1537.67	768.8425	2	27.3	3.98E+05	✓

The major enzymatic incubation products have peak areas shaded with dark blue. Details on the disulphide bonds in metabolite sequences are provided with the [*] symbol, chains are separated with / symbol.

Major insulin catabolites found in enzymatic incubation (1h)

Insulin sequence: 51 AA
Avg. Mol Wt: 5807.7, Accurate mass: 5803.67
Formula: C257H383N65O7756

Insulin: Sites of Enzymatic Cleavage for Trypsin & Chymotrypsin
Trypsin: C-Term Cleavage at K,R; Chymotrypsin: C-Terminal Cleavage at L,F,W



Conclusions

- Detection and structural characterization of CP metabolites in biological matrix represent great analytical challenges
- Recently released MetabolitePilot™ 2.0 software employs multiple search mechanisms and automated peptide spectral interpretation tools
- Results from detection and identification of insulin metabolites formed in incubation with Trypsin and Chymotrypsin demonstrate that the software is capable of rapidly identifying metabolites of cyclic peptides
- The application MetabolitePilot™ software in identifying in-vivo metabolites of CP will be studied further

Acknowledgement

Bristol-Myers Squibb
Asoka Ransinghe
William Humphreys
SCIEX
MetabolitePilot™ software Development Team