

Per- and polyfluorinated alkyl compounds (PFAS) analysis of textiles and food contact paper products

Craig M. Butt¹ and Mikyanny Reyes¹

¹ SCIEX, USA

Introduction

This poster describes a targeted method to improve the characterization of PFAS in textiles and food contact materials (FCMs) to ultimately improve our knowledge of PFAS exposure sources. The enhanced sensitivity of the Triple Quad version of SCIEX 7500 system allows for ultra-low detection limits in these commercial products. PFAS in textiles can be absorbed through the skin, or washed off during cleaning to subsequently enter water systems. PFAS in FCMs can be transferred to food items and ingested (Figure 1).

PFAS are used extensively in commercial products, primarily for their stain-resistant properties to repel both water and grease. Although these properties are advantageous for consumers, their potential to act as PFAS exposure sources to humans is concerning to scientists and regulators.

Methods


- Various textiles (tablecloth, pillow protector and several shirts) were purchased from commercial sources. All products identified as containing some PFAS compounds or having stain repellency properties.
- FCMs were obtained from local eateries without knowledge of potential PFAS content.
- A 10 x 10 cm piece of textile or FCM was cut into small pieces, placed in a 15 mL plastic centrifuge tube, spiked with mass-labelled PFAS standards, and 10 mL of methanol added. The tube was sonicated for 30 min, centrifuged and a 1 mL aliquot was transferred to a polypropylene vial for analysis
- Chromatography performed using an Exion AD system modified to replace the fluoropolymer tubing with PEEK. The delay and analytical columns were the Phenomenex Luna Omega PS C18. Mobile phases were water and methanol, both modified with 10mM ammonium acetate. A 12 min gradient was used with a 10 µl injection volume.
- Samples were analyzed using the Triple Quad version of the SCIEX 7500 system in negative electrospray ionization mode. Data were acquired using the Scheduled MRM with compound-specific optimized parameters.




SCIEX 7500 system

Detection of PFAS in commercial textiles

A) PFHxA



B) 6:2 FTS



C) 6:2 diPAPs




Figure 2. Selected PFAS compounds detected in the textile samples. Extracted ion chromatograms (XICs) for PFHxA (top panel), 6:2 FTS (middle panel) and 6:2 diPAPs (bottom panel) in the extraction blank, tablecloth, pillow protector and three different t-shirts.

Textiles key findings:

- C4-C11 perfluorinated carboxylic acids (PFCAs) were detected but their profiles varied by sample
- C5-C11 PFCAs were detected in the pillow protector and C4-C6 PFCAs were detected in t-shirt #2
- 6:2 fluorotelomer sulfonic acid (FTS) and 6:2 diPAPs were detected in t-shirt #3 (Figure 2)
- Perfluorinated sulfonic acids were not detected in any sample

Food contact materials




Figure 3. Detection of diPAPs compounds in a food wrapper. Suite of 6:2 diPAPs, 6:2/8:2 diPAPs and 8:2 diPAPs were detected in 4 of 5 FCM samples.

FCM key findings:

- PFAS profiles differed by sample
- 6:2 diPAPs, 6:2/8:2, diPAPs and 8:2 diPAPs detected in 4 of 5 samples (Figure 3)
- C4-C8 PFCAs were detected in only 2 of 5 samples
- 6:2 FTS detected in 2 samples that did not contain PFCAs

TRADEMARKS/LICENSING

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use, Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to www.sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures. Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see www.sciex.com/trademarks). The images shown may be for illustration purposes only and may not be an exact representation of the product and/or the technology. © 2022 DH Tech. Dev. Pte. Ltd. MKT-26241-A



Figure 1. The entire PFAS life cycle influences human exposure, including textiles and food contact materials